

ACCEPTED

SEP 1 1981

General Insecticide,
Fungicide and Rodenticide Act,
registered under the pesticide
EPA Reg. No. 36736-1

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS & DOMESTIC ANIMALS

DANGER

KEEP OUT OF THE REACH OF CHILDREN

DANGER - Breathing of vapors harmful. Exposure to eyes and skin may cause severe irritation. Avoid contact with eyes, skin or clothing

IN CASE OF CONTACT - Remove all clothing and shoes - flush skin or eyes with plenty of water for at least 15 minutes - **GET MEDICAL ATTENTION IMMEDIATELY.**

PHYSICAL OR CHEMICAL HAZARDS

Contents Under Pressure:
Do not use or store near heat or flame.

Environmental Hazards

LIQUIFIED GAS DANGEROUSLY REACTIVE

STORAGE & DISPOSAL

Store in cool, well ventilated area. Avoid exposure to heat or direct sunlight. When empty return to manufacturer only.

ETHYLENE OXIDE

Active Ingredient
ETHYLENE OXIDE 100%

This product contains 7.25 lbs. of Ethylene Oxide per gallon (at 20° C).

KEEP OUT OF REACH OF CHILDREN

DANGER

STATEMENT OF PRACTICAL TREATMENT

IF INHALED - Remove to fresh air, if breathing has stopped administer oxygen.

IF ON SKIN - Remove contaminated clothing, shower affected areas with water for 15 minutes.

IF IN EYES - Irrigate with copious amounts of water for 15 minutes.

ARC

Distributed by:
ARC CHEMICAL Div. of **Batchem CORPORATION**
P.O. Box 180, Slate Hill, N.Y. 10973

EPA Registration No. 36736-1
EPA Est. No. 36736-NY-01

Net Contents. 400 lbs.

DIRECTIONS FOR USE

GENERAL CLASSIFICATION

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

The product is limited to use by medical professionals or appropriately trained technical personnel in medical and industrial use areas.

FOR THE FOLLOWING INTENDED USES.

1. As a sterilizing gas for industrial and/or medical uses in 100% Ethylene Oxide Commercial Gas Sterilizers Only. Follow the sterilization procedures specified in the gas sterilizers' manufacturers manual and our booklet entitled, "Principles & Practice of Ethylene Oxide Sterilization."

2. Repacking of a sterilizing gas. Refer to our technical brochure, Drummed Ethylene Oxide for chemical and physical properties.

3. Formulation of sterilizing gas mixtures. Refer to our technical brochure Drummed Ethylene Oxide for chemical and physical properties.

Use within 60 days of fill date