UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

JAN 2 9 2007

Ms. Joanne Cashin ARC Specialty Products Balchem Corporation P.O. Box 600 New Hampton, NY 10958

Subject:

36736-3

Sterilizing Gas 3 EPA Registration No.: 36736-3 Submission Date: October 26, 2006 EPA Receipt Date: October 31, 2006

Dear Ms. Cashin,

The following amendment submitted in connection with registration under FIFRA, as amended, is acceptable.

• Deletion of Basil

 Compliance with Tolerance Reassessment and Risk Management Decision for Ethylene Oxide

General Comments

A stamped copy of the accepted labeling is enclosed. Submit one (1) copy of your final printed labeling before distributing or selling the product bearing the revised labeling.

Should you have any questions concerning this letter, please contact Jacqueline McFarlane at (703) 308-6416.

Sincerely,

Velma Noble Product Manager (31) Regulatory Management Branch I Antimicrobials Division (7510P)

·									
		CONCURRENCES					·		
SYMBOL							1		
SURNAME									
DATE		4			*************				
EPA Form	n 1320-1A (1/90)			Printed on Recycled	Paper	"U.S. Go	OFFICI errorent Parales Office	AL FILE COPY	

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

DANGER! CAUSES EYE AND SKIN BURNS, HARMFUL IF INHALED, MAY CAUSE NERVOUS SYSTEM DAMAGE.

DANGER! CANCER HAZARD AND REPRODUCTIVE HAZARD.

EFFECTS OF OVEREXPOSURE: May be fatal if inhaled in high concentrations. May cause irritation of the respiratory tract, chest tightness, headache, nausea, vomiting, diarrhea, light-headed feeling, dizziness, weakness, drowsiness, cyanosis, loss of coordination, convulsions, coma, delayed lung injury (fluid in the lungs), immediate or delayed skin irritation or blisters, allergic skin reaction

OTHER POSSIBLE DELAYED HEALTH EFFECTS: May cause nervous system injury, cataracts, adverse reproductive effects, chromosomal and mutagenic changes, and cancer.

I PPM TWA (as per the Ethylene Oxide Standard 29 CFR 1910.1047). PEL

FI · 5 PPM - excursion limit 15 minutes ODOR:

Ether-like at high concentrations. Exposure to toxic levels may occur without warning to or detection by the user.

PRECAUTIONS- Do not breathe vapor. Do not swallow. Do not get in eyes, on skin, or on clothing. Store and use with adequate ventilation in accordance with the Ethylene Oxide Standard (29 CFR 1910.1047).

PHYSICAL OR CHEMICAL HAZARDS: DANGER! - HIGHLY FLAMMABLE LIQUID AND GAS UNDER PRESSURE.

Contents under pressure. Do not use near flame, sparks, hot surfaces, or allow sources of ignition near the sterilization/fumigation area. Sterilizing Gas 3 is extremely flammable and reactive. Ground all

equipment (including this container) to prevent sparks. LEAK - In case of leak evacuate area and keep personnel upwind. Shut off all sources of ignition. Use self-contained breathing apparatus and protective clothing, and shut off leak if without risk FIRE - In case of fire move container away from fire if without risk.

Use water spray or fog nozzle to keep container cool.


ENVIRONMENTAL HAZARDS:

Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans or other waters unless in accordance with the requirements of a National Pollutant Discharge Effluent System (NPDES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office of the EPA.

It is imperative that users of this material be familiar with ARC Specialty Products' Material Safety Data Sheet for Sterilizing Gas 3, he label and valve tag attached to this cylinder.

EPA Registration No. 36736-3 EPA Establishment No. D 36736-NY-01 D 36736-SC-01 G 36736-MO-01

801139 (Rev. 10/06 - supercedes 12/03)



It is a violation of Federal Law to use this product in a manner

and method of monitoring to be used) refer to the ethylene oxide gas

1. AS A STERILANT AND FUMIGANT GAS:

Keep Out of Reach of Children

DANGER PELIGRO PRECAUTION AL USARIO: Si usted no lee Ingles, no use este producto hasta que la etiqueta le hava sido explicada ampliamente.

Users must follow requirements of the OSHA Occupational Exposure Standard for ethylene oxide (29 CFR 1910,1047).

FIRST AID

IN ALL CASES OF OVEREXPOSURE GET MEDICAL ATTENTION IMMEDIATELY. TAKE PERSON TO A DOCTOR OR EMERGENCY TREATMENT FACILITY AT ONCE.

IF INHALED - Move person to fresh air. Keep warm. If person is Sterilizer manufacturers' Operators Manuals.

immediately.

IF INTALED - Move person to item and the give artificial statisticals, and aseptic packaging, (see 21 CFR 201.1(d)(5)), or to 3. protocol becaution, preferably by mouth-to-mouth, if possible. If breathing is difficult, give oxygen. Call a poison control center or doctor for further treatment advice, even if there are no symptoms. Keep under motional load on cosmetics, whole and ground spices or other motional observation expensions may be delayed. This product may be used only to sterilize medical or laboratory items, 2.

clearance

medical observation - symptoms may be delayed. IF ON SKIN OR CLOTHING- Take off contaminated clothing and shoes. Rinse skin immediately with plenty of water for 15-20 This product may not be used on or in any form of basil. minutes. Call a poison control center or doctor for treatment advice

Aerate, wash or clean contaminated clothing and discard leather After August 1, 2008, this product may only be applied to or

goods. IF SWALLOWED - Call poison control center or doctor sterilization method that uses a single sterilization chamber immediately for treatment advice. Have person sip two glasses of to pre-condition and aerate with an alternating vacuum and water if able to swallow. Do not induce vomiting. Do not give aeration purging procedure. If you wish to employ an anything by mouth to an unconscious person. IF IN EYES - Hold eye open and rinse slowly and gently with water alternative method to that described below, you must contact for 15-20 minutes. Call a poison control center or doctor for the Environmental Protection Agency Office of Pesticide treatment advice Programs for instruction on how to receive authorization.

HOTLINE NUMBER - Have the product container or label with you when calling a poison control center or doctor, or going for

treatment. You may also contact 1-800-424-9300 for emergency "Place spices in the treatment chamber. Assure that the medical treatment information. mixture of ethylene oxide and air is compatible with the NOTE TO PHYSICIAN - Skin exposure to Sterilizing Gas 3 will chamber design, then, introduce into the chamber a commonly result in skin irritation with extensive blister formation. concentration of Ethylene Oxide not to exceed 500 mg/L, with At high concentrations, severe conjunctivitis can occur. Irritation of the respiratory tract may occur, but without acute lung edema. Symp- a dwell time not to exceed 6 hours. Then evacuate the gas toms of systemic intoxication are headache, nausea, vomiting, from the chamber using a sequence of not less than 21 steam incoordination, and cardiac irregularities. Treatment is symptomatic. washes (injections and evacuations) between 1.5 PSIA (27" Refer to ARC Specialty Products' Sterilizing Gas 3 MSDS. If unable Hg) and 5.0 PSIA (20"Hg) while maintaining a minimum to locate MSDS for this product, please call ARC Specialty Products chamber temperature of 115° F." at the telephone number below and request that one be sent

Net Contents: lbs.

B. Sterilization/fumigation with Sterilizing Gas 3 must be performed only in vacuum or gas tight chambers designed for use with Sterilizing Gas 3.

C. Sterilizing Gas 3 cycle parameters depend on several sterilizing/fumigating variable factors: pre-conditioning (if any); exposure time; chantoer air pressure; gas concentration; types and quantities of items to be sterilized/fumigated, packaging, load configuration in the chamber, microbial challenge method; desired degree of disinfection; and the desired performance of the sterilized/fumigated product and package.

D. The sterilization/fumigation cycle parameters should be those prescribed by the equipment manufacturer. If other cycle parameters are used, the safety and efficacy of the alternate cycle parameters must be validated and are the responsibility of the user. NEVER USE PARAMETERS WHICH ALLOW FLAMMABLE MIXTURES OF STERILIZING GAS 3 AND AIR TO ENTER THE CHAMBER.

inconsistent with its labeling. It is the employer's responsibility to follow the requirements of the Ethylene Oxide Standard (29 CFR 1910.1047). 2. INSTALLATION AND OPERATION OF THE CONTAINER -This product may be used only in facilities that meet the requirements of Follow the directions on the value tag attaching this container.

the Ethylene Oxide Standard (29 CFR 1910.1047). This product may be used only in non-portable (commercial) vacuum or gas-tight chambers STORAGE A DISPOSA

designed for use with 80% ethylene oxide. 20% carbon dioxide. This Do not contaminate food, feed, or water by storage and disposal. product may be used only by persons who have been trained in **PESTICIDE STORAGE** Use in accordance with tag attached to valve. accordance with the Ethylene Oxide Standard (29 CFR 1910.1047), Store in cool, well-ventilated area. Avoid exposure to heat or direct When used to sterilize health care items, this product must be used in sunlight as may cause polymerization.

non-portable (commercial) ethylene oxide gas sterilizers that have FDA PESTICIDE DISPOSAL Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray, or mixture of rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste representative at A. For complete use directions (including type of surfaces, objects, or the nearest EPA Regional Office for guidance.

items/products recommended for treatment, pre-cleaning instructions, CONTAINER DISPOSAL When empty, return to supplier only. concentration of gas per unit volume of closed space to be treated, Before returning container to supplier: exposure time/temperature, relative humidity, ventilation/aeration time, 1.

Cylinder pressure should be below 250 psig when empty. Check that container valve is closed.

Replace valve plug tightly in valve outlet.

Check container valve and plug for leaks prior to shipment.

USE THIS PRODUCT WITHIN 120 DAYS OF ITS PURCHASE.



In case of accident or emergency, call CHEMTREC at 800-424-9300

> Distributed By: **ARC Specialty Products Balchem Corporation** P.O. Box 600 · New Hampton, NY 10958 Tel: 845-326-5611 Fax: 845-326-5706 www.arcspecialtyproducts.com