

UNITED STATES ENVIRONMENTAL PROTECTION AGENCY Washington, D.C. 20460

OFFICE OF PREVENTION, PESTICIDES AND TOXIC SUBSTANCES

SEP 15 2008

Joanne K. Cashin Registration Compliance Advisor Arc Specialty Products P.O. Box 600 New Hampton, New York 10958

Subject:

Sterilizing Gas 8

EPA Registration 36736-7

Amendment Dated: June 16, 2008 EPA Received Date: June 17, 2008

Dear Joanne K. Cashin

The following amendment submitted in connection with registration under the Federal Insecticide, Fungicide, and Rodenticide Act (FIFRA) as amended, is acceptable.

Proposed Amendment

Revised Label Adding Language for Ethylene Oxide

General Comments

A stamped copy of the accepted label is enclosed for your records. Submit (3) copies of your final printed label before distributing or selling the product bearing the revised labeling.

Should you have any questions or comments concerning this letter, please contact Velma Noble at (703)308-6233.

Sincérely

CONCURRENCES

YMBOL

TO INCLUDING THE PRODUCT Manager (31)

Regulatory Management Branch1

Artificionial Division (7310)

EPA Form 1320-1A (1/90)

Printed on Recycled Pares

OFFICIAL FUE

ACCEPTED SEP 1 5 2008 Under the Federal Insecticide, Fungicide, and Riggenticida Act as amended, for the pesicido registerad undor EPA Ruy, No 36736-7

PRECAUTIONARY STATEMENTS

......

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

DANGERI CAUSES EYE AND SKIN BURNS, HARMFUL IF INHALED, MAY CAUSE NERVOUS SYSTEM

DANGER! CANCER HAZARD AND REPRODUCTIVE HAZARD.

EFFECTS OF OVEREXPOSURE: May be fatal if infralled in high concentrations. May cause irritation of the respiratory tract, chest lightness, headsche, nauses, voretring, diamhes, lightness def feeling, dizziness, weakness, droveshess, cyanosis, loss of coordination, convadiance, come, delayed harp july; (Bata in the lungs), immediate or delayed this initiation or bidsets, along such resction. OWER PASSAILE SELAYED HEALTH EFFECTS: convictions a flects.

May cause nervous system injury, cataracts, adverse reproductive effects, chromosomal and mutagenic changes, and cancer.

PEL: 1 PPM TWA (es per the Ethylene Oxide Standard 29 CPR 1910, 1047). EL: 5 PPM - accuration 11mrt 15 minuths. ODOR: Ether-libs at high consentrations. Exposure to toxic levels may occur without warming to ar detection by the user.

PRECAUTIONS- On not breathe vapor. Do not swallow, Do not get in eyes, on skin, or on clothing. Store and use with adequate ventilation in accordance with the Ethylene Oxide Standard (29 CFR 1910.1047).

PERSONAL PROTECTION EQUIPMENT (PPE):
A material that is chemical-registant to this product is butyl rubber

A material wat is chemical-resistant to the product is ourly rupper.

A material must weer at a minimum:
ong-deeved shirt and long plants,
flose plus social and long plants,
flose plus social must be a made of the plants.

Demical-resistant glowes, and

The morphoyer shalt provide a respirator that is adequate to protect
the health of the employee and ensure compliance with all other
OSHA statutory and regulatory requirements (see 29 CFR
Hein hands could have yet or soin contact with ETO or ETO solutions,
coils all the country of the plants o

- Follow the respirator manufacturer's user's instructions for
- results from the confirmation of the confirmat
- Respirator users must be trained using a program that conforms to OSHA's requirements (see 29CFR Part 1910,134).

In case of accident or emergency, call CHEMTREC at 1-800-424-9300.

DO NOT REMOVE THIS LABEL USE THIS PRODUCT WITHIN 120 DAYS OF ITS PURCHASE Respirator users must be examined by a qualified medical practitioner to ensure physical ability to safely wear the style of respirator to be worn. A qualified medical practitioner is a physical or of the forested health can be of the safely safely respirator or consists of a questionnaire that safely respirator consists of a questionnaire that safe about medical conditions (such as a heart condition) that would be problematic or respirator use. If concerning are identified, then additional evaluations, such as a physical exam, might be necessary. The initial evaluation must be often before respirator use begins, it does not need to be repeated use to begin a first of respirator use conditions change (see 250.FR Part 1910, 13-5).

Follow manufacturer's instructions for cleaning/maintaining PPE. If no such instructions for washables exist, use detargent, and hot water. Keep and wash PPE separately from other laundry.

Users should wash hands before eating, drinking, chewing gum, using tobacco, or using the tollet.

Users should remove clothing/PPE immediately if posticide gets name. Then wash thoroughly and out on clean ciothing.

Users should remove PPE immediately after handling this product. Wash the outside of gloves before removing. As soon as possible, wash thoroughly and change into clean clothing.

PHYSICAL OR CHEMICAL HAZARDS:
DANGERI - NONF-LAMMABLE LIQUID AND GAS UNDER PRESSURE.
Containts under pressure. Do not use near flame, sparks, hot surfaces, or
allow sources of syntion near the sitezation/flumigation ares. Ground sill
exaginment (including this content) to prevent parks.
LEAK - in case of leak evacuate area and teep personnel upwind. Set
of all sources of lignifion. Use self-contained breathing apparatus and
protective clothing, and struct off leak if without risk.
FIRE - In case of the move container away from fire if without risk. Use

FIRE - In case of the move container away from the if webout risk. Use water spray or for noze to lose pomeliar eroot. ENVIRONMENTAL HAZARDS: Do not discharge effluent containing this product into takes, steams, ponds, estabries, occars or other waters unless in accordance with the requirements of a National Politant Olscharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing the product to severe systems without previously notifying the local sewage treatment parts authority or guidance confact your State Weter Doard or Regional it is imperative that users of this material be timidar with ARC Specially Productal Machinel Selety Data Sheet for Starkzing Gas 8, the label and valve tag attached to this cylinder.

EPA Registration No. 36736-7 EPA Establishment No.

☐ 36736-NY-01

T 36736-MO-01

GAS 8

STERILANT/FUMIGANT ACTIVE INGREDIENT: ETHYLENE ONIDE... 8.5%
(CAS NO. 75-21-8) INERT INGREDIENT: CARBON DIONIDE. ..91.5%

STERILIZING

Keep Out of Reach of Children

DANGER PELIGRO

PRECAUTION AL USARIO: Si usted no lee Ingles, no use este producto hasta que la etiqueta le hava sido explicada ampliamente

Users must follow requirements of the OSHA Occupational Exposure Standard for ethylene oxide (29 CFR 1910.1047).

FIRST AID

IN ALL CASES OF OVEREXPOSURE GET MEDICAL ATTENTION IMMEDIATELY, TAKE PERSON TO A DOCTOR

IF INHALED - Move person to fresh air. Keep warm, if person is not breathing, call 91 for an embutance, then give artificial respiration, preserveby by mouth-to-mouth, of possible. If breathing is difficult, give oxygen. Call a poison control center or doctor for further treatment advice, even if there are no symptoms. Keep under medical observat

across, even it there are no symptoms. Keep under medical observation—symptoms may be delayed.

IF ON SKIN OR CLOTHING—Text contaminated clothing and sold since skin immediately with plenty of water for 15-20 minutes. Call a poisson control center or dector for treatment advice. Aserts, wesh or clean contrainated clothing and discard lesther goods.

IF SWALLOWED—Call poisson control center or dector immediately for treatment advice. Have person spit two glasses of water if able to swallow. Do not induce voretting. Do not give enything by mouth to an unconscious person.

unconscious person.

If IN EYES – Hold eye open and rinse slowly and gently with water for 15-20 minutes. Call a poison control center or doctor for treatment.

advice.

HOTLINE NUMBER — Have the product container or label with you when calling a poison control center or doctor, or going for treatment. You may also contact 1-800-424-9300 for emergency medical treatment.

Total may also contact in County and the state of the sta

Refer to ARC Specialty Products' Stankzing Gas 8 MSDS, if unable to locate MSDS for this product, please call ARC Specialty Products at the telephone number below and request that one be sent immediately.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistant with its labeling. Employers in feedless that use ETO must comply with all of the negatiments for ETO use specified in 22 GPR 1910.1047. This product has required to the constraint of th

In contract sterilization facalities, including facilities treating medical equipment and supplies, musical instruments, liberar/museum aritacts, cosmetics, and spices the following requirements must be followed:

Sterilization/furrigation with ETO must be performed only in vacuum or gas tight chambers designed for use wiff ETO of the sterilization of the state of the sterilization of the sterilization of all employees including office staff. Information and training must be provided to all employees in the facility at the time of initial assignment and annually therester. The safety training must include, at a minimum, the following information:

- The most recent monitored ambient levels of ETO in the facility;
 The potential health effects from the levels of ETO in the facility;
 The emergency response plan and how to respond in an emergency;
 The availability of the Material Safety Data Sheet and other materials related to the health hazards of exposure to ETO.

In order to reduce ambient levels of ethylene oxide, lengthy facility aeration is encouraged, it can reduce potential long-term risks to employees not directly involved in the ethylene oxide applications.

Air monitoring should included the entire facility including office space, break areas, and loading/unloading areas.

1. AS A STERILANT AND FUMRGANT GAS: A. For complete use directions (including type of surfaces, objects, or immerploature, recommended for treatment, pre-clashing instructions, concentration of gas per unit volume of closed space to be treated, exposure transfaragement, enlative hurseity, vertilation/verariation ren, and method of manufacturers. Operation Manuals.
This product may be used only to startize medical or laboratory flems, pharmaceuticals, and saeptic packaging, (see 21 CFR 201.1(d)(5)) or to reduce microbial load on commence, whose and ground spices or other assauding materials (see 40 CFR 180.151), and artifacts, archival material or items of being sold products.

Distributed By: **ARC Specialty Products Balchern Corporation** P.O. Box 600 · New Hampton, NY 10958 Tel: 845-326-5611 Fax: 845-326-5706 www.arcspecialtyproducts.com

Net Contents:

This product may not be used on or in any form of heaft.

After August 1, 2008, this product may only be applied to or on spices. dried vegetables or seasonings utilizing an ETO sterilization method that uses a single sterilization chamber to pre-condition and serata with an afternating vacuum and seration purping procedure. If you with no afternating vacuum and seration purping procedure, if you wish to employ an afternative method to that described below, you must contact the Environmental Protection Agency Office of Pesticide Programs for instruction on how to receive authorization.

Place spices in the treatment chamber. Assure that the mixture of ethylene oxide and air is compatible with the chember design, then, introduce into the chamber a concentration of Ethylene Oxide not to exceed 500 mg/L, with a dwell time not to exceed 6 hours. Then evacuate the gas from the chamber using a sequence of not less than 21 strem washes (injections and evacuations) between 1.5 PSIA (27" Hg) and 5.0 PSIA (20"Hg) while maintaining a minimum chamber

amperature of 110-7.

B. Shriftzion/furnigation with Stariiting Gas 8 must be performed only in vacuum or gas light chambers designed for use with Starizing Gas 8.

C. Stariiting Gas 8. cycle parameters depend on several stariiting/furnigating variable factors; pre-conditioning (if any); exposure time; chamber as pressure; gas concentration, types and quantities of items to be stariitized/furnigated; packageling, load configuration in the chamber; so the stariitized/furnigated product end package.

D. The stariitization/furnigated product end package.

D. The stariitization/furnigated or product end package.

D. The stariitization/furnigated product end package.

D. The stariitization/furnigated product end package are used the stariity and efficacy of the altamate cycle parameters must be validated and are the responsibility of the user. NEVER USE PARAMETERS WHICH ELOY IL MANDER.

. INSTALLATION AND OPERATION OF THE CONTAINER - Follow the

STORAGE AND DISPOSAL storage and disposal.

DO not contaminate food, feed, or water by PESTICIDE STORAGE

PESTICIDE STORAGE
Use in accordance with tag attached to valve. Store in cool, well-ventilated usen.
Avoid exposure to heat or durect consight us may cause polymorated.
Pesticide Wastes are acutely hazardous. Improper disposand of excess pesticide,
stays, or metture of finasks is a violation of Federal I aw. If these wastes cannot be disposed of by the according to label instructions, contact your State Pesticide or
Environmental Control / Agency, or the Hazardous Waste representative at the
nearest EPA Regional Office for guidance.

- When empty, return to supplier only. Refore returning container to supplier:

 1. Cylinder pressue should be below 600 prig when empty. Check that container valve is closed.
- Replace valve plug tightly in valve outlet
 Check container valve and plug for leaks prior to shipment.

801468 (Rev. 06/08 - supersedes 10/06)

