


107
(FRONT PANEL)

CAVICIDE®

HOSPITAL DISINFECTANT/ DECONTAMINANT CLEANER
SALON DISINFECTANT/CLEANER
VETERINARY DISINFECTANT/CLEANER
CLIPPER BLADE DISINFECTANT/CLEANER
DISINFECTANT/CLEANER
DENTAL DECONTAMINANT/CLEANER
ONE-STEP DISINFECTANT/CLEANER
INSTITUTIONAL DISINFECTANT/CLEANER
MEDICAL DECONTAMINANT/CLEANER
SURGICAL DECONTAMINANT/CLEANER
SURFACE/INSTRUMENT DISINFECTANT/CLEANER
SURFACE DISINFECTANT/DECONTAMINANT CLEANER
LABORATORY SURFACES DECONTAMINANT CLEANER
SPRAY-ON SURFACE DISINFECTANT



FOR PROFESSIONAL USE

IMMERSION SOLUTION • SURFACE CLEANER/DISINFECTANT • ULTRASONIC SOLUTION
SONIC SOLUTION • SOAKING SOLUTION • PRESOAK SOLUTION • INSTRUMENT SOLUTION
SURFACE SOLUTION • NON-POROUS SURFACE SOLUTION • CLEANER • DISINFECTANT

BACTERICIDAL • VIRUCIDAL • FUNGICIDAL • TUBERCULOCIDAL** • PSEUDOMONICIDAL
STAPHYLOCIDAL •

READY TO USE

CONTAINS BIODEGRADABLE DETERGENT

ACTIVE INGREDIENTS:

Diisobutylphenoxyethoxyethyl dimethyl benzyl ammonium chloride	0.25%
Isopropanol	15.30%
INERT INGREDIENTS	84.45%
TOTAL	100.00%

KEEP OUT OF REACH OF CHILDREN

CAUTION

PRECAUTIONARY STATEMENTS: Harmful to Humans and Domestic Animals.

AVOID CONTAMINATION OF FOOD. MAY CAUSE EYE IRRITATION. AVOID DIRECT CONTACT WITH EYES. IN
CASE OF DIRECT EYE CONTACT, IMMEDIATELY FLUSH EYES WITH PLENTY OF WATER FOR AT LEAST 15
MINUTES. IF IRRITATION PERSISTS, SEEK MEDICAL ATTENTION.

NET CONTENTS: 1 U.S. GALLON / 3.785 Liters (NET WEIGHT 8.22 lbs / 3.73 kg)

MICRO-ASEPTIC PRODUCTS, INC.

887 E. WILMETTE ROAD • PALATINE, IL 60067 USA

EPA REG. NO. 38526-1 EPA EST. NO. 39234-WV-001 REORDER NO: CO4-128

(LEFT PANEL)

CAVICIDE

- NON STAINING
- NON CORROSIVE
- NON IRRITATING
- NO DILUTION
- READY TO USE
- NO TOXIC FUMES

EFFECTIVE AGAINST:

- Staphylococcus aureus
- Pseudomonas aeruginosa
- Salmonella choleraesuis
- Mycobacterium tuberculosis var: bovis (BCG)**
- Trichophyton mentagrophytes
- Aspergillus niger
- Herpes simplex virus type 1 and 2 *
- Poliovirus type 1 and 2 *
- Coxsackievirus *
- Human Immunodeficiency Virus (HIV-1) (AIDS virus) *
- Mold and Mildew

(* on inanimate surfaces)

(** in ten minutes at room temperature (20°C))

DIRECTIONS FOR USE: It is a violation of U.S. Federal law to use this product in a manner inconsistent with its labeling.

DESCRIPTION:

Cavicide is a multi-purpose, broad spectrum, ready to use, highly effective cleaner and disinfectant for use on the surfaces of inanimate objects. It is especially useful in hospital operating rooms, emergency departments, isolation areas, neonatal units, dental operatories, surgical suites, animal care facilities, beauty salons, salon settings, manicure salons, skin care salons, barber shops, bathrooms, tanning salons, out-patient surgical centers, daycare centers, schools, ambulances, police and fire vehicles, jails, prisons, morgues, patient care areas, laboratories, food preparation areas, storage areas, health club facilities, and other critical care areas where environmental infection control and the control of cross contamination is important.

Safe for cleaning/decontamination of delicate medical/dental/surgical/salon/environmental/equipment/implements and instrumentation. Cavicide will effectively clean and disinfect, when used as directed, such items as: infant incubators and bassinets, infant care cribs and warmers, infant care equipment surfaces, oxygen hoods, anesthesia and respiratory therapy equipment surfaces, operating room tables and lights, laboratory equipment and surfaces, physical therapy (PT) equipment surfaces, neck brace appliances and cervical collars, whirlpool tanks, hydro-therapy equipment and tanks, stretchers, mayo stands, countertops, toilets, sinks, refrigerator units, floors, walls, handrails, door knobs, bed railings, bathing units, bath tubs, shower stalls, cabinets, shampoo bowls, manicure tables, chairs, workstations, tanning beds, hair dryers, telephones, diaper changing stations, baby cribs, hair clippers, shears, razors, hair cutting implements, clipper blades, salon surfaces, scissors, combs, brushes, manicure implements, washable nail files, hair rollers, animal cages, veterinary care surfaces, dental operatory surfaces, dental countertops, dental chairs, unit stools, light lense covers, curing lights, Panoramic bite blocks, pliers, impression materials (hard rubber, condensation silicone and vinyl polysilicones), dentures and other inanimate surfaces, including those made of plastics (such as: polycarbonate, polyvinylchloride, polypropylene and polystyrene), non-porous vinyl and upholstery, stainless steel, painted surfaces, plexiglas, glass, and other hard non-porous surfaces.

APPLICATIONS:

SURFACES: (Where appropriate, follow Universal Precautions.)

For disinfecting non-critical devices/medical equipment and other surfaces:

- Spray/apply Cavicide directly to surface, thoroughly wetting area to be disinfected. (**** Visibly soiled surfaces should be pre-cleaned.) Allow surface to remain wet for 2 minutes. (**FOR TUBERCULOCIDAL ACTIVITY:** Allow surface to remain wet for 10 minutes at room temperature (20°C).) Follow by wiping surface with a fresh, clean, paper or cloth towel; or rinse and either allow surface to air dry or wipe rinsed surface dry using a fresh, clean, paper or cloth towel. Discard towel.

- Cavicide completely inactivates the HIV-1 (AIDS virus) on hard, non-porous surfaces in the presence of a moderate amount of organic soil (5% blood serum) with a contact time of 2 minutes at room temperature (20-25°C).

****** For pre-cleaning visibly soiled medical equipment and other surfaces prior to disinfection:**

- Apply Cavicide directly to surface. Allow to remain wet for about 30 seconds. Wipe surface clean using a clean paper or cloth towel or rinse surface and either wipe dry or allow to air dry. Discard dirty towel.

INSTRUMENT/IMPLEMENT CLEANING INSTRUCTIONS: (Where appropriate, Universal Precautions.)

For use as a pre-cleaning instrument decontaminant solution:

- Fill appropriate size container with a sufficient amount of undiluted Cavicide so allow for complete submersion of instruments/objects. Place objects into Cavicide solution, cover and allow to soak prior to transporting and/or cleaning and terminal sterilization / high level disinfection. Rinse thoroughly. Change solution daily or more often as needed if the solution becomes diluted or visibly soiled.

For use as an instrument pre-clean/pre-transport spray:

- Place instruments onto or into a suitable container. Without handling instruments, thoroughly spray Cavicide solution onto instruments so as to thoroughly drench all surfaces. Cover instruments and transport to appropriate cleaning area. Rinse instruments and follow with the appropriate cleaning and terminal sterilization/high level disinfection process.

For use as an instrument ultrasonic cleaning solution:

- Thoroughly rinse objects under running water to remove visible gross debris. Add 1 ounce Cavicide per liter of water in ultrasonic unit, immerse instruments into mixture of solution and activate ultrasonic unit for 5 minutes or longer if necessary. Remove objects and rinse thoroughly. Change solution as needed. Critical and semi-critical devices should be followed by appropriate terminal sterilization/high level disinfection process.

For use as a manual cleaner:

- Rinse dirty objects under running tap water to remove visible gross debris. Place pre-rinsed objects into a solution of 1 ounce Cavicide per liter of ordinary tap water. Scrub objects using a stiff bristle brush. (Objects should be submerged as scrubbing.) Rinse thoroughly. Change solution as needed. Critical and semi-critical devices should be followed by appropriate terminal sterilization/high level disinfection process.

For use on hair clippers, shears, barber implements:

- While the clipper is running, hold it in the downward position and spray undiluted Cavicide directly onto the blades two or three times so as to thoroughly wet the blades. (Avoid getting the spray on the clipper case or allowing it to run into the inside of the clipper housing.) Allow to remain wet for 2 minutes before wiping dry with a clean cloth. Lubricate as per clipper manufacturer's instructions.

For cleaning salon implements:

- First, spray object using a pump spray bottle of Cavicide. Scrub/wipe away visible debris using a soft bristle brush or soft cloth. Immerse implements into an undiluted solution of Cavicide for 10 minutes at room temperature (69°F). Remove and wipe dry. No rinsing is necessary. Change solution weekly or more often if solution becomes visibly soiled.

1076

INSTRUMENT/IMPLEMENT/SMALL OBJECT DISINFECTION INSTRUCTIONS:

(Where appropriate, follow Universal Precautions)

***** For disinfection of non-critical, pre-cleaned instruments/objects:**

- Instruments/objects must be thoroughly pre-cleaned to remove excess organic debris, rinsed and then rough dried. (Clean and rinse the lumens of hollow instruments before filling with solution or before immersion.) Immerse instruments/objects into undiluted Cavicide solution and allow to remain submerged for 10 minutes at room temperature (20°C) for tuberculocidal kill. Remove and rinse or wipe dry prior to use. Change solution daily or more often as needed if the solution becomes diluted or visibly soiled.

***** For disinfection of non-critical, pre-cleaned instruments/objects in ultrasonic units:**

- Instruments/objects must be thoroughly pre-cleaned to remove excess organic debris, rinsed and rough dried. (Clean and rinse the lumens of hollow instruments before filling with solution or before immersion.) Immerse pre-cleaned instruments/objects into undiluted Cavicide solution and activate ultrasonic unit for 10 minutes. Remove and rinse or wipe dry prior to use. Change solution daily or more often as needed if the solution becomes diluted or visibly soiled.

******* This product is not to be used as a terminal sterilant / high level disinfectant on any surface or instrument that (1.) is introduced directly into or in contact with the blood-stream or normally sterile areas of the body, or (2.) that contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization / high level disinfection.

MOLD AND MILDEW: To control mold and mildew on clean, hard surfaces, apply so as to wet entire surface thoroughly with Cavicide. Allow to air dry after application. Repeat application in seven days or as necessary to maintain control.

CAVICIDE EFFECTIVELY KILLS HIV ON PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS IN HEALTH-CARE SETTINGS OR OTHER SETTINGS IN WHICH THERE IS AN EXPECTED LIKELIHOOD OF SOILING OF INANIMATE SURFACES/OBJECTS WITH BLOOD/BODY FLUIDS, AND IN WHICH THE SURFACES/OBJECTS CAN BE ASSOCIATED WITH THE POTENTIAL FOR TRANSMISSION OF HUMAN IMMUNODEFICIENCY VIRUS TYPE 1 (HIV-1) (ASSOCIATED WITH AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 (HUMAN IMMUNODEFICIENCY VIRUS OR AIDS VIRUS) OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUID:

Personal Protection: Wear appropriate barrier protection such as latex gloves, gowns, masks or eye coverings.

Cleaning Procedure: Blood and other bodily fluids must be thoroughly cleaned from surfaces and objects before disinfection with Cavicide.

Contact Time: While the HIV-1 virus is inactivated in 2 minutes, use the recommended contact time for the disinfection of other organisms listed on this label.

Infectious Materials Disposal: Cleaning materials used that may contain blood or other bodily fluids should be autoclaved and/or disposed of in accordance with local regulations for infectious materials disposal.

For product information, please contact our technical service department at 1-800-536-4129 (TOLL FREE).

STORAGE: Store in a cool place.

PESTICIDE DISPOSAL: Dilute with water. Dispose of in ordinary sanitary sewer.

CONTAINER DISPOSAL: Do not reuse empty container. Wrap empty container and place into ordinary trash receptacle.

Cavicide spray bottles are refillable.

Manufactured For:

Micro-Aseptic Products, Inc. • 887 E. Wilmette Rd. • Palatine, IL 60067 USA

Revised 8/30/94