

Reg no: 38526-2

PM-31

1 of 2

AMENDED LABEL COPY

TO REFLECT EPA

REQUIREMENTS

MANIKAID No. 7 SURFACE DISINFECTANT  
(Peppermint fragrance)

2-19-86

(AWN)

Staphylocidal, Salmonellacidal, Pseudomonocidal  
after 10 minutes contact time and Virucidal\* on hard nonporous  
surfaces after 30 seconds contact time

FOR HOSPITAL/GENERAL DISINFECTANT ~~PARAMEDIC~~ & MMR/CPR USE

\*Kills Herpes Simplex Virus Types 1 and 2 on Nonporous,  
inanimate surfaces after 30 SECONDS Contact Time

ACTIVE INGREDIENTS:

Diisobutylphenoxyethoxyethyl dimethyl- benzyl ammonium chloride .....	0.05%
Ethyl alcohol .....	18.00%
INERT INGREDIENTS .....	81.95%

KEEP OUT OF REACH OF CHILDREN

CAUTION

Professional Size  
Net Contents: 6 fl.ozs. (180 ml.)  
EPA Reg. No. 38526-2  
EPA Est. No. 26213 MD-2

~~MICROECOLOG~~ PRODUCTS CORPORATION  
7002 Golden Ring Rd.  
Baltimore, MD. 21237  
ARC PO# WAD 0000345

ACCEPTED

NOV 23 1987

Under the Fungicide, Insecticide, and Rodenticide Act, as amended, for the pesticide registered under EPA Reg. No. 38526-2

MANIKAID CONTROL NO. 7

A DISINFECTANT PRODUCT FOR  
MANIKIN USE. Meets EPA Requirements  
for use as a hospital disinfectant.

This product is used as a hospital/general disinfectant by MMR/CPR students, paramedic and training personnel. When tested as a hospital disinfectant, it killed Staphylococcus aureus, Salmonella choleraesuis and Pseudomonas aeruginosa in 10 minutes when mixed with 5% calf serum and dried as a film onto hard, nonporous inanimate surfaces. Manikaid Control No. 7 killed Herpes Simplex Virus Type 1 and 2 when mixed with a 5% calf serum and dried onto hard nonporous, inanimate surfaces.

DIRECTIONS:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.


- 1. Remove cap to expose pump sprayer valve.
- ~~2. When used on any hard surface item as a hospital/general disinfectant allow to set for 30 seconds. If wiped, use a sterile aseptic wipe.~~
- 2. Initially disinfect mouth and nose parts of manikin by spraying until thoroughly wet. Allow to remain wet for at least 30 seconds for virucidal effectiveness and/or 10 minutes for bactericidal results, and then wipe dry with gauze.
- 3. Repeat application after each student use.
- 4. Upon completion of the instruction period, disinfect manikin nose and mouth parts as in #3 above.

Storage and Disposal

Wrap and dispose of empty container in regular trash collection.

Distributed By:

MICROECOLOGY PRODUCTS CORPORATION  
7002 Golden Ring Road  
Baltimore, MD 21237  
301/574-1344



BEST AVAILABLE COPY