



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
WASHINGTON, DC 20460

OFFICE OF CHEMICAL SAFETY
AND POLLUTION PREVENTION

March 24, 2021

Rebecca Lei
Senior Regulatory Affairs Specialist
Metrex Research
1717 West Collins Avenue
Orange, CA 92867

Subject: PRIA Label and CSF Amendment –Adding new organisms and efficacy claims to the label and revising Basic and Alternate CSFs 1 and 4
Product Name: Caviwipes
EPA Registration Number: 46781-8
Application Date: 10/01/2020
Action Case Number: 00217530

Dear Ms. Lei:

The amended label and CSFs referred to above, submitted in connection with registration under the Federal Insecticide, Fungicide and Rodenticide Act, as amended, are acceptable. This approval does not affect any conditions that were previously imposed on this registration. You continue to be subject to existing conditions on your registration and any deadlines connected with them.

A stamped copy of your labeling is enclosed for your records. This labeling supersedes all previously accepted labeling. Pursuant to 40 CFR 156.10(a)(6), you must submit one copy of the final printed labeling before you release the product for shipment with the new labeling. In accordance with 40 CFR 152.130(c), you may distribute or sell this product under the previously approved labeling for 18 months from the date of this letter. After 18 months, you may only distribute or sell this product if it bears this new revised labeling or subsequently approved labeling. "To distribute or sell" is defined under FIFRA section 2(gg) and its implementing regulation at 40 CFR 152.3.

Please note that the record for this product currently contains the following CSFs:

- Basic CSF dated 03/02/2021
- Alternate CSF 4 dated 03/24/2021

Should you wish to add/retain a reference to the company's website on your label, then please be aware that the website becomes labeling under FIFRA and is subject to review by the Agency. See FIFRA section 2(p)(2). If the website is false or misleading, the product would be misbranded and unlawful to sell or distribute under FIFRA section 12(a)(1)(E). 40 CFR 156.10(a)(5) lists examples of statements EPA may consider false or misleading. In addition, regardless of whether a website is referenced on your product's label, claims made on the

website may not substantially differ from those claims approved through the registration process, FIFRA section 12(a)(1)(B). Therefore, should the Agency find or if it is brought to our attention that a website contains false or misleading statements or claims substantially differing from the EPA approved registration, the website will be referred to the EPA's Office of Enforcement and Assurance.

Your release for shipment of the product constitutes acceptance of these conditions. If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA section 6. If you have any questions, please contact Jake McFarley by phone at (703) 347-0123, or via email at McFarley.Jake@epa.gov.



Eric Miederhoff
Product Manager 31
Regulatory Management Branch I
Antimicrobials Division (7510P)
Office of Pesticide Programs

Enclosure: Master Label

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

ACCEPTED

03/24/2021

Under the Federal Insecticide, Fungicide and Rodenticide Act as amended, for the pesticide registered under EPA Reg. No. 46781-8

CaviWipes™

[TUBERCULOCIDAL**] [BACTERICIDAL] [VIRUCIDAL*] [FUNGICIDAL]

<p>Active Ingredients: Diisobutylphenoxyethoxyethyl dimethyl benzyl ammonium chloride.0.28% Isopropanol 17.20% Inert Ingredients..... 82.52% Total..... 100.00%</p>	<p>KEEP OUT OF REACH OF CHILDREN CAUTION [See other panel for additional precautionary statements.]</p>
<p>EPA Reg. No. 46781-8 EPA Est. [56952-WI-001] [•71979-SC-001] [•75272-WI-001] [•51435-WI-004] [•46781-MI-1] [insert EPA Est. No. from CSF]</p> <p>[Beginning of batch code indicates Est No.] [End of batch code indicates Est No.]</p> <p>Made in USA</p>	<p>Reorder No. x-x Net Contents: #[[lbs.] [oz]] # Pre-Saturated Towelette # x # in</p> <p>[Packed For:] [Distributed By:] [Sold By:] [Manufactured By [or for]:] [Manufacturer]</p> <p>Metrex Research 28210 Wick Road Romulus, MI 48174</p> <p>[To reorder] [,] [For product or technical information] contact Metrex at 800-841-1428 or visit our web site at www.metrex.com [For SDS information please visit www.metrex.com]</p>

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRECAUTIONARY STATEMENTS
Hazards to Humans and Domestic Animals

CAUTION

Harmful if absorbed through skin. Causes moderate eye irritation. Avoid contact with eyes, skin or clothing. Wash hands before eating, drinking, chewing gum, using tobacco, or using the toilet. Remove contaminated clothing and wash clothing before reuse.

FIRST AID

If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.

If on skin: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.

[For Chemical Emergency ONLY (spill, leak, fire, exposure or accident), call CHEMTREC at 1-800-424-9300]

For [ALL] [other] inquiries about this product, call Metrex Customer Care at 1-800-841-1428, Monday through Friday, between 6 a.m. and 4 p.m. Pacific Time.

Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

Packaging Statements and Wipe Dispensing Instructions:
[Use A or B or C as applicable depending on container type:]

[A - Refillable, Non-refillable canister, buckets, tubs or pails]

[To dispense towelettes] [container] [canister] [pail] [bucket] [tubs] [Dispensing Directions]: [To start feed, remove cover and discard seal. From center of towelette roll, pull up a towelette corner, twist it into a point and thread it through the hole in the container cover. Pull through about one inch. Replace cover. Pull out first towelette and tear off at a *[(any angle)] [45°-90° angle]*. Remaining towelettes feed automatically, ready for next use. When not in use, keep center cap closed to prevent solution loss.]

[B - Flat Pack]

[To dispense towelettes] [from flat pack:] [from Flat Pack wipes container:]

[Grasp the white tab that states "PULL HERE". Gently swing the lid open and grasp a towelette. Remove towelette from package and then swing lid closed. Gently tap the lid to ensure seal is made between it and the towelette package.]

[C - Single use packet]

Packet Directions:

[Tear open packet. Unfold towelette and use as directed. Discard as directed.]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.
[Use A or B or C as applicable depending on container type]

[A – Nonrefillable container, canister, buckets or pails] STORAGE AND DISPOSAL:

Product Storage: [Store at room temperature.] [When not in use, keep center cap of lid closed to prevent moisture loss.]

Towelette Disposal: Do not re-use towelette. Dispose of used towelette in trash. [Do not flush in toilet.]

Container Handling: Nonrefillable container. Do not reuse or refill this container. Offer for recycling, if available. If recycling is not available, discard in trash.

[B – Flat pack]

STORAGE AND DISPOSAL

Product Storage: [Store at room temperature.] [When not in use, keep [flap] [tab] closed to prevent moisture loss.]

Towelette Disposal: Do not reuse towelette. Dispose of used towelette and empty packet in trash. [Do not flush in toilet.]

[Flat Pack] [or Dispenser] Handling: Non-refillable container. Do not reuse or refill this flat pack. If recycling is not available, discard in trash.

[C – Packet]

STORAGE AND DISPOSAL

Product Storage: [Store at room temperature.]

Towelette Disposal: Do not reuse towelette. Dispose of used towelette and empty packet in trash. [Do not flush in toilet.]

[Packet] [or Dispenser] Handling: Non-refillable container. Do not reuse or refill this packet. Discard packet in trash.

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operator(ies) in the use site tables are interchangeable as applicable.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Pre-Cleaning Instructions:

Use one CaviWipes towelette to completely pre-clean surface of all visible debris.

For use as a disinfectant:

CaviWipes is a two-step disinfectant product. One wipe is required to pre-clean then a second wipe is required to disinfect. Repeated use of the product may be required to ensure that the surface remains visibly wet for 3 minutes at room temperature (68°F/20°C) for Tuberculocidal activity and effectiveness against *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Trichophyton interdigitale* and *Salmonella enterica*. [Use 2-minute contact time to disinfect against SARS-CoV-2*^{††} [and Human Coronavirus*] on hard, non-porous surfaces]. [To kill SARS-CoV-2*^{††} [and Human Coronavirus*] on hard, non-porous surfaces, use 2-minute contact time.] [Pre-cleaning [step] is not required for use against SARS-CoV-2*^{††} [and Human Coronavirus*] on hard, non-porous surfaces].

For use as a Virucide* against HIV-1, HBV, HCV, SARS-COV-2*^{††} [and Human Coronavirus*] applications and effectiveness against Methicillin Resistant *Staphylococcus aureus* (MRSA), Vancomycin Resistant *Enterococcus faecalis* (VRE) and *Staphylococcus aureus* with reduced susceptibility to vancomycin: Use a second CaviWipes towelette to thoroughly wet the surface. Repeated use of the product may be required to ensure that the surface remains visibly wet for 2 minutes at room temperature (68°F/20°C). [Pre-cleaning [step] is not required for use against SARS-CoV-2*^{††} [and Human Coronavirus*] on hard, non-porous surfaces].

[The following direction is to be included when hard, non-porous exterior surfaces of ultrasound transducers or probes is listed as a use site.]

[To disinfect external surfaces of ultrasound transducers: Pre-clean ultrasound gel from external transducer to be disinfected [by using a clean wash cloth]. Then thoroughly wet the probe with a Caviwipes towelette and allow the pre-cleaned surface. Repeated use of the product may be required to ensure that the surface remains visibly wet for 3 minutes at 68°F/20°C.] [Use 2-minute contact time to disinfect against SARS-CoV-2*^{††} [and Human Coronavirus*] on hard, non-porous surfaces]. [To kill SARS-CoV-2*^{††} [and Human Coronavirus*] on hard, non-porous surfaces, use 2-minute contact time.] [Pre-cleaning [step] is not required for use against SARS-CoV-2*^{††} [and Human Coronavirus*] on hard, non-porous surfaces].

[The following statement is to be included on production labels when food contact surfaces are listed.]

[+] A thorough rinse with potable water is required for surfaces in direct contact with food.]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Directions for Use *Continued*

[The following HIV-1 | HBV | HCV instructions are to be included on production labels when claims against HIV-1 or HBV or HCV are listed.]

CaviWipes effectively kills [Human Immunodeficiency Virus Type 1 (HIV-1),] [Hepatitis B Virus (HBV),] [[and] Hepatitis C Virus (HCV)] on precleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects can be associated with the potential for transmission of [HIV-1,] [HBV,] [[and] [HCV].

Special Instructions for Cleaning and Decontamination Against [HIV-1,] [HBV,] [[and] HCV] of Surfaces or Objects Soiled With Blood or Body Fluids:

Personal Protection: When using this disinfecting towelette, wear barrier protection [:] disposable protective gloves, protective gown, face mask and eye covering.

Cleaning Procedure: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of this disinfecting towelette. A **CaviWipes** towelette can be used for this purpose.

Contact Time: Thoroughly wet surface with a new **CaviWipes** towelette. Repeated use of the product may be required to ensure that the surface remains visibly wet for 2 minutes. This contact time will not be sufficient for some of the organisms listed on the label. Refer to label for contact times against other organisms.

Disposal of Infectious Materials: Cleaning materials used that may contain blood or other body fluids must be autoclaved and disposed of in accordance with Federal, State and local regulations for infectious materials disposal.

[The following paragraph will be included in its entirety when the production label references health care use and applications on medical devices.]

This product is not to be used as a terminal sterilant/ high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization/high-level disinfection.

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

DISINFECTION

[CaviWipes is a multi-purpose [disinfectant] [cleaning product] [cleaner] for use on hard, non-porous inanimate surfaces.]

[CaviWipes is effective against the following microorganisms on hard, non-porous surfaces when used according to disinfection directions:]

[This product has demonstrated effectiveness against the following organisms with a 3 minute contact time:]

[This product has demonstrated effectiveness against SARS-CoV-2^{††} [and Human Coronavirus*] on hard, non-porous surfaces with a 2-minute contact time [in the presence of 5% soil load]].

[Disinfects (*insert organisms*) on hard non-porous environmental surfaces listed on this label]

[When used according to disinfection directions, CaviWipes provides effective Tuberculocidal / Virucidal* / Bactericidal / Fungicidal activity against the following organisms:]

Pathogen:	Bacteria	Viruses*	Bloodborne Pathogens**	TB	Fungi	SARS-CoV-2 ^{††}
Contact Time:	3 min[utes]	2 min[utes]	2 min[utes]	3 min[utes]	3 min[utes]	2 min[ute]

ORGANISMS	
	[Strain and/or ATCC number]
<i>Bacteria (3 Minute Contact)</i>	
Staphylococcus aureus	[ATCC 6538]
Salmonella enterica [^^]	[ATCC 10708]
Pseudomonas aeruginosa [^]	[ATCC 15442]
<i>Multidrug-Resistant Bacteria (2 Minute Contact)</i>	
Methicillin Resistant <i>Staphylococcus aureus</i> (MRSA)	[ATCC 33591]
Vancomycin Resistant Enterococcus faecalis (VRE)	[ATCC 51299]
Staphylococcus aureus with reduced susceptibility to Vancomycin	[ATCC 700699]
<i>Mycobacterium (3 Minute Contact)</i>	
<i>Mycobacterium tuberculosis var: bovis (BCG)**</i>	[ATCC 35743]
<i>Mold, Fungi (3 Minute Contact)</i>	
<i>Trichophyton interdigitale</i>	[ATCC 9533]
<i>Viruses [Enveloped]: (2 Minute Contact)*</i>	
Herpes simplex virus [type 1]*	[ATCC VR-260]
Herpes simplex virus [type 2]*	[ATCC VR-734]
[Duck Hepatitis Virus Surrogate for] Hepatitis B Virus (HBV) *, **	[Grimaud, HepadnaVirus Testing]
[Bovine Viral Diarrhea Virus Surrogate for] Hepatitis C Virus (HCV) *, **	[NADL, ATCC VR-1422]
Human Immunodeficiency Virus (HIV-1)*, **	[HIV-1, Zeptomatrix]
Influenza A2 Virus*	[H3N2, A/Hong Kong/8/68 Charles River Laboratories]
Severe Acute Respiratory Syndrome-related Coronavirus 2* -or- SARS-CoV-2 ^{††}	[Strain: USA-WA1/2020] [BEI Resources, NR-52281]
Human Coronavirus [(not associated with Severe Acute Respiratory Syndrome or SARS)]*	[ATCC VR-740] [[Strain] 229E]

* Virus

†† COVID-19 Virus, qualifier will only be added next to SARS-CoV-2 when COVID-19 virus is mentioned on Product Labels and/or Marketing Materials.

^ Pseudomonocidal [Note to Reviewer: the qualifier ^ will only be used when the term Pseudomonocidal is used]

^^ Staphylocidal [Note to Reviewer: the qualifier ^^ will only be used when the term Staphylocidal is used]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operator(ies) in the use site tables are interchangeable as applicable.

[The following graphics are optional for canister containers only. They will be affixed to the label/canister so that it does not interfere with any label text. Due to size constraints, they will not appear on individually wrapped packets.]



THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS

▲ to represent any use site(s) listed in below as applicable.

◆ to represent any surface(s) listed in below as applicable.

- [# pre-saturated towelettes] [X" x Y"]
- [3 minute contact time]
- [2 minute contact -or- kill time against SARS-CoV-2*^{††} [and Human Coronavirus*] [on visibly clean surfaces] on hard, non-porous surfaces]
- [99.9% -or- 3 [three] log reduction against SARS-CoV-2*^{††} [and Human Coronavirus*] on hard, non-porous surfaces [in the presence of 5% soil load] [in [just] [2] [two] minutes] on hard, non-porous surfaces]
- [Kills 99.9% -or- 3 [three] log of SARS-CoV-2*^{††} [and Human Coronavirus*] on hard, non-porous surfaces [in the presence of 5% soil load] [in [just] [2] [two] minutes] on hard, non-porous surfaces]
- [A Metrex Research® brand]
- [A Metrex® brand]
- [A CaviWipes™ Brand]
- [Absorbs] [Absorbent] [wipe]
- [[A] cleaner/disinfectant against SARS-CoV-2*^{††} [and Human Coronavirus*] [in the presence of 5% soil load] [in [just] [2] [two] minutes] on hard, non-porous surfaces]
- [A disinfectant cloth for use in ▲ on ◆]
- [[Aids] [Helps] in reducing -or- in the reduction of cross-contamination between treated surfaces]
- [An effective method to clean against odor-causing organisms on hard, non-porous surfaces listed on this label].
- [A disinfectant towelette for use on ◆ and other ◆]
- [After use, just toss wipe(s) away]
- [Always ready-to-use]
- [An effective method to deodorize and disinfect against odor-causing organisms on hard nonporous surfaces]
- [An effective way to disinfect hard, non-porous environmental surfaces [against SARS-CoV-2*^{††} [and Human Coronavirus*]] [in the presence of 5% soil load] [in [just] [2] [two] minutes] on hard, non-porous surfaces]
- [An effective way to disinfect hard, non-porous environmental surfaces listed on this label.]
- [A non-woven disposable towelettes for use in ▲ and other critical care areas where control between treated surfaces is required]
- [Antibacterial]
- [Antimicrobial]
- [A solution impregnated in a wiping cloth]
- [[Appropriate] [for] Cleaning Patient Discharge]
- [Appropriate] for [Occupied Room] [and] [Terminal Room] [Discharge] Cleaning]
- [Assures proper dilution strength, product effectiveness and standardizes technique]
- [Bacteria -and/or- virus* -and/or- germ killing [formula]]
- [Bactericidal] [Virucidal*] [Fungicidal] [Tuberculocidal] [in 3 minutes]
- [Bactericide -or- Bactericidal] [(in 3 minutes)]
- [Bathroom disinfectant wipes]
- [Big [textured] wipes] [3 minute kill time] [2 minutes [kill/contact time] against/for SARS-CoV-2*^{††} [and Human Coronavirus*] on hard, non-porous surfaces]
- [Broad spectrum disinfectant]
- [Built-in [tough] scrubbing fibers to help clean effectively]
- [CaviWipes have been tested for surface safety on common ▲ surfaces [and equipment] to ensure compatibility with surfaces including: [Stainless Steel] [Glazed Porcelain] [Glass] [Plastic] [Glazed Ceramic]
- [CaviWipes Disinfecting Towelettes are non-woven disposable towelettes pre-saturated with CaviCide.]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [CaviWipes kills germs that cause odors]
- [CaviWipes is effective against odor causing germs]
- [CaviWipes will effectively clean and disinfect when used according to disinfection directions]
- [Can help reduce cross contamination on treated surfaces]
- [Cleaner and disinfectant [towelette] [wipe]
- [Cleaning Formula]
- [Cleans] [◆] [in ▲]]
- [[Cleans and] disinfects ◆]
- [[Cleans] [and] [disinfects] [on the go]]
- [[Cleans and] disinfects [hard, nonporous surfaces]]
- [Cleans, deodorizes and disinfects]
- [Clean Home]
- [Cleans the toughest soils]
- [Cleaning Made Easy]
- [Clinical Clean]
- [[Cleans [and /or Disinfects] [grocery cart handles or shopping cart handles or grocery cart child seats or shopping] cart child seats]
- [Cleans grocery carts and shopping carts]
- [[Cloth-like texture] [Textured] [strong] [for] [with] scrubbing power]
- [Convenient]
- [[Convenient] [Portable] way to [disinfect] [clean] hard non-porous surfaces [of ◆]]
- [Cuts Cleaning time]
- [Deep cleaning] [formula]]
- [Delivers a clean you can see and smell]
- [Designed for Home Healthcare -and/or- patients -and/or caregivers]
- [Designed -or- developed for Home Caregivers]
- [Designed -or- developed for use in [▲] Home Healthcare -or- Home Care]
- [Designed for killing SARS-CoV-2*^{††} [and Human Coronavirus*] on hard, nonporous surfaces [of ◆] [in ▲] [in the presence of 5% soil load] [in [just] [2] [two] minutes]]
- [[Destroys] [Deactivates] [Disinfects] [[Helps] Eliminate 99.9%] SARS-CoV-2*^{††} [and Human Coronavirus*] [of ◆] [in ▲] [in the presence of 5% soil load] [in [just] [2] [two] minutes] on hard, non-porous surfaces]
- [Delivers consistent [amount of solution] concentration every time [!]]
- [Deodorizes, and kills common household germs]
- [Deodorizes by killing the germs that cause odors]
- [Deodorize hard, non-porous surfaces] [◆] [in ▲]]
- [Designed for the ▲ [[industry]]]
- [Designed] [for] [with] [the] Environmental Services Professional[s][in] [mind][!]
- [Disinfectant]
- [Disinfect]
- [Disinfect external surfaces of ◆]
- [Disinfect surfaces [◆]]
- [Disposable]
- [Disinfectant cleaner [that cleans, disinfects and deodorizes]]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Disinfecting formula]
- [Disinfects] [♦] [in ▲]]
- [[Disinfects][Disinfectant][Cleaner] on the go]
- [Disinfects [and deodorizes] ♦]
- [Disinfecting towelettes]
- [Disinfects and deodorizes in 3 minutes]
- [Disinfects hard, non-porous surfaces in 3 minutes]
- [Disinfects [and deodorizes] hard, non-porous surfaces [of ♦] [in ▲] against SARS-CoV-2*^{††} [and Human Coronavirus*] [in [just] 2 minutes] [in the presence of 5% soil load]]
- [[Disinfects] hard surfaces [that people touch]]
- [Disinfects ♦]
- [Disinfects [in [just] three [3] minutes]
- [[Disinfects and deodorizes by killing germs and their odors]
- [[Disinfects hard, nonporous surfaces [throughout the ▲]
- [Disinfect [the] precleaned surface]
- [Dispenses one wipe at a time]
- [Does not contain dyes]
- [DO NOT FLUSH WIPE DOWN TOILET]
- [Double sided-scrubbing wipes]
- [Durable, lint-free wipe designed for sensitive equipment surfaces]
- [Easy]
- [Easy and convenient to use]
- [[Easy] [and convenient] way to clean [and] [disinfect]
- [Easy] pop-up-[tub][canister][pail][lid][tab][dispenser][bucket]
- [Easy [and convenient] way to disinfect [bathrooms][office][nursing stations][high touch] hard, non-porous surfaces everyday]
- [Easy] [convenient] [dispensing] lid [closure][cap][package]
- [Easy to use]
- [Easy to handle]
- [Easy to refill]
- [Easy to use package [canister][container] keeps [your] wipes moist]]
- [Easy way to [clean] [and] [disinfect] [bathroom] [or exterior toilet surfaces] [hard, non- porous surfaces]]
- [Easy and convenient disinfecting [throughout the ▲]]
- [Easy [and convenient] way to [touch-up] [bathroom] [office] everyday]
- [Easy start [tab][lid] [(insert package type)]
- [Easy] [and convenient] [way to clean] [and disinfect]
- [Easy way to clean [and disinfect] [bathroom] [or exterior toilet surfaces] [hard, non- porous surfaces]
- [Economical]
- [Economy]
- [Effective against [Insert any organism from list of organisms]]
- [Effective against Influenza A [Virus]]
- [Effective -or- For daily use against SARS-CoV-2*^{††} [and Human Coronavirus*] [in [just] 2 minutes] [in the presence of 5% soil load] on hard, non-porous surfaces]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Effective against bathroom bacteria and viruses* [use at least two bacteria and two viruses* from this list:] [Salmonella enterica] [Staphylococcus aureus] [Herpes Simplex Virus Type 2] [Hepatitis B Virus (HBV)*] [Hepatitis C Virus (HCV)*]]
- [Effective against -or- Kills a range of bacteria including Staphylococcus aureus, Salmonella enterica, Pseudomonas aeruginosa [in the presence of 5% organic soil [(5% blood serum)]]
- [Effectively disinfects against MRSA]
- [[Effective against -or- Kills (insert any organism(s) from Organisms table) [in the presence of 5% organic soil [(5% blood serum)]]]
- [Effective in the presence of 5% blood serum contamination]
- [Effective against flu virus, Influenza A2 Virus*]
- [[Effectively] disinfects [hard, nonporous] [environmental] surfaces]
- [Effectively removes everyday soils]
- [Efficacy tests have demonstrated that [Product Name] [this product] [are] [is] an effective bactericide, virucide*, tuberculocide, and fungicide in the presence of organic soil (5% blood serum)]
- [Efficacy tests have demonstrated that [Product Name] [this product] [is] effective against SARS-CoV-2*† [in [just] 2 minutes] in the presence of 5% soil load] on hard, non-porous surfaces]
- [Eliminates odors]
- [Eliminates odors from fecal matter and urine]
- [Eliminates [daily] mixing and measuring]
- [Eliminates [lingering] [food+] odors [on surfaces]]
- [[Eliminates odors at their source [caused by germs -or- bacteria]]
- [[Eliminates -or- Reduces odors [caused by germs -or- bacteria]]
- [EVS [Environmental Services] [use]]
- [Exhibits effective virucidal* activity]
- [Extra large [disinfecting] [Y" x Z" compared to our A" x B" size] wipes]
- [Extra large] [Disinfecting Towelettes]
- [Convenient] [portable] way to disinfect hard non-porous surfaces]
- [Kill[s] -and/or- Effective against [99.9% [of]] Salmonella enterica [(Salmonella)]]
- [Kill[s] -and/or- Effective against flu [virus*] [on environmental -and/or hard, nonporous surfaces]]
- [Kill[s] -and/or- Effective against (insert any organism from EPA approved list)]
- [Stop[s] -and/or- help[s] reduce[s] the spread of [99.9% of] bacteria-and/or- viruses* -and/or- germs on hard, nonporous surfaces [in your ▲] [this cold -and/or- flu season]]
- [Kill[s] -and/or- Effective against SARS-CoV-2*† [and Human Coronavirus*] on hard, non-porous surfaces [of ♦] [in ▲] [in [just] [as short as] [2] [two]minutes] [in the presence of 5% soil load] [during this [COVID-19] pandemic -or- outbreak]]
- [For [daily] [general] use [especially] [during this [COVID-19] pandemic -or- outbreak]]
- [For daily use]
- [For] general use
- [For more information visit us at www.metrex.com]
- [For use on bathroom surfaces]
- [For use on nursery surfaces]
- [For use on ♦ surfaces]
- [For use in critical care areas [where control of cross contamination between hard non-porous surfaces is of prime importance]]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [For] Surface cleaning
- [For Use by Home Caregivers]
- [For Use in Home Healthcare]
- [For use in ▲]
- [Formulated -or- Designed to] remove[s] biological odors including Urine -and/or- Feces - and/or-Vomit]
- [Fragrance Free]
- [Freshens]
- [Free of fragrance]
- [Fungicide]
- [Fungicidal]
- [Fungicidal in 3 minutes]
- [Germicidal disposable wipes]
- [Germicide -or- Germicidal]
- [Gets rid of [or eliminates] odors]
- [Great for everyday use]
- [[Heavy] cloth-like towel]
- [[Helps] [kill][s] [eliminate[s] 99.9%] odor-causing bacteria]
- [Helps reduce cross contamination between treated surfaces]
- [Institutional [size]]
- [To clean, [Just] [Simply] wipe [exterior toilet surface] [bathroom germs] [away]]
- [To clean, [Just [pull] wipe and toss]]
- [Kills [bacteria] [viruses*] [fungi]]
- [Kills bacteria, viruses*]
- [Kills [bacteria] [viruses*] [fungi] [on precleaned, hard, nonporous surfaces]
- [Kill germs and bacteria instead of paper towels, sponges and rags [Use this product to]]
- [Kills odor causing bacteria]
- [Kills -or- Effective against bacteria]
- [Kills -or- Effective against germs]
- [Kills -or- Effective against viruses*]
- [Kills [or eliminates 99.9%] of bacteria]
- [Kills [or eliminates 99.9%] of viruses*]
- [Kills [or eliminates 99.9%] of fungus]
- [Kills odor-causing bacteria -or- germs]
- [Kills -or- Effective against [bacteria][and][germs][and][viruses*]
- [[Kills multi drug resistant organisms [(list organisms)]
- [Kills multi drug resistant bacteria [(list organisms)] [in [3] [three] minutes]]
- [Kills] [Effective against] [Insert any organism from list of organisms] [in] [just] [3 minutes] [Kill bacteria in 3 minutes] [Kill viruses* in 2 minutes] [Kill TB in 3 minutes] [Kill [Effective against] SARS-CoV-2*^{††} [and Human Coronavirus*] in 2 minutes on hard, non-porous surfaces]
- [Kills -or- Disinfects -or- Effective against -or- Inactivates SARS-CoV-2*^{††} [and Human Coronavirus*] on hard, non-porous surfaces [of ♦] [in ▲] [in [just] [as short as] [2] [two]minutes] [in the presence of 5% soil load]]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Kitchen+ and Bath Wipes]
- [Large[r] Surface[s]]
- [Large[r] [wipe][s]]
- [Larger wipe covers [more surface area] larger area [with less effort][!]]
- [Lead Free]
- [Leaves bathroom fresh and clean in [just [3 minutes]]]
- [Leaves bathroom fresh and clean in [just] minutes]
- [Leaves [bathroom] [office] clean [and disinfected] [and smelling fresh]]
- [Less] [Reduced] [packaging] [less][to throw away]
- [Lint free] [wipe]
- [Low linting [wipe]]
- [Making home healthcare [manageable] [easy] [simple]]
- [Maximum cleaning, minimum effort]
- [Meets [the disinfection recommendations of] OSHA[’s] Bloodborne Pathogen Guidelines -or- Standards]
- [Moist wipes]
- [Moisture lock lid]
- [Multi-pack]
- [Multi-purpose disinfectant]
- [Multi-surface]
- [Next wipe pops up automatically [everytime] [!]]
- [Neutralizes odors [doesn’t [just] mask them]]
- [No canister needed]
- [No [odors] [or] [fumes]]
- [No [pinched fingers] [rips] [tears]]
- [No sponges, no rags, no buckets, no mess]
- [No threading [required]]
- [No fragrance added]
- [No lint wipe]
- [No measuring] [required]
- [No mixing] [required]
- [No sponges, no rags, no buckets, no mess]
- [Non-scratching [wipe]]
- [Non-woven disposable cloth]
- [[Occupied Room] [and] [Terminal Room] [Discharge] Cleaning Wipe[s] [Towel][s]]]
- [Odor [elimination] [neutralization]]
- [On the go disinfectant]
- [On-the-go disinfectant [Package] [insert type of package] made from [x%] recycled materials]
- [One [wipe] at a time dispensing]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Pre-cleaning [step] is not required for use against SARS-CoV-2*^{††} [and Human Coronavirus*] on hard, non-porous surfaces]
- [[Passes] [Meets] [Satisfies] requirements for killing SARS-CoV-2*^{††} [and Human Coronavirus*] on hard, nonporous surfaces [in [2] [two]minutes] [in the presence of 5% soil load]]
- [Package] [Insert type of package] made from [x%] recycled materials]
- [Picks up [particles] [dust] [dirt] [hair] [crumbs]]
- [Pop-up [dispensing] [top] [tab] [lid]]
- [Portable]
- [Pouch[es]]
- [Powers through tough messes]
- [Pull n/ Go [pack]]
- [Pre-measured]
- [Pre-moistened [towelettes]
- [Premixed], [premeasured], [premoistened]
- [Premoistened [wet] household cleaning wipe]
- [Pre-measured to have the recommended amount of germicidal solution, wipe after wipe]
- [Pre-moistened] [Easy] [and convenient] way to disinfect]
- [Pre-moistened with specially formulated cleaners to handle daily clean-ups without bottles, sponges or buckets]
- [Prevents formulation errors]
- [[Product Name] are useful in a variety of [healthcare], [nursery/daycare] [and] [▲] applications listed on the [this] label.]
- [Pseudomonocidal[^]]
- [Ready to use]
- [[Ready-to-use] disinfectant [for use in ▲]
- [Ready-to-use [2-minute] disinfectant [disinfection] for [killing] SARS-CoV-2*^{††} [and Human Coronavirus*] [in the presence of 5% soil load] on hard, non-porous surfaces]
- [Recommended for cleaning [and disinfecting] ♦.]
- [Recommended for disinfecting exterior toilet surfaces - or - around the exterior toilet surface]
- [Recyclable] [package]
- [Reduce]
- [Reduces] [cross contamination between treated hard non-porous surfaces]
- [Reduces] [packaging] waste [by X%] [when compared to regular canister or bucket or pail]
- [Reduces cross-contamination to (insert any organism from EPA approved list) [in 3 minutes on treated surfaces.]
- [Reduces cross-contamination [from (insert any organism from EPA approved list) [in 3 minutes) on treated (insert non-porous surface)]]].
- [Reduces cross-contamination of SARS-CoV-2*^{††} [and Human Coronavirus*] on treated hard, non-porous surfaces [in the presence of 5% soil load]]
- [Renders SARS-CoV-2*^{††} [and Human Coronavirus*] inactive on hard, nonporous surfaces [in [just] [as short as] [2] [two] minutes] [in the presence of 5% soil load]]
- [Requires no mixing or measuring and are ready-to-use]
- [Re-sealable [stay-moist] pouch [package] [container] inside]
- [Sample]
- [Scrubbing power of a sponge in a disinfecting [towelettes]]
- [Shelf-stable formula]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [So] [surfaces are more than [just] [clean] [they're] [disinfected]
- [Specifically -or- Specially designed for the Medical/Dental/ Veterinary/Health Care [▲] industry]
- [Staphylocidal ^^]
- [Starts cleaning on contact]
- [Strong, durable low-lint wipe [designed for sensitive equipment surfaces]]
- [[Strong,] [Durable,] cloth-like wipes]
- [Strong, durable textured wipe [cleans] [and] [deodorizes]]
- [[Strong,] [Durable,] non-woven disposable cloth]
- [[Strong,] [Durable,] non-woven wipe does not leave lint on surfaces]
- [Suitable for [▲]]
- [Suitable for use on [◆]]
- [Suitable for use on [medical grade] stainless steel]
- [Suitable for use on hard, nonporous inanimate environmental surfaces]
- [Strong cleaning power]
- [Sturdy][Container][Dispenser][Canister][Pail][Bucket]
- [Tear at perforation]
- [Textured for better [cleaning] [scrubbing] power]
- [Textured design allows for] [scrubbing] [cleaning] [removing] [wiping] [tough] [stains] [soils]
- [messes]
- [Textured wipe for effective cleaning]
- [Textured] wipe is low-linting [designed] [for sensitive equipment surfaces]
- [Textured design [cleans] [deodorizes] hard, non-porous surfaces] [helps remove soil from surfaces]]
- [Textured] [strong] [towelettes] [for] [with] [scrubbing power]
- [The 3-minute [germicidal wipe] [disinfectant wipe]
- [The 2-minute [germicidal wipe] [disinfectant wipe] against SARS-CoV-2*^{††} [and Human Coronavirus*] [in the presence of 5% soil load] on hard, non-porous surfaces]
- [The answer to your [disinfecting] [cleaning] needs]
- [[The] smell of clean]
- [[The] easy -and/or- convenient way to disinfect] [with CaviWipesM]
- [This product] controls the spread of germs on most hard, nonporous treated surfaces
- [This product is -or- has been tested according to AOAC efficacy testing requirements -or standards for hospital disinfection]
- [This product has demonstrated effectiveness against SARS-CoV-2*^{††} [and Human Coronavirus*] in 2 minutes [in the presence of 5% soil load] on hard, non-porous surfaces]
- [[Thick], [durable] wipes are bactericidal, tuberculocidal, and virucidal* in 3 minutes]
- [Thick, large wipe [cleans] [disinfects]]
- [This] [outer] [selling] [package] [container] [canister] contains X% post-consumer [material] [plastic] [resin]
- [Time to reorder]
- [Trial]
- [Tuberculocidal (3 minutes)]
- [Tub with pop-up action]
- [Virucide* -or- Virucidal* against SARS-CoV-2*^{††} [and Human Coronavirus*] [in [just] 2 minutes] [in the presence of 5% soil load] on hard, non-porous surfaces]
- [[Wipes away -or- out 99.9%] [Kills] SARS-CoV-2*^{††} [and Human Coronavirus*] [in the presence of 5% soil load] [in [just] 2 minutes] on hard, non-porous surfaces]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

CLEANING AND STAIN REMOVAL

To clean -and/or- remove stains: Clean surface with product; [rinse and] let air dry. [If streaking is observed, wipe with a clean, damp [cloth] [or] [paper towel].] [[Works on] [effective against] [removes the following stains:]]

Claims: <ul style="list-style-type: none">• Cleans• Removes stains• Removes (<i>insert stain(s)/soil(s) shown to the right:</i>)	Stains Soils <ul style="list-style-type: none">• Blood• Body oils• Bodily fluids• Common soils• Dead skin• Dirt• Fingerprints• Food+ residue[s] -or- soil[s]• Food stains• Grease• Grime• Laboratory stains	Stains Soils continued <ul style="list-style-type: none">• Mucus• [Other] common soils -and/or- stains• [Other] organic matter• Phlegm• Protein-based matter• Saliva• Spit• Spittle• Sputum• Stool• Throw-up• Urine• Vomit
---	---	---

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 1 Medical Use Sites:
CaviWipes is recommended for use on hard non-porous surfaces in [▲]

3D Imaging Services Centers	Eye Care [Facilities][Centers]	Optometric [Facilities] [Offices]
Acute Care Facilities	Eye Surgical Centers	Oral Health Care Facilities
Alternate Care Facilities	Forensic Labs	Orthodontic Offices
Ambulatory Surgery Centers (ASC)	Full-Care Nursing Facility	Orthopedics
Animal Research Facilities	Gastrointestinal (GI) [Labs] [Laboratories]	Outpatient Care Facilities
Assisted Living Facility	Gastrointestinal (GI) Centers	Outpatient [Surgery][Facilities] [Centers]
Autopsy Rooms	Gastrointestinal (GI) Suites	Out-Patient Surgical Centers (OPSC)
Bathroom	Home[health] Care Services	Pathology Rooms
Bioskill Laboratories	Hospice	Patient Care Areas
Blood Drawing Centers	Hospitals	Patient Rooms
Blood Apheresis Facilities	Hospital Operating Rooms	Pediatric Exam Rooms
Blood Units	Hydrotherapy [Facilities][Rooms]	Pediatric Intensive Care Units (PICU)
Burn Center	In vitro Fertilizer Center (IVF)	Pediatric Offices
Cardiac Intensive Care Unit (CICU)	Instrument Reprocessing [Facilities] [Departments] [Center]	Personal Care Service Facilities
CCU	Instrument Reprocessing [Rooms]	Pharmacy Clean Rooms
CAT Scan [Lab] [Laboratories]	Intensive Care Units (ICU)	Phototherapy Rooms
Cath Labs	Isolation Areas	Physical Therapy Facilities
Central Service Areas	Isolation Rooms	Physicians' Offices
Central Services Department	Labor Units	Plasma Apheresis [Centers][Facilities][Units]
Central Supply Areas	[Lab] [Laboratories]	Post [Ops][Operatory] Rooms
Central Supply Rooms	Lasik Eye Centers	Pre[Ops] [Operatory] Areas
Clean Rooms	Lasik Units	Psychiatric Care Facility
Clinical Areas	Laundry Rooms	Radiology Department
Clinics	Life Science Laboratories	Recovery Rooms
Critical Care Units	Long-Term Acute Center (LTAC)	Rehabilitation Facilities
Computer Science Labs	Long-Term Care [Facilities]	Respiratory Therapy Facilities
CS [Labs] [Laboratories]	Medical Clinics	Skilled Nursing [Facility] [Homes]
Dialysis Centers	Medical Facilities	Environmental Surfaces of Sterile Processing Departments
Dialysis Clinics	Medical Related Facilities	Environmental Surfaces of Sterile Storage Areas
Dissection [Labs] [Laboratories]	Medivans	Surgery Centers
Doctor's Offices	Medical Intensive Care Unit (MICUs)	Surgery Rooms
Donation Centers for [blood] [plasma] semen milk aphaeresis	Magnetic Resonance Imaging (MRI) [facilities] [laboratories] [labs]	Surgical Centers

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 1 Medical Use Sites Continued:

Drug Preparation [Centers] [Facilities] [Rooms]	ENT Facility	Out-patient surgical centers (OPSC)
Drug Preparation Rooms	Hospitals	Operating Rooms
Electrolysis Facilities	Neonatal Intensive Care Units (NICU)	Ophthalmic [Facilities] [Offices]
Emergency Departments	Neonatal Nursery Areas	Surgical Suites
Emergency Medical Services (EMS) Units	Newborn Nursery Center	Teaching Hospitals
Emergency Medical Settings	Nuclear Medicine Facilities	Tissue Banks
Emergency Rooms (ER)	Nurses' Stations	Transplant Center
Emergency Vehicles	Nursing Homes	Trauma Center
Endo Suites	Nursing Stations	University Research Lab
Endoscopy [Centers] [Facilities] [Departments]	Occupational Therapy [Centers] [Facilities]	Wound Center
Endoscopy Centers	Oncology Areas	X-Ray Rooms

Table 2 Institutional Use Sites:

CaviWipes is recommended for use on hard non-porous surfaces in [▲]

Adult Schools	Early Learning Centers	Middle Schools
Air Transport Facilities	Elementary Schools	Military Installations
Airplane Interiors	[Fire] [Rescue Departments]	Nursery Schools
Airports	Fire Vehicles	Police Vehicles
Bathrooms	Forestry Services	Pre-K
Casinos	Grammar Schools	Preparatory Schools
Child Care Facilities	High Schools	Preschools
Churches	Hospitality Facilities	Primary Schools
Colleges	Hotels	Prisoner Detention Facilities
Community Colleges	Industrial Facilities	Prisoner Lockups
Conference Centers	Infant Care	Prisons
Conference Rooms	Intermediate Schools	School Buses
Convention Facilities	Jails	School Kitchens+
Correctional Facilities	Junior High Schools	Schools
Daycare Centers	Kindergartens	[University] [Universities]
[Dormitories] [Dorms]	Learning Institutions	

Table 3 Dental Use Sites:

CaviWipes is recommended for use on hard non-porous surfaces in [▲]

Dental [Laboratories] [Lab]	Dental Office	Dental [Operatory] [Operatories]
-----------------------------	---------------	-------------------------------------

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 4 Vet Use Sites:

CaviWipes is recommended for use on hard non-porous surfaces in [▲]

Animal Care Facilities	Avian Units	Pet Homes
Animal Control Vehicles	Bathroom	Pet Motels
Animal Dental Operatory	Doggy Day Care Centers	Pet Spas
Animal Detention	Feline Facilities	Pet Stores
Animal Research Facilities	Feline Units	Small Animal Facilities
Animal Shelters	Kennel	Veterinary
Avian [Care] Facilities	Pet Areas	Veterinary Facilities

Table 5 Food Prep+ Use Sites:

CaviWipes is recommended for use on hard non-porous surfaces in [▲]

Bakeries+	Food Processing Facilities+	Restaurants+
Banquet Hall Facilities+	Grocery Stores +	School Kitchens+
Butcher Shops+	Kitchen+	Supermarkets+
Food Preparation Areas+	Produce Areas+	

[+ A thorough rinse with potable water is required for surfaces in direct contact with food.]

Table 6 Miscellaneous Use Sites:

CaviWipes is recommended for use on hard non-porous surfaces in [▲]

Airplanes	Factories	Money Handling Facilities
Arboretums	Foots Spas	Morgues
Arcades	Foreign Current Exchange Facilities	Mortuaries
[Athletic] [Training] Facilities	Funeral Homes	Motels
Auto Facilities	Gaming Facilities	Movie Theaters
Banks	Garden Centers	Outdoor Playground
Barber Shops	General Public Facilities	Personal Care Service Facilities
[Bathroom(s) or public Restrooms]	Golf Courses	Public Areas
Beauty Salons	Grounds Keeping Services	Public Facilities
Body Piercing Salons	Gymnasium	Public Transportation Vehicles
Body Service Facilities	Health Club Facilities	Recreational [center] [facilities]
Body Storage Facilities	Helicopter Interiors	Recycling Centers
Botanical Gardens	Household Bathrooms	Salon Settings
Bowling Alleys	Housekeeping Rooms	Shelters
Buses	Indoor Playground	Shower Rooms
Cabs	Laundromats	Skincare Salons
Cadaver Processing areas	Laundry Rooms	Storage [Areas] [Rooms]
Cars	Lawn Services	Subway
Centers	Limousine	Tanning Salons
Condominiums	Living Rooms	Tattoo Parlors
Cruise Ships	Locker Rooms	Taxi
Currency Exchange	Lockup Facilities	Theme Parks
Day Spa	Manicure Salons	Transport Vehicles
Delivery Units	Manufacturing [Facilities] [Plant]	
Facilities	Massage Facility	

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces:
CaviWipes is recommended for use on hard non-porous surfaces of [♦]

Ambulances [Equipment Surfaces]	Bathroom Fixtures	Carts
Anesthesia Equipment	Bathtubs	Cash Register
Anesthesia Machines	Bed Pans	CAT Scan Machines
Animal Cages	Bed Railings [Bedside Rails]	Ceilings
Animal Equipment	Bedside Commode	Cell Phones
Aphaeresis Machines	Behind and Under Sinks	Center Venous Lines
Appliances+	Bidet	Chairs
Arm Rests	Biohazard Container	Changing Tables
ATM Surfaces	Bladder Scan	Chemotherapy Hoods
Autoclave Surfaces	Blood Glucose Meters	Child Seats+
Automatic Feeders	Blood Pressure Cuff	Chrome
Automatic Teller Machines (ATM)	Booster Chairs+	Clipper Blades
Baby Cribs	Booster Seats+	Closet Handles
Baby Warmers	BP Monitors	Coated Mattresses
Baked Enamel	Brushes	Coated Pillows
Bars	Cabinet Handles	Combs
Basins	Cabinets	Compression Devices
Baskets	Cafeteria Trays+	Computer [Keyboards] [Mouse]
Bassinet	Cages	Computers
Bathing Units	Cardiac Gym Equipment	Copy Machines
<ul style="list-style-type: none"> • [For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [▲surfaces] [and other critical care areas where control of cross contamination between hard nonporous treated surfaces is important. [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.] • [Surfaces which are made of plastic [(i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene)], and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.] • + A thorough rinse with potable water is required for surfaces in direct contact with food. • + Any washable (food and non-food contact) surface where disinfection is required 		

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**Table 7 Surfaces Continued:
CaviWipes is recommended for use on hard non-porous surfaces of [♦]**

Counters	Diaper Changing Stations	Examination Table
Countertops+	Diaper Pails	Exterior Surface of Digital X-Ray Sensors
Cribs+	Dictating Equipment	Exterior Surfaces of Amalgamators
Crutches	Digital Phosphor Plates	Exterior Surfaces of Appliances such as Refrigerators+μ, Microwaves+μ, Food Processors+ and Stoves+μ
Defibrillators	Dining Tables+	Exterior surfaces of CICU Equipment
Dental Care Facilities Hand Pieces	Dish Racks+	Exterior Surfaces of Dental Operatories
Dental Chairs	Display Cabinets	Exterior Surfaces of Dental Equipment
Dental Countertops	Diving Equipment	Exterior Surfaces of Diagnostic Equipment
Exterior Surfaces of Dental Equipment Lines	Diving Masks+	Exterior Surfaces of Dialysis Machines
Dental Equipment Tubing	Diving Snorkels+	Exterior surfaces of EKG Machines
Dental Impression Trays	Doorknobs	Exterior surfaces of Electronic Thermometers
Dental Operatory Surfaces	Dopplers	Exterior surfaces of Electrosurgical Generators
Dental Unit Instrument Trays	Drain Boards+	Exterior surfaces of Emergency Resuscitative Equipment
Dentist Chair	Drained Whirlpool Tank Surfaces	Exterior Surfaces of Endodontic Equipment (Apex Locators, Pulp Testers, Motors)
Desk Tops	Emergency Vehicles	Exterior surfaces of Endodontic equipment such as apex locators
Desks	Endoscope Cabinet	Exterior surfaces of ICU Equipment
Dialogistic Equipment	Exam Tables	Exterior surfaces of Incubation Equipment
Dental Curing Lights¥		
<ul style="list-style-type: none"> • [For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [▲ surfaces] [and other critical care areas where control of cross contamination between hard nonporous treated surfaces is important. [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.] • [Surfaces which are made of plastic [(i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene)], and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.] • + A thorough rinse with potable water is required for surfaces in direct contact with food. • + Any washable (food and non-food contact) surface where disinfection is required. • μ Allow surfaces to come to room temperature prior to treatment. • ¥ For surfaces that are not inserted into the human mouth. 		

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**Table 7 Surfaces Continued:
CaviWipes is recommended for use on hard non-porous surfaces of [♦]**

Exterior surfaces of Incubators	Exterior surfaces of Patient Care Monitor	Exterior surfaces of Ultrasound Monitors
Exterior surfaces of Infusion Pumps	Exterior surfaces of Pneumatic Tourniquets	Exterior Surfaces of Vanity Tops
Exterior surfaces of Isolettes	Exterior surfaces of Probes	Exterior Toilet Surfaces
Exterior surfaces of IV Pumps	Exterior surfaces of Pulp Tester Motors	External Surfaces of Glucometers
Exterior surfaces of Kendrick Extrication Device (KED)	Exterior surfaces of Pulse Oximeter	External Surfaces of Ultrasound Transducers
Exterior surfaces of Life Support Equipment	Exterior surfaces of Radiology Equipment	External Surfaces of Urinals
Exterior surfaces of Medical Equipment	Exterior Surfaces of Refrigerator Units+	Exterior Surfaces of Veterinary Equipment
Exterior surfaces of Medical Warmers	Exterior Surfaces of Respiratory Therapy Equipment	Eye Glasses (do not use on contact lenses)
Exterior surfaces of Microwaves+	Exterior surfaces of Respiratory Ventilators	Faucets+
Exterior surfaces of Monitors	Exterior surfaces of Resuscitators	Fax Machines
Exterior surfaces of Monitors Electrodes	Exterior surfaces of Sequential Compression Devices	Feed Racks+
Exterior surfaces of MRI Machines	Exterior or Exposed Interior Surfaces of Medical Equipment	Filing Cabinets
Exterior surfaces of OLED Monitors	Exterior surfaces of Sterilizers	Filing [Doors][Handles
Exterior surfaces of Ophthalmoscopes	Exterior surfaces of Tympanic Thermometers	Fire [Vehicles] [Trucks]
Exterior surfaces of Oscopes	Exterior surfaces of Ultrasound Machines	Floors
<ul style="list-style-type: none"> • [For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [▲surfaces] [and other critical care areas where control of cross contamination between hard nonporous treated surfaces is important. [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.] • [Surfaces which are made of plastic [(i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene)], and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.] • + A thorough rinse with potable water is required for surfaces in direct contact with food. • + Any washable (food and non-food contact) surface where disinfection is required. 		

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces Continued:
CaviWipes is recommended for use on hard non-porous surfaces of [♦]

Food Cases+	Hair Cutting Implements	IV Stands
Food Service Countertops	Hair Dryers	Jacuzzi
Food Serving Trays+	Hair Rollers	Kennels
For use as a cleaner and disinfectant on cured, set and hardened, hard nonporous [dental impressions] [dental impression materials]	Hampers	Keyboards
For use as a cleaner and disinfectant on hard nonporous [dental impressions] [dental impression materials] prior to sending to a dental lab.	Handles	Knee Walkers
Formica	Handrails	Knobs
Fountains+	Hard non-porous surfaces of vehicles including steering wheel	[Lab]oratory Equipment
Freezers+ μ	Hard Plastic or Metal Flooring	[Lab]oratory Surfaces
Garbage Cans	Hard Surface Weight Lifting Surfaces	Laminated Surfaces
Glass	Hard, non-porous environmental surfaces	Lamps
Glazed Ceramic Tile	Headsets	Light Lens Covers
Glazed Porcelain	High Chairs+	Lights
Glazed Tile Walls	High Touch Objects	Linoleum
Glazed Tiles	Hoods+	Loupes
Golf Cart	Hospital Bed (Railings, Frames)	Manicure Implements
Grab Bar	Infant Incubators (interior and exterior surfaces)	Manicure Tables
Grocery Carts	Infant/Child Care Equipment Surfaces	Marlite
Gurneys	Inside of Drawers	Mayo tray
Gym Equipment	Instrument Stands	Mayo tray table
Hair Clippers	IV Poles	Nail/Hair Care Implements
<ul style="list-style-type: none"> • [For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [▲surfaces] [and other critical care areas where control of cross contamination between hard nonporous treated surfaces is important. [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.] • [Surfaces which are made of plastic [(i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene)], and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.] • + A thorough rinse with potable water is required for surfaces in direct contact with food. • +Any washable (food and non-food contact) surface where disinfection is required. • μ Allow surfaces to come to room temperature prior to treatment. 		

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces Continued:
CaviWipes is recommended for use on hard non-porous surfaces of [♦]

Non-food contact surfaces+	Police [Vehicles][Cars]	Silk Peel™ with Dermallnfusion™
Nurse Call [Button] [Device] [Cord]	Porcelain Enamel	Sinks+
Operating Room Furniture	Plastic Posey Vest	Slit Lamps
Operating Room Lights	Public Restrooms	Smart Phone Chargers
Operating Room Tables	Razors	Spine Back Boards
Operating Tables	Reception Areas	Stainless Steel
Operatory Light Switches	Reception Counters	Stalls
Outdoor Furniture+ (excluding wood frames and upholstery)	Refrigerator Carts+μ	Stand
Ovens+μ	Refrigerators+μ	Stethoscopes
Overbed table	Rubber	Stools
Oxygen Hoods	Safety Straps	Storage unit
Patient Delivery Equipment	Salad Bar Sneeze Guards+	Stoves+μ
Patient Lift Equipment	Salon Surfaces	Stovetops+μ
Patient Monitoring Equipment	School Buses	Stretchers
Patient Scales	Scissors	Suction Pumps
Patient Support Equipment	Seats	Support Bars
Pens	Shampoo Bowls	Surgery Microscopes
Phone Chargers	Shears	Surgery Robots
Phone Cradle	Shelves	Surgical Headlights
Physical Therapy (PT) equipment surfaces	Shopping Carts	Surgical Light Sources
Pipes	Shower Chair	Synthetic Marble
Plastic [Laminate]	Shower Fixtures	Tables+
Plastic Mattress Covers	Shower Stalls	Tank
Playpens	Showers	Tanning Beds
<ul style="list-style-type: none"> • [For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [▲surfaces] [and other critical care areas where control of cross contamination between hard nonporous treated surfaces is important. [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.] • [Surfaces which are made of plastic [(i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene)], and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.] • + A thorough rinse with potable water is required for surfaces in direct contact with food. • +Any washable (food and non-food contact) surface where disinfection is required. • μ Allow surfaces to come to room temperature prior to treatment. 		

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces Continued:
CaviWipes is recommended for use on hard non-porous surfaces of [♦]

Telemetry Equipment	Trash Cans	Walls
Telephones	Tray Tables	Warmers
Tiles+	Troughs	Wash Basin
Toilet (Rim, Seat, Tank, Handle,)	Tubing	Washable Nail Files
Touch Screen Monitors	Vanity Tops+	Watering Appliances
Towel Dispensers	Vaporizer	Weight-lifting Surfaces
Toys (Plastic or hard nonporous)	Veterinary Care Surfaces	Wheelchairs
Transfer Bench	Vinyl	Whitening Lamps
Workstations	Walkers	
<ul style="list-style-type: none"> • [For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [▲ surfaces] [and other critical care areas where control of cross contamination between hard nonporous treated surfaces is important. [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.] • [Surfaces which are made of plastic [(i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene)], and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.] • + A thorough rinse with potable water is required for surfaces in direct contact with food. • +Any washable (food and non-food contact) surface where disinfection is required. 		