

46781-11

11/28/2001

Page 134

[Outside container]

NOV 28 2001

**PREMIDYNE**  
Liquid Medical Waste Sanitizer  
For the Disposal of Suction Canister Contents

46781-11

The efficacy of this product has been demonstrated in the presence of 10% to 50% whole blood.

**ACTIVE INGREDIENT:**

Sodium Dichloro-s-Triazinetrione Dihydrate.....99%

**OTHER INGREDIENTS:**.....1%

**Total**.....100%

Available Chlorine: 55.5%

**KEEP OUT OF REACH OF CHILDREN  
DANGER**

<b>FIRST AID</b>	
<b>If in eyes</b>	<ul style="list-style-type: none"> <li>• Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.</li> <li>• Call a poison control center or doctor for treatment advice.</li> </ul>
<b>If inhaled</b>	<ul style="list-style-type: none"> <li>• Move person to fresh air.</li> <li>• If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably mouth-to-mouth if possible.</li> <li>• Call a poison control center or doctor for further treatment advice.</li> </ul>
<b>If on skin or clothing</b>	<ul style="list-style-type: none"> <li>• Take off contaminated clothing.</li> <li>• Rinse skin immediately with plenty of water for 15-20 minutes.</li> <li>• Call a poison control center or doctor for treatment advice.</li> </ul>
<b>If swallowed</b>	<ul style="list-style-type: none"> <li>• Call a poison control center or doctor immediately for treatment advice.</li> <li>• Have person sip a glass of water if able to swallow.</li> <li>• Do not induce vomiting unless told to do so by a poison control center or doctor.</li> <li>• Do not give anything by mouth to an unconscious person.</li> </ul>
For emergency information on product, call Metrex at 1-800-841-1428, Monday through Friday, between 6 a.m. and 4 p.m. Pacific Time. After 4 p.m., call Chemtrec 24-Hour Emergency Service at 1-800-424-9300	
Have the product container or label with you when calling a poison control center or doctor, or going for treatment.	
<b>NOTE TO PHYSICIAN:</b> Probable mucosal damage may contraindicate the use of gastric lavage.	

See side panel for additional precautionary statements

**PRECAUTIONARY STATEMENTS**

Hazards to Humans and Domestic Animals

**DANGER**

Corrosive. Causes irreversible eye damage. May be fatal if inhaled. Harmful if swallowed or absorbed through the skin. Do not get in eyes, on skin or on clothing. Do not breathe dust, vapor, or spray mist. Wear goggles, face shield, or safety glasses. Wash thoroughly with soap and water after handling. Remove and wash contaminated clothing before reuse.

**DIRECTIONS FOR USE**

It is a violation of federal law to use this product in a manner inconsistent with its labeling.

**The efficacy of this product has been demonstrated in the presence of 10% to 50% whole blood.**

Premidyne can be used in all areas of medical facilities where the suctioning of fluids may occur, *e.g.*, operating rooms, emergency rooms and acute care areas, for the disposal of suction canisters contents.

Premidyne should be introduced 30 minutes before disposal of the suction canister contents.

Before the medical procedure begins, open the access port and insert the water soluble unit of Premidyne into the suction canister. Re-cap the suction canister. Begin the medical procedure. Activation of suction will provide thorough agitation. Following the procedure, transport the suction canister to a hopper and dispose of the contents according to institutional procedures and Federal, State and local regulations. Dispose of the empty suction canister according to institutional procedures and Federal, State and local regulations.

This package contains 35 water soluble units of the size checked below.

- 4 g (0.14 oz.) each unit sanitizes 500 ml of suction canister waste
- 9.6 g (0.34 oz.) each unit sanitizes 1200 ml of suction canister waste
- 12 g (0.42 oz.) each unit sanitizes 1500 ml of suction canister waste
- 16 g (0.56 oz.) each unit sanitizes 2000 ml of suction canister waste
- 24 g (0.85 oz.) each unit sanitizes 3000 ml of suction canister waste

**STORAGE AND DISPOSAL**

Do not contaminate water, food or feed by storage or disposal.

**Pesticide Storage:** Store in a cool, dry, well-ventilated area.

**Pesticide Disposal:** Dispose of waste according to Federal, State and local regulations.

**Container Disposal:** Do not reuse container. Triple rinse. Then offer for recycling or reconditioning, or puncture and dispose of in a sanitary landfill, or incineration or, if allowed by state and local authorities, by burning. If burned, stay out of smoke.

3 7 4

For product or technical information, contact Metrex Customer Care at 800-841-1428 or visit our website at [www.metrex.com](http://www.metrex.com).

EPA Reg. No. 46781-11

EPA Est. No. 71441-OH-002

Packed for:  
Metrex Research Corporation  
28210 Wick Road  
Romulus, MI 48174

Net Contents:

This package contains 35 water soluble units of the size checked below.

- Reorder No. xxxx  
Net Wt. 4 g (0.14 oz.)
- Reorder No. xxxx  
Net Wt. 9.6 g (0.34 oz.)
- Reorder No. xxxx  
Net Wt. 12 g (0.42 oz.)
- Reorder No. xxxx  
Net Wt. 16 g (0.56 oz.)
- Reorder No. xxxx  
Net Wt. 24 g (0.85 oz.)

[ Water soluble unit packaging label for the 4 g, 9.6 g, 12 g, 16 g and 24 g packets]

**PREMIDYNE**

The efficacy of this product has been demonstrated in the presence of 10% to 50% whole blood.

**ACTIVE INGREDIENT:**

Sodium Dichloro-s-Triazinetrione Dihydrate.....99%

**OTHER INGREDIENTS:**.....1%

**Total**.....100%

Available Chlorine: 55.5%

KEEP OUT OF REACH OF CHILDREN

**DANGER**

EPA Reg. No. 46781-11

EPA Est. No. 71441-OH-002

See outside container for additional information.

Net Wt. xx g (xx oz.) sanitizes xxxx ml