



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
WASHINGTON, DC 20460

OFFICE OF CHEMICAL SAFETY
AND POLLUTION PREVENTION

November 10, 2021

Rebecca Lei
Agent
Metrex Research
1717 West Collins Avenue
Orange, CA 92867

Subject: Label and CSF Amendment – Addition of Pathogens and Revision of CSFs
Product Name: CaviCide
EPA Registration Number: 46781-12
Received Date: 5/24/2021
Action Case Number: 00303011

Dear Rebecca Lei:

The amended label and CSF(s) referred to above, submitted in connection with registration under the Federal Insecticide, Fungicide and Rodenticide Act (FIFRA), as amended, are acceptable. This approval does not affect any conditions that were previously imposed on this registration. You continue to be subject to existing conditions on your registration and any deadlines connected with them.

A stamped copy of your labeling is enclosed for your records. This labeling supersedes all previously accepted labeling. Pursuant to 40 CFR 156.10(a)(6) you must submit one copy of the final printed labeling before you release the product for shipment with the new labeling. In accordance with 40 CFR 152.130(c), you may distribute or sell this product under the previously approved labeling for 18 months from the date of this letter. After 18 months, you may only distribute or sell this product if it bears this new revised labeling or subsequently approved labeling. “To distribute or sell” is defined under FIFRA section 2(gg) and its implementing regulation at 40 CFR 152.3.

Please note that the record for this product currently contains the following CSF(s):

- Basic CSF dated 10/07/2021
- Alternate CSF 1 dated 10/07/2021


Should you wish to add/retain a reference to the company’s website on your label, then please be aware that the website becomes labeling under FIFRA and is subject to review by the Agency. See FIFRA section 2(p)(2) If the website is false or misleading, the product would be misbranded and unlawful to sell or distribute under FIFRA section 12(a)(1)(E). 40 CFR 156.10(a)(5) lists examples of statements EPA may consider false or misleading. In addition, regardless of whether a website is referenced on your product’s label, claims made on the website may not substantially differ from those claims approved through the registration process,

Page 2 of 2
EPA Reg. No. 46781-12
Action Case No. 00303011

FIFRA section 12(a)(1)(B). Therefore, should the Agency find or if it is brought to our attention that a website contains false or misleading statements or claims substantially differing from the EPA approved registration, the website will be referred to the EPA's Office of Enforcement and Assurance.

Your release for shipment of the product constitutes acceptance of these conditions. If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA section 6. If you have any questions, please contact Perri Moeller by phone at (202) 566-0813, or via email at moeller.perri@epa.gov.

Sincerely,

A handwritten signature in black ink that reads "Steven Snyderman". The signature is written in a cursive style with a horizontal line underneath the name.

Steven Snyderman, Product Manager 33
Regulatory Management Branch II
Antimicrobials Division (7510P)
Office of Pesticide Programs

Enclosure: Annotated Label

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

CaviCide 1™

[TUBERCULOCIDAL] [BACTERICIDAL] [VIRUCIDAL*^[^]] [FUNGICIDAL]

ACTIVE INGREDIENTS

Didecyldimethylammonium chloride0.76%
Ethanol 7.50%
Isopropanol.....15.00%
[Inert] Other Ingredients†.....76.74%
Total.....100.00%

KEEP OUT OF REACH OF CHILDREN

WARNING

[See other panel for additional precautionary statements.]

EPA Reg. No. 46781-12

EPA Est. No. [insert EPA Est. No of Producer List on Box 2 of CSF]
[56952-WI-001] [75272-WI-001]
[71979-SC-001] [51435-WI-004]
[56952-CHN-001] [73027-WI-001]
[46781-MI-001] [93725-TN-001]
[90532-WI-001] [90532-WI-002]
[90532-WI-003] [0312-WI-3]

[Beginning of batch code indicates Est No.]

[End of batch code indicates Est No.]

Made in USA

[CaviCide 1 is available in the following packaging configurations]:

- 4 gallons/CS
- 48 – 2 oz bottles/CS
- 12 – 8 oz bottles/CS
- 12 – 24 oz bottles/CS
- 2 – 2.5 gallons/CS
- 1 – 55 gallon drum/CS
- 1 - 350 gallon transport tote]

[1 minute]

[1 Step Germicidal Spray]

[Fragrance Free]

[Ready-To-Use Surface Disinfectant Cleaner]

Reorder No. XXXXXX

Net Contents:

X [[Gallon][gal]] [[Ounce][oz]] [[Fluid Ounce][fl. oz]]/ X [[Liter][L]] [[Milliliters][mL]]
(Net Weight x lbs./ x kg.)
(Net Weight x lbs/ x g)

[Packed For:]

[Distributed By:]

[Sold By:]

[Manufactured By]

Metrex Research
28210 Wick Road
Romulus, MI 48174

[To reorder] [,] [For product or technical information] contact Metrex at 800-841-1428 or visit our web site at

www.metrex.com

[For SDS information please visit

www.metrex.com

ACCEPTED

11/10/2021

Under the Federal Insecticide, Fungicide and Rodenticide Act as amended, for the pesticide registered under EPA Reg. No. 46781-12

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals

WARNING Causes substantial but temporary eye injury. Do not get in eyes or on clothing. Wear [appropriate protective eyewear such as] goggles, face shield, or safety glasses. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse. [See other panel for additional precautionary statements.]

Physical or Chemical Hazards: Combustible. Do not use or store near heat or open flame.

FIRST AID

If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, and then continue rinsing. Call a poison control center or doctor for treatment advice.

For Chemical Emergency ONLY (spill, leak, fire, exposure or accident), call CHEMTREC at [insert phone number]. For [ALL] [other]] inquiries about this product, call Metrex Customer [Service] [Care] at [insert phone number], Monday through Friday, between 6 [or any time] a.m. and 4 [or any time] p.m. Pacific Time.

Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

[The following Environmental Hazards statement is to be included only on labeling of containers that are 5 gallons and larger.]

Environmental Hazards: This pesticide is toxic to fish. Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office of the EPA.

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

EMERGING VIRAL PATHOGEN CLAIMS

This product qualifies for emerging viral pathogen claims per the EPA’s ‘Guidance to Registrants: Process for Making Claims Against Emerging Viral Pathogens not on EPA-Registered Disinfectant Labels’ when used in accordance with the appropriate use directions indicated below.

This product meets the criteria to make claims against certain emerging viral pathogens from the following viral category:

- Enveloped Viruses
- Large Non-enveloped
- Small Non-enveloped

| | |
|---|---|
| <i>For an emerging viral pathogen that is a/an...</i> | <i>...follow the directions for use for the following organisms on the label:</i> |
| Enveloped virus | Rotavirus |
| Large, non-enveloped virus | Feline Calicivirus |

[Product name] has demonstrated effectiveness against viruses similar to [name of emerging enveloped/small and large non-enveloped virus] on hard, non-porous surfaces. Therefore, [product name] can be used against [name of emerging enveloped/small and large non-enveloped virus] when used in accordance with the directions for use against (cited above) on hard, non-porous surfaces. Refer to the [CDC or OIE] website at [pathogen-specific website address] for additional information.

[Name of illness/outbreak] is caused by [name of emerging enveloped/small and large non-enveloped virus]. [Product name] kills similar viruses and therefore can be used against [name of emerging enveloped/small and large non-enveloped virus] when used in accordance with the directions for use against (cited above) on hard non-porous surfaces. Refer to the [CDC or OIE] website at [website address] for additional information.

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.

[Note to Reviewer: The following directions apply to both Residential and Non-Residential Use Sites mentioned on the Production Labels. Non-Residential Use Sites Disposal Directions may be used when Product is offered for both Residential/ Household use sites and Non-Residential (Professional/ Institutional and/or Industrial) use sites.]

[The following section applies to NON-REFILLABLE containers only]

PESTICIDE STORAGE: [Store at room temperature] -or- [Store product in a cool [dry] area, [away from direct sunlight and heat] to avoid deterioration.]

[For Spills and Leaks: Stop leak. Move containers from spill area. Absorb with a liquid absorbent material and place in an appropriate waste disposal container. Dispose of via a licensed waste disposal contractor. Refer to product SDS for additional information.]

Pesticide Disposal: Wastes resulting from the use of this product must be disposed of on site or at an approved waste disposal facility.

Container Disposal: Nonrefillable container. Do not reuse or refill this container. [Offer for recycling, if available.] [Triple rinse. Then offer for recycling or reconditioning, or puncture and dispose of in a sanitary landfill, or incineration or, if allowed by state and local authorities, by burning. If burned, keep out of smoke.]

[The following section applies to REFILLABLE containers only]

PESTICIDE STORAGE: [Store at room temperature] -or- [Store product in a cool [dry] area, [away from direct sunlight and heat] to avoid deterioration.]

[For Spills and Leaks: Stop leak. Move containers from spill area. Absorb with a liquid absorbent material and place in an appropriate waste disposal container. Dispose of via a licensed waste disposal contractor. Refer to product SDS for additional information.]

Pesticide Disposal: Wastes resulting from the use of this product may be disposed of according to applicable Federal, State or local procedures.

Container Handling: Refillable container [to be refilled by users]. Refill this container with CaviCide 1 [this product solution] only. Do not reuse this container for any other purpose. Cleaning the container before final disposal is the responsibility of the person disposing of the container. [A 50 gallon rinse with DI water is required to clean a 350 gallon tote.] Cleaning before refilling is the responsibility of the refiller.

[The following recycling statements apply to both REFILLABLE and NON-REFILLABLE containers. Choose one of the following recycling statements.]

[Offer for recycling, if available.]

[Triple rinse. Then offer for recycling or reconditioning, or puncture and dispose of in a sanitary landfill, or incineration or, if allowed by state and local authorities, by burning. If burned, keep out of smoke.]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

[Packaging Statements and Operating Instructions]

[Refill only with this product] [to be refilled by users]

[Use A – C as applicable depending on container type:]

A.

Operating Trigger

To Operate:

1. Turn -or- Rotate nozzle [counterclockwise] to “spray” -or- “stream” -or- ON [setting]. -or- Adjust [nozzle] to desired pattern.
2. When almost empty tilt bottle downward -or- at a downward angle [while spraying].

B.

To Refill [Spray Bottles] [to be refilled by users]:

1. Remove trigger sprayer -or- cap [from empty bottle]
2. Unscrew cap on refill and pour contents directly into empty bottle.
3. Replace trigger sprayer -or- cap [and use as you normally would].

C.

To Refill [to be refilled by users]: [Remove:] Remove cap from empty bottle. [Pour:] Unscrew cap on refill and pour contents directly into empty bottle. Use a funnel to minimize product spills -or- waste. [Use:] Replace cap and use as you normally would.

[CaviCide 1 is available in following packaging [configurations]:]

- [48 – 2oz bottles/CS]
- [12 – 24oz bottles/CS]
- [4 – 1 gallons/CS]
- [2 – 2.5 gallons/CS]
- [1 – 55 gallon drum/CS]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

CLEANING/DECONTAMINANT/DEODORIZING directions for use:

FOR CLEANING NON-INSTRUMENT SURFACES:

Remove visible soil loads. Spray CaviCide 1 directly onto surface. Allow to remain visibly wet for 30 seconds. Wipe surface using a towel.

FOR CLEANING INSTRUMENTS PRIOR TO DISINFECTION:

Remove visible soil loads.

[Listing of at least one of the following cleaning instructions, A – F, is required.]

[A]

[For use as a cleaning and pre-decontaminant immersion solution: Place instruments in a container with a sufficient amount of CaviCide 1 so as to allow for complete submersion, cover and immerse for 5 minutes at room temperature (68 °F/20 °C). Remove and rinse instruments. Follow with appropriate disinfection process. Discard solution after each use.]

[B]

[For use as a cleaning and pre-decontaminant spray: Place instruments in a container. Spray CaviCide 1 so as to visibly drench all instrument surfaces and allow to remain visibly wet for 30 seconds. Remove and rinse instruments. Follow with appropriate disinfection process.]

[C]

[For use as an ultrasonic cleaning solution: Visibly rinse instruments under running water to remove visible solids. Add 1 ounce CaviCide 1 per liter of water to ultrasonic unit. Immerse instruments in solution and activate ultrasonic unit for 5 minutes or longer, if necessary. Remove and rinse instruments. Follow with appropriate disinfection process. Discard solution after each use.]

[D]

[For use as a manual cleaner: Visibly rinse instruments under running water to remove visible solids. Place instruments in a container with 1 ounce CaviCide 1 per liter of water. Scrub instruments while submerged for 30 seconds using a stiff bristle brush until visibly clean. Remove and rinse instruments. Follow with appropriate disinfection process. Discard solution after each use.]

[E]

[For cleaning hair clippers and electric shears: While clipper/electric shear is running, hold it in the downward position and spray CaviCide 1 to thoroughly wet the blades. Allow to remain visibly wet for 30 seconds. Avoid spraying on the clipper case or dripping into the clipper housing. Rinse clipper/electric shears and follow with appropriate lubricating and disinfection process.]

[F]

[For cleaning salon implements and shears: Spray with CaviCide 1 to visibly wet the implement/shear. Allow to remain visibly wet for 30 seconds. Wipe away visible debris using a soft bristle brush or cloth. Rinse implement/shear and follow with appropriate lubricating and disinfection process.]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

[DISINFECTION] DIRECTIONS FOR USE

[**Note to Reviewer:** Contact time can be presented in 1 [one] min[ute] or 60 sec[onds] on Production Label.]

FOR USE AS A DISINFECTANT ON NON-INSTRUMENT SURFACES:

Spray CaviCide 1 6-8 inches from the surface, visibly wetting area to be disinfected. Allow surface to remain visibly wet for 1 minute [^3 minutes for Adenovirus II] at room temperature (68 °F/20 °C). Wipe surface using a clean towel or allow to air dry. Visible soil loads must be removed prior to disinfecting. [Use 1-minute contact time to disinfect against SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces.] [To kill SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces, use 1-minute contact time.]

FOR USE AS A DISINFECTANT ON NON-CRITICAL MEDICAL DEVICES¹, INSTRUMENTS AND IMPLEMENTS:

Instruments must be thoroughly cleaned to remove excess organic debris, rinsed and dried. Thoroughly clean and rinse lumens of hollow instruments. Spray CaviCide 1 6-8 inches from the surface, visibly wetting area to be disinfected. Allow surface to remain visibly wet for 1 minute [^3 minutes for Adenovirus II] at room temperature (68 °F/20 °C). Remove and rinse instruments. Wipe dry prior to use. [Use 1-minute contact time to disinfect against SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces.] [To kill SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces, use 1-minute contact time.]

¹ Non-critical medical devices are items that come in contact only with intact skin.

[Note: Critical and semi-critical devices must be followed by appropriate terminal sterilization/high-level disinfection process.]

[The following instruction is to be included when hard, non-porous exterior surfaces of ultrasound transducers or probes are listed as a use site.]

[To disinfect external surfaces of ultrasound transducers: Pre-clean ultrasound gel from external transducer to be disinfected [by using a clean wash cloth]. Spray CaviCide 1 directly onto surface, visibly wetting area to be disinfected. Allow surface to remain visibly wet for 1 minute [^3 minutes for Adenovirus II] at 68 °F/20 °C.] [Use 1-minute contact time to disinfect against SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces.] [To kill SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces, use 1-minute contact time.]

[The following paragraph will be included in its entirety when the production label references health care use.]

[This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization/high-level disinfection.]

[Always refer to manufacturer's care instructions before using on equipment -or- devices.]

[The following HIV-1 / HBV / HCV instructions are to be included on production labels when claims against HIV-1, HBV or HCV are listed.]

CaviCide 1 effectively kills [Human Immunodeficiency Virus Type 1 (HIV-1),] [Hepatitis B Virus (HBV),] [[and] Hepatitis C Virus (HCV)] on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects can be associated with the potential for transmission of [HIV-1,] [HBV,] [[and] [HCV].

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST [HIV-1,] [HBV,] [[AND] HCV] OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUID:

Personal Protection: Wear barrier protection [such as]: protective gloves, gowns, masks and eye coverings. **Cleaning Procedure:** Blood and other bodily fluids must be thoroughly cleaned from surfaces and objects before disinfection with CaviCide 1.

Contact Time: A contact time of 1 minute is required for HIV-1, HBV and HCV. [This contact time will not be sufficient for some of the organisms listed on the label. Refer to label for contact times against other organisms.]

Disposal Infectious Materials: Dispose of material in accordance with Federal, State or local regulations for infectious materials disposal.

‡SPECIAL LABEL INSTRUCTIONS FOR CLEANING PRIOR TO DISINFECTION AGAINST *CANDIDA AURIS*:

Personal Protection: Wear appropriate barrier protection [such as]: gloves, gowns, masks, or eye covering.

Cleaning Procedure: Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with the product. Pre-cleaning is to include vigorous wiping and/or scrubbing and all visible soil is removed. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading the organism. Restrooms are to be cleaned last. Do not reuse soiled cloths.

Infectious Waste Disposal: Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

The following statement is to be included on production labels when food contact areas or surfaces are listed.

+ A thorough rinse with potable water is required for surfaces in direct contact with food.

CLEANING AND STAIN REMOVAL DIRECTIONS FOR USE:

To clean -and/or- remove stains: Spray surface and; [rinse and] let air dry. [If streaking is observed, wipe with a clean, damp [cloth] [or] [paper towel].] [[Works on] [effective against] [removes the following stains:]] [CLEANING] [and] [DEODORIZING]: [Spray surface and wipe with a paper towel until visibly clean. Let air dry.]

| | | |
|---|--|--|
| <p>Claims:</p> <ul style="list-style-type: none"> • Cleans • Removes stains • Removes (<i>insert stain(s) /soil(s) shown to the right:</i>) | <p>Stains/ Soils</p> <ul style="list-style-type: none"> • Blood • Body oils • Bodily fluids • Coffee [stains] • Common soils • Dead skin • Dirt • Fingerprints • Food+ residue[s] -or soil[s] • Food stains • Fruit [stains] • Grease • Grime • Laboratory stains | <p>Stains/ Soils continued</p> <ul style="list-style-type: none"> • Mucus • [Other] common soils -and/or- stains • [Other] organic matter • Phlegm • Protein-based matter • Saliva • Spit • Spittle • Sputum • Throw-up • Tea [stains] • Urine • Vomit |
|---|--|--|

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

ALLERGEN DESTRUCTION

To **[clean and] destroy specified allergens:** Spray CaviCide 1 6-8 inches from the surface. [Visibly wetting area to be disinfected.] Allow surface to remain visibly wet for 1 minute [Rinse.] Let air dry. [If streaking is observed, wipe with a clean, damp [cloth or] paper towel after contact time has expired.]

Claims:

- Breaks Down Many Allergens[!]
- Destroys -or- deactivates dust mite matter
- [Penetrates and] breaks down dust mite matter -or- particles, cockroach matter -or- particles, and pet dander allergens [found in dust]

The following statement is to be included on production labels when animal housing facilities are listed (Including Veterinary Clinics, Farms and Poultry Houses)

For cleaning and disinfecting the following hard nonporous surfaces: feeding and watering equipment, utensils, instruments, cages, kennels, stables, catteries, etc.

1. Remove all animals and feeds from premises, animal transportation vehicles, crates, cages etc.
2. Remove all litter, droppings and manure from floors, walls and surfaces of facilities (including pens, chutes, and barns) occupied or traversed by animals.
3. Empty all troughs, racks, and other feeding and watering appliances.
4. Thoroughly clean all surfaces with soap or detergent and rinse with potable water. Saturate surfaces by spraying with product for a period of 1 minute. [3 minutes for Adenovirus Type II] [1 minute for SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces]
5. Spray all handling and restraining equipment such as leashes, muzzles, halters or ropes, as well as forks, shovels, and scrapers used for removing litter and manure.
6. Ventilate buildings and other closed spaces.
7. Do not house animals or employ equipment until treatment has been absorbed, set or dried.
8. Thoroughly scrub treated feed racks, troughs, mangers, automatic feeders, fountains, and waterers with soap or detergent, and rinse with potable water before reuse.

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

DISINFECTION [Organism List]

[DESCRIPTION] [CaviCide 1 is a multi-purpose [disinfectant]/ [cleaning product] [cleaner] for use on hard, non-porous inanimate surfaces.]

[CaviCide 1 is effective against the following microorganisms on hard, non-porous surfaces when used as directed:]

[This product has demonstrated effectiveness against the following organisms with a 1-minute contact time: **(insert organisms)**] [3 minutes for Adenovirus II]

[Disinfects **(insert organisms)** on hard non-porous environmental surfaces listed on this label] [When used as directed, CaviCide 1 provides effective Tuberculocidal / Virucidal*^[^] / Bactericidal / Fungicidal activity against the following organisms:]

[Note to Reviewer: the following table is optional on Production Labels.]

| Pathogen: | Bacteria, Enveloped Viruses*, Bloodborne Pathogens**, TB, Fungicidal, SARS-CoV-2*, Measles, Candida auris‡ | Non-enveloped Virus* ^[^] |
|---------------|--|--|
| Contact Time: | 1 min[ute] | [1 min[ute]] ^[^] 3 min[utes] |

| ORGANISMS: | [Strain and/or ATCC number] |
|--|------------------------------------|
| <i>Mycobacterium Bacteria (1 Minute Contact)</i> | |
| [<i>Mycobacterium tuberculosis</i> var: bovis (BCG) (TB)] | [Strain: Organon Teknika, Corp.] |
| <i>Bacteria (1 Minute Contact)</i> | |
| [<i>Acinetobacter baumannii</i>] | [ATCC 19606] |
| [<i>Burkholderia cepacia</i>] | [ATCC 25416] |
| [<i>Enterobacter cloacae</i>] | [ATCC 13047] |
| [<i>Klebsiella pneumoniae</i>] | [ATCC 4352] |
| [<i>Pseudomonas aeruginosa</i>] | [ATCC 15442] |
| [<i>Salmonella enterica</i>] | [ATCC 10708] |
| [<i>Staphylococcus aureus</i>] | [ATCC 6538] |
| <i>[Pathogenic] Fungi (1 Minute Contact)</i> | |
| [<i>Candida albicans</i>] | [ATCC 10231] |
| [<i>Candida auris</i> ‡] | [AR-BANK#0381] [from CDC] |
| [<i>Trichophyton interdigitale</i>] [(previously known as <i>Trichophyton mentagrophytes</i>)] | [ATCC 9533] |
| <i>Drug-Resistant Bacteria (1 Minute Contact)</i> | |
| [Carbapenem-Resistant <i>Klebsiella pneumoniae</i> [(CRKP)]] [CRE] | [ATCC BAA-2470] |
| [ESBL <i>Escherichia coli</i>] [ESBL <i>E. coli</i>] | [ATCC BAA-196] |
| [Methicillin Resistant <i>Staphylococcus aureus</i> (MRSA)] | [ATCC 33592] |
| [Methicillin Resistant <i>Staphylococcus epidermidis</i> (MRSE)] | [ATCC 51625] |
| [Multi-drug resistant [MDR] <i>Acinetobacter baumannii</i>] | [ATCC BAA-1605] |
| [Vancomycin Resistant <i>Enterococcus faecalis</i> (VRE)] | [ATCC 51299] |
| [Vancomycin Intermediate <i>Staphylococcus aureus</i>] [VISA] | [ATCC 700699] |

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

| <i>ORGANISMS continued</i> | [Strain and/or ATCC number] |
|---|---|
| <i>Viruses* [Non-enveloped], [Virucidal*^[^]] [†:] (1 Minute Contact) (^3 minute Contact)</i> | |
| [Adenovirus II*, [^]] (3 minutes) | [ATCC VR-846] |
| [Feline Calicivirus (surrogate for Norwalk and Norwalk-like viruses)*] [Norovirus*] [Norwalk Virus*] [Norwalk-like virus*] [Feline Calicivirus*] (1 minute) | [ATCC VR-782] |
| [Rotavirus*] (1 minute) | [ATCC VR-2018] |
| <i>Viruses* [Enveloped], [Virucidal*^[^]] [†:] (1 Minute Contact)</i> | |
| [Duck Hepatitis Virus Surrogate for [Hepatitis B Virus (HBV)]*] [**] | [Strain Grimaud] |
| [Bovine Viral Diarrhea Virus Surrogate for [Hepatitis C Virus (HCV)]*] [**] | [NADL] |
| [Herpes Simplex Virus Type 1*] | [ATCC VR-260] |
| [Herpes Simplex Virus Type 2*] | [ATCC VR-734] |
| [Human Coronavirus (not associated with Severe Acute Respiratory Syndrome or SARS)*] | [ATCC VR-740] |
| [Human Immunodeficiency Virus (HIV-1)*] [**] | [ZeptoMetrix] |
| [Influenza A Virus, Strain A/ [Hong Kong 8/68] [H3N2 Virus]]* | [Charles River Labs] |
| Severe Acute Respiratory Syndrome-related Coronavirus 2* -or- SARS-CoV-2* -or- COVID-19 Virus* | [Strain: USA-WA1/2020] [BEI Resources, NR-52281] |
| Measles* | [Strain: Edmonston] [ATCC VR-24] |

* Virus

[** Bloodborne pathogens]

[^ Virucidal – One minute contact time will not be sufficient for Adenovirus II.] [**Note to Reviewer:** this definition will only be included when Adenovirus Type II is listed Production Labels.]

† [One-Step Cleaning/ Disinfection] when Disinfection Directions are followed.

‡ [Candida auris –] Follow [the] Special Instructions for Cleaning required prior to Disinfection [listed on label].

[**Note to Reviewer:** Contact time can be presented in 1 [one] min[ute] or 60 sec[onds] on Production Label.]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS

▲ to represent any use site(s) listed in below as applicable.

◆ to represent any surface(s) listed in below as applicable.

[Note to Reviewer: Contact time can be presented in 1 [one] min[ute] or 60 sec[onds] on Production Label.]

- [99.9% -or 3 [three] log reduction against/for SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces [in [just] 1 minute] [in the presence of [5%] visible soil] [in one step]]
- [Kills 99.9% -or- 3 [three] log of SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces [in [just] 1 minute] [in one step] [in the presence of [5%] visible soil]]
- [99.999% -or 5 [five] log reduction for [*Candida auris*†] -or- [*C. auris*†] [in [just] 1 minute -or- 60 seconds]
- [Kills a minimum of 99.999% -or- 5 [five] log of *Candida auris*† -or- *C. auris*† on hard, non-porous surfaces] [in [just] 1 minute -or- 60 seconds]
- [1 Step Germicidal Spray]
- [1 Step Germicidal Spray against/ for SARS-CoV-2* or- COVID-19 Virus* on hard, nonporous surfaces]
- [1 minute] [3 minutes for Adenovirus II]
- [1 minute kill time] [3 minutes for Adenovirus II]
- [1 minute contact time] [3 minutes for Adenovirus II]
- [1 minute contact/kill time against/for SARS-CoV-2* or- COVID-19 Virus* on hard, nonporous surfaces [in the presence of [5%] soil load] [in one-step]]
- [A Metrex Research® brand]
- [A Metrex® brand]
- [A CaviCide Brand]
- [Advanced Cleaning Formula]
- [[Aids] [Helps] in reducing -or- in the reduction of cross-contamination between treated hard, non- porous surfaces]
- [[Aids] [Helps] in reducing -or- in the reduction of cross-contamination caused by SARS-CoV-2* or- COVID-19 Virus* between treated hard, non- porous surfaces]
- [An effective method to clean against odor-causing organisms on hard, non-porous surfaces listed on this label].
- [[A] one-step† cleaner/disinfectant]
- [[A] one-step† cleaner/disinfectant against/for SARS-CoV-2* or- COVID-19 Virus on hard, nonporous surfaces]
- [[A] One Step† disinfectant against [*Candida albicans*] (and/or) [*Candida auris*†]]
- [A disinfectant product for use [on ◆] [and other ◆] [in ▲]]
- [Always ready-to-use]
- [An effective method to deodorize and disinfect against odor-causing organisms on hard non- porous surfaces]
- [An effective way to disinfect hard, non-porous environmental surfaces]
- [An effective way to disinfect hard, non-porous environmental surfaces listed on this label.]
- [A disinfectant for use [in ▲] [on ◆ [surfaces]] and other critical care areas where control between hard non-porous treated surfaces is required]
- [Antibacterial]
- [Antimicrobial]
- [Fungicidal]
- [Appropriate] [Perfect] [Ideal] for [Occupied Room] [and] [Terminal Room] Cleaning]
- [Assures product effectiveness and standardizes technique]
- [Bacteria -and/or- virus* -and/or- germ killing [formula]]
- [Bactericidal] [Virucidal*[^3 minutes for Adenovirus II]] [Fungicidal] [Tuberculocidal] [in 1 minute]
- [Bactericide -or- Bactericidal] [(in 1 minute)]
- [Bathroom disinfectant]
- [Barber Disinfectant/Cleaner]
- [Broad spectrum disinfectant]
- [CaviCide 1 kills germs that cause odors]
- [CaviCide 1 is effective against odor causing germs]
- [Can help reduce cross contamination between treated, hard non-porous surfaces]
- [Cleaner and disinfectant] [in One]
- [[Cleans and] disinfects [◆] [in ▲]]
- [Clean and disinfect hard, non-porous surfaces easily]
- [[Cleans] [and] [disinfects] [on the go]]
- [Cleans and/or deodorizes and/or disinfects [hard, nonporous surfaces]]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Cleans, deodorizes and kills common household germs]
- [Cleans] [disinfects] [Cleaner and disinfectant] [in one] [all in one [easy] step†]
- [Clean Home]
- [Cleans the toughest soils]
- [Cleaning Made Easy]
- [[Cleans [and /or Disinfects] [grocery cart handles or shopping cart handles or grocery cart child seats or shopping] cart child seats]
- [Cleans grocery carts and shopping carts]
- [Clinical Clean]
- [Clipper Blade Disinfectant/Cleaner]
- [Convenient]
- [Convenient and easy cleaning all around [the office] [business] [workplace]]
- [[Convenient] [Portable] way to [disinfect] [clean] hard non-porous surfaces [of ♦] [in ▲]]
- [Cuts Cleaning time]
- [Deep cleaning] [formula]] [Delivers a clean you can see and smell]
- [Dental Office Disinfectant/Cleaner]
- [Demonstrates 99.999% -or- 5 [five] log reduction against for [*Candida auris*‡] -or- [*C. auris*‡] [in [just] 1 minute]
- [Demonstrates 99.9% -or- 3 [three] log reduction against SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces [in [just] [1] [one] minute] [in one step] [in the presence of [5%] visible soil]]
- [Designed for killing SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces [of ♦] [in ▲] [in [just] [1] [one] minute] [in one step] [in the presence of [5%] visible soil]]
- [Destroys] [Deactivates] [Disinfects] [[Helps] Eliminate] [99.9%] [tough-to-kill] SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces [of ♦] [in ▲] [in [just] [1] [one] minute] [in one step†] [in the presence of [5%] visible soil]
- [Designed for Home Healthcare -and/or- patients -and/or caregivers]
- [Designed -or- developed for Home Caregivers]
- [Designed -or- developed for use [on ♦ [surfaces]] in [▲]] [**Note to Reviewer:** this claim will not be listed on End-Use Product Labels.]
- [Delivers consistent [amount of solution] every time [!]]
- [Delivers consistent concentration every time [!]]
- [Deodorizes, and kills common household germs]
- [Deodorizes by killing the germs that cause odors]
- [Deodorize hard, non-porous surfaces] [♦] [in ▲]]
- [Designed for the [▲] [industry]]
- [Designed] [for] [with] [the] Environmental Services Professional[s] [in] [mind][!]
- [Disinfectant]
- [Disinfect]
- [Disinfect [external] hard, non-porous surfaces of ♦ [in ▲]]
- [Disinfectant cleaner [that cleans, disinfects and deodorizes] in one step†]
- [Disinfecting formula]
- [[Disinfects] [Disinfectant] [Cleaner] on the go]
- [Disinfects [and deodorizes] [♦] [in one step†] [in ▲]]
- [Disinfects [bacteria] [viruses*] [[multi] drug-resistant [bacteria] [pathogens]] [bloodborne pathogens**] [[pathogenic] fungi] [*Candida auris*‡ -or- *C. auris*‡] [in [just] 1 minute]
- [Disinfects and deodorizes in 1 minute] [3 minutes for Adenovirus II]
- [Disinfects hard, non-porous surfaces in 1 minute] [3 minutes for Adenovirus II]
- [Disinfects hard non-porous environmental surfaces listed on this label [Use to]]
- [[Disinfects] hard surfaces [that people touch]]
- [Disinfects [in [just] one [1] minute] [3 minutes for Adenovirus II]]
- [[Disinfects and deodorizes by killing germs and their odors]
- [[Disinfects hard, nonporous surfaces [♦] [throughout the ▲]
- [Disinfects [as it cleans] [and scrubs in one easy step†]]
- [Does not contain dyes]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Easy† [and convenient] to use]
 - [[Easy†] [and convenient] [way to clean] [and] [disinfect]]
 - [Easy One Step† Cleaning [and disinfecting] (except when precleaning is required for visible soil and blood and body fluids)]
 - [Easy† [and convenient] way to disinfect [bathrooms] [office] [nursing stations] [high touch] hard, non-porous surfaces everyday]
 - [Easy to handle]
 - [Easy† way to [clean] [and] [disinfect] [bathroom] [or exterior toilet surfaces] hard, non-porous surfaces]
 - [Easy† and convenient disinfecting [♦ [surfaces]] [throughout the ▲]]
 - [Easy, one-step† cleaning and disinfecting]
 - [Easy, one-step† [cleaning and] disinfecting against SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces]
 - [Easy† [and convenient] way to [touch-up] [bathroom] [office] everyday]
 - [Economical]
 - [Economy]
 - [Effective -or- Activity -or- For daily use against (*Insert any organism(s) from Organism List*) [and] *Candida auris*† - or- *C. auris*†] [in [just] [1] [one] minute]
 - [Effective -or- Activity -or- For daily use -or- Disinfects against SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces [in [just] [1] [one] minute] [in one step†] [in the presence of [5%] visible soil]]
 - [Effective against Influenza A and H1N1 [Influenza A Virus]
 - [Effective against (*Insert any organism from list of organisms*)
 - [Effective against bathroom bacteria and viruses* [use at least two bacteria and two viruses from this list:]: [ESBL *E. coli*] [*Salmonella enterica*] [*Staphylococcus aureus*] [Herpes Simplex Virus Type 2*] [Hepatitis B Virus (HBV)*] [Hepatitis C Virus (HCV)*] [Adenovirus II] [Influenza A, H3N2 Virus*] [*Enterobacter cloacae*] [*Burkholderia cepacia*] [*Carbapenem-Resistant Klebsiella pneumoniae* (CRKP)(CRE)]
 - [Effective against -or- Kills a range of bacteria including *Staphylococcus aureus*, *Salmonella enterica*, *Pseudomonas aeruginosa* [in the presence of organic soil [(5% blood serum)]]
 - [Effectively disinfects against MRSA]
 - [[Effective against -or- Kills (*insert any organism(s) from organisms list*) [in the presence of organic soil [(5% blood serum)]]
 - [Effective in the presence of 5% blood serum contamination]
 - [Effective against flu virus (Influenza A)]
 - [Effective against flu virus (Influenza A, H3N2 Virus)]
 - [Effective daily cleaner]
 - [Effectively disinfects against MRSA]
 - [Effectively disinfects against Norovirus]
 - [[Effectively] disinfects hard, nonporous [environmental] surfaces]
 - [Effectively removes everyday soils]
 - [Efficacy tests have demonstrated that [Product Name] [this product] [are] [is] an effective bactericide, virucide*, tuberculocide, and fungicide in the presence of organic soil (5% blood serum)]
 - [Efficacy tests have demonstrated that [Product Name] [this product] [is] an effective virucide* against/for SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces [in the presence of [5%] visible soil] [in [1] [one] minute]
 - [Eliminates odors]
 - [Eliminates odors from fecal matter and urine]
 - [Eliminates airborne particles associated with sprays]
 - [Eliminates [daily] mixing and measuring]
 - [Eliminates [lingering] [food+] odors [on surfaces]]
 - [[Eliminates odors at their source [caused by germs -or- bacteria-and/or- mildew]]
 - [[Eliminates -or- Reduces odors [caused by germs -or- bacteria -and/or- mildew]]
- [Note to Reviewer:** Claims above with “Eliminates” are NOT associated with pathogens but rather with odor or mixing / measuring; as such, these claims do NOT need to be qualified with “99.9%”.]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [EVS [Environmental Services] [use]]
- [Exhibits effective virucidal* activity]
- [Kill[s] -and/or- Effective against 99.9% [of]] [ESBL *Eschericia coli*] [ESBL *E. coli*]]
- [Kill[s] -and/or- Effective against 99.9% [of]] [*Salmonella enterica*] [(*Salmonella*)]
- [Kill[s] -and/or- Effective against flu (Influenza A) [virus*] [on environmental - and/or- hard, nonporous surfaces]]
- [Kill[s] -and/or- Effective against the flu (Influenza A) [virus*[es]]]
- [Kill[s] -and/or- Effective against (*insert any organism from organisms list*)]
- [Kill[s] -and/or- Help[s] reduce[s] the spread of -and/or- Effective against SARS-CoV-2* -or- COVID-19 Virus*] [in [just] [1] [one] minute] [in one step] [in the presence of [5%] visible soil] on hard, nonporous surfaces]
- [For daily use]
- [For] general use
- [For [daily] [general] [professional] use [especially] [during this [COVID-19] pandemic -or- outbreak]]
- [For use on hard, non-porous bathroom surfaces]
- [For use on hard, non-porous nursery surfaces]
- [For use on [♦ [surfaces]] [in ▲]]
- [For use in critical care areas between treated hard, non-porous surfaces [where control of cross contamination is of prime importance]
- [For] Surface cleaning
- [For Use by Home Caregivers]
- [For Use in Home Healthcare]
- [Formulated -or- Designed to] remove[s] biological odors including Urine -and/or- Feces -and/or- Vomit]
- [Fragrance Free]
- [Freshens]
- [Free of fragrance]
- [Fresh smell] [Helps reduce cross contamination between treated hard, non-porous surfaces]
- [Fungicide]
- [Fungicidal [disinfectant]] [against *T. interdigitale* -and/or- *C. auris*† -and-or *Candida albicans*]]
- [Fungicidal in 1 minute]
- [Germicide -or- Germicidal]
- [Gets rid of [or eliminates] odors] [**Note to Reviewer:** this “eliminates” claim is NOT associated with pathogens but “odors” only; as such, it does not need to be qualified with “99.9%”]
- [Great for everyday use]
- [[Helps] [kill][s] [eliminate][s] odor-causing bacteria]
- [Helps reduce cross contamination between treated hard, non-porous surfaces]
- [Healthcare Disinfectant/Cleaner]
- [Hospital [disinfectant] [disinfection] [use only when hospital use sites on label]
- [Institutional [size]]
- [Intermediate Level Disinfectant]
- [Just] [Simply] [spray and wipe] [exterior toilet surface] [bathroom germs] [away]
- [Kills Adenovirus II]
- [Kills [the] [Pandemic] SARS-CoV-2* -or- COVID-19 Virus* [in [just] [1] [one] minute] on hard, nonporous surfaces [in one step] [in the presence of [5%] visible soil]]
- [Kills [bacteria] [viruses*] [[pathogenic] fungi]
- [Kills bacteria in 1 minute]
- [Kills viruses* in 1 minute [(Adenovirus type 2 in 3 minutes)]
- [Kill [Effective against] [[pathogenic] fungi] [including] [*Candida albicans* -or- *C. albicans*] [*Candida auris*† -or- *C. auris*†] in [just] 1 minute]
- [Kills -or- Disinfects -or- Eliminates -or- [Is] Effective against Measles virus* on hard, nonporous surfaces [in [just] [1] [one] minute] [in one step] [in the presence of [5%] soil load]]
- [Kills -or- Disinfects -or- Eliminates -or- [Is] Effective against *Candida auris*† -or- *C. auris*† [in [just] 1 minute] on hard, non-porous surfaces]

[**Note to Reviewer:** “Eliminates” will be qualified with 99.999% if used.]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Kills [Norovirus] [Norwalk virus]]
- [Kills HIV-1]
- [Kills HCV]
- [Kills HBV]
- [Kills Rotavirus]
- [Kills Adenovirus II]
- [Kills bacteria, viruses*]
- [Kills [bacteria] [viruses*] [[pathogenic] fungi] on hard, nonporous surfaces]
- [Kills] [germs] [or bacteria] [while it cleans] [that cause odor]
- [Kills odor causing bacteria]
- [Kills [Disinfects] Norovirus [Feline Calicivirus] [in] 1 [one] minute]
- [Kills -or- Effective against bacteria]
- [Kills -or- Effective against germs]
- [Kills -or- Effective against viruses*]
- [Kills [or eliminates 99.9%] of bacteria]
- [Kills [or eliminates 99.9%] of viruses*]
- [Kills [or eliminates 99.9%] of fungus]
- [Kills odor-causing bacteria -or- germs]
- [[Kills multi drug resistant organisms [(list organisms)]]
- [Kills multi drug resistant bacteria [(list organisms)] [in [1] [one] minute]]
- [[Kills] [germs] [or bacteria] [while it cleans] [that cause odor]
- [Kills] [common] [household] [bacteria] [viruses*] [[pathogenic] fungi] [germs] [(insert any organism from organism list) [in 1 minute] [3 minutes for Adenovirus II]
- [Kills] [Eliminates 99.9%] [(list organisms) in [just] [1] minute] [(insert any organism from organism list) [in 1 minute] [3 minutes for Adenovirus II]
- [Kills -or- Destroys -or- Disinfects -or- Eliminates 99.9%-or- Effective against SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces [in [just] [1] [one] minute] [in one step†] [in the presence of [5%] soil load]]
- [Kills] [the] flu virus (Influenza A, H3N2 Virus)]
- [Kills] [viruses* and bacteria [including the flu virus (Influenza A, H3N2 Virus)]]
- [Kills TB in 1 Minute]
- [Kitchen+ and Bath]
- [Kitchen+ and Bath Solution for Hard, Non-porous Surfaces]
- [Large[r] Surface[s]]
- [Lead Free] [**Note to Reviewer:** this claim will not be listed on End-Use Product Labels.]
- [Leaves bathroom fresh and clean in [just] [1] minute]]
- [Leaves [bathroom] [office] clean [and disinfected] [and smelling fresh]
- [Leaves nothing behind]
- [Lightly scented -or- quick-dissipating scent]
- [Making home healthcare [more manageable] [easier] [simpler]]
- [Maximum cleaning, minimum effort]
- [Meets surface disinfection recommendations from OSHA[’s] Bloodborne Pathogen Guidelines -or- Standards]
- [Mercury free] [**Note to Reviewer:** this claim will not be listed on End-Use Product Labels.]
- [Multi-purpose disinfectant [for Health Care Facilities, Dental Facilities and Residential purposes]]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Multi-surface]
- [Neutralizes odors [doesn't [just] mask them]
- [No canister needed]
- [No fragrance added]
- [No measuring] [required]
- [No mixing] [required]
- [Non-abrasive formula]
- [[Occupied Room] [and] [Terminal Room] Cleaning Product]
- [Odor [elimination] [neutralization]
- [On the go disinfectant]
- [One-step† cleaner and disinfectant]
- [One step† [cleaning] [disinfecting] [product]] [liquid]]
- [One-step† cleaning and disinfecting [saves time]]
- [One step† cleaning [and/or disinfecting] [and/or deodorizing] disinfect and/or touch-up clean]
- [One step† cleaning [and/or disinfecting] [and/or deodorizing] disinfect and/or touch-up clean] [cleans and/or disinfects and/or deodorizes] [anytime] [daily] [every day] [in one easy step] on hard non-porous surfaces [of ♦] [in ▲] [bathroom and/or kitchen surfaces] [on household surfaces] [on ♦ [surfaces]]
- [One step† disinfectant cleaner that is effective against a broad spectrum of bacteria]
- [One step germicide]
- [One-step† disinfectant cleaner designed for general cleaning and disinfecting hard, nonporous environmental surfaces in ▲ [on ♦ [surfaces]]]
- [Perfect for cleaning [and disinfecting] [♦ [surfaces]] [in ▲]]
- [Perfect for disinfecting exterior toilet surfaces – or – around the exterior toilet surface]
- [Picks up [particles] [dust] [dirt] [hair] [crumbs]]
- [Portable]
- [Powers through tough messes]
- [Ready to Use [RTU] prevents formulation errors] [Premixed] [Premeasured] [Easy] [and convenient] way to disinfect]
- [[Product Name] are useful in a variety of [healthcare], [nursery/daycare] [and] [▲ [on ♦ [surfaces]]] applications listed on the [this] label.]
- [Easy† way to clean [and disinfect] [bathroom] [or exterior toilet surfaces] hard, non-porous surfaces]
- [Ready to use]
- [[Ready-to-use] disinfectant [for use in ▲ [on ♦ [surfaces]]]
- [Ready-to-use [one-step†] disinfectant [disinfection] for [killing] SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces]
- [Ready-to-use [one-step†] disinfectant [disinfection] for [killing] (***insert any organism from organism list***) [*Candida auris*† -or- *C. auris*†]]
- [Recommended for use in ▲ [on ♦ [surfaces]]]
- [Recyclable] [package]
- [Reduces cross contamination between treated hard, non-porous surfaces]
- [Reduces cross-contamination between treated hard, non-porous surfaces [from *Candida auris*† -or- *C. auris*†] [from (***insert any organism from organism list***) [in 1 minute)] [3 minutes for Adenovirus II].
- [Reduce [cross-contamination of] SARS-CoV-2* -or- COVID-19 Virus* between treated, hard, non-porous surfaces [of ♦] [in ▲]]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Reduces] [packaging] waste [by X%] [when compared to regular canister or bucket or pail]
- [Requires no mixing or measuring and is ready-to-use]
- [Salon Disinfectant/Cleaner]
- [Sample]
- [Saves Time]
- [School Disinfectant/Cleaner]
- [Shelf-stable formula]
- [So] [surfaces are more than [just] [clean] [they're] [disinfected]
- [Specifically -or- Specially designed for the Medical/Dental/ Veterinary/Health Care [▲] industry]
- [Starts cleaning on contact]
- [Suitable for use [on ♦ [surfaces]] [in ▲]]
- [Suitable for use on [titanium-coated] [medical grade] stainless steel]
- [Suitable for use on hard, nonporous inanimate environmental surfaces]
- [Surface Disinfectant/Cleaner]
- [Strong cleaning power]
- [The 1-minute [germicide] [disinfectant] [3 minutes for Adenovirus II]]
- [The 1-minute germicide/ disinfectant against *Candida auris*† -or- *C. auris*† (*insert any organism(s) from Organism List*)]
- [The answer to your [disinfecting] [cleaning] needs]
- [The] smell of clean]
- [[The] easy† -and/or- convenient way to disinfect] [with CaviCide™ 1/ this product]
- [This product was tested according to AOAC efficacy testing methods -or- standards for hospital disinfection]
- [This product has demonstrated effectiveness against Influenza A Virus]
- [This product has demonstrated effectiveness against SARS-CoV-2* -or- COVID-19 Virus* [and is expected to inactivate] SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces [in 1 minute] [in the presence of [5%] visible soil] [in one-step]]
- [This product has demonstrated effectiveness against Adenovirus II with a 3 minute contact time]
- [Time to reorder]
- [Trial]
- [Tuberculocidal (1 minute)]
- [Twin pack]
- [Two in one action [cleans and disinfects in one easy step†]]
- [Tuberculocide -or- Tuberculocidal]
- [To kill (*insert any organism(s) from Organism List*)], let stand [undisturbed] for 1 minute]
- [To kill *Candida auris*† -or- *C. auris*† [,] [on treated surfaces especially during outbreaks] [of *Candida auris*† -or- *C. auris*†] [let stand [undisturbed] for 1 minute]]
- [To kill SARS-CoV-2* -or- COVID-19 Virus* on treated hard, non-porous surfaces [especially during this/ the [COVID-19] outbreak] [,] [let stand [undisturbed] for 1 minute]]
- [To reduce *Candida auris*† -or- *C. auris*† on hard, non-porous treated surfaces [of ♦] [in ▲]]
- [Usable on variety of soils]
- [Use in [public] [common places] [in ▲] [on ♦ [surfaces]] where [bacteria -and/or- viruses*-and/or- germs may be of concern] [on hard, nonporous surfaces]
- [Use to] [Deodorize hard, non-porous surfaces]
- [Use to pre-clean critical or semi-critical medical devices prior to sterilization / high-level disinfection]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Use this product to] [Kill bacteria in 1 minute] [Kill viruses*^[^] in 1 minute] [3 minute contact time for Adenovirus II] [Kill TB in 1 minute] [Kill TB/HIV-1 in 1 minute] [Kill *Candida auris*† -or- *C. auris*† in 1 minute] [Kill SARS-CoV-2* -or- COVID-19 Virus* in 1 minute on hard, nonporous surfaces]
- [Use where cross-contamination between treated hard, non-porous surfaces is of prime importance]
- [Use on hard, non-porous surfaces only: [stainless steel] [Formica®] [Porcelain] [Plastic] [Glazed Ceramic] [glass] [glass tables] [glazed tile] [Test on small inconspicuous area first.]
- [Value]
- [Versatile]
- [Veterinary Disinfectant/ Cleaner]
- [[Virucide*^[^] -or- Virucidal*^[^]] [in 1 min] [3-minute contact time for Adenovirus II]
- [Virucide* -or- Virucidal* against SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces [in [just] 1 minute] [in the presence of [5%] visible soil] [in one-step†]]
- [Where cross contamination between treated surfaces is of prime importance]
- [Disinfects] [Kills] [common] [household] [bacteria] [viruses*] [germs]
- [Disinfects -or- Eliminates 99.9%] [Kills] [list organisms] in [just] [1 minute] [3 minutes for Adenovirus II]
- [Disinfects] [Kills] [the] flu virus* (Influenza A, H3N2 Virus)
- [Disinfects] [Kills] [viruses* and bacteria] [including the flu virus (Influenza A, H3N2 Virus)]
- [[Disinfects] [Kills] SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces [in one step] [in the presence of [5%] visible soil] [in [just] 1 minute]]
- [[X% reduction in packaging [material]]]

[Note: The following claims can be used only for the first 6 months of the product on shelf]

- [Better] [Improved] [New] [!]
- [Better] [Improved] [New] [New and improved] [formula] [packaging]
- [New] Brought to you by the makers of CaviCide [and CaviWipes]
- [New Claims]
- [New Claims and Contact Times]
- [New Claims – Kills TB in 1 Minute]
- [New and improved]
- [New] [RTU] [Ready To Use] [Solution] [Disinfectant]
- [New] [Surface Disinfectant]
- [New Kill Claims]
- [New – Saves Time]
- [New] – [Kills TB in 1 Minute]

[Note – the chart shown below shall be considered optional – e.g. inclusive in a bracket]

[When used as directed, CaviCide 1 provides effective Tuberculocidal** / Virucidal*^[^] / Fungicidal activity against the following organisms:]

| Tuberculocidal | Associated with: |
|---|-------------------------|
| <i>Mycobacterium tuberculosis</i> var: bovis (BCG) [(Organon Teknika, Corp.)] [(TB)] | Tuberculosis (BCG) |
| *Viruses | |
| Human Immunodeficiency Virus (HIV-1)* | AIDS |
| Influenza A, H3N2 Virus* | Common flu |
| Severe Acute Respiratory Syndrome-related Coronavirus 2* -or- SARS-CoV-2* -or- COVID-19 Virus* | COVID-19 [Virus] |
| Measles* virus | Measles |
| Fungus/Yeast | |
| <i>Trichophyton interdigitale</i> | Athletes Foot Fungus |
| <i>Candida auris</i> † | |

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**Table 1 Medical Use Sites:
CaviCide 1 is recommended for use on hard non-porous surfaces in [▲]
[in health care settings] [such as (insert use site)]:**

| | | |
|---|---|---|
| 3D Imaging Service Centers | Endoscopy [Centers] [Facilities] [Departments] | Ophthalmic [Facilities] [Offices] |
| Acute Care Facilities | Endoscopy Centers | Optometric [Facilities] [Offices] |
| Alternate Care Facilities | ENT Facility | Oral Health Care Facilities |
| Ambulances | Eye Care [Facilities] [Centers] | Orthodontic Offices |
| Ambulatory Surgery Centers (ASC) | Eye Surgical Centers | Orthopedic Offices |
| Animal Research Facilities | Forensic Labs | Outpatient Care Facilities |
| Assisted Living Facility | Full-Care Nursing Facility | Outpatient [Surgery] [Facilities] [Centers] |
| Autopsy Rooms | Gastrointestinal (GI) [Labs] [Laboratories] | Out-Patient Surgical Centers (OPSC) |
| Bathroom | Gastrointestinal (GI) Centers | Pathology Rooms |
| Bioskill Laboratories | Gastrointestinal (GI) Suites | Patient Care Areas |
| Blood Drawing Centers | Home [Health] Care [Services] | Patient Care Rooms |
| Blood Apheresis Facilities | Hospice | Patient Rooms |
| Blood Units | Hospital Kitchens | Pediatric Exam Rooms |
| Burn Center | Hydrotherapy [Facilities] [Rooms] | Pediatric Intensive Care Units (PICU) |
| Cardiac Intensive Care Unit (CICU) | In vitro Fertilizer Center (IVF) | Pediatric Offices |
| CAT Scan [Lab] [Laboratories] | Instrument Reprocessing [Facilities] [Departments] [Center] | Personal Care Service Facilities |
| Cath Labs | Intensive Care Units (ICU) | Pharmacy Clean Rooms |
| Central Service Areas | Isolation Areas | Phototherapy Rooms |
| Central Services Department | Isolation Rooms | Physical Therapy Facilities |
| Central Supply Areas | Labor Units | Physicians' Offices |
| Central Supply Rooms | [Lab] [Laboratories] | Plasma Apheresis [Centers] [Facilities] [Units] |
| Clean Rooms | Lasik Eye Centers | Post [Ops] [Operatory] Rooms |
| Clinical Areas | Lasik Units | Pre [Ops] [Operatory] Areas |
| Clinics | Laundry Rooms | Psychiatric Care Facility |
| Critical Care Units | Life Science Laboratories | Radiology Department |
| CS [Labs] [Laboratories] | Long-Term Acute Center (LTAC) | Recovery Rooms |
| Dialysis Centers | Long-Term Care [Facilities] | Rehabilitation Facilities |
| Dialysis Clinics | Medical Facilities [Clinics] | Respiratory Therapy Facilities |
| Dissection [Labs] [Laboratories] | Medical Related Facilities | Skilled Nursing [Facility] [Homes] |
| Doctor's Offices | Medivans | Sterile Processing Departments |
| Donation Centers [blood] [plasma] semen milk aphaeresis | Medical Intensive Care Unit (MICUs) | Sterile Storage Areas |
| Drug Preparation [Centers] [Facilities] [Rooms] | Magnetic Resonance Imaging (MRI) [facilities] [laboratories] [labs] | Surgery Centers |
| Doctor's Offices | Neonatal Intensive Care Units (NICU) | Surgery Rooms |
| Donation Centers [blood] [plasma] semen milk aphaeresis | Neonatal Nursery Areas | Surgical Centers |
| Drug Preparation [Centers] [Facilities] [Rooms] | Neonatal Units | Surgical Suites |
| Drug Preparation Rooms | Newborn Nursery | Teaching Hospitals |
| Electrolysis Facilities | Nuclear Medicine Facilities | Tissue Banks |
| Emergency Departments | Nurses' Stations | Transplant Center |
| Emergency Medical Services (EMS) Units | Nursing Homes | Trauma Center |

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

| | | |
|----------------------------|--|-------------------------|
| Emergency Medical Settings | Nursing Stations | University Research Lab |
| Emergency Rooms (ER) | Occupational Therapy [Centers] [Facilities] | Wound Center |
| Emergency Vehicles | Oncology Areas | X-Ray Rooms |
| Endo Suites | Operating Rooms | |

Table 2 Institutional Use Sites:

CaviCide 1 is recommended for use on hard non-porous surfaces in [▲]

| | | |
|------------------------------|-----------------------------|-------------------------------|
| Adult School | Early Learning Centers | Middle School |
| Air Transport Facilities | Elementary Schools | Military Installations |
| Airplane Interiors | [Fire] [Rescue Departments] | Nursery School |
| Airports | Fire Vehicles | Police Vehicles |
| Bathroom | Forestry Services | Pre-K School [Room] |
| Casino | Grammar School | Preparatory School |
| Child Care Facilities | High School | Preschool |
| Churches | Hospitality Facilities | Primary School |
| Colleges | Hotels | Prisoner Detention Facilities |
| Community Colleges | Industrial Facilities | Prisoner Lockups |
| Conference [Centers] [Rooms] | Infant Care Centers | Prisons |
| Conference Rooms Fixtures | Intermediate School | School Buses |
| Convention Facilities | Jails | School Kitchens+ |
| Correctional Facilities | Junior High School | Schools |
| Daycare Centers | Kindergarten | [University] [Universities] |
| [Dormitories] [Dorm] | Learning Institutions | |

Table 3 Dental Use Sites:

CaviCide 1 is recommended for use on hard non-porous surfaces in [▲]

| | | |
|-----------------------------|---------------|-------------------------------------|
| Dental [Laboratories] [Lab] | Dental Office | Dental [Operatory] [Operatories] |
|-----------------------------|---------------|-------------------------------------|

Table 4 Vet Use Sites:

CaviCide 1 is recommended for use on hard non-porous surfaces in [▲]

| | | |
|----------------------------|------------------------|-------------------------|
| Animal Care Facilities | Avian Units [Hatchery] | Pet Homes |
| Animal Control Vehicles | Bathroom | Pet Motels |
| Animal Dental Operatory | Doggy Day Care Centers | Pet Spas |
| Animal Detention Room | Feline Facilities | Pet Stores |
| Animal Research Facilities | Feline Units | Small Animal Facilities |
| Animal Shelters | Kennels | Veterinary Offices |
| Avian [Care] Facilities | Pet Areas | Veterinary Facilities |

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 5 Food Prep+ Use Sites:

CaviCide 1 is recommended for use on hard non-porous surfaces in [▲]

| | | |
|--------------------------|-----------------------------|------------------|
| Bakeries+ | Food Processing Facilities+ | Restaurants+ |
| Banquet Hall Facilities+ | Grocery Stores + | School Kitchens+ |
| Butcher Shops+ | Kitchen+ | Supermarkets+ |
| Food Preparation Areas+ | Produce Areas+ | |

[+ A thorough rinse with potable water is required for surfaces in direct contact with food.]

Table 6 Miscellaneous Use Sites:

CaviCide 1 is recommended for use on hard non-porous surfaces in [▲]

| | | |
|----------------------------------|-------------------------------------|------------------------------------|
| Airplanes | Factories | Massage Facility |
| Arboretums | Foots Spas | Money Handling Facilities |
| Arcades | Foreign Current Exchange Facilities | Morgues |
| [Athletic] [Training] Facilities | Funeral Homes | Mortuaries |
| Auto Facilities | Gaming Facilities | Motels |
| Banks | Garden Centers | Movie Theaters |
| Barber Shops | Golf Courses | Outdoor Playground |
| Bathrooms | Grounds Keeping Services | Personal Care Service Facilities |
| Beauty Salons | Gymnasium | Public Transportation Vehicles |
| Body Piercing Salons | Health Club Facilities | Recreational [center] [facilities] |
| Body Service Facilities | Helicopter Interiors | Recycling Centers |
| Body Storage Facilities | Household Bathrooms | Salon Settings |
| Botanical Gardens | Housekeeping Rooms | Weather Shelters |
| Bowling Alleys | Indoor Playground | Shower Rooms |
| Buses | Laundromats | Skincare Salons |
| Cabs | Laundry Rooms | Storage [Areas] [Rooms] |
| Cadaver Processing areas | Lawn Services | Subway |
| Cars | Limousine | Tanning Salons |
| Condominiums | Living Rooms | Tattoo Parlors |
| Cruise Ships | Locker Rooms | Taxi |
| Currency Exchange | Lockup Facilities | Theme Parks |
| Day Spa | Manicure Salons | Transport Vehicles |
| Delivery Units | Manufacturing [Facilities] [Plant] | |

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces:

CaviCide 1 is recommended for use on hard non-porous surfaces [of] [such as] [♦]:

| | | |
|--|-----------------------|---|
| Ambulance Equipment Surfaces | Biohazard Container | Combs |
| Animal Cages | Blood Pressure Cuff | Computer |
| Animal Equipment | Booster Chairs+ | Computer Keyboards |
| Appliances+ | Booster Seats+ | Computer Mouse |
| Arm Rests | Brushes | Copy Machines |
| ATM Surfaces | Cabinet Handles | Counter Tops+ |
| Autoclave Surfaces | Cabinets | Counters |
| Automatic Feeders | Cafeteria Trays+ | Cribs+ |
| Automatic Teller Machines (ATM) | Cages | Crutches |
| Baby Warmers | Cardiac Gym Equipment | Dental Chairs |
| Baby Cribs | Carts | Dental Countertops |
| Baked Enamel | Cash Register | Dental Curing Lights |
| Bars | Ceilings | For use as a cleaner and disinfectant on cured, set and hardened, hard non-porous [dental impressions] [dental impression materials] |
| Basins | Cell Phones | Dental Impression Trays |
| Baskets | Chairs | For use as a cleaner and disinfectant on hard nonporous [dental impressions] [dental impression materials] prior to sending to a dental lab |
| Bathing Units | Changing Tables | Dental Operatory Surfaces |
| Bathroom Fixtures | Chemotherapy Hoods | Dental Unit Instrument Trays |
| Bathtubs | Child Seats+ | Dentist Chair |
| Bed Pans | Chrome | Desk Tops |
| Bed Railings [Bedside Rails] | Clipper Blades | Desks |
| Bedside Commode | Closet Handles | Diaper Changing Stations |
| Behind and Under Sinks | Coated Mattresses | Diaper Pails |
| Bidet | Coated Pillows | Digital Phosphor Plates |
| <ul style="list-style-type: none"> • For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [♦ [surfaces] [in ▲]] [and other critical care areas where cross-contamination between treated hard, non-porous surfaces is important] [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.] • [Surfaces which are made of plastic (i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene), and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.] • + A thorough rinse with potable water is required for surfaces in direct contact with food. • + Any washable (food and non-food contact) surface where disinfection is required. | | |

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces continued:

| | | |
|---|--|---|
| Dining Tables+ | Exterior Surfaces of Toilet | Exterior Surfaces of Endodontic Equipment (Apex Locators, Pulp Testers, Motors) |
| Dish Racks+ | Exterior or Exposed Interior Surfaces of Medical Equipment | Exterior Surfaces of Radiology Equipment |
| Display Cabinets | Exterior Surfaces of Dental Equipment | Exterior Surfaces of Refrigerator Units+ |
| Diving Equipment | Exterior Surfaces of Amalgamators | Exterior Surfaces of Respiratory Therapy Equipment |
| Diving Masks+ | Exterior Surfaces of Anesthesia Equipment | Exterior Surfaces of Urinals |
| Diving Snorkels+ | Exterior Surfaces of Aphaeresis Machines | Exterior Surfaces of Vanity Tops |
| Door Knobs | Exterior Surfaces of Appliances such as Refrigerators+, Microwaves+, Food Processors+ and Stoves+π | External Surfaces Emergency Resuscitative Equipment |
| Dopplers | Exterior Surfaces of Center Venous Lines | External Surfaces of Bladder Scanners |
| Drain Boards+ | Exterior Surfaces of CICU Equipment | External Surfaces of Blood Glucose Meters |
| Drained Whirlpool Tank Surfaces | Exterior Surfaces of Defibrillators | External Surfaces of BP Monitors |
| Enamel | Exterior Surfaces of Dental Care Facilities Hand Pieces | External Surfaces of CAT Scan Machines |
| Endoscope Cabinet | Exterior Surfaces of Dental Equipment Lines | External Surfaces of Compression Devices |
| Equipment Surfaces | Exterior Surfaces of Dental Equipment Tubing | External Surfaces of Dialogistic Equipment |
| Exam Tables | Exterior Surfaces of Dental Operatories | External Surfaces of Dictating Equipment |
| Examination Table | Exterior Surfaces of Diagnostic Equipment | External Surfaces of Digital X-Ray Sensors |
| External Surfaces of Kendrick Extrication Device (KED) | Exterior Surfaces of Dialysis Machines | External Surfaces of EKG Machines |
| <ul style="list-style-type: none"> • For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [♦ [surfaces] [in ▲]] [and other critical care areas where cross-contamination between treated hard, non-porous surfaces is important] [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.] • [Surfaces which are made of plastic (i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene), and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.] • + A thorough rinse with potable water is required for surfaces in direct contact with food. • + Any washable (food and non-food contact) surface where disinfection is required. • π Allow surfaces to come to room temperature prior to disinfection. | | |

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces continued:

| | | |
|--|---|---|
| External Surfaces of Electronic Thermometers | External Surfaces of Oscopes | External Surfaces of Ultrasound Transducers |
| External Surfaces of Electrosurgical Generators | External Surfaces of Patient Care Monitor Unit | External Surfaces of Vaporizer |
| External Surfaces of Glucometers | External Surfaces of Patient Monitoring Equipment | External Surfaces of Veterinary Equipment |
| External Surfaces of Headsets | External Surfaces of Probes | External Surfaces of Vital Sign Monitors |
| External Surfaces of ICU Equipment | External Surfaces of Pulp Tester Motors | External Surfaces of X-Ray Machines |
| External Surfaces of Incubators | External Surfaces of Pulse Oximeter | Eye Glasses (do not use on contact lenses) |
| External Surfaces of Infusion Pumps | External Surfaces of Sequential Compression Devices | Faucets |
| External Surfaces of Intubation Equipment | External Surfaces of Respiratory Ventilators | Fax Machines |
| External Surfaces of Isolettes | External Surfaces of Resuscitators | Feed Racks |
| External Surfaces of IV Pumps | External Surfaces of Sterilizers | Filing Cabinets |
| External Surfaces of Life Support Equipment | External Surfaces of Suction Pumps | Filing Doors |
| External Surfaces of Loupes | External Surfaces of Surgery Microscopes | Filing Handles |
| External Surfaces of Medical Equipment | External Surfaces of Surgery Robots | Fire Trucks |
| External Surfaces of Medical Warmers | External Surfaces of Surgical Headlights | Floors |
| External Surfaces of Monitors Electrodes | External Surfaces of Surgical Light Sources | Food Cases+ |
| External Surfaces of MRI Machines | External Surfaces of Tympanic Thermometers | Food Serving Trays+ |
| External Surfaces of OLED Monitors | External Surfaces of Ultrasound Machines | Formica |
| External Surfaces of Ophthalmoscopes | External Surfaces of Ultrasound Monitors | Fountains+ |
| <ul style="list-style-type: none"> • For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [♦ [surfaces] [in ▲]] [and other critical care areas where cross-contamination between treated hard, non-porous surfaces is important] [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.] • [Surfaces which are made of plastic (i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene), and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.] • + A thorough rinse with potable water is required for surfaces in direct contact with food. • + Any washable (food and non-food contact) surface where disinfection is required. | | |

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces continued:

| | | |
|---|-------------------------------------|---|
| Exterior Surfaces of Freezers+ | IV Poles | Outdoor Furniture+ (excluding wood frames and upholstery) |
| Garbage Cans | IV Stands | Ovens+π |
| Glass | Jacuzzi | Overbed table |
| Glazed Ceramic Tile | Kennels | Oxygen Hoods |
| Glazed Porcelain | Knee Walkers | Patient Chairs |
| Glazed Tiled Walls | Knobs | Patient Delivery Equipment |
| Glazed Tiles | Laboratories | Patient Lift Equipment |
| Golf Cart | Laboratory Equipment | Patient Scales |
| Grab Bar | Laboratory Surfaces | Patient Support Equipment |
| Grocery Cart | Laminated Surfaces | Pens |
| Gurneys | Lamps | Phone Chargers |
| Gym Equipment | Light Lens Covers | Phone Cradle |
| Hair Clippers | Lights | Physical Therapy (PT) Equipment surfaces |
| Hair Cutting Implements | Linoleum | Pipes |
| Hair Dryers | Manicure Implements | Plastic Laminate |
| Hair Rollers | Manicure Tables | Plastic Mattress Covers |
| Hampers | Marlite | Playpens |
| Handles | Plastic Mattress Covers | Pneumatic Tourniquets |
| Handrails | Mayo tray stands | Police Cars |
| Hard Plastic or Metal Flooring | Mayo tray | Porcelain Enamel |
| Hard, non-porous environmental surfaces | Monitors | Plastic Posey Vest |
| High Chairs+ | Nail/Hair Care Implements | Public Restrooms |
| High Touch Objects | Non-food contact surfaces+ | Razors |
| Hoods+ | Nurse Call [Button] [Device] [Cord] | Reception Areas |
| Hospital Bed (Railings, Frames) | Operating Room Furniture | Reception Counters |
| Hospital Patient Bed | Operating Room Lights | Exterior Surfaces of Refrigerator Carts |
| Infant Incubators (interior and exterior surfaces) | Operating Room Tables | Exterior Surfaces of Refrigerators+ |
| Inside of Drawers | Operating Tables | Safety Straps |
| Instrument Stands | Operatory Light Switches | Salad Bar Sneeze Guards+ |
| <ul style="list-style-type: none"> • For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [♦ [surfaces] [in ▲]] [and other critical care areas where cross-contamination between treated hard, non-porous surfaces is important] [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.] • [Surfaces which are made of plastic (i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene), and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.] • + A thorough rinse with potable water is required for surfaces in direct contact with food. • + Any washable (food and non-food contact) surface where disinfection is required. • π Allow surfaces to come to room temperature prior to disinfection. | | |

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces continued:

| | | |
|---|-----------------------------------|--------------------------|
| Salon Surfaces | Stethoscopes | Transport Vehicles |
| School Buses | Stools | Trash Cans |
| Scissors | Storage | Tray Tables |
| Shampoo Bowls | Stoves+π | Troughs |
| Shears | Stovetops+π | Tubing |
| Shelves | Stretchers | Vanity Tops+ |
| Shopping Carts | Support Bars | Veterinary Care Surfaces |
| Shower Chair | Synthetic Marble | Vinyl |
| Shower Fixtures | Tables+ | Walkers |
| Shower Stalls | Tanning Beds | Wash Basin |
| Showers | Telemetry Equipment | Washable Nail Files |
| SilkPeel™ with DermalInfusion™ | Telephones | Watering Appliances |
| Sinks+ | Tiles+ | Weight-lifting Surfaces |
| Slit Lamps | Toilet above the water line | Wheelchairs |
| Smart Phone Chargers | Toilet (Rim, Seat, Tank, Handle,) | Whitening Lamps |
| Spine Back Boards | Touch Screen Monitors | Workstations |
| Stainless Steel | Towel Dispensers | Transfer Bench |
| Stalls | Plastic Toys | Stands |
| <ul style="list-style-type: none"> • For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [♦ [surfaces] [in ▲]] [and other critical care areas where cross-contamination between treated hard, non-porous surfaces is important] [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.] • [Surfaces which are made of plastic (i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene), and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.] • + A thorough rinse with potable water is required for surfaces in direct contact with food. • + Any washable (food and non-food contact) surface where disinfection is required. • π Allow surfaces to come to room temperature prior to disinfection. | | |