



## OFFICE OF CHEMICAL SAFETY AND POLLUTION PREVENTION

WASHINGTON, D.C. 20460

February 26, 2024

Rebecca Lei  
Principal Specialist  
Regulatory Affairs  
Metrex Research  
Electronic Transmittal: [Rebecca.Lei@envistaco.com](mailto:Rebecca.Lei@envistaco.com)

Subject: Notification per PRN 98-10 – Label Notification for the Addition of Optional Graphics for the subject product per PR Notice 98-10 Section II. H.  
Product Name: CaviCide 1  
EPA Registration Number: 46781-12  
Received Date: November 1, 2023  
Action Case Number: 00489755

Dear Rebecca Lei:

The Agency is in receipt of your Application for Pesticide Notification under Pesticide Registration Notice (PRN) 98-10 for the above referenced product. The Antimicrobials Division (AD) has conducted a review of this request for its applicability under PRN 98-10 and finds that the action requested falls within the scope of PRN 98-10.

The label submitted with the application has been stamped “Notification” and will be placed in our records.

Should you wish to add/retain a reference to the company’s website on your label, then please be aware that the website becomes labeling under the Federal Insecticide Fungicide and Rodenticide Act and is subject to review by the Agency. If the website is false or misleading, the product would be misbranded and unlawful to sell or distribute under FIFRA section 12(a)(1)(E). 40 CFR 156.10(a)(5) lists examples of statements EPA may consider false or misleading. In addition, regardless of whether a website is referenced on your product’s label, claims made on the website may not substantially differ from those claims approved through the registration process. Therefore, should the Agency find or if it is brought to our attention that a website contains false or misleading statements or claims substantially differing from the EPA approved registration, the website will be referred to the EPA’s Office of Enforcement and Assurance.

Page 2 of 2  
EPA Reg. No. 46781-12  
Action Case Number: 00489755

If you have any questions, please contact Emily Saunders at 202-566-0624 or by email at [saunders.emily@epa.gov](mailto:saunders.emily@epa.gov).

Sincerely,

A handwritten signature in black ink that reads "Steven Snyderman". The signature is written in a cursive style with a horizontal line underneath the name.

Steven Snyderman, Product Manager 33  
Regulatory Management Branch II  
Antimicrobials Division (7510M)  
Office of Pesticide Programs

Enclosure: Stamped Final Label

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

# CaviCide 1™

[TUBERCULOCIDAL] [BACTERICIDAL] [VIRUCIDAL\*<sup>[^]</sup>] [FUNGICIDAL]

**ACTIVE INGREDIENTS**

Didecyldimethylammonium chloride .....0.76%  
Ethanol ..... 7.50%  
Isopropanol.....15.00%  
**[Inert] Other Ingredients†.....76.74%**  
**Total.....100.00%**

KEEP OUT OF REACH OF CHILDREN

**WARNING**

[See other panel for additional precautionary statements.]

EPA Reg. No. 46781-12  
EPA Est. No. [insert EPA Est. No of Producer List on Box 2 of CSF]  
[56952-WI-001] [75272-WI-001]  
[71979-SC-001] [51435-WI-004]  
[56952-CHN-001] [73027-WI-001]  
[46781-MI-001] [93725-TN-001]  
[90532-WI-001] [90532-WI-002]  
[90532-WI-003] [0312-WI-3]

[Beginning of batch code indicates Est No.]  
[End of batch code indicates Est No.]  
Made in USA

[CaviCide 1 is available in the following packaging configurations]:

- 4 gallons/CS
- 48 – 2 oz bottles/CS
- 12 – 8 oz bottles/CS
- 12 – 24 oz bottles/CS
- 2 – 2.5 gallons/CS
- 1 – 55 gallon drum/CS
- 1 - 350 gallon transport tote]

[1 minute]  
[1 Step Germicidal Spray]  
[Fragrance Free]  
[Ready-To-Use Surface Disinfectant Cleaner]

Reorder No. XXXXXX  
Net Contents:  
X [[Gallon][gal]] [[Ounce][oz]] [[Fluid Ounce][fl. oz]]/ X [[Liter][L]] [[Milliliters][mL]]  
(Net Weight x lbs./ x kg.)  
(Net Weight x lbs/ x g)

[Packed For:]  
[Distributed By:]  
[Sold By:]  
[Manufactured By]

Metrex Research  
28210 Wick Road  
Romulus, MI 48174

[To reorder] [,] [For product or technical information] contact Metrex at 800-841-1428 or visit our web site at

[www.metrex.com](http://www.metrex.com)

[For SDS information please visit

[www.metrex.com](http://www.metrex.com)

**NOTIFICATION**

46781-12

The applicant has certified that no changes, other than those reported to the Agency have been made to the labeling. The Agency acknowledges this notification by letter dated:

02/26/2024

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

## PRECAUTIONARY STATEMENTS

### Hazards to Humans and Domestic Animals

**WARNING** Causes substantial but temporary eye injury. Do not get in eyes or on clothing. Wear [appropriate protective eyewear such as] goggles, face shield, or safety glasses. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse. [See other panel for additional precautionary statements.]

**Physical or Chemical Hazards: Combustible.** Do not use or store near heat or open flame.

#### FIRST AID

**If in eyes:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, and then continue rinsing. Call a poison control center or doctor for treatment advice.

**For Chemical Emergency ONLY** (spill, leak, fire, exposure or accident), call CHEMTREC at [insert phone number]. For [ALL] [other]] inquiries about this product, call Metrex Customer [Service] [Care] at [insert phone number], Monday through Friday, between 6 [or any time] a.m. and 4 [or any time] p.m. Pacific Time.

Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

***[The following Environmental Hazards statement is to be included only on labeling of containers that are 5 gallons and larger.]***

Environmental Hazards: This pesticide is toxic to fish. Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office of the EPA.

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

### EMERGING VIRAL PATHOGEN CLAIMS

This product qualifies for emerging viral pathogen claims per the EPA’s ‘Guidance to Registrants: Process for Making Claims Against Emerging Viral Pathogens not on EPA-Registered Disinfectant Labels’ when used in accordance with the appropriate use directions indicated below.

This product meets the criteria to make claims against certain emerging viral pathogens from the following viral category:

- Enveloped Viruses
- Large Non-enveloped
- Small Non-enveloped

<i>For an emerging viral pathogen that is a/an...</i>	<i>...follow the directions for use for the following organisms on the label:</i>
<b>Enveloped virus</b>	<b>Rotavirus</b>
<b>Large, non-enveloped virus</b>	<b>Feline Calicivirus</b>

*[Product name] has demonstrated effectiveness against viruses similar to [name of emerging enveloped/small and large non-enveloped virus] on hard, non-porous surfaces. Therefore, [product name] can be used against [name of emerging enveloped/small and large non-enveloped virus] when used in accordance with the directions for use against (cited above) on hard, non-porous surfaces. Refer to the [CDC or OIE] website at [pathogen-specific website address] for additional information.*

*[Name of illness/outbreak] is caused by [name of emerging enveloped/small and large non-enveloped virus]. [Product name] kills similar viruses and therefore can be used against [name of emerging enveloped/small and large non-enveloped virus] when used in accordance with the directions for use against (cited above) on hard non-porous surfaces. Refer to the [CDC or OIE] website at [website address] for additional information.*

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

## STORAGE AND DISPOSAL

**Do not contaminate water, food, or feed by storage and disposal.**

[Note to Reviewer: The following directions apply to both Residential and Non-Residential Use Sites mentioned on the Production Labels. Non-Residential Use Sites Disposal Directions may be used when Product is offered for both Residential/ Household use sites and Non-Residential (Professional/ Institutional and/or Industrial) use sites.]

### [The following section applies to NON-REFILLABLE containers only]

**PESTICIDE STORAGE:** [Store at room temperature] -or- [Store product in a cool [dry] area, [away from direct sunlight and heat] to avoid deterioration.]

[For Spills and Leaks: Stop leak. Move containers from spill area. Absorb with a liquid absorbent material and place in an appropriate waste disposal container. Dispose of via a licensed waste disposal contractor. Refer to product SDS for additional information.]

**Pesticide Disposal:** Wastes resulting from the use of this product must be disposed of on site or at an approved waste disposal facility.

**Container Disposal:** Nonrefillable container. Do not reuse or refill this container. [Offer for recycling, if available.] [Triple rinse. Then offer for recycling or reconditioning, or puncture and dispose of in a sanitary landfill, or incineration or, if allowed by state and local authorities, by burning. If burned, keep out of smoke.]

### [The following section applies to REFILLABLE containers only]

**PESTICIDE STORAGE:** [Store at room temperature] -or- [Store product in a cool [dry] area, [away from direct sunlight and heat] to avoid deterioration.]

[For Spills and Leaks: Stop leak. Move containers from spill area. Absorb with a liquid absorbent material and place in an appropriate waste disposal container. Dispose of via a licensed waste disposal contractor. Refer to product SDS for additional information.]

**Pesticide Disposal:** Wastes resulting from the use of this product may be disposed of according to applicable Federal, State or local procedures.

**Container Handling:** Refillable container [to be refilled by users]. Refill this container with CaviCide 1 [this product solution] only. Do not reuse this container for any other purpose. Cleaning the container before final disposal is the responsibility of the person disposing of the container. [A 50 gallon rinse with DI water is required to clean a 350 gallon tote.] Cleaning before refilling is the responsibility of the refiller.

**[The following recycling statements apply to both REFILLABLE and NON-REFILLABLE containers. Choose one of the following recycling statements.]**

[Offer for recycling, if available.]

[Triple rinse. Then offer for recycling or reconditioning, or puncture and dispose of in a sanitary landfill, or incineration or, if allowed by state and local authorities, by burning. If burned, keep out of smoke.]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

## [Packaging Statements and Operating Instructions]

[Refill only with this product] [to be refilled by users]

***[Use A – C as applicable depending on container type:]***

### **A.**

Operating Trigger

To Operate:

1. Turn -or- Rotate nozzle [counterclockwise] to “spray” -or- “stream” -or- ON [setting]. -or- Adjust [nozzle] to desired pattern.
2. When almost empty tilt bottle downward -or- at a downward angle [while spraying].

### **B.**

To Refill [Spray Bottles] [to be refilled by users]:

1. Remove trigger sprayer -or- cap [from empty bottle]
2. Unscrew cap on refill and pour contents directly into empty bottle.
3. Replace trigger sprayer -or- cap [and use as you normally would].

### **C.**

To Refill [to be refilled by users]: [Remove:] Remove cap from empty bottle. [Pour:] Unscrew cap on refill and pour contents directly into empty bottle. Use a funnel to minimize product spills -or- waste. [Use:] Replace cap and use as you normally would.

**[CaviCide 1 is available in following packaging [configurations]:]**

- [48 – 2oz bottles/CS]
- [12 – 24oz bottles/CS]
- [4 – 1 gallons/CS]
- [2 – 2.5 gallons/CS]
- [1 – 55 gallon drum/CS]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

## DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

***CLEANING/DECONTAMINANT/DEODORIZING directions for use:***

### **FOR CLEANING NON-INSTRUMENT SURFACES:**

Remove visible soil loads. Spray CaviCide 1 directly onto surface. Allow to remain visibly wet for 30 seconds. Wipe surface using a towel.

### **FOR CLEANING INSTRUMENTS PRIOR TO DISINFECTION:**

Remove visible soil loads.

***[Listing of at least one of the following cleaning instructions, A – F, is required.]***

#### ***[A]***

**[For use as a cleaning and pre-decontaminant immersion solution:** Place instruments in a container with a sufficient amount of CaviCide 1 so as to allow for complete submersion, cover and immerse for 5 minutes at room temperature (68 °F/20 °C). Remove and rinse instruments. Follow with appropriate disinfection process. Discard solution after each use.]

#### ***[B]***

**[For use as a cleaning and pre-decontaminant spray:** Place instruments in a container. Spray CaviCide 1 so as to visibly drench all instrument surfaces and allow to remain visibly wet for 30 seconds. Remove and rinse instruments. Follow with appropriate disinfection process.]

#### ***[C]***

**[For use as an ultrasonic cleaning solution:** Visibly rinse instruments under running water to remove visible solids. Add 1 ounce CaviCide 1 per liter of water to ultrasonic unit. Immerse instruments in solution and activate ultrasonic unit for 5 minutes or longer, if necessary. Remove and rinse instruments. Follow with appropriate disinfection process. Discard solution after each use.]

#### ***[D]***

**[For use as a manual cleaner:** Visibly rinse instruments under running water to remove visible solids. Place instruments in a container with 1 ounce CaviCide 1 per liter of water. Scrub instruments while submerged for 30 seconds using a stiff bristle brush until visibly clean. Remove and rinse instruments. Follow with appropriate disinfection process. Discard solution after each use.]

#### ***[E]***

**[For cleaning hair clippers and electric shears:** While clipper/electric shear is running, hold it in the downward position and spray CaviCide 1 to thoroughly wet the blades. Allow to remain visibly wet for 30 seconds. Avoid spraying on the clipper case or dripping into the clipper housing. Rinse clipper/electric shears and follow with appropriate lubricating and disinfection process.]

#### ***[F]***

**[For cleaning salon implements and shears:** Spray with CaviCide 1 to visibly wet the implement/shear. Allow to remain visibly wet for 30 seconds. Wipe away visible debris using a soft bristle brush or cloth. Rinse implement/shear and follow with appropriate lubricating and disinfection process.]



THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

## [DISINFECTION] DIRECTIONS FOR USE

[**Note to Reviewer:** Contact time can be presented in 1 [one] min[ute] or 60 sec[onds] on Production Label.]

### **FOR USE AS A DISINFECTANT ON NON-INSTRUMENT SURFACES:**

Spray CaviCide 1 6-8 inches from the surface, visibly wetting area to be disinfected. Allow surface to remain visibly wet for 1 minute [<sup>^</sup>3 minutes for Adenovirus II] at room temperature (68 °F/20 °C). Wipe surface using a clean towel or allow to air dry. Visible soil loads must be removed prior to disinfecting. [Use 1-minute contact time to disinfect against SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces.] [To kill SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces, use 1-minute contact time.]

### **FOR USE AS A DISINFECTANT ON NON-CRITICAL MEDICAL DEVICES<sup>1</sup>, INSTRUMENTS AND IMPLEMENTS:**

Instruments must be thoroughly cleaned to remove excess organic debris, rinsed and dried. Thoroughly clean and rinse lumens of hollow instruments. Spray CaviCide 1 6-8 inches from the surface, visibly wetting area to be disinfected. Allow surface to remain visibly wet for 1 minute [<sup>^</sup>3 minutes for Adenovirus II] at room temperature (68 °F/20 °C). Remove and rinse instruments. Wipe dry prior to use. [Use 1-minute contact time to disinfect against SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces.] [To kill SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces, use 1-minute contact time.]

<sup>1</sup> Non-critical medical devices are items that come in contact only with intact skin.

[Note: Critical and semi-critical devices must be followed by appropriate terminal sterilization/high-level disinfection process.]

***[The following instruction is to be included when hard, non-porous exterior surfaces of ultrasound transducers or probes are listed as a use site.]***

[To disinfect external surfaces of ultrasound transducers: Pre-clean ultrasound gel from external transducer to be disinfected [by using a clean wash cloth]. Spray CaviCide 1 directly onto surface, visibly wetting area to be disinfected. Allow surface to remain visibly wet for 1 minute [<sup>^</sup>3 minutes for Adenovirus II] at 68 °F/20 °C.] [Use 1-minute contact time to disinfect against SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces.] [To kill SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces, use 1-minute contact time.]

***[The following paragraph will be included in its entirety when the production label references health care use.]***

[This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization/high-level disinfection.]

[Always refer to manufacturer's care instructions before using on equipment -or- devices.]

***[The following HIV-1 / HBV / HCV instructions are to be included on production labels when claims against HIV-1, HBV or HCV are listed.]***

**CaviCide 1** effectively kills [Human Immunodeficiency Virus Type 1 (HIV-1),] [Hepatitis B Virus (HBV),] [[and] Hepatitis C Virus (HCV)] on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects can be associated with the potential for transmission of [HIV-1,] [HBV,] [[and] [HCV].

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST [HIV-1,] [HBV,] [[AND] HCV] OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUID:**

**Personal Protection:** Wear barrier protection [such as]: protective gloves, gowns, masks and eye coverings. **Cleaning Procedure:** Blood and other bodily fluids must be thoroughly cleaned from surfaces and objects before disinfection with CaviCide 1.

**Contact Time:** A contact time of 1 minute is required for HIV-1, HBV and HCV. [This contact time will not be sufficient for some of the organisms listed on the label. Refer to label for contact times against other organisms.]

**Disposal Infectious Materials:** Dispose of material in accordance with Federal, State or local regulations for infectious materials disposal.

**‡SPECIAL LABEL INSTRUCTIONS FOR CLEANING PRIOR TO DISINFECTION AGAINST *CANDIDA AURIS*:**

**Personal Protection:** Wear appropriate barrier protection [such as]: gloves, gowns, masks, or eye covering.

**Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with the product. Pre-cleaning is to include vigorous wiping and/or scrubbing and all visible soil is removed. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading the organism. Restrooms are to be cleaned last. Do not reuse soiled cloths.

**Infectious Waste Disposal:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

***The following statement is to be included on production labels when food contact areas or surfaces are listed.***

+ A thorough rinse with potable water is required for surfaces in direct contact with food.

**CLEANING AND STAIN REMOVAL DIRECTIONS FOR USE:**

**To clean -and/or- remove stains:** Spray surface and; [rinse and] let air dry. [If streaking is observed, wipe with a clean, damp [cloth] [or] [paper towel].] [[Works on] [effective against] [removes the following stains:]] [CLEANING] [and] [DEODORIZING]: [Spray surface and wipe with a paper towel until visibly clean. Let air dry.]

<p><b>Claims:</b></p> <ul style="list-style-type: none"> <li>• Cleans</li> <li>• Removes stains</li> <li>• Removes (<i>insert stain(s) /soil(s) shown to the right:</i>)</li> </ul>	<p><b>Stains/ Soils</b></p> <ul style="list-style-type: none"> <li>• Blood</li> <li>• Body oils</li> <li>• Bodily fluids</li> <li>• Coffee [stains]</li> <li>• Common soils</li> <li>• Dead skin</li> <li>• Dirt</li> <li>• Fingerprints</li> <li>• Food+ residue[s] -or soil[s]</li> <li>• Food stains</li> <li>• Fruit [stains]</li> <li>• Grease</li> <li>• Grime</li> <li>• Laboratory stains</li> </ul>	<p><b>Stains/ Soils continued</b></p> <ul style="list-style-type: none"> <li>• Mucus</li> <li>• [Other] common soils -and/or- stains</li> <li>• [Other] organic matter</li> <li>• Phlegm</li> <li>• Protein-based matter</li> <li>• Saliva</li> <li>• Spit</li> <li>• Spittle</li> <li>• Sputum</li> <li>• Throw-up</li> <li>• Tea [stains]</li> <li>• Urine</li> <li>• Vomit</li> </ul>
---	--	--

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

### ALLERGEN DESTRUCTION

To **[clean and] destroy specified allergens:** Spray CaviCide 1 6-8 inches from the surface. [Visibly wetting area to be disinfected.] Allow surface to remain visibly wet for 1 minute [Rinse.] Let air dry. [If streaking is observed, wipe with a clean, damp [cloth or] paper towel after contact time has expired.]

**Claims:**

- Breaks Down Many Allergens[!]
- Destroys -or- deactivates dust mite matter
- [Penetrates and] breaks down dust mite matter -or- particles, cockroach matter -or- particles, and pet dander allergens [found in dust]

***The following statement is to be included on production labels when animal housing facilities are listed (Including Veterinary Clinics, Farms and Poultry Houses)***

For cleaning and disinfecting the following hard nonporous surfaces: feeding and watering equipment, utensils, instruments, cages, kennels, stables, catteries, etc.

1. Remove all animals and feeds from premises, animal transportation vehicles, crates, cages etc.
2. Remove all litter, droppings and manure from floors, walls and surfaces of facilities (including pens, chutes, and barns) occupied or traversed by animals.
3. Empty all troughs, racks, and other feeding and watering appliances.
4. Thoroughly clean all surfaces with soap or detergent and rinse with potable water. Saturate surfaces by spraying with product for a period of 1 minute. [3 minutes for Adenovirus Type II] [1 minute for SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces]
5. Spray all handling and restraining equipment such as leashes, muzzles, halters or ropes, as well as forks, shovels, and scrapers used for removing litter and manure.
6. Ventilate buildings and other closed spaces.
7. Do not house animals or employ equipment until treatment has been absorbed, set or dried.
8. Thoroughly scrub treated feed racks, troughs, mangers, automatic feeders, fountains, and waterers with soap or detergent, and rinse with potable water before reuse.

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

## DISINFECTION [Organism List]

**[DESCRIPTION]** [CaviCide 1 is a multi-purpose [disinfectant]/ [cleaning product] [cleaner] for use on hard, non-porous inanimate surfaces.]

[CaviCide 1 is effective against the following microorganisms on hard, non-porous surfaces when used as directed:]

[This product has demonstrated effectiveness against the following organisms with a 1-minute contact time: **(insert organisms)**] [3 minutes for Adenovirus II]

[Disinfects **(insert organisms)** on hard non-porous environmental surfaces listed on this label] [When used as directed, CaviCide 1 provides effective Tuberculocidal / Virucidal\*<sup>[^]</sup> / Bactericidal / Fungicidal activity against the following organisms:]

**[Note to Reviewer:** the following table is optional on Production Labels.]

<b>Pathogen:</b>	<b>Bacteria, Enveloped Viruses*, Bloodborne Pathogens**, TB, Fungicidal, SARS-CoV-2*, Measles virus*, Candida auris‡</b>	<b>Non-enveloped Virus*<sup>[^]</sup></b>
Contact Time:	1 min[ute]	[1 min[ute] <sup>[^]</sup> 3 min[utes]

<b>ORGANISMS:</b>	<b>[Strain and/or ATCC number]</b>
<b><i>Mycobacterium Bacteria (1 Minute Contact)</i></b>	
[ <i>Mycobacterium tuberculosis</i> var: bovis (BCG) (TB)]	[Strain: Organon Teknika, Corp.]
<b><i>Bacteria (1 Minute Contact)</i></b>	
[ <i>Acinetobacter baumannii</i> ]	[ATCC 19606]
[ <i>Burkholderia cepacia</i> ]	[ATCC 25416]
[ <i>Enterobacter cloacae</i> ]	[ATCC 13047]
[ <i>Klebsiella pneumoniae</i> ]	[ATCC 4352]
[ <i>Pseudomonas aeruginosa</i> ]	[ATCC 15442]
[ <i>Salmonella enterica</i> ]	[ATCC 10708]
[ <i>Staphylococcus aureus</i> ]	[ATCC 6538]
<b><i>[Pathogenic] Fungi (1 Minute Contact)</i></b>	
[ <i>Candida albicans</i> ]	[ATCC 10231]
[ <i>Candida auris</i> ‡]	[AR-BANK#0381] [from CDC]
[ <i>Trichophyton interdigitale</i> ] [(previously known as <i>Trichophyton mentagrophytes</i> )]	[ATCC 9533]
<b><i>Drug-Resistant Bacteria (1 Minute Contact)</i></b>	
[Carbapenem-Resistant <i>Klebsiella pneumoniae</i> [(CRKP)]] [CRE]	[ATCC BAA-2470]
[ESBL <i>Escherichia coli</i> ] [ESBL <i>E. coli</i> ]	[ATCC BAA-196]
[Methicillin Resistant <i>Staphylococcus aureus</i> (MRSA)]	[ATCC 33592]
[Methicillin Resistant <i>Staphylococcus epidermidis</i> (MRSE)]	[ATCC 51625]
[Multi-drug resistant [MDR] <i>Acinetobacter baumannii</i> ]	[ATCC BAA-1605]
[Vancomycin Resistant <i>Enterococcus faecalis</i> (VRE)]	[ATCC 51299]
[Vancomycin Intermediate <i>Staphylococcus aureus</i> ] [VISA]	[ATCC 700699]

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

<b><i>ORGANISMS continued</i></b>	<b>[Strain and/or ATCC number]</b>
<b><i>Viruses* [Non-enveloped], [Virucidal*<sup>[^]</sup>] [†:] (1 Minute Contact) (^3 minute Contact)</i></b>	
[Adenovirus II* <sup>[^]</sup> ] (3 minutes)	[ATCC VR-846]
[Feline Calicivirus (surrogate for Norwalk and Norwalk-like viruses)*] [Norovirus*] [Norwalk Virus*] [Norwalk-like virus*] [Feline Calicivirus*] (1 minute)	[ATCC VR-782]
[Rotavirus*] (1 minute)	[ATCC VR-2018]
<b><i>Viruses* [Enveloped], [Virucidal*<sup>[^]</sup>] [†:] (1 Minute Contact)</i></b>	
[Duck Hepatitis Virus Surrogate for [Hepatitis B Virus (HBV)]*] [**]	[Strain Grimaud]
[Bovine Viral Diarrhea Virus Surrogate for [Hepatitis C Virus (HCV)]*] [**]	[NADL]
[Herpes Simplex Virus Type 1*]	[ATCC VR-260]
[Herpes Simplex Virus Type 2*]	[ATCC VR-734]
[Human Coronavirus (not associated with Severe Acute Respiratory Syndrome or SARS)*]	[ATCC VR-740]
[Human Immunodeficiency Virus (HIV-1)*] [**]	[ZeptoMetrix]
[Influenza A Virus, Strain A/ [Hong Kong 8/68] [H3N2 Virus]]*	[Charles River Labs]
Severe Acute Respiratory Syndrome-related Coronavirus 2* -or- SARS-CoV-2* -or- COVID-19 Virus*	[Strain: USA-WA1/2020] [BEI Resources, NR-52281]
Measles virus*	[Strain: Edmonston] [ATCC VR-24]

\* Virus

[\*\* Bloodborne pathogens]

[^ Virucidal – One minute contact time will not be sufficient for Adenovirus II.] [**Note to Reviewer:** this definition will only be included when Adenovirus Type II is listed Production Labels.]

† [One-Step Cleaning/ Disinfection] when Disinfection Directions are followed.

‡ [*Candida auris* –] Follow [the] Special Instructions for Cleaning required prior to Disinfection [listed on label].

[**Note to Reviewer:** Contact time can be presented in 1 [one] min[ute] or 60 sec[onds] on Production Label.]

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**PRODUCT CLAIMS**

▲ to represent any use site(s) listed in below as applicable.

◆ to represent any surface(s) listed in below as applicable.

[**Note to Reviewer:** Contact time can be presented in 1 [one] min[ute] or 60 sec[onds] on Production Label.]

- [99.9% -or 3 [three] log reduction against/for SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces [in [just] 1 minute] [in the presence of [5%] visible soil] [in one step]]
- [Kills 99.9% -or- 3 [three] log of SARS-CoV-2\* -or- COVID-19 Virus\* on hard, non-porous surfaces [in [just] 1 minute] [in one step] [in the presence of [5%] visible soil]]
- [99.999% -or 5 [five] log reduction for [*Candida auris*†] -or- [*C. auris*†] [in [just] 1 minute -or- 60 seconds]
- [Kills a minimum of 99.999% -or- 5 [five] log of *Candida auris*† -or- *C. auris*† on hard, non-porous surfaces] [in [just] 1 minute -or- 60 seconds]
- [1 Step Germicidal Spray]
- [1 Step Germicidal Spray against/ for SARS-CoV-2\* or- COVID-19 Virus\* on hard, nonporous surfaces]
- [1 minute] [3 minutes for Adenovirus II]
- [1 minute kill time] [3 minutes for Adenovirus II]
- [1 minute contact time] [3 minutes for Adenovirus II]
- [1 minute contact/kill time against/for SARS-CoV-2\* or- COVID-19 Virus\* on hard, nonporous surfaces [in the presence of [5%] soil load] [in one-step]]
- [A Metrex Research® brand]
- [A Metrex® brand]
- [A CaviCide Brand]
- [Advanced Cleaning Formula]
- [[Aids] [Helps] in reducing -or- in the reduction of cross-contamination between treated hard, non- porous surfaces]
- [[Aids] [Helps] in reducing -or- in the reduction of cross-contamination caused by SARS-CoV-2\* or- COVID-19 Virus\* between treated hard, non- porous surfaces]
- [An effective method to clean against odor-causing organisms on hard, non-porous surfaces listed on this label].
- [[A] one-step† cleaner/disinfectant]
- [[A] one-step† cleaner/disinfectant against/for SARS-CoV-2\* or- COVID-19 Virus on hard, nonporous surfaces]
- [[A] One Step† disinfectant against [*Candida albicans*] (and/or) [*Candida auris*†]]
- [A disinfectant product for use [on ◆] [and other ◆] [in ▲]]
- [Always ready-to-use]
- [An effective method to deodorize and disinfect against odor-causing organisms on hard non- porous surfaces]
- [An effective way to disinfect hard, non-porous environmental surfaces]
- [An effective way to disinfect hard, non-porous environmental surfaces listed on this label.]
- [A disinfectant for use [in ▲] [on ◆ [surfaces]] and other critical care areas where control between hard non-porous treated surfaces is required]
- [Antibacterial]
- [Antimicrobial]
- [Fungicidal]
- [Appropriate] [Perfect] [Ideal] for [Occupied Room] [and] [Terminal Room] Cleaning]
- [Assures product effectiveness and standardizes technique]
- [Bacteria -and/or- virus\* -and/or- germ killing [formula]]
- [Bactericidal] [Virucidal\*[^3 minutes for Adenovirus II]] [Fungicidal] [Tuberculocidal] [in 1 minute]
- [Bactericide -or- Bactericidal] [(in 1 minute)]
- [Bathroom disinfectant]
- [Barber Disinfectant/Cleaner]
- [Broad spectrum disinfectant]
- [CaviCide 1 kills germs that cause odors]
- [CaviCide 1 is effective against odor causing germs]
- [Can help reduce cross contamination between treated, hard non-porous surfaces]
- [Cleaner and disinfectant] [in One]
- [[Cleans and] disinfects [◆] [in ▲]]
- [Clean and disinfect hard, non-porous surfaces easily]
- [[Cleans] [and] [disinfects] [on the go]]
- [Cleans and/or deodorizes and/or disinfects [hard, nonporous surfaces]]

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**PRODUCT CLAIMS continued**

- [Cleans, deodorizes and kills common household germs]
- [Cleans] [disinfects] [Cleaner and disinfectant] [in one] [all in one [easy] step†]
- [Clean Home]
- [Cleans the toughest soils]
- [Cleaning Made Easy]
- [[Cleans [and /or Disinfects] [grocery cart handles or shopping cart handles or grocery cart child seats or shopping] cart child seats]
- [Cleans grocery carts and shopping carts]
- [Clinical Clean]
- [Clipper Blade Disinfectant/Cleaner]
- [Convenient]
- [Convenient and easy cleaning all around [the office] [business] [workplace]]
- [[Convenient] [Portable] way to [disinfect] [clean] hard non-porous surfaces [of ♦] [in ▲]]
- [Cuts Cleaning time]
- [Deep cleaning] [formula]] [Delivers a clean you can see and smell]
- [Dental Office Disinfectant/Cleaner]
- [Demonstrates 99.999% -or- 5 [five] log reduction against for [*Candida auris*‡] -or- [*C. auris*‡] [in [just] 1 minute]
- [Demonstrates 99.9% -or- 3 [three] log reduction against SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces [in [just] [1] [one] minute] [in one step] [in the presence of [5%] visible soil]]
- [Designed for killing SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces [of ♦] [in ▲] [in [just] [1] [one] minute] [in one step] [in the presence of [5%] visible soil]]
- [Destroys] [Deactivates] [Disinfects] [[Helps] Eliminate] [99.9%] [tough-to-kill] SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces [of ♦] [in ▲] [in [just] [1] [one] minute] [in one step†] [in the presence of [5%] visible soil]
- [Designed for Home Healthcare -and/or- patients -and/or caregivers]
- [Designed -or- developed for Home Caregivers]
- [Designed -or- developed for use [on ♦ [surfaces]] [in [▲]] **[Note to Reviewer:** this claim will not be listed on End-Use Product Labels.]
- [Delivers consistent [amount of solution] every time [!]]
- [Delivers consistent concentration every time [!]]
- [Deodorizes, and kills common household germs]
- [Deodorizes by killing the germs that cause odors]
- [Deodorize hard, non-porous surfaces] [♦] [in ▲]]
- [Designed for the [▲] [industry]]
- [Designed] [for] [with] [the] Environmental Services Professional[s] [in] [mind]!!]
- [Disinfectant]
- [Disinfect]
- [Disinfect [external] hard, non-porous surfaces of ♦ [in ▲]]
- [Disinfectant cleaner [that cleans, disinfects and deodorizes] in one step†]
- [Disinfecting formula]
- [[Disinfects] [Disinfectant] [Cleaner] on the go]
- [Disinfects [and deodorizes] [♦] [in one step†] [in ▲]]
- [Disinfects [bacteria] [viruses\*] [[multi] drug-resistant [bacteria] [pathogens]] [bloodborne pathogens\*\*] [[pathogenic] fungi] [*Candida auris*‡ -or- *C. auris*‡] [in [just] 1 minute]
- [Disinfects and deodorizes in 1 minute] [3 minutes for Adenovirus II]
- [Disinfects hard, non-porous surfaces in 1 minute] [3 minutes for Adenovirus II]
- [Disinfects hard non-porous environmental surfaces listed on this label [Use to]]
- [[Disinfects] hard surfaces [that people touch]]
- [Disinfects [in [just] one [1] minute] [3 minutes for Adenovirus II]]
- [[Disinfects and deodorizes by killing germs and their odors]
- [[Disinfects hard, nonporous surfaces [♦] [throughout the ▲]]
- [Disinfects [as it cleans] [and scrubs in one easy step†]]
- [Does not contain dyes]

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**PRODUCT CLAIMS continued**

- [Easy† [and convenient] to use]
  - [[Easy†] [and convenient] [way to clean] [and] [disinfect]]
  - [Easy One Step† Cleaning [and disinfecting] (except when precleaning is required for visible soil and blood and body fluids)]
  - [Easy† [and convenient] way to disinfect [bathrooms] [office] [nursing stations] [high touch] hard, non-porous surfaces everyday]
  - [Easy to handle]
  - [Easy† way to [clean] [and] [disinfect] [bathroom] [or exterior toilet surfaces] hard, non-porous surfaces]
  - [Easy† and convenient disinfecting [♦ [surfaces]] [throughout the ▲]]
  - [Easy, one-step† cleaning and disinfecting]
  - [Easy, one-step† [cleaning and] disinfecting against SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces]
  - [Easy† [and convenient] way to [touch-up] [bathroom] [office] everyday]
  - [Economical]
  - [Economy]
  - [Effective -or- Activity -or- For daily use against (*Insert any organism(s) from Organism List*) [and] *Candida auris*† -or- *C. auris*†] [in [just] [1] [one] minute]
  - [Effective -or- Activity -or- For daily use -or- Disinfects against SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces [in [just] [1] [one] minute] [in one step†] [in the presence of [5%] visible soil]]
  - [Effective against Influenza A and H1N1 [Influenza A Virus]
  - [Effective against (*Insert any organism from list of organisms*)
  - [Effective against bathroom bacteria and viruses\* [use at least two bacteria and two viruses from this list:]] [ESBL *E. coli*] [*Salmonella enterica*] [*Staphylococcus aureus*] [Herpes Simplex Virus Type 2\*] [Hepatitis B Virus (HBV)\*] [Hepatitis C Virus (HCV)\*] [Adenovirus II] [Influenza A, H3N2 Virus\*] [*Enterobacter cloacae*] [*Burkholderia cepacia*] [*Carbapenem-Resistant Klebsiella pneumoniae* (CRKP)(CRE)]
  - [Effective against -or- Kills a range of bacteria including *Staphylococcus aureus*, *Salmonella enterica*, *Pseudomonas aeruginosa* [in the presence of organic soil [(5% blood serum)]]
  - [Effectively disinfects against MRSA]
  - [[Effective against -or- Kills (*insert any organism(s) from organisms list*) [in the presence of organic soil [(5% blood serum)]]
  - [Effective in the presence of 5% blood serum contamination]
  - [Effective against flu virus (Influenza A)]
  - [Effective against flu virus (Influenza A, H3N2 Virus)]
  - [Effective daily cleaner]
  - [Effectively disinfects against MRSA]
  - [Effectively disinfects against Norovirus]
  - [[Effectively] disinfects hard, nonporous [environmental] surfaces]
  - [Effectively removes everyday soils]
  - [Efficacy tests have demonstrated that [Product Name] [this product] [are] [is] an effective bactericide, virucide\*, tuberculocide, and fungicide in the presence of organic soil (5% blood serum)]
  - [Efficacy tests have demonstrated that [Product Name] [this product] [is] an effective virucide\* against/for SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces [in the presence of [5%] visible soil] [in [1] [one] minute]
  - [Eliminates odors]
  - [Eliminates odors from fecal matter and urine]
  - [Eliminates airborne particles associated with sprays]
  - [Eliminates [daily] mixing and measuring]
  - [Eliminates [lingering] [food+] odors [on surfaces]]
  - [[Eliminates odors at their source [caused by germs -or- bacteria-and/or- mildew]]
  - [[Eliminates -or- Reduces odors [caused by germs -or- bacteria -and/or- mildew]]
- [Note to Reviewer:** Claims above with “Eliminates” are NOT associated with pathogens but rather with odor or mixing / measuring; as such, these claims do NOT need to be qualified with “99.9%”.]



**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**PRODUCT CLAIMS continued**

- [EVS [Environmental Services] [use]]
- [Exhibits effective virucidal\* activity]
- [Kill[s] -and/or- Effective against 99.9% [of]] [ESBL *Eschericia coli*] [ESBL *E. coli*]]
- [Kill[s] -and/or- Effective against 99.9% [of]] [*Salmonella enterica*] [(*Salmonella*)]]
- [Kill[s] -and/or- Effective against flu (Influenza A) [virus\*] [on environmental - and/or- hard, nonporous surfaces]]
- [Kill[s] -and/or- Effective against the flu (Influenza A) [virus\*[es]]]
- [Kill[s] -and/or- Effective against (*insert any organism from organisms list*)]
- [Kill[s] -and/or- Help[s] reduce[s] the spread of -and/or- Effective against SARS-CoV-2\* -or- COVID-19 Virus\*] [in [just] [1] [one] minute] [in one step] [in the presence of [5%] visible soil] on hard, nonporous surfaces]
- [For daily use]
- [For] general use
- [For [daily] [general] [professional] use [especially] [during this [COVID-19] pandemic -or- outbreak]]
- [For use on hard, non-porous bathroom surfaces]
- [For use on hard, non-porous nursery surfaces]
- [For use on [♦ [surfaces]] [in ▲]]
- [For use in critical care areas between treated hard, non-porous surfaces [where control of cross contamination is of prime importance]
- [For] Surface cleaning
- [For Use by Home Caregivers]
- [For Use in Home Healthcare]
- [Formulated -or- Designed to] remove[s] biological odors including Urine -and/or- Feces -and/or- Vomit]
- [Fragrance Free]
- [Freshens]
- [Free of fragrance]
- [Fresh smell] [Helps reduce cross contamination between treated hard, non-porous surfaces]
- [Fungicide]
- [Fungicidal [disinfectant]] [against *T. interdigitale* -and/or- *C. auris*† -and-or *Candida albicans*]]
- [Fungicidal in 1 minute]
- [Germicide -or- Germicidal]
- [Gets rid of [or eliminates] odors] [**Note to Reviewer:** this “eliminates” claim is NOT associated with pathogens but “odors” only; as such, it does not need to be qualified with “99.9%”]
- [Great for everyday use]
- [[Helps] [kill][s] [eliminate][s] odor-causing bacteria]
- [Helps reduce cross contamination between treated hard, non-porous surfaces]
- [Healthcare Disinfectant/Cleaner]
- [Hospital [disinfectant] [disinfection] [use only when hospital use sites on label]
- [Institutional [size]]
- [Intermediate Level Disinfectant]
- [Just] [Simply] [spray and wipe] [exterior toilet surface] [bathroom germs] [away]
- [Kills Adenovirus II]
- [Kills [the] [Pandemic] SARS-CoV-2\* -or- COVID-19 Virus\* [in [just] [1] [one] minute] on hard, nonporous surfaces [in one step] [in the presence of [5%] visible soil]]
- [Kills [bacteria] [viruses\*] [[pathogenic] fungi]]
- [Kills bacteria in 1 minute]
- [Kills viruses\* in 1 minute [(Adenovirus type 2 in 3 minutes)]]
- [Kill [Effective against] [[pathogenic] fungi] [including] [*Candida albicans* -or- *C. albicans*] [*Candida auris*† -or- *C. auris*†] in [just] 1 minute]
- [Kills -or- Disinfects -or- Eliminates -or- [Is] Effective against Measles virus\* on hard, nonporous surfaces [in [just] [1] [one] minute] [in one step] [in the presence of [5%] soil load]]
- [Kills -or- Disinfects -or- Eliminates -or- [Is] Effective against *Candida auris*† -or- *C. auris*† [in [just] 1 minute] on hard, non-porous surfaces]

[**Note to Reviewer:** “Eliminates” will be qualified with 99.999% if used.]

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**PRODUCT CLAIMS continued**

- [Kills [Norovirus] [Norwalk virus]]
- [Kills HIV-1]
- [Kills HCV]
- [Kills HBV]
- [Kills Rotavirus]
- [Kills Adenovirus II]
- [Kills bacteria, viruses\*]
- [Kills [bacteria] [viruses\*] [[pathogenic] fungi] on hard, nonporous surfaces]
- [Kills] [germs] [or bacteria] [while it cleans] [that cause odor]
- [Kills odor causing bacteria]
- [Kills [Disinfects] Norovirus [Feline Calicivirus] [in] 1 [one] minute]
- [Kills -or- Effective against bacteria]
- [Kills -or- Effective against germs]
- [Kills -or- Effective against viruses\*]
- [Kills [or eliminates 99.9%] of bacteria]
- [Kills [or eliminates 99.9%] of viruses\*]
- [Kills [or eliminates 99.9%] of fungus]
- [Kills odor-causing bacteria -or- germs]
- [[Kills multi drug resistant organisms [(list organisms)]]
- [Kills multi drug resistant bacteria [(list organisms)] [in [1] [one] minute]]
- [[Kills] [germs] [or bacteria] [while it cleans] [that cause odor]
- [Kills] [common] [household] [bacteria] [viruses\*] [[pathogenic] fungi] [germs] **[(insert any organism from organism list)]** [in 1 minute] [3 minutes for Adenovirus II]
- [Kills] [Eliminates 99.9%] [(list organisms)] [in [just] [1] minute] **[(insert any organism from organism list)]** [in 1 minute] [3 minutes for Adenovirus II]
- [Kills -or- Destroys -or- Disinfects -or- Eliminates 99.9%-or- Effective against SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces [in [just] [1] [one] minute] [in one step†] [in the presence of [5%] soil load]]
- [Kills] [the] flu virus (Influenza A, H3N2 Virus)]
- [Kills] [viruses\* and bacteria [including the flu virus (Influenza A, H3N2 Virus)]]
- [Kills TB in 1 Minute]
- [Kitchen+ and Bath]
- [Kitchen+ and Bath Solution for Hard, Non-porous Surfaces]
- [Large[r] Surface[s]]
- [Lead Free] **[Note to Reviewer: this claim will not be listed on End-Use Product Labels.]**
- [Leaves bathroom fresh and clean in [just] [1] minute]]
- [Leaves [bathroom] [office] clean [and disinfected] [and smelling fresh]
- [Leaves nothing behind]
- [Lightly scented -or- quick-dissipating scent]
- [Making home healthcare [more manageable] [easier] [simpler]]
- [Maximum cleaning, minimum effort]
- [Meets surface disinfection recommendations from OSHA['s] Bloodborne Pathogen Guidelines -or- Standards]
- [Mercury free] **[Note to Reviewer: this claim will not be listed on End-Use Product Labels.]**
- [Multi-purpose disinfectant [for Health Care Facilities, Dental Facilities and Residential purposes]]

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**PRODUCT CLAIMS continued**

- [Multi-surface]
- [Neutralizes odors [doesn't [just] mask them]
- [No canister needed]
- [No fragrance added]
- [No measuring] [required]
- [No mixing] [required]
- [Non-abrasive formula]
- [[Occupied Room] [and] [Terminal Room] Cleaning Product]
- [Odor [elimination] [neutralization]
- [On the go disinfectant]
- [One-step† cleaner and disinfectant]
- [One step† [cleaning] [disinfecting] [product]] [liquid]]
- [One-step† cleaning and disinfecting [saves time]]
- [One step† cleaning [and/or disinfecting] [and/or deodorizing] disinfect and/or touch-up clean]
- [One step† cleaning [and/or disinfecting] [and/or deodorizing] disinfect and/or touch-up clean] [cleans and/or disinfects and/or deodorizes] [anytime] [daily] [every day] [in one easy step] on hard non-porous surfaces [of ♦] [in ▲] [bathroom and/or kitchen surfaces] [on household surfaces] [on ♦ [surfaces]]
- [One step† disinfectant cleaner that is effective against a broad spectrum of bacteria]
- [One step germicide]
- [One-step† disinfectant cleaner designed for general cleaning and disinfecting hard, nonporous environmental surfaces in ▲ [on ♦ [surfaces]]]
- [Perfect for cleaning [and disinfecting] [♦ [surfaces]] [in ▲]]
- [Perfect for disinfecting exterior toilet surfaces – or – around the exterior toilet surface]
- [Picks up [particles] [dust] [dirt] [hair] [crumbs]]
- [Portable]
- [Powers through tough messes]
- [Ready to Use [RTU] prevents formulation errors] [Premixed] [Premeasured] [Easy] [and convenient] way to disinfect]
- [[Product Name] are useful in a variety of [healthcare], [nursery/daycare] [and] [▲ [on ♦ [surfaces]]] applications listed on the [this] label.]
- [Easy† way to clean [and disinfect] [bathroom] [or exterior toilet surfaces] hard, non-porous surfaces]
- [Ready to use]
- [[Ready-to-use] disinfectant [for use in ▲ [on ♦ [surfaces]]]
- [Ready-to-use [one-step†] disinfectant [disinfection] for [killing] SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces]
- [Ready-to-use [one-step†] disinfectant [disinfection] for [killing] (***insert any organism from organism list***) [*Candida auris*† -or- *C. auris*†]]
- [Recommended for use in ▲ [on ♦ [surfaces]]]
- [Recyclable] [package]
- [Reduces cross contamination between treated hard, non-porous surfaces]
- [Reduces cross-contamination between treated hard, non-porous surfaces [from *Candida auris*† -or- *C. auris*†] [from (***insert any organism from organism list***) [in 1 minute)] [3 minutes for Adenovirus II].
- [Reduce [cross-contamination of] SARS-CoV-2\* -or- COVID-19 Virus\* between treated, hard, non-porous surfaces [of ♦] [in ▲]]

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**PRODUCT CLAIMS continued**

- [Reduces] [packaging] waste [by X%] [when compared to regular canister or bucket or pail]
- [Requires no mixing or measuring and is ready-to-use]
- [Salon Disinfectant/Cleaner]
- [Sample]
- [Saves Time]
- [School Disinfectant/Cleaner]
- [Shelf-stable formula]
- [So] [surfaces are more than [just] [clean] [they're] [disinfected]
- [Specifically -or- Specially designed for the Medical/Dental/ Veterinary/Health Care [▲] industry]
- [Starts cleaning on contact]
- [Suitable for use [on ♦ [surfaces]] [in ▲]]
- [Suitable for use on [titanium-coated] [medical grade] stainless steel]
- [Suitable for use on hard, nonporous inanimate environmental surfaces]
- [Surface Disinfectant/Cleaner]
- [Strong cleaning power]
- [The 1-minute [germicide] [disinfectant] [3 minutes for Adenovirus II]]
- [The 1-minute germicide/ disinfectant against *Candida auris*† -or- *C. auris*† (*insert any organism(s) from Organism List*)]
- [The answer to your [disinfecting] [cleaning] needs]
- [The] smell of clean]
- [[The] easy† -and/or- convenient way to disinfect] [with CaviCide™ 1/ this product]
- [This product was tested according to AOAC efficacy testing methods -or- standards for hospital disinfection]
- [This product has demonstrated effectiveness against Influenza A Virus]
- [This product has demonstrated effectiveness against SARS-CoV-2\* -or- COVID-19 Virus\* [and is expected to inactivate] SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces [in 1 minute] [in the presence of [5%] visible soil] [in one-step]]
- [This product has demonstrated effectiveness against Adenovirus II with a 3 minute contact time]
- [Time to reorder]
- [Trial]
- [Tuberculocidal (1 minute)]
- [Twin pack]
- [Two in one action [cleans and disinfects in one easy step†]]
- [Tuberculocide -or- Tuberculocidal]
- [To kill (*insert any organism(s) from Organism List*)], let stand [undisturbed] for 1 minute]
- [To kill *Candida auris*† -or- *C. auris*† [,] [on treated surfaces especially during outbreaks] [of *Candida auris*† -or- *C. auris*†] [let stand [undisturbed] for 1 minute]]
- [To kill SARS-CoV-2\* -or- COVID-19 Virus\* on treated hard, non-porous surfaces [especially during this/ the [COVID-19] outbreak] [,] [let stand [undisturbed] for 1 minute]]
- [To reduce *Candida auris*† -or- *C. auris*† on hard, non-porous treated surfaces [of ♦] [in ▲]]
- [Usable on variety of soils]
- [Use in [public] [common places] [in ▲] [on ♦ [surfaces]] where [bacteria -and/or- viruses\*-and/or- germs may be of concern] [on hard, nonporous surfaces]
- [Use to] [Deodorize hard, non-porous surfaces]
- [Use to pre-clean critical or semi-critical medical devices prior to sterilization / high-level disinfection]

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**PRODUCT CLAIMS continued**

- [Use this product to] [Kill bacteria in 1 minute] [Kill viruses\*[^] in 1 minute] [3 minute contact time for Adenovirus II] [Kill TB in 1 minute] [Kill TB/HIV-1 in 1 minute] [Kill *Candida auris*† -or- *C. auris*† in 1 minute] [Kill SARS-CoV-2\* -or- COVID-19 Virus\* in 1 minute on hard, nonporous surfaces]
- [Use where cross-contamination between treated hard, non-porous surfaces is of prime importance]
- [Use on hard, non-porous surfaces only: [stainless steel] [Formica®] [Porcelain] [Plastic] [Glazed Ceramic] [glass] [glass tables] [glazed tile] [Test on small inconspicuous area first.]
- [Value]
- [Versatile]
- [Veterinary Disinfectant/ Cleaner]
- [[Virucide\*[^] -or- Virucidal\*[^]] [in 1 min] [3-minute contact time for Adenovirus II]
- [Virucide\* -or- Virucidal\* against SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces [in [just] 1 minute] [in the presence of [5%] visible soil] [in one-step†]]
- [Where cross contamination between treated surfaces is of prime importance]
- [Disinfects] [Kills] [common] [household] [bacteria] [viruses\*] [germs]
- [Disinfects -or- Eliminates 99.9%] [Kills] [list organisms] in [just] [1 minute] [3 minutes for Adenovirus II]
- [Disinfects] [Kills] [the] flu virus\* (Influenza A, H3N2 Virus)
- [Disinfects] [Kills] [viruses\* and bacteria] [including the flu virus (Influenza A, H3N2 Virus)]
- [[Disinfects] [Kills] SARS-CoV-2\* -or- COVID-19 Virus\* on hard, nonporous surfaces [in one step] [in the presence of [5%] visible soil] [in [just] 1 minute]]
- [[X% reduction in packaging [material]]]

**[Note: The following claims can be used only for the first 6 months of the product on shelf]**

- [Better] [Improved] [New] [!]
- [Better] [Improved] [New] [New and improved] [formula] [packaging]
- [New] Brought to you by the makers of CaviCide [and CaviWipes]
- [New Claims]
- [New Claims and Contact Times]
- [New Claims – Kills TB in 1 Minute]
- [New and improved]
- [New] [RTU] [Ready To Use] [Solution] [Disinfectant]
- [New] [Surface Disinfectant]
- [New Kill Claims]
- [New – Saves Time]
- [New] – [Kills TB in 1 Minute]

**[Note – the chart shown below shall be considered optional – e.g. inclusive in a bracket]**

[When used as directed, CaviCide 1 provides effective Tuberculocidal\*\* / Virucidal\*[^] / Fungicidal activity against the following organisms:]

<b>Tuberculocidal</b>	<b>Associated with:</b>
<i>Mycobacterium tuberculosis</i> var: bovis (BCG) [(Organon Teknika, Corp.)] [(TB)]	Tuberculosis (BCG)
<b>*Viruses</b>	
Human Immunodeficiency Virus (HIV-1)*	AIDS
Influenza A, H3N2 Virus*	Common flu
Severe Acute Respiratory Syndrome-related Coronavirus 2* -or- SARS-CoV-2* -or- COVID-19 Virus*	COVID-19 [Virus]
Measles* virus	Measles
<b>Fungus/Yeast</b>	
<i>Trichophyton interdigitale</i>	Athletes Foot Fungus
<i>Candida auris</i> †	

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**Table 1 Medical Use Sites:  
CaviCide 1 is recommended for use on hard non-porous surfaces in [▲]  
[in health care settings] [such as (insert use site)]:**

3D Imaging Service Centers	Endoscopy [Centers] [Facilities] [Departments]	Ophthalmic [Facilities] [Offices]
Acute Care Facilities	Endoscopy Centers	Optometric [Facilities] [Offices]
Alternate Care Facilities	ENT Facility	Oral Health Care Facilities
Ambulances	Eye Care [Facilities] [Centers]	Orthodontic Offices
Ambulatory Surgery Centers (ASC)	Eye Surgical Centers	Orthopedic Offices
Animal Research Facilities	Forensic Labs	Outpatient Care Facilities
Assisted Living Facility	Full-Care Nursing Facility	Outpatient [Surgery] [Facilities] [Centers]
Autopsy Rooms	Gastrointestinal (GI) [Labs] [Laboratories]	Out-Patient Surgical Centers (OPSC)
Bathroom	Gastrointestinal (GI) Centers	Pathology Rooms
Bioskill Laboratories	Gastrointestinal (GI) Suites	Patient Care Areas
Blood Drawing Centers	Home [Health] Care [Services]	Patient Care Rooms
Blood Apheresis Facilities	Hospice	Patient Rooms
Blood Units	Hospital Kitchens	Pediatric Exam Rooms
Burn Center	Hydrotherapy [Facilities] [Rooms]	Pediatric Intensive Care Units (PICU)
Cardiac Intensive Care Unit (CICU)	In vitro Fertilizer Center (IVF)	Pediatric Offices
CAT Scan [Lab] [Laboratories]	Instrument Reprocessing [Facilities] [Departments] [Center]	Personal Care Service Facilities
Cath Labs	Intensive Care Units (ICU)	Pharmacy Clean Rooms
Central Service Areas	Isolation Areas	Phototherapy Rooms
Central Services Department	Isolation Rooms	Physical Therapy Facilities
Central Supply Areas	Labor Units	Physicians' Offices
Central Supply Rooms	[Lab] [Laboratories]	Plasma Apheresis [Centers] [Facilities] [Units]
Clean Rooms	Lasik Eye Centers	Post [Ops] [Operatory] Rooms
Clinical Areas	Lasik Units	Pre [Ops] [Operatory] Areas
Clinics	Laundry Rooms	Psychiatric Care Facility
Critical Care Units	Life Science Laboratories	Radiology Department
CS [Labs] [Laboratories]	Long-Term Acute Center (LTAC)	Recovery Rooms
Dialysis Centers	Long-Term Care [Facilities]	Rehabilitation Facilities
Dialysis Clinics	Medical Facilities [Clinics]	Respiratory Therapy Facilities
Dissection [Labs] [Laboratories]	Medical Related Facilities	Skilled Nursing [Facility] [Homes]
Doctor's Offices	Medivans	Sterile Processing Departments
Donation Centers [blood] [plasma] semen milk aphaeresis	Medical Intensive Care Unit (MICUs)	Sterile Storage Areas
Drug Preparation [Centers] [Facilities] [Rooms]	Magnetic Resonance Imaging (MRI) [facilities] [laboratories] [labs]	Surgery Centers
Doctor's Offices	Neonatal Intensive Care Units (NICU)	Surgery Rooms
Donation Centers [blood] [plasma] semen milk aphaeresis	Neonatal Nursery Areas	Surgical Centers
Drug Preparation [Centers] [Facilities] [Rooms]	Neonatal Units	Surgical Suites
Drug Preparation Rooms	Newborn Nursery	Teaching Hospitals
Electrolysis Facilities	Nuclear Medicine Facilities	Tissue Banks
Emergency Departments	Nurses' Stations	Transplant Center
Emergency Medical Services (EMS) Units	Nursing Homes	Trauma Center

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Emergency Medical Settings	Nursing Stations	University Research Lab
Emergency Rooms (ER)	Occupational Therapy [Centers] [Facilities]	Wound Center
Emergency Vehicles	Oncology Areas	X-Ray Rooms
Endo Suites	Operating Rooms	

**Table 2 Institutional Use Sites:**  
***CaviCide 1 is recommended for use on hard non-porous surfaces in [▲]***

Adult School	Early Learning Centers	Middle School
Air Transport Facilities	Elementary Schools	Military Installations
Airplane Interiors	[Fire] [Rescue Departments]	Nursery School
Airports	Fire Vehicles	Police Vehicles
Bathroom	Forestry Services	Pre-K School [Room]
Casino	Grammar School	Preparatory School
Child Care Facilities	High School	Preschool
Churches	Hospitality Facilities	Primary School
Colleges	Hotels	Prisoner Detention Facilities
Community Colleges	Industrial Facilities	Prisoner Lockups
Conference [Centers] [Rooms]	Infant Care Centers	Prisons
Conference Rooms Fixtures	Intermediate School	School Buses
Convention Facilities	Jails	School Kitchens+
Correctional Facilities	Junior High School	Schools
Daycare Centers	Kindergarten	[University] [Universities]
[Dormitories] [Dorm]	Learning Institutions	

**Table 3 Dental Use Sites:**  
***CaviCide 1 is recommended for use on hard non-porous surfaces in [▲]***

Dental [Laboratories] [Lab]	Dental Office	Dental [Operatory] [Operatories]
-----------------------------	---------------	----------------------------------

**Table 4 Vet Use Sites:**  
***CaviCide 1 is recommended for use on hard non-porous surfaces in [▲]***

Animal Care Facilities	Avian Units [Hatchery]	Pet Homes
Animal Control Vehicles	Bathroom	Pet Motels
Animal Dental Operatory	Doggy Day Care Centers	Pet Spas
Animal Detention Room	Feline Facilities	Pet Stores
Animal Research Facilities	Feline Units	Small Animal Facilities
Animal Shelters	Kennels	Veterinary Offices
Avian [Care] Facilities	Pet Areas	Veterinary Facilities

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**Table 5 Food Prep+ Use Sites:**

***CaviCide 1 is recommended for use on hard non-porous surfaces in [▲]***

Bakeries+	Food Processing Facilities+	Restaurants+
Banquet Hall Facilities+	Grocery Stores +	School Kitchens+
Butcher Shops+	Kitchen+	Supermarkets+
Food Preparation Areas+	Produce Areas+	

[+ A thorough rinse with potable water is required for surfaces in direct contact with food.]

**Table 6 Miscellaneous Use Sites:**

***CaviCide 1 is recommended for use on hard non-porous surfaces in [▲]***

Airplanes	Factories	Massage Facility
Arboretums	Foots Spas	Money Handling Facilities
Arcades	Foreign Current Exchange Facilities	Morgues
[Athletic] [Training] Facilities	Funeral Homes	Mortuaries
Auto Facilities	Gaming Facilities	Motels
Banks	Garden Centers	Movie Theaters
Barber Shops	Golf Courses	Outdoor Playground
Bathrooms	Grounds Keeping Services	Personal Care Service Facilities
Beauty Salons	Gymnasium	Public Transportation Vehicles
Body Piercing Salons	Health Club Facilities	Recreational [center] [facilities]
Body Service Facilities	Helicopter Interiors	Recycling Centers
Body Storage Facilities	Household Bathrooms	Salon Settings
Botanical Gardens	Housekeeping Rooms	Weather Shelters
Bowling Alleys	Indoor Playground	Shower Rooms
Buses	Laundromats	Skincare Salons
Cabs	Laundry Rooms	Storage [Areas] [Rooms]
Cadaver Processing areas	Lawn Services	Subway
Cars	Limousine	Tanning Salons
Condominiums	Living Rooms	Tattoo Parlors
Cruise Ships	Locker Rooms	Taxi
Currency Exchange	Lockup Facilities	Theme Parks
Day Spa	Manicure Salons	Transport Vehicles
Delivery Units	Manufacturing [Facilities] [Plant]	



**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**Table 7 Surfaces:**

***CaviCide 1 is recommended for use on hard non-porous surfaces [of] [such as] [♦]:***

Ambulance Equipment Surfaces	Biohazard Container	Combs
Animal Cages	Blood Pressure Cuff	Computer
Animal Equipment	Booster Chairs+	Computer Keyboards
Appliances+	Booster Seats+	Computer Mouse
Arm Rests	Brushes	Copy Machines
ATM Surfaces	Cabinet Handles	Counter Tops+
Autoclave Surfaces	Cabinets	Counters
Automatic Feeders	Cafeteria Trays+	Cribs+
Automatic Teller Machines (ATM)	Cages	Crutches
Baby Warmers	Cardiac Gym Equipment	Dental Chairs
Baby Cribs	Carts	Dental Countertops
Baked Enamel	Cash Register	Dental Curing Lights
Bars	Ceilings	For use as a cleaner and disinfectant on cured, set and hardened, hard non-porous [dental impressions] [dental impression materials]
Basins	Cell Phones	Dental Impression Trays
Baskets	Chairs	For use as a cleaner and disinfectant on hard nonporous [dental impressions] [dental impression materials] prior to sending to a dental lab
Bathing Units	Changing Tables	Dental Operatory Surfaces
Bathroom Fixtures	Chemotherapy Hoods	Dental Unit Instrument Trays
Bathtubs	Child Seats+	Dentist Chair
Bed Pans	Chrome	Desk Tops
Bed Railings [Bedside Rails]	Clipper Blades	Desks
Bedside Commode	Closet Handles	Diaper Changing Stations
Behind and Under Sinks	Coated Mattresses	Diaper Pails
Bidet	Coated Pillows	Digital Phosphor Plates
<ul style="list-style-type: none"> <li>• For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [♦ [surfaces] [in ▲]] [and other critical care areas where cross-contamination between treated hard, non-porous surfaces is important] [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.]</li> <li>• [Surfaces which are made of plastic (i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene), and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.]</li> <li>• + A thorough rinse with potable water is required for surfaces in direct contact with food.</li> <li>• + Any washable (food and non-food contact) surface where disinfection is required.</li> </ul>		

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**Table 7 Surfaces continued:**

Dining Tables+	Exterior Surfaces of Toilet	Exterior Surfaces of Endodontic Equipment (Apex Locators, Pulp Testers, Motors)
Dish Racks+	Exterior or Exposed Interior Surfaces of Medical Equipment	Exterior Surfaces of Radiology Equipment
Display Cabinets	Exterior Surfaces of Dental Equipment	Exterior Surfaces of Refrigerator Units+
Diving Equipment	Exterior Surfaces of Amalgamators	Exterior Surfaces of Respiratory Therapy Equipment
Diving Masks+	Exterior Surfaces of Anesthesia Equipment	Exterior Surfaces of Urinals
Diving Snorkels+	Exterior Surfaces of Aphaeresis Machines	Exterior Surfaces of Vanity Tops
Door Knobs	Exterior Surfaces of Appliances such as Refrigerators+, Microwaves+, Food Processors+ and Stoves+π	External Surfaces Emergency Resuscitative Equipment
Dopplers	Exterior Surfaces of Center Venous Lines	External Surfaces of Bladder Scanners
Drain Boards+	Exterior Surfaces of CICU Equipment	External Surfaces of Blood Glucose Meters
Drained Whirlpool Tank Surfaces	Exterior Surfaces of Defibrillators	External Surfaces of BP Monitors
Enamel	Exterior Surfaces of Dental Care Facilities Hand Pieces	External Surfaces of CAT Scan Machines
Endoscope Cabinet	Exterior Surfaces of Dental Equipment Lines	External Surfaces of Compression Devices
Equipment Surfaces	Exterior Surfaces of Dental Equipment Tubing	External Surfaces of Dialogistic Equipment
Exam Tables	Exterior Surfaces of Dental Operatories	External Surfaces of Dictating Equipment
Examination Table	Exterior Surfaces of Diagnostic Equipment	External Surfaces of Digital X-Ray Sensors
External Surfaces of Kendrick Extrication Device (KED)	Exterior Surfaces of Dialysis Machines	External Surfaces of EKG Machines
<ul style="list-style-type: none"> <li>• For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [♦ [surfaces] [in ▲]] [and other critical care areas where cross-contamination between treated hard, non-porous surfaces is important] [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.]</li> <li>• [Surfaces which are made of plastic (i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene), and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.]</li> <li>• + A thorough rinse with potable water is required for surfaces in direct contact with food.</li> <li>• + Any washable (food and non-food contact) surface where disinfection is required.</li> <li>• π Allow surfaces to come to room temperature prior to disinfection.</li> </ul>		

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**Table 7 Surfaces continued:**

External Surfaces of Electronic Thermometers	External Surfaces of Oscopes	External Surfaces of Ultrasound Transducers
External Surfaces of Electrosurgical Generators	External Surfaces of Patient Care Monitor Unit	External Surfaces of Vaporizer
External Surfaces of Glucometers	External Surfaces of Patient Monitoring Equipment	External Surfaces of Veterinary Equipment
External Surfaces of Headsets	External Surfaces of Probes	External Surfaces of Vital Sign Monitors
External Surfaces of ICU Equipment	External Surfaces of Pulp Tester Motors	External Surfaces of X-Ray Machines
External Surfaces of Incubators	External Surfaces of Pulse Oximeter	Eye Glasses (do not use on contact lenses)
External Surfaces of Infusion Pumps	External Surfaces of Sequential Compression Devices	Faucets
External Surfaces of Intubation Equipment	External Surfaces of Respiratory Ventilators	Fax Machines
External Surfaces of Isolettes	External Surfaces of Resuscitators	Feed Racks
External Surfaces of IV Pumps	External Surfaces of Sterilizers	Filing Cabinets
External Surfaces of Life Support Equipment	External Surfaces of Suction Pumps	Filing Doors
External Surfaces of Loupes	External Surfaces of Surgery Microscopes	Filing Handles
External Surfaces of Medical Equipment	External Surfaces of Surgery Robots	Fire Trucks
External Surfaces of Medical Warmers	External Surfaces of Surgical Headlights	Floors
External Surfaces of Monitors Electrodes	External Surfaces of Surgical Light Sources	Food Cases+
External Surfaces of MRI Machines	External Surfaces of Tympanic Thermometers	Food Serving Trays+
External Surfaces of OLED Monitors	External Surfaces of Ultrasound Machines	Formica
External Surfaces of Ophthalmoscopes	External Surfaces of Ultrasound Monitors	Fountains+
<ul style="list-style-type: none"> <li>• For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [♦ [surfaces] [in ▲]] [and other critical care areas where cross-contamination between treated hard, non-porous surfaces is important] [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.]</li> <li>• [Surfaces which are made of plastic (i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene), and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.]</li> <li>• + A thorough rinse with potable water is required for surfaces in direct contact with food.</li> <li>• + Any washable (food and non-food contact) surface where disinfection is required.</li> </ul>		

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**Table 7 Surfaces continued:**

Exterior Surfaces of Freezers+	IV Poles	Outdoor Furniture+ (excluding wood frames and upholstery)
Garbage Cans	IV Stands	Ovens+π
Glass	Jacuzzi	Overbed table
Glazed Ceramic Tile	Kennels	Oxygen Hoods
Glazed Porcelain	Knee Walkers	Patient Chairs
Glazed Tiled Walls	Knobs	Patient Delivery Equipment
Glazed Tiles	Laboratories	Patient Lift Equipment
Golf Cart	Laboratory Equipment	Patient Scales
Grab Bar	Laboratory Surfaces	Patient Support Equipment
Grocery Cart	Laminated Surfaces	Pens
Gurneys	Lamps	Phone Chargers
Gym Equipment	Light Lens Covers	Phone Cradle
Hair Clippers	Lights	Physical Therapy (PT) Equipment surfaces
Hair Cutting Implements	Linoleum	Pipes
Hair Dryers	Manicure Implements	Plastic Laminate
Hair Rollers	Manicure Tables	Plastic Mattress Covers
Hampers	Marlite	Playpens
Handles	Plastic Mattress Covers	Pneumatic Tourniquets
Handrails	Mayo tray stands	Police Cars
Hard Plastic or Metal Flooring	Mayo tray	Porcelain Enamel
Hard, non-porous environmental surfaces	Monitors	Plastic Posey Vest
High Chairs+	Nail/Hair Care Implements	Public Restrooms
High Touch Objects	Non-food contact surfaces+	Razors
Hoods+	Nurse Call [Button] [Device] [Cord]	Reception Areas
Hospital Bed (Railings, Frames)	Operating Room Furniture	Reception Counters
Hospital Patient Bed	Operating Room Lights	Exterior Surfaces of Refrigerator Carts
Infant Incubators (interior and exterior surfaces)	Operating Room Tables	Exterior Surfaces of Refrigerators+
Inside of Drawers	Operating Tables	Safety Straps
Instrument Stands	Operatory Light Switches	Salad Bar Sneeze Guards+
<ul style="list-style-type: none"> <li>• For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [♦ [surfaces] [in ▲]] [and other critical care areas where cross-contamination between treated hard, non-porous surfaces is important] [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.]</li> <li>• [Surfaces which are made of plastic (i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene), and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.]</li> <li>• + A thorough rinse with potable water is required for surfaces in direct contact with food.</li> <li>• + Any washable (food and non-food contact) surface where disinfection is required.</li> <li>• π Allow surfaces to come to room temperature prior to disinfection.</li> </ul>		

**THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL**

Statements in brackets, [ ], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**Table 7 Surfaces continued:**

Salon Surfaces	Stethoscopes	Transport Vehicles
School Buses	Stools	Trash Cans
Scissors	Storage	Tray Tables
Shampoo Bowls	Stoves+π	Troughs
Shears	Stovetops+π	Tubing
Shelves	Stretchers	Vanity Tops+
Shopping Carts	Support Bars	Veterinary Care Surfaces
Shower Chair	Synthetic Marble	Vinyl
Shower Fixtures	Tables+	Walkers
Shower Stalls	Tanning Beds	Wash Basin
Showers	Telemetry Equipment	Washable Nail Files
SilkPeel™ with DermalInfusion™	Telephones	Watering Appliances
Sinks+	Tiles+	Weight-lifting Surfaces
Slit Lamps	Toilet above the water line	Wheelchairs
Smart Phone Chargers	Toilet (Rim, Seat, Tank, Handle,)	Whitening Lamps
Spine Back Boards	Touch Screen Monitors	Workstations
Stainless Steel	Towel Dispensers	Transfer Bench
Stalls	Plastic Toys	Stands
<ul style="list-style-type: none"> <li>• For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [♦ [surfaces] [in ▲]] [and other critical care areas where cross-contamination between treated hard, non-porous surfaces is important] [Use on hard non-porous surfaces and equipment listed on this label made of stainless steel, plastic, Formica, and glass.]</li> <li>• [Surfaces which are made of plastic (i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene), and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.]</li> <li>• + A thorough rinse with potable water is required for surfaces in direct contact with food.</li> <li>• + Any washable (food and non-food contact) surface where disinfection is required.</li> <li>• π Allow surfaces to come to room temperature prior to disinfection.</li> </ul>		

**[OPTIONAL GRAPHICS]**

**[The following graphics are optional for all containers. They will be affixed to the label/ container such that they do not interfere with any label text. Color of graphics may be changed. The following graphics will ONLY be included on Product Labels when Adenovirus II\* is NOT listed on Product Labels.]**



Disinfection



Disinfection



Disinfection