



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
WASHINGTON, DC 20460

OFFICE OF CHEMICAL SAFETY
AND POLLUTION PREVENTION

March 12, 2020

Rebecca Lei
Senior Regulatory Affairs Specialist
Metrex Research
1717 West Collins Avenue
Orange, CA 92867

Subject: PRIA Label and CSF Amendment – Add Two Organisms, First Aid Statement Revise
Poliovirus Contact Time, and Basic and Alternate CSF #2
Product Name: CaviWipes Bleach
EPA Registration Number: 46781-14
Application Date: October 18, 2019
Decision Number: 556620

Dear Ms. Lei:

The amended label and CSF(s) referred to above, submitted in connection with registration under the Federal Insecticide, Fungicide and Rodenticide Act, as amended, are acceptable. This approval does not affect any conditions that were previously imposed on this registration. You continue to be subject to existing conditions on your registration and any deadlines connected with them.

A stamped copy of your labeling is enclosed for your records. This labeling supersedes all previously accepted labeling. You must submit one copy of the final printed labeling before you release the product for shipment with the new labeling. In accordance with 40 CFR 152.130(c), you may distribute or sell this product under the previously approved labeling for 18 months from the date of this letter. After 18 months, you may only distribute or sell this product if it bears this new revised labeling or subsequently approved labeling. "To distribute or sell" is defined under FIFRA section 2(gg) and its implementing regulation at 40 CFR 152.3.

Please note that the record for this product currently contains the following CSF(s):

- Basic CSF dated 10/18/2019
- Alternate CSF 1 dated 05/16/2016
- Alternate CSF 2 dated 10/18/2019

Should you wish to add/retain a reference to the company's website on your label, then please be aware that the website becomes labeling under the Federal Insecticide Fungicide and Rodenticide Act and is subject to review by the Agency. If the website is false or misleading, the product would be misbranded and unlawful to sell or distribute under FIFRA section 12(a)(1)(E). 40 CFR 156.10(a)(5) list examples of statements EPA may consider false or misleading. In addition, regardless of whether a website is referenced on your product's label, claims made on the website may not substantially differ from those claims approved through the registration process. Therefore, should the Agency find or if it is brought to our attention that a website contains false or misleading statements or claims substantially differing from the EPA approved registration, the website will be referred to the EPA's Office of Enforcement and Compliance.

Page 2 of 2
EPA Reg. No. 46781-14
Decision No. 556620

Your release for shipment of the product constitutes acceptance of these conditions. If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA section 6. If you have any questions, please contact Wanda Henson by phone at (703) 308-6345 or via email at henson.wanda@epa.gov

Sincerely,

A handwritten signature in black ink, appearing to read "D. Fuller", with a long horizontal line extending to the right.

Demson Fuller, Product Manager 32
Regulatory Management Branch I
Antimicrobials Division (7510P)
Office of Pesticide Programs

Enclosure

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

CaviWipes Bleach

[TUBERCULOCIDAL] [BACTERICIDAL] [VIRUCIDAL*] [FUNGICIDAL]

<p>Active Ingredient: Sodium Hypochlorite..... 0.91% Other Ingredients†..... 99.09% Total..... 100.00%</p> <p>[†Does not include the weight of the towelette]</p> <p>[Available [in] Chlorine.0.87%] [9100 PPM [Parts Per Million] Available Sodium Hypochlorite] [8700 PPM [Parts Per Million] Available Chlorine] [Pre-moistened wipe that meets [CDC recommended] 1:10 bleach dilution]</p>	<p>KEEP OUT OF REACH OF CHILDREN</p> <div data-bbox="982 724 1377 1012" style="border: 2px solid black; padding: 10px; text-align: center;"><p>A C C E P T E D</p><p>03/12/2020</p><p>Under the Federal Insecticide, Fungicide and Rodenticide Act as amended, for the pesticide registered under EPA Reg. No. 46781-14</p></div>
<p>EPA Reg. No. 46781-14 EPA Est. No. 8251-WI-003 [(A1)], 8251-WI-004 [(A2)], 8251-WI-005 [(A3)], 8251-OH-001 [(A4)], 46781-MI-00 [(A5)], 56952-WI-001 [A6]), 2686-WI-008 [(A7)]</p> <p>[Beginning of batch code indicates Est No.] [End of batch code indicates Est No.]</p> <p>Alternate Brand Name (ABN): CaviBleach Wipes</p> <p>Made in USA</p>	<p>Reorder No. Net Contents: #[[lbs.] [oz] # Pre-Saturated Towelette # x # in</p> <p>[Packed For:] [Distributed By:] [Sold By:] [Manufacturer By [or for]:]</p> <p>Metrex Research 28210 Wick Road Romulus, MI 48174</p> <p>[To reorder] [,] [For product or technical information] contact Metrex at 800-841-1428 or visit our web site at www.metrex.com [For SDS information please visit www.metrex.com]</p>

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PHYSICAL AND CHEMICAL HAZARDS: This product contains bleach. Do not use this product with other chemicals such as ammonia, toilet bowl cleaners, rust removers, or acid, as this releases hazardous gases.

[FIRST AID]

[If on skin or clothing: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes.]

[For Chemical Emergency ONLY (spill, leak, fire, exposure or accident), call CHEMTREC at 1-800-424-9300].

For [ALL] [other] inquiries about this product, call Metrex Customer [Service] [Care] at 1-800-841-1428, Monday through Friday, between 6 a.m. and 4 p.m. Pacific Time.

Note: This product meets EPA's Toxicity Category IV by all routes of exposure. Therefore a Signal Word, Precautionary Language and First Aid Statement is not required in accordance to EPA Label Review Manual: Chapter 7. Registrant voluntarily put in the above First Aid statements.

Packaging Statements and Wipe Dispensing Instructions:

[Use A, B, or C as applicable depending on the wipe container/dispensertype.]

[A – Refillable, Non-refillable canister, buckets, tubs or pails]

[To dispense towelettes] [container] [canister][pail] [bucket] [tubs] [Dispensing Directions]:

[To start feed, remove cover and discard seal. From center of towelette roll, pull up a towelette corner, twist it into a point and thread it through the hole in the container cover. Pull through about one inch. Replace cover. Pull out first towelette and tear off at [(any angle)] [45°-90° angle]. Remaining towelettes feed automatically, ready for next use. When not in use, keep center cap closed to prevent solution loss.]

[To refill] [canister] [bucket] [pail]: Place refill [pouch] [bag] in bucket [canister]. Remove the lid from the [pail] [bucket] [canister]. [Use scissors] [to] [Cut] [Rip] [Tear] [Pull] [Open] refill pouch [bag]. Place [put] [insert] refill pouch into [bucket] [pail] [canister]. DO NOT remove wipe roll from refill pouch [bag]. From center of towelette roll, pull up a towelette corner, twist it into a point and thread it through the hole in the [bucket] [pail] [container] cover. Pull through about one inch. Replace cover. Pull out first towelette and tear off at [(any angle)] [45°-90° angle]. Remaining towelettes feed automatically, ready for next use. When not in use, keep center cap closed to prevent solution loss.]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Packaging Statements and Wipe Dispensing Instructions *continued*:
[Use A, B, or C as applicable depending on the wipe container/dispenser type.]

[B – Flat Pack]

[To dispense towelettes] [from flat pack:] [from Flat Pack wipes container:]

[Grasp the white tab that states "PULL HERE". Gently swing the lid open and grasp a towelette. Remove towelette from package and then swing lid closed. Gently tap the lid to ensure seal is made between it and the towelette package.]

[C – Single use packet]:

[Tear open packet. Unfold towelette and use as directed. Discard as directed.]

This product qualifies for emerging pathogens claims against enveloped, large non-enveloped, and small non-enveloped when used according to the direction for use against Feline Calicivirus.

[Product name] has demonstrated effectiveness against viruses similar to **[name of emerging virus, which must be a small non-enveloped, large non-enveloped, or enveloped virus]** on **[eligible surface type(s)]**. Therefore, **[product name]** can be used against **[name of emerging virus, which must be a small non-enveloped, large non-enveloped, and/or enveloped virus]** when used in accordance with the directions for use against **[name of supporting virus(es)]** on **[hard, non-porous surfaces]**. Refer to the **[CDC or OIE]** website at **[pathogen-specific website address]** for additional information.

[Name of illness/outbreak] is caused by **[name of emerging virus, which must be a small non-enveloped, large non-enveloped, and/or enveloped virus]**. **[Product name]** kills similar viruses and therefore can be used against **[name of emerging virus, which must be a small non-enveloped, large non-enveloped, and/or enveloped virus]** when used in accordance with the directions for use against **[name of supporting virus(es)]** on **[hard, non-porous surfaces]**. Refer to the **[CDC or OIE]** website at **[website address]** for additional information.

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.

[Use A, B, C, or D as applicable depending on the wipe container/dispenser type.]

[A – Nonrefillable container, canister, buckets or pails]

STORAGE AND DISPOSAL:

Product Storage: [Store at room temperature.] -or- [Store this product in a cool [dry] area, [away from direct sunlight and heat] to avoid deterioration.] [When not in use, keep center cap of lid closed to prevent moisture loss.]

Towelette Disposal: Do not re-use towelette. Dispose of used towelette in trash. [Do not flush in toilet.]

Container Handling: Nonrefillable container. Do not reuse or refill this container. Offer for recycling, if available. If recycling is not available, discard in trash

[B – Refillable container canister, buckets, tubs, or pails]

STORAGE AND DISPOSAL:

Product Storage: [Store at room temperature.] -or- [Store this product in a cool [dry] area, [away from direct sunlight and heat] to avoid deterioration.] [When not in use, keep center cap of lid closed to prevent moisture loss.]

Towelette Disposal: Do not re-use towelette. Dispose of used towelette in trash. [Do not flush in toilet.]

[Refillable] Container Handling: Refillable Container. Do not reuse or refill except as described in direction for use. Refill only with this product. Triple rinse cracked or broken containers and offer for recycling. If recycling is not available, discard in trash.

[C – Flat pack]

STORAGE AND DISPOSAL

Product Storage: [Store at room temperature.] -or- [Store this product in a cool [dry] area, [away from direct sunlight and heat] to avoid deterioration.] [When not in use, keep [flap] [tab] closed to prevent moisture loss.]

Towelette Disposal: Do not reuse towelette. Dispose of used towelette and empty packet in trash. [Do not flush in toilet.]

[Flat Pack] [or Dispenser] Handling: Non-refillable container. Do not reuse or refill this flat pack. If recycling is not available, discard in trash.

[D – Packet]

STORAGE AND DISPOSAL

Product Storage: [Store at room temperature.] -or- [Store this product in a cool [dry] area, [away from direct sunlight and heat] to avoid deterioration.] [When not in use, keep [flap] [tab] closed to prevent moisture loss.]

Towelette Disposal: Do not reuse towelette. Dispose of used towelette and empty packet in trash. [Do not flush in toilet.]

[Packet][or Dispenser] Handling: Non-refillable container. Do not reuse or refill this packet. Discard packet in trash.

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

To [clean and] disinfect [and deodorize] hard, nonporous surfaces:

Wipe hard, nonporous surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for 3 minutes. [Rinse.] Let air dry. Visible soil loads must be removed prior to disinfecting. [If streaking is observed, wipe with a clean, damp [cloth or] paper towel after 3 minute contact time has expired.]

(#)(Y)Special Instructions for Cleaning Prior to Disinfection against *Clostridium difficile* spores.

Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks and eye covering.

Cleaning Procedure: Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with [Product Name] [this product saturated with product name]. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.

Infectious Materials Disposal: Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

(#)(Y)Special Label Instructions for Cleaning Prior to Disinfection against *Candida auris*.

Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks, or eye covering.

Cleaning Procedure: Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with [Product Name] [this product saturated with product name]. Pre-cleaning is to include vigorous wiping and/or scrubbing and all visible soil is removed. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading the organism. Restrooms are to be cleaned last. Do not reuse soiled cloths.

Infectious Waste Disposal: Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

The following paragraph will be included in its entirety when the production label references medical devices:

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream, or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

DIRECTIONS FOR USE *continued*

The following HIV-1 / HBV / HCV instructions are to be included on production labels when claims against HIV-1 or HBV or HCV are listed.

(Y) CaviWipes Bleach effectively kills [Human Immunodeficiency Virus* Type 1 (HIV-1),] [Hepatitis B Virus* (HBV),] [[and] Hepatitis C Virus* (HCV)] on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects can be associated with the potential for transmission of [HIV-1,] [HBV,] [[and] [HCV].

Special Instructions for Cleaning and Decontamination Against [HIV-1,] [HBV,] [[and] HCV] of Surfaces or Objects Soiled With Blood or Body Fluids:

Personal Protection: When handling items soiled with blood or body fluids, wear appropriate barrier protection: disposable protective gloves, gowns, face masks and eye coverings.

Cleaning Procedure: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of this disinfecting towelette. A CaviWipes Bleach towelette can be used for this purpose.

Contact Time: Thoroughly wet surface with a new CaviWipes Bleach towelette. Repeated use of the product may be required to ensure that the surface remains visibly wet for 3 minutes.

Disposal of Infectious Materials: Dispose of used wipe in accordance with Federal, State or local regulations for infectious materials disposal.

The following direction is to be included when hard, non-porous exterior surfaces of ultrasound transducers or probes are listed as a use site.

To disinfect external surfaces of ultrasound transducers: Pre-clean ultrasound gel from external transducer to be disinfected [by using a clean washcloth]. Then thoroughly wet the probe with a CaviWipes Bleach towelette. Repeated use of the product may be required to ensure that the surface remains visibly wet for 3 minutes.

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

DIRECTIONS FOR USE *continued*

The following statement is to be included on production labels when animal housing facilities are listed. (Including Veterinary Clinics, Farms and Poultry Houses)]

For cleaning and disinfecting the following hard nonporous surfaces: feeding and watering equipment, utensils, instruments, cages, kennels, stables, catteries, etc.

1. Remove all animals and feeds from premises, animal transportation vehicles, crates, cages etc.
2. Remove all litter, droppings and manure from floors, walls and surfaces of facilities (including pens, chutes, and barns) occupied or traversed by animals.
3. Empty all troughs, racks, and other feeding and watering appliances.
4. Thoroughly clean all surfaces with soap or detergent and rinse with potable water. Use enough wipes for treated surface to remain visibly wet for 3 minutes.
5. Wipe all handling and restraining equipment such as leashes, muzzles, halters or ropes, as well as forks, shovels, and scrapers used for removing litter and manure.
6. Ventilate buildings and other closed spaces.
7. Do not house animals or employ equipment until treatment has been absorbed, set or dried.
8. Thoroughly scrub treated feed racks, troughs, mangers, automatic feeders, fountains, and waterers with soap or detergent, and rinse with potable water before reuse.

The following statement is to be included on production labels when food contact areas or surfaces are listed.

+ A thorough rinse with potable water is required for surfaces in direct contact with food.

[DO NOT USE AS A BABY WIPE.] [Not For Use On Skin] [DO NOT FLUSH WIPE]

[Not for cleaning or sanitizing skin.] [Do not use as a diaper wipe or for personal cleansing.] [APPLY CAREFULLY – This product contains bleach.] [Refill only with this product.] [Always refer to manufacturer's care instructions before using on equipment -or- devices.]

[The following graphics are optional for canister containers only. They will be affixed to the label/canister so that it does not interfere with any label text. Due to size constraints, they will not appear on individually wrapped packets.]



THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

DISINFECTION

[CaviWipes Bleach Disinfecting Towelettes are non-woven disposable towelettes pre-saturated with [CaviCide Bleach] [CaviWipes Bleach Solution] [CaviWipes Bleach [Solution].]

[CaviWipes Bleach Disinfecting Towelettes are multi-purpose [disinfectant]/ [cleaning] wipes for use on hard, non-porous inanimate surfaces.]

[CaviWipes Bleach are effective against the following microorganisms on hard, non-porous surfaces when used as directed:]

[This product has demonstrated effectiveness against the following organisms with a 3 minute contact time:]

[Disinfects (*insert organisms*) on hard non-porous environmental surfaces listed on this label]

[When used as directed, CaviWipes Bleach provides effective Tuberculocidal / Virucidal* / Bactericidal / Fungicidal activity against the following organisms:]

Pathogen:	Bacteria	Viruses*	Bloodborne Pathogens**	TB	Parvo viruses*	<i>Clostridium difficile</i> ‡	Mold/Mildew /Fungi	<i>Candida auris</i> #	Poliovirus*
Contact Time:	3 min[utes]	3 min[utes]	3 min[utes]	3 min[utes]	3 min[utes]	3 min[utes]	3 min[utes]	3 min[utes]	1 min[ute]

ORGANISMS

	[Strain and/or ATCC number]
Bacteria (3 Minute Contact)	
<i>Clostridium difficile</i> spores ‡	[ATCC 43598]
<i>Staphylococcus aureus</i>	[ATCC 6538]
<i>Pseudomonas aeruginosa</i>	[ATCC 15442]
<i>Salmonella enterica</i>	[ATCC 10708]
<i>Burkholderia cepacia</i>	[ATCC 25416]
<i>Escherichia coli</i> O157:H7	[ATCC 35150]
<i>Klebsiella pneumoniae</i>	[ATCC 4352]
<i>Legionella pneumophila</i>	[ATCC 33153]
<i>Listeria monocytogenes</i>	[ATCC 19117]
<i>Serratia marcescens</i>	[ATCC 14756]
<i>Streptococcus pyogenes</i>	[ATCC 12344]
<i>Neisseria gonorrhoeae</i>	[ATCC 43069]
Drug-Resistant Bacteria (3 Minute Contact)	
<i>Enterobacter cloacae</i>	[ATCC BAA-2468]
Multidrug Resistant <i>Acinetobacter baumannii</i> (MRAB) (Resistant to Ceftazidime and Gentamicin)	[ATCC BAA-1605]
Carbapenem-resistant <i>Klebsiella pneumoniae</i> (CRKP) (CRE)	[ATCC BAA-1705]
Methicillin Resistant <i>Staphylococcus aureus</i> (MRSA)	[ATCC 33592]
Penicillin Resistant <i>Streptococcus pneumoniae</i>	[ATCC 700677]
Vancomycin Resistant <i>Staphylococcus aureus</i> (VISA)	[NARSA VRS1]
Vancomycin Resistant <i>Enterococcus faecalis</i> (VRE)	[ATCC 51575]
Extended spectrum β -lactamase (ESBL) <i>Escherichia coli</i> , (Resistant to Ceftazidime and Penicillin)	[ATCC BAA-196]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

ORGANISMS <i>continued</i>	
	[Strain and/or ATCC number]
<i>Mycobacterium</i> (3 Minute Contact)	
<i>Mycobacterium tuberculosis</i> var: <i>bovis</i> (BCG) (TB)	
Mold, Mildew, Fungi (3 Minute Contact)	
<i>Candida albicans</i>	[ATCC 10231]
<i>Candida auris</i> #	[AR-BANK#0381] [from CDC]
<i>Aspergillus brasiliensis</i>	[ATCC 16404]
<i>Trichophyton interdigitale</i> (formerly <i>Trichophyton mentagrophytes</i>)	[ATCC 9533]
Viruses * [Non-enveloped] [:] (3 Minute Contact)	
Adenovirus [type 2] [(Adenoid 6)] *	[ATCC VR-846]
Hepatitis A [Virus] [Human [(HAV)] *	[Strain HM-175]
Feline Calicivirus Surrogate for Norovirus [Feline Calicivirus]*	[ATCC VR-782]
Rhinovirus [type 37] *	[ATCC VR-1147] [(Strain 151-1)]
Enterovirus [type D68] (EV-D68) *	[ATCC VR-561]
Human Rotavirus [strain WA] *	[ATCC VR-2018]
Poliovirus [type 1], Chat Strain * [Non-enveloped] [(3 Minute Contact)] [(1 Minute Contact)]	[ATCC VR-1562]
Viruses * [Enveloped] [:] (3 Minute Contact)	
Herpes simplex virus [type 1] *	[ATCC VR-260] [Strain HF]
Herpes simplex virus [type 2] *	[ATCC VR-734] [Strain G]
Human Coronavirus *	[ATCC VR-740] [Strain 229E]
Influenza A virus Strain A/Hong Kong 8/68 [(H3N2)] *	[Charles River Labs]
Influenza B virus Strain B/Hong Kong/5/72 *	[ATCC VR-823]
Duck Hepatitis B virus (HBV) [surrogate for human Hepatitis B virus] *, **	[Strain Grimaud]
Bovine-Viral Diarrhea Virus Hepatitis C virus (HCV) [BVDV surrogate for Hepatitis C virus] *, **	[Strain NADL] [ATCC VR-1422]
HIV-1 [(Human Immunodeficiency Virus type 1)] *, **	[Strain III B (B)] [ZeptoMetrix]
Respiratory syncytial virus [(RSV)] *	[ATCC VR-26] [(Strain Long)]
Measles virus *	[Strain: Edmonston, ATCC VR-24]
Parvoviruses * [Non-enveloped] (3 Minute Contact)	
Canine Parvovirus *	[ATCC VR-2017] [(Cornell Strain)]

¥ *Clostridium difficile* spores

Follow [the] Special Instructions for Cleaning required prior to Disinfection [listed on label]

*Candida auris*

Follow [the] Special Instructions for Cleaning required prior to Disinfection [listed on label]

*** Virus**

**** Bloodborne pathogens**

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS

▲ to represent any use site(s) listed in below as applicable.

◆ to represent any surface(s) listed in below as applicable.

- [___# pre-saturated towelettes] [X" x Y"]
- [3 minute contact time]
- [[99.9999% [6 [six] log] reduction for [*Clostridium difficile*¥] -or- [*C. difficile*¥] -or- [*C. diff*¥] [(*C. diff*)¥] [spores]
- [[99.999% [5 [five] log] reduction for [*Candida auris*#] -or- [*C. auris*#]
- [Kills a minimum of [[99.999% [5 [five] log] of [*Candida auris*#] -or- [*C. auris*#] on hard, non-porous surfaces.]
- [A 1:10 sodium hypochlorite solution similar to that recommended by the CDC for effective disinfecting against *Clostridium difficile*¥ -or- *C. difficile*¥ -or- *C. diff*¥ [(*C. diff*)¥] [spores]]
- [A 1:10 sodium hypochlorite solution recommended by the CDC for effective environmental disinfection against [*Candida auris*#] -or- [*C. auris*#]
- [A Metrex Research® brand]
- [A Metrex® brand]
- [A CaviWipes™ Brand]
- [Absorbs] [Absorbent] [wipe]
- [A disinfectant cloth for use in ▲ on ◆]
- [[Aids] [Helps] in reducing -or- in the reduction of cross-contamination between treated surfaces]
- [An effective method to clean against odor-causing organisms on hard, non-porous surfaces listed on this label].
- [[A] one-step ¥ cleaner/disinfectant]
- [A disinfectant towelette for use on ◆ and other ◆]
- [After use, just toss wipe(s) away]
- [Always ready-to-use][
- [An effective method to deodorize and disinfect against odor-causing organisms on hard non-porous surfaces, when used [according to] [following] disinfection [directions] [instructions]]
- [An effective way to disinfect hard, non-porous environmental surfaces]
- [An effective way to disinfect hard, non-porous environmental surfaces listed on this label.]
- [A non-woven disposable towelettes for use in ▲ and other critical care areas where control between treated surfaces is required]
- [Antibacterial]
- [Antimicrobial]
- [Antifungal]
- [A One Step ¥ Cleaner [and disinfectant] [except when precleaning is required for [*C. diff*¥,] [*C. auris*#,] blood and body fluids]]
- [A solution impregnated in a wiping cloth]
- [[Appropriate] [for] Cleaning Patient Discharge]
- [Appropriate] for [Occupied Room] [and] [Terminal Room] [Discharge] Cleaning]
- [Assures proper dilution strength, product effectiveness and standardizes technique]
- [Attacks and kills [*Clostridium difficile*¥] [*C. difficile*¥] [*C. diff*¥] [(*C. diff*)¥] [spores] [on precleaned, hard, nonporous surfaces]
- [Attacks and/or Kills and/or Inactivates [*Candida auris*#] [*C. auris*#] [on precleaned, hard, nonporous surfaces]
- [Bacteria -and/or- virus* -and/or- germ fighting [formula]]
- [Bactericidal] [Virucidal*] [Fungicidal] [Sporicidal] [Tuberculocidal] [in 3 minutes]
- [Bactericide -or- Bactericidal] [(in 3 minutes)]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Bathroom disinfectant wipes]
- [Big [textured] wipes] [3 minute kill time]
- [Broad spectrum disinfectant]
- [Built-in [tough] scrubbing fibers to help clean effectively]
- [CaviWipes Bleach Wipes have been tested for surface safety on common ▲ surfaces [and equipment] to ensure compatibility with [hospital-grade] surfaces including: [Stainless Steel] [Porcelain] [Glass] [Plastic] [Glazed Ceramic]
- [CaviWipes Bleach kills germs that cause odors]
- [CaviWipes Bleach is effective against odor causing germs]
- [CaviWipes Bleach will effectively clean and disinfect when used as directed]
- [Can help reduce cross contamination on treated surfaces]
- [Cleaner and disinfectant [towelette] [wipe]]
- [Cleaning Formula]
- [Cleans] [◆] [in ▲]]
- [[Cleans and] disinfects ◆, when used [according to] [following] disinfection [directions] [instructions]]
- [[Cleans and] disinfects [hard, nonporous surfaces] [on the go], when used according to disinfection [directions] [instructions]]
- [Cleans, deodorizes and disinfects, when used [according to] [following] disinfection [directions] [instructions]]
- [Cleans] [disinfects] [Cleaner and disinfectant] [in one] [towelette] [wipe] [all in one [easy] step], when used [according to] [following] disinfection [directions] [instructions]]
- [Clean Home]
- [Cleans the toughest soils]
- [Cleaning Wipe[s] [Towel][s] refill[s] [able]]
- [Cleaning Made Easy]
- [Cleaner] [on the go]
- [Clinical Clean]
- [[Cleans [and /or Disinfects] [grocery cart handles or shopping cart handles or grocery cart child seats or shopping] cart child seats, when used [according to] [following] disinfection [directions] [instructions]]
- [Cleans grocery carts and shopping carts]
- [[Cloth-like texture] [Textured] [strong] [for] [with] scrubbing power]
- [Concentration of hypochlorite in CaviWipes Bleach meets the disinfection requirements of a 1:10 dilution of bleach]
- [Contains [number] canisters [plus][and][number refills]
- [Contains no abrasive harsh acids]
- [Contains 0.91% sodium hypochlorite]
- [Contains (name of active ingredient), the [same] active ingredient used in [Brand-name] This product is not manufactured or distributed by ® /or TM/]
- [Contains a (1:10 sodium hypochlorite solution similar to the CDC recommendation for effective disinfecting against [*Clostridium difficile* ¥] [*C. difficile*¥] [*C. diff*¥] spores)].
- [Contains a (1:10 sodium hypochlorite solution similar to the CDC recommendation for effective disinfecting against [*Candida auris*#] [*C. auris*#]).
- [Contains a (1:10 sodium hypochlorite solution similar to that recommended for effective disinfecting against [*Clostridium difficile* ¥] [*C. difficile*¥] [*C. diff*¥] spores) following CDC recommendation].

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Contains a (1:10 sodium hypochlorite solution similar to that recommended for effective disinfecting against [*Candida auris*#] [*C. auris*#] following CDC recommendation].
- [Contains a (1:10 sodium hypochlorite solution similar to that recommended for effective disinfecting hard surfaces against [*Clostridium difficile*¥] [*C. difficile*¥] [*C. diff*¥] spores)]
- [Contains a (1:10 sodium hypochlorite solution similar to that recommended for effective disinfecting hard surfaces against [*Candida auris*#] [*C. auris*#]]
- [Convenient]
- [[Convenient] [Portable] way to [disinfect] [clean] hard non-porous surfaces [of ♦], when used [according to] [following] disinfection [directions] [instructions]]
- [Cuts Cleaning time]
- [Deep cleaning] [formula]]
- [Delivers a clean you can see and smell]
- [Demonstrates [99.9999%] [6 [six] log] reduction against for [*Clostridium difficile*¥] -or- [*C. difficile*¥] -or- [*C. diff*¥] [(*C. diff*)¥] [spores]
- [Demonstrates [99.999%] [5 [five] log] reduction against for [*Candida auris*#] -or- [*C. auris*#]
- [Designed for Home Healthcare -and/or- patients -and/or caregivers]
- [Designed -or- developed for Home Caregivers]
- [Designed -or- developed for use in [▲] Home Healthcare -or- Home Care]
- [Demonstrates [99.9999%] [6 [six] log] 6 [six] log reduction against *Clostridium difficile*¥ [*C. difficile*¥] [*C. diff*¥] spores [in [3] [three] minutes]]
- [Demonstrates [99.999%] [5 [five] log] 5 [five] log reduction against *Candida auris*# [*C. auris*#] [[in [3] [three] minutes]]
- [Delivers consistent [amount of solution] concentration every time [!]]
- [Deodorizes, and kills common household germs]
- [Deodorizes by killing the germs that cause odors]
- [Deodorize hard, non-porous surfaces] [♦] [in ▲]]
- [Designed for the ▲ [[industry]]]
- [Designed for killing *Clostridium difficile*¥ -or- *C. difficile*¥ -or- *C. diff*¥ [(*C. diff*)¥] [spores] on surfaces in ▲]
- [Designed for killing *Candida auris*# -or- *C. auris*# on surfaces in ▲]
- [Designed for killing [*Clostridium difficile*¥] [*C. difficile*¥] [*C. diff*¥] [(*C. diff*)¥] [spores] [*] [on precleaned, hard, nonporous surfaces] in ▲.]
- [Designed for killing [*Candida auris*#] [*C. auris*#] [*] [on precleaned, hard, nonporous surfaces] in ▲.]
- [Designed] [for] [with] [the] Environmental Services Professional[s] [in] [mind][!]]
- [Destroys [Deactivates] [[Helps] Eliminate [99.9 –] 99.9999%] [tough-to-kill] *Clostridium difficile*¥ -or- *C. difficile*¥-or- *C. diff*¥ [(*C. diff*)¥] [spores] [on ▲]]
- [Destroys [Deactivates] [[Helps] Eliminate [99.9 –] 99.999%] [tough-to-kill] *Candida auris*# -or- *C. auris*# [on ▲]]
- [Disinfectant]
- [Disinfect]
- [Disinfect external surfaces of ♦]
- [Disinfect surfaces [♦]]
- [Disposable]
- [Disinfectant cleaner [that cleans, disinfects and deodorizes] in one step¥, when used [according to] [following] disinfection [directions] [instructions]]
- [Disinfecting formula]
- [Disinfects] [♦] [in ▲]]
- [[Disinfects] [Disinfectant] on the go]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Disinfects [and deodorizes] ♦ [in one step], when used [according to] [following] disinfection [directions] [instructions]]
- [Disinfecting towelettes]
- [Disinfects and deodorizes in 3 minutes, when used [according to] [following] disinfection [directions] [instructions]]
- [Disinfects hard, non-porous surfaces in 3 minutes]
- [[Disinfects] hard surfaces [that people touch]]
- [Disinfects [in [just] three [3] minutes]]
- [[Disinfects and deodorizes by killing germs and their odors, when used [according to] [following] disinfection [directions] [instructions]]
- [Disinfects hard, nonporous surfaces [throughout the ▲]]
- [Disinfects [as it cleans] [and scrubs in one easy step], when used [according to] [following] disinfection [directions] [instructions]]
- [Disinfect [the] precleaned surface]
- [Dispenses one wipe at a time]
- [Does not contain VOCs]
- [Does not contain dyes]
- [DO NOT FLUSH WIPE DOWN TOILET]
- [Double sided-scrubbing wipes]
- [Durable, lint-free wipe designed for sensitive equipment surfaces]
- [Easy]
- [Easy and convenient to use]
- [[Easy] [and convenient] way to clean [and] [disinfect], when used [according to] [following] disinfection [directions] [instructions]]
- [Easy One Step Cleaning [and disinfecting] (except when precleaning is required for visible soil and blood and body fluids)]
- [[Easy] pop-up-[tub] [canister] [pail] [lid] [tab] [dispenser] [bucket]]
- [Easy [and convenient] way to disinfect [bathrooms] [office] [nursing stations] [high touch] hard, non-porous surfaces everyday]
- [Easy] [convenient] [dispensing] lid [closure][cap][package]
- [Easy drop-in refills]
- [Easy to use]
- [Easy to handle]
- [Easy to refill]
- [Easy to use package [canister] [container] [refill] keeps [your] wipes moist]
- [Easy way to [clean] [and] [disinfect] [bathroom] [or exterior toilet surfaces] [hard, non-porous surfaces], when used [according to] [following] disinfection [directions] [instructions]]
- [Easy and convenient disinfecting [throughout the ▲]]
- [Easy, one-step cleaning and disinfecting]
- [Easy [and convenient] way to [touch-up] [bathroom] [office] everyday]
- [Easy start [tab][lid] [(insert package type)]]
- [[Easy] [and convenient] [way to clean] [and disinfect], when used [according to] [following] disinfection [directions] [instructions]]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Easy way to clean [and disinfect] [bathroom] [or exterior toilet surfaces] [hard, non-porous surfaces], when used [according to] [following] disinfection [directions] [instructions]]
- [Economical]
- [Economy]
- [Economical] Refill[s] [one] [1] [two] [2] [three] [3] [four] [4] [various numbers] pack]]
- [Economical] Refills [now] available]
- [Effective against ***[Insert any organism from list of organisms]*** [and] *Clostridium difficile* -or- *C. difficile* -or- *C. diff* [(*C. diff*)] [spores] [and] *Candida auris* -or- *C. auris*]
- [Effective against Avian Influenza A* and H1N1* [Pandemic] [Influenza A Virus]
- [Effective against Measles* -or- Measles Virus*]
- [Effective against Poliovirus* [in [just] [as short as] [1] [one] minute]]
- [Effective against bathroom bacteria and viruses* [use at least two bacteria and two viruses* from this list:] [ESBL *E. coli*] [*Salmonella enterica*] [*Staphylococcus aureus*] [Herpes Simplex Virus Type 2] [Hepatitis B Virus (HBV)*] [Hepatitis C Virus (HCV)*] [Adenovirus II] [Influenza A, H3N2 Virus*] [Measles Virus*]
- [Effective against -or- Kills a range of bacteria including *Staphylococcus aureus*, *Salmonella enterica*, *Pseudomonas aeruginosa* [in the presence of organic soil [(5% blood serum)]]
- [Effectively disinfects against MRSA]
- [[Effective against -or- Kills (***insert any organism(s) from Organisms table***) [in the presence of organic soil [(5% blood serum)]]
- [Effective in the presence of 5% blood serum contamination]
- [Effective against flu virus (Influenza A/Hong Kong, H3N2 Virus*)]
- [Effective against flu virus Measles* -or- Measles virus*]
- [[Effective -or- For daily use against *Clostridium difficile* -or- *C. difficile* -or- *C. diff* [(*C. diff*)] [spores]]
- [[Effective -or- For daily use against *Candida auris* -or- *Candida auris*]]
- [[Effectively] disinfects [hard, nonporous] [environmental] surfaces]
- [Effectively removes everyday soils]
- [Efficacy tests have demonstrated that [Product Name] [this product] [are] [is] an effective bactericide, virucide*, tuberculocide, and fungicide in the presence of organic soil (5% blood serum)]
- [Eliminates odors]
- [Eliminates odors from fecal matter and urine]
- [Eliminates airborne particles associated with sprays]
- [Eliminates [daily] mixing and measuring]
- [Eliminates [lingering] [food+] odors [on surfaces]]
- [[Eliminates odors at their source [caused by germs -or- bacteria-and/or- mildew]]
- [[Eliminates -or- Reduces odors [caused by germs -or- bacteria -and/or- mildew]]
- [EVS [Environmental Services] [use]]
- [Exhibits effective virucidal* activity]
- [Extra large [disinfecting] [Y" x Z" compared to our A" x B" size] wipes]
- [Extra large] [Disinfecting Towelettes]
- [Convenient] [portable] way to disinfect hard non-porous surfaces]
- [Fight[s] -and/or- kill[s] -and/or- Effective against 99.9% [of]] *Escherichia coli* [(*E. coli*)]]
- [Fight[s] -and/or- kill[s] -and/or- Effective against 99.9% [of]] *Salmonella enterica* [(*Salmonella*)]]
- [Fight[s] -and/or- kill[s] -and/or- Effective against Avian flu [virus*] [on environmental -and/or- hard, nonporous surfaces]]
- [Fight[s] -and/or- kill[s] -and/or- Effective against cold -and/or- flu [virus[es*]]
- [Fight[s] -and/or- kill[s] -and/or- Effective against Measles* - or- Measles virus*]
- [Fight[s] -and/or- kill[s] -and/or- Effective against Poliovirus*] [in [just] [as short as] [1] [one] minute]]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Fight[s] -and/or- Kill[s] -and/or- Effective against (*insert any organism from EPA approved list*)]
- [Fight[s] -and/or- reduce[s] the spread of [99.9% of] [illness-causing] bacteria-and/or- viruses* -and/or- germs on hard, nonporous surfaces [in your ▲] [this cold -and/or- flu season]]
- [For daily use]
- [For] general use
- [For more information visit us at [www.\(insert website\).com](http://www.(insert website).com)]
- [For use on bathroom surfaces]
- [For use on nursery surfaces]
- [For use on ♦ surfaces]
- [For use in critical care areas [where control of [the hazards of] cross contamination from hard non-porous surfaces is of prime importance]
- [For] Surface cleaning
- [For Use by Home Caregivers]
- [For Use in Home Healthcare]
- [For use in ▲]
- [Formulated -or- Designed to] remove[s] biological odors including Urine -and/or- Feces -and/or- Vomit]
- [Formulated to kill *Clostridium difficile*¥ -or- *C. difficile*¥ -or- *C. diff*¥ on [(*C. diff*)¥] [spores] [precleaned, hard, nonporous surfaces] [when used according to instructions]]
- [Formulated to kill *Candida auris*# -or- *C. auris*# on [precleaned, hard, nonporous surfaces] [when used according to instructions]]
- [Fragrance Free]
- [Freshens]
- [Free of fragrance]
- [Fresh [lingering] bleach scent]
- [Fulfills CDC's 1:10 dilution of bleach] [1:10 dilution of household]
- [Fungicide]
- [Fungicidal]
- [Fungicidal in 3 minutes]
- [Germicidal disposable wipes]
- [Germicide -or- Germicidal]
- [Gets rid of [or eliminates] odors]
- [Great for everyday use]
- [[Heavy] cloth-like towel]
- [[Helps] [kill][s] [eliminate [s] 99.9%] odor-causing bacteria]
- [Reduces cross contamination on treated surfaces]
- [Includes ___count refill [pack] [contains]___wipes]
- [Institutional [size]]
- [Just] [Simply] wipe [exterior toilet surface] [bathroom germs] [away]
- [Just [pull] wipe and toss]
- [Kills and/or inactivates spores of *Clostridium difficile*¥ on hard, non-porous surfaces]
- [Kills and/or inactivates *Candida auris*# -or- *C. auris*# on hard, non-porous surfaces]
- [Kills Adenovirus II*]
- [Kills Measles* -or- Measles virus*]
- [Kills Poliovirus*] [in [just] [as short as] [1] [one] minute]]
- [Kills Pandemic 2009 H1N1 Influenza A Virus* [(formerly called swine flu)]]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Kills [bacteria] [viruses*] [fungi] *Clostridium difficile* -or- *C. difficile* -or- *C. diff* [(*C. diff*)] [spores] [*Candida auris*# -or- *C. auris*#]
- [Kills bacteria, viruses*, mold, mildew]
- [Kills [bacteria] [viruses*] [fungi] [*C. difficile* spores] [*Candida auris*# -or- *C. auris*#] [on precleaned, hard, nonporous surfaces]
- [Kills] [Is effective against] 99.9999% of [*Clostridium difficile*] [*C. difficile*] spores in 3 [minutes] [min] [on precleaned, hard, nonporous surfaces]]
- [Kills] [Is effective against] 99.999% of [*Candida auris*#] [*C. auris*#] in 3 [minutes] [min] [on precleaned, hard, nonporous surfaces]]
- [Kill germs and bacteria instead of germ paper towels, sponges and rags [Use this product to]]
- [[Kills] [Disinfects] *Clostridium difficile* -or- *C. difficile* -or- *C. diff* [(*C. diff*)] [spores] in 3 min[utes] [on hard, nonporous surfaces]]
- [[Kills] [Disinfects] *Candida auris*# -or- *C. auris*# in 3 min[utes] [on hard, nonporous surfaces]]
- [Kills odor causing bacteria]
- [Kills -or- Destroys mold and mildew]
- [Kills [Disinfects] Norovirus* [Feline Calicivirus] [in] 3 [three] minutes]
- [Kills [Disinfects] Measles* -or- Measles virus* [in] 3 [three] minutes]
- [Kills [Disinfects] Poliovirus* [in] [just] [as short as] 1 [one] minute]
- [Kills -or- Effective against bacteria]
- [Kills -or- Effective against germs]
- [Kills -or- Effective against viruses*]
- [Kills [or eliminates 99.9%] of bacteria]
- [Kills [or eliminates 99.9%] of viruses*]
- [Kills [or eliminates 99.9%] of fungus]
- [[Kills -or- Is effective against [99.9999%] *Clostridium difficile* -or- [(*C. difficile* -or- *C. diff* [(*C. diff*)])] spores] [in 3 min[utes]] [on hard, nonporous surfaces]]
- [[Kills -or- Is effective against [99.999%] *Candida auris*# -or- *C. auris*# [in 3 min[utes]] [on hard, nonporous surfaces]]
- [Kills [Effective against] 99.9999% of *Clostridium difficile* [*C. difficile*] [*C. diff*] spores [in [3] [three] minutes]
- [Kills [Effective against] 99.999% of *Candida auris*# -or- *C. auris*# [in [3] [three] minutes]
- [[Kills the spore form of *Clostridium difficile* -or- *C. difficile* -or- *C. diff* [(*C. diff*)] [spores]]
- [Kills odor-causing bacteria -or- germs]
- [Kills -or- Effective against [bacteria] [and] [germs] [and] [viruses*]]
- [Kills the spore form of *Clostridium difficile* [or *C. difficile*]]
- [[Kills multi drug resistant organisms [(list organisms)]]
- [Kills multi drug resistant bacteria [(list organisms)] [in [3] [three] minutes]]
- [Kills] [Effective against] [**Insert any organism from list of organisms**] [in] [just] [3 minutes] [Kill bacteria in 3 minutes] [Kill viruses* in 3 minutes] [Kill Poliovirus* in 1 minute] [Kill Measles* -or- Measles virus* in 3 minutes] [Kill TB in 3 minutes] [Kill TB/HIV-1 in 3 minutes]
- [Kill *Clostridium difficile* -or- *C. difficile* -or- *C. diff* [(*C. diff*)] [spores] in 3 minutes]
- [Kill *Candida auris*# -or- *C. auris*# in 3 minutes]
- [Kitchen+ and Bath Wipes]
- [Large[r] Surface[s]]
- [Large[r] [wipe][s]]
- [Larger wipe covers [more surface area] larger area [with less effort][!]]
- [Lead Free]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Leaves a fresh bleach scent]
- [Leaves bathroom fresh and clean in [just [3 minutes]]]
- [Leaves bathroom fresh and clean in [just] minutes]
- [Leaves [bathroom] [office] clean [and disinfected] [and smelling fresh]]
- [Less] [Reduced] [packaging] [less] [to throw away]
- [Lightly scented -or- quick-dissipating scent]
- [Lint free] [wipe]
- [Low linting [wipe]]
- [Making home healthcare [manageable] [easy] [simple]]
- [Maximum cleaning, minimum effort]
- [Meets [the disinfection requirements of] OSHA['s] Bloodborne Pathogen Guidelines -or- Standards]
- [Mercury free]
- [Mildewcide -or- Mildewcidal]
- [Moist wipes]
- [Moisture lock lid]
- [Multi-pack]
- [Multi-purpose disinfectant]
- [Multi [-refill] [pack]
- [Multi-surface]
- [Next wipe pops up automatically [everytime] [!]]
- [Neutralizes odors [doesn't [just] mask them]
- [Non-abrasive formula]
- [No need for gloves]
- [No need for [gloves], [mask], [goggles], [protective eye] -or- [face mask] needed]
- [No canister needed]
- [No harsh chemical [odors] [or] [fumes]]
- [No [pinched fingers] [rips] [tears]]
- [No sponges, no rags, no buckets, no mess]
- [No threading [required]]
- [No fragrance added]
- [No lint wipe]
- [No need for gloves]
- [No need for [gloves], [mask], [goggles], [protective eye], [and] [or] [face mask]
- [No measuring] [required]
- [No mixing] [required]
- [No sponges, no rags, no buckets, no mess]
- [Non-scratching [wipe]]
- [Non-woven disposable cloth]
- [Non-abrasive formula]
- [Nonylphenol ethoxylate (NPE) -and/or- alkylphenol ethoxylate (APE)Free]
- [[Occupied Room] [and] [Terminal Room] [Discharge] Cleaning Wipe[s] [Towel][s] Refill[s][able]]
- [Odor [elimination] [neutralization]
- [On the go disinfectant]
- [One-step✓ cleaner and disinfectant [wipe] [saves time]]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [One step cleaning [and/or disinfecting] [and/or deodorizing] [disinfect and/or touch-up clean], when used [according to] [following] disinfection [directions] [instructions]]
- [One step cleaning [and/or disinfecting] [and/or deodorizing] [disinfect and/or touch-up clean] [cleans and/or disinfects and/or deodorizes] [anytime] [daily] [every day] [in one easy step] [on hard non-porous] [on ▲] [bathroom and/or kitchen+] surfaces [on household surfaces] [on ♦ surfaces], when used [according to] [following] disinfection [directions] [instructions]]
- [One step disinfectant cleaner that is effective against a broad spectrum of bacteria when used as directed]
- [One step germicide]
- [One-step disinfectant cleaner designed for general cleaning and disinfecting hard, nonporous environmental surfaces in ▲]
- [On-the-go disinfectant [Package] [insert type of package] made from [x%] recycled materials]
- [One [wipe] at a time dispensing]
- [Package] [Insert type of package] made from [x%] recycled materials]
- [Parvovirucide* -or- Parvovirucidal*]
- [Passes] [Meets] [Satisfies] [EPA] requirements for killing [*Clostridium difficile* spores] [*C. difficile*] spores [on precleaned, hard, nonporous surfaces].
- [Passes] [Meets] [Satisfies] [EPA] requirements for killing [*Candida auris*] [*C. auris*] [on precleaned, hard, nonporous surfaces].
- [Picks up [particles] [dust] [dirt] [hair] [crumbs]]
- [Pop-up [dispensing] [top] [tab] [lid]]
- [Portable]
- [Pouch(es)]
- [Powers through tough messes]
- [Pull n/ Go [pack]]
- [Pre-measured]
- [Pre-moistened [towelettes]]
- [Premixed], [premeasured], [premoistened]
- [Premoistened [wet] household cleaning wipe]
- [Pre-measured to have the recommended amount of germicidal solution, wipe after wipe]
- [Pre-moistened] [Easy] [and convenient] way to disinfect]
- [Pre-moistened with specially formulated cleaners to handle daily clean-ups without bottles, sponges or buckets]
- [Prevents formulation errors]
- [[Product Name] are pre-moistened with a stabilized bleach and detergent solution and are ready-to-use.]
- [[Product Name] are useful in a variety of [healthcare], [nursery/daycare] [and] [▲] applications listed on the [this] label.]
- [Pseudomonocidal]
- [Ready to use]
- [Ready-to-use disinfectant [disinfection] for [killing] *Clostridium difficile* -or- *C. difficile* -or- *C. diff* [(*C. diff*)] [spores]]
- [Ready-to-use disinfectant [disinfection] for [killing] *Candida auris* -or- *C. auris*]
- [[Ready-to-use] disinfectant [for use in ▲]]
- [Recommended for cleaning [and disinfecting] ♦]
- [Recommended for disinfecting exterior toilet surfaces - or - around the exterior toilet surface]
- [Recyclable] [package]
- [Reduce]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continue

- [Reduces cross-contamination of [*Clostridium difficile*☞] [*Clostridium difficile*☞ (*C. diff*☞)] [*C. difficile*☞] [*C. diff*☞] on treated surfaces]
- [Reduces cross-contamination of [*Candida auris*#] [*C. auris*#] on treated surfaces]
- [Reduces] [risk of cross contamination] [on treated hard non-porous surfaces]
- [Reduces] [packaging] waste [by X%] [when compared to regular canister or bucket or pail]
- [Reduces cross-contamination to (***insert any organism from EPA approved list***) [in 3 minutes on treated surfaces.]
- [Reduces cross-contamination [from (***insert any organism from EPA approved list***) [in 3 minutes) on treated (***insert non-porous surface***)]].
- [Renders *Clostridium difficile*☞ -or- *C. difficile*☞ -or- *C. diff*☞ [(*C. diff*)☞] [spores] inactive [on precleaned, hard, nonporous surfaces]]
- [Renders *Candida auris*# -or- *C. auris*# inactive [on precleaned, hard, nonporous surfaces]]
- [Refill package contains [one] [1] [two] [2] [three] [3] [four] [4] [various numbers] [wet] [pre- moistened] [wipe] [refill(s)] [roll(s)]
- [Refill [size] [pack] [value pack]
- [Refillable [with____# wipes] [dispenser] refill [pack] [for dispenser]
- [Respiratory illnesses attributable to Pandemic 2009 H1N1 are caused by Influenza A virus*. This product is a broad-spectrum hard surface disinfectant that has been shown to be effective against (Influenza A tested and listed on this label) and is expected to inactivate all Influenza A viruses* including Pandemic 2009 H1N1 (formerly called swing flu).]
- [Requires no mixing or measuring and are ready-to-use]
- [Re-sealable [stay-moist] pouch [package] [container] inside]
- [Re-use your [empty] canisters [again [and again] [!]
- [Reuse[able] [Container][Dispenser][Canister][Pail][Bucket]]
- [Sample]
- [Scrubbing power of a sponge in a disinfecting [towelettes]]
- [Shelf-stable formula]
- [So] [surfaces are more than [just] [clean] [they're] [disinfected]]
- [Specifically -or- Specially designed for the Medical/Dental/ Veterinary/Health Care [▲] industry]
- [Staphylocidal]
- [Starts cleaning on contact]
- [Strong, durable low-lint wipe [designed for sensitive equipment surfaces]]
- [[Strong,] [Durable,] cloth-like wipes]
- [Strong, durable textured wipe [cleans] [and] [deodorizes]]
- [[Strong,] [Durable,] non-woven disposable cloth]
- [[Strong,] [Durable,] non-woven wipe does not leave lint on surfaces]
- [Suitable for [▲]]
- [Suitable for use on [◆]]
- [Suitable for use on equipment requiring bleach-based products]
- [Suitable for use on [titanium-coated] [medical grade] stainless steel]
- [Suitable for use on hard, nonporous inanimate environmental surfaces]
- [Surfactant free]
- [Strong cleaning power]
- [Sturdy][Container][Dispenser][Canister][Pail][Bucket]
- [Tear at perforation]
- [Textured for better [cleaning] [scrubbing] power]
- [Textured design allows for] [scrubbing] [cleaning] [removing] [wiping] [tough] [stains] [soils] [messes]
- [Textured wipe for effective cleaning]
- [Textured] wipe is low-linting [designed] [for sensitive equipment surfaces]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Textured design [cleans] [deodorizes] hard, non-porous surfaces] [helps remove soil from surfaces]]
- [Textured] [strong] [towelettes] [for] [with] [scrubbing power]
- [The 3-minute [germicidal wipe] [disinfectant wipe]]
- [The answer to your [disinfecting] [cleaning] needs]
- [The] smell of clean
- [[The] easy -and/or- convenient way to disinfect] [with CaviWipes™ Bleach]
- [This product] controls the spread of germs on most hard, nonporous treated surfaces
- [This product has demonstrated effectiveness against Influenza A Virus* and is expected to inactivate all Influenza A Viruses* including Pandemic 2009 H1N1 Influenza A Virus*]
- [This product has demonstrated effectiveness against Measles* -or- Measles virus*]
- [This product has demonstrated effectiveness against Poliovirus* in 1 minute]
- [This product meets AOAC efficacy testing requirements -or- standards for hospital disinfection]
- [[Thick], [durable] [and] [heavy duty] wipes are bactericidal, tuberculocidal, and virucidal* in 3 minutes]
- [Thick, large wipe that [cleans] [disinfects]]
- [This] [outer] [selling] [package] [container] [canister] contains X% post-consumer [material] [plastic] [resin]
- [Time to reorder]
- [Trial]
- [Tuberculocidal (3 minutes)]
- [Tub with pop-up action]
- [To break down and kill [inactivate] [*Clostridium difficile*¥] [*C. difficile*¥] [*C. diff*¥] [(*C. diff*)¥] [spores] [on surfaces] [on precleaned, hard, nonporous surfaces].
- [To break down and kill [inactivate] [*Candida auris*#] [*C. auris*#] [on surfaces] [on precleaned, hard, nonporous surfaces].
- [To kill [inactivate] *Clostridium difficile*¥ -or- *C. difficile*¥ -or- *C. diff*¥ [(*C. diff*)¥] [spores] [on surfaces]
- [To kill [inactivate] *Candida auris*# -or- *C. auris*# [on surfaces]
- To fight [*Clostridium difficile*] [*C. difficile*] [*C. diff*] [(*C. diff*)] [spores] [on precleaned, hard, nonporous [surfaces].
- To fight [*Candida auris*#] [*C. auris*#] [on precleaned, hard, nonporous [surfaces].
- [To kill *Clostridium difficile*¥ -or- *C. difficile*¥ -or- *C. diff*¥ [(*C. diff*)¥] [spores], let stand for 3 min[utes]
- [To kill *Candida auris*# -or- *C. auris*#, let stand for 3 min[utes]
- [To kill *Clostridium difficile*¥ -or- *C. difficile*¥ -or- *C. diff*¥ [(*C. diff*)¥] [spores][,] on treated surfaces especially during outbreaks [of *Clostridium difficile*¥ -or- *C. difficile*¥ -or- *C. diff*¥ [(*C. diff*)¥] [let stand for three -or- 3 min[utes]]
- To kill *Candida auris*# -or- *C. auris*# on treated surfaces especially during outbreaks [of *Candida auris*# -or- *C. auris*#] [let stand for three -or- 3 min[utes]]
- [To reduce *Clostridium difficile*¥ -or- *C. difficile*¥ -or- *C. diff*¥ [(*C. diff*)¥] [spores] on treated surfaces in [▲]
- [To reduce *Candida auris*# -or- *C. auris*# on treated surfaces in [▲]
- [To reduce *Clostridium difficile*¥ -or- *C. difficile*¥ -or- *C. diff*¥ [(*C. diff*)¥] [spores] on precleaned, hard, nonporous treated surfaces [of ♦][in ▲]]
- [To reduce *Candida auris*# -or- *C. auris*# on precleaned, hard, nonporous treated surfaces [of ♦] [in ▲]]
- [To reduce cross-contamination from hard non-porous dental impressions for both dentist office personnel and dental laboratory personnel.]
- [[Tub] [Pail] [Bucket] refill] [available]
- [Twin pack]
- [Two in one action [cleans and disinfects] in one easy step¥]
- [Tuberculocide -or- Tuberculocidal]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Usable on variety of soils]
- [Use in [public] [common places] [▲] where [bacteria -and/or- viruses*-and/or- germs may be of concern] [on hard, nonporous surfaces]
- [Use to] [Deodorize hard, non-porous surfaces]
- [Use to] refill[s] [count] [wipe] [canisters] [tubs] [bucket] [tub] [pail]
- [Use to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization / high-level disinfection]
- [Use wipes to kill germs and bacteria instead of germ paper towels, sponges and rags]
- [Use this product to] [Kill bacteria in 3 minutes] [Kill viruses* in 3 minutes] [Kill Poliovirus* in 1 minute] [Kill Measles* -or- Measles virus* in 3 minutes] [Kill TB in 3 minutes] [Kill TB/HIV-1 in 3 minutes] [Kill *Clostridium difficile* -or- *C. difficile* -or- *C. diff*][(*C. diff*)] [spores] in 3 minutes] [Kill *Candida auris*# -or- *C. auris*# in 3 minutes]
- [Use where control of the hazards of cross-contamination between treated surfaces is of prime importance]
- [Use on hard, non-porous surfaces only: [stainless steel] [Formica®] [Porcelain] [Plastic] [Glazed Ceramic] [glass] [glass tables] [glazed tile] [Test wipe on small inconspicuous area first.]
- [Value]
- [Versatile]
- [Virtually [Lint Free] [streak free] [clean] [shine]
- [Virucide* -or- Virucidal*] [in 3 mins]
- [Virtually lint free]
- [VOC free]
- [Wet wipes]
- [When used according to label directions, [this product] reduces cross-contamination of [*Clostridium difficile*][*Clostridium difficile* (*C. diff*)] [*C. difficile*][*C. diff*] on treated surfaces]
- [When used according to label directions, [this product] reduces cross-contamination of [*Candida auris*#] [*C. auris*#] on treated surfaces]
- [Where control of the hazards of cross contamination between treated surfaces is of prime importance]
- [Wipe away odors]
- [Wipes out 99.9%] [Kills] [germs] [or bacteria] [while it cleans] [that cause odor]
- [Wipe] [material] made from [x %] [sustainable] [ingredients] [fibers] [may only be used when cellulose wipe material used]
- [Wipes out 99.9%] [Kills] [common] [household] [bacteria] [viruses*] [fungi] [clinically relevant organisms] [germs] [***[(insert any organism from EPA approved list)***] in 3 minutes]
- [Wipes out 99.9%] [Kills] [Eliminates 99.9%] [list organisms] in [just] [3 minute] [***[(insert any organism from EPA approved list)***] in 3 minutes]
- [Wipes out 99.9%] [Kills] [the] flu virus* (Influenza A/Hong Kong, H3N2 Virus*)]
- [Wipes out 99.9%] [Kills] [the] Measles* -or- Measles virus*]
- [Wipes out 99.9%] [Kills] [the] Poliovirus* in [just] [1 minute]
- [Wipes out 99.9%] [Kills] [viruses* and bacteria] [including the flu virus* (Influenza A/Hong Kong, H3N2 Virus*)]
- [Wipes out 99.9%] [Kills] [viruses* and bacteria] [including the Measles* -or- Measles virus*]
- [X" x Y" Wipes]
- [X% [more] [wipes] [free] [vs. (*insert a Metrex disinfectant wipes*) count] [canisters] [packs]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

[Note: The following claims can be used only for the first 6 months of the product on shelf]

- [New] [Strong] [Thick] Wipe [!]
- [New] [&] [dispensing] lid [closure] [cap] [package] [dispensing] [better [than before]]
- [New]
- [New] [Package] [packaging]
- [New thick and larger embossed wipe cleans more surface area]
- [Now] [Reinforced] [With] [Quilted] [strength]
- [Now with [%] more towel]
- [Now with [#] more wipes]
- [Scrubs] (&) [cleans]
- [New Strong durable [wipe]
- [New Thick [wipe] [with] cloth-like texture]
- [New Large] [quilted] wipes]
- [New Thick wipes]
- [Multi-purpose [cleaner] [disinfectant] [Multi [-refill] pack] New [large] thick [textured] wipe]
- [New] [Big] [Strong] [Thick] [wipe] [!]
- [New] [dispensing] [lid] [closure] [cap] [package] [dispenses]

OPTIONAL GRAPHICS

[THE FOLLOWING GRAPHIS ARE OPTIONAL FOR ALL CONTAINERS. THEY WILL BE ARRIXED TO THE LABEL SUCH THAT THEY DO NOT INTERFERE WITH ANY LABEL TEXT.]



THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

CLEANING AND STAIN REMOVAL

To clean -and/or- remove stains: Wipe surface with a wipe; [rinse and] let air dry. [If streaking is observed, wipe with a clean, damp [cloth] [or] [paper towel].] [[Works on] [effective against] [removes the following stains:]]

Claims:	Stains/Soils	Stains/Soils continued
<ul style="list-style-type: none">• Cleans• Removes stains• Removes (<i>insert stain(s)/soil(s) shown to the right:</i>)	<ul style="list-style-type: none">• Blood• Body oils• Bodily fluids• Coffee [stains]• Common soils• Dead skin• Defecation• Dirt• Fecal matter• Feces• Fingerprints• Food+ residue[s] -or- soil[s]• Food stains• Fruit [stains]• Grease• Grime• Laboratory stains	<ul style="list-style-type: none">• Mildew stains• Mold stains• Mucus• [Other] common soils -and/or- stains• [Other] organic matter• Phlegm• Protein-based matter• Saliva• Spit• Spittle• Sputum• Stool• Throw-up• Tea [stains]• Urine• Vomit

ALLERGEN DESTRUCTION

To [clean and] destroy specified allergens: Wipe surface. Use enough wipes for treated surfaces to remain visibly wet for 3 minute. [Rinse.] Let air dry. [If streaking is observed, wipe with a clean, damp [cloth or] paper towel after contact time has expired.]

Claims:

- Breaks Down Many Allergens[!]
- Destroys -or- deactivates dust mite matter
- [Penetrates and] breaks down dust mite matter -or- particles, cockroach matter -or- particles, and pet dander allergens [found in dust]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 1 Medical Use Sites:
CaviWipes Bleach are recommended for use on hard non-porous surfaces in [▲]

3D Imaging Service Centers	Endoscopy Centers	Ophthalmic [Facilities] [Offices]
Acute Care Facilities	ENT Facility	Optometric [Facilities] [Offices]
Alternate Care Facilities	Eye Care [Facilities][Centers]	Oral Health Care Facilities
Ambulances	Eye Surgical Centers	Orthodontic Offices
Ambulatory Surgery Centers (ASC)	Forensic Labs	Orthopedic Offices
Animal Research Facilities	Full-Care Nursing Facility	Outpatient Care Facilities
Assisted Living Facility	Gastrointestinal (GI) [Labs] [Laboratories]	Outpatient [Surgery] [Facilities] [Centers]
Autopsy Rooms	Gastrointestinal (GI) Centers	Out-Patient Surgical Centers (OPSC)
Bathroom	Gastrointestinal (GI) Suites	Pathology Rooms
Bioskill Laboratories	Home Care Services	Patient Care Areas
Blood Drawing Centers	Hospice	Patient Care Rooms
Blood Apheresis Facilities	Hospital [Operating Rooms] [Kitchens+]	Patient Rooms
Blood Units	Hydrotherapy [Facilities] [Rooms]	Pediatric Exam Rooms
Burn Center	In vitro Fertilizer Center (IVF)	Pediatric Intensive Care Units (PICU)
Cardiac Intensive Care Unit (CICU)	Instrument Reprocessing [Facilities] [Departments] [Center]	Pediatric Offices
CAT Scan [Lab] [Laboratories]	Intensive Care Units (ICU)	Personal Care Service Facilities
Cath Labs	Isolation Areas	Pharmacy Clean Rooms
Central Service Areas	Isolation Rooms	Phototherapy Rooms
Central Services Department	Labor Units	Physical Therapy Facilities
Central Supply Areas	[Lab] [Laboratories]	Physicians' Offices
Central Supply Rooms	Lasik Eye Centers	Plasma Apheresis [Centers] [Facilities] [Units]
Clean Rooms	Lasik Units	Post [Ops] [Operatory] Rooms
Clinical Areas	Laundry Rooms	Pre [Ops] [Operatory] Areas
Clinics	Life Science Laboratories	Psychiatric Care Facility
Critical Care Units	Long-Term Acute Center (LTAC)	Radiology Department
	Long-Term Care [Facilities]	Recovery Rooms
CS [Labs] [Laboratories]	Medical Facilities [Clinics]	Rehabilitation Facilities
Dialysis Centers	Medical Related Facilities	Respiratory Therapy Facilities
Dialysis Clinics	Medivans	Skilled Nursing [Facility] [Homes]
Dissection [Labs] [Laboratories]	Medical Intensive Care Unit (MICU)	Sterile Processing Departments
Doctor's Offices	Magnetic Resonance Imaging (MRI) [facilities] [laboratories] [labs]	Sterile Storage Areas
Donation Centers [blood] [plasma] semen milk aphaeresis	Neonatal Intensive Care Units (NICU)	Surgery Centers

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Drug Preparation [Centers] [Facilities] [Rooms]	Neonatal Nursery Areas	Surgery Rooms
Drug Preparation Rooms	Neonatal Units	Surgical Centers
Electrolysis Facilities	Newborn Nursery	Surgical Suites
Emergency Departments	Nuclear Medicine Facilities	Teaching Hospitals
Emergency Medical Services (EMS) Units	Nurses' Stations	Tissue Banks
Emergency Medical Settings	Nursing Homes	Transplant Center
Emergency Rooms (ER)	Nursing Stations	Trauma Center
Emergency Vehicles	Occupational Therapy [Centers] [Facilities]	University Research Lab
Endo Suites	Oncology Areas	Wound Center
Endoscopy [Centers] [Facilities] [Departments]	Operating Rooms	X-Ray Rooms

Table 2 Institutional Use Sites:

CaviWipes Bleach are recommended for use on hard non-porous surfaces in [▲]

Adult School	Early Learning Centers	Middle Schools
Air Transport Facilities	Elementary Schools	Military Installations
Airplane Interiors	[Fire][Rescue Departments]	Nursery Schools
Airports	Fire Vehicles	Police Vehicles
Bathroom	Forestry Services	Pre-K School [Room]
Casinos	Grammar School	Preparatory School
Child Care Facilities	High School	Preschool
Churches	Hospitality Facilities	Primary School
Colleges	Hotels	Prisoner Detention Facilities
Community Colleges	Industrial Facilities	Prisoner Lockups
Conference [Centers] [Rooms]	Infant Care Centers	Prisons
Conference Rooms Fixtures	Intermediate School	School Buses
Convention Facilities	Jails	School Kitchens+
Correctional Facilities	Junior High School	Schools
Daycare Centers	Kindergarten	[University] [Universities]
[Dormitories] [Dorm]	Learning Institutions	

Table 3 Dental Use Sites:

CaviWipes Bleach are recommended for use on hard non-porous surfaces in [▲]

Dental [Laboratories] [Lab]	Dental Office	Dental [Operatory] [Operatories]
-----------------------------	---------------	-------------------------------------

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 4 Vet Use Sites:

CaviWipes Bleach are recommended for use on hard non-porous surfaces in [▲]

Animal Care Facilities	Avian Units	Pet Homes
Animal Control Vehicles	Bathroom	Pet Motels
Animal Dental Operatory	Doggy Day Care Centers	Pet Spas
Animal Detention Room	Feline Facilities	Pet Stores
Animal Research Facilities	Feline Units	Small Animal Facilities
Animal Shelters	Kennel	Veterinary Offices
Avian [Care] Facilities	Pet Areas	Veterinary Facilities

Table 5 Food Prep+ Use Sites:

CaviWipes Bleach are recommended for use on hard non-porous surfaces in [▲]

Bakeries+	Food Processing Facilities+	Restaurants+
Banquet Hall Facilities+	Grocery Stores +	School Kitchens+
Butcher Shops+	Kitchen+	Supermarkets+
Food Preparation Areas+	Produce Areas+	

[+ A thorough rinse with potable water is required for surfaces in direct contact with food.]

Table 6 Miscellaneous Use Sites:

CaviWipes Bleach are recommended for use on hard non-porous surfaces in [▲]

Airplanes	Facilities	Manufacturing [Facilities] [Plant]
Arboretums	Factories	Money Handling Facilities
Arcades	Foots Spas	Morgues
[Athletic] [Training] Facilities	Foreign Current Exchange Facilities	Mortuaries
Auto Facilities	Funeral Homes	Motels
Banks	Gaming Facilities	Movie Theaters
Barber Shops	Garden Centers	Outdoor Playground
Bathroom	General Public Facilities	Personal Care Service Facilities
Beauty Salons	Golf Courses	Public Areas
Body Piercing Salons	Grounds Keeping Services	Public Facilities
Body Service Facilities	Gymnasium	Public Transportation Vehicles
Body Storage Facilities	Health Club Facilities	Recreational [center] [facilities]
Botanical Gardens	Helicopter Interiors	Recycling Centers
Bowling Alleys	Household Bathrooms	Salon Settings
Buses	Housekeeping Rooms	Shelters
Cabs	Indoor Playground	Shower Rooms
Cadaver Processing areas	Laundromats	Skincare Salons
Cars	Laundry Rooms	Storage [Areas] [Rooms]
Centers	Lawn Services	Subway
Condominiums	Limousine	Tanning Salons
Cruise Ships	Living Rooms	Tattoo Parlors
Currency Exchange	Locker Rooms	Taxi
Day Spa	Lockup Facilities	Theme Parks
Delivery Units	Manicure Salons	Transport Vehicles
Dissection Work Tables	Massage Facility	

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

**Table 7 Surfaces:
CaviWipes Bleach are recommended for use on hard non-porous surfaces of [♦]**

Ambulance Equipment Surfaces	Cafeteria Trays+	Dental Impression Trays
Animal Cages	Cages	For use as a cleaner and disinfectant on hard nonporous [dental impressions] [dental impression materials] prior to sending to a dental lab.
Animal Equipment	Cardiac Gym Equipment	Dental Operatory Surfaces
Appliances+	Carts	Dental Unit Instrument Trays
Arm Rests	Cash Register	Dentist Chair
ATM Surfaces	Ceilings	Desk Tops
Autoclave Surfaces	Cell Phones	Desks
Automatic Feeders	Chairs	Diaper Changing Stations
Automatic Teller Machines (ATM)	Changing Tables	Diaper Pails
Baby Warmers	Chemotherapy Hoods	Diaper Pails
Baked Enamel	Child Seats+	Digital Phosphor Plates
Bars	Chrome	Dining Tables+
Basins	Clipper Blades	Dish Racks+
Baskets	Closet Handles	Display Cabinets
Bathing Units	Coated Mattresses	Diving Equipment
Bathroom Fixtures	Coated Pillows	Diving Masks
Bathtubs	Combs	Diving Snorkels
Bed Pans	Computer	Door Knobs
Bed Railings [Bedside Rails]	Computer Keyboards	Dopplers
Bedside Commode	Computer Mouse	Drain Boards+
Behind and Under Sinks	Copy Machines	Drained Whirlpool Tank Surfaces
Bidet	Counter Tops+	Enamel
Biohazard Container	Counters	Endoscope Cabinet
Blood Pressure Cuff	Cribs+	Equipment Surfaces
Booster Chairs+	Crutches	Exam Tables
Booster Seats+	Dental Chairs	Examination Table
Brushes	Dental Countertops	External Surfaces of Kendrick Extrication Device (KED)
Cabinets Handles	Dental Curing Lights	Exterior Surfaces of Toilet
Cabinets	For use as a cleaner and disinfectant on hard nonporous [dental impressions] [dental impression materials] prior to sending to a dental lab.	Exterior or Exposed Interior Surfaces of Medical Equipment

- + A thorough rinse with potable water is required for surfaces in direct contact with food.
- ***For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [▲ surfaces]***
- ***Any washable (food and non-food contact) surface where disinfection is required+***

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces continued:

Exterior Surfaces of Dental Equipment	Exterior Surfaces of Vanity Tops	External Surfaces of Intubation Equipment
Exterior Surfaces of Amalgamators	External Surfaces Emergency Resuscitative Equipment	External Surfaces of Isolettes
Exterior Surfaces of Anesthesia Equipment	External Surfaces of Bladder Scanners	External Surfaces of IV Pumps
Exterior Surfaces of Aphaeresis Machines	External Surfaces of Blood Glucose Meters	External Surfaces of Life Support Equipment
Exterior Surfaces of Appliances such as Refrigerators+, Microwaves+, Food Processors+ and Stoves+	External Surfaces of BP Monitors	External Surfaces of Loupes
Exterior Surfaces of Center Venous Lines	External Surfaces of CAT Scan Machines	External Surfaces of Medical Equipment
Exterior Surfaces of CICU Equipment	External Surfaces of Compression Devices	External Surfaces of Medical Warmers
Exterior Surfaces of Defibrillators	External Surfaces of Dialogistic Equipment	External Surfaces of Monitors Electrodes
Exterior Surfaces of Dental Care Facilities Hand Pieces	External Surfaces of Dictating Equipment	External Surfaces of MRI Machines
Exterior Surfaces of Dental Equipment Lines	External Surfaces of Digital X-Ray Sensors	External Surfaces of OLED Monitors
Exterior Surfaces of Dental Equipment Tubing	External Surfaces of EKG Machines	External Surfaces of Ophthalmoscopes
Exterior Surfaces of Dental Operatories	External Surfaces of Electronic Thermometers	External Surfaces of Otopscopes
Exterior Surfaces of Diagnostic Equipment	External Surfaces of Electrosurgical Generators	External Surfaces of Patient Care Monitor Unit
Exterior Surfaces of Dialysis Machines	External Surfaces of Equipment	External Surfaces of Patient Monitoring Equipment
Exterior Surfaces of Endodontic Equipment (Apex Locators, Pulp Testers, Motors)	External Surfaces of Glucometers	External Surfaces of Probes
Exterior Surfaces of Radiology Equipment	External Surfaces of Headsets	External Surfaces of Pulp Tester Motors
Exterior Surfaces of Refrigerator Units+	External Surfaces of ICU Equipment	External Surfaces of Pulse Oximeter
Exterior Surfaces of Respiratory Therapy Equipment	External Surfaces of Incubators	External Surfaces of Sequential Compression Devices
Exterior Surfaces of Urinals	External Surfaces of Infusion Pumps	External Surfaces of Respiratory Ventilators
<ul style="list-style-type: none"> • + A thorough rinse with potable water is required for surfaces in direct contact with food. • <i>For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [▲ surfaces]</i> • <i>Any washable (food and non-food contact) surface where disinfection is required+</i> 		

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces continued:

External Surfaces of Resuscitators	Fire Trucks	Handles
External Surfaces of Sterilizers	Floors	Handrails
External Surfaces of Suction Pumps	Food Cases+	Hard Plastic or Metal Flooring
External Surfaces of Surgery Microscopes	Food Serving Trays+	Hard, non-porous environmental surfaces
External Surfaces of Surgery Robots	Formica	High Chairs+
External Surfaces of Surgical Headlights	Fountains	High Touch Objects
External Surfaces of Surgical Light Sources	Freezers+	Hoods+
External Surfaces of Tympanic Thermometers	Garbage Cans	Hospital Bed (Railings, Frames)
External Surfaces of Ultrasound Machines	Glass	Hospital Patient Bed
External Surfaces of Ultrasound Monitors	Glazed Ceramic Tile	Infant Incubators (interior and exterior surfaces)
External Surfaces of Ultrasound Transducers	Glazed Porcelain	Inside of Drawers
External Surfaces of Vaporizer	Glazed Tiled Walls	Instrument Stands
External Surfaces of Veterinary Equipment	Glazed Tiles	IV Poles
External Surfaces of Vital Sign Monitors	Golf Cart	IV Stands
External Surfaces of X-Ray Machines	Grab Bar	Jacuzzi
Eye Glasses (do not use on contact lenses)	Grocery Cart	Kennels
Facilities [Surfaces]	Gurneys	Knee Walkers
Faucets	Gym Equipment	Knobs
Fax Machines	Hair Clippers	Laboratories
Feed Racks	Hair Cutting Implements	Laboratory Equipment
Filing Cabinets	Hair Dryers	Laboratory Surfaces
Filing Doors	Hair Rollers	Laminated Surfaces
Filing Handles	Hampers	Lamps
<ul style="list-style-type: none"> • + A thorough rinse with potable water is required for surfaces in direct contact with food. • <i>For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [▲ surfaces]</i> • <i>Any washable (food and non-food contact) surface where disinfection is required+</i> 		

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses, use sites and surfaces may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces continued:

Light Lens Covers	Pipes	Spine Back Boards
Lights	Plastic Laminate	Stainless Steel
Linoleum	Plastic Mattress Covers	Stalls
Manicure Implements	Playpens	Stands
Manicure Tables	Pneumatic Tourniquets	Stethoscopes
Marlite	Police Cars	Stools
Mattress Covers	Porcelain Enamel	Storage
Mayo	Posey Vest	Stoves+
Mayo Tray	Public Restrooms	Stovetops+
Monitors	Razors	Stretchers
Nail/Hair Care Implements	Reception Areas	Stretches
Non-food contact surfaces+	Reception Counters	Support Bars
Nurse Call [Button] [Device] [Cord]	Refrigerator Carts	Synthetic Marble
Operating Room Furniture	[Exterior Surfaces of] [Outside of] Refrigerators+	Tables+
Operating Room Lights	Safety Straps	Tanning Beds
Operating Room Tables	Salad Bar Sneeze Guards+	Telemetry Equipment
Operating Tables	Salon Surfaces	Telephones
Operatory Light Switches	School Buses	Tiles+
Outdoor Furniture+ (excluding wood frames and upholstery)	Scissors	Toilet
Ovens+	Shampoo Bowls	Toilet (Rim, Seat, Tank, Handle,)
Overbed table	Shears	Touch Screen Monitors
Oxygen Hoods	Shelves	Towel Dispensers
Patient Chairs	Shopping Carts	Toys
Patient Delivery Equipment	Shower Chair	Transfer Bench
Patient Lift Equipment	Shower Fixtures	Transport Vehicles
Patient Scales	Shower Stalls	Trash Cans
Patient Support Equipment	Showers	Tray Tables
Pens	Silk Peel™ with DermalInfusion™	Troughs
Phone Chargers	Sinks+	Tubing
Phone Cradle	Slit Lamps	Vanity Tops+
Physical Therapy (PT) Equipment surfaces	Smart Phone Chargers	Veterinary Care Surfaces
Vinyl	Washable Nail Files	Wheelchairs
Walkers	Watering Appliances	Whitening Lamps
Wash Basin	Weight-lifting Surfaces	Workstations

- + A thorough rinse with potable water is required for surfaces in direct contact with food.
- ***For use on Hard Non-Porous [Food Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [▲ surfaces]***
- ***Any washable (food and non-food contact) surface where disinfection is required+***