



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
WASHINGTON, DC 20460

OFFICE OF CHEMICAL SAFETY
AND POLLUTION PREVENTION

February 8, 2021

Rebeca Lei
Senior Regulatory Specialist
Metrex Research
1717 West Collins Avenue
Orange, CA 92867

Subject: Label Amendment – Add efficacy claim for SARS-CoV-2 to the label
Product Name: CWN-07-W
EPA Registration Number: 46781 -17
Application Date: September 28, 2020
Decision Number: 566548

Dear Ms. Lei:

The amended label referred to above, submitted in connection with registration under the Federal Insecticide, Fungicide and Rodenticide Act, as amended, is acceptable. This approval does not affect any conditions that were previously imposed on this registration. You continue to be subject to existing conditions on your registration and any deadlines connected with them.

Because you have opted to add statements pertaining to emerging viral pathogens to your label as described in the [August 19, 2016, Guidance To Registrants: Process For Making Claims Against Emerging Viral Pathogens Not On EPA-Registered Disinfectant Labels \(“Guidance”\)](https://www.epa.gov/sites/production/files/2016-09/documents/emerging_viral_pathogen_program_guidance_final_8_19_16_001_0.pdf), https://www.epa.gov/sites/production/files/2016-09/documents/emerging_viral_pathogen_program_guidance_final_8_19_16_001_0.pdf, you are subject to the following additional terms of registration:

1. You may make statements pertaining to emerging viral pathogens only through the following communications outlets: technical literature distributed exclusively to health care facilities, physicians, nurses and public health officials, "1-800" consumer information services, social media sites and company websites (non-label related). These statements shall not appear on marketed (final print) product labels.
2. Your statements pertaining to emerging viral pathogens must adhere to the format approved on the Agency-accepted master label.
3. You may make statements pertaining to emerging viral pathogens only upon a disease outbreak that meets all the following criteria:
 - a. The causative organism must be a virus that causes an infectious disease that has appeared in a human or animal population in the U.S. for the first time, or that may have existed previously but is rapidly increasing in incidence or geographic range.
 - i. For human disease, the outbreak is listed in one of the following Centers for Disease Control (CDC) publications:

- A. CDC Current Outbreak List for “U.S. Based Outbreaks” (www.cdc.gov/outbreaks),
 - B. CDC Current Outbreak List for “Outbreaks Affecting International Travelers” with an “Alert” or “Advisory” classification (www.cdc.gov/outbreaks) (also released through the CDC’s Health Alert Network (HAN) notification process)
 - C. Healthcare-Associated Infections (HAIs) Outbreaks and Patient Notifications page (www.cdc.gov/hai/outbreaks)
- ii. For animal disease, the outbreak is identified as an infectious disease outbreak in animals within the U.S. on the World Organization for Animal Health (OIE) Weekly Disease Information page (www.oie.int/wahis_2/public/wahid.php/Diseaseinformation/WI).
- b. The CDC or OIE has identified the taxonomy, including the viral family and/or species, of the pathogen and provides notice to the public of the identity of the emerging virus that is responsible for an infectious disease outbreak. Based on the taxonomy of the outbreak pathogen identified by the CDC or OEI, the pathogen's viral subgroup is enveloped.
 - c. The virus can be transmitted via environmental surfaces (non-vector transmission), and environmental surface disinfection has been recommended by the CDC, OIE or EPA to control the spread of the pathogen.
4. You may begin communicating statements pertaining to emerging viral pathogens only upon CDC or OIE’s publication per term B.3.a. of an outbreak of an emerging viral pathogen meeting all of the criteria of term B.3. You must cease and remove all such non-label communications intended for consumers no later than 24 months after the original publication of the outbreak per term B.3.a., unless the Agency issue written guidance to the contrary due to continued public health concerns. The emerging pathogen claim language may remain on the master label.
 5. Terms B.1 through B.4 above shall become immediately void and ineffective if registration for use against SARS-CoV-2 is suspended or cancelled or no longer meets the criteria for a disinfectant claim (see EPA Product Performance Test Guideline 810.2200). In addition, terms B.1 through B.4 above shall become immediately void and ineffective upon your receipt of evidence of ineffectiveness against any pathogen in a less-resistant Spaulding category.

A stamped copy of your labeling is enclosed for your records. This labeling supersedes all previously accepted labeling. You must submit one copy of the final printed labeling before you release the product for shipment with the new labeling. In accordance with 40 CFR 152.130(c), you may distribute or sell this product under the previously approved labeling for 18 months from the date of this letter. After 18 months, you may only distribute or sell this product if it bears this new revised labeling or subsequently approved labeling. “To distribute or sell” is defined under FIFRA section 2(gg) and its implementing regulation at 40 CFR 152.3.

Should you wish to add/retain a reference to the company's website on your label, then please be aware that the website becomes labeling under the Federal Insecticide Fungicide and Rodenticide Act and is subject to review by the Agency. If the website is false or misleading, the product would be misbranded and unlawful to sell or distribute under FIFRA section 12(a)(1)(E). 40 CFR 156.10(a)(5) list examples of statements EPA may consider false or misleading. In addition, regardless of whether a website is referenced on your product's label, claims made on the website may not substantially differ from those claims approved through the registration process. Therefore, should the Agency find or if it is brought to our attention that a website contains false or misleading statements or claims substantially differing from the EPA approved registration, the website will be referred to the EPA's Office of Enforcement and Compliance.

Your release for shipment of the product constitutes acceptance of these conditions. If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA section 6. If you have any questions, please contact Mohammad Alavi by phone at (703) 347-0522, or via email at Alavi.mohammad@epa.gov.

Sincerely,



Jacqueline Hardy Manager 34
Regulatory Management Branch II
Antimicrobials Division (7510P)
Office of Pesticide Programs

Enclosure Stamped Label

CWN-07-W

[TUBERCULOCIDAL] [BACTERICIDAL] [VIRUCIDAL*]
[Fungicidal[†] [Pathogenic Fungi[†]] against *C. auris*†]

<p>Active Ingredient:</p> <p>Isopropyl Alcohol..... 17.20%</p> <p>Didecyl dimethyl ammonium chloride0.17%</p> <p>Alkyl (C14-50%, C12-40%, C16-10%) dimethyl benzyl ammonium chloride0.11%</p> <p>Other Ingredients† 82.52%</p> <p>Total 100.00%</p> <p>[†Does not include the weight of the towelette]</p>	<p>KEEP OUT OF REACH OF CHILDREN</p> <p>CAUTION [See other panel for additional Precautionary and First Aid statements.]</p> <div data-bbox="967 669 1360 957" style="border: 2px solid black; padding: 5px; text-align: center;"><p>A C C E P T E D</p><p>2/08/2021</p><p>Under the Federal Insecticide, Fungicide and Rodenticide Act as amended, for the pesticide registered under EPA Reg. No. 46781-17</p></div>
<p>EPA Reg. No. 46781-17</p> <p>EPA Est. No. 8251-WI-003 [(A1)], 8251-WI-004 [(A2)], 8251-WI-005 [(A3)], 8251-OH-001 [(A4)], 46781-MI-001 [(A5)], 56952-WI-001 [(A6)] 2686-WI-008 [(A7)], 71979-SC-001 [(A8)] 75272-WI-001 [(A9)], 93725-TN-001 [(A10)] 0312-WI-3 [(A11)]</p> <p>[Beginning of batch code indicates Est No.] [End of batch code indicates Est No.]</p> <p>Made in USA</p>	<p>Reorder No: Net Contents: #[[lbs.] [oz]] # Pre-Saturated Towelette # x # in.</p> <p>[Packed For:] [Distributed By:] [Sold By:] [Manufactured By [or for]:] [Manufacturer:]</p> <p>Metrex Research 28210 Wick Road Romulus, MI 48174</p> <p>[To reorder] [,] [For product or technical information] contact Metrex at 800-841-1428 or visit our web site at www.metrex.com</p> <p>[For SDS information, visit www.metrex.com]</p>

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRECAUTIONARY STATEMENTS Hazards to Humans and Domestic Animals

CAUTION

Causes moderate eye irritation. Avoid contact with eyes or clothing. Wear protective eyewear. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

PHYSICAL AND CHEMICAL HAZARDS: Combustible. Do not use or store near heat or open flame.

FIRST AID

If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.

[For Chemical Emergency ONLY (spill, leak, fire, exposure or accident), call CHEMTREC at 1-800-424-9300]

For [ALL] [other] inquiries about this product, call Metrex Customer Care at 1-800-841-1428, Monday through Friday, between 6 a.m. and 4 p.m. Pacific Time.

Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Packaging Statements and Wipe Dispensing Instructions:

[Use A, B, or C as applicable depending on the wipe container/dispenser type.]

[A – Refillable, Non-refillable canister, buckets, tubs or pails]

[To dispense towelettes] [container] [canister] [pail] [bucket] [tubs] [Dispensing Directions]:

[To start feed, remove cover and discard seal. From center of towelette roll, pull up a towelette corner, twist it into a point and thread it through the hole in the container cover. Pull through about one inch. Replace cover. Pull out first towelette and tear off at [(any angle)] [45°-90° angle]. Remaining towelettes feed automatically, ready for next use. [When not in use, keep center cap closed to prevent solution loss.]]

[To refill] [canister] [bucket] [pail]: Place refill [pouch] [bag] in bucket [canister]. Remove the lid from the [pail] [bucket] [canister]. [Use scissors] [to] [Cut] [Rip] [Tear] [Pull] [Open] refill pouch [bag]. Place [put] [insert] refill pouch into [bucket] [pail] [canister]. DO NOT remove wipe roll from refill pouch [bag]. From center of towelette roll, pull up a towelette corner, twist it into a point and thread it through the hole in the [bucket] [pail] [container] cover. Pull through about one inch. Replace cover. Pull out first towelette and tear off at [(any angle)] [45°-90° angle]. Remaining towelettes feed automatically, ready for next use. When not in use, keep center cap closed to prevent solution loss.]

[B – Flat Pack]

[To dispense towelettes] [from flat pack:] [from Flat Pack wipes container:]

[Grasp the white tab that states "PULL HERE". Gently swing the lid open [peel back lid] and grasp a towelette. Remove towelette from package and then swing lid closed [close lid]. Gently tap the lid to ensure seal is made between it and the towelette package.]

[C – Single Use Packet]:

[Tear open packet. Unfold towelette and use as directed. Discard as directed.]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

EMERGING VIRAL PATHOGEN CLAIMS

This product qualifies for emerging viral pathogen claims per the EPA's 'Guidance to Registrants: Process for Making Claims Against Emerging Viral Pathogens not on EPA-Registered Disinfectant Labels' when used in accordance with the appropriate use directions indicated below. This product meets the criteria to make claims against certain emerging viral pathogens from the following viral category:

- Enveloped Viruses
- Large Non-enveloped
- Small Non-enveloped

<i>For an emerging viral pathogen that is a/an...</i>	<i>...follow the directions for use for the following organisms on the label:</i>
Enveloped virus	Adenovirus Type 2
Large, non-enveloped virus	Feline Calicivirus
Small, non-enveloped virus	Feline Calicivirus Enterovirus Type 71

[Product name] has demonstrated effectiveness against viruses similar to [name of emerging enveloped/small and large non-enveloped virus] on hard, non-porous surfaces. Therefore, [product name] can be used against [name of emerging enveloped/small and large non-enveloped virus] when used in accordance with the directions for use against (cited above) on hard, non-porous surfaces. Refer to the [CDC or OIE] website at [pathogen-specific website address] for additional information.

[Name of illness/outbreak] is caused by [name of emerging enveloped/small and large non-enveloped virus]. [Product name] kills similar viruses and therefore can be used against [name of emerging enveloped/small and large non-enveloped virus] when used in accordance with the directions for use against (cited above) on hard non-porous surfaces. Refer to the [CDC or OIE] website at [website address] for additional information.

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.

Note to Reviewer: Instructions for Residential and/or Non-Residential Use Sites will be used depending on Use Sites mentioned on the Production Labels.

[Use A, B or C below as applicable depending on the wipe container/ dispenser type.]

- **[A – Nonrefillable Containers, Canisters, Buckets or Pails]**
- **[B – Refillable Containers, Canisters, Buckets or Pails [to be refilled by the user]]**
- **[C – Flat Packs or Packets]**

[Note to Reviewer: the following Pesticide Storage directions apply to Residential and Non-Residential Use Sites for A, B or C container types above]

PESTICIDE STORAGE: [Store at room temperature] -or- [Store product in a cool [dry] area, [away from direct sunlight and heat] to avoid deterioration.] [When not in use, keep [center cap of lid] [flap] [tab] closed to prevent moisture loss.]

Residential Use Sites:

[B – Refillable Containers, Canisters, Buckets or Pails [to be refilled by the user]]

TOWELETTE DISPOSAL: Dispose of used towelette in trash after use. [Do not re-use towelette.] [Do not flush in toilet.]

[Refillable] Container Handling: Refillable Container. Do not reuse or refill except as described in direction for use. Refill only with this product. Offer for recycling. If recycling is not available, discard in trash.

[C – Flat Packs or Packets]

TOWELETTE DISPOSAL: Dispose of used towelette in trash after use. [Do not re-use towelette.] [Do not flush in toilet.]

Container Handling: Nonrefillable container. Do not reuse or refill this container. Wrap container and put in trash or offer for recycling if available.

Non-Residential Use Sites:

[A – Nonrefillable Containers, Canisters, Buckets or Pails]

[C – Flat Packs or Packets]

PESTICIDE DISPOSAL: To avoid waste, use all material in this container by application according to label directions. If wastes cannot be avoided, offer remaining product to a waste disposal facility or pesticide disposal program. [Do not re-use towelette.]

Container Handling: Nonrefillable container. Do not reuse or refill this container. Offer for recycling if available or dispose of in a sanitary landfill or by other procedures approved by state and local authorities.

[B – Refillable Containers, Canisters, Buckets or Pails [to be refilled by the user]]

PESTICIDE DISPOSAL: To avoid waste, use all material in this container by application according to label directions. If wastes cannot be avoided, offer remaining product to a waste disposal facility or pesticide disposal program. [Do not re-use towelette.]

Container Handling: Refillable Container. Do not reuse or refill except as described in direction for use. Refill only with this product. Offer for recycling if available or dispose of in a sanitary landfill or by other procedures approved by state and local authorities.

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

To [Clean and] Disinfect [and Deodorize] Hard, Non-Food Contact, Non-Porous Surfaces:

Remove Visible filth and soil [loads] [prior to disinfecting]. To disinfect hard, non-porous surfaces. Use [enough] [Product Name] wipes to treat the surface, allow surface to remain visibly wet for 2 minutes [contact time] [dwell time]. [Rinse.] [Let air dry]. [If streaking is observed, wipe with a clean, damp [cloth or] paper towel after 2-minute contact time [dwell time] has expired.] [Use 2-minute contact time [dwell time] for (Tuberculosis) [-or-] (TB).] [Use 2-minute contact time to disinfect against SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces.] [To kill SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces, use 2-minute contact time.]

(#) Special Label Instructions for Cleaning Prior to Disinfection against *Candida auris*:

Personal Protection: Wear [appropriate] barrier protection [such as] gloves, gowns, masks, or eye covering.

Cleaning Procedure: Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with [Product Name] [this product saturated with product name]. Pre-cleaning is to include vigorous wiping and/or scrubbing and all visible soil is removed. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading the organism. Restrooms are to be cleaned last. Do not reuse soiled cloths.

Infectious Waste Disposal: Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

[The following HIV-1 / HBV / HCV instructions are to be included on production labels when claims against HIV-1 or HBV or HCV are listed:]

(Y) This product effectively kills [Human Immunodeficiency Virus* Type 1 (HIV-1),] [Hepatitis B Virus* (HBV),] [[and] Hepatitis C Virus* (HCV)] on pre-cleaned hard, nonporous environmental surfaces/objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/ body fluids can be associated with the potential for transmission of [HIV-1,] [HBV,] [[and] [HCV].

Special Instructions for Cleaning and Disinfection Against [HIV-1,] [HBV,] [[and] HCV] of Surfaces or Objects Soiled with Blood or Body Fluids:

Personal Protection: When handling items soiled with blood or body fluids, wear [appropriate] barrier protection: [such as] disposable protective gloves, gowns, face masks and eye coverings.

Cleaning Procedure: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of this disinfecting towelette. A [Product Name] towelette can be used for this purpose.

Contact Time: Thoroughly wet surface with a new [Product Name] towelette. Ensure surface remains visibly wet for 2 minutes.

Disposal of Infectious Materials: Blood and other body fluids must be autoclaved and disposed of in addition to used wipe in accordance with Federal, State or local regulations for infectious materials disposal.

For use as a Disinfectant on Hard, Non-Porous Surfaces of Non-Critical Medical Devices, Instruments and Implements:

Instruments must be thoroughly cleaned to remove excess organic debris, rinsed and dried. Thoroughly clean and rinse lumens of hollow instruments. Then thoroughly treat [wet] the instruments with the [Product Name] product towelette. Ensure surface remains visibly wet for 2 minutes. [Use 2-minute contact time [dwell time] for (Tuberculosis) [-or-] (TB).] [Use 2-minute contact time to disinfect against SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces.] [To kill SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces, use 2-minute contact time.]

Remove and rinse instruments. Wipe dry prior to use. Discard wipe after each use.

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

DIRECTIONS FOR USE *continued*

[The following direction is to be included when hard, non-porous exterior surfaces of ultrasound transducers or probes are listed as a use site:]

To disinfect external surfaces of ultrasound transducers that do not go into the human body only: Pre-clean ultrasound gel from external transducer to be disinfected [by using a clean washcloth]. Then thoroughly [treat] wet the probe with the product [Product Name] towelette. Ensure surface remains visibly wet for 2 minutes. [Use 2-minute contact time [dwell time] for (Tuberculosis) [-or-] (TB).] [Use 2-minute contact time to disinfect against SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces.] [To kill SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces, use 2-minute contact time.]

[The following paragraph will be included in its entirety when the production label references medical devices:]

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream, or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean critical or semi-critical medical devices prior to sterilization or high level disinfection.

This product may be used to disinfect non-critical medical devices, which come in contact only with intact skin.

[The following statement is to be included on production labels when animal housing facilities are listed: (Including Veterinary Clinics, Farms and Poultry Houses)]

For cleaning and disinfecting the following hard, non-food contact, non-porous surfaces: feeding and watering equipment, utensils, instruments, cages, kennels, stables, catteries, etc.

1. Remove all animals and feeds from premises, animal transportation vehicles, crates, cages etc.
2. Remove all litter, droppings and manure from floors, walls and surfaces of facilities (including pens, chutes, and barns) occupied or traversed by animals.
3. Empty all troughs, racks, and other feeding and watering appliances.
4. Thoroughly clean all surfaces with soap or detergent and rinse with potable water. Ensure surface remains visibly wet for 2 minutes.
5. Wipe all handling and restraining equipment such as leashes, muzzles, halters or ropes, as well as forks, shovels, and scrapers used for removing litter and manure.
6. Ventilate buildings and other closed spaces.
7. Do not house animals or employ equipment until treatment has been absorbed, set or dried.
8. Thoroughly scrub treated feed racks, troughs, mangers, automatic feeders, fountains, and waterers with soap or detergent, and rinse with potable water before reuse.

+ For non-food contact surfaces only.

[DO NOT USE AS A BABY WIPE.] [Not for Use On Skin] [DO NOT FLUSH WIPE] [DO NOT RE-USE WIPE] [SINGLE USE ONLY]

[Not for cleaning skin.] [Do not use as a diaper wipe or for personal cleansing.] [APPLY CAREFULLY] [Refill only with this product.] [Always refer to manufacturer's care instructions before using on equipment -or- devices.]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

DISINFECTION [Organism List]

[[Product Name] [Disinfecting Towelettes] are non-woven disposable towelettes pre-saturated with the [Product Name] solution.]

[[Product Name] [Disinfecting Towelettes] are multi-purpose [disinfectant]/ [cleaning] wipes for use on hard, non-food contact, non-porous inanimate surfaces.]

[[Product Name] [This product] is effective against the following microorganisms on hard, non-food contact, non-porous surfaces when used as directed:]

[[Product Name] [This product] has demonstrated effectiveness against the following organisms with a 2-minute contact time [dwell time]. Use 2-minute contact time [dwell time] for Tuberculosis -or- TB.] [Use 2-minute contact time [dwell time] for SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces.]

[Disinfects (*insert organisms from Organism List*) on hard, non-food contact, non-porous environmental surfaces listed on this label]

[When used as directed, [this product] [Product Name] provides effective Bactericidal / Virucidal* / Pathogenic Fungi[†] / Tuberculocidal activity against the following organisms:]

Pathogen:	Bacteria	Viruses*	Bloodborne Pathogens**	ESKAPE Pathogens***	TB	<i>Candida auris</i> ‡	SARS-CoV-2*
Contact Time:	2 min[utes]	2 min[utes]	2 min[utes]	2 min[utes]	2 min[utes]	2 min[utes]	2 min[utes]

ORGANISMS	
	[Strain and/or ATCC number]
Bacteria (2 Minute Contact)	
<i>Staphylococcus aureus</i> ***, [^^]	[ATCC 6538]
<i>Pseudomonas aeruginosa</i> ***, [^]	[ATCC 15442]
<i>Salmonella enterica</i>	[ATCC 10708]
<i>Legionella pneumophila</i>	[ATCC 33153]
<i>Bordetella pertussis</i>	[ATCC 12743]
<i>Moraxella catarrhalis</i>	[ATCC 25240]
<i>Enterococcus faecium</i> ***	[ATCC 51559]
<i>Klebsiella pneumoniae</i> ***	[ATCC 4352]
<i>Acinetobacter baumannii</i> ***	[ATCC 19606]
<i>Escherichia coli</i> O157:H7	[ATCC 35150]
<i>Micrococcus luteus</i>	[ATCC 14408]
<i>Streptococcus mutans</i>	[ATCC 25175]
<i>Burkholderia cepacia</i>	[ATCC 25416]
<i>Stenotrophomonas maltophilia</i>	[ATCC 25596]
<i>Streptococcus pyogenes</i>	[ATCC 12344]
[Multi] Drug-Resistant Bacteria (2 Minute Contact)	
Penicillin-Resistant <i>Streptococcus pneumoniae</i> (PRSP)	[ATCC 700677]
Vancomycin-Resistant <i>Enterococcus faecalis</i> (VRE) ***	[ATCC 51575]
Methicillin-Resistant <i>Staphylococcus aureus</i> (MRSA) ***, ^^	[ATCC 33592]
Carbapenem-resistant <i>Klebsiella pneumoniae</i> (CRKP) ***	[ATCC BAA-1705]
Multidrug-Resistant <i>Acinetobacter baumannii</i> (MRAB) *** (Resistant to Ceftazidime and Gentamicin)	[ATCC BAA-1605]
NDM-1 (New Delhi metallo-beta-lactamase 1) [-Positive] <i>Enterobacter cloacae</i> ***	[AR-0038] [from CDC]
Vancomycin-Intermediate <i>Staphylococcus aureus</i> (VISA) ***, ^^	[ATCC 700699]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

ORGANISMS <i>continued</i>	
	[Strain and/or ATCC number]
<i>Mycobacterium</i> (2 Minute Contact)	
<i>Mycobacterium tuberculosis</i> var: <i>bovis</i> (BCG) (TB)	[ATCC 35743]
Other Organisms [Pathogenic Fungi^π] (2 Minute Contact)	
<i>Candida auris</i> ‡	[AR-BANK#0381] [from CDC]
Viruses * [Non-enveloped][:] (2 Minute Contact)	
Adenovirus [[Type] 2] *	[ATCC VR-846] [Adenoid 6]
Feline Calicivirus surrogate for Human Norovirus [Feline Calicivirus] * -or- [Human] Norovirus – Feline Calicivirus as surrogate *	[ATCC VR-782]
Feline Calicivirus *	[ATCC VR-782]
Human Rotavirus *	[ATCC VR-2018] [strain WA]
Enterovirus [[Type] 71] *	[ATCC VR-1432] [Strain H]
Echovirus [[Type] 11] *	[ATCC VR-41] [Strain Gregory]
Viruses * [Enveloped][:] (2 Minute Contact)	
Herpes simplex virus [type 1] *	[ATCC VR-733] [Strain F(1)]
Herpes simplex virus [type 2] *	[ATCC VR-734] [Strain G]
Influenza A virus [Strain A/Hong Kong 8/68] [(H3N2)] *	[ATCC VR-544] [Strain A/Hong Kong 8/68]
Influenza B virus [Strain B/Hong Kong/5/72] *	[ATCC VR-823] [Strain B/Hong Kong/5/72]
Parainfluenza virus [type 3] *	[ATCC VR-93] [Strain C243]
Human orthopneumovirus [(Respiratory Syncytial Virus)] *	[ATCC VR-26] [Strain Long]
SARS-associated Coronavirus *	[CDC 200300592] ZeptoMetrix]
Measles virus *	[ATCC VR-24] [Strain Edmonston]
Duck Hepatitis B virus (DHBV) [surrogate for Human Hepatitis B virus (HBV)] *, ** -or- Hepatitis B virus (HBV) - Duck Hepatitis B virus as surrogate *, **	[Strain Grimaud] [HepadnaVirus Testing, Inc.]
Bovine Viral Diarrhea virus (BVDV) [surrogate for Human Hepatitis C virus (HCV)] *, ** -or- Human Hepatitis C virus (HCV) - Bovine Viral Diarrhea virus as surrogate *, **	[Strain NADL] [American BioResearch Lab]
HIV-1 [(Human Immunodeficiency virus type 1)] *, **	[Strain IIIB] [ZeptoMetrix]
Severe Acute Respiratory Syndrome-related Coronavirus 2* -or- SARS-CoV-2* -or- COVID-19 Virus*	[Strain: USA-WA1/2020] [BEI Resources, NR-52281]

* Virus

** Bloodborne pathogens

*** ESKAPE pathogens, Not for California.

‡ *Candida auris*

Follow [the] Special Instructions for Cleaning required prior to Disinfection [listed on label]

^ Pseudomonacidal [Note to Reviewer: the qualifier ^ will only be used when the term Pseudomonacidal is used]

^^ Staphylocidal [Note to Reviewer: the qualifier ^^ will only be used when the term Staphylocidal is used]

π Pathogenic Fungi / Fungicidal against *Candida auris* ‡ only

Y HIV-1 / HBV / HCV

Follow [the] Special Instructions prior to Disinfection [listed on label]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS

▲ **to represent any use site(s) listed in below as applicable.**

◆ **to represent any surface(s) listed in below as applicable.**

- [# pre-saturated towelettes] [X" x Y"]
- [2-minute contact time [dwell time] for bacteria [,] [or] viruses* [,] [or] [multi] drug-resistant bacteria [pathogens] [,] [or] bloodborne pathogens** [,] [or] ESKAPE pathogens*** [,] [or] Tuberculosis -or- TB [,] [or] *Candida auris*† -or- *C. auris*† [,] [or] SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces]
- [2-minute contact time [dwell time]]
- [2-minute contact time [dwell time] against or for TB -or- tuberculosis]
- [2-minute contact time [dwell time] against/for SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces [in the presence of [5%] soil load]]
- [99.999% -or 5 [five] log reduction for [*Candida auris*†] -or- [*C. auris*†] [in [just] 2 minutes]
- [99.9% -or 3 [three] log reduction against/for SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces [in [just] 2 minutes] [in the presence of [5%] soil load] [in one-stepY]]
- [Kills a minimum of 99.999% -or- 5 [five] log of *Candida auris*† -or- *C. auris*† on hard, non-porous surfaces] [in [just] 2 minutes]
- [Kills 99.9% -or- 3 [three] log of SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces [in [just] 2 minutes] [in one stepY] [in the presence of [5%] soil load]]
- [A Metrex Research® brand]
- [A Metrex® brand]
- [A CaviWipes™ Brand]
- [A disinfectant cloth for use in ▲ on ◆]
- [[Aids] [Helps] in reducing -or- in the reduction of cross-contamination between treated hard non-porous surfaces]
- [An effective method to clean against odor-causing organisms on hard, non-porous surfaces listed on this label]
- [[A] one-stepY cleaner/disinfectant]
- [[A] one-stepY [cleaner and] disinfectant [towelette] against SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces]
- [A disinfectant towelette for use on ▲ and ◆ surfaces]
- [After use, just toss wipe(s) away]
- [Always ready-to-use]
- [An effective method to deodorize and disinfect against odor-causing organisms on hard, non-porous surfaces]
- [An effective way to disinfect hard, non-porous environmental surfaces] [listed on this label]
- [A non-woven disposable towelettes for use in ▲ and other critical care areas where reduction of cross-contamination between treated hard non-porous surfaces is required]
- [Antibacterial]
- [Antimicrobial]
- [A One-StepY Cleaner [and disinfectant] [except when precleaning is required for [*C. auris*†,] blood and body fluids]]
- [A solution impregnated in a wiping cloth]
- [[Appropriate] [for] Cleaning Patient Discharge]
- [[Appropriate] for Occupied Room -and/or Terminal Room Cleaning -and/or- Room Discharge Cleaning] (**Note to Reviewer: this claim is for cleaning only.**)
- [[Appropriate] [for] Cleaning and disinfecting dental operatories]
- [[Appropriate] [for] Cleaning and disinfecting medical sites]
- [[Appropriate] [for] Cleaning and disinfecting home]
- [Assures proper dilution strength, product effectiveness and standardizes technique]
- [Kills and/or Inactivates (*insert any organism(s) from Organism List*)] [on hard, non-porous surfaces]]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Bacteria -and/or- virus* -and/or- germ killing [formula]]
- [Bactericide -or- Bactericidal] [in [just] 2 minutes]
- [Bathroom disinfectant wipes]
- [Big [textured] virucidal* [,] [or] bactericidal wipes [,] [or] tuberculocidal [,] [or] [pathogenic fungi™] in 2-minute kill time]
- [Bleach-free cleaning -or- cleaner] [formula]
- [Bonus]
- [Broad spectrum disinfectant]
- [Built-in [tough] scrubbing fibers to help clean effectively]
- [[Product Name] has been tested for surface safety [,] [or] compatibility on common hard, non-porous surfaces [and equipment] to ensure compatibility with surfaces commonly used in hospital/ medical setting including: [acrylic] [brass] [polycarbonate] [titanium-coated] [medical grade] [stainless steel] [Formica®] [Glazed Porcelain] [Plastic] [Glazed Ceramic] [glass] [glass tables] [glazed tile] [Vinyl]]
- [[Product Name] has great compatibility with [medical devices] [dental equipment] [♦] [in ▲] when used to clean -or- pre-clean]
- [[Product Name] has great compatibility with non-critical, non-semi-critical [medical devices] [dental equipment] [♦] [in ▲] when used to disinfect]
- [[Product Name] kills germs that cause odors] [is effective against odor causing germs]
- [[Product Name] will effectively clean and disinfect when used as directed]
- [Pathogenic Fungi™] [in [just] 2 minutes]
- [Can help reduce cross contamination between treated hard non-porous surfaces]
- [Cleaner and disinfectant [towelette] [wipe]]
- [Cleaning Formula]
- [[Cleans] [and] [disinfects] [♦] [in ▲]]
- [[Cleans] [and] [disinfects] [on the go]]
- [[Cleans and] disinfects [hard, non-porous surfaces]] [on [high-touch] [frequently used] areas]
- [Cleans, deodorizes and disinfects]
- [Cleans] [disinfects] [Cleaner and disinfectant] [in one] [all in one [easy] step]
- [Clean Home]
- [Cleans the toughest soils]
- [Cleaning Wipe[s] [Towel] [s] refill[s][able]]
- [Cleaning Made Easy]
- [[Cleans [and /or Disinfects] [grocery cart handles or shopping cart handles or grocery cart child seats or shopping] cart child seats]
- [Cleans grocery carts and shopping carts]
- [[Cloth-like texture] [Textured] [strong] [for] [with] scrubbing power]
- [Compatible for use on hard, nonporous surfaces commonly found in [dental] [healthcare] [medical] [settings] such as [♦] [in ▲]]
- [Contains [number] canisters [plus] [and] [number refills]]
- [Contains (name of active ingredient), the [same] active ingredient used in [Brand name]. This product is not manufactured or distributed by © /or M/]
- [Convenient]
- [[Convenient] [Portable] way to [disinfect] [clean] hard, non-porous surfaces [of ♦] [in ▲]]
- [Cuts Cleaning time]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Deep cleaning [formula]]
- [Delivers a clean [▲] you can see and smell]
- [Demonstrates 99.999% -or- 5 [five] log reduction against for [*Candida auris*†] -or- [*C. auris*†] [in [just] [2] [two] minutes]]
- [Demonstrates 99.9% -or- 3 [three] log] reduction against SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces [in [just] [2] [two] minutes] [in one stepY] [in the presence of [5%] soil load]]
- [Designed for Home Healthcare -and/or- patients -and/or caregivers]
- [Designed for professional dental settings]
- [Designed -or- Developed for use in Dental [Offices -or- Labs - or Laboratories -or- Operatory -or- Operatories]
- [Designed -or- developed for Home Caregivers]
- [Designed -or- developed for use in [▲] in Home Healthcare -or- Home Care]
- [Delivers consistent [amount of solution] concentration every time [!]]
- [Deodorizes, and kills household germs]
- [Deodorizes by killing the germs that cause odors]
- [Deodorize hard, non-porous surfaces] [of ♦] [in ▲]]
- [Designed for the ▲ [industry]]
- [Designed for killing (*insert any organism(s) from Organism List*) [on hard, non-porous surfaces [of ♦] [in ▲]]]
- [Designed for killing SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces [of ♦] [in ▲] [in one stepY] [in the presence of [5%] soil load] [in [just] [2] [two] minutes]]
- [[Designed] [for] [with] [the] Environmental Services [in] [mind]!]]
- [Destroys [Deactivates] [[Helps] Eliminate 99.9%] (*insert any organism(s) from Organism List*) [of ♦] [in ▲]]
- [Destroys -or- Deactivates -or- Disinfects -or- [Helps] Eliminate 99.9% SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces [of ♦] [in ▲] [in one stepY] [in the presence of [5%] soil load] [in [just] [2] [two] minutes]]
- [Disinfectant]
- [Disinfect]
- [Disinfect external surfaces [of ♦] [in ▲]]
- [Disinfect [♦] surfaces]
- [Dispense one wipe at a time]
- [Disposable]
- [Disinfectant cleaner [that cleans, disinfects and deodorizes] in one [easy] stepY]
- [Disinfecting formula]
- [Disinfects] [♦] [in ▲]]
- [[Disinfects] [Disinfectant] [Cleaner] on the go]
- [Disinfects [and deodorizes] ♦ [in one [easy] stepY]]
- [Disinfecting towelettes]
- [Disinfects and deodorizes in [just] 2 minutes] (*insert any organism(s) from Organism List*)
- [Disinfects hard, non-porous surfaces in [just] 2 minutes] (*insert any organism(s) from Organism List*)
- [[Disinfects] hard, non-food contact, non-porous surfaces [that people [frequently] touch -or- use]]
- [Disinfects [in [just] two [2] minutes] (*insert any organism(s) from Organism List*)]
- [Disinfects [in [just] two [2] minutes]
- [Disinfects [bacteria] [viruses*] [[multi] drug-resistant [bacteria] [pathogens]] [bloodborne pathogens**] [ESKAPE pathogens***] [*Candida auris*† -or- *C. auris*†] [in [just] two [2] minutes]]
- [Disinfects [Tuberculosis -or- TB] [in [just] two [2] minutes]]
- [[Disinfects and deodorizes by killing germs and their odors]
- [[Disinfects hard, non-porous surfaces [throughout the ▲]]
- [[Disinfects [as it cleans] [and scrubs] in one [easy] stepY]]
- [Disinfect [the] surface] [in one [easy] stepY]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Does not contain dyes]
- [DO NOT FLUSH WIPE DOWN TOILET]
- [Double sided-scrubbing wipes]
- [Durable, lint-free wipe designed for sensitive equipment surfaces]
- [Easy]
- [Easy and convenient to use]
- [[Easy] [and convenient] way to clean [and] [disinfect]
- [Easy One-Step] Cleaning [and disinfecting] (except when precleaning is required for *Candida auris*† -or- *C. auris*†, visible soil, blood and body fluids)]
- [Easy, one-step] Cleaning [and] disinfecting against SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces]
- [Easy] pop-up-[tub] [canister] [pail] [lid] [tab] [dispenser] [bucket]]
- [Easy [and convenient] way to disinfect hard, non-porous surfaces [bathrooms] [office] [nursing stations] [high touch areas] [frequently used areas] everyday]
- [[Easy] [convenient] [dispensing] lid [closure] [cap] [package]]
- [Easy drop-in refills]
- [Easy to use]
- [Easy to handle]
- [Easy to refill]
- [Easy to use package [canister] [container] [refill] keeps [your] wipes moist]]
- [Easy way to [clean] [and] [disinfect] [bathroom] [or exterior toilet surfaces] [hard, non-porous surfaces]]
- [Easy and convenient disinfecting [throughout the ▲]]
- [Easy [and convenient] one-step] cleaning and disinfecting]
- [Easy [and convenient] way to [touch-up] [bathroom] [office] everyday]
- [Easy start [tab] [lid] [(insert package type)]
- [[Easy] [and convenient] way to [clean] [disinfect] [use]]
- [Easy way to clean [and disinfect] [bathroom] [or exterior toilet surfaces] [hard, non-porous surfaces]]
- [Easy to read label]
- [Economical]
- [Economy]
- [Economical] Refill[s] [one] [1] [two] [2] [three] [3] [four] [4] [various numbers] pack]]
- [[Economical] Refills [now] available]]
- [Effective -or- For daily use against (***Insert any organism(s) from Organism List***) [and] *Candida auris*† -or- *C. auris*†]
- [Effective -or- For daily use against Measles virus*, Avian Influenza A* and H1N1* [Pandemic] [Influenza A Virus*]]
- [Effective -or- For daily use against SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces [in one step] [in the presence of [5%] soil load]]
- [Effective against bathroom bacteria and viruses* [use at least two bacteria and two viruses* from this list:] [*E. coli* [O157:H7]] [*Salmonella enterica*] [*Staphylococcus aureus*] [Herpes Simplex Virus* [Type 2]] [Hepatitis B Virus* (HBV)] [Hepatitis C Virus* (HCV)] [Adenovirus* [Type II]] [Influenza A, H3N2 Virus*] [in [just] 2 minutes] [in the presence of organic soil [(5% blood serum)]
- [Effective against -or- Kills a range of bacteria including *Staphylococcus aureus*, *Salmonella enterica*, *Pseudomonas aeruginosa*]
- [Effectively disinfects [Effective] against (***Insert any organism(s) from Organism List***) [in [just] 2 minutes]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [[Effective against -or- Kills [in the presence of [organic soil load] [(5% blood serum)]]
- [[Effectively] disinfects [hard, non-porous] [environmental] surfaces]
- [Effectively removes everyday soils]
- [[Efficacy tests have demonstrated that] [Product Name] [this product] [are] [is] an effective bactericide, virucide*, tuberculocide in the presence of organic soil load [5% blood serum]]
- [Eliminates odors]
- [Eliminates odors from fecal matter and urine]
- [Eliminates airborne particles associated with sprays]
- [Eliminates [daily] mixing and measuring]
- [Eliminates [lingering] [food] odors on hard, non-porous surfaces]
- [[Eliminates -or- Reduces odors [at their source] [caused by germs -or- bacteria-and/or- mildew]]
- [EVS [Environmental Services] [use]]
- [Exhibits effective [bactericidal] [virucidal*] [tuberculocidal] activity]
- [Extra large] [disinfecting] [Y" x Z" compared to our A" x B" size] [wipes] [Disinfecting Towelettes]
- [Kill[s] -and/or- Effective against [99.9%] (*insert any organism(s) from Organism List*) [on environmental -and/or- hard, non-porous surfaces]]
- [Kill[s] -and/or- Effective against [99.9%] cold -and/or- flu [virus(es)*, Avian flu [virus*] in [just] 2 minutes]
- [Help[s] reduce[s] the spread of [99.9% of] [illness-causing] [clinically relevant] [bacteria] [pathogens] [organisms] -and/or- viruses* -and/or- germs between treated hard, non-porous surfaces [in your ▲] [in this cold -and/or- flu season]]
- [Kill[s] -and/or- Help[s] reduce[s] the spread of -and/or- Effective against SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces] [in [just] [2] [two] minutes] [in one step] [in the presence of [5%] soil load]]
- [For [daily] [general] use]
- [For [daily] [general] [professional] use [especially] [during this [COVID-19] pandemic -or- outbreak]]
- [For more information visit us at www.metrex.com -or- www.kerrdental.com]
- [For use on bathroom surfaces]
- [For use on nursery surfaces]
- [For use [in ▲] [on ♦ surfaces]]
- [For use in critical care areas where reduction of cross contamination between hard non-porous surfaces is of prime importance]]
- [For surface cleaning -and/or- deodorizing -and/or- disinfecting]
- [For Use by Home [Caregivers] [Healthcare]]
- [Formulated -or- Designed to remove[s] biological odors including Urine -and/or- Feces -and/or- Vomit]
- [Formulated to kill (*insert any organism(s) from Organism List*) [on hard, non-porous surfaces] [when used according to instructions]]
- [Fragrance Free] [Free of fragrance]
- [Freshens]
- [[Fungicidal™ disinfectant] against *C. auris*‡]
- [Furniture Safe]
- [Germs are no match with [Product Name]]
- [Germicidal disposable wipes]
- [Germicide -or- Germicidal]
- [Gets rid of [eliminates] odors]
- [Great for [everyday] [daily] [general] use]
- [[Heavy] cloth-like towel]
- [[Helps] [kill][s] [eliminate][s] odor-causing bacteria]
- [Helps prevent cross contamination between treated hard non-porous surfaces]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Improved canister and lid]
- [Improved dispensing system]
- [Includes ___count refill [pack] [contains]___wipes]
- [Institutional [size]]
- [[Just] [Simply] wipe [exterior toilet surface] [to Kills -or- eliminates 99.9% bathroom germs]]
- [Just [pull] wipe and toss]
- [Kills -and/or- inactivates [(**Insert any organism(s) with from Organism List**)] [in [just] 2 minutes] [on hard, non-porous surfaces]]
- [Kills Pandemic 2009 H1N1 Influenza A Virus* [(formerly called swine flu)]] [in [just] 2 minutes]
- [Kills [the] [Pandemic] SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces [in [just] [2] [two] minutes] [in one step∨] [in the presence of [5%] soil load]]
- [Kills [bacteria] [viruses*] [germs] [pathogenic fungi™] [*Candida auris*† -or- *C. auris*†] [mold] [mildew] [clinically relevant organisms [pathogens]] [on hard, non-porous surfaces]]
- [Kill germs and bacteria instead of paper towels, sponges and rags [Use this product to]]
- [Kills -or- Destroys [mold and mildew] [odor causing bacteria]]
- [Kills -or- Disinfects [Human] Norovirus* [Feline Calicivirus] [Measles virus]] [in [just] 2 [two] minutes]
- [Kills -or- eliminates 99.9% -or- Effective against [bacteria] [germs] [viruses*] [pathogenic fungi™] [*Candida auris*† -or- *C. auris*†] [clinically relevant organisms [pathogens]]
- [Kills -or- Disinfects -or- Eliminates -or- [Is] Effective against [99.999%] *Candida auris*† -or- *C. auris*† [in [just] 2 min[utes]] [on hard, non-porous surfaces]]
- [Kills -or- Disinfects -or- Effective against [bacteria] [and] [germs] [and] [viruses*] [and] [multi] drug resistant bacteria [pathogens] [and] [ESKAPE pathogens***] [and] [pathogenic fungi™] [and] [*Candida auris*† -or- *C. auris*†] [and] [clinically relevant organisms [pathogens]] [in [just] [2] [two] minutes] [on hard, non-porous surfaces]]
- [Kills -or- Disinfects -or- Eliminates -or- [Is] Effective against [Feline Calicivirus] [[Human] Norovirus*] [Measles virus] [Tuberculosis -or- TB] [(**Insert any organism(s) from Organism List**)] [in [just] 2 minutes] [on hard, non-porous surfaces]]
- [Kill [Effective against] bacteria in 2 minutes] [Kill [Effective against] viruses* in 2 minutes] [Kill [Effective against] HIV-1 in 2 minutes] [Kill bloodborne pathogens** in 2 minutes] [Kill ESKAPE pathogens*** in 2 minutes] [Kill [Effective against] SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces in 2 minutes]
- [[Kills] [Effective against] Tuberculosis -or- TB in [just] 2 minutes]
- [Kill [Effective against] [pathogenic fungi™] [*Candida auris*† -or- *C. auris*†] in [just] 2 minutes]
- [Kills -or- Destroys -or- Disinfects -or- Eliminates 99.9%-or- Effective against SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces [in [just] [2] [two] minutes] [in one step∨] [in the presence of [5%] soil load]]
- [Kitchen+ and Bath Wipes]
- [Large[r] Surface[s]]
- [Large[r] [Wipe][s]]
- [Larger wipe covers [more surface area] larger area [with less effort] [!]]
- [Lead Free]
- [Leaves [bathroom] [office] [▲] [♦] fresh and clean [and disinfected] [and smelling fresh] in [just] minutes]
- [Less] [Reduced] [packaging] [less] [to throw away]
- [Lightly scented -or- quick-dissipating scent]
- [Low linting [Lint free] [No lint] [wipe]]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Made for Healthcare Workers]
- [Made to reduce the risk of cross contamination between treated hard non-porous surfaces]
- [Making home healthcare [manageable] [easy] [simple]]
- [Maximum cleaning, minimum effort]
- [Meets [the disinfection requirements of] OSHA['s] Bloodborne Pathogen Guidelines -or- Standards]
- [Mercury free]
- [Mildewcide -or- Mildewcidal]
- [Moist wipes]
- [Moisture lock lid]
- [Multi-pack]
- [Multi-purpose disinfectant]
- [Multi [-refill] [pack]]
- [Multi-surface]
- [Next wipe pops up automatically [everytime] [!]]
- [Neutralizes odors [doesn't [just] mask them]]
- [No canister needed]
- [No [pinched fingers] [rips] [tears]]
- [No sponges, no rags, no buckets, no mess]
- [No threading [required]]
- [No fragrance added]
- [No measuring] [No mixing] [required]
- [Non-scratching [wipe]]
- [Non-woven disposable cloth]
- [Nonylphenol ethoxylate (NPE) -and/or- alkylphenol ethoxylate (APE) Free]
- [[Occupied Room] [and] [Terminal Room] [Discharge] Cleaning -and/or- Disinfecting Wipe[s] [Towel][s] Refill[s][able]]
- [Odor [elimination] [neutralization]
- [On the go disinfectant]
- [On -or- Included in EPA's List N [Disinfectants for Use Against SARS-CoV-2 (COVID-19)]]
- [One step∇ [cleaning] [disinfecting] [wipe]] [in one easy step∇]
- [One step∇ cleaning and disinfecting [saves time]]
- [One step∇ cleaning [and/or disinfecting] [and/or deodorizing] disinfect and/or touch-up clean] [cleans and/or disinfects and/or deodorizes] [anytime] [daily] [every day] [in one easy step∇] [on hard, non-porous] [on ▲] [bathroom and/or kitchen+] surfaces [on household surfaces] [on ♦ surfaces]]
- [One step∇ disinfectant cleaner [that is] effective against a broad spectrum of bacteria -and/or- viruses* when used as directed]
- [One step∇ [germicide] [bactericide] [virucide*] [against Pathogenic Fungi[™]]]
- [One step∇ disinfectant cleaner designed for [general] cleaning and disinfecting on hard, non-porous environmental surfaces in ▲]
- [On-the-go disinfectant [Package] [insert type of package] made from [x%] recycled materials]
- [One [wipe] at a time dispensing]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Package] [Insert type of package] made from [x%] recycled materials]
- [[Passes] [Meets] [Satisfies] requirements for killing *Candida auris*† -or- *C. auris*† [on hard, non-porous surfaces]]
- [Passes -or- Meets -or Satisfies [efficacy] [testing] requirements for killing SARS-CoV-2* -or- COVID-19 Virus* on hard, nonporous surfaces [in the presence of [5%] soil load] [in [2] [two] minutes]
- [Picks up [particles] [dust] [dirt] [hair] [crumbs]]
- [Pop-up [dispensing] [top] [tab] [lid]]
- [Portable]
- [Pouch[es]]
- [Powers through tough messes]
- [Pull n/ Go [pack]]
- [Pull, Clean, and Disinfect [in] [only] [just] [2 minutes]] [as a [virucide] [bactericide] [Tuberculocide] [against pathogenic fungi™]]
- [Premixed], [premeasured], [premoistened] [towelettes]
- [Premoistened [wet] household cleaning wipe]
- [Pre-measured to have the recommended amount of germicidal solution, wipe after wipe]
- [Pre-moistened [Easy] [and convenient] way to disinfect]
- [Pre-moistened with specially formulated cleaners to handle daily clean-ups without bottles, sponges or buckets]
- [Prevents formulation errors]
- [[Product Name] are pre-moistened with a stabilized [Product Name] [disinfectant] solution and detergent solution and are ready- to-use.]
- [[Product Name] are useful in a variety of [healthcare], [medical], [dental], [nursery/daycare] [and] [▲] applications listed on the [this] label.]
- [Pseudomonacidal^] [in [just] 2 minutes]
- [Ready-to-use disinfectant [disinfection] for [killing] *Candida auris*† -or- *C. auris*†]
- [Ready-to-use [one-step √] disinfectant [disinfection] for [killing] SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces]
- [[Ready-to-use] disinfectant [for use in ▲]
- [Recommended for cleaning [and disinfecting] ◆]
- [Recommended for disinfecting exterior toilet surfaces -or- around the exterior toilet surface]
- [Recyclable] [package]
- [Reduces [packaging] waste [by X%] [when compared to regular canister or bucket or pail]]
- [Reduces [risk of] cross-contamination [from (*insert any organism(s) from Organism List*)] [in [just] 2 minutes] between treated hard non-porous surfaces]
- [Renders (*insert any organism(s) from Organism List*) inactive [on hard, non-porous surfaces]] [in [just] 2 minutes]
- [Refill package contains [one] [1] [two] [2] [three] [3] [four] [4] [various numbers] [wet] [pre- moistened] [wipe] [refill(s)] [roll(s)]]
- [Refill [size] [pack] [value pack]]
- [Refillable [with ____# wipes] [dispenser] refill [pack] [for dispenser]]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Requires no mixing or measuring and are ready-to-use]
- [Re-sealable [stay-moist] pouch [package] [container] inside]
- [Re-use your [empty] canisters [again [and again] [!]]]
- [Reuse[able]] [Container] [Dispenser] [Canister] [Pail] [Bucket]
- [Sample]
- [Scrubbing power of a sponge in a disinfecting [towelettes]]
- [Shelf-stable] [formula]
- [So] [surfaces are more than [just] [clean] [they're] [disinfected]]
- [Specifically -or- Specially designed for the Medical/ Dental/ Veterinary/ Health Care [▲] industry]
- [Specifically -or- Specially designed [for] [with] [dentists] [hygienists] [patients] [in mind]]
- [Staphylocidal^^] [in [just] 2 minutes]
- [Starts cleaning on contact]
- [Strong, durable, low-lint [lint free] [cloth-like] wipe [designed for sensitive equipment surfaces]]
- [Strong, durable textured wipe [cleans] [and] [deodorizes] [and] [disinfects]]
- [[Strong,] [Durable,] non-woven disposable cloth [wipe]] [does not leave lint on surfaces]
- [[Suitable] [Safe] for [▲]]
- [[Suitable] [Safe] for use on [◆]]
- [Suitable for use on hard, non-porous surfaces such as [acrylic] [brass] [polycarbonate] [titanium-coated] [medical grade] [stainless steel] [Formica®] [Glazed Porcelain] [Plastic] [Glazed Ceramic] [glass] [glass tables] [glazed tile] [Vinyl]]
- [[Suitable] [Safe] for use on hard, non-porous inanimate environmental surfaces]
- [Strong cleaning power]
- [[Strong] [Sturdy] canister lid's hinge]
- [Sturdy] [Container] [Dispenser] [Canister] [Pail] [Bucket]
- [Tear at perforation]
- [Textured for better [cleaning power] [scrubbing power] [disinfecting]]
- [Textured design allows for] [scrubbing] [cleaning] [removing] [wiping] [tough] [stains] [soils] [messes]
- [Textured wipe for effective cleaning]
- [[Textured] wipe is low-linting [lint-free] [designed] [for sensitive equipment surfaces]]
- [[Textured design [cleans] [deodorizes] hard, non-porous surfaces] [helps remove soil from surfaces]]
- [Textured] [strong] [towelettes] [for] [with] [scrubbing power]
- [The 2-minute [germicidal wipe] [disinfectant wipe] [bactericidal wipe] [virucidal* wipe]]
- [The 2-minute wipe against *Candida auris*† -or- *C. auris*† (***insert any organism(s) from Organism List***)]
- [The answer to your [disinfecting] [cleaning] needs]
- [[The] smell of clean]
- [[The] easy -and/or- convenient way to disinfect] [with [Product Name]]
- [[This product] reduces cross-contamination [of germs] between hard, non-porous treated surfaces]
- [This product has demonstrated effectiveness against Influenza A Virus* and is expected to inactivate all Influenza A Viruses* including Pandemic 2009 H1N1 Influenza A Virus*] [in [just] 2 minutes]
- [This product has demonstrated effectiveness against SARS-CoV-2* -or- COVID-19 Virus* in 2 minutes on hard, non-porous surfaces [in the presence of [5%] soil load] [in one-step √]]
- [This product meets AOAC efficacy testing requirements -or standards for hospital disinfection]
- [[Thick], [durable] [and] [heavy duty] wipes [that are] [bactericidal -and/or- virucidal* -and/or- in 2 minutes] [tuberculocidal in 2 minutes]]
- [Thick, large wipe [that] [cleans] [disinfects] [deodorizes]]
- [[This] [outer] [selling] [package] [container] [canister] contains X% post-consumer [material] [plastic] [resin]]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Time to reorder]
- [Trial] [Pack] [Canister]
- [Tuberculocidal] [in [just] 2 minutes]
- [Tub with pop-up action]
- [To break down -and/or kill -and/or- inactivate (*insert any organism(s) from Organism List*)] [on surfaces] [on hard, non-porous surfaces].
- [To kill (*insert any organism(s) from Organism List*)], let stand [undisturbed] for 2 min[utes]
- [To kill *Candida auris*† -or- *C. auris*† [,] on treated surfaces especially during outbreaks] [of *Candida auris*† -or- *C. auris*†] [let stand [undisturbed] for two -or- 2 min[utes]]
- [To kill SARS-CoV-2* -or- COVID-19 Virus* on treated hard, non-porous surfaces [especially during this [COVID-19] outbreak] [,] [let stand [undisturbed] for two -or- 2 min[utes]]
- [To reduce *Candida auris*† -or- *C. auris*† on hard, non-porous treated surfaces [of ♦] [in ▲]]
- [To reduce [cross-contamination of] SARS-CoV-2* -or- COVID-19 Virus* between treated, hard, non-porous surfaces [of ♦] [in ▲]]
- [To reduce cross-contamination between treated, hard non-porous surfaces of dental impression trays for both dentist office personnel and dental laboratory personnel.] (**Note to Reviewer: Dental impression trays do not go into the human mouth.**)
- [[Tub] [Pail] [Bucket] refill] [available]]
- [Twin pack]
- [Two in one action] [cleans and disinfects in one easy step∇]
- [Tuberculocide -or- Tuberculocidal] [in [just] 2 minutes]
- [Usable on variety of soils]
- [Use in [public] [common places] [▲] where [bacteria -and/or- viruses*-and/or- germs may be of concern]]
- [Use to] [Clean -and/or Disinfect -and/or- Deodorize hard, non-porous surfaces]
- [[Use to] refill[s] [count] [wipe] [canisters] [tubs] [bucket] [tub] [pail]]
- [Use to pre-clean critical or semi-critical medical devices prior to sterilization / high-level disinfection]
- [Use wipes to kill germs -and/or- bacteria -and/or- and viruses* instead of paper towels, sponges and rags]
- [Use this product to] [Kill bacteria in 2 minutes] [Kill viruses* in 2 minutes] [Kill TB in 2 minutes] [Kill HIV-1 in 2 minutes] [Kill *Candida auris*† -or- *C. auris*† in 2 minutes] [Kill clinically relevant organisms [pathogens]] [Kill ESKAPE pathogens*** in 2 minutes] [Kills bloodborne pathogens** in 2 minutes] [Kill SARS-CoV-2* -or- COVID-19 Virus* in 2 minutes on hard, non-porous surfaces]
- [Use where reduction of cross-contamination between treated hard non-porous surfaces is of prime importance]
- [Use on hard, non-porous surfaces only: [acrylic] [brass] [polycarbonate] [titanium-coated] [medical grade] [stainless steel] [Formica®] [Glazed Porcelain] [Plastic] [Glazed Ceramic] [glass] [glass tables] [glazed tile] [Vinyl]] [Test wipe on small inconspicuous area first]
- [Value]
- [Versatile]
- [Virtually [Lint Free] [residue free] [streak free] [clean] [shine]]
- [Virucide* -or- Virucidal*] [in [just] 2 minutes]
- [Virucide* -or- Virucidal* against SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces [in [just] 2 minutes] [in the presence of [5%] soil load] [in one-step ∇]]
- [VOC free]
- [Wet wipes]
- [When used according to label directions, [this product] reduces cross-contamination of (*insert any organism(s) from Organism List*) between treated hard non-porous surfaces]
- [Where control of cross contamination between treated hard non-porous surfaces is of prime importance]
- [Quick Cleaning action -or- activity]
- [Wipe away odors]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

PRODUCT CLAIMS continued

- [Wipes out 99.9%] [Kills] [germs] [bacteria] [viruses*] [while it cleans] [that cause odor]
- [[Wipe] [material] made from [x %] [sustainable] [ingredients] [fibers]] [may only be used when cellulose wipe material used]
- [Wipes away -or- out 99.9%] [Kills] [Eliminates 99.9%] [household] [bacteria] [viruses*] [pathogenic fungi™] [*Candida auris*† -or- *C. auris*†] [clinically relevant organisms [pathogens]] [germs] [(*insert any organism(s) from Organism List*)] [in [just] 2 minutes]
- [Wipes away -or- out 99.9%] [Kills] [Eliminates 99.9%] [household] [(*insert any organism(s) from Organism List*)] [in [just] 2 minutes]
- [Wipes away -or- out 99.9%] [Kills] [Eliminates 99.9%] [viruses* and bacteria [including (*insert any organism(s) from Organism List*)] [in [just] 2 minutes]
- [[Wipes away -or- out 99.9%] [Kills] SARS-CoV-2* -or- COVID-19 Virus* on hard, non-porous surfaces [in one step]] [in the presence of [5%] soil load] [in [just] 2 minutes]]
- [X" x Y" Wipes]
- [X% [more] [wipes] [free] [vs. (*insert a Metrex disinfectant wipes*) count] [canisters] [packs]]

[Note: The following claims can be used only for the first 6 months of the product on shelf]

- [New] [Strong] [Thick] Wipe [!]
- [New] [&] [dispensing] lid [closure] [cap] [package] [dispensing] [better [than before]]
- [New]
- [New] [Product Name] Same [trusted] [great] [reliable] brand
- [New] [Package] [Packaging]
- [New thick and larger embossed wipe cleans more surface area]
- [Now] [Reinforced] [With] [Quilted] [strength]
- [Now with [%] more towel]
- [Now with [#] more wipes]
- [Scrubs] (&) [cleans]
- [New Strong durable [wipe]
- [New Thick [wipe] [with] cloth-like texture]
- [New Large] [quilted] wipes]
- [New Thick wipes]
- [Multi-purpose [cleaner] [disinfectant] [Multi [-refill] pack] New [large] thick [textured] wipe]]
- [New] [Big] [Strong] [Thick] [wipe] [!]
- [New] [dispensing] [lid] [closure] [cap] [package] [dispenses]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Alternate Brand Names

- CaviSurf
- CaviQuat
- CaviWipes 3D
- CaviWipes DT
- CaviWipes JS
- CaviWipes Blue
- CaviWipes +
- CaviWipes Plus
- CaviWipes 360
- CaviWipes Pur
- CaviWipes U
- CaviWipes Gold
- Metrewipe 2M
- Cavi[Wipes] X2
- Cavi[Wipes] IP
- CaviWipes Medical 2M
- CaviWipes Medical M2
- CaviWipes M2
- CaviWipes Formula 7
- CaviWipes LA
- CaviWipes Z
- CaviWipes Xtra
- CaviWipes Avant
- CaviWipes Alcazar
- Cavi Avant
- Cavi Alcazar
- CaviWipes No. 7
- Cavi-7
- CaviWipes Number 7
- CaviWipes #7
- Metrewipe Medical 2M
- Cavi[Wipes] NexGen
- Cavi[Wipes] 2.0
- Cavi[Wipes] II
- CaviWipes Medical LA
- CaviWipes V2

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

OPTIONAL GRAPHICS

[The following graphics are optional for canister containers only. They will be affixed to the label/canister such that they do not interfere with any label text. Due to size constraints, they will not appear on individually wrapped packets.]



[The following graphics are optional for all containers. They will be affixed to the label such that they do not interfere with any label text.]



THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

[Note to Reviewer: following section is optional.]

[CLEANING AND STAIN REMOVAL]

[To clean -and/or- remove stains: Wipe surface with a wipe; [rinse and] let air dry. [If streaking is observed, wipe with a clean, damp [cloth] [or] [paper towel].] **[[Works on] [effective against] [removes the following stains:]]]**

<i>Claims:</i> <ul style="list-style-type: none">• Cleans• Removes stains• Removes (<i>insert stain(s) / soil(s) shown to the right:</i>)	<i>Stains/Soils</i> <ul style="list-style-type: none">• Blood• Body oils• Bodily fluids• Coffee [stains]• Common soils• Dead skin• Defecation• Dirt• Fecal matter• Feces• Fingerprints• Food+ residue[s] -or- soil[s]• Food [stains]• Fruit [stains]• Grease• Grime• Laboratory stains	<i>Stains/Soils continued</i> <ul style="list-style-type: none">• Mildew stains• Mold stains• Mucus• [Other] common soils -and / or- stains• [Other] organic matter• Phlegm• Protein-based matter• Saliva• Spit• Spittle• Sputum• Stool• Throw-up• Tea [stains]• Urine• Vomit
--	---	---

[Note to Reviewer: following section is optional.]

[ALLERGEN DESTRUCTION]

[To [clean and] destroy specified allergens: Wipe surface. Use enough wipes to treat surfaces for 2 minutes. [Use 2-minute contact time [dwell time].] [Rinse.] Let air dry. [If streaking is observed, wipe with a clean, damp [cloth or] paper towel after contact time [dwell time] has expired.]

[Claims:]

- [Breaks Down Allergens[!]]
- [Destroys -or- deactivates dust mite matter]
- [[Penetrates and] breaks down dust mite matter -or- particles, cockroach matter -or- particles, and pet dander allergens [found in dust]]

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 1 Medical Use Sites:
[Product Name] are recommended for use on hard, non-porous surfaces in [▲]

3D Imaging Service Centers	Endoscopy Centers	Ophthalmic Facilities -or- Offices
Acute Care Facilities	ENT [(Ear, Nose, Throat)] Facility	Optometric Facilities -or- Offices
Alternate Care Facilities	Eye Care [Facilities] [Centers]	Oral Health Care Facilities
Ambulances	Eye Surgical Centers	Orthodontic Offices
Ambulatory Surgery Centers (ASC)	Forensic Labs	Orthopedic Offices
Animal Research Facilities	Full-Care Nursing Facility	Outpatient Care Facilities
Assisted Living Facility	Gastrointestinal [(GI)] [Labs] [Laboratories]	Outpatient [Surgery] [Facilities] [Centers]
Autopsy Rooms	Gastrointestinal [(GI)] [Centers] [Suites]	Out-Patient Surgical Centers (OPSC)
Bathroom	Hallways	Pathology Rooms
Bioskill Laboratories	Home Care Services	Patient Care Areas
Blood Drawing Centers	Hospice	Patient Care Rooms
Blood Apheresis Facilities	Hospital [Operating Rooms] [Kitchens+]	Patient Rooms
Blood Units	Hydrotherapy [Facilities] [Rooms]	Pediatric Exam Rooms
Burn Center	<i>In vitro</i> Fertilizer Center (IVF)	Pediatric Intensive Care Units (PICU)
Cardiac Intensive Care Unit (CICU)	Instrument Reprocessing [Facilities] [Departments] [Center]	Pediatric Offices
CAT Scan [Lab] [Laboratories]	Intensive Care Units (ICU)	Personal Care Service Facilities
Cath Labs	Isolation Areas	Pharmacy Clean Rooms
Central Service Areas	Isolation Rooms	Phototherapy Rooms
Central Services Department	Labor Units	Physical Therapy Facilities
Central Supply Areas	[Lab] [Laboratories]	Physicians' Offices
Central Supply Rooms	Lasik Eye Centers	Plasma Apheresis Centers -or- Facilities -or- Units
Clean Rooms	Lasik Units	Post [Ops] [Operatory] Rooms
Clinical Areas	Laundry Rooms	Pre [Ops] [Operatory] Areas
Clinics	Life Science Laboratories	Psychiatric Care Facility
Critical Care Units	Long-Term Acute Center (LTAC)	Radiology Department
CS [(Blood Culture)] [Labs] [Laboratories]	Long-Term Care [Facilities]	Recovery Rooms
Dialysis Centers	Medical [Facilities] [Clinics] [Sites] [Areas]	Rehabilitation Facilities
Dialysis Clinics	Medical Related Facilities	Respiratory Therapy Facilities
Dissection [Labs] [Laboratories]	Medivans	Skilled Nursing Facilities -or- Homes
Doctor's Offices	Medical Intensive Care Unit (MICU)	Environmental Surfaces of Sterile Processing Departments
Donation Centers [blood] [plasma] [semen] [milk] [aphaeresis]	Magnetic Resonance Imaging [(MRI)] [facilities] [laboratories] [labs]	Environmental Surfaces of Sterile Storage Areas

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Drug Preparation Centers -or- Facilities -or- Rooms	Neonatal Intensive Care Units (NICU)	Surgery Centers Surgery Rooms
Electrolysis Facilities	Neonatal Nursery Areas	Surgical Centers
Emergency Departments	Neonatal Units	Surgical Suites
Emergency Medical Services (EMS) Units	Newborn Nursery	Teaching Hospitals
Emergency Medical Settings	Nuclear Medicine Facilities	Tissue Banks
Emergency Rooms (ER)	Nurses' Stations	Transplant Center
Emergency Vehicles	Nursing [Homes] [Stations]	Trauma Center
Endo Suites	Occupational Therapy [Centers] [Facilities]	University Research Labs
Endoscopy [Centers] [Facilities] [Departments]	Oncology Areas	Waiting Rooms
X-Ray Rooms	Operating Rooms	Wound Center
Extended Care Facilities	Skilled Nursing Facilities	

Table 2 Institutional Use Sites:

[Product Name] are recommended for use on hard, non-porous surfaces in [▲]

Adult School	Early Learning Centers	Middle Schools
Air Transport Facilities	Elementary Schools	Military Installations
Airplane Interiors	[Fire] [Rescue Departments]	Nursery Schools
Airports	Fire Vehicles	Police Vehicles
Bathroom	Forestry Services	Pre-K School [Room]
Casinos	Grammar School	Preparatory School
Child Care Facilities	High School	Preschool
Churches	Hospitality Facilities	Primary School
Colleges	Hotels	Prisoner Detention Facilities
Community Colleges	Industrial Facilities	Prisoner Lockups
Conference [Centers] [Rooms]	Infant Care Centers	Prisons
Conference Rooms Fixtures	Intermediate School	School Buses
Convention Facilities	Jails	School Kitchens+
Correctional Facilities	Junior High School	Schools
Daycare Centers	Kindergarten	[University] [Universities]
[Dormitories] [Dorm]	Learning Institutions	

Table 3 Dental Use Sites:

[Product Name] are recommended for use on hard, non-porous surfaces in [▲]

Dental Laboratories -or- Labs	Dental Office	Dental Operatory -or- Operatories
-------------------------------	---------------	-----------------------------------

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 4 Vet Use Sites:

[Product Name] are recommended for use on hard, non-porous surfaces in [▲]

Animal Care Facilities	Avian Units	Pet Homes
Animal Control Vehicles	Bathroom	Pet Motels
Animal Dental Operatory	Doggy Day Care Centers	Pet Spas
Animal Detention Room	Feline Facilities	Pet Stores
Animal Research Facilities	Feline Units	Small Animal Facilities
Animal Shelters	Kennel	Veterinary Offices
Avian [Care] Facilities	Pet Areas	Veterinary Facilities

Table 5 Non-Food Surfaces of Food Prep+ Use Sites:

[Product Name] are recommended for use on hard, non-porous surfaces in [▲]

Bakeries+	Food Processing Facilities+	Restaurants+
Banquet Hall Facilities+	Grocery Stores +	School Kitchens+
Butcher Shops+	Kitchen+	Supermarkets+
Food Preparation Areas+	Produce Areas+	

+ For non-food contact surfaces only.

Table 6 Miscellaneous Use Sites:

[Product Name] are recommended for use on hard, non-porous surfaces in [▲]

Airplanes	Facilities	Manufacturing [Facilities] [Plant]
Arboretums	Factories	Money Handling Facilities
Arcades	Foots Spas	Morgues
Athletic -or- Training Facilities	Foreign Current Exchange Facilities	Mortuaries
Auto Facilities	Funeral Homes	Motels
Banks	Gaming Facilities	Movie Theaters
Barber Shops	Garden Centers	Outdoor Playground
Bathroom	General Public Facilities	Personal Care Service Facilities
Beauty Salons	Golf Courses, not for use on grass	Public Areas
Body Piercing Salons	Grounds Keeping Services	Public Facilities
Body Service Facilities	Gymnasium [Gym]	Public Transportation Vehicles
Body Storage Facilities	Health Club Facilities	Public Restrooms
Botanical Gardens, not for use on plants/ grass	Helicopter Interiors	Recreational Centers -or- Facilities
Bowling Alleys	Household Bathrooms	Recycling Centers
Buses	Housekeeping Rooms	Salon Settings
Cabs	Indoor Playground	Shelters
Cadaver Processing areas	Laundromats	Shower Rooms
Cars	Laundry Rooms	Skincare Salons
Centers	Lawn Services	Storage Areas -or- Rooms
Condominiums	Limousine	Subway
Cruise Ships	Living Rooms	Tanning Salons
Currency Exchange	Locker Rooms	Tattoo Parlors
Day Spa	Lockup Facilities	Taxis
Delivery Units	Manicure Salons	Theme Parks
Dissection Work Tables	Massage Facility	Transport Vehicles
Universities	Trade Schools	

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces:

[Product Name] is recommended for use on hard, non-porous surfaces of [♦].

Ambulance -and/or- Ambulance Equipment Surfaces	Dental impressions [materials] [prior to sending to a dental lab]	
Animal Cages	Cages	Dental Operatory Surfaces
Animal Equipment	Cardiac Gym Equipment	Dental Unit Instrument Trays
Non-Food Contact Surfaces of Appliances+	Carts	Dentist Chair
Arm Rests	Cash Register	Desk Tops
ATM Surfaces	Ceilings	Desks
Autoclave Surfaces	Cell Phones	Diaper Changing Stations
Automatic Feeders	Chairs	Diaper Pails
Automatic Teller Machines (ATM)	Changing Tables	Digital Sensors
Baby Warmers	Chemotherapy Hoods	Digital Phosphor Plates
Baked Enamel	Child Seats+	Non-Food Contact Surfaces of Tables+
Bars	Chrome	Non-Food Contact Drying Racks+
Basins	Clipper Blades	Display Cabinets
Baskets	Closet Handles	Diving Equipment
Bathing Units	Coated Mattresses	Diving Masks
Bathroom Fixtures	Coated Pillows	Diving Snorkels
Bathtubs	Combs	Door Knobs
Bed Pans	Computer	Dopplers
Bed Railings [Bedside Rails]	Computer Display	Non-Food Contact Drain Boards+
Bedside Commode	Computer [Keyboards] [Mouse]	Drained Whirlpool® Tank Surfaces
Behind and Under Sinks	Copy Machines	Enamel
Bidet	Non-Food Contact Counter Tops+	Endoscope Cabinet
Biohazard Container	Non-Food Contact Counters	Equipment
Blood Pressure Cuff	[Baby] Cribs	Equipment Surfaces
Booster Chairs+	Crutches	Examination -or- Exam Table
Booster Seats+	Dental Chairs	External Surfaces of Kendrick Extrication Device (KED)
Brushes	Dental Countertops	Exterior Surfaces of Toilet
Cabinets Handles	Dental Curing Lights¥	Exterior or Exposed Interior Surfaces of Medical Equipment
Cabinets	Dental Impression Trays	Bassinets

- **+ For non-food contact surfaces only.**
- ¥ For surfaces that are not inserted into the human mouth and body.
- µ Allow surfaces to come to room temperature prior to treatment.
- ***For use on Hard, Non-Porous [Non-Food Contact Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [▲ surfaces].***
- ***Any hard, non-porous, non-food contact surface where disinfection is required+.***

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces continued:

Exterior Surfaces of Dental Equipment [‡]	Exterior Surfaces of Non-Food Contact Vanity Tops ⁺	External Surfaces of Intubation Equipment
Exterior Surfaces of Amalgamators [‡]	External Surfaces of Emergency Resuscitative Equipment	External Surfaces of Isolettes
Exterior Surfaces of Anesthesia Equipment -or- Machines [‡]	External Surfaces of Bladder Scan[ners]	External Surfaces of IV Pumps
Exterior Surfaces of Aphaeresis [Machines]	External Surfaces of Blood Glucose Meters	External Surfaces of Life Support Equipment
Exterior Surfaces of Appliances such as Refrigerators ⁺ , Microwaves ⁺ , Food Processors ⁺ , Stoves ⁺	External Surfaces of Blood Pressure (BP) Monitors	External Surfaces of Loupes
Exterior Surfaces of Center Venous Lines	External Surfaces of CAT Scan Machines	External Surfaces of Medical Equipment
Exterior Surfaces of CICU Equipment	External Surfaces of Compression Devices	External Surfaces of Medical Warmers
Exterior Surfaces of Defibrillators	External Surfaces of Dialogistic Equipment	External Surfaces of Monitors Electrodes
Exterior Surfaces of Dental Care Facilities Hand Pieces	External Surfaces of Dictating Equipment	External Surfaces of MRI Machines
Exterior Surfaces of Dental Equipment Lines	External Surfaces of Digital X-Ray Sensors	External Surfaces of OLED, LED Monitors
Exterior Surfaces of Dental Equipment Tubing	External Surfaces of EKG (Electrocardiogram) Machines	External Surfaces of Ophthalmoscopes
Exterior Surfaces of Dental Operatories	External Surfaces of Electronic Thermometers	External Surfaces of Oscopes
Exterior Surfaces of Diagnostic Equipment	External Surfaces of Electrosurgical Generators	External Surfaces of Patient Care Monitor [Unit]
Exterior Surfaces of Dialysis Machines	External Surfaces of Equipment	External Surfaces of Patient Monitoring Equipment
Exterior Surfaces of Endodontic Equipment (Apex Locators, Pulp Testers, Motors)	External Surfaces of Glucometers	External Surfaces of Probes
Exterior Surfaces of Radiology Equipment	External Surfaces of Headsets	External Surfaces of Pulp Tester Motors
Exterior Surfaces of Refrigerator Units	External Surfaces of ICU Equipment	External Surfaces of Pulse Oximeter
Exterior Surfaces of Respiratory Therapy Equipment	External Surfaces of Incubators	External Surfaces of Sequential Compression Devices
Exterior Surfaces of Urinals	External Surfaces of Infusion Pumps	External Surfaces of Respiratory Ventilators
<ul style="list-style-type: none"> • + For non-food contact surfaces only. • ‡ For surfaces that are not inserted into the human mouth and body. • μ Allow surfaces to come to room temperature prior to treatment. • <i>For use on Hard, Non-Porous [Non-Food Contact Prep Surfaces⁺] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [▲ surfaces].</i> • <i>Any hard, non-porous, non-food contact surface where disinfection is required⁺.</i> 		

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces continued:

External Surfaces of Non-Critical Medical Devices	Fire Trucks -or- Vehicles	Handles
External Surfaces of Resuscitators	Floors	Handrails
External Surfaces of Instrument Reprocessors (<i>Note to Reviewer: Product can be used to disinfect exterior/ external surfaces of Instrument Reprocessors that are not to be sterilized.</i>)	Non-Food Contact Surfaces of Food Cases+	Hard Plastic or Metal Flooring
External Surfaces of Suction Pumps	Non-Food Contact Surfaces of Food Serving Trays+, Countertops+	Hard, non-porous environmental surfaces
External Surfaces of Surgery Microscopes	Furniture	Hard, non-porous environmental surfaces of Vehicles including steering wheel
External Surfaces of Surgery Robots	Formica®	High Chairs+
External Surfaces of Surgical Headlights	Exterior Surfaces of Drinking/ Water Fountains+	High Touch/ Frequently Used Areas/ Objects
External Surfaces of Surgical Light Sources	Exterior Surfaces of Freezers	Hoods+
External Surfaces of Tympanic Thermometers	Garbage Cans	Hospital Bed (Railings, Frames)
External Surfaces of Ultrasound Machines	Glass	Hospital Patient Bed
External Surfaces of Ultrasound Monitors	Glazed Ceramic Tile	Infant Incubators (interior and exterior surfaces)
External Surfaces of Ultrasound Transducers	Glazed Porcelain	Non-Food Contact Surfaces of Infant/ Child Care Equipment+
External Surfaces of Vaporizer	Glazed Tile Walls	Inside, Outside of Drawers
External Surfaces of Veterinary Equipment	Glazed Tiles	Instrument Stands
External Surfaces of Vital Sign Monitors	Golf Cart	IV [Stands] [Poles]
External Surfaces of X-Ray Machines	Grab Bar	Jacuzzi
Eye Glasses (do not use on contact lenses)	Grocery Cart	Kennels
Facilities [Surfaces]	Gurneys	Knee Walkers
Faucets	Gym Equipment	Knobs
Fax Machines	Hair Clippers	Laboratories
Non-Food Contact Surfaces of Feed Racks+	Hair Cutting Implements	Laboratory Equipment
Filing Cabinets	Hair Dryers	Laboratory Surfaces
Filing Doors	Hair Rollers	Laminated Surfaces
Filing Handles	Hampers	Lamps
<ul style="list-style-type: none"> • + For non-food contact surfaces only. • ≠ For surfaces that are not inserted into the human mouth and body. • μ Allow surfaces to come to room temperature prior to treatment. • <i>For use on Hard, Non-Porous [Non-Food Contact Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [▲ surfaces].</i> • <i>Any hard, non-porous, non-food contact surface where disinfection is required+.</i> 		

THE INFORMATION IN THIS BOX IS NOT PART OF THE LABEL

Statements in brackets, [], are optional or instructional. Bold, italicized text and italicized words within parentheses are information for the reader and are not included on production labels. Sequence of statements is optional unless specified in 40 CFR 156.10. Variables such as capitalizing, pluralizing, singularizing, bolding, italicizing and abbreviating are variable unless specified in 40 CFR. Lists of uses and use sites may be used partially or in their entirety. The term 'towelette (s)' is interchangeable with wipe(s), towel(s), and cloth(s) and the term 'container' is interchangeable with dispenser(s), bucket(s), pail(s), tub(s) throughout the label. The term center(s), facility(ies), area(s), unit(s), clinic(s), room(s), lab(s), department(s), suites, office(s), operatory(ies) in the use site tables are interchangeable as applicable.

Table 7 Surfaces continued:

Light Lens Covers	Non-Food Contact Pipes	Stainless Steel
Lights	Plastic Laminate	Stalls
Linoleum	Plastic Mattress Covers	Stands
Manicure Implements	Playpens	Stethoscopes
Manicure Tables	Pneumatic Tourniquets	Stools
Marlite	Police Cars -or- Vehicles	Storage [Units]
Mattress Covers	Porcelain Enamel	Non-Food Contact Surfaces of Stoves+, μ
Mayo Tray -or- Table	[Plastic] Posey Vest	Non-Food Contact Surfaces of Stovetops+, μ
Monitors	Public Restrooms	Stretchers
Nail/Hair Care Implements	Razors	Stretches
Non-food contact surfaces+	Reception [Areas] [Counters]	Support Bars
Nurse Call Button -or- Device -or- Cord	Exterior Surfaces of Refrigerator Carts	Synthetic Marble
NaugaHyde®	Exterior Surfaces of Refrigerators	Tables+
Operating Room [Furniture] [Light] [Light Switches] [Tables]	Seats	Tanning Beds
Operating Tables	Safety Straps	Telemetry Equipment
Outdoor Furniture+ (excluding wood frames and upholstery)	Salad Bar Sneeze Guards+	Telephones
Exterior Surfaces of Ovens	Salon Surfaces	Tiles+
Overbed table	School Buses	Toilet (Rim, Seat, Tank, Handle)
Oxygen Hoods	Scissors	Touch Screen Monitors
Patient Monitoring Equipment	Shampoo Bowls	Towel Dispensers
Patient Chairs	Shears	Toys, not for infant/ toddler toys
Patient Delivery Equipment	Shelves	Transfer Bench
Patient Lift Equipment	Shopping Carts	Transport Vehicles
Patient Scales	Shower [Fixtures] [Chair] [Stalls]	Trash Cans
Patient Support Equipment	Synthetic Marble	Tray Tables
Pens	Silk Peel™ with DermalInfusion™	Troughs
Phone Chargers	Sinks+	Tubing
Phone Cradle	Slit Lamps	Vanity Tops+
Physical Therapy (PT) Equipment surfaces	Smart Phone Chargers	Veterinary Care Surfaces
Wash Basin	Spine Back Boards	Vinyl
Washable Nail Files	Tanks	Walkers
Workstations	Walls	Wheelchairs
Rubber [material]	Watering Appliances	Whitening Lamps
Titanium-Coated [material]	Weight-lifting Surfaces	Titanium-Coated [material]
Brass [material]	Acrylic [material]	Polycarbonate [material]
Medical Grade [material]		
<ul style="list-style-type: none"> • + For non-food contact surfaces only. • ¥ For surfaces that are not inserted into the human mouth and body. • μ Allow surfaces to come to room temperature prior to treatment. • <i>For use on Hard, Non-Porous [Non-Food Contact Prep Surfaces+] [Surfaces] [Medical Surfaces] [Veterinary Surfaces] [Dental Surfaces] [Institutional Surfaces] [▲ surfaces].</i> • <i>Any hard, non-porous, non-food contact surface where disinfection is required+.</i> 		