

PM32

46851-1

753
FEB 18 1986

ADM Medical Division, Inc.
dba Omnitec Medical Corporation
25283 Cabot Road, Suite 109
Laguna Hills, CA 92653

Attention: Palmer B. Ford

Gentlemen:

Subject: Omni II
EPA Registration No. 46851-1
Your Submission Dated January 2, 1986

The revised draft label referred to above, submitted in connection with registration under section 3(c)(7)(A) of the Federal Insecticide, Fungicide, and Rodenticide Act (FIFRA), is acceptable provided that you comply with the following conditions:

1. Make the labeling changes listed below on the finished printed labels.
 - a. Delete the ADA seal. It is misleading labeling in that it indirectly implies recommendation of the pesticide by this agency (§162.10(a)(5)(v) of 40 CFR Ch. 1, July 7, 1985).
 - b. Change the precautionary labeling to conform to the regulations, 40 CFR Ch. 1 (July 7, 1985) §162.10, p. 97. Your accepted label of August 17, 1983 was correct in this regard.
 - c. The finished printed labels must comply with the general labeling requirements of §162.10 of Title 40 of the Code of Federal Regulations.
2. Submit five (5) copies of your final printed labeling, as defined in the A-79 Enclosure, before you release the product for shipment, to comply with FR Notice 82-2 of June 18, 1982.

86907:Lucke:C-2:KENCO:2/10/86:3/7/86:ab:vo/ek

SYMBOL								
SURNAME	✓							
DATE	FEB 18 1986							

2

If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA section 6(e). Your release for shipment of the product bearing the amended labeling constitutes acceptance of these conditions.

A stamped copy of the label is enclosed for your records.

Sincerely yours,



A. E. Castillo
Product Manager (32)
Disinfectants Branch
Registration Division (TS-767C)

Enclosures

OMNI-II™

INSTRUMENT & EQUIPMENT DISINFECTANT

Active Ingredients
o-phenylphenol
o-Benzyl-p-chlorophenol
Inert Ingredients

9.0%
1.0%
90.0%
100.0%

- Bacteriocidal
- Fungicidal --- destroys pathogenic fungi
- Virucidal*
- Tuberculocidal

* Viruses including Herpes virus types 1 & 2, Polio virus types 1 & 2, Vaccinia virus, Influenza virus type A1 (New Jersey) and Influenza B (Maryland)

EPA Reg No 46851-1
Est No 7368-CA-3



DIRECTIONS FOR USE:

Omni II should be diluted 1:32 with water for disinfection of medical and dental instruments and devices. Rinse thoroughly after immersion. Discard the used solution daily, or sooner if it becomes diluted or soiled. It is in violation of federal law to use this product in a manner inconsistent with its labeling.

Clean all excess organic debris from items to be disinfected and rough dry. Immerse in 1:32 concentration of Omni II for 10 minutes to destroy bacteria including *Escherichia coli* O157, *Staphylococcus aureus*, *Salmonella typhosa*, *Pseudomonas aeruginosa*, *Mycobacterium tuberculosis* III. Tuberculosis contamination is suspected, increase immersion time to 20 minutes. For added safety, fungi including *Trichophyton interdigitale*, viruses including Herpes virus types 1 & 2, Polio virus types 1 & 2, Vaccinia virus, Influenza virus

To prepare 1:32 use dilution	
WATER	OMNI II
1 pint	1/2 oz. (15 ml)
1 quart	1 oz. (30 ml)
1/2 gallon	2 oz. (60 ml)
1 gallon	4 oz. (120 ml)

type A1 (New Jersey) and Influenza B (Maryland) destroy the etiologic agent for demonstration.

For use on metals, plastics, rubber, repeated exposure. Don't store. When diluted as directed, this product is based on previously unorganic material or water introduced. The product effectiveness.

Omni II should not be used for items with deep and/or narrow

Do not reuse empty container.

PRECAUTIONS

DANGER: Corrosive. Causes eye or on clothing. Wear goggles or face shield. Harmful or fatal if swallowed.

FIRST AID: In case of contact, immediately flush with water for at least 15 minutes. For eyes, call a physician before reuse. If swallowed, egg whites, gelatin solution, or milk of water. Avoid alcohol. Call

NOTE TO PHYSICIAN: Probable mucous membrane irritation.

KEEP OUT OF REACH OF CHILDREN

DANGER:

Manufactured & Distributed by:
ADM MEDICAL DIVISION

Contents: 1/2 gallon (64 fluid oz.)

25283 Cabot, Laguna Hills, Ca. 92653

DIRECTIONS FOR USE:

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING.

OMNI II should be diluted 1:32 with water for disinfection of medical and dental instruments and devices. after immersion. Other hard surface areas such as countertops should be mopped, sponged or thoroughly wet the surface with OMNI II. Disinfect daily, or sooner if it becomes diluted or soiled.

ACCEPTED
with CONCURRENCE
in EPA Letter #

FEB 13 1985

Under the
EPA

46851-1

BEST AVAILABLE COPY

LABEL CHANGES REQUESTED:

✓ 1) PARAGRAPH #1 IN DIRECTIONS FOR USE TO READ

✗ 2) PRECAUTIONS CHANGED TO

✗ 3) ADDITION OF ADA LABEL

✗ 4) ADDRESS CHANGE TO:
LAGUNA HILLS, CA

PRECAUTIONS: Harmful if swallowed. Avoid contact with mucous membranes. Wear rubber gloves recommended. Thoroughly shield eyes. If contact occurs, flush with water. Seek medical attention. Avoid container thoroughly with