

# Omni II

[Alternate Brand Names: "Vital Defense-D" and "Prophene Plus"]

## Concentrated Surface Disinfectant and Holding Solution

Tuberculocidal, Virucidal;  
Bactericidal, Fungicidal

**Active Ingredients**

o-phenylphenol	9.0%
o-benzyl-p-chlorophenol	1.0%
Inert Ingredients	90.0%
Total	100.0%

Contents: 1 Quart  
(32 oz. - 946 ml)

Reorder No. OMN232

E.P.A. Reg. No. 46851-  
E.P.A. Est. No. 46851-CD-3

**DANGER: Keep Out of Reach of Children**  
See Side Panel for Additional Precautionary Statements

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV ON SURFACES FOR OBJECTS SOILED WITH BLOOD OR BODY FLUIDS**

**PERSONAL PROTECTION:** When handling items soiled with blood or body fluids or when conducting procedures which might lead to contact with blood or body fluids, always wear items such as heavy duty nitrile gloves, gowns, masks and eye coverings.

**CLEANING PROCEDURES:** Blood and body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant.

**DISPOSAL OF INFECTIOUS MATERIAL:** Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

**PRECAUTIONARY STATEMENTS**  
Hazard to Humans

**DANGER: Corrosive:** Causes eye and skin damage. Do not get in eyes, on skin or on clothing. Wear protective eyewear and rubber gloves when handling. Harmful if swallowed.

**STATEMENT OF PRACTICAL TREATMENT**  
In case of contact, immediately flood eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash contaminated clothing before reuse. If swallowed, drink promptly a large quantity of water. Avoid alcohol. Call a physician immediately.

**NOTE TO PHYSICIAN:** Probable mucosal damage may contraindicate the use of gastric lavage.

**STORAGE & DISPOSAL:** Keep Omni II covered in storage or use. Do not contaminate water, food or feed by storage or disposal. Do not reuse empty container. Triple rinse container, puncture and dispose of in a sanitary landfill or incinerator. Used disinfectant solutions are toxic. Store, use and dispose of according to label instructions.

**DIRECTIONS FOR USE**

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING. Omni II concentrated solution and must be diluted with water prior to instrument soak and/or decontamination of dental and medical instruments or disinfection of hard environmental surfaces. This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

**Soak/Decontamination** - Immerse items in diluted Omni II solution for instrument soak and/or decontamination of medical and dental instruments. Rinse articles thoroughly following soak/decontamination. Organic materials or water introduced into the diluted product may affect product's effectiveness. Discard the used solution daily or sooner if it becomes diluted or soiled.

**Ultrasonic Cleaning** - Place instruments into diluted Omni II and activate ultrasonic unit. When instruments are clean, remove and rinse thoroughly. Change solution daily or sooner if it becomes diluted or soiled.

**Cleaning** - Apply Omni II directly to surface using mop or spray bottle. Wipe surface to remove debris and solution using a clean paper or cloth towel and discard.

**Hard Surface Disinfection** - diluted Omni II solution should be mopped, sponged, sprayed or wiped onto hard environmental surfaces such as

walls, floors and countertops in a manner to thoroughly wet the surface. A properly prepared solution of Omni II intended for use as a hard surface disinfectant, has a shelf life of 60 days (use life) when stored in a closed container such as a spray bottle.

A contact time of 10 minutes at 20°C (68°F) with a 1:32 dilution of Omni II has been shown to be effective on bacteria including: Escherichia Coli 0113, Staphylococcus aureus, Salmonella typhosa, Pseudomonas aeruginosa, Mycobacterium tuberculosis; fungi including: Trichophyton mentagrophytes; viruses\* including influenza virus type A1 (New Jersey) and Influenza B (Maryland), Herpes virus types 1 and 2, Vaccinia virus, and Polio virus type 1.

Omni II when diluted 1:32 has been shown to be effective in destroying the HIV-1 (AIDS virus) in 60 seconds at 20°C (68°F). Although efficacy in 60 second contact time has been shown to be adequate against HIV-1, this time would not be sufficient for other organisms. Use 10 minute contact time for disinfection against all the organisms claimed. Effective in waters up to 400 ppm hardness calculated as calcium carbonate.

For use on metals, plastics, rubber and glass. Glass may etch over prolonged or repeated exposure to concentrated solutions. Do not store lensed articles in concentrated Omni II.

MADE IN U.S.A. BY  
COTTRELL LTD., A DIVISION OF  
CERTOL INTERNATIONAL LLC  
5430 JOLIET STREET  
DENVER, CO 80239 USA  
1 308 THE EDGE

TO PREPARE 1:32 USE-DILUTION	
WATER	OMNI II
To Make 1 Pint	
15 1/4 fl. oz.	1/2 fl. oz. (15 ml.)
To Make 1 Quart	
31 fl. oz.	1 fl. oz. (30 ml.)
To Make 1/2 Gallon	
62 fl. oz.	2 fl. oz. (60 ml.)
To Make 1 Gallon	
124 fl. oz.	4 fl. oz. (120 ml.)

**ACCEPTED**  
OCT 30 2002  
Under the Federal Insecticide, Fungicide, and Rodenticide Act as amended, for the pesticide, registered under EPA Reg. No. 46851-1

Note to EPA reviewers: This label (proposed draft for 46851-1) has been recently reviewed and approved (Sept. 17, 2002) for use with Cottrell's "re-pack" product 46851-11, which has an identical formulation and use sites as this product. This "minor" label amendment request will allow Cottrell to use this brand name and label for its product bearing the number 46851-1. Alternate brand names Vital Defense-D and Prophene Plus also will be used for this product, when redesignated by this amendment.

46851-1

10-30-2002

4/