

Omni II™

Concentrated Surface Disinfectant and Holding Solution

Tuberculocidal, Virucidal, Bactericidal, Fungicidal

Active Ingredients
 9.9% Chloroxydine
 1.0% Benzalkonium Chloride
 90.0% Water
 100.0% Total

EPA Reg. No. 46851-1
EPA Est. No. 46851-CO-033

Contents: 1 Quart (32 oz. - 946 ml)
Reorder No. OMN232

KEEP OUT OF REACH OF CHILDREN

DANGER: See Side Panel for Additional Precautionary Statements

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV ON SURFACES OR OBJECTS SOILED WITH BLOOD OR BODY FLUIDS
 When handling items soiled with blood or body fluids or when cleaning procedures which might lead to contact with blood or body fluids, always wear items such as disposable latex gloves, gowns, masks and eye coverings.
CLEANING PROCEDURES: Blood and body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant.
CONTACT TIME: Omni II when diluted 1:32 has been shown to be effective in destroying the HIV-1 (AIDS virus) in 60 seconds at 20°C. This contact time will not control other common types of viruses and bacteria.
DISPOSAL OF INFECTIOUS MATERIAL: Blood and other body fluids should be decontaminated and disposed of according to local regulations for infectious waste disposal.
PRECAUTIONARY STATEMENT: Hazards to Humans and Domestic Animals
DANGER: Corrosive. Causes irreversible eye damage and skin burns. Harmful if inhaled, swallowed or absorbed through skin. Do not get in eyes, on skin or on clothing. Avoid breathing spray mist. Wear goggles or face shield (safety glasses), protective clothing (long-sleeved shirt and long pants, sock, flat shoes and gloves) and avoid contact with such as water runoff (spray).
USERS SAFETY REQUIREMENTS: Federal requirements: For use on surfaces soiled with blood or body fluids, use of PPE is required. For use on surfaces soiled with blood or body fluids, use of PPE is required. For use on surfaces soiled with blood or body fluids, use of PPE is required.


USERS SAFETY RECOMMENDATIONS: Users must wash hands before eating, drinking, chewing gum, using tobacco, or using the toilet. Users must remove clothing immediately if possible gaiters. Then wash thoroughly and put on clean clothing. Users must remove PPE immediately after handling the product. Wash the outside of gloves before removing. As soon as possible, wash thoroughly and clean clothing.
FIRST AID
IF IN EYES: Hold eye open and flush slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue flushing. Call a poison control center or doctor for treatment advice.
IF ON SKIN OR CLOTHING: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.
IF SWALLOWED: Call a Poison Control Center or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told by a Poison Control Center or doctor. Do not give anything by mouth to an unconscious person.
IF INHALED: Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration (preferably mouth-to-mouth) if possible. Call a Poison Control Center or doctor for further treatment advice.
 In case of an emergency for additional information call toll free 800-535-5653. Have the product container or label with you when calling a poison control center or doctor, or going for treatment.
NOTE TO PHYSICIAN
 Possible mucous membrane damage may contraindicate the use of gastric lavage.

DIRECTIONS FOR USE
 IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING.
 Omni II is designed for use on hard non-porous surfaces and/or in holding solutions in medical, dental, laboratory and veterinary facilities. Thoroughly clean surface of items to be disinfected and rough dry. Omni II is a concentrated solution and must be diluted 1:32 with water prior to use. Use the following directions for dental and medical instruments or disinfectant and hold environmental surfaces. This product is not to be used as a terminal sterilant/high-level disinfectant on any surgical instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.
Soak/Decontamination: Immerse items in diluted Omni II solution for treatment soak and/or decontamination of medical and dental instruments. Rinse articles thoroughly following soak/decontamination. Organic materials or water introduced into the diluent may affect product effectiveness. Discard the used solution daily or sooner if it becomes diluted or soiled.
Ultrasonic Cleaning: Place instruments into diluted Omni II and activate ultrasonic unit. When instruments are clean, remove and rinse thoroughly. Change solution daily or sooner if it becomes diluted or soiled.

Cleaning: Apply Omni II directly to surface using mop or spray bottle. Wipe surface to remove debris and solution using a clean paper or cloth towel and discard.
Hard Surface Disinfection: diluted Omni II solution should be mopped, sprayed, or wiped onto hard environmental surfaces such as walls, floors and countertops in a manner to thoroughly wet the surface. A properly prepared solution of Omni II intended for use as a hard surface disinfectant, has a shelf life of 60 days (up to 90) when stored in a closed container such as a spray bottle.
 A contact time of 10 minutes at 20°C (68°F) of a 1:32 dilution of Omni II has been shown to be effective on bacteria including Escherichia Coli 0111, Staphylococcus aureus, Salmonella typhosa, Pseudomonas aeruginosa, Mycobacterium tuberculosis, phages, viruses including influenza virus (type A) (New Jersey) and Influenza B (Maryland), herpes virus types 1 and 2, Vaccinia virus, and Polio virus 1.
 Omni II when diluted 1:32 has been shown to be effective in destroying the HIV-1 (AIDS virus) in 60 seconds at 20°C (68°F). Although efficacy in 60 second contact time has been shown to be adequate against HIV-1, this time would not be sufficient for other organisms. Use 10 minute contact time for disinfection against all the organisms claimed. Effective in waters up to 400 ppm hardness calculated as calcium carbonate.
 For use on metals, plastics, rubber and glass. Glass may etch over prolonged or repeated exposure to concentrated solutions. Do not store lensed articles in concentrated Omni II.

TO PREPARE 1:32 USE-DILUTION WATER - OMI II	
To Make: 1 Quart	1 fl. oz. (30 ml)
To Make: 1/2 Gallon	2 fl. oz. (60 ml)
To Make: 1 Gallon	4 fl. oz. (120 ml)

STORAGE AND DISPOSAL
 Do not contaminate water, food or feed by storage or disposal.
PESTICIDE STORAGE: Keep Omni II covered in storage or in use.
PESTICIDE DISPOSAL: Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture, or empty container is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste Representative at the nearest EPA Regional Office for guidance.
CONTAINER DISPOSAL: Do not reuse empty container. Triple rinse container, puncture and dispose of in a sanitary landfill, or incinerate.



ACCEPTED
 SEP - 7 2004

Under the Federal Insecticide, Fungicide, and Rodenticide Act as amended, for the pesticide, registered under EPA Reg. No. 46851-1

MANUFACTURED BY:
CERTOL INTERNATIONAL, LLC
 6120 E. 58TH AVE.
 COMMERCE CITY, CO 80022 USA
 1 800 THE EDGE
 (1 800 843 3343)

OMN232B.US