

TIME-MIST

METERED AIR SANITIZER

TIME-MIST Air Sanitizer is engineered for highly effective use in the **TIME-MIST** Automatic Aerosol Dispenser.

The **TIME-MIST** Dispenser automatically actuates the special metered valve every fifteen minutes, treating the air against odors and providing effective air sanitization. Within 24 hours of continued use **TIME-MIST** Air Sanitizer will prove highly effective in reducing strep, staph, and other airborne bacteria virus particles. One **TIME-MIST** Dispenser is designed to treat a room up to 30' x 20' x 10'.

OPERATION: Remove protective cap from top of can and insert can in **TIME-MIST** Dispenser. Suspend Dispenser 7" from floor.

EPA REG NO 7710-1
EPA Est 10007-GA-1

ACCEPTE
10 JAN 1981

IMPORTANT

When used with the **TIME-MIST** Automatic Dispenser, one unit should be installed for each 6000 cubic foot area to be treated. It is usually advisable to mount the dispenser 7" from the floor. Do not use in areas where the temperature is below 50°F. Do not use in areas where the humidity is above 70%.

DISPOSAL: Do not reuse empty container. Wrap container and put in trash. Do not use in areas where the temperature is below 50°F. Do not use in areas where the humidity is above 70%.

CAUTION: Do not use in nurseries or rooms where infants or aged persons are confined. Do not use in kitchens, dining areas or other areas where food is prepared or served. Avoid contact with eyes.

Manufactured for: **TIME-MIST INC.**
135 South Main Street
Thomaston, Connecticut 06787

TIME-MIST

SPECIAL FORMULA AIR SANITIZER

3000 SPRAYS

METERED FOR EFFICIENCY AND ECONOMY

ACTIVE INGREDIENTS
Triethylene Glycol 9.15%
Dipropylene Glycol 3.43%
Aryl dimethyl Benzyl Ammonium Chloride 17.1%
Essential Oils 50.0%
Ethanol 31.93%

INERT INGREDIENTS 54.40%
*C₁₀ 61% C₁₂ 23%
C₁₄ 11% C₁₆ 0% and C₁₈ 5%

1. Sanitizes, deodorizes and refreshes the air. 2. Reduces airborne bacteria and many virus particles. 3. Minimizes danger of spreading germs and viruses through the air.

CAUTION: Keep Out of Reach of Children. Contents under pressure. See back panel for other cautions.

THIS TIME-MIST AEROSOL IS OF THE ENVIRONMENTAL FORMULA AND ENGINEERED FOR HIGHLY EFFECTIVE USE IN THE TIME-MIST AUTOMATIC AEROSOL DISPENSER.
NET WT. 7 oz.

READ INSTRUCTIONS
BEFORE USING THIS PRODUCT
AND KEEP WITH ITS LABELING.

Rinse articles thoroughly with sterile water after disinfection or sterilization.

For Complete Sterilization — to destroy resistant pathogenic spores including — *Bacillus subtilis* and *Clostridium sporogenes* — immerse articles for minimum of 10 hours. Remove instruments using sterile technique and rinse thoroughly in sterile water.

All instruments and devices must be thoroughly cleaned to remove excess organic debris, rinsed and rough dried before immersion in Omnicide. Hypodermic needles cannot be effectively disinfected or sterilized with Omnicide.

All labeling claims of effectiveness are based on previously unused batches of product.

Dating Record
OMNICIDE HAS A 28-DAY
SHELF LIFE AFTER ACTIVATION.

ACTIVATED _____

EXPIRATION _____

Omnicide[★]

STERILIZING and DISINFECTING SOLUTION

For use in hospitals, dental offices, nursing homes, and health care institutions.

Cat No. U8998

28 Day Activated Dialdehyde

Qty. 1 Gallon
(3.79 liters)

- Sporicidal
- Virucidal†
- Fungicidal
- Tuberculocidal
- Bactericidal
- Pseudomonacidal

Non-Rusting Complete Sterilant Pleasant Fragrance

† Influenza virus types A1 (New Jersey) and B (Maryland), Herpes simplex virus types 1 and 2, Polio virus types 1 and 2 and Vaccinia virus.

E.P.A. Establishment No. 7368-CA-3
E.P.A. Registration No. 46851-2

Active Ingredient	Inert Ingredients	Total
Glutaraldehyde		100.0%

KEEP OUT OF REACH OF CHILDREN
DANGER:
See side panel for first aid instructions.

ACCEPTED

AUG 11 1983

Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under 46851-2
EPA Reg. No.

Recommended for Disinfection or Sterilization of: MEDICAL and DENTAL INSTRUMENTS, CARBON STEEL CUTTING INSTRUMENTS, BURS, ALUMINUM, GLASS, PLASTIC or RUBBER ARTICLES, THERMOMETERS, MIRRORS, RESPIRATORY EQUIPMENT.

STORAGE AND DISPOSAL—STORE IN A COOL PLACE

PROHIBITIONS:

Do not contaminate water, food or feed by storage or disposal. Open dumping is prohibited. Do not reuse empty container.

PESTICIDE DISPOSAL

Pesticide, spray mixture, or rinse water that cannot be used according to label instructions must be disposed of according to Federal or approved state procedures under Subtitle C of the Resource Conservation and Recovery Act.

CONTAINER DISPOSAL

Triple rinse (or equivalent). Then offer for recycling or reconditioning, or dispose of in a sanitary landfill, or by incineration if allowed by state and local authorities.

GENERAL

Consult Federal, state or local disposal authorities for approved alternative procedures.

PRECAUTIONARY STATEMENTS:

HAZARDS TO HUMANS AND DOMESTIC ANIMALS.

DANGER. Corrosive. Causes burns. Causes eye damage and skin irritation. Do not get in eyes, on skin or on clothing. May be absorbed through skin. Wear goggles or face shield and rubber gloves when handling.

HARMFUL IF SWALLOWED. Avoid contact with food. Use in ventilated area. Use containers which can be tightly closed during use or storage of solution.

STATEMENT OF PRACTICAL TREATMENT

Eyes — in case of skin contact, immediately flush thoroughly with water. Skin — in case of skin contact, immediately flush thoroughly with water. Ingestion: If swallowed, drink milk, egg whites, gelatin solution or, if these are not available, large quantities of water. Avoid alcohol. Call a physician.

Note to Physician: Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed.

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH THE LABEL DIRECTIONS.

DIRECTIONS FOR ACTIVATION AND USE:

To activate solution — add activator contents to this one gallon container — recap — shake. Activated solution changes color to blue. Do not dilute. Discard daily or sooner if used solution becomes soiled or diluted.


KEEP OMNICIDE COVERED IN STORAGE OR USE & STORE IN A COOL PLACE.

Always clean away excess organic debris, rinse and rough dry instruments or equipment before immersion.

FOR DISINFECTION:

Immerse articles for a minimum of 10 minutes to destroy vegetative pathogens including — *Escherichia coli* 0113, *Staphylococcus aureus*, *Salmonella typhosa*, *Pseudomonas aeruginosa*, *Mycobacterium tuberculosis* (if tuberculosis contamination is suspected, increase immersion time to 20 minutes); fungi including — *Trichophyton interdigitale*, viruses including — Herpes simplex virus types 1 & 2, Polio virus types 1 & 2, Vaccinia virus, Influenza virus types A1 (New Jersey) and B (Maryland). The ability of this product to destroy the etiologic agent for viral hepatitis has not been established or demonstrated.

APB



FRONT PANEL

Omnicide*

Sterilizing and Disinfecting
Solution Activator

Contents: 57 grams

ADM Medical Equipment

BACK PANEL

Activator

To be used in conjunction with
Omnicide

Directions For Use:

1. Add contents of container to 1 liter of water.
2. Use for disinfecting and sterilizing.

DANGER:

1. Do not inhale dust or fumes.
2. Do not get in eyes or on skin.
3. If in eyes, flush with water.
4. If on skin, wash with water.
5. Keep out of reach of children.
6. Store in a cool, dry place.

ACCEPTED

AUG 11 1983

Under the Federal Insecticide,
Fungicide, and Rodenticide Act,
as amended, for the pesticide
registered under
EPA Reg. No. 46857-2