

PM-31

Reg # 46851-2

INTENDED USE

Sterilization: This solution should be used for the sterilization of heat sensitive medical equipment for which alternative methods of sterilization are not suitable. Medical equipment which should always be sterilized is that which is categorized as critical (e.g., used in procedures in which contact will be made with tissue that is normally considered sterile).

High Level Disinfection: This solution should be used for the high level disinfection of heat sensitive medical equipment for which sterilization is not practical. Medical equipment which should always be subjected to high level disinfection or sterilization is that which will be used in procedures categorized as semi-critical (e.g., used in procedures in which contact will be made with mucous membranes or other body surfaces which are not normally considered sterile).

Intermediate Level Disinfection: This solution should be used for disinfection of medical equipment for which a risk of cross contamination exists. Medical equipment which should always be disinfected is that which is to be used in procedures categorized as non-critical (e.g., used in procedures in which contact will only be made with intact skin).

PRECAUTIONARY STATEMENTS**HAZARDS TO HUMANS AND DOMESTIC ANIMALS**

DANGER: Keep Out of Reach of Children.

Direct contact may cause eye damage and skin irritation/damage. Do not get into eyes, on skin or on clothing. Wear goggles or face shield and chemical resistant gloves when handling or pouring. Avoid contamination of food. Use in a well ventilated area in closed containers.

STATEMENT OF PRACTICAL TREATMENT: In case of contact, immediately flush eyes or skin with copious amounts of water for at least

15 minutes. For eyes, get medical attention. Harmful if swallowed. Drink large quantities of water and call a physician immediately.

NOTE TO PHYSICIAN: Probable mucosal damage from oral exposure may contraindicate the use of gastric lavage.

STORAGE AND DISPOSAL

Store at controlled room temperature 15°C-30°C (59°-86°F).

Pesticide Disposal: Discard residual solution in drain. Flush thoroughly with water.

Container Disposal: Do not reuse empty container. Wrap container and put in trash.

NFPA HAZARDOUS CHEMICAL CHART

Health Hazard	Flammability
1 (Slightly Hazardous)	1 (Slightly Hazardous)
2 (Moderately Hazardous)	2 (Moderately Hazardous)
3 (Severely Hazardous)	3 (Severely Hazardous)
4 (Extremely Hazardous)	4 (Extremely Hazardous)
5 (Fatal)	5 (Fatal)
6 (Lethal)	6 (Lethal)
7 (Corrosive)	7 (Corrosive)
8 (Very Corrosive)	8 (Very Corrosive)
9 (Extremely Corrosive)	9 (Extremely Corrosive)
10 (Fatal)	10 (Fatal)

Product Name: OmniCide
Hazardous Chemicals: Glutaraldehyde

Route of Entry:

- ☒ Inhalation
- ☒ Ingestion
- ☒ Skin/Eye Absorption

Personal Protection: Chemical Resistant Gloves, Goggles or Face Shield

COTTRELL, LTD.
7399 SOUTH TUCSON WAY
ENGLEWOOD, CO 80112 USA
1 800 THE EDGE



H8160M322805

OmniCide™

Long Life Reusable Activated Dialdehyde

Sterilizing and Disinfecting Solution

Sporicidal, Virucidal[†], Fungicidal,
Pseudomonacidal, Bactericidal,
Tuberculocidal

Recommended for use in hospitals, medical/dental clinics, nursing homes, healthcare institutions and veterinary clinics.

Active Ingredient

Glutaraldehyde	2.4%
Inert Ingredient	97.6%
TOTAL	100.0%

E.P.A. Registration No. 46851-2

E.P.A. Establishment No. 37265-CA-01

E.P.A. Establishment No. 39754-WI-1

(See shoulder of container for E.P.A. Est. No.)

Note: Contents of attached vial must be added to solution before this product is effective. See "DIRECTIONS FOR USE - ACTIVATE."

[†] See package insert for specific viruses claimed.

Contents: 1 Quart
(946 ml)
Reorder No.
OM32/28

DIRECTIONS FOR USE

IT IS VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING

Activation: Activate the OmniCide solution by adding the entire contents of the Activator Bottle (which is attached to the OmniCide solution container) to the container. Place cap on container and shake well. On activation the solution immediately changes color to blue, thereby indicating solution is ready to use. Record the date of activation on the indicated space below, in a log book or a label affixed to any secondary container used for activated solution. See package insert for additional instructions and information regarding activated solution.

Cleaning/Decontamination: Blood and other body fluids must be thoroughly cleaned from the surfaces, lumens, and objects before application of the disinfectant or sterilant. Blood and other body fluids should be autoclaved and disposed of according to all applicable Federal, State and local regulations for infectious waste disposal.

For complete disinfection or sterilization of medical instruments and equipment, thoroughly clean, rinse and rough dry objects before immersing in OmniCide solution. Cleanse and rinse the lumens of hollow instruments before filling with OmniCide solution. See package insert for additional cleaning/decontamination instructions.

Sterilization: Immerse medical equipment/devices completely in OmniCide solution for a minimum of ten hours at 68°F (20°C) to eliminate all microorganisms, including *Clostridium sporogenes* and *Bacillus*

subtilis spores. Remove solution using sterile technique. High level disinfection for complete information on sterilization.

High Level Disinfection: Immerse medical equipment in OmniCide solution 45 minutes at 68°F to eliminate large numbers of pathogens, but including *Mycobacterium tuberculosis*.

Remove devices from solution and rinse in water or potable water used is depend on use of the instrument. See package insert for complete information on high level disinfection.


Intermediate Level Disinfection: Immerse medical instruments in OmniCide for a minimum of 10 minutes at 68°F (20°C) to disinfect bacteria, all pathogen viruses and 99.99% strains. (Quantitative disinfection represented by Davis for complete instructions on intermediate level disinfection.)

Refer to package insert for usage/product data.

Note Dates Here:
Activation Date:

Check solution upon use with ProChek Concentration Indicator. Discard after 2 ProChek G Glutaraldehyde Indicator in Control Number and Unactivated Sol.

DANGER: Keep Out of Reach of Children
See Side Panel for Additional Precautionary Statements



OmniCide™

Long Life Reusable
Activated Dialdehyde

Sterilizing and Disinfecting Solution

COTTRELL, LTD.
7399 SOUTH TULSON WAY
ENGLEWOOD, CO 80112 USA
1 800 THE EDGE

A. INTENDED USE/DIRECTIONS FOR USE

OmniCide long life activated dialdehyde solution is a liquid chemical sterilant and a high or intermediate level disinfectant when used according to the DIRECTIONS FOR USE.

1. **Germicide Level of Activity:** OmniCide can be used at the following germicide levels of activity:

Sterilant: OmniCide is a sterilant when used or reused, according to the Directions For Use, at or above its Minimum Effective Concentration (MEC) as determined by the ProCheck G Glutaraldehyde Concentration Indicator Test Strip at 68°F (20°C) with an immersion time of ten hours for a use period not to exceed 28 days.

High Level Disinfectant: OmniCide is a high level disinfectant when used or reused, according to the DIRECTIONS FOR USE, at or above its Minimum Effective Concentration (MEC) as determined by the ProCheck G Glutaraldehyde Concentration Indicator Test Strip at 68°F (20°C) with an immersion time of 45 minutes for a use period not to exceed 28 days.

Intermediate Level Disinfectant: OmniCide is an intermediate level disinfectant when used or reused, at or above its Minimum Effective Concentration (MEC) as determined by the ProCheck G Glutaraldehyde Concentration Indicator Test Strip at 68°F (20°C) with an immersion time of ten minutes according to the DIRECTIONS FOR USE in SECTION E for a use period not to exceed 28 days.

A ten minute immersion at 68°F (20°C) will destroy vegetative bacteria including *S. aureus*, *S. choleraesuis*, *P. aeruginosa*, *E. coli*, all pathogenic fungi, representative viruses and 99.99% of *Mycobacterium* strains (Quantitative TB Method) as represented by *bovis*. ** HIV-1 (AIDS virus) is destroyed in 60 seconds with full strength solution.

2. **Reuse Period:** OmniCide has demonstrated efficacy in the presence of 5 percent organic soil contamination and a simulated amount of microbiological burden during reuse. OMNICODE SOLUTION MAY BE REUSED ONLY WHILE THE MINIMUM EFFECTIVE CONCENTRATION (MEC) AS DETERMINED BY THE PROCHECK G GLUTARALDEHYDE CONCENTRATION INDICATOR TEST STRIP, PH AND TEMPERATURE MEET THE REQUIREMENTS BASED UPON MONITORING AS DESCRIBED IN THE PROCHECK G DIRECTIONS FOR USE. Efficacy of OmniCide solution during its use-life must be verified by the ProCheck G Glutaraldehyde Concentration Indicator Test Strip to determine that at least the MEC as determined by the ProCheck G Glutaraldehyde Concentration Indicator Test Strip is present. This product must be discarded

after 28 days, even if the ProCheck G Concentration Indicator indicates pass.

3. **General Information:** Choose a germicide with the that is appropriate for the reusable medical device. Follow the reusable device labeling and standard in absence of complete instructions, use the following:

First, for patient contacting devices, determine whether the device to be processed is a critical, semi-critical, or non-critical device.

- A critical device routinely penetrates the skin during use or are otherwise used in normally sterile areas.
- A semi-critical device makes contact with mucous membranes but does not ordinarily penetrate sterile areas.
- A non-critical device contacts only intact skin.

Second, determine if sterilization, high level or intermediate level disinfection is required.

Critical device:

Sterilization is required.

Semi-critical device:

Although sterilization is recommended whenever possible, high level disinfection is acceptable (e.g., GI endoscope) if the device is used in the airway, diaphragm-fitting rings.

Non-critical device:

Medical Equipment Surfaces: Intermediate level disinfection is recommended.

Third, determine the time required to achieve the level of sterilization required for the specified medical device.

4. The germicidal activity of OmniCide was demonstrated in performance, clinical and simulated use on the following organisms:

ACCEPTED

MAR 4 1996

Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the pesticide registered under EPA Reg. No. 46851-2



Organisms - Omnicide Sterilization/Disinfection Times

Spores - 10 hours

- Bacillus subtilis*
- Clostridium sporogenes*

Vegetative Organisms - 10 minutes

- Staphylococcus aureus*
- Salmonella choleraesuis*
- Pseudomonas aeruginosa*
- Escherichia coli*

Vegetative Organisms - 45 minutes

- Mycobacterium bovis*

Fungi - 10 minutes

- Trichophyton interdigitale*
- Trichophyton mentagrophytes*

Non-lipid Small Virus* - 10 minutes

- Polio 2

Lipid Medium Virus†

- Herpes simplex (α, β, γ)*
- HIV-1 (AIDS virus)* (60 sec.)

* Testing was performed using Omnicide solution which had been diluted to 1.5 percent using 5 percent bovine calf serum

5. **Material Compatibility:** Omnicide solution is recommended for usage with medical devices made from the materials shown below.

- | | |
|---------------------|----------------------------|
| • polypropylene | • vinyl and Tygon tubing** |
| • ABS | • Nickel plating* |
| • polyethylene* | • acrylic bar* |
| • polycarbonate | • polyethylene tubing* |
| • black oxide steel | • PVC* |
| • Mylar* | |

* represents four weeks of exposure ** represents 10 days of exposure

Following sterilization or disinfection, the sterilized or disinfected medical device should be rinsed according to the **DIRECTIONS FOR USE, RINSING** (SECTION E.4), and dried according to manufacturer's instructions.

6. **Pre-cleaning / gent Compatibility:** Omnicide is compatible with enzymatic detergents which are neutral in pH, low foaming and easily rinsed from equipment. Detergents that are either highly acidic or alkaline are contraindicated as pre-cleaning agents since improper rinsing could effect the efficacy of the Omnicide solution by altering its pH.

3

B. CONTRAINDICATIONS

1. **Sterilant Usage:** Routine biological monitoring is not possible with Omnicide solution and therefore Omnicide solution should NOT be used to sterilize reusable medical devices that are compatible with other sterilization processes that can be biologically monitored.

Omnicide solution should not be used for sterilization of critical devices intended for single use (e.g., catheters).

2. **High Level Disinfectant Usage:** Omnicide solution should not be used to high level disinfect a semi-critical device when sterilization is practical.

3. **Endoscope Usage:** Omnicide solution is not the method of choice for sterilization of rigid endoscopes which the device manufacturer indicates are compatible with steam sterilization.

C. WARNINGS

Omnicide ACTIVATED DIALDEHYDE SOLUTION IS HAZARDOUS TO HUMANS AND DOMESTIC ANIMALS.

DANGER: Keep Out of Reach of Children.

CONTAINS: Glutaraldehyde

1. Direct contact is corrosive to exposed tissue, causing eye damage and skin irritation/damage. Do not get into eyes, on skin or on clothing.
2. Avoid contamination of food.
3. Use in well ventilated area in closed containers.

In case of contact, immediately flush eyes or skin with copious amounts of water for at least 15 minutes. For eyes, get medical attention.

Harmful if swallowed. Drink large quantities of water and call a physician immediately.

Probable mucosal damage from oral exposure may contraindicate the use of gastric lavage.

Emergency, safety, or technical information about Omnicide solution can be obtained from Cottrell, Ltd. Hotline at 1-800-THE EDGE (843-3343), or by contacting Intolac at 1-800-535-5053 or by contacting your Cottrell, Ltd. representative.

4

D. PRECAUTIONS

1. Appropriate hand, eye and face protection, as well as should be worn when cleaning and sterilizing/disinfecting and equipment.

More detailed information regarding the handling of the compatible materials is included in the MSDS sheet.

2. Contaminated, reusable devices MUST BE THOROUGH disinfection or sterilization, since residual contamination effectiveness of the germicide.

3. The user MUST adhere to the DIRECTIONS FOR USE will affect the safety and effectiveness of the germicide.

4. The reusable device manufacturer should provide the reprocessing procedure for that device using Omnicide.

E. DIRECTIONS FOR USE


1. **Activation**

Activate the Omnicide solution by adding the entire bottle (which is attached to the Omnicide solution or Place cap on container and shake well. On activation, changes color to blue, thereby indicating solution is ready should then be immediately tested with the PreCheck Concentration Indicator Test Strip upon activation and assure glutaraldehyde concentration is above the MEI intended for use in manual (bucket and tray) systems. ABS, polyethylene, glass-filled polypropylene or special plastics. Record the date of activation (mixing date) on the Omnicide solution container label in the space on log book or a label affixed to any secondary container solution. Omnicide must be discarded after 28 days. Glutaraldehyde Concentration Indicator Test Strip indicates.

2. **Cleaning/Decontamination**

Blood and other body fluids must be thoroughly cleared lumens, and objects before application of the disinfectant and other body fluids should be autoclaved and disposed applicable Federal, State and local regulations for infection.

5



For complete disinfection or sterilization of medical instruments and equipment, thoroughly clean, rinse and rough dry objects before immersing in OmniCide solution. Cleanse and rinse the lumens of hollow instruments before filling with OmniCide solution. Refer to the reusable device manufacturer's labeling for instructions on disassembly, decontamination, cleaning and leak testing of their equipment.

3. Usage

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING.

a. Sterilization (Bucket/Tray Manual System)

Prior to immersing medical equipment/devices, test the solution with a ProCheck G Glutaraldehyde Concentration Indicator Test Strip to assure that the glutaraldehyde concentration is above its MEC.

Immerse medical equipment/devices completely in OmniCide solution for a minimum of **ten hours at 68°F (20°C)** to eliminate all microorganisms including *Clostridium sporogenes* and *Bacillus subtilis* spores. Remove equipment from the solution using sterile technique and rinse thoroughly with sterile water following the rinsing instructions below.

b. High Level Disinfection (Bucket/Tray Manual System)

Prior to immersing medical equipment/devices, test the solution with a ProCheck G Glutaraldehyde Concentration Indicator Test Strip to assure that the glutaraldehyde concentration is above its MEC.

Immerse medical equipment/devices completely in OmniCide solution for a minimum of **45 minutes at 68°F (20°C)** to destroy all pathogenic microorganisms, except for large numbers of bacterial endospores, but including *Mycobacterium bovis* (Quantitative TB Method). Remove devices and equipment from the solution and rinse thoroughly following the rinsing instructions below.

c. Intermediate Level Disinfectant

Prior to immersing medical equipment/devices, test the solution with a ProCheck G Glutaraldehyde Concentration Indicator Test Strip to assure that the glutaraldehyde concentration is above its MEC.

OmniCide is an intermediate level disinfectant when used or reused, at full strength for a maximum of 28 days or above its Minimum Effective Concentration (MEC) as determined by the ProCheck G Glutaraldehyde Concentration Indicator Test Strip at **68°F (20°C) with an immersion time of ten minutes** for a use period not to exceed 28 days.

4. Rinsing Instructions

Following immersion in OmniCide solution, thoroughly rinse the equipment or medical device by immersing in two gallons of water. Repeat this procedure a second time with a fresh two gallon volume of water.

For endoscopic instruments with lumens, a minimum of 500 ml of water should be flushed through lumens during each separate rinse unless otherwise noted by the device or equipment manufacturer. Use fresh volumes of water for each rinse. Discard the water following each rinse. Do not reuse the water for rinsing or any other purpose as it will become contaminated with glutaraldehyde.

Refer to the reusable device/equipment manufacturer's labeling for rinsing instructions.

Sterile Water Rinse: Critical devices which are sterilized with OmniCide must be rinsed with sterile water.

Potable Water Rinse: A sterile water rinse is recommended when practical, for all devices. Alternatively, a high quality potable water (one that meets Federal Clean Water Standards at point of use) may be used.

The use of potable water for rinsing, increases the risk of contaminating the device or medical equipment with *Pseudomonas* and atypical (fast growing) mycobacteria that are often present in potable water supplies. The devices (e.g., colonoscope) need to be completely dried, because any moisture remaining provides an ideal situation for rapid colonization of bacteria. Additionally, mycobacteria are highly resistant to drying, therefore, rapid drying will avoid possible colonization but may not result in a device free from atypical mycobacteria. A final rinse using a 70 percent isopropyl alcohol solution should be used to speed the drying process and reduce the numbers of any organism present as a result of rinsing with potable water.

F. REUSE

OmniCide solution has demonstrated efficacy in the presence of soil contamination and a simulated amount of microbiological contamination. The glutaraldehyde concentration of this product during its use is determined by the ProCheck G Glutaraldehyde Concentration Indicator Test Strip. The solution is above the Minimum Effective Concentration (MEC) if the ProCheck G Glutaraldehyde Concentration Indicator Test Strip indicates a positive result. OmniCide must be discarded after 28 days, even if the ProCheck G Glutaraldehyde Concentration Indicator Test Strip indicates a positive result.

G. MONITORING OF GERMICIDE TO ENSURE SPECIFICATION


It is recommended that the OmniCide solution be tested with a ProCheck G Glutaraldehyde Concentration Indicator Test Strip prior to use to insure that the appropriate concentration of glutaraldehyde is maintained against a dilution which may lower the concentration of the MEC. During the use of OmniCide as a high level disinfectant, it is also highly recommended that a thermometer and timer be used to ensure minimum conditions are met. The pH of the activated solution should be checked to verify that the pH of the solution is between 8 and 10.

H. POST-PROCESSING HANDLING AND STORAGE OF REUSABLE DEVICES

Sterilized or disinfected reusable devices are either to be used or stored in a manner to minimize contamination. Refer to the manufacturer's labeling for additional storage and/or handling instructions.

I. STORAGE CONDITIONS AND EXPIRATION DATE

1. Prior to activation, OmniCide solution should be stored in a sealed container at controlled room temperature 59°-86°F (15°-30°C).
2. The expiration date of the unactivated OmniCide solution is indicated on the activator. The activator will be found on the side of the immediate container.
3. The use period for activated OmniCide is for no longer than 28 days following activation or as indicated by ProCheck G Glutaraldehyde Concentration Indicator Test Strip. Once activated, the solution should be further diluted prior to its usage.



J. SAFETY INFORMATION

Emergency, safety, or technical information about OmniCide solution can be obtained from Cottrell, Ltd. at 1-800-THE EDGE (843-3343), Infotrac at 1-800-535-5053, or by contacting your Cottrell, Ltd. representative.

K. USER TRAINING

The user should be adequately trained in the decontamination and disinfection or sterilization of medical devices and the handling of liquid chemical germicides. Additional information about OmniCide solution can be obtained from Cottrell, Ltd. at 1-800-843-3343, or by contacting your local Cottrell, Ltd. representative.

L. DISPOSAL INFORMATION

0.946 L (1 quart), 1.892 L (1/2 gallon) and 3.785 L (1 gallon) size container must be triple rinsed and disposed of in accordance with Federal, State or local regulations.

M. REORDER INFORMATION

Order	Description	Case Contains
OM32/28	1 quart (0.946 L)	4 qts/case
OM128/28	1 gallon (3.785 L)	4 gals/case
PCG660	ProChek G Concentration Indicator Strips	60 strips/canister (6 canisters/case)
PCG615	ProChek G Concentration Indicator Strips	15 strips/canister (16 canisters/case)

ENDOSCOPE REPROCESSING

A. INTENDED USE/DIRECTIONS FOR USE

OmniCide is a liquid chemical sterilant and high level disinfectant for flexible endoscopes when used according to the DIRECTIONS FOR USE (SEE SECTION E).

1. Germicide Level of Activity: OmniCide can be used at the following germicide levels of activity:

Flexible endoscopes, when expected to penetrate the skin or mucous membranes or are used in otherwise normally sterile tissues of the body during use are critical devices and therefore, are required to be sterile.

Sterilant: OmniCide is a sterilant for flexible endoscopes when used or reused, according to the DIRECTIONS FOR USE, at full strength for a maximum of 28 days or above its Minimum Effective Concentration (MEC) as determined by the ProChek G Glutaraldehyde Concentration Indicator Test Strip at 68°F (20°C) with an immersion time of 10 hours for a use period not to exceed 28 days.

Flexible endoscopes when expected to come in contact without penetration of mucous membranes are semi-critical devices and therefore may be high level disinfected.

High Level Disinfectant: OmniCide is a high level disinfectant for flexible endoscopes when used or reused, according to the DIRECTIONS FOR USE, at or above its Minimum Effective Concentration (MEC) as determined by the ProChek G Glutaraldehyde Concentration Indicator Test Strip at 68°F (20°C) with an immersion time of 45 minutes according to the DIRECTIONS FOR USE in SECTION E for a period not to exceed 28 days.

B. GENERAL PROCEDURE FOR HIGH LEVEL DISINFECTION OF FLEXIBLE ENDOSCOPES

(This procedure is recommended in the absence of specific directions from the device manufacturer.)

1. Trained Personnel

- a. Personnel involved in the reprocessing of endoscopes should have the ability to read, understand, and implement instructions from manufacturer's and regulatory agencies as they relate to endoscopic disinfection.
- b. The person(s) to whom the job of reprocessing endoscopes is given should have the opportunity to become completely familiar with the mechanical aspects of the endoscopic equipment.

- c. Training should include familiarization with OSHA's in-house policies on how to appropriately and safely use chemical germicides.

- d. Training should also include information on the safe use of instruments that are contaminated with body fluids should include familiarization with universal precautions.

2. Cleaning of flexible endoscopes

a. Cleaning at the Examination Room

Reflex of body fluids from the patient may occur in channels. Cleaning of endoscopes and accessories promptly after removing the endoscope from the patient and drying of secretions.

1. Personnel should don all personnel protective equipment.

2. Prepare an enzyme detergent (e.g., ProEZ) or one recommended by the manufacturer.

3. Gently wipe all debris from the insertion tube with a soft cloth.

4. Place the distal end of the flexible endoscope into the detergent solution and aspirate through the biopsy channel 5-10 seconds or until the solution is visibly clear. Rinse the detergent solution and air several times.

5. Flush or blow out air and water channels in accordance with endoscope manufacturer's instructions.


6. Transport the endoscope to the reprocessing area.

b. Cleaning at the Reprocessing Area

1. Attach any necessary water-tight caps to the electric of the umbilicus.

2. Before proceeding with any further cleaning steps, the endoscope should be leak tested. (Refer to manufacturer's instructions). Follow the manufacturer's instructions if the endoscope appears damaged.

3. Fill a sink or basin with a freshly-made enzyme (e.g., ProEZ) detergent solution.



4. Immerse the endoscope. All channels should be irrigated with copious amounts of detergent and tap water to soften, moisten, and dilute the organic debris. All detachable parts (e.g., hoods and suction valves) should be removed and soaked in the detergent solution. The insertion tube should be washed with detergent solution and rinsed.
 5. Use a small soft brush to scrub all detachable parts.
 6. Use a brush to clean under the suction valve, air/water valve and biopsy port openings.
 7. Brush the entire suction/biopsy system including the body, the insertion tube, and the umbilicus of the endoscope in accordance with the manufacturer's instructions.
 8. Accessible channel(s) should be brushed to remove particulate matter, and the detergent solution must be suctioned or pumped through all channels to remove dislodged material. Channel irrigators may be useful for this step. Fill all channels with detergent solution and soak 2-5 minutes in accordance to the ProEZ instructions.
3. **Rinse after cleaning**
 - a. Rinse the endoscope and all detachable parts in clean water.
 - b. Rinse all channels well with water to remove debris and detergent.
 - c. Purge water from all channels and wipe dry the exterior of the endoscope with a soft clean cloth to prevent dilution of the OmniCide disinfectant used in subsequent steps.
 4. **Manual Disinfection**
 - a. Activate the OmniCide solution by adding the entire content of the Activator bottle which is attached to the OmniCide solution container. Recap the container and shake well. On activation, the solution immediately changes color to blue, thereby indicating the solution is ready. Use ProChek G Glutaraldehyde Concentration Indicator Test Strips after activation to determine that the solution is above the Minimum Effective Concentration (MEC) before use.
Test the activated OmniCide solution using the ProChek G Glutaraldehyde Concentration Indicator at 68°F (20°C), before each use.
 - b. Attach channel irrigators/adapters and cover the biopsy port in accordance with the manufacturer's instructions.


- c. Pour the activated OmniCide into an appropriate sized basin.
 - d. Completely immerse the endoscope in the basin of OmniCide.
NOTE: In order to prevent damage to the endoscope, DO NOT soak any other accessory equipment with the endoscope.
 - e. Inject the OmniCide solution into all channels of the endoscope until it can be seen exiting the opposite end of each channel. Assure that all channels are filled with disinfectant and that no air pockets remain within the channels.
 - f. Cover the disinfectant soaking basin with a tight fitting lid to minimize chemical vapor exposure.
 - g. Soak the endoscope for 45 minutes. Use a timer to ensure adequate soaking time.
 - h. Before completely removing the endoscope from the disinfectant, flush all channels with air to remove disinfectant.
5. **Rinse After Manual Disinfection**
 - a. **Rinse 1:** Fill a basin with two gallons of water (preferably sterile water). Place the endoscope into the basin and thoroughly rinse the exterior of the scope. Attach channel irrigators/adapters to the endoscope and flush with 500 ml of water through the channel irrigator. Empty basin.
 - b. **Rinse 2:** Fill a basin with two gallons of water (preferably sterile water). Place the endoscope into the basin and thoroughly rinse the exterior of the scope and flush with 500 ml of water through the channel irrigator.
 - c. Purge all channels with air.
 - d. Flush all channels with 70% alcohol until the alcohol can be seen exiting the opposite end of each channel.
 - e. Purge all channels with air.
 - f. Remove all adapters and devices.

6. Storage

- a. Dry the exterior of the endoscope with a soft cloth (pre). DO NOT attach detachable parts to the endoscope. Storage of endoscopes with the removable parts at risk of trapping liquid inside the instrument and for drying of the channels and channel openings. To prevent waterborne organisms, the endoscope and all distal parts should be thoroughly dried prior to storage.
- b. Hang the endoscope vertically with the distal tip in a ventilated, dust-free cabinet.

REFERENCES:

1. ASTM F 1518-94, Standard for Cleaning and Disinfection of Flexible Fiberoptic for the Examination of the Hollow Vessels. Current edition approved May 15, 1994. Ph.
2. Martin, M.A., MD, Reichelderfer, M., APIC guidelines for infection prevention and Association for Professionals in Infection Control and Epidemiology, Inc., A.I.C. A Control, 1994; 22: 19-30.
3. Vesley, D. et al. Significant factors in the disinfection and sterilization of flexible. Journal of Infection Control, December 1992, pp. 297.
4. Axon, A.T.R., Bond, B., Botwell, P.M., Cowen, A.E., Fleisher, D.E. and Tandon, R.K. disinfection. Working Party Reports, Blackwell Scientific Publications, 1990, 46-5.



ULTRASOUND TRANSDUCER REPROCESSING
(Endocavity, Endovaginal, Endorectal, Transesophageal, etc.)
(Refer to the device manufacturer's instruction for specific procedures)

A. INTENDED USE/DIRECTIONS FOR USE

OmniCide activated dialdehyde solution is a liquid chemical sterilant, high level disinfectant and intermediate level disinfectant for ultrasound transducers when used according to the DIRECTIONS FOR USE.

1. **Germicide Level of Activity:** OmniCide can be used at the following levels of activity:

Sterilant: OmniCide is a sterilant for ultrasound transducers when used or reused according to the DIRECTIONS FOR USE, at or above its Minimum Effective Concentration (MEC) as determined by the ProChek G Glutaraldehyde Concentration Indicator Test Strip at 68°F (20°C) with an immersion time of ten hours for a use period not to exceed 28 days.

High Level Disinfectant: OmniCide is a high level disinfectant for ultrasound transducers when used or reused, according to the DIRECTIONS FOR USE, at or above its Minimum Effective Concentration (MEC) as determined by ProChek G Glutaraldehyde Concentration Indicator Test Strip at 68°F (20°C) with an immersion time of 45 minutes according to the DIRECTIONS FOR USE in SECTION E for a use period not to exceed 28 days.

Intermediate Level Disinfectant: OmniCide is an intermediate level disinfectant when used or reused according to the DIRECTIONS FOR USE, at or above its Minimum Effective Concentration (MEC) as determined by the ProChek G Glutaraldehyde Concentration Indicator Strip at 68°F (20°C) with an immersion time of 15 minutes according to the DIRECTIONS FOR USE in SECTION E for a use period not to exceed 28 days.

B. GENERAL PROCEDURE FOR HIGH LEVEL DISINFECTION OF ULTRASOUND TRANSDUCERS (Endocavity, Endovaginal, Endorectal, Transesophageal, etc.)

(Refer to the device manufacturer's instruction for specific procedures)

1. Trained Personnel

- a. Personnel involved in the reprocessing of ultrasound transducers should have the ability to read, understand, and implement instructions from manufacturers and regulatory agencies as they relate to transducer disinfection.

15

- b. The person(s) to whom the job of reprocessing ultrasound transducers is given should have the opportunity to become completely familiar with the mechanical aspects of the ultrasound equipment.
- c. Training should include familiarization with OSHA regulations and in-house policies on how to appropriately and safely handle liquid chemical germicides.
- d. Training should also include information on the safe handling of instruments that are contaminated with body fluids after use. This should include familiarization with universal precautions.

2. Cleaning of Transducers

Cleaning of transducers and accessories should be performed promptly after patient use to prevent drying of secretions.

- a. Personnel should don all personnel protective equipment.
- b. Prepare an enzyme detergent (e.g., Pro EZ) or one recommended by the transducer manufacturer.
- c. Gently wipe all debris from transducer surfaces with a moistened gauze or the like.
- d. Immerse the distal end of the transducer into the water and enzyme detergent solution for 2-5 minutes or as recommended by the device manufacturer.

NOTE: Refer to manufacturer's instructions for recommended immersion depth.

- e. Rinse all surfaces with water and dry by wiping with a soft clean cloth.

3. Manual Disinfection

Caution: OmniCide may discolor the exterior of some transducers, however, the acoustic or scanning performance is not impaired. Check with the device manufacturer for specific recommendations before proceeding.

- a. Activate the OmniCide long life activated dialdehyde by adding the entire content of the Activator bottle which is attached to the OmniCide solution container. Recap the container and shake well. On activation, the solution immediately changes color to blue, thereby indicating the solution is ready. Use ProChek G Glutaraldehyde Concentration Indicator Test Strips after activation to determine that the solution is above the Minimum Effective Concentration (MEC) before use.

16

Test the activated OmniCide solution using the ProChek G Glutaraldehyde Concentration Indicator at 68°F (20°C) each use.

- b. Pour the activated OmniCide into an appropriate size container.
- c. Immerse the distal end of the transducer into the OmniCide solution.

Note: Refer to manufacturer's instructions for recommended immersion depth.

- d. Cover the disinfectant soaking basin to minimize chemical vapor exposure.
- e. Soak the transducer for 45 minutes. Use a timer to ensure adequate soaking time.

4. Rinse After Manual Disinfection

- a. **Rinse 1:** Fill a basin with a minimum of 1 gallon of water (preferably sterile water). Place the transducer into the recommended depth and allow to soak for 5 minutes.
- b. **Rinse 2:** Holding the transducer over the basin, flush with water (preferably sterile water) for one minute.
- c. Dry the transducer by wiping with a soft (preferably sterile) cloth.

5. Storage

Refer to manufacturer's recommendation.

ACCEPTED
MAR 4 1996

Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the pest registered under EPA Reg. No. 46851-2