

FRONT PANEL



Omnicide 14*

STERILIZING and DISINFECTING SOLUTION

For use in hospitals, dental offices, nursing homes, and health care institutions.

14 Day Activated Dialdehyde

Cat. No. U9006

- Sporicidal
- Virucidal†
- Fungicidal
- Tuberculocidal
- Bactericidal
- Pseudomonacidal

Qty. 1 Gallon
(3.79 liters)

Non-Rusting

Complete Sterilant


Pleasant Fragrance

E.P.A. Establishment No. 7388-CA-3

E.P.A. Registration No. 46851-4

KEEP OUT OF REACH OF CHILDREN
DANGER

ADM MEDICAL DIVISION
25283 CABOT Laguna Hills, Ca 92653



Clear Material (Reverse)

46851-4

PM 31

1063

LEFT PANEL

DIRECTIONS FOR ACTIVATION AND USE:

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT ITS WITH LABELING.

To activate solution — add activator contents to this one gallon container — recap — shake. Activated solution changes color to green. Do not filter. Always clean away excess organic debris, rinse and rough dry instruments or equipment before immersion. Omnicide 14™ may be USED AND REUSED in a manual system as a disinfectant for up to 34 days after activation. For Tuberculocidal activity use the sterilization procedure or employ a previously unused solution for the disinfection procedure.

FOR DISINFECTION:

Immersion articles for a minimum of 10 minutes to destroy vegetative pathogens including — *Escherichia coli* 0113, *Staphylococcus aureus*, *Salmonella typhosa*, *Pseudomonas aeruginosa*, *Mycobacterium tuberculosis*, fungi including — *Trichophyton interdigitale*, viruses including — Herpes simplex virus types 1 & 2, Polio virus types 1 & 2, Vaccinia virus, Influenza virus types A1 (New Jersey) and B (Maryland). The ability of this product to destroy the etiologic agent for viral hepatitis has not been established or demonstrated.

Rinse articles thoroughly with sterile water after disinfection or sterilization.

For Complete Sterilization — to destroy resistant pathogenic spores including — *Bacillus subtilis* and *Clostridium sporogenes* — immerse articles for a minimum of 10 hours. Remove instruments using aseptic technique and rinse thoroughly in sterile water.

All instruments and devices must be thoroughly cleaned to remove excess organic debris, rinsed and rough dried before immersion in Omnicide 14. Hypodermic needles cannot be effectively disinfected or sterilized with Omnicide 14.

REUSABLE

Dating Record
OMNICIDE 14
HAS A 14-DAY SHELF
LIFE AFTER
ACTIVATION.



ACTIVATED _____

EXPIRATION _____

2.5

ACCEPTED

JAN 1

10.5 —

AMERICAN DENTAL
Art #: 86752
Sleeve Size: 10.5 x 4.5
Art Size: 10.25 x 4.125

46851-4

23

Recommended for Disinfection or Sterilization of: MEDICAL INSTRUMENTS, CARBON STEEL CUTTING INSTRUMENTS, BURS, ALUMINUM, GLASS, PLASTIC or RUBBER ARTICLES, THERMOMETERS, MIRRORS, RESPIRATORY EQUIPMENT.

STORAGE. Keep Omnicide 14 covered in storage or use. Store in a cool place. Do not contaminate water, food or feed by storage or disposal.

DISPOSAL. Triple rinse container (or equivalent) then puncture and dispose of in a sanitary landfill or incinerate. Used disinfectant/sterilant solutions are acutely hazardous. Store, use, and dispose of according to label instructions.

8-85

PRECAUTIONARY STATEMENTS:

HAZARDS TO HUMANS AND DOMESTIC ANIMALS. DANGER. Corrosive to tissue. Causes burns. Causes eye damage and skin irritation. Do not get in eyes, on skin or

on clothing. May be absorbed through skin. Wear goggles or face shield and rubber gloves when handling.

HARMFUL IF SWALLOWED. Avoid contact with food. Use in ventilated area. Use containers which can be tightly closed during use or storage of solution.

STATEMENT OF PRACTICAL TREATMENT

Eyes — in case of eye contact, immediately flush thoroughly with water. **Skin** — in case of skin contact, immediately flush thoroughly with water. **Ingestion:** If swallowed, drink milk, egg whites, gelatin solution or, if these are not available, large quantities of water. **Avoid alcohol.** Call a physician.

Note to Physician: Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed.

4.5

10.5

Printed)

ACCEPTED

JAN 1

USEPA Reg No. 4685-10

Give
Print
GR