US ENVIRONMENTAL PROTECTION AGENCY OFFICE OF PESTICIDES PROGRAMS REGISTRATION DIVISION (75-767) WASHINGTON, DC 20480	EPA REGISTRATION NO. 46851-5	16 301 1967
	TERM OF ISSUANCE	
NOTICE OF PESTICIL : L REGISTRATION	NAME OF PESTICIDE PRODUCT	
(Under the Federal Ins. icide, rungicide,	OMET-II	
end Rodenticide Act, as amended)		
NAME AND ADDRESS OF REGISTRANT (Include ZIP code)		
r		
ADM Medical pivision, Inc.		
25283 Cabot Road, Suite 106 Laguna Hills, CA 92653		
1		
-	J	
	,	
NOTE: Changes in labeling formula differing in substance from that accepted in connection with this registration must be submitted to and accepted by the Registration Division prior to use of the label in commerce, in any correspondence on this product always refer to the above U.S. EPA registration number.		
On the basis of information furnished by the registrant, the above named peaticide is hereby Registered/Reregistered under the Federal Insecticide, Fungicide, and Rodenticide Act.		
A copy of the labeling accepted in connection with this Registration/Reregistration is returned herewith.		
Registration is in no way to be construed as an indorsement or approval of this product by this Agency. In order to protect health and the environment, the Administrator, on his motion, may at any time suspend or cancel the registration of a pesticide in accordance with the Act. The acceptance of any name in connection with the registration of a product under this Act is not to be construed as giving the registrant a right to exclusive use of the name or to its use if it has been covered		
by others. This product is conditionally registered in accordance with FirkA		
section 3(c)(7)(A) provided that you:  1. Submit/cite all data required for registration/reregistration		
of your product under FIFRA section 3(c)(5) when the Agency requires all		
registrants of similar products to submit such data.		
<ol> <li>The results of 12 months Storage and Stability data must be received by the Agency within 2 years after the granting of this</li> </ol>		
registration.		
3. Add the phrase "LPA Registration No. 46851-5" to your label . before you release the product for shipment.		
4. Submit two (2) copies of your final printed labeling before you		
release the product for shipment, kefer to the A-79 enclosure for a further description of final printed labeling.		
A copy of your proposed label with our acceptance stamp is included		
with this registration with the understanding that you will comply with the conditions in the attachment to the registration notice.		
The registration notice.		
HIM I		
Jery Kempter		
Product Manager (32)		
Disinfectants branch Registration Division (TS-767C)		
Łnclosures	-	
	DESI A	VAILABLE COPY
ATTACHMENT IS APPLICABLE		
SIGNATURE OF APPROVING OFFICIAL		DATE
	·	

92153:I:Pringle:K-9:KENCO:7/13/87:7/22/87:cdb:vo:dag:cdb

# INSTRUMENT & EQUIPMENT DISINFECTANT

## Non-Rusting

### Pleasant Fragrance

Active Ingredients

100.0

- \* Bactericidal
- \* Fungicidal destroys pathogenic fungi
- \* Virucidal\*
- \* Tuberculocidal
- Viruses including Herpes virus types 1 & 2, Polio virus types 1 & 2, Vaccinia virus, Influenza virus type A1 (New Jersey) and influenza type 8 (Maryland).

# KEEP OUT OF REACH OF CHILDREN CAUTION:

See handling precautions and practical treatments on back panel.

Contents: 1/2 gallon (64 fluid oz.)

Manufactured & Distributed by:
ADM MEDICAL DIVISION
25283 Cabot Road Laguna Hills, CA 92653

## DIRECTIONS FOR USE:

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING

OMC II is ready to use for disinfection of medical and dental instruments and devices. Rinse thoroughly after immersion. Other hard surface areas such as walls, floors and countertops should be mopped, sponged, sprayed or wiped in a manner to thoroughly wet the surface with OMC II. Discard the used solution daily, or sooner if it becomes diluted or soiled.

Crean mll excess organic debris from items to be disinfected and rough dry, immerse in OMC II for 10 minutes at 20°C to destroy bacteria including Escherichia coli 0113, Staphylococcus aureus, Salmonella typhosa, Pseudomonas meruginosa, Mycobacterium tuberculosis (if tuberculosis contamination is suspected, increase immersion time to 20 minutes, for added safety); fungi including Trichophyton interdigitale; viruses including Herpes virus types 1 & 2. Polio virus types 1 & 2. Vaccinia virus, influenza virus type Al (New Jersey) and Influenza B (Maryland).

for use on metals, plastics, rubber and glass. Glass may etch over prolonged or repeated exposure. The control of the control

JUL 1 6 1987

8006-080 Hother the Federal Insecticite, Progress, and Radenticide Act Progress, and Radenticide Act Progress, and Radenticide Act Progress, and Radenticide as accorded, for the Penticide

BEST AVAILABLE COPY

When diluted as directed, OMC II is stable. All label claims of effectiveness are based on previously unused batches of OMC II diluted as directed.

OMC II should not be used for disinfection of hypodermic needles or any items with deep and/or narrow crevices.

### DISPOSAL

Do not reuse empty container. Rinse container and discard in trash collection. Waste from the use of this product may be disposed of on site or on approved waste disposal facility.

PRECAUTIONARY STATEMENTS:

CAUTION: Hay cause eye irritation. Avoid contact with ......

STATEMENT OF PRACTICAL TREATMENT:

If in eyes: Flush with plenty of water. Get medical attention if irritation persists.

NOTE TO PHYSICIAN:

Probable aucosal damage may contraindicate the uce of gastric lavage.