

LIFEEX SURFACTANT IODINE GERMICIDE

\*\*\*\*\*  
\*\*  
\*\*  
\*\*  
\*\* For Disinfecting Instruments and Devices \*\*  
\*\*  
\*\*\*\*\*

.....  
• Contains .....  
• 100 ppm .....  
• Titratable .....  
• Iodine .....  
.....

READY TO USE

FOR DISINFECTING AND CLEANING  
HARD SURFACE INSTRUMENTS AND DEVICES

NON-STAINING and ODORLESS

RECORDED  
17 SEP 1982  
47923-4  
FOR THE FEDERAL BUREAU OF INVESTIGATION  
AND OTHER AGENCIES FOR ALL  
FEDERAL BUREAU OF INVESTIGATION  
RECORDS

Active Ingredients:  
Nonylphenoxypoly(ethyleneoxy)ethanol-iodine complex... 0.163%  
Inert Ingredients..... 99.837%

\* Providing 100 ppm titratable iodine.  
EPA REG. NO. 47923-U

Keep Out Of Reach Of Children  
DANGER

See bottom panel for additional precautions and first aid.

For use as a disinfectant of non porous hard surfaces and articles in:

- \* HOSPITALS, CLINICS, MEDICAL, AND DENTAL OFFICES.
- \* NURSING HOMES AND OTHER HEALTH CARE FACILITIES.
- \* FOOD AND DRUG PROCESSING FACILITIES.
- \* RESEARCH FACILITIES.
- \* INSTITUTIONS.
- \* VETERINARY, LIVESTOCK, AND ANIMAL HANDLING FACILITIES.

WIZ INC.  
U.S.A.

FOR INDUSTRY AND PROFESSIONAL USE  
Only for Sale to, Use, and Storage by Servicepersons.

Net Contents 1 Gallon

4149 LAB. DATA LIFEX 2, EPA SYMBOL 4747-1, Page 1  
Parsons Paper Label Text - September 8, 1982  
L I F E X

For all hard surfaces of hard and animal contact, instruments, and devices, not for medical, dental, and veterinary instruments, test and laboratory equipment, electronic equipment, contacts, probes, syringes, thermometers, needles and pliers, and all other hard surface containers, dispensers, and applicators of hormonal drugs, fluids, and all other research and equipment containers and all other hard surfaces.

Directions For Use

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.  
Use frequently as instruments decontaminant. FOR SMALL INSTRUMENTS: (1) Rinse all instruments such as syringes, (2) Pre-clean the instruments in detergent and water and remove all adhering water or (such as blood, heavy dirt, and contaminants) by wiping or scrubbing. Then, rinse with water. (3) Soak instruments in LIFEX 2 bath for 10 minutes or until next use. (4) Rinse instruments with sterile or potable water. The LIFEX 2 bath container should remain covered, and the solution should be changed daily or more often if needed. Change the solution immediately if it becomes cloudy, dirty, or a different color. FOR LARGE INSTRUMENTS: (1) Apply LIFEX 2 by spray or other suitable applicator so as to thoroughly wet the surface. (2) Wipe surface with clean sponge, paper towel, or disposable wipe to disperse the LIFEX 2 and to dislodge any surface dust and dirt. Do not reuse sponges, paper towels, or disposable wipes with potable water - discard rinse water when dirty. (3) Rubby LIFEX 2 and saturate surface so that it remains wet for a minimum of 10 minutes. (4) Rinse with sterile or potable water. For dirty or extra dirty surfaces, pre-clean the surface with water and surface detergent, and rinse. Then, proceed as above, beginning with step 1. LIFEX 2 is light-moderate tea color, and the solution should not be used if there is a noticeable change in color. Keep container cap tightly closed. Avoid contact with starch composition materials, as the iodine will turn blue - which is removable with soap and water.

Effectiveness

Iodophors yielding free iodine are established germicides of viruses, bacteria, and fungi and, the LIFEX 2 composition has been determined effective on: Salmonella, Tuberculosis, Escherichia coli, Streptococcus, Pseudomonas aeruginosa, Trichophyton meningitidis, Influenza A2 (Hong Kong), and Herpes Simplex.

Disposal

Do not reuse container except to refill it with LIFEX 2. Crush or cut empty container and discard in trash collection; or, rinse container thoroughly with water before discarding.

PRECAUTIONARY STATEMENTS

Hazard to Humans and Domestic Animals

Danger. Corrosive. Causes eye damage. Causes skin irritation. Protect eyes and skin when handling. Do not get in eyes, on skin, or on clothing. Harmful if swallowed. Avoid contamination of food.  
FIRST AID: In case of contact with eyes and skin, flush immediately with plenty of water. For eyes, call a physician. Remove and wash contaminated clothing before reuse. If swallowed, drink promptly a large quantity of milk, egg whites, gelatin solution, or if these are not available, drink large quantities of water. Avoid alcohol. Call a physician immediately.

- USE LIFEX 1 - For Disinfecting Rooms and Equipment (150 PPM titratable iodine).  
LIFEX 2 - For Disinfecting Instruments and Devices (1000 PPM titratable iodine).  
LIFEX 3 - For Sanitizing Food and Drug Equipment (25 PPM titratable iodine).

LIFEX 2 is manufactured in the UNITED STATES OF AMERICA by WIZ CHEMICAL CORPORATION. Copyright 1982, WIZ CHEMICAL CORPORATION, P.O. BOX 687, BALA CYNWYD, PA. 19380. LIFEX 2 is a trademark owned by WIZ CHEMICAL CORPORATION.