

UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

OCT 25 1989

Inhalation Plastics, Inc.
7790 North Merimac Avenue
Niles, IL 60648


Attention: James Lekkas
Manager Q.A.

Gentlemen:

Subject: Solution II Concentrate
EPA Registration No. 48026-2
Your Submission Dated October 4, 1989

The labeling referred to above, submitted in connection with registration under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, is acceptable. A stamped copy is enclosed for your records.

Sincerely yours,



John H. Lee
Product Manager (31)
Antimicrobial Program Branch
Registration Division (H7505C)

003317E/CIA-ATL-89-16;KEMCO:10/20/89:12/02/89:dg:sw:vo:ek:dg

CONCURRENCES

SYMBOL							
SUBJECT							

252

SOLUTION II CONCENTRATE when used as recommended, is effective

as a disinfectant against *Staphylococcus aureus*, *Salmonella choleraesuis* and *Pseudomonas aeruginosa*.

Virucidal Activity* This product when used on environmental, inanimate hard surfaces at 1 ounce per 1 gallon of water is effective against influenza A; herpes simplex; Adenovirus Type 2; and Vaccinia Virus.

For disinfection of previously cleaned surgical instruments, respiratory devices and dental equipment such as scissors, scalpels, mouthpieces, plastic tubing, simply submerge instruments into solution containing 1 ounce of **SOLUTION II CONCENTRATE** per 1 gallon of water for 10 minutes.

This product has not been tested for effectiveness against *Mycobacterium tuberculosis* and must not be relied upon when a tuberculocidal product is desired.

STORAGE AND DISPOSAL

DO NOT CONTAMINATE WATER, FOOD OR FEED BY STORAGE OR DISPOSAL
OPEN DUMPING IS PROHIBITED
DO NOT REUSE EMPTY CONTAINER

PESTICIDE DISPOSAL - Pesticides that cannot be used or chemically reprocessed should be disposed of in a landfill approved for pesticides

CONTAINER DISPOSAL - Triple rinse (or equivalent); and dispose in an incinerator or landfill approved for pesticide containers.

IP
SOLUTION II
CONCENTRATE
GERMICIDE

Cat. No. 0212

DISINFECTANT DEODORIZER VIRUCIDE*

ACTIVE INGREDIENTS:

n-Alkyl (60% C₁₂, 30% C₁₄, 5% C₁₇, 5% C₁₈) dimethyl benzyl ammonium chloride 5%
n-Alkyl (68% C₁₂, 32% C₁₄) dimethyl ethylbenzyl ammonium chloride 5%

INERT INGREDIENTS: 90%
TOTAL 100%

EPA REG NO 48026-2
EPA EST. NO. 48026-IL-1

**ONLY FOR SALE TO, USE, AND STORAGE BY SERVICE PERSONS
KEEP OUT OF REACH OF CHILDREN**

DANGER

See side panel for additional precaution statements

**CONCENTRATED
DILUTE BEFORE USE**

Manufactured by:
INHALATION PLASTICS, INC.
7790 N. MERRIMAC AVE.
NILES, IL 60648

NET CONTENTS 8 FL. OZ.

DIRECTIONS FOR GENERAL CLASSIFICATION

It is a violation of Federal product in a manner including labeling.

SOLUTION II CONCENTRATE designed specifically for processing plants, dairies, and other institutions where

is of prime importance

DISINFECTION - To disinfect cleaned hard inanimate surfaces (walls, floors, table tops), apply **CONCENTRATE** with mop or mechanical sprayer so as to allow to remain wet for 10 minutes and air dry.

DILUTION RATES: General Hospital Disinfection - add 1 ounce of water.

For disinfection, a pre-rinse is required and recommended solution for each use.

PRECAUTIONARY STATEMENTS

Hazard to Human and Domestic Animals
DANGER

Keep Out of Reach of Children. Avoid eye and skin damage. Do not use on clothing. Harmful if swallowed. Wear goggles or face shield when handling. Avoid contact with skin. Remove and wash contaminated clothing.

STATEMENT OF PRACTICE

In case of contact, immediately flush skin with plenty of water for at least 15 minutes. For eyes, call a physician. For all contaminated clothing, remove and wash thoroughly. If swallowed, drink milk, egg solution or if these are not available, large quantities of water. Call a physician.

NOTE TO PHYSICIAN

Probable mucosal damage may result from the use of gastric lavage. In case of circulatory shock, respiration, or convulsion may be needed.

OCT 23 1982

48026-2