

PM 31

48920-1

1071

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS Causes severe eye damage and skin irritation. Do not get in eyes, on skin or clothing. Wear glasses when handling. Harmful if swallowed. ENVIRONMENTAL HAZARDS Do not discharge concentrate into lakes or streams. PHYSICAL HAZARD Do not use near heat or open flames.

STATEMENT OF PRACTICAL TREATMENT In case of contact immediately flush eyes with plenty of water for at least 15 minutes. If irritation persists see a physician. Wash from skin with water. Wash contaminated clothing before reuse. If swallowed, drink promptly a large quantity of gelatin solution or, if not readily available, drink a large quantity of water. Avoid alcohol. Call a physician immediately. Wash thoroughly with soap and water after handling.

STORAGE AND DISPOSAL: Do not contaminate water, food or feed by storage or disposal. Store only in a tightly closed original container, in an area inaccessible to children and away from heat and open flames. Do not reuse container. Triple rinse container. Then offer for recycling, or puncture and dispose of in a sanitary landfill, or incineration, or, if allowed by state and local authorities, by burning. If burned stay out of smoke.

MFG BY IABCO MANUFACTURING INC P O BOX 519 LAKE CHARLES, LA 70602

NET CONTENTS

HERPICIDE SPRAY

HOSPITAL DISINFECTANT / CLEANER

BACTERICIDAL\* VIRUCIDAL\* FUNGICIDAL\* TUBERCULOCIDAL

Ready to use

Kills Herpes simplex type 1 and type 2 viruses on hard, inanimate environmental surfaces in 1 minute. It provides 24 hours of residual virucidal activity against Herpes simplex type 1 and type 2, in the presence of moisture, if not removed from the surface. Kills human immunodeficiency virus Type 1 (HIV-1)(AIDS virus) in one minute in the presence of 5% blood serum.

Active ingredients:

Alkyl (C <sub>14</sub> 58%, C <sub>16</sub> 28%, C <sub>12</sub> 14%) dimethyl benzyl ammonium chloride	2.0%
Isopropanol	40.0%
Inert ingredients	58.0%
Total	100.0%

When applied as directed to inanimate environmental surfaces and left on the surface for 10 minutes, HERPICIDE SPRAY is also Virucidal against Influenza A2/Japan, Influenza B, Adeno virus type 2, Polio type 1, and Vaccinia V506. Bactericidal against Pseudomonas aeruginosa, Staphylococcus aureus, Salmonella choleraesuis. Tuberculocidal against Mycobacterium tuberculosis bovis BGG on pre-cleaned surfaces. Fungicidal against Trichophyton mentagrophytes. Also effectively controls mold and deodorizes hard surfaces.

CONTAINS ONLY BIODEGRADABLE MATERIALS AND SURFACTANTS TO AID CLEANING.

CAUTION: KEEP OUT OF REACH OF CHILDREN

SEE SIDE PANEL FOR ADDITIONAL PRECAUTIONARY STATEMENTS.

EPA Reg. No.- 48920-1. EPA Est. No - 48920-LA-001

PATENTED COMPOSITION - US 4,797,420

MULTIPURPOSE BROAD SPECTRUM HOSPITAL CLEANER / DISINFECTANT FOR IMMERSION, ULTRASONIC AND SPRAY APPLICATION ON INANIMATE OBJECTS. RECOMMENDED FOR USE IN HOSPITALS, MEDICAL & DENTAL OFFICES, NURSING HOMES, VETERINARY CLINICS, INSTITUTIONS, HOMES, PUBLIC RESTROOMS AND ALL OTHER TREATED AREAS/SURFACES WHERE THE CONTROL OF CROSS-CONTAMINATION IS IMPORTANT.


APPLICATIONS: For use in hospital operating rooms, emergency rooms, laboratories, therapy units, isolation areas, neonatal care areas, infant incubators, transport equipment, oxygen hoods, anesthesia and therapy equipment, tables and high risk areas where cleaning and decontamination are important. Use to clean and decontaminate surgical and dental equipment and instrumentation prior to sterilization. level disinfection.

DIRECTIONS FOR USE: It is a violation of the federal law to use this product in a manner inconsistent with its labeling. This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

SPRAY APPLICATION: For Hospital Disinfection of large equipment and surfaces. Remove gross filth from the surface prior to application by spraying and wiping with a paper towel. Discard towel. Spray surface until thoroughly wet and maintain a wet surface for 10 minutes at room temperature then dry. Tuberculocidal Disinfection, pre-clean the surface by spraying and wiping clean with a paper towel. Discard towel. Re-spray, maintain a wet surface for 10 minutes at room temperature and allow to dry. For only Herpes simplex type 1 and type 2 viruses, spray the surface until thoroughly wet, wait 1 minute then air dry. Follow hospital directions for mold control and to deodorize hard surfaces.

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV OF SURFACES /OBJECTS SOILED WITH BLOOD/BODY FLUIDS - Personal Protection - Wear disposable latex gloves, gowns, mask and eye goggles when handling HIV contaminated fluids. Cleaning Procedure - Blood and other body fluids must be thoroughly cleaned from the surfaces and objects before application of the disinfectant. Thoroughly spray the surfaces and objects with the disinfectant and keep for 1 minute. Contact Time - This contact time will not control other viruses and bacteria. Maintain surface for 10 minutes for broader control. Disposal of Infectious Materials - Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

ULTRASONIC AND IMMERSION SOLUTION: To thoroughly clean and decontaminate instruments prior to high level disinfection/sterilization, clean and rinse instruments under cool running water to remove gross debris. Allow to drain and rough dry. Fill a container with the disinfectant that will allow complete submersion of the instruments. Immerse instruments into the solution and allow to remain for 10 minutes at room temperature or activate ultrasonic unit for 10 minutes. Remove and rinse prior to sterilization/high level disinfection. Change solution daily or more often if the solution becomes dirty or diluted.



181