

Actril™

CAT. NO. 78270

COLD STERILANT

FOR USE IN HEALTH CARE ENVIRONMENTS

Actril™ Cold Sterilant is intended for use in Health Care Environments. When used as described in the enclosed package insert, it can perform as a Sterilant, broad spectrum disinfectant, hospital disinfectant, cleaner-sanitizer, germicidal spray disinfectant. As a sterilant Actril™ is reusable, sporicidal, virucidal, fungicidal, tuberculocidal, and will kill Hepatitis B and AIDS viruses.

ACTIVE Ingredients

Hydrogen Peroxide 0.60%

Peroxyacetic Acid 0.06%

INERT Ingredients 99.14%

Total 100.00%

READ PACKAGE INSERT BEFORE USE:

NET CONTENTS:

1 U.S. Gallon (3.78 LITERS)

DANGER

KEEP OUT OF REACH OF CHILDREN

EPA Reg. No. 52252-7

EPA Est. No. 52252-MN-01

WARNING: Hazard to humans and domestic animals. Corrosive. Causes eye damage. Do not get in eyes. Wear goggles or face shield and rubber gloves when handling. Wash thoroughly after handling.

TREATMENT: In case of contact with eyes, flush with large amounts of water for at least 15 minutes. Get prompt medical attention. Note to Physician: Toxicologic damage may contraindicate the use of gastric lavage.

STORAGE: Store in original container. Do not expose to direct sunlight. Maintain temperature below 75°F (24°C). Avoid contact with combustible materials. Store in original closed container. For chemical emergency, spill, leak, fire, exposure and accident, call Chemtrec, day or night (800) 424-9300, (262) 455-7616.

DISPOSAL: Do not contaminate water, food or feed by storage or disposal. Wastes resulting from the use of this product may be disposed of on-site in a sanitary sewer or at an approved waste disposal facility. Triple rinse empty container with water. Then offer for recycling or reconditioning, or puncture and dispose of in a sanitary landfill, incinerate, or if allowed by state and local authorities dispose of by burning. If burned, stay clear of the smoke.



Minntech Corporation
14905 28th Avenue North
Mpls. Minnesota 55447 U.S.A.
Telephone (612) 553-3300
(800) 328-3345
Telex 29-0825-RSAMPLS PLDM
FAX 612 553-3387

87187-000/C

NOT RECOMMENDED FOR USE

10/1/90

Actril™

COLD STERILANT

For Use In Health Care Environments

PACKAGE INSERT

Active Ingredients

Hydrogen Peroxide	0.80%
Peroxyacetic Acid	0.06%

Inert Ingredients 99.14%

Total 100.00%

KEEP OUT OF REACH OF CHILDREN

DANGER

EPA Reg. No.: 52252-7

EPA Est. No.: 52252-MN-01

PRECAUTIONARY STATEMENTS

Hazard to humans and domestic animals.

Eyes: Corrosive, causes eye damage. Do not get in eyes. Wear goggles or face shield and rubber gloves when handling. Wash thoroughly after handling.

Statement of Practical Treatment

In case of contact with eyes, flush with large amounts of water for at least fifteen (15) minutes. Get prompt medical attention.

Note to Physician: Probable mucosal damage may contraindicate the use of gastric lavage.

PRODUCT CLAIMS

- Sterilant * undiluted 11 hrs. 20°C
- Sporicidal * undiluted 11 hrs. 20°C
- Bactericidal * undiluted 11 hrs. 20°C
- Virucidal * undiluted 11 hrs. 20°C
 - Hepatitis B * undiluted 11 hrs. 20°C
 - HIV * undiluted 11 hrs. 20°C
- Fungicidal * undiluted 11 hrs. 20°C
- Tuberculocidal * undiluted 10 min. 20°C
- Broad Spectrum Disinfectant * undiluted 10 min. 20°C
- Hospital Disinfectant and Pseudomonacidal * diluted 35X 30 min. 20°C
 - Use Fresh Solution
 - Do Not Reuse
- Cleaner-Sanitizer (non-food contact surfaces) * 50X dilution 10 min. 20°C
 - Do Not Reuse
- Germicidal Spray Disinfectant * undiluted 30 sec. 20°C
- Non-Staining
- Pre-Activated/No Activation Needed
- Fast-Acting
- Reusable
 - As a sterilant for 30 days
 - As a broad spectrum disinfectant for 30 days (when used undiluted)

DIRECTIONS FOR USE

1. It is a violation of Federal (U.S.A.) Law to use this product in a manner inconsistent with its labeling.
2. Do not use after Expiration Date.
3. Do not allow Actril™ Cold Sterilant to mix with alkaline substances such as bleach (Sodium hypochlorite) or other oxidizing agents.
4. Use purified water (e.g. deionized) for making dilutions, cleaning and rinsing.
5. Some materials may be incompatible with Actril™. Test material prior to use.
6. Reuse of diluted Actril™ is not recommended.

DIRECTIONS

1. Remove any material from sterilized, polished by detergent clean rinse.
2. Immerse the item in sufficient volume of Cold Sterilant for all passages. Hold in the minimum of 11 (68°F).
3. Remove items with sterile water level of 10 ppm Actril™ Residue.
4. The solution is stable up to 30 days hours immersion.

DIRECTIONS BROAD SPECTRUM

1. For broad spectrum such as dental equipment, e.g. remove any material from disinfected, polished by detergent (e.g. deionized) by detergent clean rinse.
2. Immerse the item in sufficient volume of Cold Sterilant and/or where passages require item in contact solution for 20°C temperature.
3. Remove items with purified water acceptable level of Actril™ Residue.
4. The solution is stable up to 30 days minutes immersion.
5. For tuberculous with unused solution for 10 minute thoroughly.

NOT REVIEWED
In accordance with FR Notice 62-2.
Based on Draft Labeling Dated

1. Remove any obvious debris or organic material from the surface to be disinfected. This can often be accomplished by rinsing with purified water or by detergent cleaning followed by a water rinse.
2. Immerse the item to be disinfected in sufficient volume of freshly diluted 35X Actril™ Cold Sterilant using purified water, to cover the item and fill passages requiring disinfection. Hold in the disinfecting solution for a minimum of 30 minutes at 20°C temperature (68°F).
3. Remove items after 30 minutes and rinse with purified water until rinse water shows levels of 10 ppm or less when tested with Actril™ Residual Test Strips.
4. The solution cannot be reused as a hospital disinfectant.

DIRECTIONS FOR USE AS A CLEANER-SANITIZER

1. Using water or mechanical action to remove heavy soil or gross filth from hard surfaces such as formica, stainless steel or vinyl surfaces.
2. Apply a freshly made 50X dilution of Actril™, made using purified water, to the surface or immerse items to be sanitized in the solution. Allow 10 minutes of contact time. Let air dry or rinse with purified water; drain excess if possible and allow to air dry.
3. Product may be applied by cloth, mop, or sponge.
4. May not be reused as a cleaner-sanitizer.

DIRECTIONS FOR USE AS A GERMICIDAL SPRAY DISINFECTANT

1. Spray undiluted Actril™ onto cleaned surface using a plastic spray bottle.
2. Allow to remain on surface for 30 seconds at 20°C temperature (68°F).
3. Let air dry or rinse with purified water; drain off excess water if possible and allow to dry.

Storage

1. Store in shipping carton.
2. Do not expose to direct sunlight.
3. Maintain temperature below 75°F (24°C).
4. Avoid contact with combustible materials.
5. Store in original closed container.
6. For chemical emergency, spill, leak, fire, exposure and accident, call Chemtec, day or night (800) 424-9300, (202) 483-7616.

Pesticide Disposal

Do not contaminate water, food or feed by storage or disposal.

Wastes resulting from the use of this product may be disposed of on-site in a sanitary sewer or at an approved waste disposal facility.

Container Disposal

Triple rinse empty container with water. Then offer for recycling or reconditioning; or puncture and dispose of in a sanitary landfill; incinerate; or if allowed by state and local authorities dispose of by burning. If burned, stay clear of the smoke.

Primarily Intended for sterilization, disinfection or sanitization of:

Stainless Steel
Plastic Items and Surfaces
Anesthesia Equipment
Medical Equipment
Surgical Instruments
Respiratory Equipment
Hard Surfaces
Surgical Equipment
Floors
Walls
Bathroom Fixtures
Glass
Dental Instruments
Endotracheal Tubes
Formica
Vinyl



MINNTECH

Minntech Corporation
14905 28th Avenue North
Minneapolis, Minnesota 55447 U.S.A.
Telephone (612) 553-3300 • (800) 328-3346
Telex 29-0625 ASAMPLS PLOH
FAX 612-553-3387

renal systems

Division of Minntech Corporation
14905 28th Avenue North
Minneapolis, Minnesota 55447 U.S.A.
Telephone (612) 553-3300 • (800) 328-3346
Telex 29-0625 ASAMPLS PLOH
FAX 612-553-3387

NOT REVIEWED
In Accordance with FDA Regulation 82-2.

20