

US ENVIRONMENTAL PROTECTION AGENCY
OFFICE OF PESTICIDES PROGRAMS
REGISTRATION DIVISION (TS-767)
WASHINGTON, DC 20460

EPA REGISTRATION NO.
55195-4
TERM OF ISSUANCE

DATE OF ISSUANCE

NOTICE OF PESTICIDE: REGISTRATION
 REREGISTRATION

(Under the Federal Insecticide, Fungicide
and Rodenticide Act, as amended)

NAME OF PESTICIDE PRODUCT

ColdCide-0.25%
Disinfecting Wipes

NAME AND ADDRESS OF REGISTRANT (Include ZIP code)

Metrex Research Corporation
10270 S. Progress Way
Parker, CO 80134

NOTE: Changes in labeling formula differing in substance from that accepted in connection with this registration must be submitted to and accepted by the Registration Division prior to use of the label in commerce. In any correspondence on this product always refer to the above U.S. EPA registration number.

On the basis of information furnished by the registrant, the above named pesticide is hereby Registered/Reregistered under the Federal Insecticide, Fungicide, and Rodenticide Act.

A copy of the labeling accepted in connection with this Registration/Reregistration is returned herewith.

Registration is in no way to be construed as an indorsement or approval of this product by this Agency. In order to protect health and the environment, the Administrator, on his motion, may at any time suspend or cancel the registration of a pesticide in accordance with the Act. The acceptance of any name in connection with the registration of a product under this Act is not to be construed as giving the registrant a right to exclusive use of the name or to its use if it has been covered by others.

This product is conditionally registered in accordance with FIFRA sec. 3(c)(7)(A) provided that you:

1. Submit/cite all data required for registration/reregistration of your product under FIFRA section 3(c)(5) when the Agency requires all registrants of similar products to submit such data; and submit acceptable responses required for reregistration of your product under FIFRA section 4.

2. Make the labeling changes listed below before you release the product for shipment:

a. Add the phrase "EPA Registration No. 55195-4."

Submit five (5) copies of your final printed labeling before you release the product for shipment. Refer to the A-79 enclosure for a further description of final printed labeling.


If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA sec. 6(e). Your release for shipment of the product constitutes acceptance of these conditions.

ATTACHMENT IS APPLICABLE

SIGNATURE OF APPROVING OFFICIAL

A stamped copy of the label is enclosed for your records.

Sincerely,



Walter Francis
Acting Product Manager (31)
Antimicrobial Program Branch
Registration Division (7505)

Enclosures

)

)

BEST AVAILABLE COPY

New Light Mint Scent

ColdCide

0.25% DISINFECTING WIPES

A stable, low pH formula for use in Hospitals, Dental Offices and Veterinary Facilities

- Bactericidal
- Virucidal*
- Pseudomonacidal
- Fungicidal
- Tuberculocidal

Active Ingredients:
 Acetic Acid
 Glutaraldehyde
 O-Phenylphenol
 P-Tertiary Amine
 Inert Ingredients
 TOTAL
 ** Do not include

Contents: 40 Wipes Saturated with Disinfectant

CAUTION
KEEP OUT OF REACH OF CHILDREN
See side panel for additional precautionary statements

This product is not to be used as a terminal sterilant/high level disinfectant on (1) instruments that are introduced directly into the human body, either into or in contact with the blood, the body, or (2) contacts intact mucous membranes but which does not ordinarily otherwise enter normally sterile areas of the body. This product may be used to disinfect or semi-critical medical devices prior to sterilization or high level disinfection.

ACCEPTED WITH COMMENTS
in EPA Form 7260-10
Date: Jan 30 1995

Underwritten by...
Fungicide and...
See... for the...
registered under EPA Reg. No. 55195-4

New Light Mint Scent

ColdCide®

0.25% DISINFECTING WIPES

A disinfectant and deodorant for hospitals, medical and dental offices, and veterinary facilities.
PREPARATION OF DISPENSER:
 Open top of tub, remove and discard foil covering. Remove wipes one at a time as needed. Close dispenser tightly after use to retain moisture.

Laser Graphics — 629-7725 — 05/18/93
 Laser Job No. 3050003G — Customer P.O. 00000



Laser Graphics — 629-7725 — 05/19/93
 Laser Job No. 3050003G — Customer P.O. 00000
 Metrex Research Corp. — ColdCide 0.25% Wipes
 OP: KMS DISK: 77 FONT(S): 6,128,131 FILE(S): 316

BEST AVAILABLE COPY

PRECAUTIONARY STATEMENTS:

Hazard to Humans and Domestic Animals.

CAUTION: MAY BE IRRITATING TO EYES. Do not get in eyes, on skin or on clothing. Wear safety glasses and rubber gloves when handling. Harmful if swallowed. Avoid contamination of food. Avoid reuse of unwashed, disinfectant-soaked gloves. Avoid prolonged inhalation.


STATEMENT OF PRACTICAL TREATMENT:

If in eyes, immediately flush with water for at least 15 minutes. Get medical attention. If on skin, wash thoroughly with water. Wash contaminated clothing before reuse.

WASH THOROUGHLY AFTER HANDLING WIPES.

STORAGE AND DISPOSAL:

Store in cool area. Do not reuse wipes
 Dispose of used wipe and empty dispenser in trash.



Laser Graphics — 629-7725 — 05/18/93
 Laser Job No. 3050003G — Customer P.O. 00000
 Metrex Research Corp. — ColdCide 0.25% Wipes
 OP: KMS DISK: 77 FONT(S): 6,128,131 FILE(S): 314

THE ACTIVITY OF COLDCIDE-0.25% DISINFECTING WIPES AGAINST THE HIV-1 (AIDS VIRUS)

ColdCide-0.25% Disinfecting Wipes completely inactivate the HIV-1 (AIDS Virus) on hard, non-porous surfaces when used at a contact time of 30 seconds at room temperature (20-25°C). Although efficacy at a 30-second contact time has been shown to be adequate against HIV-1, this time would not be sufficient for other organisms. Use a 20-minute contact time for disinfection against all of the organisms claimed on the label.

ColdCide-0.25% Disinfecting Wipes kill HIV-1 on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids, in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus type-1 (HIV-1) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

Personal Protection: Specific barrier protection items to be used when handling items soiled with blood or body fluids are disposable latex gloves, gowns, masks or eye coverings.

Cleaning Procedure: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of ColdCide-0.25% Disinfecting Wipe.

Contact Time: While the HIV-1 virus is inactivated in 30 seconds, use the recommended contact time for disinfection of other organisms on this label.

Disposal of Infectious Materials: Materials that may contain blood and other body fluids should be autoclaved and disposed of according to federal, state and local regulations for infectious waste disposal.

ASSOCIATED
 WITH COMMENTS
 In EPA Letter Dated:

July 1995

Under the provisions of the title,
 Sample 1 and 2 are being tested
 for HIV-1 and HIV-2. Sample
 1 and 2 are EPA Reg. No.

53195-4

New Light Mint Scent



A disinfectant in
 facilities.
PREPARATION C
 Open top of tub,
 as needed. Close