



HOME USE: For disinfection of Respiratory Therapy Equipment in the home refer to technical Bulletin "CONTROL III DISINFECTANT FOR HOME CARE".

NOTE: CONTROL III DISINFECTANT, GERMICIDE has not been tested for effectiveness against Mycobacterium Tuberculosis and must not be relied upon when a Tuberculocidal product is desired.

USED FOR:

Disinfection of walls, equipment, heat sensitive instruments and tubing in operating rooms, recovery rooms and intensive care rooms.

**PRECAUTIONARY STATEMENTS:
HAZARDS TO HUMANS AND
DOMESTIC ANIMALS**

DANGER

Keep out of reach of children. Corrosive. Causes severe eye and skin damage. Do not get in eyes, on skin, or on clothing. Wear goggles or face shield and rubber gloves when handling. Harmful if swallowed. Avoid contamination of food. Remove and wash contaminated clothing.

STORAGE AND DISPOSAL

- Do not contaminate water, food or feed by storage or disposal
- Do not store on side
- Do not reuse empty container

PESTICIDE DISPOSAL

Pesticide wastes are acutely hazardous. Improper disposal of pesticide spray or mixture of rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Agency, or the Hazardous Waste representative at the nearest EPA Regional Office of guidance.

CONTAINER DISPOSAL

Do not reuse container. Rinse thoroughly before discarding in trash. Triple rinse (or equivalent). Then offer for recycling or reconditioning, or puncture and dispose of in a sanitary landfill, or incineration, or if allowed by State and local authorities by burning, if burned, stay out of smoke.

Manufactured by Maril Products, Inc. (714) 844-7718

**CONTROL III
disinfectant
GERMICIDE**

Rapid Action — Long Life

ACTIVE INGREDIENTS

- n-alkyl (80% C₁₄, 20% C₁₈, 5% C₈, 5% C₁₂) dimethyl benzyl ammonium chloride _____ 10%
- n-alkyl (88% C₁₂, 32% C₁₄) dimethyl ethyl benzyl ammonium chloride _____ 10%

INERT INGREDIENTS _____ 80%
100%

**ON INANIMATE ENVIRONMENTAL
SURFACES EFFECTIVE AGAINST**

Pseudomonas aeruginosa • Adenovirus type 1 • Staphylococcus aureus • Vaccinia virus • Salmonella choleraesuis • Influenza A₂ virus • Escherichia coli • Herpes simplex virus

**KEEP OUT OF REACH OF CHILDREN
DANGER**

STATEMENT OF PRACTICAL TREATMENT

In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash contaminated clothing before reuse.

If swallowed, drink one or two large glasses of milk, egg whites, gelatin solution, or—if these are not available—drink one or two large glasses of water. Avoid alcohol. Call a physician immediately.

NOTE TO PHYSICIAN: Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiration depression, and convulsion may be needed.

This product is a concentrate and must be diluted before use.

See left panel for additional precautionary statements.

NET 16 FL. OZ. (1 PT.) EPA Est. #40873 CA-01
EPA Reg. #55384-3



MARIL PRODUCTS, INC.

820 South B Street, Suite A, Torrance, California 92660 U.S.A.

DIRECTIONS FOR USE:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

General:

1. Remove heavy soil as required.
2. Apply CONTROL III DISINFECTANT as directed.
3. Will not harm surfaces which are safe in water.

OPERATING ROOM GERMICIDE:

Mix 1/2 ounce CONTROL III DISINFECTANT in one gallon of water. Resulting mixture contains 780 ppm active ingredients. Apply to surfaces of furniture and equipment in operating room. Allow to air dry.

COLD DISINFECTANT:

Clean instruments with CONTROL III Detergent. Disinfect as follows: Mix one ounce CONTROL III Disinfectant in one gallon of water. Immerse instruments for 10 minutes or more. Active concentration is 1560 ppm.

INSTRUMENT STORAGE:

Mix one ounce CONTROL III DISINFECTANT in one gallon of water. Add 6 gms. sodium carbonate per gallon of solution. Solution can be retained for one month.

**FOOD PROCESSING EQUIPMENT
AND DISHES**

Remove all gross food soils with a pre-rinse or by scraping. Then wash all surfaces with a good liquid detergent such as CONTROL III PRE-CLEANER. After washing, rinse all surfaces thoroughly with potable water. To sanitize apply a solution of 1 ounce CONTROL III DISINFECTANT in 8 gallons of water (200 ppm active). Surfaces should remain wet for one minute followed by adequate draining and air drying. Fresh solution should be prepared daily or when use solution becomes visibly dirty. For mechanical application, use solution may not be reused for sanitizing applications.

Apply to sink tops, counter tops and other stationary surfaces by cloth or brush. Dishes, silverware, glasses and other similar size food processing equipment can be sanitized by immersion in a 1/2 ounce per 4 gallons of water dilution of CONTROL III DISINFECTANT GERMICIDE. No potable rinse is required.

LAUNDRY:



Use 1/2 ounce CONTROL III per 25 lbs. dry weight of laundry in the final rinse.

REVISED 6/22/88

- . HALF THE COST OF VINEGAR
- . GOOD FOR TWO WEEKS
- . TIME TESTED SINCE 1974

EPA #55364-3

MARIL PRODUCTS, INC.
620 South "B" Street, Ste. A
TUSTIN, CA 92680



WHAT IS CONTROL III?



CONTROL III is a time proven germicide. It is an excellent disinfectant for respiratory therapy equipment in the home. CONTROL III has been used by Visiting Nurse Associations, Home Care Organizations and Hospital Home Rehabilitation programs since 1974.

CONTROL III, when used in an overall cleaning program, will help control infection or reinfection due to contaminated nebulizers, tubing or reservoirs used for respiratory therapy. 1.

HOW TO PREPARE CONTROL III.

Prepare a disinfecting solution of CONTROL III DISINFECTANT GERMICIDE BY ADDING 1 ounce (one tablespoon) CONTROL III to 2 quarts of water in a good grade plastic container. This solution may be reused for up to 14 days. Check solution daily with test strips to ensure adequate disinfection level of CONTROL III. If the solution becomes visibly dirty or falls below 1,000 ppm dump the used solution and make a fresh solution of CONTROL III. Keep disinfecting solution covered when not in use. 2.

1. CONTROL III DISINFECTANT GERMICIDE has not been tested for effectiveness against Mycobacterium tuberculosis and must not be relied upon when a tuberculocidal product is desired.
2. At this use dilution CONTROL III is neither a skin irritant nor sensitization agent.




DISINFECTION OF RESPIRATORY EQUIPMENT WITH CONTROL III

1. Disassemble equipment according to manufacturer's recommendations.
2. Wash nebulizer, flex tubing, mouthpiece and manifold valve thoroughly in a mild detergent (such as CONTROL III PRE-CLEANER) and warm tap water. DO NOT USE SOAP. Rinse parts with tap water.
3. Immerse equipment parts in solution of CONTROL III DISINFECTANT GERMICIDE and soak for 10 minutes.
4. Rinse all equipment parts with hot tap water (allow water to run for 3 minutes before rinsing parts) and allow to air dry.
5. After parts are dry cover them with a clean towel until equipment is reassembled.

ACCEPTED

OCT 17 1983

Under the provisions of the
Federal Food, Drug, and Cosmetic Act, as amended,
periodically tested under
EPA Reg No.



PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

DANGER

Keep out of reach of children. Corrosive. Causes severe eye and skin damage. Do not get in eyes, on skin or on clothing. Wear goggles or face shield and rubber gloves when handling. Harmful if swallowed. Avoid contamination of food. Remove and wash contaminated clothing.

STORAGE AND DISPOSAL

Storage: Store in original container placed in locked storage area in a dry place no lower than 50°F or higher than 120°F.

Disposal: Do not reuse container. Securely wrap original container in several layers of newspaper and discard in trash.

ACCEPTED
OCT 17 1998
This material has been inspected and found to be in accordance with the requirements of the...
LEAFING

