

UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

JUL 20 1989

Wet Products, Inc.
620 South F Street, Suite A
Castroville, CA 92680

Attention: James C. Swims

Gentlemen:

Subject: Control III Laboratory Germicide
EPA Registration No. 55364-4
Your Application Dated April 17, 1989

The labeling referred to above, submitted in connection with registration under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, is acceptable. A stamped copy is enclosed for your records.

However, at the next printing of your product labeling:

1. Include the statement "Hazards to Humans and Domestic Animals" under the "Precautionary Statements" heading.
2. Revise the "First Aid" heading to read "Statement of Practical Treatment."
3. Revise the referral statement to read "See left panel for additional precautionary statements."

If you have any questions concerning this letter, please contact Martha DeLaney at (703) 557-6982.

Sincerely yours,



John H. Lee
Product Manager (31)
Antimicrobial Program Branch
Registration Division (H7505C)

Enclosure

54771:A-4:I:DeLaney:L-3:KENCO:7/12/89:8/21/89:KA:GJ:VO:DE:KA

CONCURRENCE

SYMBOL

GERMICIDE

RAPID ACTION — READY TO USE

Bactericidal • Fungicidal • * Virucidal

Requires No Activators • Contents Are Ready For Use

ACTIVE INGREDIENTS — n-alkyl (60% C₁₂, 30% C₁₄, 5% C₁₆, 5% C₁₂)
dimethyl benzyl ammonium chloride 0.0781%
n-alkyl (68% C₁₂, 32% C₁₄) dimethyl ethyl
benzyl ammonium chloride 0.0781%

INERT INGREDIENTS 99.8438%

100.0000%

KEEP OUT OF REACH OF CHILDREN

WARNING

SEE FIRST AID and other precautions on left panel

One Gallon (128 fl. oz.)

EPA Est. #40873 CA-01 EPA Reg. #55304-4



MARIL PRODUCTS, INC.

620 South B Street, Suite A, Tustin, California 92680 U.S.A.

BEST AVAILABLE COPY

2. Will not harm surfaces.

3. Before disinfecting, clean off all contaminants. Use of CONTROL III Detergent Germicide is recommended. Rinse well before applying Laboratory Disinfectant/Germicide solution.

4. CONTROL III Laboratory Disinfectant Germicide does not affect cements.

5. Solution can remain in this container indefinitely without loss of effectiveness.

Operating Room Germicide:

Apply CONTROL III Laboratory Disinfectant Germicide solution to a clean sponge or cloth. Wipe down surgery table, instrument stands, and other equipment. Allow to dry in air.

Surgical Instruments:

Pre-clean instruments with CONTROL III Detergent Germicide (one ounce per gallon of water). Rinse thoroughly in running water. Place instruments in a sterile auto-clavable plastic tray or other suitable tray. Fill tray with CONTROL III Laboratory Disinfectant/Germicide solution to fully cover all instruments. Soak for ten minutes or longer for disinfection.

Surgical Instrument Storage:

Place instruments in clean storage tray container. Fill container with CONTROL III Laboratory Disinfectant/Germicide solution to fully immerse instruments. Cover container. In normal use storage solution can be retained for thirty days. Replace solution more often if frequent instrument removals and additions are made.

BEST AVAILABLE COPY

Effective Against:

Pseudomonas aeruginosa (PRD-10); *Staphylococcus aureus* (ATCC 6538); *Salmonella choleraesuis* (ATCC 10708); Influenza A₂ virus; Herpes simplex; Adenovirus type #5; Vaccinia viruses; Trichophyton interdigitale, on inanimate hard surfaces.

Note: Bacteria → killed within ten minutes.

PRECAUTIONARY STATEMENTS

WARNING:

Keep Out of Reach of Children. Avoid getting in eyes. Avoid contamination of food. May be harmful if swallowed. Rinse empty container thoroughly with water and discard it.

FIRST AID:

Immediately flush eyes with plenty of water. If irritation persists, obtain medical attention. If swallowed, drink large quantities of water.

STORAGE AND DISPOSAL:

Do not contaminate water, food or feed by storage or disposal.

CONTAINER DISPOSAL:

Do not reuse container (bottle, can, jar). Rinse thoroughly before discarding in trash.

STORAGE:

Store in a dry place no lower in temperature than 50°F or higher than 120°F.

Manufacturer of:
MARIL PRODUCTS, INC.
(714) 544-719

HIV-1: See Technical Bulletin "Efficacy of CONTROL III Products For the Control Of The Human Immunodeficiency Virus Type 1 (AIDS VIRUS)"

BEST AVAILABLE COPY

ACCEPTED
COMMENTS:
A Letter Dated:

20 1989

Federal Insecticide,
Fungicide, and Rodenticide Act
amended, for the pesticide
registered under EPA Reg.

55364--



MARIL PRODUCTS, INC.

Efficacy of Control III Products for the control of the Human Immunodeficiency Virus Type 1 (AIDS Virus).

Both Control III Disinfectant Germicide and Control III Laboratory Germicide kill the HIV-1 (AIDS Virus) on hard, non-porous surfaces in the presence of a moderate amount of organic soil (10% bovine serum) when diluted in hard water (400 ppm CA CO₃) at a level of 781 ppm active quaternary for a contact time of 10 minutes at room temperature.

Preparation of 781 ppm active quaternary: Control III Disinfectant Germicide diluted at 1/2 ounce per gallon of water equals 781 ppm active quaternary. Control III Laboratory Germicide diluted one part Laboratory Germicide in one part water equals 781 ppm active quaternary.

Kills Human Immunodeficiency Virus Type 1 (HIV-1) (AIDS Virus) on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids.

Special instructions for cleaning and decontamination against HIV-1 of surfaces/objects soiled with blood/body fluids:

PERSONAL PROTECTION: When handling items soiled with blood or body fluids use disposable latex gloves, gowns, masks or eye coverings.

CLEANING PROCEDURES: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of Control III Disinfectant Germicide, and Control III Laboratory Germicide.

CONTACT TIME: Allow Surface to remain wet for 10 minutes.

DISPOSAL OF INFECTIOUS MATERIALS: Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

BEST AVAILABLE COPY