**PROTECTION AGENCY** INITED STA

APR 8 2010

Mr. David Swain Agent for Maril Products, Inc. 320 West 6<sup>th</sup> Street Tustin, CA 92780

Subject: **Control III Laboratory Germicide** EPA Registration Number 55364-4 Notification Application Dated March 10, 2010 EPA Received Date March 12, 2010

Dear Mr. Swain:

This will knowledge receipt of your notification, submitted under the provisions of FIFRA section 3c 9. Based on a review of the submitted material the following comment apply.

#### **Proposed Amendment:**

Revise the Storage and Disposal Per PR Notice 2007-4

#### General Comments:

Based on a review of the material submitted, the following comment apply:

The notification is acceptable. A copy has been inserted in your file for future reference.

Should you have any questions or comments concerning this letter, please contact Drusilla Copeland at (703) 308-6224.

Sincerely,

Ducilla Copelind

Velma Noble

Product Manager (31) **Regulatory Management Branch I** Antimicrobials Division (7510P)

#### CONCURRENCES

SYMBOL SURNAME DATE OFFICIAL FILE COPY Printed on Recycled Paper EPA Form 1320-1A (1/90)

3)								2011	
	Washington, DC 20460			Form Approve		Registration Amendment		D. Approval expires 2-28-9 OPP Identifier Number	
		Application	for Pesticide	- Section	on I				
1. Company/Product Number 55364-4			2. EPA Product Manager 3. Velma Noble				. Proposed Classification		
4. Company/Product (Name) Control III Laboratory Germicide			<b>PM#</b> 31	PM#				Restricte	
5. Name and Address of A Maril Products, Inc. 320 West 6th Street Tustin, CA 92780	pplicant <i>(Include ZIP Cc</i>	ode)	(b)(i), my p to: EPA Reg.	roduct is	w. In accordanc similar or identica	al in cor	mposition an	d labeling	
Chack II II			Product I	Name					
			Section - II						
Amendment - Explain Resubmission in rev Notification - Explain Explanation: Use additi		Ag "M Oti	Final printed labels in repsonse to Agency letter dated "Me Too" Application. Other - Explain below.						
Notification of updated stor regulations at 40 CFR 152 a violation of 18 U.S.C. Se Notice 98-10 and 40 CFR of FIFRA. Please contact D	.46, and no other change c, 1001 to willfully make a 152.46, this product may	s have been made any false statement be in violation of FI and provide an ele	to the labeling or con to EPA. I further und IFRA and I may be su ectronic response to o	ifidential sta lerstand that ubject to en	tement of formula of t if this notification is orcement action and	f this pro- not con: penaltie	duct. I understa sistent with the es under sectio	and that it i terms of F	
			Section - III						
1. Material This Product Will Be Packaged In:         Child-Resistant Packaging         Yes         ✓         No         * Certification must		No. per I	Water Soluble Packaging Yes V No If "Yes" No. per Package wgt containe		2. Type of Container Metal Plastic Glass Paper Other (Specify)				
be submitted		1							
3. Location of Net Contents Information       4. Size(s) I         I       I         Container       I			Container 1 gal	5	5. Location of Label Directions				
6. Manner in Which Label		Lithograph Paper glu Stenciled	h od	Other	Leasan				
			Section - IV			.n.e.e	0		
1. Contact Point (Complet	te items directly below i	for identification o	f individual to be co	ntected, if					
<b>Name</b> David Swain		Tit Aç	<b>ie</b> gent for Maril Produ	ucts, inc.		260) 244	<b>e No. (Include</b> -62၇၀ ၀ ၂ ၂ ၀	Area Cod	
	tements   have made or any knowlinglly false or le law.		attachments therete		accurate and comp		6. Date Appli Beceived		
2. Signature			3. Title Agent for Maril Products, Inc.						
- 1/m h									
4. Typed Name		5. [	Date / /						
David Swain			3/10/10	)					
PA Form 8570-1 (Rev. 3-9	94) Previous editions are	e obsolete.		White	- EPA File Copy (	original)	Yellow - /	Applicant (	

A Form 8570-1 (Rev. 3-94) Previous editions are obsolete. Er/



March 10, 2010

Ms. Velma Noble, PM 31 Document Processing Desk (NOTIF) Office of Pesticide Programs (7504P) U. S. Environmental Protection Agency Room S-4900, 4th Floor One Potomac Yard 2777 S Crystal Drive Arlington, VA 22202

SUBJECT:

Notification (Container and Containment Regulations) per PR Notice 2007-4 Control III Laboratory Germicide EPA Reg. No.: 55364-4

#### Dear Ms. Noble,

On behalf of our client, Maril Products, Inc., enclosed is a notification revising the storage and disposal section to comply with the Final Pesticide Container and Containment Regulations. This is a Notification per PR Notice 2007-4; therefore, no PRIA fee is required.

The new statements, as specified in the guidance, have been included as required text in the attached labeling. The new statements are highlighted in yellow on the enclosed label. No other changes in the existing labeling have been made.

This notification, as stipulated by the guidance, is consistent with the provisions of PR Notice 98-10 and EPA regulations at 40 CFR 152.46, and no other changes have been made to the labeling or confidential statement of formula of this product. I understand that it is a violation of 18 U.S.C. Sec, 1001 to willfully make any false statement to EPA. I further understand that if this notification is not consistent with the terms of PR Notice 98-10 and 40 CFR 152.46, this product may be in violation of FIFRA and I may be subject to enforcement action, and penalties under sections 12 and 14 of FIFRA.

Please contact me at (260) 244-6270 or at dswain@srcconsultants.com if you have any questions about o vovos

Sincerely,

David Swain Agent for Maril Products, Inc. cc: A. Whitehead

0000

0000

0000

Phone: 260-244-6270 Fax: 260-244-6273

P.O. Box 1014 Columbia City, IN 46725

www.srcconsultants.com

#### BACTERICIDAL **CONTROLI FUNGICIDAL \*VIRUCIDAL REQUIRES NO ACTIVATORS** CONTENTS ARE READY FOR USE ACTIVE INGREDIENTS n-alkvl (60% C14, 30% C16, 5% C18, 5% C12) dimethyl benzyl ammonium chloride LABORATORY GERMICIDE ......**0.0781%** n-alkyl (68% C12, 32% C14) dimethyl ethyl benzyl ammonium chloride **KEEP OUT OF REACH OF CHILDREN** CAUTION SEE LEFT PANEL FOR ADDITIONAL PRECAUTIONARY STATEMENTS ONE GALLON (128 FL. OZ.) EPA EST. #40873 CA-01/GA-01 EPA REG. # 55364-4 READY TO USE " MARIL PRODUCTS, INC. 320 West 6th Street

HOPU

Tustin, CA 92780 USA 800-546-7711 • 714-544-7711

# DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

# **DISINFECTION:**

Control III Laboratory Germicide is a "Ready-To-Use" disinfectant which is designed for use as a disinfectant in hospitals, doctors' offices, dentists' offices, nursing homes, clinics and surgi-centers. It is formulated to disinfect hard, non-porous, inanimate surfaces such as counter tops, beds, instruments, fixtures and equipment.

1. Prior to disinfection with Control III Laboratory Germicide, all surfaces must be thoroughly cleaned to remove gross filth or heavy soil.

2. Use full strength. Do not dilute.

3. Apply Control III Laboratory Germicide with a clean sponge, cloth or by a trigger sprayer so as to wet all surfaces thoroughly. For disinfection, surfaces must remain wet for ten minutes. Allow to air dry.

4. Control III Laboratory Germicide

is classified as a low-level disinfectant.

# **OPERATING ROOM GERMICIDE:**

Apply CONTROL III Laboratory Germicide solution to a clean sponge or cloth. Wipe down surgery tables, instrument stands, and other equipment so as to wet all surfaces thoroughly. Allow to air dry. (See Disinfection directions)

# SURGICAL INSTRUMENT STOR-AGE:

Place instruments in clean storage tray container. Fill container with CONTROL III Laboratory Germicide solution to fully immerse instruments. Cover container. In normal use storage solution can be retained for thirty days. Replace solution more often if frequent instrument removals and additions are made. Instruments should be autoclaved after soaking in Control III prior to reuse. (See Disinfection directions)

(((C) () (C) () (C) c,c,c,c,c

#### **Effective Against:**

Pseudomonas aeruginosa (PRD-10); Staphylococcus aureus (ATCC 6538); Salmonella choleraesuis (ATCC 10708); Influenza A2 virus\*; Herpes simplex type 1\*; Adenovirus type #5; Vaccinia viruses\*; Trichophyton mentagrophytes, on hard, non-porous, inanimate surfaces.

Note: Control III Laboratory Germicide has not been tested for effectiveness against Mycobacterium Tuberculosis and must not be relied upon when a Tuberculocidal product is desired.

HIV-1: See Technical Bulletin "Efficacy of Control III Products For the Control of the Human Immunodeficiency Virus Type 1 (AIDS VIRUS).

#### KILLS HIV-1 ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY

**FLUIDS** in health care settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS).

#### SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINA-TION AGAMST MIV- OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS. PERSONAL PROTECTION:

Clean up should always be done wearing protective crex ploves, masks, and eye protection.

# **CLEANING PROCEDURES:**

Blood and other body fluids

containing HIV must be thoroughly cleaned from surfaces and objects before application of this product.

#### CONTACT TIME:

Allow surface to remain wet for 10 minutes.

# DISPOSAL OF INFECTIOUS MATERIALS:

Bloody body fluids, cleaning materials and clothing should be autoclaved and disposed of according to local regulations for infectious waste disposal.

#### PRECAUTIONARY STATEMENTS: HAZARDS TO HUMANS AND DOMESTIC ANIMALS CAUTION:

Causes moderate eye irritation. Avoid contact with eyes, skin, or clothing. Wash thoroughly with soap and water after handling.

#### STORAGE AND DISPOSAL:

Do not contaminate water, food, or feed by storage and disposal.

#### **PESTICIDE STORAGE:**

Store in a dry place no lower in temperature than 50°F or higher than 120°F

# PESTICIDE DISPOSAL:

Wastes resulting from the use of this product must be disposed of on site or at an approved waste disposal facility.

# **CONTAINER HANDLING:**

Nonrefillable container. Do not reuse or refill this container. Offer for recycling, if available or wrap container and put in trash.

#### FIRST AID

If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first five minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.

**If on skin or clothing:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.

**If swallowed:** Call a poison control center or doctor immediately for treatment advice. Have a person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by the poison control center or doctor. Do not give anything by mouth to an unconscious person.

**If inhaled:** Move person to fresh air. If person is not breathing, call 911 or an ambulance then give artificial respiration, preferably by mouth-tomouth, if possible. Call a poison control center or doctor for further treatment advice.

Note to Physician: Probable mucosal damage may contraindicate the use of gastric lavage.

# HOT LINE NUMBER

For emergency medical advice, call your local poison control center (1-800-222-1222) or doctor. Have the product container or label with you when seeking medical advice or treatment.

P/N 10005 Rev. A