



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

1 of 6

PM 31

55364-7

11/23/98

 U.S. ENVIRONMENTAL PROTECTION AGENCY Office of Pesticide Programs Antimicrobial Division (7505W) 401 "M" St., S.W. Washington, D.C. 20460	EPA Reg. Number: 55364-7	Date of Issuance: NOV 23 1998
	Term of Issuance: Conditional	
	Name of Pesticide Product: CONTROL III ELITE	
NOTICE OF PESTICIDE: <input checked="" type="checkbox"/> Registration <input type="checkbox"/> Reregistration (under FIFRA, as amended)		
Name and Address of Registrant (include ZIP Code): MARIL PRODUCTS, INC 620 S. B STREET, STE A. TUSTIN, CA 92780-4318		
<p><b>Note:</b> Changes in labeling differing in substance from that accepted in connection with this registration must be submitted to and accepted by the Registration Division prior to use of the label in commerce. In any correspondence on this product always refer to the above EPA registration number.</p> <p>On the basis of information furnished by the registrant, the above named pesticide is hereby registered/reregistered under the Federal Insecticide, Fungicide and Rodenticide Act.</p> <p>Registration is in no way to be construed as an endorsement or recommendation of this product by the Agency. In order to protect health and the environment, the Administrator, on his motion, may at any time suspend or cancel the registration of a pesticide in accordance with the Act. The acceptance of any name in connection with the registration of a product under this Act is not to be construed as giving the registrant a right to exclusive use of the name or to its use if it has been covered by others.</p> <p>This product is conditionally registered in accordance with FIFRA sec. 3(c)(7)(A) provided that you:</p> <ol style="list-style-type: none"> <li>Submit and/or cite all data required for registration/reregistration of your product under FIFRA sec. 3(c)(5) when the Agency requires all registrants of similar products to submit such data; and submit acceptable responses required for reregistration of your product under FIFRA section 4.</li> <li>Make the following label changes:                     <ol style="list-style-type: none"> <li>Revise the EPA Registration Number to read, "EPA Reg. No. 5534-7.</li> <li></li> </ol> </li> </ol>		
Signature of Approving Official:  PM-31/AD/OPP	Date: NOV 23 1998	

EPA Form 8570-6

CONCURRENCES

SYMBOL							
SURNAME							
DATE							

page 2  
EPA Reg. No. 55364-7

Submit five final printed copies of labeling before you release the product for shipment.

If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA sec. 6(e). Your release for shipment of the product constitutes acceptance of these conditions.

A stamped copy of the label is enclosed for your records.

Sincerely,



Velma Noble  
Product Manager 31  
Regulatory Management Branch I  
Antimicrobial Division



**LEFT PANEL**

**KILLS HIV-1 ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings (Hospitals, Nursing Homes) or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS).

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.**

**PERSONAL PROTECTION:** When handling items soiled with blood or body fluids use disposable latex gloves, gowns, masks, or eye coverings.

**CLEANING PROCEDURES:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of this product.

**CONTACT TIME:** Allow surface to remain wet for 10 minutes.

**DISPOSAL OF INFECTIOUS MATERIALS:** Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

**PRECAUTIONARY STATEMENTS  
HAZARDS TO HUMANS AND DOMESTIC ANIMALS**

**WARNING**

**KEEP OUT OF REACH OF CHILDREN.** Causes eye and skin irritation. Do not get in eyes, on skin or on clothing. Harmful if swallowed. Avoid contamination of food. Remove contaminated clothing and wash before reuse. Wash thoroughly with soap and water after handling.

**STATEMENT OF PRACTICAL TREATMENT**

In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician.

If swallowed, drink egg whites or gelatin solution, or if these are not available, drink large quantities of water. Call a physician immediately.

**NOTE TO PHYSICIAN**

Probable mucosal damage may contraindicate the use of gastric lavage.

**STORAGE AND DISPOSAL**

**DO NOT CONTAMINATE WATER, FOOD, OR FEED BY STORAGE OR DISPOSAL**

ACCEPTED  
with COMMENTS  
in EPA Letter Dated  
NOV 23 1998  
55364-17

**LEFT PANEL**

**CONTAINER DISPOSAL** - Do not reuse container (bottle, can, or jar). Rinse thoroughly before discarding in trash.

**STORAGE** - Store in a dry place no lower in temperature than 50 F or higher than 120 F.

ACCEPTED  
with COMMENTS  
in EPA Letter Dated

NOV 23 1998

55364-7

**RIGHT PANEL**

**DIRECTIONS FOR USE**

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

Control III ELITE is specifically designed for the cleaning and disinfection of non-critical medical devices such as bed frames, wheel chairs, physical therapy equipment and other hard non-porous environmental surfaces in hospitals, nursing homes, durable medical equipment suppliers, doctors offices, schools, gymnasiums, homes and animal clinics. Control III ELITE may also be used to preclean or decontaminate instruments in hospitals or other medical facilities prior to sterilization or high level disinfection.

**DISINFECTION AND CLEANING:**

Remove gross filth or heavy soil prior to application of Control III ELITE. Apply solution as is - **DO NOT DILUTE** - to equipment with a cloth, sponge or by spraying so as to thoroughly wet all surfaces. For disinfection surfaces must remain wet for ten minutes. Allow to air dry or wipe off with a clean cloth or sponge.

**PRECLEANING AND SOAKING:**

Control III ELITE can be used for precleaning and soaking small medical devices prior to sterilization or high level disinfection. All items to be soaked must be precleaned to remove gross filth or heavy soil. Items should then be fully immersed in Control III ELITE solution for ten minutes; then thoroughly rinsed and allowed to air dry.

**TUBERCULOCIDAL ACTIVITY:**

This product exhibits disinfectant efficacy against Mycobacterium tuberculosis (BCG) in 10 minutes at 20 degrees Centigrade when used as directed on previously cleaned hard nonporous inanimate surfaces.

**BACTERICIDAL ACTIVITY:** When used as directed, this product exhibits effective disinfectant activity against the organisms: Staphylococcus aureus, Salmonella choleraesuis, Pseudomonas aeruginosa, Escherichia coli 0157:H7, and meets the requirements for hospital use.

**\*VIRUCIDAL ACTIVITY:** This product, when used on environmental, inanimate, non-porous surfaces exhibits effective virucidal activity against HIV-1 (associated with AIDS) and Canine Parvovirus.

ACCEPTED  
with COMMENTS  
In EPA Letter Dated:

NOV 23 1998

Under the Federal Insecticide,  
Fungicide, and Rodenticide Act  
amended, for the pesticide,  
registered under EPA Reg. No.

55364-7