


STORAGE AND DISPOSAL: Store this product in a cool, dry area, away from direct sunlight and heat to avoid deterioration. In case of spill, flood areas with large quantities of water. Product or rinsates that cannot be used should be diluted with water before disposal in a sanitary sewer. Refill empty container only with DISPATCH® Hospital Cleaner Disinfectant with Bleach Refill. When discarding container, place in trash collection. Do not contaminate food or feed by storage, disposal, or cleaning of equipment.

PRECAUTIONARY STATEMENTS: Hazards to Humans and Domestic Animals. **CAUTION:** Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling. **FIRST AID:** If in eyes: hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for further treatment advice. **PHYSICAL AND CHEMICAL HAZARDS:** Do not use this product with ammonia or acids such as vinegar, rust removers, or toilet bowl cleaners. This product not intended for use on fabrics.


EPA Reg. No. 56392-7
EPA Est. No. 62801-MI-1
74592-OH-1
50757-WI-1

©2005 Caltech Industries, Inc.
Midland, MI 48642
1-800-234-7700
®Trademark of
Caltech Industries, Inc.
68832-200-7A



0 = NONE
1 = SLIGHT
2 = MODERATE
3 = HIGH
4 = EXTREME
OX = OXIDIZER

REORDER NO. 68832
See Expiration Date on Bottle



DISPATCH®

**Hospital Cleaner Disinfectant
with Bleach (5500 ppm)**

- Cleans & Disinfects in One Step • Stabilized
- Ready-to-use • Bleach & Detergent Solution

Kills bacteria, viruses, & fungi on hard, non-porous surfaces:

In 30 Seconds:

Mycobacterium tuberculosis (TB)

In 1 Minute:

Staphylococcus aureus, Methicillin resistant *Staphylococcus aureus* (MRSA), *Salmonella choleraesuis*, *Pseudomonas aeruginosa*, *Streptococcus pyogenes*, Vancomycin resistant *Enterococcus faecalis* (VRE), *Clostridium difficile* (vegetative), *Trichophyton mentagrophytes*, HIV-1 (AIDS Virus)¹, Herpes Simplex Virus (HSV-2)², Rotavirus³, Poliovirus, Type 1 (Mahoney)⁴, Hepatitis B Virus (HBV), Hepatitis C Virus (HCV)

Active ingredient: Sodium Hypochlorite 0.55%
Other Ingredients: 99.45%

**KEEP OUT OF REACH OF CHILDREN
CAUTION**

READ COMPLETE DIRECTIONS AND PRECAUTIONS ON REAR PANEL
NET 32 FL. OZ. (1 Quart)

DISPATCH® Hospital Cleaner Disinfectant with Bleach

• For hospital non-porous surfaces such as stainless steel, plastic, ceramic tile, laminated plastic countertops, vinyl, fiberglass, and porcelain, which must be cleaned and disinfected • A rapid-acting cleaner disinfectant that can't scratch • For use in outpatient clinic, laboratory, clinic, emergency room, radiology, operating rooms, physical therapy, physician office, nursing home, dental clinic, hospital central service area, kitchen and food service area (non-food surfaces).

DIRECTIONS FOR USE: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. To Clean: Hold bottle 6-8" from the surface. Direct gentle spray to soiled area or apply with damp applicator. Wipe clean with a damp applicator such as sponge/brush/paper towel. Stubborn stains — allow 5 minute contact. No rinsing needed. If streaking is observed, rinse with water. Polish with dry cloth. To Clean and Disinfect in One Step: Spray surface thoroughly until completely wet. Allow 1 minute at room temperature (68-77° F) to wet surface and disinfect. Wipe clean with a damp sponge, cloth, or paper towel. No rinsing needed. To Kill TB allow 30 seconds contact at room temperature. Gross soil must be removed prior to cleaning and disinfecting.

To Pre-clean Instruments Prior to Terminal Sterilization/high level disinfection. Place instruments in a suitable container. Spray to thoroughly drench all surfaces or pour sufficient undiluted product over instruments to allow for complete submersion. Remove instruments, rinse and follow with an appropriate disinfection process. (Critical and semi-critical devices must be followed by the appropriate terminal sterilization/high level disinfection process).

Special Instructions for Cleaning and Decontamination Against HIV-1 (AIDS Virus), HBV, and HCV on Surfaces/Objects Soiled with Blood/Body Fluids. KILLS HIV, HBV, and HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV). • Personal Protection: When handling items soiled with blood/body fluids, personal protective equipment such as disposable plastic gloves, gowns, masks, or eye coverings should be used. • Cleaning Procedure: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before disinfecting. • Contact Time: Leave surfaces wet for 1 minute. • Disposal of Infectious Materials: Blood, body fluids, and cleaning materials should be autoclaved and disposed according to federal, state and local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

ACCEPTED

OCT 18 2005

Under the Federal Insecticide, Fungicide, and
Rodenticide Act as amended, for the
pesticide, registered under
EPA Reg. No. 56392-7

56392-7

10/18/2005

Page 181