

FILE:VINCLNF.COM COLOR: Composite SCPPFN: Invoice's default TIME: Fri Sep 10 13:56:10 1993

HOSPITAL DISINFECTANT/CLEANER

VIRAHOL®

U.S. PATENT 5,145,663

ACTIVE INGREDIENT - Isopropanol . . .	70%
PROPYLENE GLYCOL	10%
INERT INGREDIENTS	20%
TOTAL	100%

**KEEP OUT OF REACH OF CHILDREN
WARNING**

SEE RIGHT PANEL FOR ADDITIONAL PRECAUTIONARY STATEMENTS

© 1993 VERIDIEN CORPORATION



VERIDIEN CORPORATION
11800 28TH STREET NORTH
ST. PETERSBURG, FL 33716
(813) 572-6636
EPA Reg. No. 60142 - 1
EPA Est. No. 60142 - FL - 1


NET CONTENTS - 1 U.S. GALLON - (128 Fluid Ounces) 3.785 Liters

ACCEPTED
SEP 16 1993
Under the Federal Insecticide, Fungicide, and
Rodenticide Act, as amended, for the
pesticide registered under
EPA Reg. No. 60142-1

BEST AVAILABLE COPY

182

2 7 2



FILE:VCLM788.CDR COLOR: Composite SCREEN: Device's default TIME: Tue Sep 14 17:38:48 1993

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS & DOMESTIC ANIMALS

WARNING

Causes eye irritation. Do not get in eyes. In case of eye contact, immediately flush eyes with plenty of water. Get medical attention if irritation persists. Harmful if swallowed.

PHYSICAL OR CHEMICAL HAZARDS: Flammable, keep away from heat & open flame

STORAGE AND DISPOSAL

STORAGE

Store in a cool, dry place.

CONTAINER DISPOSAL

Triple rinse (or equivalent). Then offer for recycling or reconditioning, or puncture and dispose of in a sanitary landfill, or incineration, or, if allowed by state and local authorities, by burning. If burned, stay out of smoke.

PESTICIDE DISPOSAL

Do not contaminate water, food or feed by storage or disposal. Wastes resulting from use of this product may be disposed of on site or at an appropriate waste disposal facility.

EFFECTIVE AGAINST A BROAD RANGE OF PATHOGENS

HIV-1 (Associated with AIDS)
Staphylococcus aureus
Poliovirus Type I

Mycobacterium tuberculosis
Salmonella choleraesuis
Trichophyton mentagrophytes

Pseudomonas aeruginosa
Herpes simplex Type II

COMPATIBLE IN CONTACT WITH THESE MATERIALS

Aluminum
Epoxy
Nylon
Silicone
Titanium

Buna N
Natural Rubber
Polypropylene
Stainless Steel(304,316,440)
Tygon

Cast Bronze
Neoprene
Rigid PVC
Teflon
Viton

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Virahol® Hospital Disinfectant is designed for effective cleaning, disinfection, and deodorizing of hard inanimate surfaces such as walls, sink tops, tables, chairs, telephones and bed frames.

TO DISINFECT CONTAMINATED SURFACES AND OBJECTS:

For heavily soiled areas, or surfaces of objects soiled with blood or body fluids, a precleaning step is required. Use a clean sponge or lint-free cloth.

On precleaned surfaces, apply Virahol® Hospital Disinfectant with a mop, cloth, sponge or mechanical sprayer so as to thoroughly wet precleaned surfaces. DO NOT DILUTE. Allow surfaces to remain wet for 10 minutes.

Kills HIV-1 on precleaned environmental surfaces & objects previously soiled with blood or body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces, objects with blood or body fluids, and in which the surfaces & objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus (HIV-1) (Associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 ON HARD INANIMATE SURFACES AND/OR OBJECTS SOILED WITH BLOOD OR BODY FLUIDS:

Personal Protection: When handling items soiled with blood or body fluids use disposable latex gloves, gown, mask and or eye coverings.

Cleaning Procedures: Blood and body fluids must be thoroughly cleaned from surfaces and objects before application of Virahol® Hospital Surface Disinfectant/Cleaner.

Disposal of Infectious Materials: Blood and other body fluids should be autoclaved and disposed of according to federal, state, and local regulations for infectious waste disposal.

Contact Time: Apply Virahol® Hospital Surface Disinfectant/Cleaner with a mop, cloth, sponge or mechanical sprayer so as to thoroughly wet precleaned surfaces. DO NOT DILUTE. Allow surfaces to remain wet for 1 minute to insure total inactivation of HIV virus. This contact time will not control other common types of viruses and bacteria. For example, 10 minute contact time is required to control Mycobacterium tuberculosis.

BEST AVAILABLE COPY

ACCEPTED
SEP 16 1993
Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended, for the pesticide registered under EPA Reg. No. 62142-1

272