#### UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

APR 2 9 2002

Albina Otte Veridien Corporation 2875 MCI Drive, Suite B Pinellas Park, FL 33782

RE: "Virahol Hospital Surface Disinfectant Towelette"

EPA Reg. No. 60142-3

Your Submissions Dated 4/5/02 and 4/24/02 (Responding to EPA letter Dated 3/4/02)

The submission referred to above, submitted in connection with registration under the Federal Insecticide, Fungicide, and Rodenticide Act as amended, which responds to the Agency letter dated 3/4/02 by providing a revised label, generic data matrix and data submitters list for isopropanol, is acceptable with the following comment.

The type size on the final labels must be increased in size. The reviewed label has a typesize that is smaller than 6 point. Use larger more readable typesize.

A stamped copy of the acceptable label is enclosed for your records.

If you have any questions about the comments in this letter, please feel free to contact Tony Kish at 703-308-9443, or myself at 703-308-6341.

Sincerely,

Marshall Swindell,

Product Manager Team 33,

Regulatory Management Branch I Antimicrobials Division (7510C)

SURNAME DATE

CONCURRENCES

SYMBOL

SURNAME

OFFICIAL FILE COPY

#### **TUBERCULOCIDE • BACTERICIDE • FUNGICIDE • VIRUCIDE**

## **VIRAHOL® HOSPITAL SURFACE** DISINFECTANT TOWELETTE

FOR DISINFECTING, CLEANING AND DEODORIZING HARD, NON-POROUS **ENVIRONMENTAL SURFACES NON-STAINING** 

#### **KEEP OUT OF REACH OF CHILDREN CAUTION**

SEE SIDE PANELS FOR ADDITIONAL PRECAUTIONARY STATEMENTS.

**ACTIVE INGREDIENT:** ISOPROPANOL...... 70% OTHER INGREDIENTS\* ... 30%

TOTAL . . . . . . . . . . . 100%

\*Does not include the weight of towelette

#### EFFECTIVE AGAINST A BROAD RANGE OF PATHOGENS

TUBERCULOCIDAL: Mycobacterium bovis (BCG) BACTERICIDAL: Staphylococcus aureus [MRSA], ATCC#33591; Vancomycin-resistant Enterococcus faecium (VRE) ATCC#51559; Streptococcus pyogenes, ATCC#19615; Pseudomonas aeruginosa, ATCC#15442; Staphylococcus aureus, ATCC#6538; Salmonella choleraesuis, ATCC#10708; Escherichia coli. ATCC#11229

FUNGICIPAL: Trichophyton mentogrophytes, ATCC#9533; VIRUCIPAL: Herpes simplex type 1, ATCC#VR260; HIV-1 (AIDS Virus) Inactivates HIV-1 on precleaned environmental surfaces and objects previously soiled with blood or body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces or objects with blood or body fluids, and in which the surfaces and objects likely to be sailed with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus (HIV-1) (AIDS virus).

DIRECTIONS FOR USE: It is a violation of Federal Law to use this product in a manner inconsistent

FOR DISINFECTING, CLEANING AND DEODORIZING HARD, NON-POROUS SURFACES:

GENERAL INFORMATION: An effective disinfectant and deadorant towelette for hospitals, clinics, medical and dental offices, laboratories, industrial clean rooms, nursing homes, ambulances, veterinary offices, hotels, restaurants, schools, professional salons and health spas. This product is designed for effective cleaning. disinfecting and deadorizing of hard non-paraus, inanimate surfaces such as walls, sink tops, stainless stee surfaces, equipment, tables, chairs, telephones, tallet seats and bed frames. When using strong cleaners and disinfecting agents, it is always advisable to first test in an inconspicuous area. This product may be used in kitchen preparation areas on counters, sinks, tables, chairs, thermometers and appliances. For surfaces which may come into contact with food a potable water rinse is required. This product may also be used to deodorize areas which generally are hard to keep fresh smelling, such as garage storage areas, empty garbage bins and cons, restrooms and other areas which are prone to adors caused by micro-organisms. This product is not to be used as a terminal sterilant/high level distinfectant on any surface or instrument that [1] is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or [2] contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body.

APPLICATION INSTRUCTIONS: Pull out towelette from dispenser and wipe surface to be disinfected with towelette. Make certain to well surface thoroughly. Surface must remain visibly wet for a full 10 minutes. Use additional towelettes if needed to assure a 10 minute wet contact time. For surfaces which may come into contact with food, a potable water rinse is required. Dispose of towelette in regular trash. SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 ON HARD INANIMATE SURFACES AND/OR OBJECTS SORED WITH BLOOD OR BODY FEIRDS.

Personal Protection: When handling items soiled with blood or body fluids, use disposable latex gloves, gown, mask and/or eye coverings, General Procedures: Blood and body fluids must be thoroughly cleaned from surfaces and objects before application of this product. Contact Time: Surface must remain visibly wet for a full 10 minutes. Use additional towelettes if needed to assure a 10 minute wet contact time. Disposal of Infectious Motertals; Blood and other body fluids should be autoclaved and disposed of according to federal, state, and local regulations for infectious waste disposal. No warranty, either expressed or implied, is given to use contrary to instructions given herein

> ACCEPTED with COMMENTS in and Letter Datoi: APR 2 9 2002

Under the Rederal Intentición. Fungicide, and Poder Scide Act as amended, for the pacticide. registered mater STA flog (to.

60142-3

#### **TUBERCULOCIDE • BACTERICIDE • FUNGICIDE • VIRUCIDE**

#### NET CONTENTS: 160 Towelettes 5.5" X 7.9"

#### PRECAUTIONARY STATEMENTS

#### HAZARD TO HUMANS & DOMESTIC ANIMALS KEEP OUT OF REACH OF CHILDREN

Caution, Causes moderate eye irritation. Harmful if inhaled or absorbed through the skin. Avoid contact with eyes, skin, or clothing. Wash thoroughly with soop and water after handling.

#### FIRST AID

Have the product container or label with you when calling a poison control center or doctor or going for treatment.

If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing the eye. Call a poison control center or doctor for further treatment and advice.

If on skin or dething: Take off contaminated clothing. Rinse skin with plenty of water for 15-20 minutes. Call a poison control center or doctor for further treatment and advice. If Inhelet: Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably mouth-to-mouth, if possible. Call a poison control center or doctor for further treatment and advice.

PHYSICAL OR CHEMICAL HAZARDS: Flammable. Keep away from heat, sparks and open flame. STORAGE AND DISPOSAL: Do not contaminate water, food, or feed by storage and disposal

PESTICIDE STORAGE: Store in a cool, dry place, at room temperature.

CONTAINER DISPOSAL: Do not reuse container. Discard empty container in trash.

PESTIQUE DISPOSAL: Discard used towelettes and empty container in regular trash.

Veridien Corporation 2875 MCI Drive, Suite B Pinellas Park, FL 33782-6105 1-800-345-5444

UPC code

EPA Reg. No. 60142-3 EPA Est. No. 60142-FL-1 (A) EPA Est. No. 073573-ISR-001 (B)

The letter in () following the Est. No. corresponds to the letter suffix of the lot code and identifies which Est. produced this lot.

lot#

# VIRAHOL® HOSPITAL SURFACE DISINFECTANT TOWELETTE

FOR DISINFECTING, CLEANING AND DEODORIZING HARD, NON-POROUS ENVIRONMENTAL SURFACES NON-STAINING

### KEEP OUT OF REACH OF CHILDREN CAUTION

SEE SIDE PANELS FOR ADDITIONAL PRECAUTIONARY STATEMENTS.

ACTIVE INGREDIENT:

ISOPROPANOL..... 70%

OTHER INGREDIENTS\* ... 30%

TOTAL ...... 100%

\*Does not include the weight of towelette.

ACCEPTED with COMMENTS in EPALeHO Dated

APR 2 9 2002

Under the Vertex (1904), site.
Hemmelt is extracted as a second control of the control of the temperature of the control of th

60142-3