



062296-00001-092598

Printed: 10:43:02 Friday, 28 May, 1999 # 61 / 3181

Systems Integration Group, Inc.

PM 34

62296-1

9-25-98

10/6

Cristina Griffin
RBR Production, Inc.
c/o Delta Analytical Corp.
7910 Woodmont Avenue, Suite #1000
Bethesda, MD 20814

SEP 25 1998

Dear Ms. Griffin:

Subject: Change signal word and precautionary statements
Let's Touch
EPA Registration No. 62296-1
Amendment application of March 9, 1998

The amendment referred to above, submitted in connection with registration under FIFRA sec. 3(c)(7)(A), is acceptable provided that you:

1. Submit and/or cite all data or material required for registration/reregistration of your product under FIFRA sec. 3(c)(5) or FIFRA sec.4 when the Agency requires all registrants of similar products to submit such data.
2. Submit two copies of your final printed labeling before you release the product for shipment. Refer to the A-79 Enclosure for further description of final printed labeling.

If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA sec. 6(e). Your release for shipment of the product constitutes acceptance of these conditions.

CONCURRENCES

SYMBOL	7510W						
SURNAME	Jenkins						
DATE	9-25-98						

EPA Form 1320-1A (1/90)


Printed on Recycled Paper

OFFICIAL FILE COPY

2 of 6

A stamped copy of the label is enclosed for you records.

Sincerely your,



Adam Heyward
Product Manager (34)
Regulatory Management Branch II
Antimicrobials Division (7510W)

Enclosure

3066

LOOSEN CAP TO FILL RESERVOIR
DO NOT SQUEEZE BOTTLE WHILE POURING

ISABEL CRISTINA

LET'S TOUCH

HOSPITAL AND TUBERCULOCIDAL METAL SALON DISINFECTANT

ACCEPTED
with COMMENTS
in EPA Letter Dated:

SEP 25 1998

Under the Federal Insecticide,
Fungicide, and Rodenticide Act as
amended, for the pesticide,
registered under EPA Reg. No.

62296-1

For Disinfection and Storage
of All Quality Metal
Instruments Used In Salons

Bactericidal Fungicidal Tuberculocidal
Virucidal Staphylocidal Pseudomonacidal

Active Ingredients:

o-Benzyl-p-chlorophenol 5.25%
o-Phenylphenol 0.99%

Inert Ingredients 93.76%

TOTAL 100.00%

KEEP OUT OF REACH OF CHILDREN

WARNING

STATEMENT OF PRACTICAL TREATMENT

In case of eye contact, hold eyelids open and flush with a gentle stream of water for 15 minutes. Get prompt medical attention. For skin contact, flush with water. If swallowed drink large quantities of water. Do not induce vomiting. See back panel for additional precautionary statements and directions for use.

Net Contents:

ISABEL CRISTINA
THE PERFECT BLEND OF BEAUTY AND SCIENCE
Manufactured by: RBR Productions, Inc.
P.O. Box 3599, Teaneck, NJ 07666
800-247-4130 (outside NJ) 201-313-5700 (in NJ)

EPA Reg. No. 62296-1

EPA Est. No. 8284-PA-1

4076

LET'S TOUCH

Let's Touch is a hospital and tuberculocidal disinfectant concentrate especially formulated for the beauty care profession. It provides broad spectrum protection: bactericidal, fungicidal, virucidal, staphylocidal, pseudomonacidal, tuberculocidal, and prevents cross contamination from all metal salon instruments when used as directed. Use it to disinfect all precleaned metal salon instruments including manicurist nippers and cuticle pushers, shears, and metal skin care instruments. Let's Touch does not rust or dull quality metal instruments, is pH buffered, and is biodegradable.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Before immersion in Let's Touch, all instruments and devices must be thoroughly cleaned to remove excess dirt, rinsed and then rough dried. The removal of heavy soils and the cleaning of all surfaces prior to application of this product is required.

When using Let's Touch for the first time, thoroughly clean and rinse trays and containers in order to remove any residues. Mixing Let's Touch in a container which had another disinfectant in it can cause a brown gel like substance to form. If this happens, thoroughly clean and rinse all contaminated instruments, trays or containers and remix Let's Touch.

Use one part Let's Touch to 32 parts of water for disinfection of precleaned metal instruments and equipment. Let's Touch is distributed in concentrated form in Twin-Neck containers and in 1 oz. packets.

DIRECTIONS FOR MIXING WITH PACKETS AND TWIN-NECK CONTAINER

Water	1 Pint (16 oz.)	1 Quart (32 oz.)	1/2 Gallon (64 oz.)	1 Gallon (128 oz.)
Let's Touch	1/2 fl. oz.	1 fl. oz. or 1 packet	2 fl. oz. or 2 packets	4 fl. oz. or 4 packets

Always add 1 packet or 1 oz. of Let's Touch to every 32 fl. oz. of water. When using the Twin-Neck Container, squeeze the bottle until the easy-measure reservoir is filled to the desired amount.

DIRECTIONS FOR MIXING WITH THE MANICURIST AND HAIRSTYLIST SYSTEMS

Manicurist System

1. Fill the tall plastic mixing jar to the 32 oz. fill line with cold tap water.
2. Add 1 oz. of Let's Touch (1 packet). Cap and mix by turning upside down a few times.
3. Pour the mixed solution into the glass manicurist jar to the suggested use level line.
4. Store the remaining solution until needed.

Hairstylist System

1. Remove the plunger mechanism.
2. Fill the glass jar to the 32 oz. fill line with cold tap water.
3. Add 1 oz. Let's Touch (1 packet). Mix solution by replacing plunger and moving it up and down.

Always leave instruments in Let's Touch solution at all times and overnight for storage. Replace solution daily or earlier if cloudy or dirty.

5046

FOR INSTRUMENT DISINFECTION - After cleaning instruments, place articles in Let's Touch solution for ten minutes. NOTE: Articles must be fully submerged in solution at all times.

FOR METAL INSTRUMENT STORAGE - A solution of one part Let's Touch to 32 parts of cold water (see table above) does not rust or dull, and does not stain or otherwise attack metal instruments such as quality salon shears, manicurist and skin care implements. The solution may be used for prolonged (overnight) storage of all but non metal instruments, plastic or rubber items.

LET'S TOUCH IS A COMPLETE PRODUCT. DO NOT MIX WITH OTHER CHEMICALS. Use only as directed. If frozen, thaw and remix before use.

Let's Touch is recommended for use in salons, hospitals, dental offices or other facilities requiring disinfection of metal instruments.

FOR HOSPITAL AND DENTAL OFFICE USE - Let's Touch is a concentrated disinfectant to be used at a 1:32 dilution with water for cold decontamination/disinfection of pre-cleaned equipment. For disinfection, immerse instruments for 10 minutes.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced with the bloodstream or normal sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

The 1:32 dilution kills pathogenic bacteria and fungi such as S. aureus, S. choleraesuis, Ps. aeruginosa, M. tuberculosis (effective in 10 minutes at 20° C), and C. albicans. Let's Touch is virucidal against Influenza A2 (Hong Kong), Herpes Simplex 1 and 2, Type 2 Adenovirus and Vaccinia on hard, inanimate surfaces in 10 minutes at 20° C.

Minimum use Dilution Confirmation (A.O.A.C.-10 Min.- 20° C) S. choleraesuis 1:32, S. aureus 1:32, Ps. aeruginosa-1:32

A 1:32 dilution is effective against M. tuberculosis on hard, inanimate surfaces in 10 minutes at 20° C. Remove heavy soil or gross filth and thoroughly clean surfaces.

When tested by the EPA-approved Dilution Method the HIV-1 (AIDS) virus was completely inactivated by a 1:32 solution in 10 minutes at 25° C.

**PRECAUTIONARY STATEMENTS
HAZARD TO HUMANS AND DOMESTIC ANIMALS**

WARNING: Causes substantial but temporary eye injury. Do not get in eyes or on clothing. Wear protective eyewear (goggles, face shield or safety glasses). Harmful if swallowed. Avoid prolonged skin contact. Wash thoroughly with soap and water after handling. Remove contaminated clothing and wash clothing before reuse.

STORAGE AND DISPOSAL

Prohibitions: Do not contaminate water, food or feed by storage or disposal. Open dumping is prohibited. Do not reuse empty container.

Pesticide Disposal: This germicide, its solutions or rinsing from empty containers should be disposed of in a toilet or service sink served by a sanitary sewer or in a landfill approved for pesticides.

Container Disposal: Triple rinse and deposit in a waste container for incineration or burial in a landfill approved for pesticide containers.

General: Consult federal, state and local disposal authorities for approved alternative procedures such as limited open burning.

KILLS HIV (AIDS virus) ON PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV (AIDS VIRUS) OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.

Personal Protection: Disposable latex or vinyl gloves, gowns, masks, and/or eye coverings as appropriate must be worn during all cleaning and decontamination procedures of blood and other body fluids.

Cleaning Procedures: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before applying disinfectant.

Disinfectant Use and Contact Time: Effective against HIV-1 (AIDS virus) on hard non-porous surfaces/objects in the presence of organic soil (5% blood serum). Prepare disinfectant by mixing one part Let's Touch to 32 parts water. Leave surfaces wet for 10 minutes.

Disposal of Infectious Materials: Blood and other body fluids should be autoclaved and disposed of according to Federal, State, and local regulations for infectious waste disposal.

.....
.....
.....
.....
.....
.....