

## DIRECTIONS FOR USE:

It is a violation of federal law to use this product in a manner inconsistent with its labeling.

## DIRECTIONS FOR ACTIVATING:

To prepare the activated stock solution for 14-day reuse disinfection and sterilization, add contents of the attached vial to one gallon of the solution. Do not dilute. The shelf life of the activated unused stock solution is 14 days. THE ACTIVATED SOLUTION MAY BE REUSED CONTINUOUSLY FOR FOURTEEN (14) DAYS IN MANUAL (BUCKET OR TRAY) SYSTEMS INCLUDING USE TO KILL MYCOBACTERIUM TUBERCULOSIS.

For use in Hospitals, Nursing Homes, and Health Care Institutions. For disinfection or sterilization of medical instruments, equipment and appliances, the items must be thoroughly cleaned, rinsed, and rough-dried before immersing in the full-strength activated solution. Flush and fill cleansed lumen of any hollow instrument. Avoid entrapment of air bubbles.

**DISINFECTION:** Immerse item completely for 10 minutes at room temperature.

\*\*\*When used as directed for disinfection, Medcide 14 will kill vegetative pathogens such as *Pseudomonas aeruginosa* and viruses (Poliovirus types 1 & 2, Coxsackie B-1, B-3, Echo type 6, Adenovirus type 2, Herpes simplex type 1, Influenza A2, Hong Kong, Vaccinia) on surfaces of inanimate objects. Where a 100% tuberculocidal result is required, immerse items 10 minutes on first use of activated solution and 25 minutes on reused solution at room temperature (25°C).

1 U.S. GALLON (3.785 LITERS)

**Medsafe International, Inc.**

**STERILIZATION:** Immerse item completely for 10 hours at room temperature. Using aseptic technique, remove item from the solution and rinse thoroughly with sterile water.

\*\*When used as directed for 14-day sterilization, Medcide 14 will kill resistant pathogenic spores of *Clostridium sporogenes*, *Bacillus subtilis* and *Clostridium tetani*.

**Directions for use and activation must be followed:**

Record activation and expiration dates on label

**Formulated for the Disinfection and Sterilization of:**  
*Sharp Instruments • Thermometers Oral and Rectal • Metal and Rubber Objects • Lensed Instruments • Plastic Objects • Anesthesia Equipment • Respiratory Therapy Equipment • Flexible Fiberoptic Endoscopes.*

\*MEDCIDE 14™ WILL KILL 100% OF MYCOBACTERIUM TUBERCULOSIS IN 10 MINUTES AT ROOM TEMPERATURE (25°C) ON FIRST USE OF ACTIVATED SOLUTION. WHERE A TUBERCULOCIDAL RESULT IS DESIRED ON REUSED SOLUTION, IMMERSE OBJECTS 25 MINUTES AT ROOM TEMPERATURE (25°C)

## STORAGE AND DISPOSAL

### STORE IN A COOL PLACE

**PROHIBITIONS:** Do not contaminate water, food or feed by storage or disposal. Open dumping is prohibited. Do not reuse empty container.

**CONTAINER DISPOSAL:** Triple rinse (or equivalent), then offer for recycling or reconditioning, or dispose of in a sanitary landfill, or by incineration if allowed by state and local authorities.

**PESTICIDE DISPOSAL:** Pesticide, spray mixture, or rinse water that cannot be used according to label instructions must be disposed of according to federal or approved state procedures under Subtitle C of the Resource Conservation and Recovery Act.

**GENERAL:** Consult Federal, State or Local Disposal authorities

## PRECAUTIONARY STATEMENTS:

### Hazards to Humans and Domestic Animals

**DANGER:** Corrosive, causes burns, causes eye damage and skin irritation. Do not get in eyes, on skin or on clothing. May be absorbed through skin. Wear goggles or face shield and rubber gloves when handling.

**HARMFUL IF SWALLOWED.** Avoid contact with food. Use in ventilated area.

### STATEMENT OF PRACTICAL TREATMENT

**EYES** — In case of eye contact, immediately flush thoroughly with water and seek prompt medical attention.

**SKIN** — In case of skin contact, immediately flush thoroughly with water.

**INGESTION** — If swallowed, drink milk, egg whites, gelatin solution or, if these are not available, large quantities of water. Avoid alcohol. Call a physician.

**NOTE TO PHYSICIAN** — Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed.

Harm can be caused to carbon steel instruments immersed for over 24 hours. Medcide 14™ solution contains a rust inhibitor but will not inhibit rust caused by poor or damaged plating or improper care of instruments.

### NOTE DATES HERE

Activation Date

Expiration Date

**Reusable Disinfecting and Sterilizing Solution**

# Medcide 14

## Activated Dialdehyde Solution

**ACTIVE INGREDIENT:**

Glutaraldehyde	2.0%
INERT INGREDIENTS	98.0%
	100.0%

**DANGER:**

KEEP OUT OF REACH OF CHILDREN. SEE LEFT PANEL FOR ADDITIONAL PRECAUTIONARY STATEMENTS.

**Bactericidal**  
**Fungicidal**  
**Tuberculocidal\***

**Pseudomonacidal**  
**Sporicidal\*\***  
**Virucidal\*\*\***

## **DIRECTIONS FOR USE:**

It is a violation of federal law to use this product in a manner inconsistent with its labeling.

## **DIRECTIONS FOR ACTIVATING:**

To prepare the activated stock solution for 28-day reuse disinfection and sterilization, add contents of the attached vial to one gallon of the solution. Do not dilute. The shelf life of the activated unused stock solution is 28 days. THE ACTIVATED SOLUTION MAY BE REUSED CONTINUOUSLY FOR TWENTY-EIGHT (28) DAYS IN MANUAL (BUCKET OR TRAY) SYSTEMS INCLUDING USE TO KILL MYCOBACTERIUM TUBERCULOSIS.

For use in Hospitals, Nursing Homes, and Health Care Institutions. For disinfection or sterilization of medical instruments, equipment and appliances, the items must be thoroughly cleaned, rinsed, and rough-dried before immersing in the full-strength activated solution. Flush and fill cleansed lumen of any hollow instrument. Avoid entrapment of air bubbles.

**DISINFECTION:** Immerse item completely for 10 minutes at room temperature.

\*\*\*When used as directed for disinfection, Medcide 28™ will kill vegetative pathogens such as *Pseudomonas aeruginosa* and viruses (Poliovirus types 1 & 2, Coxsackie B-1, B-3; Echo type 6, Adenovirus type 2; Herpes simplex type 1; Influenza A2, Hong Kong; Vaccinia) on surfaces of inanimate objects. Where a 100% tuberculocidal result is required, immerse items 10 minutes on first use of activated solution and 25 minutes on reused solution at room temperature (25°C).

**1 U.S. GALLON (3.785 LITERS)**

**Medsafe International, Inc.**

SPOKANE, WASHINGTON 99201

**STERILIZATION:** Immerse item completely for 10 hours at room temperature. Using aseptic technique, remove item from the solution and rinse thoroughly with sterile water.

\*\*\*When used as directed for 28-day sterilization, Medcide 28 will kill resistant pathogenic spores of *Clostridium sporogenes* and *Bacillus subtilis*.

**Directions for use and activation must be followed:**  
Record activation and expiration dates on label.

**Formulated for the Disinfection and Sterilization of:**  
*Sharp Instruments • Thermometers • Oral and Rectal • Metal and Rubber Objects • Lensed Instruments • Plastic Objects • Anesthesia Equipment • Respiratory Therapy Equipment • Flexible Fiberoptic Endoscopes.*

\*MEDCIDE 28™ WILL KILL 100% OF MYCOBACTERIUM TUBERCULOSIS IN 10 MINUTES AT ROOM TEMPERATURE (25°C) ON FIRST USE OF ACTIVATED SOLUTION. WHERE A TUBERCULOCIDAL RESULT IS DESIRED ON REUSED SOLUTION, IMMERSIVE OBJECTS 25 MINUTES AT ROOM TEMPERATURE (25°C).

## **STORAGE AND DISPOSAL** **STORE IN A COOL PLACE**

**PROHIBITIONS:** Do not contaminate water, food or feed by storage or disposal. Open dumping is prohibited. Do not reuse empty container.

**CONTAINER DISPOSAL:** Triple rinse (or equivalent), then offer for recycling or reconditioning, or dispose of in a sanitary landfill, or by incineration if allowed by state and local authorities.

**PESTICIDE DISPOSAL:** Pesticide, spray mixture, or rinse water that cannot be used according to label instructions must be disposed of according to federal or approved state procedures under Subtitle C of the Resource Conservation and Recovery Act.

**GENERAL:** Consult Federal, State or Local Disposal authorities for approved alternatives.

**Reusable Disinfecting and Sterilizing Solution**

# Medcide 28

## Activated Dialdehyde Solution

**ACTIVE INGREDIENT:**

Glutaraldehyde	2.0%
INERT INGREDIENTS	98.0%
	<hr/> 100.0%

**DANGER:**

KEEP OUT OF REACH OF CHILDREN SEE LEFT PANEL FOR ADDITIONAL PRECAUTIONARY STATEMENTS.

**Bactericidal  
Fungicidal  
Tuberculocidal\***

**Pseudomonacidal  
Sporicidal\*\*  
Virucidal\*\*\***