

pm32

63163-1

63163-1

1 of 4

Phenosol

JUL 27 1992

K.C. Laboratories
625 North Avenue West
Missoula, MT 59801

This product is conditionally registered in accordance with FIFRA section 3(c)(7)(A) provided that you:

1. Submit/cite all data or other material required for registration/ reregistration of your product under FIFRA section 3(c)(5) and FIFRA section 4 when the Agency requires all registrants of similar products to submit such data.

2. Add the phrase "EPA Registration No. 63163-1" to your label before you release the product for shipment.

3. The accelerated storage stability data is acceptable on an interim basis. The results of 12 months storage and stability tests should be received by the Agency within 15 months after the granting of registration.

4. Under "First Aid" statement "If swallowed," add, "Drink lots of water."

Submit five (5) copies of your final printed labeling before you release the product for shipment. Refer to the A-79 enclosure for a further description of final printed labeling.

62961:1:WP50:Pringle:WP32-10:KEVRIC:07/14/92:08/10/92

Pringle 7/15/92

2 of 4

If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA section 6(e). Your release for shipment of the product constitutes acceptance of these conditions.

A stamped copy of the label is enclosed for your records.



Ruth G. Douglas
Product Manager (32)
Antimicrobial Program Branch
Registration Division (H7505C)

Enclosures



E.P.A. Registration No. _____ Establishment No. _____

PHENOSOL

INSTRUMENT & EQUIPMENT DISINFECTANT

Non-Rusting	Non-Staining	Pleasant Odor	Economical
Active Ingredients:		Virucidal*	Tuberculocidal
2-benzyl-4-chlorophenol	1.0%	Fungicidal	Bactericidal
2-phenylphenol	9.0%		
Inert ingredients	90.0%		
	100.0%		

*Viruses including Polio virus types 1 & 2, Vaccinia virus, Herpes virus types 1 & 2, Influenza virus type A1 (New Jersey) and Influenza B (Maryland) on hard, non-porous, inanimate surfaces.

KEEP OUT OF REACH OF CHILDREN

DANGER: See back panel for additional precautionary statements.

Net Contents: 1/2 gallon (64 fl ozs) — 1.875 liters

Manufactured by:
KC Laboratories
625 North Ave. West
Missoula, MT 59801

ACCEPTED
with COMMENTS
by EPA, Letter Dated

JUL 27 1982

Under the Fungicide, Insecticide and Rodenticide Act as amended, for the pesticide registered under EPA Reg. No.

63163-1

BEST AVAILABLE COPY

PROCEDURE FOR DISINFECTION:

PHENOSOL should be diluted 1:32 with water for disinfection. All surgical instruments and devices must be thoroughly cleaned to remove excess organic debris, rinsed and rough dried prior to immersion in PHENOSOL. Clean and rinse the lumens of hollow instruments before filling with solution. Immerse in 1:32 concentration of PHENOSOL for a minimum of 10 minutes to destroy the indicated bacteria: *Salmonella typhosa*, *Escherichia coli* 0113, *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Mycobacterium tuberculosis* (If tuberculosis contamination is suspected, thoroughly pre-clean surface prior to application of this product. PHENOSOL is effective in 10 minutes at 20° Celsius), fungi including *Trichophyton interdigitale*. PHENOSOL is an effective fungicide for use on shower room floors, locker room benches and bath mats. PHENOSOL is effective against the following viruses: Polio virus types 1 & 2, Influenza virus type A1 (New Jersey) and Influenza B (Maryland) on hard, non-porous, inanimate surfaces. The ability of this product to destroy the etiologic agent for viral hepatitis has not been established or demonstrated. Rinse disinfected items thoroughly after immersion. Discard the used solution daily or sooner if it becomes soiled or diluted.

For use on metal, glass, rubber and plastics including teflon, polyvinyl chloride, polyethylene, polystyrene, polypropylene and polymethylpentene. Polycarbonate plastic may be adversely affected by prolonged storage in PHENOSOL. Do not store lensed instruments in PHENOSOL.

All label claims of effectiveness are based on previously unused batches of PHENOSOL diluted as directed. Organic material or water introduced into the diluted product may affect the product effectiveness.

PHENOSOL should not be used for the disinfection of hypodermic needles or any item with deep and/or narrow crevices.

IT IS IN VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING

PREPARATION OF 1:32 USE DILUTION:

WATER

1 pint
1 quart
1/2 gallon
1 gallon

PHENOSOL

1/2 oz. (15 ml.)
1 oz. (30 ml.)
2 oz. (60 ml.)
4 oz. (120 ml.)

STORAGE AND DISPOSAL:

DO NOT reuse empty container. Triple rinse with water and dispose in trash. Do not contaminate water, food or feed by storage or disposal. Open dumping is prohibited. Wastes resulting from the use of this product may be disposed of on site or at an approved waste disposal facility. This product should be stored at room temperature, preferably in a cool dry place.

PRECAUTIONARY STATEMENT — HAZARD TO HUMANS AND DOMESTIC ANIMALS

DANGER: Causes eye and skin damage. Corrosive to tissue. Do not get in eyes, on skin or on clothing. Wear goggles or face shield and rubber gloves when handling. Harmful or fatal if swallowed.

FIRST AID:

In case of contact, immediately flush eyes or skin thoroughly with water. For eyes, call a physician. If swallowed, avoid alcohol, call a physician. If on clothing, remove contaminated clothing immediately. Wash before reuse.

Note to Physician: Probable mucosal damage may contraindicate the use of gastric lavage.

PHENOSOL is recommended for use in hospitals, dental offices, laboratories, nursing homes, public and health care institutions.

BEST AVAILABLE COPY