

Important Notice:

Before using Sterilex Ultra Liquid B Solution 1 please read the following information.

1. This product is to be used in conjunction with regular testing of dental unit water.
2. The manufacturer of the dental unit should be consulted before use of this product about compatibility of Sterilex Ultra Liquid B Solution 1 with the dental unit.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Always Use Sterilex Ultra Liquid B Solution 1 with Sterilex Ultra Liquid B Solution 2

1. Initial Start-Up Treatment: Treat the Dental Unit Water Lines (DUWL) for three consecutive nights starting on Monday.
2. Routine Treatment: Treat DUWL one night/week (Monday-Thursday recommended) after Initial Start-up treatment.
3. If two or more weekly treatments are skipped, repeat initial start-up treatment.

How To Apply Sterilex Ultra Liquid B Solutions 1 and 2:

1. Add one bottle (3 ounces) of Sterilex Ultra Liquid B Solution 1 and one bottle (3 ounces) of Sterilex Ultra Liquid B Solution 2 in a measuring cup or other empty container and stir. Pour appropriate amount of the combined mixture into an empty external dental unit water bottle. Do not use PET bottles (See note below). After Sterilex Ultra Liquid B Solution 1 and Sterilex Ultra Liquid B Solution 2 have been mixed together use within 24 hours.
2. Run Sterilex Ultra Liquid B mixture through the system until the pink mixture appears at the end of the Air/Water syringes and handpiece lines. Always remove the handpiece. It may be advisable to remove coupler depending on the type of coupler. Contact manufacturer for specific recommendations.
3. Allow the Sterilex Liquid Ultra mixture to sit in the unit overnight (Place the ends of water lines into a sink or cuspidor in case any pink mixture drips overnight).
4. At the beginning of the next workday, discard any remaining Sterilex Liquid Ultra mixture. Rinse the external water bottle with water.
5. Fill the external water bottle with water. Flush each line (Air/Water syringe, handpiece lines) until the bottle is empty.



Sterilex Ultra Liquid B Solution 1

(This is a 2 part liquid product)

Always Use With

Sterilex Ultra Liquid B Solution 2

- Removes* biofilm from dental unit water lines
- Prevents and suppresses formation of biofilm in dental unit water lines
- Kills biofilm bacteria
- Maintains dental unit water line effluent water <500 cfu/ml **

Note: Do Not Use Thin Walled Polyethylene Terephthalate (PET) Bottles. We Recommend Either A High Density Polyethylene Bottle With A Minimum Thickness Of 0.08 Inches Or A High/Low Density Polyethylene Blend With A Minimum Thickness Of 0.14 Inches.

* Laboratory studies have shown >90% removal
 ** Control of bacterial contamination in the DUWL is affected by the quality of incoming water. Sterile water is recommended by the Organization for Safety and Asepsis Procedures (OSAP).

Active Ingredient:	
Hydrogen Peroxide.....	1.05%
Other Ingredients.....	98.95%
	TOTAL 100.00%

KEEP OUT OF THE REACH OF CHILDREN

DANGER

(See side panel for additional precautionary and first aid statements)

STERILEX CORPORATION

11409 Cronhill Drive

Owings Mills, MD 21117

Phone: 1-800-511-1659 Fax: 410-581-8864

NET WEIGHT _____ EPA REG. NO. 63761-9
 LOT NO. _____ EPA EST. NO. 73516-MD-01
 Made in USA

PRECAUTIONARY STATEMENTS

Hazards To Humans and Domestic Animals.


DANGER. Corrosive. Causes irreversible eye damage. Harmful if absorbed through skin. Do not get in eyes, on skin, or on clothing. Wear protective eyewear (goggles, face shield, or safety glasses) and rubber gloves when handling. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, or using tobacco. Remove and wash contaminated clothing before reuse.

FIRST AID

Have the product container or label with you when calling a poison control center or doctor, or going for treatment.
IF IN EYES: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.
IF ON SKIN OR CLOTHING: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by improper storage or disposal.
PESTICIDE STORAGE: Keep product in closed original container when not in use. Store in a cool, dry area away from potential sources of heat, open flames, sunlight or other chemicals. DO NOT store above 30°C (or 86°F).
PESTICIDE DISPOSAL: Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide or rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste Representative at the nearest EPA Regional Office for guidance.
CONTAINER DISPOSAL: Containers less than or equal to 1 gallon: Securely wrap original container in several layers of newspaper and discard in trash. Containers greater than 1 gallon: Triple-rinse then dispose of in sanitary landfill or by other approved state or local procedures.



63761-9

07/14/2005

11