

67470-6

12/18/2008

1097



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
WASHINGTON, D.C. 20460

OFFICE OF PREVENTION, PESTICIDES AND TOXIC SUBSTANCES

DEC 18 2008

Mr. John Jones
Honeywell
101 Columbia Road
Morristown, NJ 07962-103

SUBJECT: Ethylene Oxide
EPA Registration No.: 67470-6
Notification Date: November 14, 2008
EPA Receipt Date: November 20, 2008

Dear Mr. Jones,

This letter acknowledges receipt of your notification under the provisions of FIFRA section 3(c)9 and PR Notice 91-2.

- Update Container Disposal statements per 2007-4
- Delete Non-FIFRA advisory language

General Comments

Based on a review of the submitted materials, the notification to update container disposal statements and delete non-FIFRA advisory language is acceptable and apart of the records on file.

Should you have any questions or comments concerning this letter, please contact Jacqueline McFarlane at (703) 308-6416 or Velma Noble at (703) 308-6233.

Sincerely,

Velma Noble
Velma Noble
Product Manager (31)

CONCURRENTLY Laboratory Management Branch I
Antimicrobials Division (7510P)

SYMBOL							
URL NAME							
DATE							



United States
Environmental Protection Agency
 Washington, DC 20460

Registration
Amendment
 Other

OPP Identifier Number.

Application for Pesticide - Section I

1. Company/Product Number Honeywell/67470-6	2. EPA Product Manager Velma Noble	3. Proposed Classification <input checked="" type="checkbox"/> None <input type="checkbox"/> Restricted
4. Company/Product (Name) Honeywell/Ethylene Oxide	PM# 31	
5. Name and Address of Applicant (Include Zip Code) Honeywell International, Inc. 101 Columbia Road Morristown, NJ 07962-1053 <input type="checkbox"/> Check if this is a new address	6. Expedited Review. In accordance with FIFRA Section 3(c)(3) (b)(i), my product is similar or identical in composition and labeling to: EPA Reg. No. _____ Product Name _____	

Section - II

<input type="checkbox"/> Amendment - Explain below.	<input type="checkbox"/> Final printed labels in response to Agency letter dated _____
<input type="checkbox"/> Resubmission in response to Agency letter dated _____	<input type="checkbox"/> "Me Too" Application.
<input checked="" type="checkbox"/> Notification - Explain below.	<input type="checkbox"/> Other - Explain below.

Explanation: Use additional page(s) if necessary. (For Section I and Section II.)
 This notification adds language required by PR Notice 2007-4 and deletes non-FIFRA advisory language. Time and experience show that the non-FIFRA advisory language is no longer required for equipment maintenance. This notification is consistent with the provisions of PR Notice 98-10 and EPA regulations at 40 C.F.R. Section 152.46, and no other changes have been made to the labeling or the confidential statement of formula of this product. I understand that it is a violation of 18 U.S.C. Section 1001 to willfully make any false statement to EPA. I further understand that if this notification is not consistent with the terms of PR Notice 98-10 and 40 C.F.R. Section 152.46, this product may be in violation of FIFRA and I may be subject to enforcement action and penalties under FIFRA Sections 12 and 14.

Section - III

1. Material This Product Will Be Packaged In:

Child-Resistant Packaging <input type="checkbox"/> Yes* <input type="checkbox"/> No	Unit Packaging <input type="checkbox"/> Yes <input type="checkbox"/> No	Water Soluble Packaging <input type="checkbox"/> Yes <input type="checkbox"/> No	2. Type of Container <input type="checkbox"/> Metal <input type="checkbox"/> Plastic <input type="checkbox"/> Glass <input type="checkbox"/> Paper <input type="checkbox"/> Other (Specify) _____
* Certification must be submitted			

3. Location of Net Contents Information <input type="checkbox"/> Label <input type="checkbox"/> Container	4. Size(s) Retail Container	5. Location of Label Directions <input type="checkbox"/> On Label <input type="checkbox"/> On Labeling accompanying product
6. Manner in Which Label is Affixed to Product <input type="checkbox"/> Lithograph <input type="checkbox"/> Other <input type="checkbox"/> Paper glued <input type="checkbox"/> Stenciled		

Section - IV

1. Contact Point (Complete items directly below for identification of individual to be contacted, if necessary, to process this application.)

Name John W. Jones	Title Global Leader Product Stewardship	Telephone No. (Include Area Code) (973) 455-4779
-----------------------	--	---

Certification
 I certify that the statements I have made on this form and all attachments thereto are true, accurate and complete. I acknowledge that any knowingly false or misleading statement may be punishable by fine or imprisonment or both under applicable law.

2. Signature 	3. Title Global Leader Product Stewardship	6. Date Application Received (Stamped)
4. Typed Name John W. Jones	5. Date November 14, 2008	

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS.

DANGER! EXTREMELY FLAMMABLE LIQUID AND GAS UNDER PRESSURE. MAY CAUSE EXPLOSIVE MIXTURES WITH AIR. CAUSES EYE AND SKIN BURNS. HARMFUL IF INHALED. HARMFUL IF SWALLOWED. MAY CAUSE RESPIRATORY AND NERVOUS SYSTEM DAMAGE. **DANGER!** CANCER HAZARD AND REPRODUCTIVE HAZARD.

EFFECTS OF OVEREXPOSURE: May be fatal if inhaled in high concentrations. May cause irritation of respiratory tract, chest tightness, headache, nausea, vomiting, diarrhea, lightheaded feeling, dizziness, weakness, drowsiness, cyanosis, loss of coordination, convulsions, coma, delayed lung injury (fluid in lungs), immediate or delayed skin irritation and blisters, allergic skin reaction.

OTHER POSSIBLE DELAYED HEALTH EFFECTS: May cause nervous system injury, cataracts, adverse reproductive effects, chromosomal and mutagenic changes, and cancer.

PEL: 1PPM-TWA Ethylene Oxide (OSHA-29CFR1910.1047)

EL: 5PPM-excursion limit, 15 minutes.

ETHYLENE OXIDE

STERILANT-FUMIGANT GAS.

ACTIVE INGREDIENT: ETHYLENE OXIDE (CAS 75-21-8) 100.0%

KEEP OUT OF REACH OF CHILDREN

DANGER PELIGRO

PRECAUCION AL USUARIO: Si usted no lee Ingles, no use este producto hasta que la etiqueta la haya sido explicada ampliamente.

Users must follow the requirements of the OSHA occupational exposure standard for ethylene oxide (29 CFR 1910.1047).

ODOR: Ether-like at high concentrations. Exposure to toxic levels may occur without warning or detection by the user.

PRECAUTIONS: Do not breathe vapor. Do not swallow. Do not get in eyes, on skin, on clothing. Store and use with adequate ventilation in accordance with 29 CFR 1910.1047

PHYSICAL AND CHEMICAL HAZARDS

Ethylene oxide is extremely flammable and reactive. Contents under pressure. Keep away from heat, flame, sparks or hot surfaces. Do not allow sources of ignition near the sterilization/ fumigation area. Use only in closed system. No part of the container may be exposed above 125°F (52°C). Close valve when not in use and when empty. Use in accordance with tag attached to valve. Ground all equipment, including containers, to avoid static sparks.

LEAK: Evacuate area and keep personnel upwind. Use self-contained breathing apparatus and protective clothing, and shut off leak if without risk.

FIRE: Do not extinguish burning gas if flow cannot be shut off immediately. Move container away from fire if without risk. Use water spray or fog nozzle to keep container cool.

Personal Protective Equipment

Some materials that are chemical resistant to this product are butyl rubber

All handlers must wear at a minimum:

- > Long-sleeved shirt and long pants,
- > Shoes plus socks,
- > Chemical-resistant gloves, and
- > The employer shall provide a respirator that is adequate to protect the health of the employee and ensure compliance with all other OSHA statutory and regulatory requirements (including 29CFR 1910.1047 and 29CFR 1910.134), under routine and reasonably foreseeable emergency situations.

When handlers could have eye or skin contact with ethylene oxide or ethylene oxide solutions such as during maintenance and repair, vessel cleaning, or cleaning up spills, they must wear:

- > Chemical-resistant attire, such as an apron, protective suit, or footwear that protects the area of the body that might contact ethylene oxide or ethylene oxide solutions, and
- > Face-sealing goggles, a full-face shield, or a full-face respirator.

1. Follow the respirator manufacturer's user's instructions for changing canisters.
2. Respirators must be fit-tested and fit-checked using a program that conforms to OSHA's requirements (see 29CFR Part 1910.134).

3. Respirator users must be trained using a program that conforms to OSHA's requirements (see 29CFR Part 1910.134).

4. Respirator users must be examined by a qualified medical practitioner to ensure physical ability to safely wear the style of respirator to be worn. A qualified medical practitioner is a physician or other licensed health care professional (PLHCP) who will evaluate the ability of a worker to wear a respirator. The initial evaluation consists of a questionnaire that asks about medical conditions (such as a heart condition) that would be problematic for respirator use. If concerns are identified, then additional evaluations, such as a physical exam, might be necessary. The initial evaluation must be done before respirator use begins. It does not need to be repeated unless the health status or respirator use conditions change (see 29CFR Part 1910.134).

Follow manufacturer's instructions for cleaning/maintaining PPE. If no such instructions for washables exist, use detergent and hot water. Keep and wash PPE separately from other laundry.

User Safety Recommendations

Users should wash hands before eating, drinking, chewing gum, using tobacco, or using the toilet.

Users should remove clothing/PPE immediately if pesticide gets inside. Then wash thoroughly and put on clean clothing.

Users should remove PPE immediately after handling this product. Wash the outside of gloves before removing. As soon as possible, wash thoroughly and change into clean clothing.

ENVIRONMENTAL HAZARDS

Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans, or other waters unless in accordance with the requirements of a National Pollution Discharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office of the EPA.

FIRST AID

IN ALL CASES OF OVEREXPOSURE, GET MEDICAL ATTENTION IMMEDIATELY. CALL THE POISON CONTROL CENTER OR DOCTOR FOR TREATMENT ADVICE.

IF INHALED: Remove exposed person to fresh air, keep warm. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Call the Poison Control Center or doctor for advice even if no symptoms are present. Keep under medical observation. Symptoms may be delayed.

IF IN EYES: Hold eyelids open and rinse slowly and gently with water for 15-20 minutes. Call the Poison Control Center or doctor for advice.

IF ON SKIN OR CLOTHING: Immediately rinse with plenty of water for 15-20 minutes while removing contaminated clothing and shoes. Call the Poison Control Center or doctor for advice. Aerate, wash or clean contaminated clothing and discard leather goods.

IF SWALLOWED: Call the Poison Control Center or doctor for advice. Give at least two glasses of water. Do not induce vomiting. Do not give anything by mouth to an unconscious person.

NOTE TO PHYSICIAN: REFER TO SECTION IV, FIRST AID MEASURES OF THE MSDSs FOR EACH INGREDIENT. To obtain MSDSs, call 1-800-522-8001.

Skin exposure to Ethylene Oxide will commonly result in skin irritation with extensive blister formation. At high concentrations severe conjunctivitis can occur. Irritation of the respiratory tract may occur, but without acute lung edema. Symptoms of systemic intoxication are headache, nausea, vomiting, incoordination, and cardiac irregularities. Treatment is symptomatic.

IN CASE OF EMERGENCY CALL:

1-800-498-5701. Have a copy of the label or the MSDS when calling a poison control center or doctor or going for treatment.

BEFORE USING OR HANDLING THIS PRODUCT YOU MUST ALSO READ AND UNDERSTAND THE HONEYWELL MATERIAL SAFETY DATA SHEET FOR THIS PRODUCT.

FOR INDUSTRIAL USE ONLY. NOT FOR USE IN HOSPITALS OR HEALTH CARE FACILITIES.

DOT/IMO Shipping Name: Ethylene Oxide


Hazard Class: 2.3, (2.1) Honeywell

ID Number: UN1040 INHALATION HAZARD

101 Columbia Rd., Morristown, NJ 07962-1053

EPA Registration No. 67470-6 EPA Establishment No. 67470-AZ-001

NET CONTENTS- BATCH- DO NOT REMOVE THIS LABEL MADE IN USA STB-0100 (11-14-2008)



6007

In order to reduce ambient levels of ethylene oxide, lengthy facility aeration is encouraged. It can reduce potential long-term risks to employees not directly involved in the ethylene oxide applications. Air monitoring should include the entire facility including office space, break areas, and loading/unloading areas

NOTE: It is a violation of Federal Law to use ethylene oxide sterilant/fumigant gas for the fumigation of beehives, airplanes, trains, buses, ships, trucks, trailers, warehouses, or other similar spaces. After August 1, 2008, this product may only be applied to or on spices, dried vegetables or seasonings utilizing an ETO sterilization method that uses a single sterilization chamber to pre-condition and aerate with an alternating vacuum and aeration purging procedure. If you wish to employ an alternative method to that described below, you must contact the Environmental Protection Agency Office of Pesticide Programs for instruction on how to receive authorization.

Place spices in the treatment chamber. Assure that the mixture of ETO and air is compatible with the chamber design, then, introduce into the chamber a concentration of ETO not to exceed 500 mg/L, with a dwell time not to exceed 6 hours. Then evacuate the gas from the chamber using a sequence of not less than 21 steam washes (injections and evacuations) between 1.5 PSIA (27" Hg) and 5.0 PSIA (20" Hg) while maintaining a minimum chamber temperature of 115° F.

B. Ethylene oxide cycle parameters depend on several sterilizing/fumigating variable factors: preconditioning (if any); exposure time; chamber air concentration; ethylene oxide concentration; chamber temperature; humidity level; types and quantities of items to be sterilized/fumigated; packaging; load configuration in the chamber; microbial challenge method; desired level of sterility assurance; and the desired performance of the sterilized; fumigated product and package.

C. The following is a list of ranges for the critical variables which must be in proper relationship for ethylene oxide to be an effective sterilizing/fumigating agent. This information should be considered general, and not as a replacement for detailed information issued by manufacturers.

TEMPERATURES - 70°F TO 150°F

PRE-VACUUM - typically 25 to 28 inches of mercury. Use vacuums compatible with the products and packages to be sterilized/fumigated, and such that explosive atmospheres are never present in the chamber.

MOISTURE - relative humidity of 33% to 80%

GAS CONCENTRATION - 250 mg/L to 1500 mg/L milligrams of ethylene oxide per liter of chamber volume.

EXPOSURE TIME - 45 minutes to 20 hours

POST-VACUUMS - Ethylene oxide is removed from the chamber and vented to an appropriate ethylene oxide capture or destruction device.

AERATION - aerate sterilized/fumigated materials before use. Do not allow any person to enter the chamber or aeration area if such entry will result in exposures to ethylene oxide above the levels established in 29 CFR 1910.1047.

Cycle parameters and post-cycle aeration parameters (temperature, time, air flow-rate) can affect residue levels. The user must determine that the parameters chosen result in goods which comply with applicable Federal and State residue requirements. For residual limits of ethylene oxide on drug products and medical products see 21 CFR 201.1 sub-section (d). For residual limits on agricultural commodities see 40 CFR 180.151.

D. The sterilization/fumigation cycle parameters should be those prescribed by the sterilizer equipment manufacturer. If other cycle parameters are used, the safety and efficacy of the alternate cycle parameters must be validated and are the responsibility of the user.

NEVER USE PARAMETERS WHICH ALLOW FLAMMABLE MIXTURES OF ETHYLENE OXIDE AND AIR TO ENTER THE CHAMBER.

E. Employers in facilities that use ethylene oxide must comply with all of the requirements for ethylene oxide use specified in 29 CFR 1910.1047.

Handwritten marks and symbols, possibly bleed-through from the reverse side of the page, including various small characters and lines.

7097

STT 0100 EO TAG BACK

GENERAL INSTRUCTIONS

1. Always check container valves and relief valves for leaks before moving cylinder into your facility.
2. This container is equipped with an eductor tube for liquid delivery. If ethylene oxide gas is required, use vaporizing equipment.
3. This container has been pressurized with nitrogen to a pressure of 50 psig (3.52kg/cm²) at 70°F (21.1°C). Vapor pressure will be higher if temperature is above 70°F (21.1°C); lower if temperature is below 70°F (21.1°C). Contact supplier if, upon receipt, container pressure is below 50 psig (3.52kg/cm²) at 70°F (21.1°C).
4. Container must be in an upright position when discharging. Cylinders must be secured to prevent falling over.
5. Liquid withdrawal valve (marked "Liquid") is provided with a CGA 510 connection which has left-hand threads.
6. EOX and SS.55 style cylinders and DOT 5P drums are also provided with a CGA 580 inert pressurizing valve (marked "Vent") which has right-hand threads. Do not discharge product from the CGA 580 inert pressurizing valve.
7. Remove protective valve plugs and make sure valve threads are undamaged. Do not attach an ordinary pipe fitting to these valves. The connections to the container valves should be brass CGA 510 and CGA 580 connectors. Use of other metals could cause damage to the brass container valves.
8. All other piping and fittings should be steel or stainless steel fittings and piping capable of withstanding the pressure to be encountered. Do not use rubber, plastics, or copper materials. Install relief devices where liquid can be trapped between valves.
9. Ground all equipment, including containers, to avoid static sparks.
10. Use only spark-proof tools.
11. Use only explosion-proof electrical equipment where ethylene oxide may be present.
12. Install check valves in the discharge line from this container to processing equipment to prevent backflow into the container.
13. Connect the CGA 580 inert pressurizing valve (marked "Vent") to a source of nitrogen using a line equipped with a pressure regulator, safety relief valve and check valve. The source of nitrogen should be used exclusively for ethylene oxide and for no other purpose. Nitrogen pressure must not exceed the service pressure of the container. Never use compressed air or other gases to pressurize the cylinder or drum.
14. To open container valves, turn counterclockwise. The liquid discharge valve is equipped with a handwheel. Do not use a wrench or other leverage device on handwheel. Use a "T" wrench to open the inert pressurizing valve on the EOX and SS.55 cylinders and 5P drums.
15. Use with adequate general and local ventilation. Vapors form from this product and may travel or be moved by air currents and ignited by pilot lights, other flames, smoking, sparks, heaters electrical equipment, static discharges or other ignition sources at location distant from product handling point.
16. Determine the quantity of product withdrawn from this container by using an appropriate scale.

STORAGE AND DISPOSAL

Do not contaminate food, feed, or water by storage and disposal.

PESTICIDE STORAGE

Store according to instructions provided on label and this tag. Store away from heat in an area with adequate ventilation. Do not store in direct sunlight. To minimize polymer growth, ethylene oxide must not be stored in any place where the temperature consistently exceeds 100°F. To control ethylene oxide polymer growth, use all sterilant gas on a first-in, first-out basis.

PESTICIDE DISPOSAL

Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray, or mixture of rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste representative at the nearest EPA Regional Office for guidance.

CONTAINER DISPOSAL

Refillable container. Refill this container with ethylene oxide only. Do not reuse this container for any other purpose. Return container to supplier for reuse. Before returning container to supplier:

- A. Pressure container with nitrogen to 50 psig total pressure at 70°F (21.1°C) and be sure container valves are closed.
- B. Replace valve plugs tightly in valve outlets. If valve plugs are not available, contact supplier.
- C. Check container valves and plugs for leaks prior to shipment. If leaks are detected, contact supplier.

EPA Registration No. 67470-6
101 Columbia Road
Morristown, NJ 07962-1053

STT-0100; 11/14/2008

