

67619-12

6/21/2012

1/16



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
WASHINGTON D C 20460

June 21, 2012

OFFICE OF CHEMICAL SAFETY
AND POLLUTION PREVENTION

J Evelyn Lawson
Federal Registration Specialist
Clorox Professional Product Company
c/o PS&RC, P O Box 493
Pleasanton, CA 94566-0803

Subject (PRIA D-460663)
CPPC TSUNAMI
EPA Registration No 67619-12
Application Dated February 3, 2012
Receipt Dated February 6, 2012

Dear Ms Lawson

This acknowledges the receipt of your Amendment application dated February 3, 2013 in connection with registration under the Federal Insecticide, Fungicide and Rodenticide Act (FIFRA), as amended, is acceptable with comments

Submission and Proposed Changes

To amend the label to include additional new microorganism claims, reduce contact times of microorganisms already present on the label, and adding an exposure-reduction claim to the label Submission included new efficacy data and in citation of data from CPPC Storm (EPA Reg#67619-13) in support for the current amendment

Findings and Comments

Based on the submitted materials the amended changes of the proposed label dated February 3, 2012 are acceptable, with comments as follows

1 P 2- Posting of website Please note that if you wish to retain a reference to the company's website on the label, then please be aware that such a reference transforms the website into labeling under the Federal Insecticide Fungicide and Rodenticide Act sec 2 (p) (2) and then the website is subject to review by the Agency If the website content is false or misleading, the product would be misbranded and its sale or distribution is unlawful to sell or distribute under FIFRA section 12(a)(1)(E) In addition, regardless of whether a website is referenced on your product's label, claims made on the website may not substantially differ from those claims approved through the registration process Although EPA has not yet determined the extent to which it will routinely review company websites, if the Agency finds or if it is brought to our attention that a website contains false or misleading statements or claims substantially differing

from claims approved through the registration process, the website may be referred to the EPA's Office of Enforcement and Compliance Assurance

2 P 3- Under "Special Instruction for Cleaning against *Clostridium difficile* spores" {Wearing PPE is mandatory}, this instruction must read "**Wear appropriate barrier protection such as gloves, gowns, masks and eye covering**"

3 P 4 - Under "**DISINFECTION**", in fifth sentence, change "should" to read "Gross filth must be removed"

4 P 6 - Since HCV virus is considered as a bloodborne pathogen, add a †† notation, and include "†,†† [Human] Hepatitis C [Virus][(HCV)][as bovine diarrhea virus]**" with its ATCC VR-1422] in the "**Bloodborne Pathogens**" list

5 P 9 -Under "**Special Instructions for Use Against HIV-1, HBV and HCV**" In the last sentence, change "should" to "must", to read "**Blood and other body fluid must be autoclaved**"


6 P 12- In Table 1 under **Use Site**, qualify "X-Ray", e g "**X-Ray machine surfaces**"

General Comments

A stamped copy of the label with comments is enclosed Submit one copy of your final printed labeling before distributing or selling the product bearing the revised labeling

This amendment and a copy of this letter have been inserted in your file for future reference

If you have any questions or comments concerning this letter, please contact David Liem at liem.david@epa.gov or call (703) 305-1284

Sincerely,

Monisha Harris
Product Manager - Team 32
Regulatory Management Branch II
Antimicrobials Division (7510P)

Att Accepted stamped Master Label with comments

Note **Bold italicized text** and (italicized text within parentheses) is information for the reader and is not part of the label [Bracketed information and phrase sequence are optional unless specified in 40 CFR 156.10 Variables such as CAPITALIZATION plurals or singular bold and italics are optional unless specified in 40 CFR The term wipe(s) is interchangeable with towel(s) towelette(s) and cloth(s) throughout except for the product brand name which can only be changed by notification or amendment Listed uses and use sites can be used partially or in their entirety Underlined text is new Strike through (text) means removed R0618161

3/16

CPPC Tsunami

ACTIVE INGREDIENT

Sodium Hypochlorite	0.55% [††]
Other Ingredients	99.45%
Total	100.00%

[Available Chlorine 0.52%]
[5200 PPM [Parts Per Million] Available Chlorine]

[††Does not include weight of dry wipe]

**KEEP OUT OF REACH OF CHILDREN
CAUTION:**

Liquid causes moderate eye irritation
See back label for additional precautionary statements
**Basic formula and alternates A02 A04 A06 A08 and A10
through 16** Contains no phosphorus

**Additional statement for single count wipe pouch
See carton for additional directions for use
or Read carton before use**

[Net] [Contents] _____
**Each package type and size will bear the
actual net contents in pounds and ounces**
Premoistened Wipes
[Individually Packaged]

ACCEPTED
with COMMENTS
in EPA Letter Dated

JUN 27 2012

Under the Federal Insecticide
Fungicide and Rodenticide Act as
amended for the pesticide
registered under EPA Reg No

67619-12

Note **Bold italicized text** and (italicized text within parenthesis) is information for the reader and is not part of the label [Bracketed information and phrase sequence are optional unless specified in 40 CFR 156.10 Variables such as CAPITALIZATION plurals or singular bold and italics are optional unless specified in 40 CFR The term wipe(s) is interchangeable with towel(s) towelette(s) and cloth(s) throughout except for the product brand name which can only be changed by notification or amendment Listed uses and use sites can be used partially or in their entirety Underlined text is new Strike through (text) means removed R0618161

9/16

A list of this product's ingredients is available at www.IngredientsInside.com
or A list of [this product's] [**product name**] ingredients is available at
[**insert company website**]

Mfd for & © YYYY Clorox Professional Products Company

1221 Broadway Oakland CA 94612

Made in the USA

EPA Reg No 67619 12

EPA Est No 5813 CA 3 (A8 or H8) GA 1 (A4) GA 2 (VG) IL 2 (VF) 2686

WI 8 (07) 8251 OH 1 (Q6) WI 2 (KW)

WI 3 (KV) WI 5 (26) 41925 WI 1 (TA) 56952 AR 1 (RL) WI 1 (RK) 63650

WI 1 (JZ)

Beginning of batch code indicates Est No



[[Use your phone to] go here and/or scan [[the] tag or [this] code] to learn or get or see [more] uses or information [on how to use [This Product]] or [how to] tips or [how to [use]] videos]

additional optional text

Clorox Customer Assistance 1 800 234 7700 Any Questions?

The experts at Clorox are standing by to answer any of your questions or Questions or Comments? Call [toll free] 1 800 234 7700 Satisfaction Guaranteed[!] At Clorox we have prided ourselves in making consistent quality products [for over 95 or **insert number** years] and we are dedicated to ensuring your total satisfaction with our product If you are not completely satisfied please call the number below or 1 800 234 7700 For MSDS information please go to www.cloroxprofessional.com Visit us at www.cloroxprofessional.com

To re order [call (**insert company name**) at (**insert phone number**)] and/or [go to (**insert website**)]

Commercial Solutions is a [registered] trademark of Clorox Professional Products Company

Clorox is a [registered] trademark of The Clorox Company

insert trademark is a [registered] trademark of Clorox Professional Products Company

PRECAUTIONARY STATEMENTS HAZARDS TO HUMANS & DOMESTIC ANIMALS

CAUTION Liquid causes moderate eye irritation Do not get in eyes or on clothing Avoid contact with clothing Wear protective eyewear Wash thoroughly with soap and water after handling and before eating drinking chewing gum using tobacco or using the toilet [[For sensitive skin or prolonged use] wear gloves]

FIRST AID If in eyes Hold eye open and rinse slowly and gently with water for 15 20 minutes Remove contact lenses if present after the first 5 minutes then continue rinsing Call a poison control center or doctor for treatment advice Have the product container or label with you when calling a poison control center or doctor or going for treatment For further information call 1 800 234 7700 or Questions? Call 1 800 234 7700

PHYSICAL AND CHEMICAL HAZARDS This product contains bleach Do not use this product with other chemicals such as ammonia toilet bowl cleaners rust removers or acid as this releases hazardous gases

For canisters or buckets not intended to be refilled by the consumer

STORAGE AND DISPOSAL

Store this product in a cool [dry] area [away from direct sunlight and heat] to avoid deterioration Dispose of wipe in trash after use Do not flush Nonrefillable container Do not reuse or refill empty container [Triple rinse promptly then] Recycle empty container or discard in trash

For canisters or buckets intended to be refilled by the consumer

STORAGE AND DISPOSAL

Store this product in a cool [dry] area [away from direct sunlight and heat] to avoid deterioration Dispose of wipe in trash after use Do not flush Do not reuse or refill except as described in the directions for use Refill only with this product [Triple rinse promptly then] Recycle empty container or discard in trash

For single count or multi count pouch

STORAGE AND DISPOSAL

Dispose of wipe in trash after use Do not flush Recycle empty pouch or wrapper or bag or discard in trash

Note **Bold italicized text** and (italicized text within parentheses) is information for the reader and is not part of the label [Bracketed] information and phrase sequence are optional unless specified in 40 CFR 156.10 Variables such as CAPITALIZATION plurals or singular bold and italics are optional unless specified in 40 CFR The term wipe(s) is interchangeable with towel(s) towelette(s) and cloth(s) throughout except for the product brand name which can only be changed by notification or amendment Listed uses and use sites can be used partially or in their entirety Underlined text is new Strike through (text) means removed R0618161

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling

Not for cleaning or sanitizing skin Do not use as a diaper wipe or for personal cleansing For surfaces that may come in contact with food a potable water rinse is required [Avoid contact with clothes fabric wood natural rubber and painted and paper surfaces Do not use on silver or chipped enamel] **APPLY CAREFULLY** – This product contains bleach [Refill only with this product] [Always refer to manufacturer's care instructions before using on equipment or devices]

DISPENSER DIRECTIONS Remove lid and discard seal Insert corner of center sheet through flap opening on lid **DO NOT** push finger through opening Replace lid and pull wipe at an angle The next sheet pops up automatically When finished close lid flap to retain moisture

TO OPEN PACKAGE [Open lid and] pull out wipe from center of roll [The] next sheet pops up automatically [When finished] close lid to retain moisture

DISPENSING DIRECTIONS Starter Kit [Pull down ribbed tab and lift up to] Remove lid from bucket [canister] [Use scissors] [to] [Cut] [Rip] [Tear] [Pull] [Open] refill pouch Insert corner of center sheet through flap opening on lid **DO NOT** push finger through opening Replace lid and pull wipe [towel] at an angle The next sheet [wipe] [towel] pops up automatically When finished close lid flap to retain moisture

REFILL [POUCH] Place refill [pouch] in bucket [canister] [Use scissors] [to] [Cut] [Rip] [Tear] [Pull] [Open] refill pouch Place [put] [insert] refill pouch into bucket [canister] **DO NOT** remove wipe roll from refill pouch Insert corner of center sheet through flap opening on lid **DO NOT** push finger through opening Replace lid and pull wipe [towel] at an angle The next sheet [wipe] [towel] pops up automatically When finished close lid flap to retain moisture

[Look for] Product expiration date on back of refill pouch Discard wipe after use To prevent cross contamination do not reuse wipe If bucket becomes contaminated [or] [damaged] replace immediately Replace bucket immediately if tabs below cap become damaged

Special Instructions for Cleaning Prior to ^{with} Disinfection against Clostridium difficile spores **Personal Protection** Wear appropriate barrier protection such as gloves gowns masks ~~and~~ eye covering **Cleaning Procedure** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with [**Product Name**] [**this product**] Cleaning is to include vigorous wiping and/or scrubbing until all visible soil is removed Special attention is needed for high touch surfaces Surfaces in patient rooms are to be cleaned in an appropriate manner such as from right to left or left to right on horizontal surfaces and top to bottom on vertical surfaces to minimize spreading of the spores Restrooms are to be cleaned last Do not reuse soiled cloths **Infectious Materials Disposal** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal

Note **Bold italicized text** and (italicized text within parentheses) is information for the reader and is not part of the label [Bracketed] information and phrase sequence are optional unless specified in 40 CFR 156.10 Variables such as CAPITALIZATION plurals or singular bold and italics are optional unless specified in 40 CFR The term wipe(s) is interchangeable with towel(s) towelette(s) and cloth(s) throughout except for the product brand name which can only be changed by notification or amendment Listed uses and use sites can be used partially or in their entirety Underlined text is new Strike through (text) means removed R0618161

6/16

SUPPLEMENTAL GRAPHICS/LOGOS/SYMBOLS

The following graphics/symbols/logos are optional and may be placed anywhere on the product label container package etc



DISINFECTION

To [clean and] disinfect (and deodorize) hard nonporous surfaces Wipe [hard nonporous] surface to be disinfected Use enough wipes for treated surface to remain visibly wet for the contact time listed below or on label [Rinse] Let air dry Gross filth should be removed prior to disinfecting [If streaking is observed wipe with a clean damp [cloth or] paper towel after appropriate contact time has expired] [For] Killing Clostridium difficile [spores] Pre clean hard nonporous surfaces by removing gross filth [loose dirt debris blood/bodily fluids etc] Wipe [hard nonporous] surface to be disinfected Use enough wipes for treated surface to remain visibly wet for 3 minutes [Rinse] Let air dry

→ must

Pathogen	Bacteria	‡ Viruses	‡‡ Bloodborne Pathogens	TB	Parvoviruses	Clostridium difficile	Mold/Mildew/Fungi
Contact Time	30 sec[onds]	1 min[ute]	1 min[ute]	3 min[utes]	3 min[utes]	3 min[utes]	5 min[utes]

and/or

ORGANISMS	
Bacteria	[Strain and/or ATCC number]
Acinetobacter baumannii*	[ATCC 15308]
Bordetella pertussis*	[ATCC 12743]
Campylobacter jejuni*	[ATCC 29428]
Carbapenem resistant Klebsiella pneumoniae [(CRKP)]*	[ATCC BAA 1705]
Clostridium difficile or C difficile or C diff [(C diff)] [spores]‡ ***	[ATCC 43598]
Community Acquired Methicillin resistant Staphylococcus aureus [CA MRSA]*	[NARSA NRS123] [Genotype USA400]
Community Acquired Methicillin resistant Staphylococcus aureus [CA MRSA]*	[NARSA NRS384] [Genotype USA300]
Enterobacter aerogenes*	[ATCC 13048]
Enterobacter cloacae*	[ATCC 13047]
Enterococcus faecalis vancomycin resistant or Vancomycin resistant Enterococcus faecalis [(VRE)]	[ATCC 51575]
ESBL (Extended Spectrum Beta Lactamase) producing Escherichia coli (ESBL producing E coli)*	[ATCC BAA 196]
ESBL (Extended Spectrum Beta Lactamase) producing Klebsiella pneumoniae (ESBL producing K pneumoniae)*	[CDC 700603]
Escherichia coli O157 H7*	[ATCC 35150]
Escherichia coli [(E coli)]*	[ATCC 11229]
Escherichia coli New Delhi Metallo Beta Lactamase 1 [(NDM 1)] [NDM 1 E coli]*	[CDC 1001728]
Klebsiella pneumoniae*	[ATCC 4352]
Klebsiella pneumoniae New Delhi Metallo Beta Lactamase 1 [(NDM 1)] [NDM 1 K pneumoniae]*	[CDC 1000527]
Legionella pneumophila [(Legionella)]*	[ATCC 33153]
Linezolid resistant Staphylococcus aureus [(LRSA)]*	[NRS 119]

- *30 sec[ond] contact time
- **1 min[ute] contact time
- ***3 min[ute] contact time

‡ Follow [the] Special Instructions for Cleaning Prior to Disinfection [listed on label]

Note ***Bold italicized text*** and (*italicized text within parenthesis*), is information for the reader and is not part of the label [Bracketed] information and phrase sequence are optional unless specified in 40 CFR 156.10 Variables such as CAPITALIZATION plurals or singular bold and italics are optional unless specified in 40 CFR The term wipe(s) is interchangeable with towel(s) towelette(s) and cloth(s) throughout except for the product brand name which can only be changed by notification or amendment Listed uses and use sites can be used partially or in their entirety Underlined text is new Strike through (~~text~~) means removed R0618161

ORGANISMS continued	
Bacteria	[Strain and/or ATCC number]
Listeria monocytogenes [(Listeria)]*	[ATCC 19117]
Methicillin resistant Staphylococcus aureus [(MRSA)]*	[ATCC 33592]
[Y] Multi drug resistant Enterococcus faecium [(MDR E faecium)]*	[ATCC 51559]
Proteus mirabilis*	[ATCC 9240]
Pseudomonas aeruginosa [(Pseudomonas)]	[ATCC 15442]
Salmonella enterica tested as Salmonella choleraesuis*	[ATCC 10708]
Serratia marcescens*	[ATCC 14756]
Shigella dysenteriae*	[ATCC 11835]
Staphylococcus aureus*	[ATCC 6538]
Streptococcus pneumoniae [(Strep)]*	[ATCC 6305]
Streptococcus pyogenes [(Strep)]*	[ATCC 12344]
Vancomycin intermediate resistant Staphylococcus aureus [(VISA)]*	[CDC HIP 5836]
Vancomycin resistant Staphylococcus aureus [(VRS)]*	[NARSA VRS1]
Mycobacterium	
Mycobacterium bovis [(BCG)] (Tuberculosis or TB)***	
Mold, Mildew Fungi	
Aspergillus brasiliensis tested as Aspergillus niger [(mildew)]****	[ATCC 16404]
Candida albicans***	[ATCC 10231]
Trichophyton mentagrophytes [(Athlete s Foot Fungus)]***	[ATCC 9533]

*30 sec[ond] contact time
 **1 min[ute] contact time
 *** 3 min[ute] contact time
 ****5 min[ute] contact time

YFollow [the] Special Instructions for Cleaning Prior to Disinfection [listed on label]

[Y] ***This strain of Enterococcus faecium has shown resistance to the following antibiotics***
 Amikacin
 Ampicillin
 Gentamicin
 Penicillin
 Tobramycin
 Vancomycin]

Note **Bold italicized text** and (italicized text within parenthesis) is information for the reader and is not part of the label [Bracketed] information and phrase sequence are optional unless specified in 40 CFR 156.10 Variables such as CAPITALIZATION plurals or singular bold and italics are optional unless specified in 40 CFR The term wipe(s) is interchangeable with towel(s) towelette(s) and cloth(s) throughout except for the product brand name which can only be changed by notification or amendment Listed uses and use sites can be used partially or in their entirety Underlined text is new Strike through (text) means removed R0618161

8/16

ORGANISMS Continued	
Viruses [Non enveloped][]	[Strain and/or ATCC number]
‡Adenovirus [type 2] [(Strain Adenoid 6)] [(common cause of respiratory infections)]**	[ATCC VR 846]
‡[Human] Hepatitis A [Virus] [(HAV)] **	[(Strain HM 175)]
‡Norovirus or Norwalk virus [(as Feline Calicivirus)]**	[ATCC VR 782]
‡Poliovirus [type 1]**	[ATCC VR 1000] [(Strain Brunhilde)]
‡Rhinovirus [type 37] [(a leading) causative agent of the common cold]**	[ATCC VR 1147] [(Strain 151 1)]
‡Rotavirus [(causative agent of viral diarrhea)]**	[(Strain WA)]
Viruses [Enveloped][]	
‡Avian Influenza A** §	[ATCC VR 2072]
‡ ††HIV [type 1]*	[(Strain HTLV IIIB)]
‡[Human] Coronavirus [(additional) causative agent of the common cold]**	[ATCC VR 740] [(Strain 229E)]
‡ ††[Human] Hepatitis B [Virus] [(HBV)] [(as duck HBV)]	
→ ††[Human] Hepatitis C [Virus] [(HCV)] [(as bovine diarrhea virus)]**	[ATCC VR 1422]
‡[Human] Herpesvirus [type 2] or Herpes simplex virus [type 2] [(causative agent of genital herpes)]**	[ATCC VR 734] [(Strain G)]
‡Influenza A [(representative of the common flu virus)]**	[ATCC VR 544] [(Strain Hong Kong)]
‡Respiratory syncytial virus [(RSV)] [(cause of respiratory infections in infants)]**	[ATCC VR 26] [(Strain Long)]
Bloodborne Pathogens	
‡ ††HIV type 1 [Human Immunodeficiency Virus]*	[(Strain HTLV IIIB)]
‡ ††[Human] Hepatitis B [Virus] [(HBV)] [(as duck HBV)]**	
→ †, †† [HUMAN] HEPATITIS E [VIRUS] [HEV]	[ATCC VR-1422]
Parvoviruses Non enveloped	
Canine parvovirus***	[ATCC VR 2017] [(Strain Cornell)]
Feline parvovirus [(Feline panleukopenia virus)]***	[ATCC VR 648] [(Strain Philips Roxane)]

*30 sec[ond] contact time
**1 min[ute] contact time
***3 min[ute] contact time

§Kills Avian Influenza A virus on precleaned hard nonporous surfaces
‡ Viruses
†† Bloodborne Pathogens

9/16

DISINFECTION continued

Now and/or New [!] and **or & Improved [!]** to be used as a claim descriptor only for the first 6 months of product on shelf

Claims

- [5 minute] broad spectrum no rinse disinfectant for non food contact surfaces
- A 1 10 sodium hypochlorite solution similar to that recommended by the CDC for effective disinfecting against Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores]
- A disinfectant cloth for use on (*insert use site(s) from table 1 5*)
- [A] one step cleaner/disinfectant
- Aids in reducing or in the reduction of cross contamination between treated surfaces
- Assures proper dilution strength product effectiveness and standardizes technique
- Bacteria and/or virus and/or germ fighting [formula]
- [Bactericidal] [Virucidal] [Fungicidal] [Sporicidal] [Tuberculocidal] [Disinfectant]
- Bactericide or Bactericidal
- Bathroom disinfectant wipes
- Broad spectrum disinfectant [with Clorox® Bleach]
- [Cleans and] disinfects (*insert use site(s) from table 1 5*)
- [Cleans and] disinfects [hard nonporous surfaces]
- Cleans deodorizes and disinfects
- Deodorizes by killing the germs that cause odors
- Designed for killing Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores] on surfaces in hospitals
- Destroys [Deactivates] [(Helps) Eliminate] [tough to kill] Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores] [on hospital surfaces]
- Disinfecting formula
- Disinfects [and deodorizes]
- Disinfects and deodorizes by killing germs and their odors
- Disinfects hard nonporous surfaces [throughout the (*insert use site(s) from table 1 5*)]
- Easy and convenient disinfecting [throughout the (*insert use site(s) from table 1 5*)]
- Easy one step cleaning and disinfecting
- Effective against [*Insert any organism from list of organisms*] [and] Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores]
- Effective against or Kills a wide range of bacteria including Staphylococcus aureus Salmonella enterica Pseudomonas aeruginosa [in the presence of organic soil [(5% blood serum)]]
- Effective against or Kills (*insert any organism(s) from table above*) [in the presence of organic soil [(5% blood serum)]]
- Effective or For daily use against Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores]
- [Effectively] disinfects [hard nonporous] [environmental] surfaces
- Eliminates odors at their source [caused by germs or bacteria and/or mildew]
- Eliminates or Reduces odors [caused by germs or bacteria and/or mildew]
- Exhibits effective virucidal activity
- Fight[s] and/or kill[s] and/or Effective against [99 9% [of]] Eschericia coli [(E coli)]
- Fight[s] and/or kill[s] and/or Effective against [99 9% [of]] Salmonella enterica [(Salmonella)]
- Fight[s] and/or kill[s] and/or Effective against Avian flu [virus] [on environmental and/or hard nonporous surfaces]
- Fight[s] and/or kill[s] and/or Effective against cold and/or flu [virus[es]]
- Fight[s] and/or kill[s] and/or Effective against [Human] Hepatitis A [(HAV)]
- Fight[s] and/or Kill[s] and/or Effective against (*insert any organism from EPA approved list*) [in (*insert appropriate contact time for respective organism*)]
- Fight[s] and/or stop[s] and/or help[s] prevent[s] and/or reduce[s] the spread of [99 9% of] [illness causing] bacteria and/or viruses and/or germs on hard nonporous surfaces [in your (*list any use site*)] [this cold and/or flu season] Formulated to kill Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores]
- Fungicide or Fungicidal
- Germicidal disposable wipes
- Germicide or Germicidal
- Just wipe to clean and disinfect
- Kills [bacteria] [viruses] [fungi] Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores]
- Kills bacteria viruses mold mildew
- Kills Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores] in 3 min[utes] [on hard nonporous surfaces]
- Kills or Is effective against [99 9999%] Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores] [in 3 min[utes]] [on hard nonporous surfaces]
- Kills major forms of microbes including the spore form of Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores]
- Kills odor causing bacteria or germs
- Kills or Destroys mold and mildew
- Kills or Effective against bacteria
- Kills or Effective against germs
- Kills or Effective against viruses
- Meets [the disinfection requirements of] OSHA[s] Bloodborne Pathogen Guidelines or Standards
- Mildewcide or Mildewcidal
- Multi purpose disinfectant
- One step cleaner and disinfectant
- One step cleaning and disinfecting [saves time]

Use Sites See Tables 1 5

10/16

DISINFECTION continued

Now and/or New [!] and **or & Improved [!]** to be used as a claim descriptor only for the first 6 months of product on shelf

- One step disinfectant cleaner designed for general cleaning and disinfecting hard nonporous environmental surfaces in health care facilities or (*insert use site(s) from table 1*)
- Parvovirucide or Parvovirucidal
- Prevents or controls mold and mildew growth
- Pseudomonocidal
- Ready to use disinfectant [disinfection] for [killing] Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores]
- [Ready to use] hospital [grade] disinfectant
- Reduces exposure to [Clostridium difficile^y] [Clostridium difficile^y (C diff)] [C difficile^y] [C diff^y] from treated surfaces
- Reduces exposure to Methicillin resistant Staphylococcus aureus [(MRSA)] from treated surfaces
- Reduces the risk of infection and cross contamination
- Renders Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores] inactive
- Staphylocidal
- The answer to your disinfecting needs
- The easy convenient way to clean and disinfect with Clorox® Bleach
- [The] easy and/or convenient way to disinfect
- [This product] controls the spread of germs on most hard nonporous surfaces
- This product meets AOAC efficacy testing requirements or standards for hospital disinfection
- To kill [inactivate] [[help] prevent the spread of] Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores] [on surfaces]
- To decontaminate Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores] laden surfaces
- To fight Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores] outbreaks
- To kill Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores] let stand for five or 3 min[utes]
- To kill Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores][] especially during outbreaks [of Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y] [let stand for three or 3 min[utes]]
- To reduce Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores] on surfaces in [*insert use site from tables 1 5*]
- To reduce transmission of Clostridium difficile^y or C difficile^y or C diff^y [(C diff)]^y [spores]
- Tuberculocide or Tuberculocidal
- Use in public or common places where [bacteria and/or viruses and/or germs may be of concern] [you are concerned about bacteria and/or viruses and/or germs] [on hard nonporous surfaces]
- Use where control of the hazards of cross contamination between treated surfaces is of prime importance
- Virucide‡ or Virucidal‡
- When used according to label directions [this product] reduces exposure to [Clostridium difficile^y] [Clostridium difficile^y (C diff^y)] [C difficile^y] [C diff^y] from contact

Note *Bold italicized text* and (italicized text within parentheses), is information for the reader and is not part of the label. [Brackets], information and phrase sequence are optional unless specified in 40 CFR 156.10. Variables such as CAPITALIZATION, plurals or singular, bold and italics are optional unless specified in 40 CFR. The term wipe(s) is interchangeable with towel(s), towelette(s) and cloth(s) throughout except for the product brand name which can only be changed by notification or amendment. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. Strike through (text) means removed. R0618161

11/16

DISINFECTION continued

The following is only for products with medical devices listed on the label

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body either into or in contact with the bloodstream or normally sterile areas of the body or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

Special Instructions for Use Against HIV 1, HBV and HCV

This product kills HIV 1, HBV and HCV on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus (HIV 1) (associated with AIDS), Human Hepatitis B Virus (HBV) and Human Hepatitis C Virus (HCV).

Special instructions for using this product to clean and decontaminate against HIV 1, HBV and HCV on surfaces/objects soiled with blood/body fluids

Personal Protection When handling items soiled with blood or body fluids, use disposable impervious gloves, gowns, masks and eye coverings.

Cleaning Procedure Blood and other body fluids must be thoroughly cleaned from surfaces and other objects before applying this product.

Contact Time Allow surfaces to remain wet for 1 minute [ute] [rinse and] let air dry. For all other organisms, see directions for contact times.

Disposal of Infectious Materials Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids ~~should~~ be autoclaved and disposed of according to local regulations for infectious waste disposal.

must

12/16

Pandemic 2009 H1N1 Influenza A virus^{###}

Claims

- Respiratory illnesses attributable to Pandemic 2009 H1N1 are caused by influenza A virus This product (**Product Name**) is a broad spectrum hard surface disinfectant that has been shown to be effective against influenza A virus and is expected to inactivate all influenza A viruses including Pandemic 2009 H1N1 (formerly called swine flu)
- This product has demonstrated effectiveness against influenza A virus and is expected to inactivate all influenza A viruses including Pandemic 2009 H1N1 influenza A virus
- This product has demonstrated effectiveness against influenza A virus and is expected to inactivate all influenza A viruses including Pandemic 2009 H1N1 influenza A virus
- This product has demonstrated effectiveness against influenza A virus and is expected to inactivate all influenza A viruses including Pandemic 2009 H1N1 influenza A virus
- Kills Pandemic 2009 H1N1 influenza A virus (formerly called swine flu)
- Kills Pandemic 2009 H1N1 influenza A virus

per Agency guidance dated October 22 2009 Guidance for Testing and Labeling Claims against Pandemic 2009 H1N1 Influenza A Virus (formerly called Swine Flu) www.epa.gov/oppad001/h1n1_guide.html

FARM PREMISE AND POULTRY HOUSE DISINFECTANT [INSTRUCTIONS]

[Directions For Use]

- 1 Remove all poultry or animals and feeds from premises trucks vehicles coops crates and enclosures
- 2 Remove all litter and manure or droppings from floors walls and surfaces of barns pens stalls chutes and other facilities and fixtures occupied or traversed by animals or poultry
- 3 Empty all troughs racks and other feeding and watering appliances
- 4 Thoroughly clean all surfaces with soap or detergent and rinse with water
- 5 Wipe all surfaces using enough wipes to keep surfaces visibly wet for 5 min[utes]
- 6 Wipe all halters ropes and other types of equipment used in handling and restraining animals as well as forks shovels and scrapers used for removing litter and manure
- 7 Ventilate buildings cars boats coops and other closed spaces Do not house livestock or poultry or employ equipment until treatment has been absorbed set or dried
- 8 Thoroughly scrub all treated feed racks mangers troughs automatic feeders fountains and waterers with soap or detergent and rinse with potable water before reuse

CLEANING AND STAIN REMOVAL

To clean and/or remove stains Wipe surface with a wipe [rinse and] let air dry [If streaking is observed wipe with a clean damp [cloth or] paper towel]

Claims

- Cleans
- Removes (*insert stain(s)/soil(s) from table below*)
- Removes stains

Stains/Soils

- | | | | |
|-------------------|------------------------------|---------------------|--------------------------------------|
| • Blood | • Fecal matter | • Grime | • [Other] common soils and/or stains |
| • Body oils | • Fingerprints | • Laboratory stains | • [Other] organic matter |
| • Coffee [stains] | • Food residue[s] or soil[s] | • Mildew stains | • Tea [stains] |
| • Dead skin | • Fruit [stains] | • Mold stains | • Urine |
| • Dirt | • Grease | | |

13/16

ALLERGEN DESTRUCTION

To [clean and] destroy specified allergens Wipe surface Use enough wipes for treated surfaces to remain visibly wet for 1 minute [Rinse] Let air dry [If streaking is observed wipe with a clean damp [cloth or] paper towel after contact time has expired]

Allergens

- Cockroach matter or particles
- Dust mite matter or particles
- Pet dander or dog and cat dander

Claims

- Breaks Down Many Allergens[!]
- Destroys Many Allergens[!]
- Destroys or deactivates dust mite matter or particles cockroach matter or particles and pet dander allergens [found in dust]
- [Penetrates and] breaks down dust mite matter or particles cockroach matter or particles and pet dander allergens [found in dust]
- Reduces dust mite matter or particles cockroach matter or particles and pet dander allergens [found in dust]
- Reduces Many Allergens[!]
- This product kills germs and destroys many allergens while cleaning hard nonporous surfaces

GENERAL CLAIMS

- A solution impregnated in a wiping cloth
- Appropriate or Perfect or Ideal for Terminal Cleaning
- Convenient
- Delivers consistent [amount of solution] concentration every time[!]
- [Designed] [for] [with] [the] Environmental Services Professional[s] [in] [mind][!]
- Disposable
- Does not contain mercury
- Does not contain nonylphenol ethoxylates (NPEs) and/or alkylphenol ethoxylates (APEs)
- Does not contain VOCs
- [Contains] the power of Clorox® Bleach
- Easy and convenient to use
- Easy to handle
- Eliminates airborne particles associated with sprays
- Eliminates [daily] mixing and measuring
- EVS [Environmental Services] [use]
- Expiration Date [*insert date or length of time*]
- [For] general use
- For use on bathroom surfaces
- For use on nursery surfaces
- Fresh [lingering] bleach scent
- [Heavy] cloth like towel
- Large[r] Surface[s]
- Larger wipe covers [more surface area] larger area [with less effort][!]
- Leaves a fresh bleach scent
- Mercury free
- Moist wipes
- No measuring [required]
- No mixing
- Nonylphenol ethoxylate (NPE) and/or alkylphenol ethoxylate (APE) free
- Pre measured
- Pre moistened [towelettes]
- Prevents formulation errors
- Ready to use
- [Refer to] see [Clorox] [Terminal] [Cleaning] protocol for additional [cleaning] information
Reuse[able] container or bucket
Pouch[es]
[Sturdy] Bucket[s]
System
Protocol[s]
12x12 in[ch]
110 count [ct]
Cloth[s]
Towel[s]
Large[r] [wipe][s]
Big[ger]
- Shelf stable formula
- Specifically or Specially designed for the Medical/Dental/Veterinary/Health Care industry
- [Strong] [Durable] cloth like wipes
- [Strong] [Durable] non woven disposable cloth
- [Strong] [Durable] non woven wipe does not leave lint on surfaces [following use]
- Suitable for hospital use
- Suitable for use on *insert surface material(s) from Table 5*
- Suitable for use on [titanium coated] [medical grade] stainless steel
- Suitable for use on hard nonporous inanimate environmental surfaces
- Terminal Cleaning Wipe[s] [Towel][s] Refill[s][able]
- [The] smell of clean
- VOC free
- VOCs below laboratory detection limits

19/16

TABLE 1 Medical

USE SITES

Ambulances or [Emergency Medical] Transport Vehicles	Headboards	Outpatient Clinics
Anesthesia Rooms or Areas	Health Care Settings or Facilities	Outpatient Surgical Centers [(OPSC)]
[Assisted Living or Full Care] Nursing Homes	Home Health Care [Settings]	Patient Care Areas
[Bed] Head [Foot] Board	Hospices	Patient Restrooms
(blood) (plasma) (semen) (milk) (apheresis) Donation centers	[Hospital] Kitchens	Patient Rooms
BP Monitors	Hospitals	[Pediatric] Examination Rooms or Areas
Carts	Intensive Care Units or ICUs	Pediatric Intensive Care Units [(PICU)]
CAT Lab[oratories]	Isolation Areas	Pharmacies
Central Service Areas	IV Stands	Physicians Offices
Central Supply Rooms or Areas	Laboratories	Physical Therapy Rooms or Areas
Chemotherapy Hoods	[Mayo] [Instrument] Stands	Radiology or X Ray Rooms or Areas
Critical Care Units or CCUs	[Medical] Clinics	Recovery Rooms
Dialysis Clinics	Medical Facilities	Rehabilitation
Dialysis Machines	Medical or Physicians or Doctors Offices	Respiratory Therapy Rooms or Areas
Emergency Rooms or ERs	Neonatal Intensive Care Units [(NICU)]	Sill[s] Ledge[s]
[Eye] Surgical Centers	Newborn or Neonatal Nurseries	Surgery Rooms or Operating Rooms or ORs
Footboards	Nurse Call Device [and Cord]	Transport Vehicles
Glucometers or Blood Glucose Monitors	Nurse Call [Device] [Button] [and Cord]	[TV] Remote[s] [Control][s]
	Nursing or Nurses Stations	X Ray <i>machine surfaces</i> ←
	Ophthalmic Offices	Ultrasound Transducers [and Probes]
	Orthopedics	

SURFACES

bed pans	external surfaces of [medical] equipment or [medical] equipment surfaces	phone cradle
behind and under sinks ^P	external [surfaces of] ultrasound transducers	pipes
cabinet handles	gurneys	plastic mattress covers
cabinets	hard nonporous [environmental] hospital or medical surfaces	reception counters or desks or areas
closet handles	high touch surfaces	shower fixtures
coated mattresses	[hospital or patient] bed railings or linings or frames	showers
coated pillows	[inner] [inside of] drawers	sinks
computer keyboards	IV poles	stretchers
counters	mattress covers	support bars
desk tops	patient chairs	tables
doorknobs		telephones
[exam or examination] tables		wash basins
exterior surfaces of air vents or air vent exteriors		wheelchairs

^Phard nonporous surfaces of

15716

TABLE 2 Dental

USE SITES

Dental Operatories
Dental or Dentists Offices

SURFACES

dental countertops
dental operator surfaces
dentist or dental chairs
hard nonporous [environmental] dental surfaces
light lens covers
reception counters or desks or areas

TABLE 3 Veterinary

USE SITES

Animal Housing Facilities
Animal Life Science Laboratories
[Animal or Pet] Grooming Facilities
Equine Farms
Farms
Kennels
Livestock and/or Swine and/or Poultry Facilities
Pet Areas
Pet Shops or Stores
Small Animal Facilities
Veterinary Clinics or Facilities
Veterinary Offices
Veterinary or Animal Hospitals

SURFACES

animal equipment
automatic feeders
cages
external surfaces of [veterinary] equipment
feed racks
fountains
hard nonporous [environmental] veterinary surfaces
pens
reception counters or desks or areas
stalls
troughs
veterinary care surfaces
watering appliances

TABLE 4 Food Service

USE SITES

Bars
Cafeterias
Commercial or Institutional Kitchens
Delis
Fast Food Chains or Restaurants
Food Preparation and Processing Areas
Food [Service or Processing] Establishments
Food Serving Areas
Other Food Service Establishments
Restaurants
School Kitchens

SURFACES

any washable (food and non food contact) surface where disinfection is required
appliances
dish racks
drain boards
food cases
food trays
freezers
hoods
microwaves
outdoor furniture (excluding wood frames and upholstery)
ovens
refrigerators
salad bar sneeze guards
stoves or stovetops

A potable water rinse is required for food contact surfaces
Do not use on glassware utensils or dishes

16/16

TABLE 5 Miscellaneous/General

USE SITES

Acute Care Facilities
Airplanes
Blood Banks
Bowling Alleys
Butcher Shops
Churches
Colleges
Correctional Facilities
Day Care Centers
Dormitories
Factories
Funeral Homes
Grocery Stores
Gymnasiums or Gyms
Health Club Facilities
Hotels
Industrial Facilities
Laundromats
Laundry Rooms
Locker Rooms
Long Term Care Facilities
Manufacturing Plants or Facilities
Military Installations
Motels
Nursing Homes
Produce Areas
Public Areas
Recreational Centers or Facilities
Restrooms or Restroom Areas
School Buses
Schools
Shelters
Shower Rooms
Storage Rooms or Areas
Supermarkets
Trains
Universities

SURFACES

[bathroom] fixtures
[bath]tubs
behind and under counters^P
behind and under sinks^P
booster chairs
cabinets
ceilings
cell[ular] or wireless or mobile or digital phones
chairs
[computer] keyboards
[computer] monitors
counters or countertops
cribs
desks
[diaper or infant] changing tables
diaper pails
dictating equipment [surfaces]
doorknobs
exterior or external toilet surfaces
exterior or external urinal surfaces
faucets
floors
garbage or trash cans
[grocery store or supermarket] carts
hampers
hand railings
headsets
highchairs
lamps
linoleum
other telecommunications equipment [surfaces]
playpens
shelves
showers or shower stalls
sinks
stall doors
tables
telephones
[tiled] walls
toilet rims
toilet seats
towel dispensers
toys
vanity tops or vanities

SURFACE MATERIALS

[baked] enamel
chrome
[common] hard nonporous [household or environmental] surfaces
Formica
glass
glazed ceramic [tile]
glazed porcelain
laminated surfaces
Marlite
plastic [laminate]
porcelain enamel
stainless steel
synthetic marble
vinyl [tile]
similar hard nonporous surfaces except for those excluded by the label

Not Recommended For Use On
or Avoid contact with
brass
chipped enamel
clear plastic
clothes
copper
fabric
natural marble
painted surfaces
paper surfaces
rubber
sealed granite
silver
unfinished wood
wood

^Phard nonporous surfaces of