



## OFFICE OF CHEMICAL SAFETY AND POLLUTION PREVENTION

WASHINGTON, D.C. 20460

February 23, 2026

### **SENT BY EMAIL**

Heather Bjornson  
heather.bjornson@sterile.com  
VELTEK ASSOCIATES INC.

Subject: Labeling Notification per Pesticide Registration Notice (PRN) 98-10 - Add the statement, "Processed, Filled, Packaged, Sterilized and Tested in the USA. Dry Wipes made in South Korea."  
Product Name: STERI-PEROX 6% WIPE  
Admin Number: 68959-15  
EPA Receipt Date: 09/03/2025  
Action Case Number: 00666404

Dear Heather Bjornson:

The U.S. Environmental Protection Agency is in receipt of your application for notification under Pesticide Registration Notice 98-10 for the above referenced product. The EPA has conducted a review of this request for its applicability under PRN 98-10 and finds that the action requested falls within the scope of PRN 98-10.

The labeling submitted with this application has been stamped "Notification" and will be placed in our records.

Should you wish to add/retain a reference to your company's website on your label, then please be aware that the website becomes labeling under the Federal Insecticide, Fungicide, and Rodenticide Act and is subject to review by the EPA. If the website is false or misleading, the product will be considered to be misbranded and sale or distribution of the product is unlawful under FIFRA section 12(a)(1)(E). 40 CFR § 156.10(a)(5) lists examples of statements the EPA may consider false or misleading. In addition, regardless of whether a website is referenced on your product's label, claims made on the website may not substantially differ from those claims approved through the registration process. Therefore, should the EPA find or if it is brought to our attention that a website contains statements or claims substantially differing from statements or claims made in connection with obtaining a FIFRA section 3 registration, the website will be referred to the EPA's Office of Enforcement and Compliance Assurance.

If you have questions, please contact Srinivas Gowda by telephone at (202) 565-0078 or via email at [gowda.srinivas@epa.gov](mailto:gowda.srinivas@epa.gov).

Sincerely,

*Steven Snyderman*

Steven Snyderman, Product Manager, Team 33  
RMB 2, AD  
Office of Pesticide Programs

[FRONT PANEL]

**NOTIFICATION**

68959-15

The applicant has certified that no changes, other than those reported to the Agency have been made to the labeling. The Agency acknowledges this notification by letter dated:

02/23/2026

**STERI-PEROX® 6% WIPE**

**WFI Formula**

**For Use in Cleanrooms and Controlled Areas**

**Hard Surface Disinfectant**

[Formulated with USP Water for Injection Quality Water (0.25 EU/mL)]

(Aseptically Filtered at 0.2 µm) (Sterile)

**Active Ingredients:**

Hydrogen Peroxide CAS# 7722-84-1.....6.4%

**Other Ingredients:**

\*Water .....93.6%

**Total:**.....100.0%

\*USP Water for Injection Quality Water

**KEEP OUT OF REACH OF CHILDREN**

**CAUTION**

*{Note to reviewer: The following reference statements available for various packaging scenarios.}*

[See [side] [other] [back] [case label] [panel] [and/or] [package insert] for [first aid] [precautionary statements] [directions for use] [additional precautions] [and] additional product information].]

[Refer to the primary case label for [additional] precautionary statements, directions for use, and additional product information]

[Peel [here] [back] for [Spanish and] additional information]

[See container for lot number and storage expiration date]

	{Lot Label Sticker}	
Reorder {xxxxxxxxxxxx}		{QR Code}
Lot: {xxxxxxx}		
EXP. {dd-MM-YYYY}		
Net Contents: {XXX (Individually Packaged) Wipes, Size, Weight}		

**EPA Reg. No.:** 68959-15

**EPA Est. No.:** For EPA Est. No. see the final character of the lot #: L = 68959 –PA-01, C = 68959-PA-03, T = 68959-PA-02

Patents: [www.sterile.com/patents](http://www.sterile.com/patents)

**Manufactured by:**  
Veltek Associates, Inc.  
15 Lee Blvd.  
Malvern, PA 19355-1234  
Tel. (610) 644-8335

Red dot indicates Sterile {will not appear on non-sterile packages} Rev. MM#/YY ●{code}

**(Made in USA)**      **VAI** {corporate logo}

<b>FIRST AID</b>	
<b>If in eyes</b>	<ul style="list-style-type: none"><li>• Hold eye open and rinse slowly and gently with water for 15-20 minutes.</li><li>• Remove contact lenses, if present, after the first five minutes, then continue rinsing eye.</li><li>• Call a poison control center or doctor for treatment advice.</li></ul>
<b>If on skin or clothing</b>	<ul style="list-style-type: none"><li>• Take off contaminated clothing.</li><li>• Rinse skin immediately with plenty of water for 15 – 20 minutes.</li><li>• Call a poison control center or doctor for treatment advice.</li></ul>
<b>If swallowed</b>	<ul style="list-style-type: none"><li>• Call a poison control center or doctor immediately for treatment advice.</li><li>• Have person sip a glass of water if able to swallow.</li><li>• Move person to fresh air.</li><li>• Do not induce vomiting unless told to by a poison control center or doctor.</li><li>• Do not give anything to an unconscious person.</li></ul>
Have the product container or label with you when calling the Poison Control Center or doctor or going to treatment. For general information on product use, etc., call the National Pesticides Information Center (NPIC) at 1-800-858-7378, Monday – Friday, 8:00 am – 12:00 pm Pacific Time; email: <a href="mailto:npic@ace.orst.edu">npic@ace.orst.edu</a> ; or website: <a href="http://npic.orst.edu/">http://npic.orst.edu/</a> . You may also contact the Poison Control Center at 1-800-222-1222 for emergency medical information.	

**EMERGENCIES:** In the case of emergency call: 1-866-928-0789 (toll free) or 1-215-207-0061 (Americas), +52 55 5004 8763 (México)

**[BACK / SIDE PANELS]**

### **PRECAUTIONARY STATEMENTS**

#### **HAZARDS TO HUMANS AND DOMESTIC ANIMALS**

**CAUTION.** Causes moderate eye irritation. Harmful if swallowed or absorbed through skin. Avoid contact with skin, eyes or clothing. Wear appropriate protective eyewear such as goggles or face shield, protective clothing and chemical resistant gloves. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, or using tobacco or using the toilet. Remove and wash clothing contaminated before reuse. Prolonged or frequently repeated skin contact may cause allergic reactions in some individuals.

**PHYSICAL AND CHEMICAL HAZARDS:** Strong oxidizing agent. Keep away from combustible materials. Do not mix this product with other chemicals, organic substances, or other cleaners. Incompatible with acids, alcohols, alkali, powdered metals, metal salts, reactive chemicals, organic materials and reducing agents.

## DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

*This product is NOT to be used as a terminal sterilant / high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.*

**LOCATIONS OF USE:** Cleanrooms and controlled areas such as those in health care institutions, biopharmaceutical, pharmaceutical, medical device and diagnostic manufacturing facilities. Use on hard non-porous, inanimate, surfaces in aseptic filling and gowning rooms, general manufacturing areas and laboratories or on: machinery, tools, tables, counters, laminar-flow benches, floors, walls, carts, shelves, made of plastic, glass, Corian™, vinyl, glazed porcelain, laminates, glazed tile, and stainless steel.

STERI-PEROX 6% WIPE are an effective hard surface disinfectant against *Pseudomonas aeruginosa* (ATCC# 15442), *Staphylococcus aureus* (ATCC#6538) in 10 minutes.

### To Disinfect Hard Nonporous Surfaces:

Pre-clean all surfaces or pre-clean with STERI-PEROX 6% WIPE and dispose.

Using fresh STERI-PEROX 6% WIPE, thoroughly wet surface to disinfect. Use as many wipes as needed to ensure the treated surface remains visibly wet for 10 minutes. Allow surface to air dry, or after 10 minutes, wipe dry if needed.

### Multi Pack

Open Package. Pull out wipes as needed and reseal. When not in use keep sealed to prevent moisture loss. (Peel (and) (or) Reseal Here).

### Individual Wrapped Pack

Tear open packet. Unfold.

### Storage and Disposal

Do not contaminate water, food, or feed by storage and disposal.

**Pesticide Storage:** Store in original container in a cool, well-ventilated area. Do not store near heat or open flame.

**Pesticide Disposal:** Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide or rinsate from container is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your state Pesticide or Environmental Control Agency, or the Hazardous Waste representative at the nearest EPA Regional Office for guidance.

#### Container Disposal:

Nonrefillable container. Do not reuse or refill this container. Offer for recycling, if available.

#### Wipes Disposal:

Nonreusable one-time use container. Discard used and unused wipes in trash. Do not flush wipes down toilet.

**Optional Label Claims:**

Disinfectant  
Hard Surface Disinfectant  
Removes Residues (from) (disinfection)  
Leaves (Low) (No) Residue  
Ready to Use  
Ideal for cleaning and disinfecting hard, non-porous surfaces  
May be used daily  
Meets USP<797> standards

Processed, Filled, Packaged, Sterilized and Tested in the USA. Dry wipes made in South Korea.

{Note: To appear as label end panel for customer use on containers: convenience item.}

<b>User Information</b>
Date Opened/Operator: <input type="text"/>
Expiration Date: <input type="text"/>
Area: <input type="text"/>