JUN - 9 2008

Dear Mr. Fenn,

The following amendment submitted in connection with registration under FIFRA, as amended, is acceptable.

- Increase active ingredient concentration

General Comments

The Confidential Statement of Formula for alternate formulation #1 dated May 28, 2008 is acceptable. It is in compliance with PR Notice 91-2 and in agreement with the label.

A stamped copy of the accepted label is enclosed. Submit three (3) copies of your final printed labeling before distributing or selling the product bearing the revised labeling.

Should you have any questions regarding this letter, please contact Jacqueline McFarlane at (703) 308-6416.

Sincerely,

Sincerely,

CONCURRENCES	
SYMBOL	Valma Noble
SURNAME	Product Manager (31)
DATE	Regulatory Management Branch I
	Antimicrobials Division (7510C)

ACCEPTED
with COMMENTS
in EPA Letter Dated:

JUN - 9 2008

Under the Federal Insecticide,
Fungicide, and Rodenticide Act as
amended, for the pesticide,
registered under EPA Reg. No.

69340-1

ONE BOX CONTENTS: 25 AMPOULES, 26 LINER BAGS, 26 CLOSURES

AN 71

ANPROLENE® STERILIZING GAS AMPOULES

Active ingredient: Ethylene oxide 89.4%
Inert ingredients: 10.6%
Total: 100.0%
Net contents 25 ampoules, 0.20 av. oz. (5.59 g) each ampoule

CE 0413

increased active ingredient
in closed total net content

PN671_VN022508.ai 2/26/2008 9:44:11 AM

KEEP OUT OF REACH OF CHILDREN • DANGER

PRECAUTIONARY STATEMENTS

Hazard to Humans and Domestic Animals.

DANGER

CAUSES IRREVERSIBLE EYE DAMAGE AND SKIN BURNS. HARMFUL IF INHALED. DO NOT BREATH VAPOR. DO NOT GET ON EYES, SKIN OR CLOTHING. DO NOT SWALLOW. CANCER HAZARD AND REPRODUCTIVE HAZARD. MAY CAUSE NERVOUS SYSTEM DAMAGE. STORE AND USE WITH ADEQUATE VENTILATION IN ACCORDANCE WITH 29 CFR 1910.1047

See user's manual for additional
precautionary statements.

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IF IN EYES: Hold eyes open and rinse slowly and gently with water for 15 - 20 minutes. Remove contact lenses, if present, after the first 5 minutes then continue rinsing eyes. Get medical treatment.

IF INHALED: Get exposed person to fresh air. Keep warm. If not breathing call 911 or an ambulance, then give artificial respiration, preferably by mouth - to - mouth, if possible. If breathing is difficult, give oxygen. Call a physician even if no symptoms are present. Keep under medical observation. Symptoms may be delayed.

IF ON SKIN: Immediately wash skin with plenty of water, while removing contaminated clothing and shoes. Call a physician. Aerate, wash, or clean contaminated clothing. Discard leather goods.

IF SWALLOWED: Drink at least two glasses of water. Do not induce vomiting. Do not give anything by mouth to an unconscious victim. Call a physician.

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ENVIRONMENTAL HAZARDS

Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans, or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office of the EPA.

PHYSICAL OR CHEMICAL HAZARDS DANGER FLAMMABLE LIQUID AND GAS UNDER PRESSURE.

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It is a violation of Federal law to use this product in a manner inconsistent with its labeling. It is the responsibility of the employer of any person engaged in the handling or application of this product to follow the requirements of 29 CFR 1910.1047. Read package insert for complete sterilization instructions and additional precautions.

IMPORTANT: For Professional Use Only

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FOR USE IN ANPROLENE STERILIZERS ONLY

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STORAGE: Store in a cool place and out of direct sunlight.

PESTICIDE DISPOSAL: Unwanted or expired ampoules should be returned to the manufacturer for disposal. Contact Andersen Products, Inc. (tel. no. 336-376-3000) for instructions. If unwanted or expired ampoules cannot be disposed of according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste Representative at the nearest EPA Regional Office for guidance.

The AN71 Refill Kit contains 25 AN71 ampoules.

CONTAINER DISPOSAL: Do not reuse empty box, empty ampoules, or liner bags. Wrap box, ampoules and bags in paper and discard in the trash.

Manufactured and Distributed by:

Andersen Sterilizers, Inc.
3154 Caroline Drive
Haw River, NC 27258 USA
336-376-8622

Distributed by:

Andersen Products, Inc.
3202 Caroline Drive
Haw River, NC 27258 USA
1-800-523-1276 • 336-376-3000

European Distributor:

H.W. Andersen Products, Ltd.
Davy Road • Clacton-on-Sea
Essex C015 4XA UK

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ENGLISH

USE ONLY IN ANPROLENE STERILIZERS

IMPORTANT: Anprolene gas is flammable. Do not use near fire, heated surfaces, or flame. Do not open before the 14 hour cycle time has elapsed. Aerate gas absorbing materials outside of Liner Bag according to enclosed instructions. FAILURE TO ADEQUATELY AERATE MAY LEAD TO CONTACT CHEMICAL BURNS.

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ENGLISH

DO NOT OPEN THIS BAG

UNROLL BAG AND LAY FLAT BEFORE BREAKING AMPOULE
CONTENTS FLAMMABLE

ANPROLENE® STERILIZING GAS

Read insert instructions for additional precautions

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**AN
73**

STERILIZING GAS AMPOULES
KEEP OUT OF REACH OF CHILDREN
DANGER

IMPORTANT: For Professional Use Only.

Users must follow the requirements of the OSHA occupational exposure standard for ethylene oxide (29 CFR 1910.1047).

Active ingredient: Ethylene oxide 89.4%
 Inert ingredients: 10.6%
 Total: 100.0%

Contents 60 ampoules, 0.20 av. oz. (5.59 g) each ampoule

EPA REGISTRATION NO. 69340-1 • EPA ESTABLISHMENT NO. 69340-NC-001

Precautionary Statements
 Hazard to Humans and Domestic Animals

DANGER

CAUSES IRREVERSIBLE EYE DAMAGE AND SKIN BURNS. HARMFUL IF INHALED. DO NOT BREATHE VAPOR. DO NOT GET ON EYES, SKIN, OR CLOTHING. DO NOT SWALLOW. CANCER HAZARD AND REPRODUCTIVE HAZARD. MAY CAUSE NERVOUS SYSTEM DAMAGE. STORE AND USE WITH ADEQUATE VENTILATION IN ACCORDANCE WITH 29 CFR 1910.1047

See user's manual for additional precautionary statements.
 Emergency Contact Number: 1-800-255-3924

FIRST AID STATEMENTS

In all cases of exposure, have the product container or label with you when you are going to treatment facility.
IF IN EYES: Hold eyes open and rinse slowly and gently with water for 15 - 20 minutes. Remove contact lenses, if present, after the first 5 minutes then continue rinsing eye. Get medical treatment.
IF INHALED: Get exposed person to fresh air. Keep warm. If not breathing call 911 or an ambulance, then give artificial respiration, preferably by mouth - to - mouth, if possible. If breathing is difficult, give oxygen. Call a physician even if no symptoms are present. Keep under medical observation. Symptoms may be delayed.
IF ON SKIN: Immediately wash skin with plenty of water, while removing contaminated clothing and shoes. Call a physician. Aerate, wash, or clean contaminated clothing and discard leather goods.

IF SWALLOWED: Drink at least 2 glasses of water. Do not induce vomiting. Do not give anything by mouth to an unconscious victim. Call a physician.

NOTE TO PHYSICIAN

Ethylene oxide is a gas. Skin exposure from contact with fabric, rubber or plastic containing residual ethylene oxide will commonly result in skin burns with extensive blister formation. At high concentrations severe conjunctivitis can occur. Respiratory tract irritation may occur, but without acute lung edema. Symptoms of systemic intoxication are headache, nausea, vomiting, incoordination, and cardiac irregularities. Treatment is symptomatic.

ENVIRONMENTAL HAZARDS

Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans, or other waters unless in accordance with the requirements of a National Pollutant Discharge Eliminations System (NPDES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office of the EPA.

PHYSICAL OR CHEMICAL HAZARDS

DANGER. FLAMMABLE LIQUID AND GAS UNDER PRESSURE.

Ethylene oxide gas is extremely flammable. Do not use near flame, electrical sparks, or hot surfaces. Ground all equipment to prevent static sparks.

AN71 • AN73

ANPROLENE®

ENGLISH:

**BRAND GASEOUS STERILANT FOR
ROOM TEMPERATURE,
ATMOSPHERIC PRESSURE STERILIZATION**

Instructions for using Anprolene gas ampoules

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die Seiten 26-31

Manufactured by:
Andersen Sterilizers, Inc.
Health Science park
3154 Caroline Drive
Haw River, NC 27258 U.S.A.

Distributed by:
Andersen Products, Inc.
Health Science Park
3202 Caroline Drive
Haw River, NC 27258 U.S.A.

H.W. Andersen Products, Ltd.
Davy Road • Gorse Lane Industrial Estates
Clacton-on-Sea • Essex C015 4XA U.K.

Distributed in Asia by:
Andersen Medical (HK), Ltd.
Unit 2, 18/F Westley Square
40 Hoi Yuen Road
Kwun Tong, Hong Kong
Tel:2851-1655

AN71 • AN73

ANPROLENE®

Brand Gaseous Sterilant for Room Temperature,
Atmospheric Pressure Sterilization
Use in Ventilated Anprolene Sterilizers Only

DIRECTIONS FOR USE

To be used only in Anprolene sterilizers, and only for hospital, medical, and veterinary sterilization.

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. It is the responsibility of the employer of any person engaged in the handling or application of this product to follow the requirements of 29 CFR 1910.1047.

Do not open the plastic bag in which each individual Anprolene ampoule is sealed.

Each ampoule of Anprolene is surrounded by a plastic break shield. The ampoule and shield are sealed in a plastic bag. The plastic bag is a gas diffusion membrane of known permeability whose function is to contain the gas given off by the ampoule and to release it at a controlled rate during the sterilization cycle. The plastic shield around the ampoule prevents the broken glass of the opened ampoule from puncturing the gas-release bag.

Typical products which may be conveniently processed in Anprolene sterilizers

Respirators, corrugated tubing
Bronchoscopes, gastroscopes, fiberscopes of all kinds
Procedure trays
Catheters - plastic, rubber, metal, glass, cloth
Anesthesia equipment - endotracheal tubes, masks
Rubber tubing
Adhesive tape
Bandages, dressing sets (reuse plastic forceps)
Syringes - plastic, rubber, glass, bulb syringes
Gloves - rubber, plastic, cloth
Surgical instruments - steel, chrome plate, brass, plastic
Optical instruments - scopes, cameras, lenses, mirrors
Electrical equipment - whether autoclavable or not
Painted equipment - metal, wood
High-speed steel - drills, burrs, chisels
Airways - plastic, rubber, metal
Fabric - cloth, rubber, plastic
Electric wire - whether autoclavable or not
Dry-cell batteries, battery cases, bulbs

Sutures - plastic, silk, cotton, stainless steel
Thermometers, applicator sticks
Rectal tubes, douche tubes - rubber plastic
Specula - plastic, metal

Preparation of materials for sterilization

Materials to be sterilized by Anprolene must first be meticulously cleaned and towel-dried. Coating of dry protein, such as dry pus, blood, or feces, protect microorganisms and slow the sterilization process. To prevent this possibility, and to be certain that your system meets the highest standard of sterilization, you must:

1. Disassemble and scrub all instruments in detergent and water to the most critical standard of cleanliness possible. If an item cannot be washed in detergent and water, then pre-sterilization humidification is required. This pretreatment must be done in a chamber having 100% relative humidity at a temperature of at least 68° F (20°C) for at least four hours. Such a chamber may be constructed by placing an AN1071 Humidichip® in a liner bag with the item and closing the bag with a twist-tie.
2. Water on instruments at the time of exposure to Anprolene may react with the gas and reduce its effectiveness. Be sure that items to be sterilized are dry before wrapping and processing. Drain or towel-dry instruments. Do not dry them in a hot-air oven.
3. Occlusive caps, plugs, or stylets must be removed from instruments so that the gas can penetrate freely. Hollow-bore needles and plastic or rubber tubing must be open at both ends and free from plugs. Syringes must be packaged disassembled, with the plunger removed from the barrel.
4. You must wrap all items individually, in cloth or paper, in a manner conventional for steam sterilization, or in Seal and Peel® Packaging. Seal and Peel is a transparent, peel-open, extended shelf-life packaging, proven to be compatible with the Anprolene Sterilizing System. Do not pack the liner bag so tightly with cloth or gauze that gas diffusion is slowed.
5. Be certain that all items have been prepared as described above. The ambient relative humidity must be at least 30% for Anprolene processing. The user must verify that this minimum exists before processing begins.

Sterilization Method using AN74i/Lx, AN2000, AN74j, AN74C, AN74D, AN74E and 74m/mx Sterilizers

1. The sterilizer must be located in a clean, well-ventilated area, away from sparks or flame. The

temperature of the room must be maintained at no less than 68°F (20°C).

2. Remove one liner bag from the AN71 or AN73 dispenser box. Prepare the bag to receive the wrapped items to be sterilized by opening it and placing it into the sterilizer container, open end out. Put the wrapped items to be sterilized into the liner bag.
3. Remove one Anprolene ampoule from the AN71 or AN73 dispenser box. Do not open the plastic bag in which it is wrapped. Rather, unroll the plastic bag. Push the ampoule gently to the center of the gas-release bag. Place the ampoule, still sealed in its plastic bag, inside the liner bag along with the material to be sterilized. For 74m/mx sterilizers, skip to paragraph 7 below.
4. Locate the flexible, purge tube that protrudes into the interior of the sterilizer. Place the purge tube into the mouth of the liner bag. Using one of the closures provided in the Anprolene kit, or the black Velcro® strap provided with newer model sterilizers and tie the liner bag around the purge tube.
5. AN74i/ix/j/2000: Press the Purge button which will draw a vacuum on the liner bag for ninety seconds. At the end of ninety seconds, break the ampoule by manipulating the pre-scored neck of the ampoule through the wall of the liner bag. Close the door and select the 12 hour sterilization cycle.
6. AN72C and AN74D/E: Grasp the ampoule through the bag and protective shield and snap the top of the ampoule. Each ampoule is pre-scored around its neck to facilitate this action. Snapping off the top of the ampoule activates it, i.e., releases the sterilant gas within the gas-release bag. Close the door and press the START button on the sterilizer control panel for the 12 sterilization cycle.
7. 74m/mx: These sterilizers do not have a purge tube or an automatic cycle timer. To run a sterilization cycle, follow steps 1 through 3 above. Press out any excess air before closing the mouth of the liner bag and securing it with the twist ties provided. Break the ampoule by manipulating the pre-scored neck of the ampoule through the wall of the liner bag. Close the door to the sterilizer and note the time the cycle started. At the end of 12 hours, open the door and remove the closure from the mouth of the bag and spread the open end of the bag to facilitate release of residual ethylene oxide from the liner bag. Close the door and allow the products to aerate for at least 2 hours before removing from the sterilizer.
8. Always lock the door and remove the key from the sterilizer to prevent unauthorized removal of goods

during the sterilization and aeration cycle.

CAUTION: FAILURE TO OBSERVE THE ABOVE INSTRUCTIONS MAY EXPOSE THE STERILIZER OPERATOR TO MORE THAN 0.5 PPM ETHYLENE OXIDE (The OSHA Established Action Level Measured as an Eight-Hour Time Weighted Average).

After the sterile materials have been removed from the liner bag, the exhausted and empty gas-release ampoule and bag along with the liner bag may be disposed of in ordinary rubbish.

Special processes

When sterilizing long lengths of plastic or rubber tubing, such as urethral catheters, heart catheters, and coils of drainage tubing, use two ampoules and double the sterilization cycle time to 24 hours. For your particular sterilizer, select the extended cycle time as follows:

AN74i/ix/j/2000: Select the 24 hour cycle instead of the 12 hour cycle.

AN72C and AN74D/E: Press START and 24 HOUR CYCLE ONLY buttons simultaneously to run a 24 hour cycle.

74m/mx: Allow 24 hours to elapse between breaking of the ampoule and the time when you open the liner bag and commence aeration.

Aerate the gas-absorbent items for at least 24 hours before use (see: Precautionary Statements).

Technical description of the system

The Anprolene Sterilizer container serves as a guard against inadvertent ignition of the contents by spark or open flame during the sterilization cycle. The cabinet also serves as a ventilation hood.

The liner bag, when properly closed with a tie or Velcro strap, has a 7 liter capacity. It is a second gas-diffusion membrane that serves to retain Anprolene long enough to sterilize the liner bag contents and releases it into the surrounding ventilated cabinet at a controlled rate. Each ampoule releases approximately 4.5 grams of Anprolene at room temperature (68°F/20°C) and sea-level atmospheric pressure.

The normal AN71/73 ampoule produces a peak concentration within the liner bag of 500 mg / L * Hr. Laboratory tests confirm that this dose will kill the most resistant spores known within the twelve-hour cycle, providing they have been re-hydrated according to our instructions (see: Preparation of Materials for Sterilization).

Storage and disposal

Do not contaminate water, food, or feed by storage or disposal.

STORAGE: AN71/73 Anprolene must be stored in a cool place out of direct sunlight. Under normal conditions its shelf life is five years. As long as the content of the ampoule is liquid at 68°F (20° C), it is sufficiently potent to use.

PESTICIDE DISPOSAL: Unwanted or expired ampoules should be returned to the manufacturer for disposal. Contact Andersen Products, Inc. (tel. no. 336-376-3000) for instructions. If unwanted or expired ampoules cannot be disposed of according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste representative at the nearest EPA Regional Office for guidance.

CONTAINER DISPOSAL: Do not reuse empty box, used ampoules, or liner bags. Wrap used box, ampoules, and liner bags in paper and discard in the trash.

Testing the efficiency of the Anprolene Sterilizer

Monitoring sterilization efficacy is extremely important. It is strongly recommended that the user establish a routine for monitoring each cycle. A color-change chemical indicator, such as the AN87 Dosimeter®, placed in the most inaccessible part of the load, will indicate whether or not the gaseous sterilant penetrated to the core of the load in adequate concentration to assure sterilization. In addition, an appropriate biological control, such as the AN80 Steritest®, should be used at least once per month to challenge the procedure. The AN80 Steritest is sensitive to the gas concentration, cycle time, and average cycle temperature.

AERATION: It is incumbent on the user to confirm that gas-absorbent items, sterilized in Anprolene, have been adequately aerated before they contact living tissue. This is a particularly important warning to those using plastic or rubber items that directly contact tissue culture preparations, ova, semen or embryos.

Nylon and polyester films are virtually impervious to Ethylene Oxide. The only other plastic waterproof wrapping material proven compatible with Anprolene is our Seal and Peel® brand roll stock. Do not use any other plastic film to wrap items to be sterilized in Anprolene.

Vacuum dehydration, chemical desiccation or prolonged exposure to ambient relative humidity below 30% has been demonstrated to produce spores highly resistant to

sterilization by Anprolene. Re-hydration of spores so changed, and hence reversion to normal sensitivity, does not seem to occur until they have been actually wetted or placed in a 100% relative humidity atmosphere. Do not attempt to sterilize materials which may be carrying dried spores without first washing the articles with water and detergents. If the nature of the material is such that the water treatments specified above are harmful, then pretreatment must be done in a chamber having 100% relative humidity at a temperature of at least 68°F (20°C) for at least four hours. Such a chamber may be made by placing an AN1071 Humidichip® in a liner bag with the item and closing the bag with a twist tie.

PHYSICAL OR CHEMICAL HAZARDS

DANGER: FLAMMABLE LIQUID AND GAS UNDER PRESSURE

Ethylene Oxide gas is extremely flammable. Do not use near flame, electrical sparks, or hot surfaces. Ground all equipment to prevent static sparks.

IMPORTANT - THE USER MUST NOT DEVIATE FROM THESE INSTRUCTIONS

Do not be fooled by the apparent simplicity of the Anprolene Sterilizing System. Its reliable use depends upon precise adherence to these instructions. AN71/73 Anprolene ampoules are designed to be used only in the 72C, 74D/E, 74i/ix, 74m/mx, 74j and AN2000 Anprolene Sterilizers with the liner bags and twist seals supplied in the refill dispenser. No other container may be used, no matter how similar it seems. Liner bags must not be reused; use only a genuine, fresh Anprolene liner bag with each load.

Precautionary Statements

Hazard to Humans and Domestic Animals

DANGER

CAUSES IRREVERSIBLE EYE DAMAGE AND SKIN BURNS. HARMFUL IF INHALED. DO NOT BREATHE VAPOR. DO NOT GET ON EYES, SKIN, OR CLOTHING. DO NOT SWALLOW. CANCER HAZARD AND REPRODUCTIVE HAZARD. MAY CAUSE NERVOUS SYSTEM DAMAGE. STORE AND USE WITH ADEQUATE VENTILATION IN ACCORDANCE WITH 29 CFR 1910.1047

See user's manual for additional precautionary statements.

Emergency Contact Number: 1-800-255-3924

Never open the sterilizer during the sterilization cycle. Avoid breathing Anprolene vapor.

Exposure to liquid or gas can cause severe eye damage. Intimate skin contact with fabric, clothing, rubber, or plastics containing residual Ethylene Oxide can cause damage to the skin. Wear goggles or face shield when handling. May be fatal if swallowed. Do not breathe vapor. Extremely hazardous liquid and gas under pressure. Inhalation may be fatal. Ethylene oxide gas can be odorless at toxic levels; therefore, exposure to toxic levels may occur without warning or capability of detection by user.

EFFECTS OF OVEREXPOSURE TO ETHYLENE OXIDE

May be fatal if inhaled in high concentrations. May cause irritation of respiratory tract, chest tightness, headache, nausea, vomiting, diarrhea, light-headed feeling, dizziness, weakness, drowsiness, cyanosis, loss of coordination, convulsions, coma, delayed lung injury (fluid in lungs), immediate or delayed skin irritation and blisters, allergic skin.

OTHER POSSIBLE DELAYED HEALTH EFFECTS

May cause nervous system injury, cataracts, adverse reproductive effects, chromosomal and mutagenic changes, and cancer.

PEL: 1 PPM-TWA Ethylene Oxide (OSHA - 29 CFR 1910.1047)

STEL: 5 PPM-excursion limit, 15 minutes

ODOR: Ether-like in high concentrations. Exposure to toxic levels may occur without warning or detection by the user.

See label for additional precautionary statements and statements of practical treatment.

Aeration: Failure to adequately air gas-absorbing materials may lead to contact chemical burns. All items that may contact living tissues must be aerated for at least 24 hours, at a minimum temperature of 68°F (20°C), prior to use. Items such as plastic instruments, foam rubber, plastic foams, vinyl tubing, rubber tubing, plastic items, rubber items, and sealed air cushioned devices (for instance, anesthesia masks), must be aerated. Instruments need not be removed from cloth or paper wrappings to obtain adequate diffusion of the residual Anprolene vapors.

Operator Exposure: It is absolutely essential that reliable forced ventilation be employed in the room where the Anprolene Sterilizer is used. This system must be capable

of ten fresh make-up air changes per hour, so that operator exposure remains within U.S. government permissible exposure limits.

FIRST AID STATEMENTS

In all cases of exposure, get medical attention immediately. Take person to a doctor or emergency treatment facility. Have the product container or label with you when going for treatment.

If in Eyes: Hold eyes open and rinse slowly and gently with water for 15 - 20 minutes. Remove contact lenses, if present, after the first 5 minutes then continue rinsing eye. Get medical treatment.

If Inhaled: Get exposed person to fresh air. Keep warm. *If not breathing call 911 or an ambulance, then give artificial respiration, preferably by mouth - to - mouth, if possible.* If breathing is difficult, give oxygen. Call physician even if no symptoms are present. Keep under medical observation. Symptoms may be delayed.

If on Skin: Immediately wash skin with plenty of water, while removing contaminated clothing and shoes. Call a physician. Aerate, wash, or clean contaminated clothing and discard leather goods.

If Swallowed Drink at least two glasses of water. Do not induce vomiting. Do not give anything by mouth to an unconscious victim. Call a physician.

See label or instructions for additional precautionary statements and statements of practical treatment.

NOTE TO PHYSICIAN

Ethylene oxide is a gas. Skin exposure from contact with fabric, rubber or plastic containing residual Ethylene Oxide will commonly result in skin burns with extensive blister formation. At high concentrations severe conjunctivitis can occur. Respiratory tract irritation may occur, but without acute lung edema. Symptoms of systemic intoxication are headache, nausea, vomiting, incoordination, and cardiac irregularities. Treatment is symptomatic.

Ampoule Fill Weights

Active ingredient: Ethylene Oxid.89.4%
Other Ingredients:10.6%
Total:100.0%
Net contents each ampoule0.20 av. oz (5.59 g)

Anprolene Modular Sterilizing System

- AN71 Anprolene 25 Ampoule Refill Dispenser:
Net contents 0.20 av.oz. (5.59g) each ampoule
- AN73 Anprolene 60-Ampoule Refill Dispenser:
Net contents 0.20 av.oz. (5.59g) each ampoule
- AN79 Anprolene 20-Ampoule Refill Dispenser:
Net contents 0.64 av.oz. (18.18g) each ampoule
- AN7916 Anprolene 16-Ampoule Refill Dispenser:
Net contents 0.64 av.oz. (18.18g) each ampoule
- AN72C Ventilated and Purged Anprolene Tray
(For use only with 5.59 gram AN71/73 ampoules)
Sterilizer with automatic cycle controls.
- AN74D/E Ventilated and Purged, High-Capacity
Anprolene Sterilizer with automatic cycle controls.
- AN74i/ix Ventilated and Purged, High-Capacity
Anprolene Sterilizer with automatic cycle controls.
- AN74j Ventilated and Purged, High-Capacity
Anprolene Sterilizer with automatic cycle controls.
- AN2000 Ventilated and Purged, High-Capacity
Anprolene Sterilizer with automatic cycle controls.
- AN74 m/mx Ventilated, High-Capacity Anprolene
Sterilizer with manual cycle controls.

Anprolene Sterilizing Accessories

- AN80 Steritest biological and chemical controls.
- AN85 Anprolene Exposure Indicator Strips.
- AN86 Anprolene Exposure Indicator Tape.
- AN87 Anprolene Dosimeter.
- AN90 Seal and Peel Electric Impulse Sealer.
- AN820 2" x 200' (5 cm x 60 m)
Seal and Peel Roll Stock.
- AN830 3" x 200' (8 cm x 60 m)
Seal and Peel Roll Stock.
- AN850 5" x 200' (13 cm x 60 m)
Seal and Peel Roll Stock.
- AN870 7" x 200' (18 cm x 60 m)
Seal and Peel Roll Stock.
- AN1071 Humidichip (2" x 2" / 5 cm x 5 cm) -
Humidity stabilizing chips.
- AN1072 Humiditube

Anprolene®, Dosimeter®, Humidichip®, Seal and Peel® and Steritest® are trademarks of Andersen Sterilizers, Inc.

Humidichip U.S. patent No. 5,082,636

Manufactured by:
Andersen Sterilizers, Inc.
Health Science Park
3154 Caroline Drive
Haw River, NC 27258 U.S.A.
Tel: 336-376-8622
Fax: 336-376-5428

Distributed by:
Andersen Products, Inc.
Health Science Park
3202 Caroline Drive
Haw River, NC 27258 U.S.A.
1-800-523-1276
www.anpro.com
mailbox@anpro.com

Distributed in Europe, the Middle East and Africa by:
H. W. Andersen Products, Ltd.
Davy Road • Clacton-on-Sea
Essex C015 4XA U.K.
Tel: 44-1-255-428-328
Fax: 44-1-255-222-987

Distributed in Asia by:
Andersen Medical (HK) Ltd.
Unit 2, 18/F Westley Square
40 Hoi Yuen Road
Kwun Tong, Hong Kong

EPA Registration No.69340-1
EPA Establishment No. 69340-NC-001

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Ampoule Fill Weights

Active ingredient: Ethylene oxide89.4%
Inert Ingredients: 10.6%
Total: 100.0%
Net contents each ampoule0.20 av. oz (5.59 g)

Anprolene Modular Sterilizing System

AN71	Anprolene 25 Ampoule Refill Dispenser: use with AN72C & D Anprolene Sterilizers.
AN73	Anprolene 60-Ampoule Refill Dispenser: use with AN72 C & D Anprolene Sterilizers.
AN79	Anprolene 20-Ampoule Refill Dispenser: use with AN74C/D/E/I/ix Anprolene Sterilizers.
AN7916	Anprolene 16-Ampoule Refill Dispenser: use with AN74C/D/E/I/ix Anprolene Sterilizers.
AN72C & D	Ventilated and Purged Anprolene Tray Sterilizer with automatic cycle controls.
AN74C/D/E	Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.
AN74i/ix	Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.
AN74J	Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.
AN2000	Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.

Anprolene Sterilizing Accessories

AN80 Steritest biological and chemical controls.
AN85 Anprolene Exposure Indicator Strips.
AN86 Anprolene Exposure Indicator Tape.
AN87 Anprolene Dosimeter.
AN90 Seal and Peel Electric Impulse Sealer.
AN820 2" x 200' (5 cm x 60 m)
Seal and Peel Roll Stock.
AN830 3" x 200' (8 cm x 60 m)
Seal and Peel Roll Stock.
AN850 5" x 200' (13 cm x 60 m)
Seal and Peel Roll Stock.
AN870 7" x 200' (18 cm x 60 m)
Seal and Peel Roll Stock.
AN1071 Humidichip (2" x 2" / 5 cm x 5 cm) humidity

stabilizing chips.
AN1072 Humiditube (5 in a pack).

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EPA Registration No. 69340-2
EPA Establishment No. 69340-NC-001

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Anprolene Modular Sterilizing System

AN71 Anprolene 25 Ampoule Refill Dispenser:
use with AN72C & D Anprolene Sterilizer
AN73 Anprolene 60-Ampoule Refill Dispenser:
use with AN72 C & D Anprolene
Sterilizers.
AN79 Anprolene 20-Ampoule Refill Dispenser:
use with AN74C/D/E/i/ix Anprolene
Sterilizers.
AN7916 Anprolene 16-Ampoule Refill Dispenser:
use with AN74C/D/E/i/ix Anprolene
Sterilizers.
AN72C & D Ventilated and Purged Anprolene Tray
Sterilizer with automatic cycle controls.
AN74C/D/E Ventilated and Purged, High-Capacity
Anprolene Sterilizer with automatic cycle
controls.
AN74i/ix Ventilated and Purged, High-Capacity
Anprolene Sterilizer with automatic cycle
controls.
AN74J Ventilated and Purged, High-Capacity
Anprolene Sterilizer with automatic cycle
controls.
AN2000 Ventilated and Purged, High-Capacity
Anprolene Sterilizer with automatic cycle
controls.

Anprolene Sterilizing Accessories

AN80 Steritest biological and chemical controls.
AN85 Anprolene Exposure Indicator Strips.
AN86 Anprolene Exposure Indicator Tape.
AN87 Anprolene Dosimeter.
AN90 Seal and Peel Electric Impulse Sealer.
AN820 2" x 200' (5 cm x 60 m)
Seal and Peel Roll Stock.
AN830 3" x 200' (8 cm x 60 m)
Seal and Peel Roll Stock.
AN850 5" x 200' (13 cm x 60 m)
Seal and Peel Roll Stock.
AN870 7" x 200' (18 cm x 60 m)

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Validation services and contract
sterilization available from:

Andersen Scientific, Inc.
3220 Caroline Drive
Haw River, NC 27258 U.S.A.
Tel: 888-933-6427
Fax: 336-376-8153

Distributed in Europe, the Middle East and Africa by:

H. W. Andersen Products, Ltd.
Davy Road • Clacton-on-Sea, Essex
Essex CO15 4XA U.K.
Tel: 44-1-255-428-328
Fax: 44-1-255-222-987

Distributed in France by:

H.W. Andersen, SARL
Batiment A - RDC - Parc d'Activites de l'Epi de Soll
55 rue Salvador Allende
59120 LOOS, France
Tel: 33-03-20-97-05-25
Fax: 33-03-20-97-61-77

Distributed in Asia by:

Andersen Medical (HK) Ltd.
Unit 2, 18/F Westley Square
40 Hoi Yuen Road
Kwun Tong, Hong Kong

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Anprolene Modular Sterilizing System

- AN71 Anprolene 25 Ampoule Refill Dispenser:
use with AN72C & D Anprolene Sterilizers.
- AN73 Anprolene 60-Ampoule Refill Dispenser:
use with AN72 C & D Anprolene
Sterilizers.
- AN79 Anprolene 20-Ampoule Refill Dispenser:
use with AN74C/D/E/I/ix Anprolene
Sterilizers.
- AN7916 Anprolene 16-Ampoule Refill Dispenser:
use with AN74C/D/E/I/ix Anprolene
Sterilizers.
- AN72C & D Ventilated and Purged Anprolene Tray
Sterilizer with automatic cycle controls.
- AN74C/D/E Ventilated and Purged, High-Capacity
Anprolene Sterilizer with automatic cycle

- controls.
- AN74i/ix Ventilated and Purged, High-Capacity
Anprolene Sterilizer with automatic cycle
controls.
- AN74J Ventilated and Purged, High-Capacity
Anprolene Sterilizer with automatic cycle
controls.
- AN2000 Ventilated and Purged, High-Capacity
Anprolene Sterilizer with automatic cycle
controls.

Anprolene Sterilizing Accessories

- AN80 Steritest biological and chemical controls.
- AN85 Anprolene Exposure Indicator Strips.
- AN86 Anprolene Exposure Indicator Tape.
- AN87 Anprolene Dosimeter.
- AN90 Seal and Peel Electric Impulse Sealer.
- AN820 2" x 200' (5 cm x 60 m)
Seal and Peel Roll Stock.
- AN830 3" x 200' (8 cm x 60 m)
Seal and Peel Roll Stock.
- AN850 5" x 200' (13 cm x 60 m)
Seal and Peel Roll Stock.
- AN870 7" x 200' (18 cm x 60 m)
Seal and Peel Roll Stock.
- AN1071 Humidichip (2" x 2" / 5 cm x 5 cm) humidity
stabilizing chips.
- AN1072 Humiditube (5 in a pack).

Anpro®, Anprolene®, Dosimeter®, Humidichip®, Seal
and Peel® and Steritest® are trademarks of Andersen
Products, Inc.

Humidichip U.S. patent No. 5,082,636

Manufactured by:
Andersen Sterilizers, Inc.
Health Science Park
3154 Caroline Drive
Haw River, NC 27258 U.S.A.
Tel: 336-376-8622
Fax: 336-376-3088

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Validation services and contract
sterilization available from:
Andersen Scientific, Inc.
3220 Caroline Drive
Haw River, NC 27258 U.S.A.
Tel: 888-933-6427
Fax: 336-376-8153

Distributed in Europe, the Middle East and Africa by:
H. W. Andersen Products, Ltd.
Davy Road • Clacton-on-Sea, Essex
Essex C015 4XA U.K.
Tel: 44-1-255-428-328
Fax: 44-1-255-222-987

Distributed in France by:
H.W. Andersen, SARL
Batiment A - RDC - Parc d'Activités de l'Epi de Soll
55 rue Salvador Allende
59120 LOOS, France
Tel: 33-03-20-97-05-25
Fax: 33-03-20-97-61-77

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Anprolene Modular Sterilizing System

AN71	Anprolene 25 Ampoule Refill Dispenser: use with AN72C & D Anprolene Steriliziers.
AN73	Anprolene 60-Ampoule Refill Dispenser: use with AN72 C & D Anprolene Steriliziers.
AN79	Anprolene 20-Ampoule Refill Dispenser: use with AN74C/D/E/i/ix Anprolene Steriliziers.
AN7916	Anprolene 16-Ampoule Refill Dispenser: use with AN74C/D/E/i/ix Anprolene Steriliziers.
AN72C & D	Ventilated and Purged Anprolene Tray Sterilizer with automatic cycle controls.
AN74C/D/E	Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle

AN74i/ix	controls. Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.
AN74J	Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.
AN2000	Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.

Anprolene Sterilizing Accessories

AN80 Steritest biological and chemical controls.
AN85 Anprolene Exposure Indicator Strips.
AN86 Anprolene Exposure Indicator Tape.
AN87 Anprolene Dosimeter.
AN90 Seal and Peel Electric Impulse Sealer.
AN820 2" x 200' (5 cm x 60 m) Seal and Peel Roll Stock.
AN830 3" x 200' (8 cm x 60 m) Seal and Peel Roll Stock.
AN850 5" x 200' (13 cm x 60 m) Seal and Peel Roll Stock.
AN870 7" x 200' (18 cm x 60 m) Seal and Peel Roll Stock.
AN1071 Humidichip (2" x 2" / 5 cm x 5 cm) humidity stabilizing chips.
AN1072 Humiditube (5 in a pack).

Anpro®, Anprolene®, Dosimeter®, Humidichip®, Seal and Peel® and Steritest® are trademarks of Andersen Products, Inc.

Humidichip U.S. patent No. 5,082,636

Manufactured by:
Andersen Steriliziers, Inc.
Health Science Park
3154 Caroline Drive
Haw River, NC 27258 U.S.A.
Tel: 336-376-8622
Fax: 336-376-3088

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Validation services and contract sterilization available from:
Andersen Scientific, Inc.
3220 Caroline Drive
Haw River, NC 27258 U.S.A.
Tel: 888-933-6427
Fax: 336-376-8153

Distributed in Europe, the Middle East and Africa by:
H. W. Andersen Products, Ltd.
Davy Road • Clacton-on-Sea, Essex
Essex CO15 4XA U.K.
Tel: 44-1-255-428-328
Fax: 44-1-255-222-987

Distributed in France by:
H.W. Andersen, SARL
Batiment A - RDC - Parc d'Activites de l'Epi de Soll
55 rue Salvador Allende
59120 LOOS, France
Tel: 33-03-20-97-05-25
Fax: 33-03-20-97-61-77

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EPA Registration No.69340-2
EPA Establishment No. 69340-NC-001