**AN
7916**

ONE BOX CONTENTS: 16 AMPOULES, 17 LINER BAGS, 16 CLOSURES

ANPROLENE®

STERILIZING GAS AMPOULES

Active ingredient: Ethylene oxide.....97.0%
Inert ingredients:.....3.0%
Total:.....100.0%
Net contents 16 ampoules, 0.64 av. oz. (18.15g) each ampoule

Manufactured and Distributed by:
Andersen Sterilizers, Inc.
3154 Caroline Drive
Haw River, NC 27258 USA
800-376-8622

European Distributor:
H.W. Andersen Products, Ltd.
Davy Road • Clacton-on-Sea
Essex C015 4XA UK

Distributed by:
Andersen Products, Inc.
3202 Caroline Drive
Haw River, NC 27258 USA
800-523-1276

EPA Registration No: 69340-2
EPA Establishment No: 69340-NC-001

KEEP OUT OF REACH OF CHILDREN • DANGER

PRECAUTIONARY STATEMENTS

Hazard to Humans and Domestic Animals.

DANGER

Causes irreversible eye damage and skin burns. Harmful if inhaled. Do not breath vapor. Do not get on eyes, skin or clothing. Do not swallow. Cancer Hazard and reproductive Hazard. May cause nervous system damage. Store and use with adequate ventilation in accordance with 29 CFR 1910.1047

See user's manual for additional precautionary statements.
Emergency Contact Number: 1-800-255-3924

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IF IN EYES: Hold eyes open and rinse slowly and gently with water for 15 - 20 minutes. Remove contact lenses, if present, after the first 5 minutes then continue rinsing eyes. Get medical treatment.

IF INHALED: Get exposed person to fresh air. Keep warm. If not breathing call 911 or an ambulance, then give artificial respiration, preferably by mouth - to - mouth, if possible. If breathing is difficult, give oxygen. Call a physician even if no symptoms are present. Keep under medical observation. Symptoms may be delayed.

IF ON SKIN: Immediately wash skin with plenty of water, while removing contaminated clothing and shoes. Call a physician. Aerate, wash, or clean contaminated clothing. Discard leather goods.

IF SWALLOWED: Drink at least two glasses of water. Do not induce vomiting. Do not give anything by mouth to an unconscious victim. Call a physician.

NOTE TO PHYSICIAN

Ethylene oxide is a gas. Skin exposure from contact with fabric, rubber or plastic containing residual ethylene oxide will commonly result in skin burns with extensive blister formation. At high concentrations severe conjunctivitis can occur. Respiratory tract irritation may occur, but without acute lung edema. Symptoms of systemic intoxication are headache, nausea, vomiting, incoordination, and cardiac irregularities. Treatment is symptomatic.

ENVIRONMENTAL HAZARDS

Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans, or other waters unless in accordance with the requirements of a National Discharge Eliminations System (NDPES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office of the EPA.

PHYSICAL OR CHEMICAL HAZARDS DANGER FLAMMABLE LIQUID AND GAS UNDER PRESSURE.

Ethylene oxide gas is extremely flammable. Do not use near flame, electrical sparks, or hot surfaces. Ground all equipment to prevent static sparks.

IMPORTANT: For Professional Use Only

Users must follow the requirements of the OSHA occupational exposure standard or ethylene oxide (29 CFR 1910.1047).

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It is a violation of Federal law to use this product in a manner inconsistent with its labeling. It is the responsibility of the employer of any person engaged in the handling or application of this product to follow the requirements of 29 CFR 1910.1047. Read package insert for complete sterilization instructions and additional precautions.

FOR USE IN ANPROLENE STERILIZERS ONLY

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STORAGE: Store in a cool place and out of direct sunlight.

PESTICIDE DISPOSAL: Unwanted or expired ampoules should be returned to the manufacturer for disposal. Contact Andersen Products, Inc. (tel. no. 336-376-3000) for instructions. If unwanted or expired ampoules cannot be disposed of according to Control Agency, contact the Hazardous Waste Representative at the nearest EPA Regional Office for guidance.

CONTAINER DISPOSAL: Do not reuse empty box, empty ampoules, or liner bags. Wrap box, ampoules and bags in paper and discard in the trash.

PN736_VM061307

69340-2

9/18/2007

Page 1918

Do not open this bag before the 14 hour cycle at or above 68°F / 20°C is complete.

EPA Reg. No. 69340-2

Shortening the mandatory 14 hour cycle may lead to unacceptable operator exposure to ethylene oxide and chemical burns to the patient. Aerate gas absorbing materials outside of the Sterilization Liner Bag according to the enclosed instructions.

The Anprolene sterilizer must be operated according to the detailed instructions available in the Owner's Manual and the Instructions enclosed in each Anprolene Refill Kit.

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ANPROLENE®

ENGLISH:

BRAND GASEOUS STERILANT FOR
ROOM TEMPERATURE,
ATMOSPHERIC PRESSURE STERILIZATION

Instructions for using Anprolene gas ampoules

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ACCEPTED

SEP 18 2007

Under the Federal Insecticide, Fungicide, and
Rodenticide Act as amended, is the
pesticide, registered under
EPA Reg. No. 69340-2

Manufactured by:

Andersen Sterilizers, Inc.
Health Science park
3154 Caroline Drive
Haw River, NC 27258 U.S.A.

CE 0413

Distributed by:

Andersen Products, Inc.
Health Science Park
3202 Caroline Drive
Haw River, NC 27258 U.S.A.

H.W. Andersen Products, Ltd.

Davy Road • Gorse Lane Industries Estates
Clacton-on-Sea • Essex C015 4XA U.K.

Distributed in Asia by:

Andersen Medical (HK), Ltd.
Unit 2, 18/F Westley Square
40 Hoi Yuen Road
Kwun Tong, Hong Kong
Tel:2851-1655

PN863_VN061307

3 8 18

AN7916

ANPROLENE®

Brand Gaseous Sterilant for Room Temperature,
Atmospheric Pressure Sterilization
Use in Ventilated Anprolene Sterilizers Only

DIRECTIONS FOR USE

To be used only in Anprolene sterilizers, and only for hospital, medical, and veterinary sterilization.

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. It is the responsibility of the employer of any person engaged in the handling or application of this product to follow the requirements of 29 CFR 1910.1047.

Do not open the plastic bag in which each individual Anprolene ampoule is sealed.

Each ampoule of Anprolene is surrounded by a plastic break shield. The ampoule and shield are sealed in a plastic bag. The plastic bag is a gas diffusion membrane of known permeability whose function is to contain the gas given off by the ampoule and to release it at a controlled rate during the sterilization cycle. The plastic shield around the ampoule prevents the broken glass of the opened ampoule from puncturing the gas-release bag.

Typical products which may be conveniently processed in Anprolene sterilizers

Respirators, corrugated tubing
Bronchoscopes, gastroscopes, fiberscopes of all kinds
Procedure trays
Catheters - plastic, rubber, metal, glass, cloth
Anesthesia equipment - endotracheal tubes, masks
Rubber tubing
Adhesive tape
Bandages, dressing sets (reuse plastic forceps)
Syringes - plastic, rubber, glass, bulb syringes
Gloves - rubber, plastic, cloth
Surgical instruments - steel, chrome plate, brass, plastic
Optical instruments - scopes, cameras, lenses, mirrors
Electrical equipment - whether autoclavable or not
Painted equipment - metal, wood
High-speed steel - drills, burrs, chisels
Airways - plastic, rubber, metal
Fabric - cloth, rubber, plastic
Electric wire - whether autoclavable or not
Dry-cell batteries, battery cases, bulbs

Sutures - plastic, silk, cotton, stainless steel
Thermometers, applicator sticks
Rectal tubes, douche tubes - rubber plastic
Specula - plastic, metal

Preparation of materials for sterilization

Materials to be sterilized by Anprolene must first be meticulously cleaned and towel-dried. Coating of dry protein, such as dry pus, blood, or feces, protect microorganisms and slow the sterilization process. To prevent this possibility, and to be certain that your system meets the highest standard of sterilization, you must:

1. Disassemble and scrub all instruments in detergent and water to the most critical standard of cleanliness possible. If an item cannot be washed in detergent and water, then pre-sterilization humidification is required. This pretreatment must be done in a chamber having 100% relative humidity at a temperature of at least 68° F (20°C) for at least four hours. Such a chamber may be constructed by placing an AN1071 Humidichip® in a liner bag with the item and closing the bag with a twist-tie.
2. Water on instruments at the time of exposure to Anprolene may react with the gas and reduce its effectiveness. Be sure that items to be sterilized are dry before wrapping and processing. Drain or towel-dry instruments. Do not dry them in a hot-air oven.
3. Occlusive caps, plugs, or stylets must be removed from instruments so that the gas can penetrate freely. Hollow-bore needles and plastic or rubber tubing must be open at both ends and free from plugs. Syringes must be packaged disassembled, with the plunger removed from the barrel.
4. You must wrap all items individually, in cloth or paper, in a manner conventional for steam sterilization, or in Anpro® Seal and Peel™ Packaging. Anpro Seal and Peel is a transparent, peel-open, extended shelf-life packaging, proven to be compatible with the Anprolene Sterilizing System. Do not pack the liner bag so tightly with cloth or gauze that gas diffusion is slowed.

Sterilization method using AN74A, AN74C/D/E, AN74V AND AN74ix

Be certain that all items have been prepared as described above. The ambient relative humidity must be at least 30% for Anprolene processing. The user must verify that this minimum exists before processing begins. Remove one liner bag from the AN7916 Anprolene dispenser box. Prepare the bag to receive the wrapped items to be

4 3 1 8

sterilized by opening it and placing it into the sterilizer container, open end out. Put the wrapped items to be sterilized into the liner bag.

Remove one ampoule from the AN7916 Anprolene dispenser box. Do not open the plastic bag in which it is wrapped. Rather, unroll the plastic bag. Push the ampoule gently to the center of the gas-release bag. Place the ampoule, still sealed in its plastic bag, inside the liner bag along with the material to be sterilized. Press out any excess air before closing the mouth of the liner bag. Locate the flexible, plastic purge tube that protrudes into the interior of the sterilizer. Place the purge tube into the mouth of the liner bag. Using one of the closures provided in the Anprolene kit, tie the liner bag around the purge tube. Grasp the ampoule through the bag and protective shield and snap the top of the ampoule. Each ampoule is pre-scored around its neck to facilitate this action. Snapping off the top of the ampoule activates it, i.e., releases the sterilant gas within the gas-release bag. Close and lock the sterilizer. The AN74 key must be removed from the sterilizer and remain in the possession of a supervisor to prevent unauthorized removal of goods during the sterilization cycle.

AN74C/D/E/i/x: Close sterilizer and lock the door. The AN74 key must be removed from the sterilizer and remain in the possession of a supervisor to prevent unauthorized removal of goods during the sterilization cycle. The sterilizer must be located in a clean, well-ventilated area, away from sparks or flame. The temperature of the room must be maintained at no less than 68°F (20°C).

At the end of the fourteen-hour cycle, unlock the sterilizer. Remove the seal from the liner bag and open it. Leave the opened liner bag and the wrapped items in the sterilizer for an additional fifteen minutes before removing them.

AN74A or AN74V: Turn the ventilation on if you are using an AN74V sterilizer, make sure that the purge pump is turned off. If you are using an AN74A sterilizer, press the CYCLE START button.

At the end of the twelve-hour sterilization cycle, the AN74A will automatically begin a two-hour purge cycle and a green light will signal when you may open the sterilizer to remove the load. If you are using an AN74V, turn the purge pump on and purge the liner bag for two hours before opening the sterilizer and removing the sterile items.

CAUTION: FAILURE TO OBSERVE THE ABOVE INSTRUCTIONS MAY EXPOSE THE STERILIZER OPERATOR TO MORE THAN 0.5 PPM ETHYLENE

OXIDE (The OSHA Established Action Level Measured as an Eight-Hour Time Weighted Average).

Aerate the gas-absorbent items for at least 24 hours before use (see: Precautionary Statements).

After the sterile materials have been removed from the liner bag, the exhausted and empty gas-release ampoule and bag along with the liner bag may be disposed of in ordinary rubbish.

Special processes

When sterilizing long lengths of plastic or rubber tubing, such as urethral catheters, heart catheters, and coils of drainage tubing, use two ampoules and double the cycle time.

If you are using an AN74 series sterilizer, increase the cycle time to 24 hours. If you are using an AN74V sterilizer, increase the cycle time to 24 hours, and purge the liner bag for two hours. With an AN74C/D/E, load the sterilizer with the items for sterilization and the two activated gas release bags, but then wait twelve hours before pushing the START button. Before loading the sterilizer, be sure that the ventilator is running. The ventilator must continue to run during the entire augmented cycle. For a 74i/ix, simply select the 24 hour cycle instead of the 12 hour cycle.

Technical description of the system

The Anprolene Sterilizer container serves as a guard against inadvertent ignition of the contents by spark or open flame during the sterilization cycle. In the case of the AN74C, AN74D, AN74E, and AN74i/ix sterilizers, the container also serves as a ventilation hood.

The liner bag, when properly closed with a tie or Velcro® strap, has a 35 liter capacity. It is a second gas-diffusion membrane that serves to retain Anprolene long enough to sterilize the liner bag contents and releases it into the surrounding vent hood at a slow enough rate to ensure that toxic levels are not reached in a properly ventilated room. Each ampoule releases approximately 17.6 grams of Anprolene at room temperature (68°F/20°C) and sea-level atmospheric pressure.

The AN7916 Anprolene ampoule produces a minimum peak concentration within the liner bag of 500 mg / 1000 cc. Tests in our laboratory confirm that this dose will kill the most resistant spores known within the twelve-hour cycle, providing they have been re-hydrated according to our instructions (see: Preparation of Materials for Sterilization).

Storage and disposal

Do not contaminate water, food, or feed by storage or disposal.

STORAGE: Anprolene must be stored in a cool place out of direct sunlight. Under normal conditions its shelf life will exceed one year. As long as the content of the ampoule is liquid at 68°F (20°C), it is sufficiently potent to use.

PESTICIDE DISPOSAL: Unwanted or expired ampoules should be returned to the manufacturer for disposal. Contact Andersen Products, Inc. (tel. no. 336-376-3000) for instructions. If unwanted or expired ampoules cannot be disposed of according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste Representative at the nearest EPA Regional Office for guidance.

CONTAINER DISPOSAL: Do not reuse empty box, used ampoules, or liner bags. Wrap used box, ampoules, and liner bags in paper and discard in the trash.

Testing the efficiency of the Anprolene Sterilizer

Monitoring sterilization efficacy is extremely important. It is strongly recommended that the user establish a routine for monitoring each cycle. A color-change chemical indicator, such as the AN87 Dosimeter®, placed in the most inaccessible part of the load, will indicate whether or not the gaseous sterilant penetrated to the core of the load in adequate concentration to assure sterilization. In addition, an appropriate biological control, such as the AN80 Steritest®, should be used at least once per month to challenge the procedure. The AN80 Steritest is sensitive to the gas concentration, cycle time, and average cycle temperature.

AERATION: It is incumbent on the user to confirm that gas-absorbent items, sterilized in AN7916 Anprolene, have been adequately aired before they contact living tissue. This is a particularly important warning to those using plastic or rubber items that directly contact tissue culture preparations, ova, semen or embryos.

Nylon and polyester films are virtually impervious to Ethylene Oxide. The only other plastic, waterproof wrapping material proven compatible with AN7916 Anprolene is our Seal and Peel® brand roll stock. Do not use any other plastic film to wrap items to be sterilized in AN7916 Anprolene.

Vacuum dehydration, chemical desiccation or prolonged exposure to ambient relative humidity below 30% has

been demonstrated to produce spores highly resistant to sterilization by AN7916 Anprolene. Re-hydration of spores so changed, and hence reversion to normal sensitivity, does not seem to occur until they have been actually wetted or placed in a 100% relative humidity atmosphere. Do not attempt to sterilize materials which may be carrying dried spores without first washing the articles with water and detergents. If the nature of the material is such that the water treatments specified above are harmful, then pretreatment must be done in a chamber having 100% relative humidity at a temperature of at least 68°F (20°C) for at least four hours. Such a chamber may be made by placing an AN1071 Humidichip® in a liner bag with the item and closing the bag with a twist tie.

PHYSICAL OR CHEMICAL HAZARDS

DANGER: FLAMMABLE LIQUID AND GAS UNDER PRESSURE

Ethylene Oxide gas is extremely flammable. Do not use near flame, electrical sparks, or hot surfaces. Ground all equipment to prevent static sparks.

IMPORTANT - THE USER MUST NOT DEVIATE FROM THESE INSTRUCTIONS

Do not be fooled by the apparent simplicity of the Anprolene Sterilizing System. Its reliable use depends upon precise adherence to these instructions. AN71 and AN73 Anprolene ampoules are designed to be used only in the AN72C, AN72D, and AN72V Anprolene Sterilizers with the liner bags and twist seals supplied in the refill dispenser. No other container may be used, no matter how similar it seems. Liner bags must not be reused; use only a genuine, fresh Anprolene liner bag with each load.

Precautionary Statements
Hazard to Humans and Domestic Animals

DANGER

CAUSES IRREVERSIBLE EYE DAMAGE AND SKIN BURNS. HARMFUL IF INHALED. DO NOT BREATHE VAPOR. DO NOT GET ON EYES, SKIN, OR CLOTHING. DO NOT SWALLOW. CANCER HAZARD AND REPRODUCTIVE HAZARD. MAY CAUSE NERVOUS SYSTEM DAMAGE. STORE AND USE WITH ADEQUATE VENTILATION IN ACCORDANCE WITH 29 CFR 1910.1047

See user's manual for additional precautionary statements.

5
8
18

Emergency Contact Number: 1-800-255-3924

Never open the sterilizer during the sterilization cycle. Avoid breathing AN7916 Anprolene vapor.

Exposure to liquid or gas can cause severe eye damage. Intimate skin contact with fabric, clothing, rubber, or plastics containing residual Ethylene oxide can cause damage to the skin. Wear goggles or face shield when handling. May be fatal if swallowed. Do not breathe vapor. Extremely hazardous liquid and gas under pressure. Inhalation may be fatal. Ethylene oxide gas can be odorless at toxic levels; therefore, exposure to toxic levels may occur without warning or capability of detection by user.

EFFECTS OF OVEREXPOSURE TO ETHYLENE OXIDE

May be fatal if inhaled in high concentrations. May cause irritation of respiratory tract, chest tightness, headache, nausea, vomiting, diarrhea, light-headed feeling, dizziness, weakness, drowsiness, cyanosis, loss of coordination, convulsions, coma, delayed lung injury (fluid in lungs), immediate or delayed skin irritation and blisters, allergic skin.

OTHER POSSIBLE DELAYED HEALTH EFFECTS

May cause nervous system injury, cataracts, adverse reproductive effects, chromosomal and mutagenic changes, and cancer.

PEL: 1PPM-TWA Ethylene Oxide (OSHA - 29 CFR 1910.1047)

STEL: 5PPM-excursion limit, 15 minutes

ODOR: Ether-like in high concentrations. Exposure to toxic levels may occur without warning or detection by the user.

See label for additional precautionary statements and statements of practical treatment.

Aeration: Failure to adequately air gas-absorbing materials may lead to contact chemical burns. All items that may contact living tissues must be aired for at least 24 hours, at a minimum temperature of 68°F (20°C), prior to use. Items such as plastic instruments, foam rubber, plastic foams, vinyl tubing, rubber tubing, plastic items, rubber items, and sealed air cushioned devices (for instance, anesthesia masks), must be aired. Instruments need not be removed from cloth or paper wrappings to obtain adequate diffusion of the residual Anprolene vapors.

Operator Exposure: It is absolutely essential that reliable

forced ventilation be employed in the room where the Anprolene Sterilizer is used. This system must be capable of ten air changes per hour, so that operator exposure remains within U.S. government permissible exposure limits.

FIRST AID STATEMENTS

In all cases of exposure, get medical attention immediately. Take person to a doctor or emergency treatment facility. Have the product container or label with you when going for treatment.

If in Eyes: Hold eyes open and rinse slowly and gently with water for 15 - 20 minutes. Remove contact lenses, if present, after the first 5 minutes then continue rinsing eye. Get medical treatment.

If Inhaled: Get exposed person to fresh air. Keep warm. *If not breathing call 911 or an ambulance, then give artificial respiration, preferably by mouth - to - mouth, if possible.* If breathing is difficult, give oxygen. Call a physician even if no symptoms are present. Keep under medical observation. Symptoms may be delayed.

If on Skin: Immediately wash skin with plenty of water, while removing contaminated clothing and shoes. Call a physician. Aerate, wash, or clean contaminated clothing and discard leather goods.

If Swallowed Drink at least two glasses of water. Do not induce vomiting. Do not give anything by mouth to an unconscious victim. Call a physician.

See label or instructions for additional precautionary statements and statements of practical treatment.

NOTE TO PHYSICIAN

Ethylene oxide is a gas. Skin exposure from contact with fabric, rubber or plastic containing residual Ethylene Oxide will commonly result in skin burns with extensive blister formation. At high concentrations severe conjunctivitis can occur. Respiratory tract irritation may occur, but without acute lung edema. Symptoms of systemic intoxication are headache, nausea, vomiting, incoordination, and cardiac irregularities. Treatment is symptomatic.

Ampoule Fill Weights

Active ingredient: Ethylene Oxide	97.0%
Other Ingredients:	3.0%
Total:	100.0%
Net contents each ampoule	0.64 av. oz (18.15 g)

Anprolene Modular Sterilizing System

- AN71 Anprolene 25 Ampoule Refill Dispenser: use with AN72C & D Anprolene Sterilizers.
- AN73 Anprolene 60-Ampoule Refill Dispenser: use with AN72 C & D Anprolene Sterilizers.
- AN79 Anprolene 20-Ampoule Refill Dispenser: use with AN74C/D/E/i/ix Anprolene Sterilizers.
- AN7916 Anprolene 16-Ampoule Refill Dispenser: use with AN74C/D/E/i/ix Anprolene Sterilizers.
- AN72C & D Ventilated and Purged Anprolene Tray Sterilizer with automatic cycle controls.
- AN74C/D/E Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.
- AN74i/ix Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.
- AN74J Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.
- AN2000 Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.

Anprolene Sterilizing Accessories

- AN80 Steritest biological and chemical controls.
- AN85 Anprolene Exposure Indicator Strips.
- AN86 Anprolene Exposure Indicator Tape.
- AN87 Anprolene Dosimeter.
- AN90 Seal and Peel Electric Impulse Sealer.
- AN820 2" x 200' (5 cm x 60 m) Seal and Peel Roll Stock.
- AN830 3" x 200' (8 cm x 60 m) Seal and Peel Roll Stock.
- AN850 5" x 200' (13 cm x 60 m) Seal and Peel Roll Stock.
- AN870 7" x 200' (18 cm x 60 m) Seal and Peel Roll Stock.
- AN1071 Humidichip (2" x 2" / 5 cm x 5 cm) humidity stabilizing chips.
- AN1072 Humiditube (5 in a pack).

Anpro®, Anprolene®, Dosimeter®, Humidichip®, Seal and Peel® and Steritest® are trademarks of Andersen Products, Inc.

Humidichip U.S. patent No. 5,082,636

Manufactured by:
Andersen Sterilizers, Inc.
Health Science Park
3154 Caroline Drive
Haw River, NC 27258 U.S.A.
Tel: 336-376-8622
Fax: 336-376-3088

Distributed by:
Andersen Products, Inc.
Health Science Park
3202 Caroline Drive
Haw River, NC 27258 U.S.A.
1-800-523-1276
www.anpro.com
mailbox@anpro.com

Distributed in Europe, the Middle East and Africa by:
H. W. Andersen Products, Ltd.
Davy Road • Clacton-on-Sea, Essex
Essex CO15 4XA U.K.
Tel: 44-1-255-428-328
Fax: 44-1-255-222-987

Distributed in Asia by:
Andersen Medical (HK) Ltd.
Unit 2, 18/F Westley Square
40 Hoi Yuen Road
Kwun Tong, Hong Kong

EPA Registration No.69340-2
EPA Establishment No. 69340-NC-001

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Ampoule Fill Weights

Active ingredient: Ethylene oxide97.0%
Inert Ingredients: 3.0%
Total:100.0%
Net contents each ampoule0.64 av. oz (18.15 g)

Anprolene Modular Sterilizing System

AN71	Anprolene 25 Ampoule Refill Dispenser: use with AN72C & D Anprolene Sterilizers.
AN73	Anprolene 60-Ampoule Refill Dispenser: use with AN72 C & D Anprolene Sterilizers..
AN79	Anprolene 20-Ampoule Refill Dispenser: use with AN74C/D/E/i/ix Anprolene Sterilizers.
AN7916	Anprolene 16-Ampoule Refill Dispenser: use with AN74C/D/E/i/ix Anprolene Sterilizers.
AN72C & D	Ventilated and Purged Anprolene Tray Sterilizer with automatic cycle controls.
AN74C/D/E	Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.
AN74i/ix	Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.
AN74J	Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.
AN2000	Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.

Anprolene Sterilizing Accessories

AN80 Steritest biological and chemical controls.
AN85 Anprolene Exposure Indicator Strips.
AN86 Anprolene Exposure Indicator Tape.
AN87 Anprolene Dosimeter.
AN90 Seal and Peel Electric Impulse Sealer.
AN820 2" x 200' (5 cm x 60 m)
Seal and Peel Roll Stock.
AN830 3" x 200' (8 cm x 60 m)
Seal and Peel Roll Stock.
AN850 5" x 200' (13 cm x 60 m)
Seal and Peel Roll Stock.
AN870 7" x 200' (18 cm x 60 m)
Seal and Peel Roll Stock.

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Distributed in Europe, the Middle East and Africa by:
H. W. Andersen Products, Ltd.
Davy Road • Clacton-on-Sea, Essex
Essex C015 4XA U.K.
Tel: 44-1-255-428-328
Fax: 44-1-255-222-987

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EPA Registration No. 69340-2
EPA Establishment No. 69340-NC-001

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Anprolene Modular Sterilizing System

- AN71 Anprolene 25 Ampoule Refill Dispenser: use with AN72C & D Anprolene Sterilizers.
- AN73 Anprolene 60-Ampoule Refill Dispenser: use with AN72 C & D Anprolene Sterilizers.
- AN79 Anprolene 20-Ampoule Refill Dispenser: use with AN74C/D/E/i/ix Anprolene Sterilizers.
- AN7916 Anprolene 16-Ampoule Refill Dispenser: use with AN74C/D/E/i/ix Anprolene Sterilizers.
- AN72C & D Ventilated and Purged Anprolene Tray Sterilizer with automatic cycle controls.
- AN74C/D/E Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.

- AN74i/ix Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.
- AN74J Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.
- AN2000 Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.

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Humidichip U.S. patent No. 5,082,636

Manufactured by:
Andersen Sterilizers, Inc.
 Health Science Park
 3154 Caroline Drive
 Haw River, NC 27258 U.S.A.
 Tel: 336-376-8622
 Fax: 336-376-3088

Distributed by:
Andersen Products, Inc.
 Health Science Park
 3202 Caroline Drive
 Haw River, NC 27258 U.S.A.
 1-800-523-1276
 www.anpro.com
 mailbox@anpro.com

Distributed in Europe, the Middle East and Africa by:
H. W. Andersen Products, Ltd.
 Davy Road • Clacton-on-Sea, Essex
 Essex CO15 4XA U.K.
 Tel: 44-1-255-428-328
 Fax: 44-1-255-222-987

Distributed in Asia by:
Andersen Medical (HK) Ltd.
 Unit 2, 18/F Westley Square
 40 Hoi Yuen Road
 Kwun Tong, Hong Kong

EPA Registration No.69340-2
 EPA Establishment No. 69340-NC-001

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Anprolene Modular Sterilizing System

- AN71 Anprolene 25 Ampoule Refill Dispenser: use with AN72C & D Anprolene Steriliziers.
- AN73 Anprolene 60-Ampoule Refill Dispenser: use with AN72 C & D Anprolene Steriliziers.
- AN79 Anprolene 20-Ampoule Refill Dispenser: use with AN74C/D/E/i/i/x Anprolene Steriliziers.
- AN7916 Anprolene 16-Ampoule Refill Dispenser: use with AN74C/D/E/i/i/x Anprolene Steriliziers.
- AN72C & D Ventilated and Purged Anprolene Tray Sterilizer with automatic cycle controls.
- AN74C/D/E Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.

- AN74i/i/x Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.
- AN74J Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.
- AN2000 Ventilated and Purged, High-Capacity Anprolene Sterilizer with automatic cycle controls.

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Humidichip U.S. patent No. 5,082,636

Manufactured by:
Andersen Steriliziers, Inc.
Health Science Park
3154 Caroline Drive
Haw River, NC 27258 U.S.A.
Tel: 336-376-8622
Fax: 336-376-3088

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