UN .. D STATES ENVIRONMENTAL PROTEC ... IN AGENCY

69687-1

12/29/2006

Page 18 7

December 29, 2006

Bob MacDonald Agent for Medtrol, Inc. c/o SRC Scientific & Regulatory Consultants, Inc. P.O. Box 1014 Columbia City, IN 46725

Subject:

Super-Chlor

EPA Registration No. 69687-1

Application Date: December 6, 2006 Receipt Date: December 14, 2006

Dear Mr. MacDonald:

This acknowledges receipt of your notification, submitted under the provision of PR Notice 98-10, FIFRA section 3(c)9.

Proposed Notification

Addition of "Total" line to ingredient statement

General Comments

Based on a review of the material submitted, the following comments apply:

The notification application is acceptable and a copy has been placed in your file for future reference.

Should you have any questions or comments concerning this letter, please contact me at (703) 308-6345.

Sincerely,

Wanda Y. Henson Product Reviewer (32) Regulatory Management Branch II Antimicrobials Division (7510P)

No. 1	CONCURRENCES	
TS10 P 7510 P		
TRHAME E. Berg & Mitchel	**************************************	
IRNAME E. Brig & Shiklu LTE 12/29/06 1/83/07	***************************************	***************************************
PA Form 1320-1A (1/90)	Printed on Recycled Paper	DEFICIAL EN T CONT

2 7 7

Please read instructions	on reverse before	oleting form.		Form App	ad. OMB No.	2070-006	O. Approval expires 2-28-9	
\$EPA	United States Environmental Protection Agency Washington, DC 20460		✓	Registr Amend Other		OPP Identifier Number		
		Applicati	on for Pesticide	- Section	n I			
1. Company/Product Number 69687-1		į.	2. EPA Product Manager 3 Emily Mitchell			oposed Classification		
4. Company/Product (Name) Super-Chlor			PM# 32					
5. Name and Address of Applicant (Include ZIP Code) Medtrol, Inc. 7157 Austin Avenue Niles, IL 60714		6. Expedited Reveiw. In accordance with FIFRA Section 3(c)(3) (b)(i), my product is similar or identical in composition and labeling to: EPA Reg. No.						
Check if	this is a new address		Product	Product Name				
			Section - II	_				
Amendment - Exp	response to Agency lett	er dated		nal printed lab jency latter d le Too" Appli her - Explain	ication.	se to		
make any false statemer product may be in violation	it to EPA. I further unders on of FIFRA and I may be	tand that if this r subject to enfor	otification is not consiste cement action and penal	nt with the ter ties under sec	ms of PR Notice tions 12 and 1	e 98-10 an 4 of FIFRA.	J.S.C. Sec, 1001 to willfully d 40 CFR 152.46, this hald@srcconsultants.com or	
1. Material This Product	Will Be Packaged In:				<u> </u>			
Child-Resistant Packagin Yes No Certification must be submitted	Unit Packaging Water Yes No If "Yes" No. per If "Yes"			Yes Me Pla No Gla		Metal Plastic Glass Paper	tai stic ss	
3. Location of Net Conte	nts Information	4. Size(s) Re	tail Container	5.1	ocation of La	bel Direction	ons	
Label								
6. Manner in Which Labe	l is Affixed to Product	Lithog Paper Stend	graph [glued iled	Other				
			Section - IV		•			
1. Contact Point (Compl	ete items directly below	for identification	on of individual to be co	ntected, if ne	cessary, to p	rocess this	application.)	
Name Bob MacDonald		Title Agent for Medtrol Inc.		Telephone No. (Include Area Code) 260-244-6270				
•	atements I have made o any knowlinglly falsa o ble law.		all attachments thereto			•	6. Date Application Received (Stamped)	
2. Signature But Man	Sould		3. Title Agent for Medtrol Inc.					
4. Typed Name Bob MacDon			5. Date 12/14/10					



Consultants, Inc.

December 6, 2006

Ms. Emily Mitchell, PM 32
Document Processing Desk (NOTIF)
Office of Pesticide Programs (7504P)
U. S. Environmental Protection Agency
Room S-4900, 4th Floor
One Potomac Yard
2777 S Crystal Drive
Arlington, VA 22202

SUBJECT:

NOTIFICATION

Super-Chlor

EPA Reg. No. 69687-1

Dear Emily,

On behalf of our client, Medtrol, Inc., enclosed is a notification to official update their ingredient statement for the above product to include the standard "Total ...100%" declaration. We are resubmitting this notification per the recent response from Wanda Henson dated 11/28/06 (copy enclosed) requesting a complete label for this notification.

The updated declaration will include the "Total" line and appear as follows:

We have enclosed two copies of the labeling. One with a simply highlight of the change and one with the change "circled" so it may be copied. Otherwise the labeling is consistent with the most recently accepted stamped label.

Thank you for your assistance on this matter. If you have any questions please contact me at (260) 244-6270 or by e-mail at bmacdonald@srcconsultants.com.

Sincerely

Bob MacDonald

Agent for Medtrol, Inc.

P.O. Box 1014 Columbia City, IN 46725

www.srcconsultants.com

Phone: 260-244-6270 Fax: 260-244-6273

SUPER-CHLOR™[Product Name]

[For wiping] [Disinfecting wipe for] hard non-porous [surfaces] [environmental surfaces] [and] [patient care equipment] [point-of-care equipment]

[Tuberculocidal] [Mycobacterium bovis – BCG (OT451C150)] [Mycobacterium bovis variant tuberculosis] [OT 451C150] [(Tuberculosis)] [(TB)] [+Virucidal] [+HIV-I], [Fungicidal] [Trichophyton mentagrophytes] [ATCC 9533], [Bactericidal] [for] [Staphylococcus aureus-MRSA] [ATCC 33592], [Enterococcus faecalis-VRE] [ATCC 51575], [Staphylococcus aureus] [ATCC 6538], [Salmonella choleraesuis] [ATCC 10708], and [Pseudomonas aeruginosa] [ATCC 15442]

Active Ingredient: Sodium hypochlorite	0.525%
Other Ingredients*:	
Total	100.000%
Other ingredients do not include the weight of to	owelette.

SUGGESTED AREAS OF USAGE:

[Medical] [,] [and] [Dental] [and] [Laboratory] [Counters] [,] [Exam Tables] [,] [Carts] [,] [Patient Care Equipment] [Point-of-Care Equipment] [,] [Telephones] [,] [Sink Tops] [,] [Toilet Seats]

KEEP OUT OF REACH OF CHILDREN CAUTION

FIRST AID: If in eyes: Hold eye open and rinse slowly and gently with water for 15 – 20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.

STORAGE AND DISPOSAL

Store product in a cool, dry area, away from heat and direct sunlight. • Do not reuse empty packet and used towelette. • Discard used towelette and empty packet in designated infectious waste container for autoclaving and disposal according to federal, state and local regulations.

(See panel on reverse for directions for use and additional precautionary statements.)

EPA Reg. No. 69687-1

EPA Est. No. 69687-1L-031

[MEDTROL logo] [©2004 Medtrol, Inc.] Net Contents [0.1][0.15][0.2][0.25]

[0.3][0.4][0.5][0.6][0.7][0.8] fl. oz. U.S. Patent No. 6,530,477,B1

[DIN: {insert # as applicable}]

{Administrative Note: Bracketed text [] is optional version of text for different packaging.}

{Administrative Note: The terms "towelette(s)", "towel(s)", "wipe(s)", "wiper(s)" and "cloth(s)" may be used interchangeably

throughout.}

4

 \searrow

Directions for use:

- It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.
- Always use personal protective equipment.
- · Gross filth and heavy soil should be removed prior to disinfection.
- Open [Super-Chlor™][Product name] packet.
- Remove premoistened [3" x 3"][3" x 4"][8" X 10"][9" x 9"][12" x 12"][14" x 14"] towelette.
- Apply towelette and wipe desired surface to be [decontaminated][disinfected].
- A 5 minute contact is required for effectiveness against organisms listed.
- Allow surface to air dry and discard used towelette and empty packet (see front panel).

{Administrative Note: If a virucidal and/or HIV claim is made the following paragraph and special instructions are required}

+KILLS HIV-1 ON PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS: in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of HIV-1(associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.

- PERSONAL PROTECTION: When handling items soiled with blood or body fluids, use disposable latex gloves, gown, masks or eye protection.
- CLEANING PROCEDURES: Blood/body fluids must be thoroughly cleaned from surfaces/objects before application of [Super-Chlor][Product Name].
- CONTACT TIME: Allow surface to remain wet for 5 minutes.
- DISPOSAL OF INFECTIOUS MATERIAL: Blood/body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

Precautionary Statements

Hazards to humans and domestic animals

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barriers or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-vritical medical devices prior to sterilization or high level disinfection.

Complete product line information available by calling:

7157 North Austin Avenue Niles, Illinois 607:4

[MEDTROL, INC. logo]

e-mail: sales@medtrol.com Website: www.medtrol.com

[BARCODE LABEL]

Ph: 800-647-7180 Fx: 800-255-3027

{Administrative Note: Bracketed text [] is optional version of text for different packaging.}

{Administrative Note: The terms "towelette(s)", "towel(s)", "wipe(s)", and "cloth(s)", may be used interchangeably

throughout.}

S

\

{OPTIONAL ADDITIONAL TEXT FOR USE ON PRODUCT LABEL OR COLLATERAL LABELING}

[[Effective against] [Kills] hospital [type] [representative] bacteria Staphylococcus aureus [ATCC 6538], Salmonella choleraesuis [ATCC 10708], and Pseudomonas aeruginosa [ATCC 15442]]

[[Effective against] [Kills] Mycobacterium bovis variant tuberculosis [OT 451C150]]

[[Effective against] [Kills] Trichophyton mentagrophytes [ATCC 9533]]

[[Effective against] [Kills] Methicillin Resistant Staphylococcus aureus (MRSA) [ATCC 33592]]

[[Effective against] [Kills] Enterococcus faecalis (VRE) [ATCC 51575]]

[[Effective against] [Kills] Staphylococcus aureus and Salmonella choleraesuis]

[A One Step Cleaner Disinfectant effective against: [Pseudomonas aeruginosa] [Staphylococcus aureus] [Salmonella choleraesuis] [Methicillin Resistant Staphylococcus aureus (MRSA)] [Vancomycin Resistant Enterococcus faecalis (VRE)] [in the presence of organic soil (5% blood serum)]

[[Effective against] [Kills] Mycobacterium bovis variant tuberculosis [OT 451C150] [in the presence of organic soil (5% blood serum)]

[Where control of the hazards of cross contamination between treated surfaces is of prime importance]

[Antibacterial]

[Aids in reducing [in the reduction of] cross contamination between treated surfaces]

[An effective method to disinfect hard, non-porous environmental surfaces]

[[Use to] [Effectively] [cleans and] disinfects common hard, non-porous surfaces]

[An effective way to clean and deodorize against odor-causing organisms and disinfect hard, non-porous surfaces]

[Cleans [,][deodorizes] and disinfects]

{OPTIONAL ADDITIONAL TEXT FOR USE ON PRODUCT LABEL OR COLLATERAL LABELING}

[[Effective against] [Kills] Mycobacterium bovis – BCG (OT 451C150) [in the presence of organic soil (5% blood serum)]

[Innovative healthcare product designed with a commitment to quality and excellence]

[For complete product line information, Call Medtrol, Inc.]

[Complete product line information available by calling:]

[[Effective against] [Kills] and the *HIV-1 (AIDS Virus) on precleaned hard non-porous environmental surfaces]

[Disinfects hard, non-porous [environmental] [healthcare] surfaces]

[Bathroom [disinfectant] wipes for hard, non-porous environmental surfaces]

[Use wipes to kill [Salmonella choleraesuis (Salmonella)] [and] [Staphylococcus aureus (Staph)] on bathroom hard, non-porous surfaces]

[[Telephone] [[cell] [cellular] [wireless] [mobile] [digital] phone[s]] disinfectant wipe]

[Disinfects [headsets] [telephones] [[cell] [cellular] [wireless] [mobile] [digital] phone[s]]

[Ready to Use]

[Contains 5,250 [parts per million] [ppm] of active sodium hypochlorite]

[[Convenient] [Premeasured] [Premoistened] [Disposable] towelette]

[No mixing] [No measuring]

[Moist inanimate surface wipes]

[Use on hard, non-porous surfaces only.]

[Test wipe on small inconspicuous area.]

[Suitable for hospital use]

[Disinfects hard, non-porous healthcare surfaces]