

FILE COPY

**UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
WASHINGTON, D.C. 20460**



Office of Pesticide Programs

February 20, 2008

Georgia Anastasiou
Lewis & Harrison
122 C Street, NW
Suite 740
Washington, DC 20001

Subject: BASF CATALYSTS LLC
Aseptrol Dental Water Treatment
EPA Registration No. 70060-25
Application Date: December 4, 2007
Receipt Date: December 5, 2007

Dear Ms. Anastasiou:

The following amendment submitted in connection with registration under the Federal Insecticide, Fungicide and Rodenticide Act (FIFRA), as amended, to revise dental delivery unit language to be more generic, is acceptable.

General Comments

A stamped copy of the accepted labeling is enclosed. Submit one copy of your final printed labeling before distributing or selling the product bearing the revised labeling.

Should you have any questions or comments concerning this letter, please call Wanda Henson at (703) 308-6345.

Sincerely,

Emily H. Mitchell
Product Manager - Team 32
Regulatory Management Branch II
Antimicrobials Division (7510P)

- All text in brackets [xxx] is optional and may or may not be included on a final label.
- All text in braces {xxx} is administrative and will not appear on a final label.

PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals

DANGER. Corrosive. Causes irreversible eye damage. Harmful if swallowed. Do not get in eyes or on clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, or using tobacco. Remove contaminated clothing before reuse.

FIRST AID

If in Eyes:	<ul style="list-style-type: none"> •Hold eye open and rinse slowly and gently with water for 15-20 minutes. •Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. •Call a poison control center or doctor for treatment advice.
If Swallowed:	<ul style="list-style-type: none"> •Call poison control center or doctor immediately for treatment advice. •Have person sip a glass of water if able to swallow. •Do not induce vomiting unless told to do so by the poison control center or doctor. •Do not give anything by mouth to an unconscious person.

Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

NOTE TO PHYSICIAN


Probable mucosal damage may contraindicate the use of gastric lavage.

ENVIRONMENTAL HAZARDS

This product is toxic to fish and other aquatic organisms.

PHYSICAL/CHEMICAL HAZARDS

Keep away from combustible organic or other readily oxidizable materials. Do not mix with any other chemicals.



**ASEPTROL®
DENTAL WATER TREATMENT**

**For the Cleaning and Control of
Microbial Contamination in Dental Unit
Waterlines**

Active Ingredients:
Sodium Chlorite 4.0%
Sodium Dichloroisocyanurate dihydrate ... 1.0%
Other Ingredients: 95.0%
Total: 100.0%

KEEP OUT OF REACH OF CHILDREN

DANGER

See [Side] {or} [Back] Panels for First Aid and Precautionary Statements

EPA Reg. No. 70060-25
EPA Est. No.

Net Weight:

BASF Catalysts, LLC
100 Campus Dr.
Florham Park, NJ 07932

Patent: 6,699,404

Sodium chlorite in this product is from Spain

[contains]



[Aseptrol and the Aseptrol logo are trademarks of Engelhard Corporation]

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

INITIAL TREATMENT:

1. Add two (2) 1.7 gram Aseptrol® Dental Water Treatment tablets for every 750 ml of potable water. Allow the tablets to dissolve.
2. Flow about half of the treated water through the dental waterlines and let the treated water sit within the lines for 5 minutes. Purge the remaining water through the lines.
3. Repeat steps 1 and 2.
4. Rinse and refill the dental unit water bottle with potable water and flush at least 500 ml (1/2 liter) of water through the lines to remove any residual chemical.
5. Flush at least 500 ml of potable water through the lines to remove any residual chemical.

ROUTINE TREATMENT:

1. Every week or as needed*, add one (1) 1.7 gram Aseptrol® Dental Water Treatment for every 750 ml of potable water. Allow the tablet to dissolve.
 2. Flow about half of the treated water through the dental waterlines and let the treated water sit within the lines for 5 minutes. Purge the remaining water through the lines.
 3. Flush at least 500 ml of potable water through the lines to remove any residual chemical.
- Note: Periodic cleaning may be done on any day of the week. *Treatment should be done if bacteria count is greater than 500 cfu/ml as recommended by the CDC.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal. Store in a cool, dry area away from sunlight and heat and in an area inaccessible to children. Keep product from contacting air prior to use. Place empty foil pouch in trash.



Under the Federal Insecticide, Fungicide, and Rodenticide Act as amended, for the pesticide, registered under EPA Reg. No. 70060-25

282