

70144-1

11/13/2008

1049

UNITED STATES ENVIRONMENTAL PROTECTION AGENCY



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY  
WASHINGTON, D.C. 20460

OFFICE OF PREVENTION, PESTICIDES AND TOXIC SUBSTANCES

NOV 13 2008

Ms. Kim Davis  
Consultant for Micro-Scientific Industries, Inc  
30856 Rocky Road  
Greely, CO 80631-9375

Subject: Opti-Cide 3  
EPA Registration No.: 70144-1  
Amendment Date: July 2, 2008  
EPA Receipt Date: July 3, 2008

Dear Ms. Davis,

The following amendment submitted in connection with registration under FIFRA section 3(c)(7)(9) is acceptable with the conditions listed below.

- Addition of microorganism, Chlamydomonas reinhardtii

**Conditions**

Revise the label as follows:

- 1) Revise the 3<sup>rd</sup>, 4<sup>th</sup>, and 5<sup>th</sup> claims on page one by deleting the term, decontamination because this product is not a liquid chemical sterilant. The term, decontamination, has been associated with liquid chemical sterilants implying that the product mitigates all forms and species of microorganisms.
- 2) Revise the 13<sup>th</sup> claim on page one by deleting the term, biodegradable because it is an implied safety claim that is false and misleading (40 CFR 156.10 (a)(5)I).
- 3) Revise the 16<sup>th</sup> claim on page one by deleting the phrase, "For Professional Use" because it is considered false and misleading since the product is not classified as a restricted use pesticide.
- 4) Delete the organism, Clostridium difficile (vegetative), from pages one and two.

The Agency has re-evaluated its acceptance of Clostridium difficile (vegetative cells). Peer reviewed scientific literature and case studies have consistently demonstrated that the C. difficile spore is the

CONCURRENCES							
SYMBOL							
SURNAME							
DATE							

2099

UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

source of Clostridium difficile-associated disease (CDAD), and is an immediate public health concern. In light of scientific guidance and supporting documentation, the Agency is certain that claims against the vegetative form of C. difficile are true statements, but are used in such a way as to give false or misleading impression to the purchaser (40 CFR 156.10(a)(5)(vii), and therefore, imply heightened efficacy. The Agency considered antimicrobial pesticides to be unique because of the critical nature of the threat to public health that may result from ineffective use of products due to obsolete or misleading labeling. As a result, any reference to claims of effectiveness against C. difficile (vegetative) or C. difficile (without supporting data against the C. difficile spores) are unacceptable. To address the growing need for products in hospital/medical settings, the Agency has determined that the use of the AOAC Sporidical Activity Test, Three Step Method, and the Quantitative Carrier Test II are acceptable for supporting C. difficile sporicide claims.

- 5) Delete the term, decontaminate, from the second paragraph, first statement on page 2.
- 6) Revise the first statement in the first paragraph on page 3 by deleting the term, infectious, and stating "Opti-Cide 3 is particularly .... chance of cross-contamination from microorganisms is of concern."
- 7) Since your product is used in animal housing facilities, farms, zoos, veterinary hospitals, the label must support the following directions in accordance to DISS/TSS 17 by stating:

**Animal Housing/Veterinary Hospitals/Animal Laboratory Disinfection Directions:**

Remove all animals and feed from premises, vehicles, and enclosures. Remove all litter and manure from floors, walls, and surfaces of barn, pens, chutes, and other facilities and fixtures occupied or traversed by animals. Empty all troughs, racks, and other feeding and watering appliances. Thoroughly clean all surfaces with soap and detergent and rinse with water. Saturate all surfaces with Opti-Cide 3. Treated surfaces must remain wet for 3 minutes (5 minutes for Canine Parvovirus). Immerse all halters, ropes, and other types of equipment used in handling and restraining animals as well as forks, shovels, scrapers used in removing litter and manure. Ventilate buildings, cars, trucks, boats, and other closed spaces. Do not house animals or employ equipment until treatment has set, dried, or absorb. Thoroughly scrub all treated feed racks, mangers, troughs, automatic feeders, fountains and water with soap or detergent and rinse with potable water before reuse.

- 8) Revise the "General Surface Cleaning" heading on page three to read "General Surface Cleaning and Deodorizing."
- 9) Revise the HCV, HBV, and HIV-1 directions on page 3 to be consistent with HIV FR Notice by stating **"KILLS HBV, HCV, AND HIV-1 ON PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS..."**
- 10) Revise the statement under "Personal Protection" of the HCV, HBV, and HIV-1 directions to read *"Specific barrier protection items to be used when handling items soiled with blood or body fluids are disposable latex gloves, gowns, masks, and eye coverings."*

CONCURRENCES							
SYMBOL							
SURNAME							
DATE							

30f9

UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

11) Revise the statement under "Infectious Materials Disposal" of the HCV, HBV, and HIV-1 directions on page four by deleting the term, should, and state "must," to bring the statement in compliance with PR Notice 2000-5, Mandatory Labeling.

12) Revise the Pre-Cleaning Instrument Direction headings by removing the term "decontamination."

13) Delete the Instrument Disinfection directions and replace them with the following:

**Surgical Instrument Presoak:** Fill container with Opti-Cide 3 to allow for complete submersion of instruments/devices. Place pre-cleaned instruments in solutions for a minimum of 3 minutes (5 minutes for Canine parvovirus), then proceed with normal sterilization process. Prepare a fresh solution at least daily or when solution becomes visibly dirty or discolored.

**Ultrasonic Bath Disinfection Directions:** Use this product to disinfect hard nonporous noncritical objects compatible with Ultrasonic cleaning units. Pour Opti-Cide 3 into bath chamber. Place pre-cleaned objects into unit and operate for a minimum of 3 minutes (5 minutes for Canine parvovirus) according to manufacturer's use directions. Remove objects and rinse with sterile water or allow to air dry. Prepare a fresh solution at least daily or when solution becomes visibly dirty or discolored.

**Barber and Beauty Salon Instrument Disinfection Directions:** Immerse pre-cleaned barber/salon tools (combs, brushes, razors, manicure/pedicure tools, clippers, trimmer blades, tweezers, and scissors in Opti-Cide 3 solution. Allow objects to remain wet for 3 minutes. Rinse thoroughly and dry before use. Prepare a fresh solution at least daily or when solution becomes visibly dirty or discolored.

**Disinfect Hair Clippers and Electric Shears:** Disconnect the clippers/shears from the electrical outlet. Hold clippers/shears in downward position and spray Opti-Cide 3 onto the blades two or three times to thoroughly wet blades. Do not spray directly into housing. Allow blades to remain wet for 3 minutes. Then wipe dry with clean cloth. Lubricate as per clipper/shear manufacturer's instructions.

14) Revise the "Storage and Disposal" statements"

(Containers equal to or less than 1 gallon)

**Storage and Disposal**

Store in a cool, well ventilated area in accessible to children and persons unfamiliar with its use. Do not reuse or refill container. Wrap container in trash or offer for recycling if available.

(Containers greater than 1 gallon)

**Storage and Disposal**

Do not contaminate food or feed by storage or disposal

**Pesticide Storage:** Store in a cool, well ventilated area in accessible to children and persons unfamiliar with its use

**Pesticide Disposal:** Waste resulting from the use of this product must be disposed of on site or at an approved waste disposal facility.

			CONCURRENCES					
SYMBOL								
SURNAME								
DATE								

4089

UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

**Container Disposal:** Do not reuse or refill container. Triple rinse empty container, then offer for recycling or reconditioning, or dispose of in a sanitary landfill. If allowed by state and local authorities, dispose by burning. If burned stay out of smoke.

**Data Summary**

Data Requirement	Means of Support	Status
Chlamydiacidal Efficacy Test	MRID No. 474674-01	Acceptable, RTU in at room temperature for 3 minutes.

**General Comments**


A stamped copy of the label accepted with conditions is enclosed. Submit three (3) copies of your final printed label before distributing or selling the product bearing the revised labeling.

Submit and/or cite all data required for registration/reregistration of your product under FIFRA section 3(c)(5) when the Agency requires all registrants of similar products to submit such data.

If the above conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA section 6(e). Your release for shipment of the product bearing the amended labeling constitutes acceptance of these conditions.

Should you have any questions or comments concerning this letter, please contact Velma Noble at (703) 308-6233 or Jacqueline McFarlane at (703) 308-6416.

Sincerely,



Velma Noble  
Product Manager (31)  
Regulatory Management Branch I  
Antimicrobials Division (7510C)

Enclosures: Stamped Label

CONCURRENCES

SYMBOL								
SURNAME								
DATE								

## Opti-Cide® 3

Veterinary/Animal/Pet Care Facility Disinfectant/Cleaner  
Barber/Salon [Clipper Blade] Disinfectant/Cleaner  
Disinfectant/~~Decontamination~~ Cleaner  
Surface Disinfectant/~~Decontamination~~ Cleaner  
Medical/Surgical/Dental Equipment Disinfectant Cleaner ~~and Decontaminant~~  
Bactericidal • Fungicidal • Virucidal • Tuberculocidal

Ready-to-Use  
Disinfectant Cleaner [and] Deodorizer  
Disinfects Cleans Deodorizes

Kills Bacteria, Viruses, Fungi, Mycobacterium Tuberculosis [in 3 minutes]  
Kills Bacteria – Viruses – Fungus on hard, non-porous surfaces such as Litter Boxes, Kennels, Cages  
And for Spot Cleaning

With [Contains] Fresh {or appropriate descriptive term} Scent  
With [Contains] Rust Inhibitor, ~~and Degreaser~~ Detergent

For [Use on] Pet Care Surfaces  
[Surface/Facility] {description taken from use directions} Disinfectant Cleaner  
~~[For Professional Use Only]~~

### Active Ingredients:

n-Alkyl (60% C <sub>14</sub> , 30% C <sub>16</sub> , 5% C <sub>12</sub> , 5% C <sub>18</sub> ) dimethyl benzyl ammonium chloride .....	0.154%
n-Alkyl (68% C <sub>12</sub> , 32% C <sub>14</sub> ) dimethyl ethylbenzyl ammonium chloride .....	0.154%

Other Ingredients .....	99.692%
Total .....	100.000%

Keep Out of Reach of Children

### CAUTION

See back [side] panel for additional precautionary statements.

ACCEPTED  
with COMMENTS  
in EPA Letter Dated:

NOV 13 2008

EPA Reg. No. 70144-1 EPA Est. 1130-IL-1

Net Contents: \_\_\_\_ oz. (\_\_\_\_ ml)

Product Order Number: \_\_\_\_

{flag graphic}

Made in the U.S.A.

Under the Federal Insecticide,  
Fungicide, and Rodenticide Act as  
amended, for the pesticide,  
registered under EPA Reg. No. 70144-1

{Back/Side Panel(s)}

Opti-Cide 3 effectively kills the following microorganisms at room temperature (69°F/20°C) with a 3-minute contact time:

**[\*Viruses:]** Hepatitis B Virus (HBV), Hepatitis C Virus (HCV), Herpes Simplex II Virus (HSV 2) G Strain, Human Immunodeficiency Virus (HIV-1) (associated with AIDS), Influenza Virus Strain A2/Hong Kong, Polio I Virus ATCC VR-58, Rhinovirus ATCC VR-286 and Rotavirus Strain WA

**[Mycobacterium:]** *Mycobacterium bovis* BCG (tuberculosis) (TB)

**[Gram Positive and Gram Negative Bacteria:]** *Acinetobacter baumannii* ATCC 19606, *Chlamydophila psittaci* ATCC VR-125, *Clostridium difficile* (negative) ATCC 9609, *Staphylococcus aureus* ATCC 6538, *Pseudomonas aeruginosa* ATCC 15442, *Salmonella enterica* ATCC 10708, Methicillin-Drug Resistant *Staphylococcus aureus* (MRSA) ATCC 33591 and *Enterococcus faecium* (VRE) ATCC 51559.

[Pathogenic Fungi:] *Trichophyton mentagrophytes* ATCC 9533

**Opti-Cide 3 effectively kills the following microorganisms with a 3-minute contact time at 69°F/20°C on hard, non-porous inanimate surfaces when used as directed:** Avian Influenza A Virus ATCC VR-2072, *Bordetella bronchiseptica* (kennel cough) ATCC 10580, \*Bovine Viral Diarrhea Virus (BVDV – Surrogate Hepatitis C Virus), Canine Distemper Virus ATCC VR-128, Canine Para-Influenza Virus ATCC VR-666, *Chlamydophila psittaci* ATCC VR-125, ~~difficile (vegetative) ATCC 9689~~, Duck Hepatitis B Virus (DHBV – Surrogate Hepatitis B Virus), Feline Leukemia Virus, Herpes Simplex II Virus (HSV 2) G Strain, Infectious Bovine *Rhinotracheitis* Virus ATCC VR-188, Influenza Virus Strain A2/Hong Kong, *Mycobacterium bovis* BCG (tuberculosis) (TB), *Pseudomonas aeruginosa* ATCC 15442, Rhinovirus ATCC VR-286, Rotavirus Strain WA, *Salmonella enterica* ATCC 10708, *Staphylococcus aureus* ATCC 6538, Methicillin-Drug Resistant *Staphylococcus aureus* (MRSA) ATCC 33591, *Streptococcus equi* ATCC 33398, *Trichophyton mentagrophytes* (ringworm) ATCC 9533 and *Vancomycin*-Resistant *Enterococcus faecium* (VRE) ATCC 51559 **and effectively kills Canine Parvovirus with a 5-minute contact time.**

### DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

#### Description:

Use where control of the hazards of cross-contamination is of concern. Aids in the reduction of cross-contamination on treated surfaces. Opti-Cide 3 is particularly effective in: animal/pet care facilities, automobile interiors, barber shops, bathrooms, buses, cadaver processing areas, child care facilities (day care, infant care, schools), food processing facilities (non-food contact surfaces), funeral homes, health clubs, hotels, laboratories, long-term care facilities, medical facilities (ambulances/EMS/paramedic units, blood banks [centers], clean rooms, dental offices, drug preparation rooms, emergency rooms, hospitals, infant/child care units, isolation units [rooms] [areas], morgues, neonatal intensive care units, nursing homes, operating rooms, patient care areas, pharmacy surfaces, plasma pheresis centers, surgical centers), police vehicles, prisoner lockups, prisons, salons (beauty, manicure/pedicure, tanning, tattoo), skin care facilities, spas, train cars and areas where the control of cross-contamination is important.

When used as directed Opti-Cide 3 effectively cleans, ~~decontaminates~~ and disinfects hard, non-porous environmental surfaces such as: ambulance patient care surfaces, anesthesia equipment surfaces, animal cages (aviary, canine, feline), athletic/training facility surfaces, baby/infant care surfaces, backboards, bathing units, bed railings, blood pressure devices, brushes, cabinets, changing tables, combs, computer surfaces, countertops, CPR training devices, CT equipment/surfaces, curing lights, dental equipment surfaces/devices, diving equipment and air packs, doorknobs, food preparation surfaces, firefighter turnout gear, floors, foot spa surfaces, gonial lenses, hair cutting/clipping implements, hair rollers, hard cervical collars/braces, helmets, hospitality and hotel surfaces, infant care equipment (bassinets, bed railings, cribs, incubators, tables, warmers [exterior surface]), keyboards, laboratory equipment/surfaces, light lens covers, lumens, manicure/pedicure instruments and tables/chairs, mannequins, medical equipment surfaces/devices, microkeratome surfaces, microscopes, MRI equipment surfaces, neonatal/nursery surfaces, obstetrical equipment surfaces, operating room lights/tables, ophthalmic equipment surfaces/devices, optical wear (excluding contact lenses), oxygen hoods, patient bed surfaces, patient stretchers, patient transport surfaces, pet care surfaces, physical therapy equipment (such as bars, chairs, empty whirlpool tanks), play-care equipment, play mats, radiology equipment surfaces, refrigerators (exterior surface), respirators, respiratory care equipment surfaces, shampoo bowls, interior ship [river, ocean, sea vessel] surfaces, stalls, shower floors/walls, skin care equipment/surfaces, surgical equipment surfaces, surgical microscopes, tables, tanning bed surfaces, tattooing equipment surfaces/devices, telephones, toilets (rims, seats and exterior surfaces), tonometer tips, toys, transportation equipment surfaces, walls, wet/dry suits,

work stations, wrestling mats, X-ray equipment/surfaces and hard non-porous use sites listed on this label made of plastics (such as polycarbonate, polyvinylchloride, polypropylene and polystyrene), acrylic, non-porous vinyl, metal, stainless steel, painted, laminated and coated surfaces, Plexiglas® and glass.

Opti-Cide 3 is particularly effective for cleaning and disinfecting animal habitats, animal isolation units, animal laboratories, animal show facilities, animal/pet care facilities, aquariums, avian habitats, grooming facilities, reptile habitats, veterinary hospitals, zoo facilities and similar animal care and house facilities where the chance of cross-contamination from ~~microorganisms~~ microorganisms is of concern. When used as directed, Opti-Cide 3 effectively cleans and disinfects hard, non-porous environmental surfaces such as: animal transport units (such as animal control vehicles and horse vans), bits, bowls, bridles, cages, clippers, examination tables, harness devices, incubators, instruments, kennels, litter boxes, pet care surfaces, sheers, troughs, veterinary equipment surfaces/devices, veterinary operating room lights and tables, warmers (exterior surface) and use sites including those made of plastics (such as polycarbonate, polyvinylchloride, polypropylene and polystyrene), acrylic, non-porous vinyl, metal, stainless steel, painted, laminated and coated surfaces, Plexiglas® and glass. For disinfection of Avian influenza A, use this product only for surfaces that are not conducive to treatment by immersion or excess liquid.

#### **Surfaces:**

**Cleaning:** Prior to surface disinfection, clean surfaces by saturating surfaces with Opti-Cide 3 and allowing surfaces to remain wet for 15 seconds; then wipe clean using a fresh paper or cloth towel. Repeat, if necessary, until surfaces are visibly clean. Make sure surfaces are free of all gross filth and heavy soil. Discard paper towels or launder cloth towels before reusing.

#### **General Surface Cleaning:**

Prior to surface disinfection, remove visible debris and gross filth such as animal excrement, visible blood or tissue and bedding materials such as blankets, litter, wood chips, hay or straw. Apply undiluted Opti-Cide 3 solution directly to surfaces by using a clean applicator such as a sponge, mop or cloth or spray directly onto soiled surface and allow surfaces to remain wet for approximately 15 seconds; then rinse and/or wipe surface until visibly clean and dry.

**General [Bactericidal, Virucidal\*, Fungicidal and Tuberculocidal] Disinfection:** Reapply Opti-Cide 3 by sponge, mop, cloth, spray or dip to precleaned surfaces. Allow surface to remain wet for 3 minutes at room temperature (69°F/20°C) then wipe dry using a clean paper or cloth towel. Discard paper towels or launder cloth towels before reusing.

#### **General Surface Disinfection:**

Apply undiluted Opti-Cide 3 to precleaned surface by sponge, mop, cloth, or immersion or spray directly onto surface. Allow surface to remain visibly wet for 3 minutes (5 minutes for canine parvovirus); then wipe surface dry or allow to air dry; or rinse surface using potable water and either wipe surface dry or allow to air dry.

~~HCV, HBV and HIV-1 Precautions:~~ Opti-Cide 3 effectively inactivates HCV, HBV and HIV-1 on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in health care or settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects can be associated with the potential for transmission of HCV, HBV and/or the human immunodeficiency virus Type 1 (HIV-1) associated with AIDS. **Special instructions for cleaning and decontamination against HCV, HBV and HIV-1 of surfaces/objects soiled with blood/body fluids: Personal Protection:** Wear appropriate barrier protection such as gloves, gowns, masks ~~and~~ eye coverings. **Cleaning Procedure:** Blood and body fluids must be thoroughly cleaned from surfaces/objects before disinfection. **Contact**

**Time:** Allow the surface to remain wet for 3 minutes prior to wiping. **Infectious Materials Disposal:** Cleaning materials used that may contain blood/body fluids ~~should~~ be autoclaved and/or disposed of in accordance with local regulations for infectious materials disposal.

This product is not to be used as a terminal sterilant or high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter sterile areas of the body. This product may be used to pre-clean or decontaminate critical and semi-critical medical devices prior to sterilization or high-level disinfection.

**Instruments:**

**For Use as a Pre-Cleaning Instrument ~~Disinfection~~ Soak Solution:** Fill appropriately sized container with a sufficient amount of undiluted Opti-Cide 3 to allow for complete submersion of instruments/devices. Place instruments/devices into solution and allow to soak at room temperature (69°F/20°C) for 3 minutes, then remove instruments/devices and rinse. Follow with the appropriate cleaning and disinfection or sterilization process. Change solution as needed, i.e., when the solution becomes diluted or visibly soiled.

**For Use as a Pre-Cleaning Instrument ~~Disinfection~~ Spray:** Place instruments into suitable container. Thoroughly spray instrument surfaces with undiluted Opti-Cide 3 to thoroughly drench surfaces. Allow instruments to remain wet for 3 minutes, thoroughly rinse, then follow with the appropriate terminal cleaning and sterilization/disinfection process.

**For Use as an Ultrasonic ~~Disinfection~~ Cleaning Solution:** Place instruments/devices into the Opti-Cide 3 pre-filled ultrasonic unit. Activate unit to operate for 3 minutes. When unit shuts off, remove instruments/devices and thoroughly rinse. For critical and semi-critical instruments/devices follow decontamination with an appropriate sterilization or high-level disinfection process as required.

**For Use as an Ultrasonic Cleaning Solution:**

Place instruments/devices into an Opti-Cide 3 pre-filled ultrasonic unit. Activate unit to operate for 3 minutes. When unit shuts off, remove instruments/devices and thoroughly rinse, then dry or allow to air dry.

**For Cleaning [To Clean] Hair Clippers and Electric Shears:** Disconnect the clipper/shear from the electrical supply. Hold it in the downward position and spray Opti-Cide 3 directly onto the blades making sure the solution does not drain into the housing. Allow to remain wet for 15 seconds and wipe dry with a clean, soft cloth. Lubricate the blades per manufacturer's instructions.

~~For Use as an Instrument/Device Disinfection Solution: Clean instrument/device to remove organic soil then rinse and rough dry. Flush lumens of hollow instruments/devices before filling with solution or immersion. Immerse instrument/device into undiluted Opti-Cide 3 solution and allow to soak for 3 minutes at room temperature (69°F/20°C). Remove and rinse or wipe instrument/device dry.~~

~~**Disinfection of Instruments, Implements and Devices:**~~

~~Spray directly on visibly clean instruments, implements or devices. Allow to remain wet for 3 minutes (5 minutes for canine parvovirus), then before use wipe surface dry or allow to air dry, or rinse surface using potable water and either wipe surface dry or allow to air dry.~~



### **Disinfection Immersion Solution:**

~~Clean instruments/devices to remove organic soil and gross filth, then rinse and rough dry. Flush lumens of hollow instruments/devices before filling with solution or immersion. Immerse instruments/device into undiluted Opti-Cide 8 solution, fill lumens and allow to soak for 3 minutes (5 minutes for canine parvovirus) at 69°F/20°C. Remove and rinse exterior and inside of lumens to remove solution then wipe instruments/devices dry or allow to air dry.~~

### **Storage and Disposal**

Do not contaminate food or feed by storage or disposal.

**Storage:** Store in a cool, well-ventilated area. **Product Disposal:** Unused product may be discarded by pouring down the drain and flushing drain with water. **Container Disposal:** Triple rinse empty container, then offer for recycling or reconditioning; or puncture and dispose of in a sanitary landfill; or incinerate; or if allowed by state and local authorities, dispose of by burning. If burned, stay out of smoke.

## **PRECAUTIONARY STATEMENTS**

### **Hazards to Humans and Domestic Animals**

**CAUTION:** Harmful if inhaled or absorbed through skin. Causes moderate eye irritation. Avoid breathing vapors and contact with eyes, skin or clothing. Wash thoroughly with soap and water after handling. Remove contaminated clothing and wash clothing before reuse.

### **First Aid**

<b>If Inhaled</b>	<ul style="list-style-type: none"><li>• Move person to fresh air.</li><li>• If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible.</li></ul>
<b>If on Skin or Clothing</b>	<ul style="list-style-type: none"><li>• Take off contaminated clothing.</li><li>• Rinse skin immediately with plenty of water for 15-20 minutes.</li></ul>
<b>If in Eyes</b>	<ul style="list-style-type: none"><li>• Hold eye open and rinse slowly and gently with water for 15-20 minutes.</li><li>• Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.</li></ul>

Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. You may also contact 1-800-255-3924 for emergency medical treatment information.

{The First Aid statements' grid format will be used if market label space permits; otherwise a paragraph format will be used.}

### **Physical or Chemical Hazards**

Do not use or store near heat or open flame.

**Manufactured in the U.S. A. for:**  
**Micro-Scientific Industries, Inc.**

[Corporate Offices:] 1225 Carnegie Street; Rolling Meadows, IL 60008 U.S.A.  
Telephone: 1-888-253-2536  
www.opticide.com

Order No: \_\_\_\_\_

[ ] Denotes alternate language

{ } Denotes language that does not appear on the market label