## UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

December 30, 2004

Timothy H. Dickens Registered Agent for Medentech Ltd. c/o Scientific Coordination, Inc. 4629 Cherry Valley Drive Rockville, MD 20853

Subject:

KLOR-KLEEN

EPA Registration No. 71847-2 Application Date: December 1, 2004 Receipt Date: December 3, 2004

Dear Mr. Dickens:

This acknowledges receipt of your notification, submitted under the provision of PR Notice 98-10, FIFRA section 3(c)9.

## Proposed Notification

addition of animal pathogens to label

## **General Comments**

Based on a review of the material submitted, the following comments apply:

- 1) The notification application is acceptable and a copy has been inserted in your file for future reference.
- 2) It is not Agency policy to stamp labels for changes submitted by notification. In order to obtain a stamped label, you must submit these changes as an Application for Amended Registration.

Should you have any questions or comments concerning this letter, please contact me at (703) 308-6345.

Sincerely,

Wanda Y. Mitchell
Product Reviewer (32)
Regulatory Management Branch II
Antimicrobials Division (7510C)

	C0	NCURRENCES		
YMBOL 17510C 1/5 C	Ċ			
TURNAME E. Bera /	( )			
PATE 12/30/04	104			
EPA Form 1320-1A (1/90)	Peinted	Printed on Paguiled Pages OFFICIAL FILE		

#### Scientific Coordination, Inc.

Helping Technology Transfer 4629 Cherry Valley Drive Rockville, MD 20853 December 1, 2004

Document Processing Desk (NOTIF) Office of Pesticide Programs (7504C) U.S. Environmental Protection Agency 1801 South Bell St. Arlington, VA 22202-4501

Subject: KLOR-KLEEN<sup>TM</sup>
EPA Reg. No. 71847-2
NOTIFICATION

#### Dear Sir or Madam:

By notification, our Client, Medentech Ltd. is adding animal pathogens to their label in accordance with Section II.B of PR Notice 98-10. Further, our Client does not believe that efficacy data is required as the proposed microorganisms listed on the enclosed draft label are solely pathogenic to animals (DIS/TSS-16).

Enclosed, please find the following:

- 1. Application for Pesticide (Other), EPA Form 8570-1
- 2. One copy of a draft label as a reference illustrating all label revisions.

H. Dickens

3. Two copies of the final label.

We would appreciate receiving a copy of the label stamped approved -notification.

If you have any questions or comments regarding this submission, please do not hesitate to contact me at (301) 570-4390 or by e-mail at dickensth@aol.com.

Sincerely,

Timothy H. Dickens, Ph.D.

Registered Agent

Encl./as listed

· NO	TIFICATION	Form Ap	prove	d.	OME	No. 2	070-00	60. Ap	oroval (	3/ expires 2-28
United States  Environmental Protection A  Washington, DC 20460		ency	Registrati Amendme				OPP	OPP Identifier Number		
	Application for	Pesticide - Sec	tion	ì						
1. Company/Product Nur 71847 - 2	nber	2. EPA Product Man Emily H. Mitchell	-				3.1	Proposed		ification
4. Company/Product (Na Klor-Kleen™	PM# 32									
Medentech Ltd., c/ 4629 Cherry Valley Rockville, MD 208	<i>?</i>	6. Expedited Rev (b)(i), my product to: EPA Reg. No Product Name	is sim	nila	r or	identi	cal in c			
·	Sec	ction - II								
Amendment - Exp	response to Agency letter July 18, 2003	Final printe Agency lett "Me Too" / Other - Exp	ter dat Applic	ted catio	on.	SONSO	to			
"Notification to add anima 98-10 and EPA regulation understand that it is a vio	itional page(s) if necessary. (For section I and Sal pathogens in accordance with Section II.B of PR ns at 40 CFR 152.46, and no other changes have blation of 18 U.S.C. Sec. 1001 to willfully make any ce 98-10 and 40 CFR 152.46, this product may be of FIFRA."	Notice 98-10." "This noting een made to the labeling alse statement to EPA. I	or the	e co	onfide inders	ential st stand th	atemen	t of form	ula of th	nis product. I not consisten

Section - III 1. Material This Product Will Be Packaged In: Unit Packaging Child-Resistent Packaging Water Soluble Packaging 2. Type of Container Yes Metal Yes Yes Plastic No Glass Paper If "Yes" If "Yes" No. per No. per Certification must Unit Packaging wgt. container Peckage wgt container Other (Specify) be submitted 3.7 kg 300 tabs 3. Location of Net Contents Information 5. Location of Label Directions 4. Size(s) Retail Container On Labeling Accompanying Product 4 Liter Capacity Container 6. Manner in Which Label is Affixed to Product Lithograph Paper glued Stenciled Other Section - IV 1. Contact Point (Complete items directly below for identification of individual to be contacted, if necessary, to process this application.) Telephone No. (Include Area Code) Timothy H. Dickens, Ph.D. Registered Agent (301) 570 4390 6 Date Application Certification Received I certify that the statements I have made on this form and all attachments thereto are true, accurate and complete. I acknowledge that any knowlingtly false or misleading statement may be punishable by fine or imprisonment or (Stamped) both under applicable law. 3. Title 2. Signature Registered Agent 5. Date 4. Typed Name December 01, 2004 Timothy H. Dickens, Ph.D.

EPA Form 8570-1 (Rev. 3-94) Previous editions are obsolete.

# KLOR-KLEEN™

Effervescent Disinfectant Tablets for Hospitals and Institutional Use.

**ACTIVE INGREDIENT:** 

Sodium Dichloroisocyanurate ........... 48.21% 

## KEEP OUT OF REACH OF CHILDREN **DANGER**

See first aid /statement of practical treatment on back panel

Each tablet will produce a 937ppm disinfectant solution when dissolved in one gallon of water.

KLOR-KLEEN<sup>™</sup> is effective against the following microorganisms on precleaned, hard, nonporous, inanimate surfaces: Salmonella choleraesuis, Staphylococcus aureus, Pseudomonas aeruginosa, Klebsiella pneumoniae, Staphylococcus epidermidis, Escherichia coli O157:H7, Staphylococcus aureus -MRSA & GRSA, Clostridium difficile, Mycobacterium bovis - BCG, Poliovirus type 1, Herpes simplex virus type 1, Hepatitis A virus, Human Immunodeficiency Virus Type 1(associated with AIDS) or (AIDS Virus), Canine Parvovirus, Newcastle Disease virus, Pseudorabies,

Underlined denotes new verbiage.

Manufactured by:

, Medentech Ltd.

Whitemill Industrial Estate, Clonard Road Wexford, Ireland Tel: 353 53 60040 Fax: 353 53 41271

E.P.A. Reg. No. 71847-2

E.P.A. Est. No. 71847-IRL-001

**NET CONTENTS** 

**BATCH NO** 

KLOR-KLEEN NOTIFICATION REFERENCE LABEL

KLOR-KLEEN™ is designed to provide effective cleaning, deodorizing and disinfection in areas where it is of prime importance to control the hazard of cross contamination on treated precleaned, hard, nonporous, inanimate surfaces.

## KLOR-KLEEN<sup>TM</sup> is for use in:

- Hospitals, nursing homes, medical and dental offices and clinics, operating rooms, isolation wards, and medical research facilities.
- Day care centers and nurseries.
- Veterinary clinics, animal life science laboratories, kennels, breeding and grooming establishments, pet animal quarters, zoos, pet shops, and other animal care facilities.

KLOR-KLEEN<sup>™</sup> is a disinfectant and cleaner that disinfects, cleans, and deodorizes precleaned, hard, nonporous, inanimate surfaces. This pre-cleaning may be accomplished with any cleaner solution including Klor-Kleen

KLOR-KLEEN™ provides effective cleaning strength that will not dull high gloss floor finishes with repeated use.

## PRECAUTIONARY STATEMENTS

#### HAZARDS TO HUMANS AND DOMESTIC ANIMALS

Danger: Corrosive. Causes irreversible eye damage. Harmful if swallowed, inhaled or absorbed through skin. Do not get in eyes, on skin or clothing. Avoid breathing dust. Wear rubber gloves and safety glasses or face shield when making up solution. Wash thoroughly with soap and water after handling. Remove contaminated clothing and wash before reuse.

FIRST AID		
If inhaled	<ul> <li>Move person to fresh air</li> <li>If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible.</li> <li>Call a poison control center or doctor for further treatment advice</li> </ul>	
lf on skin or clothing	<ul> <li>Take off contaminated clothing</li> <li>Rinse skin immediately with plenty of water for 15 – 20 minutes.</li> <li>Call a poison control center or doctor for treatment advice</li> </ul>	
If in eyes	<ul> <li>Hold eye open and rinse slowly and gently with water for 15 – 20 minutes.</li> <li>Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.</li> <li>Call a poison control center or doctor for treatment advice</li> </ul>	
If swallowed	<ul> <li>Call poison control center or doctor immediately for treatment advice.</li> <li>Have person sip a glass of water if able to swallow.</li> <li>Do not induce vomiting unless told to do so by the poison control center or doctor.</li> <li>Do not give anything by mouth to an unconscious person</li> </ul>	
	: HOT LINE NUMBER	
	container or label with you when calling a poison control center or doctor, or going for any also contact 1-800-xxx-xxx for emergency medical treatment information.  NOTE TO PHYSICIAN	
	Probable mucosal damage may contraindicate the use of gastric lavage	

## **DIRECTIONS FOR USE**

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or, (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. (This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection)

#### PREPARATION OF USE SOLUTION:

Add 1 KLOR-KLEEN<sup>™</sup> tablet to 1 gallon of water (25°C/77°F) to obtain a concentration of 937ppm (mg/Liter) of available chlorine. This solution should be freshly prepared. Use within day. Keep in a properly labelled container to protect against contamination. Unused solution should be discarded.

KLOR-KLEEN™ is a Hospital Use Disinfectant at 1 tablet per gallon.

#### Disinfection/\*Virucidal Directions:

Apply use solution to precleaned, hard, nonporous inanimate surfaces. Allow to remain wet for 10 minutes, then remove product by wiping with brush, sponge, or cloth. For heavily soiled areas, a preliminary cleaning is required. For sprayer applications use a coarse spray device. Spray 6–8 inches from surface, rub with brush, sponge, or cloth. Do not breathe spray mist.

Before using this product, food products and packaging materials must be removed from the room or carefully protected. After use, all surfaces in the area must be thoroughly rinsed with potable water.

Add 1 tablet per gallon of water to disinfect hard, non-porous surfaces. Treated surfaces must remain wet for 10 minutes. Prepare a fresh solution daily or when use dilution becomes diluted or soiled.

Effective as a tuberculocide in 10 minutes at 20°C against Mycobacterium bovis - BCG

\* Kills HIV-1 on PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS).

## SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

PERSONAL PROTECTION: Specific barrier protection items to be used when handling items soiled with blood or body fluids are disposable latex gloves, gowns, masks, or eye coverings.

CLEANING PROCEDURE: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of KLOR-KLEEN™. This pre-cleaning may be accomplished with any cleaning solution including KLOR-KLEEN™.

DISPOSAL OF INFECTIOUS MATERIALS: Blood and other body fluids should be autoclaved and disposed of according to federal, state and local regulations for infectious waste disposal. CONTACT TIME: Leave surfaces wet for 10 minutes.

#### **ANIMAL PATHOGENS:**

When used at one tablet per gallon of water, applied as outlined under Disinfection/Virucidal Directions, Klor-Kleen is effective against the following animal pathogens: Infectious Canine hepatitis, Gumboro disease, Transmissable gastroenteritis (TGE), Swine Vesicular disease, Teschen/Talfan disease, African swine fever, Avian influenza, Porcine parvovirus, Hog cholera/Classical swine fever, Avipox (fowl pox), Runting & Stunting virus (tenysynovitis), Actinobacillus pleuropheumonia, Bordetella bronchiseptica (rhinitis), Brachispira (Treponema/Serpuline) hyodysenteriae (swine dysentry), Streptococcus uberis, Streptococcus dysgalactiae, Clostridium perfringes USDA

## STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage and disposal.

**STORAGE:** Store in a cool, dry, well-ventilated area. Do not store at temperatures above 40°C/104°F.

**PESTICIDE DISPOSAL** – Do not contaminate water, food or feed by storage or disposal. Wastes resulting from the use of this product may be disposed of on-site, in a toilet, or service sink served by a sanitary sewer or an approved waste disposal facility for pesticides.

**CONTAINER DISPOSAL** – Do not reuse empty container. Rinse empty container thoroughly with water to dissolve all material before discarding. Securely wrap container in several layers of newspaper and discard in trash.

71847-2Notif.Ref.120104-01

(