

74559-3

01/21/2010

1/11



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
WASHINGTON, D.C. 20460

OFFICE OF
PREVENTION, PESTICIDES AND
TOXIC SUBSTANCES

Mr. David Swain
Regulatory Agent for,
Virox Technologies, Inc.
2770 Conventry Road
Oakville, Ontario L6H 6R1 Canada

JAN 21 2010

Mail to: Attn: David Swain
Scientific & Regulatory Consultants, Inc.
P. O. Box 1014
Columbia City, IN 46725

Subject: Accel Tb Wipes
EPA Registration Number 74559-3
Your Amendment Dated October 22nd, 2009
EPA Received Date October 23rd, 2009

The amendment referred to above, submitted in connection with under section 3(c)(7)(A) of the Federal Insecticide, Fungicide, and Rodenticide Act, FIFRA, as amended, to add the organism, "Rhinovirus Type 14" to the product labeling, is acceptable.


A stamped copy of the labeling is enclosed. Submit three (3) copies of your final printed labeling before distributing or selling the product bearing the revised labeling.

Submit and/or cite all data required for registration /reregistration of your product under FIFRA section 3(c)(5) when the Agency requires all registrants of similar products to submit such data.

If the above conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA section 6(e). Your release for shipment of the product bearing the amended labeling constitutes acceptance of these conditions.

If you have any questions concerning this letter, please contact Karen M. Leavy-Munk at (703)-308-6237.

Sincerely,

A handwritten signature in black ink, appearing to read "Karen M. Leavy-Munk". The signature is fluid and cursive, with the first name "Karen" being the most prominent.

Marshall Swindell
Product Manager 33
Regulatory Management Branch I
Antimicrobial Division(7510C)

3/11
{Administrative note: any language which appears between { } and is italicized will NOT appear on final printed labels.
Any language that appears between [] is considered interchangeable language which may appear on final printed labeling}

Accel Tb Wipes

[Ready to Use] [Easy to Use] [Towelette] [Wipe]
One-Step Disinfectant and Deodorizing Wipes
Sanitizer • Sanitizing Wipes • Cleaner • Disinfectant • Deodorizer
*Virucidal • Bactericidal • [Fungicidal] • Tuberculocidal
Meets OSHA Bloodborne Pathogen Standard for HBV, HCV & HIV-1
Kills Norovirus • Kills Rhinovirus [Type 14]
Use on non-food contact surfaces, including: [Grocery Cart Handles] [Cash Registers]
[Telephones] [Doorknobs]
Just pull, wipe & toss! • Eliminates Buckets • No need to Spray!
Pre-moistened • Pre-measured • Always Ready-to-Use
No measuring or mixing required

[Hydrogen Peroxide]

For [Commercial] [Industrial] [&] [Institutional] Use [Only]

ACTIVE INGREDIENT:

Hydrogen Peroxide..... 0.5%

INERT [OTHER] INGREDIENTS:..... 99.5%

TOTAL:..... 100.0%

KEEP OUT OF REACH OF CHILDREN

[See reference sheet [enclosed in each case] for [a complete list of organisms] [additional [features] [claims] [directions for use]] [for this product]].

Net Contents:

160 sheets (6" X 6.75") (15.24 cm X 17.14 cm)

160 sheets (8" X 9") (20.32 cm X 22.86 cm)

10 sheets (7.48" x 10") (19.0 cm x 25.4 cm) {for re-sealable package}

ACCEPTED
with COMMENTS
EPA Letter Dated:

JAN 21 2010

Under the Federal Insecticide,
Fungicide, and Rodenticide Act as
amended, for the pesticide,
registered under EPA Reg. No. 7455 9-3

4/11

{Administrative note: any language which appears between { } and is italicized will NOT appear on final printed labels.
Any language that appears between [] is considered interchangeable language which may appear on final printed labeling}

{MARKETING CLAIMS}:

A disinfectant cloth for use on telephone surfaces and other telecommunications equipment surfaces
A Hydrogen Peroxide solution impregnated in a wiping cloth
Aids in reducing [or in the reduction of] cross contamination between treated surfaces
Always ready to use
An effective counter cleaner
An effective way to clean, deodorize against odor-causing organisms and disinfect hard, non-porous surfaces
An effective way to disinfect hard, non-porous environmental surfaces
Antibacterial
Antibacterial [surface] Cleaner [cloth]
Antibacterial Cell phone Wipe
Antibacterial Surface Cleaner [Cloth]
Anti-microbial cleaner [cloth]
Anti-viral* [Anti-Virus*] [Surface] Cleaner
Assures proper dilution strength for product effectiveness
Bactericidal Cleaner
Bactericidal
Built-in Scrubbing Fibers
Clean & disinfect any room, any time
Clean & disinfect hard non-porous surfaces
Clean & disinfect hard non-porous surfaces easily
Cleans without scratching surface[s]
Cleaner
Cleaner and disinfectant towel [or cloth]
Clean[s] and deodorize[s]
Cleans and disinfects
Cleans and disinfects [cell phone surfaces] [telephone equipment] [telephone equipment surfaces] [dictating equipment surfaces] [telecommunications equipment surfaces] [devices]
Cleans and disinfects common hard, non-porous surfaces
Cleans effectively
Clean hard-to-remove stains
Cleans, deodorizes and disinfects
Cloth
Convenient
Convenient, quick & easy cleaning all around [the office] [business] [workplace]
Cuts grease
Cuts the toughest grime and grease without scratching surfaces
Disinfectant [disinfecting] surface cleaner [cloth]
Disinfects against *Staphylococcus aureus*, *Salmonella enterica*, and *Pseudomonas aeruginosa*
Disinfects environmental surfaces
Disinfects hard, non-porous environmental surfaces
Disinfects headsets and telephones
Disposable
Each pre-moistened towelette [disposable wipe] cleans, deodorizes and kills listed bacteria all around you [office] [business] [workplace]
Easy flip-top lid
Effective against flu virus [Influenza A/Hong Kong] [Hong Kong Flu]
Effective daily cleaner
Effective in 1 minute against *Staphylococcus aureus*, *Salmonella enterica*, and *Pseudomonas aeruginosa*
Effectively disinfects hard, non-porous surfaces
Effectively disinfects in the presence of organic soil [5% blood serum]
Everyday, every room disinfectant
Exhibits effective virucidal activity: kills Influenza A/Hong Kong and Herpes Simplex Type 2 [cause of Genital Herpes]
For hard, non-porous surfaces
For use in critical care areas [where control of the hazards of cross-contamination between treated surfaces is of prime importance]
For your toughest bathroom messes

ACCEPTED
with COMMENTS
EPA Letter Dated:

JAN 21 2010

Under the Federal Insecticide,
Fungicide, and Rodenticide Act as
amended, for the pesticide,
registered under EPA Reg. No. 74559-3

5/11

{Administrative note: any language which appears between { } and is italicized will NOT appear on final printed labels. Any language that appears between [] is considered interchangeable language which may appear on final printed labeling}

Grease Cutting Formula

Great for grease cutting

Great for on-the-go cleaning

Great on-the-go disinfectant [cloth]

Hard non-porous surface [cleaning &] disinfecting made easy

Heavy-duty cleaner [or cleaning] [cloth]

Kills Avian Influenza A virus on pre-cleaned environmental surfaces

Kills *Staphylococcus aureus*, *Salmonella enterica*, and *Pseudomonas aeruginosa*

Kills the flu virus Influenza A

Kills Norovirus [Norwalk virus]

Kills Human Coronavirus

Kills antibiotic-resistant organisms

Kills HIV-1 [AIDS virus] [HTLV-III_B]

Kills HCV

Kills HBV

Kills Rhinovirus [Type 14]

No mixing

No streaking

Non-woven disposable cloth

On-the-go disinfectant

Accel Tb Wipes control the spread of listed bacteria on treated hard, non-porous environmental surfaces

Accel Tb Wipes help control listed bacteria on treated surfaces while keeping your office [or workstation or workplace or business] clean and dust-free

Pre-measured

Pre-moistened [towelettes]

Pre-soaked [cloth or towelette]

Quick and easy cleaning everywhere

Quick drying

Removes tough stains

Restroom Disinfectant Wipes

Restroom Wipes [for hard, non-porous surfaces]

Scrubbing surface cleaner

Stay-moist package

Strong towelettes [or wipes]

Telephone [cleaner [and] disinfectant] wipes

These wipes are a convenient way to disinfect surfaces all around your [business] [office] [workplace]

Towel

Use on hard, non-porous surfaces only. [Test wipe on a small inconspicuous area first]

Use to disinfectant shopping [or grocery] cart handles

Versatile disinfectant [or disinfecting] cleaner

Virucidal* [or] Kills viruses* *Influenza A/Hong Kong [Hong Kong flu], Herpes Simplex Type 2, HCV, HBV, HIV-1

Virus killing cleaner*

Wipe down grocery [or shopping] cart handles

[FEATURES, CLAIMS & USES:]

[Notice to User: Accel TB may not be suitable for all items or surfaces. Prolonged exposure to copper, brass, anodized aluminum and some coatings should be avoided. If you are uncertain as to the material composition of your item, confirm with the manufacturer before proceeding [or visit the Compatibility section on www.virox.com for a listing of compatible materials.]]

For Use in Health Care, Veterinary, Animal Housing, Dental, [Life Science] Laboratory, Medical, Industrial & Institutional Facilities, Food Processing Establishments, and Hospitals

This product [These] is [are] [is a] non-woven, disposable [towelette[s]] [wipe[s]], pre-saturated with Hydrogen Peroxide [based] disinfectant cleaner.

**ACCEPTED
with COMMENTS
EPA Letter Dated:**

JAN 21 2010

Under the Federal Insecticide,
Fungicide, and Rodenticide Act as
amended, for the pesticide,
registered under EPA Reg. No. 74559-3

6/11

{Administrative note: any language which appears between { } and is italicized will NOT appear on final printed labels. Any language that appears between [] is considered interchangeable language which may appear on final printed labeling}

This product is a one-step [hospital-use] germicidal [disinfectant] cleaner and deodorant [odor-counteractant] [odor, neutralizer] designed for general cleaning, [and] disinfecting, [and] [deodorizing] [of] hard, non-porous environmental surfaces. It cleans quickly by removing dirt, grime, body oils and other common soils found in hospitals, nursing homes, [acute] [and] [long-term] care facilities, office buildings, schools and colleges [universities], [day care centers], [medical] [dental] offices, funeral homes, transportation facilities, train [depots] [stations] bus [depots] [stations], airports, casinos, cruise ships, correctional facilities, athletic facilities, athletic equipment, playgrounds, locker rooms, health clubs, exercise equipment, barber/beauty [salons] [shops], veterinary clinics, pet shops, [equine farms], animal life science laboratories, hotels, motels, public areas and restrooms, homes cleaned by [home cleaning services,] [professional cleaning services,] [medical assistance services,] and office buildings.

This product's [non-dulling] [advanced] [quick cleaning] [convenient] [easy-to-use] [highly effective] formula eliminates the time and labor normally required for rinsing. It is designed for use on [Use daily on] [the following] [hard, non-porous environmental [inanimate] surfaces]: vinyl, painted surfaces, plastic [surfaces], glass, mirrors, glazed ceramic, glazed porcelain, chrome, stainless steel, aluminum, laminated surfaces and baked enamel surfaces associated with floors, walls ceilings, tables, chairs, non-food contact countertops, telephones, fixtures, athletic equipment, exercise equipment, glazed tile, exterior surfaces of toilets, exterior surfaces of urinals, sinks, sink basins, faucets, tubs, chrome, stainless steel, shower rooms and locker room areas – any washable [non-food contact] surface where disinfection is required. This product's [non-dulling] [advanced] [quick cleaning] [convenient] [easy-to-use] [highly effective] formula eliminates the time and labor normally required for rinsing. Do not use on glasses, dishes and utensils. Not recommended for use on brass, copper or marble.

[Hospitals/Healthcare Facilities/Dental Offices/Life Science Laboratories]

This product is one-step [hospital-use] [laboratory-use] germicidal [disinfectant] cleaner and deodorant [odor-counteractant] [odor neutralizer] towelette designed for general cleaning, [and] disinfecting, [deodorizing] [of] hard, non-porous environmental surfaces. This product cleans quickly by removing dirt, grime, food residue, dead skin, blood, body soils and other organic matter commonly found in Dental Clinics, [Hospitals], [in Health Care Facilities], [laboratory settings]. It also eliminates odors, leaving [restroom] surfaces smelling clean and fresh. Use where odors are a problem.

This product cleans, disinfects and deodorizes hard, non-porous environmental [hospital] [medical] [dental] [laboratory] surfaces in one step [with no rinsing required]. Its non-abrasive formula is designed for use on [Use daily on] [the following] [hard, non-porous environmental surfaces:] [vinyl], [Naugahyde], painted surfaces, plastic [surfaces], glazed ceramic, glazed porcelain, chrome, stainless steel, laminated surfaces and baked enamel surfaces associated with floors, walls, ceilings, tables, chairs, non-food contact countertops, telephones, fixtures, glazed tile, toilet exteriors, urinal exteriors, sinks found in [health care facilities], [hospitals], [dental offices], [patient rooms], [operating rooms], [ICU areas], [shower rooms], and [locker rooms] [fume hoods and animal cages [found in laboratory settings]]. [It can be used to pre-clean and disinfect hospital items: wheelchairs, [hospital] [patient] bedrails and linings [patient bed linings] [vinyl mattress covers], wash basins, bed pans, non-critical medical equipment surfaces – almost any washable non-food contact surface [where disinfection is required].

[[Bathroom] [Restroom] Disinfectant Cleaner:]

This product is a one-step disinfectant cleaner and deodorant [odor-counteractant] [odor-neutralizer] designed for general cleaning, [and] disinfecting, [deodorizing] on hard, non-porous environmental surfaces.

This product cleans, disinfects and deodorizes surfaces by killing odor-causing microorganisms. Its [non-dulling] [advanced] [quick cleaning] [convenient] [easy-to-use] [highly effective] formula is designed for use on [Use it daily] on [hard, non-porous environmental] restroom surfaces: glazed ceramic [restroom] tiles, glazed porcelain, chrome, stainless steel and plastic surfaces associated with floors, walls, fixtures, exterior [surfaces] of toilets, exterior [surfaces] of urinals, sinks, shower rooms and locker rooms. Not recommended for use on brass, copper or marble.

[Animal Housing Facilities]

This product cleans by removing dirt, grime, blood, urine, fecal matter and other common soils found in animal housing facilities, grooming facilities, farms, kennels, pet stores, veterinary clinics, laboratories or other small animal facilities. It [also] eliminates odors leaving surfaces smelling clean and fresh.

ACCEPTED
with COMMENTS
EPA Letter Dated:
JAN 21 2010
Under the Federal Insecticide,
Fungicide, and Rodenticide Act as
amended for the pesticide,
registered under EPA Reg. No. 74559-3

7/11
{Administrative note: any language which appears between { } and is italicized will NOT appear on final printed labels. Any language that appears between [] is considered interchangeable language which may appear on final printed labeling}

This product cleans, disinfects and deodorizes hard, non-porous environmental surfaces in one step. Its [non-dulling] [advanced] [quick cleaning] [convenient] [easy-to-use] [highly effective] formula is designed for use on [Use it daily] [Use daily to clean and disinfect] [hard, non-porous surfaces:] plated or stainless steel, chrome, glazed porcelain, glazed tile, laminated surfaces, [associated with floors, walls, non-food contact countertops, cages, kennels, animal equipment] found in [barns, pens and stalls] animal housing facilities.

[CLAIMS:]

[When used as directed, this product [is] highly effective against a wide variety [broad-spectrum] of pathogenic microorganisms [including bacteria, tuberculosis, antibiotic resistant bacteria, viruses and fungi. See reference sheet [enclosed in each case] for a complete list or organism].

HOSPITAL DISINFECTANT: BACTERICIDAL in the presence of 5% serum load and 1 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

Bactericidal Activity:

Pseudomonas aeruginosa [ATCC 15442]

Staphylococcus aureus [ATCC 6538]

Escherichia coli O157:H7

Salmonella enterica [ATCC 10708]

Acinetobacter baumannii [ATCC 19606]

[formally know as *Salmonella choleraesuis*]

ACCEPTED
with COMMENTS
EPA Letter Dated:

Antibiotic-Resistant [Strains of] Bacteria [Antibiotic-Resistant Bactericidal Activity]-

Staphylococcus aureus MRSA [ATCC 33592]

Enterococcus faecalis VRE [ATCC 51575]

***VIRUCIDAL** in the presence of 5% serum load and 1 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

JAN 21 2010
Under the Federal Insecticide,
Fungicide, and Rodenticide Act as
intended for the pesticide,
registered under EPA Reg. No. 74559-3

Viruses [*Virucidal Activity]:

Herpes Simplex Virus, Type 1 [HSV-1]

[ATCC VR-733]

Influenza A/Hong Kong [ATCC VR-544]

Herpes Simplex Virus, Type 2 [HSV-2]

[ATCC VR-734]

Poliovirus Type 1, Strain Brunhilde

[ATCC VR-1000]

Human Coronavirus [ATCC VR-740]

Rhinovirus, Type 14 [ATCC VR-284]

Rotavirus WA

Kills HIV-1 [AIDS virus] [HTLV-III_B] when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

Kills HBV & HCV when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

Kills Norovirus [Feline Calicivirus as the surrogate] [VR-782] when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

[Veterinary Viruses:]

Avian Influenza A, [H3N2] [ATCC VR-2072]

Feline Calicivirus, Strain F9 [ATCC VR-782]

TUBERCULOCIDAL in the presence of 5% serum load and 5 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

Mycobacterium bovis [Tb] [OT 451C150]

[FUNGICIDAL] This product is fungicidal against the pathogenic fungi, *Trichophyton mentagrophytes* [ATCC 9533] [the fungus which causes Athlete's Foot] in the presence of 5% serum load and 10 minute contact time at 20°C (68°F) when used as directed on hard surfaces found in bathroom, shower stalls, locker rooms, or other clean, non-porous, hard surfaces commonly contacted by bare feet.]

[*Trichophyton mentagrophytes* [ATCC 9533]]

Non-Food Contact Surface Sanitizer for hard non-porous surfaces - in the presence of 5% serum load and 30 second contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

Klebsiella pneumoniae [ATCC 4352]

Staphylococcus aureus [ATCC 6538]

Escherichia coli O157:H7 [ATCC 35150]

Salmonella enterica [ATCC 10708]

[formally know as *Salmonella choleraesuis*]

Staphylococcus aureus MRSA [ATCC 33592]

Pseudomonas aeruginosa [ATCC 15442]

Enterococcus faecalis VRE [ATCC 51575]

8/11

{Administrative note: any language which appears between { } and is italicized will NOT appear on final printed labels. Any language that appears between [] is considered interchangeable language which may appear on final printed labeling}

Malodor[s] [Activity] [Counteractancy] – eliminates [destroys] odors and odor-causing bacteria in restroom areas, behind and under sinks and counters, and storage areas [and other places where bacterial growth can cause malodors].

DIRECTIONS FOR USE:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Not for personal cleansing. This is not a baby wipe. Do not flush in toilet.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

For Use as [Daily] [One-Step] Cleaner/Disinfectant:

1. Pre-clean heavily soiled areas
2. Pull towelette from dispenser [canister] and wipe hard, non-porous environmental surfaces.
3. All surfaces must remain visibly wet for 1 minute. Use a 5 minute contact time for Tuberculosis. [Use a 10 minute contact time for fungi.]
4. Allow to air dry [Wipe surface dry] [or rinse] [with potable water if necessary].

Note: All food contact surfaces must be rinsed with potable water. Do not use on glassware, utensils, or dishes.

For Use as [Daily] [Bathroom] [Restroom] [Shower Rooms] [Locker Rooms] [Athletic Equipment] [One-Step] Cleaner/Disinfectant:

1. Pre-clean heavily soiled areas
2. Pull towelette from dispenser [canister] and wipe hard, non-porous environmental surfaces.
3. To disinfect, all surfaces must remain wet for one [1] minutes. Use a 5 minute contact time for Tuberculosis. [Use a 10 minute contact time for fungi.]
4. Wipe surfaces dry [or rinse] [allow to air dry].

***KILLS HIV AND HCV AND HBV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS**

PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings [Hospitals, Nursing Homes] and other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/ body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS) Hepatitis B Virus (HBV), or Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV OR HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

Personal Protection: Disposable latex or vinyl gloves, gowns, face masks, or eye coverings as appropriate, must be worn during all cleaning of body fluids, blood, and decontamination procedures.

Cleaning Procedures: Blood/body fluids must be thoroughly cleaned from surfaces and objects before application of this product.

Contact Time: Allow surface to remain wet for 1 minute to kill HIV-1, HBV and HCV. Use a 5 minute contact time for Tuberculosis and 10 minute contact for fungi.

Disposal of Infectious Material: Blood/body fluids should be autoclaved and disposed of according to Federal, State, and local regulations for infectious waste disposal.

For Use as a Tuberculocide:

1. Pre-clean heavily soiled areas.
2. Pull towelette from dispenser [canister] and wipe hard, non-porous environmental surfaces.
3. Allow the surface to remain wet for 5 minutes. Use additional towelettes if needed to ensure the correct wet contact time.
4. Wipe surfaces dry [or rinse] [or allow to air dry].
5. Change towelette after each use [application].

ACCEPTED
with COMMENTS
EPA Letter Dated:

JAN 21 2010

Under the Federal Insecticide,
Fungicide, and Rodenticide Act as
amended for the pesticide,
registered under EPA Reg. No. 74559-3

102209
Page 6 of 9
EPA Reg. No. 74559-3

9/11

{Administrative note: any language which appears between { } and is italicized will **NOT** appear on final printed labels. Any language that appears between [] is considered interchangeable language which may appear on final printed labeling}

To Kill Fungi:

1. Pre-clean heavily soiled areas.
2. Pull towelette from dispenser [canister] and wipe hard, non-porous environmental surfaces.
3. All surfaces must remain visibly wet for 10 minutes. Use additional towelettes if needed to ensure the correct wet contact time.
4. Wipe surfaces dry [or rinse] [or allow to air dry].

For use as a Non-Food Contact Sanitizer:

1. Pre-clean heavily soiled areas.
2. Pull towelette from dispenser [canister] and wipe hard, non-porous environmental surfaces.
3. Allow the surface to remain wet for 30 seconds.
4. Wipe surfaces dry [or rinse] [or allow to air dry].
5. Change towelette after each use [application].

For Use as a Cleaner and/or Deodorizer:

Pull towelette from dispenser [canister] and wipe hard, non-porous environmental surfaces. Allow to air dry [or rinse [with potable water if necessary]].

**Animal Housing Facilities [including Veterinary Clinics, [Equine] Farms, Livestock, Swine and Poultry Houses]:
For cleaning and disinfecting the following hard, non-porous surfaces: feeding and watering equipment, utensils, instruments, cages, kennels, stables, catteries, etc.**

1. Remove all animals and feeds from premises, animal transportation vehicles, crates, etc.
2. Remove all litter, droppings and manure from floors, walls and surfaces of facilities [including pens, chutes, and barns] occupied or traversed by animals.
3. Empty all troughs, racks, and other feeding or watering appliances.
4. Thoroughly clean all surfaces with soap and detergent and rinse with potable water. Pull towelette from dispenser [canister] and wipe hard, non-porous environmental surfaces.
5. All surfaces must remain visibly wet for 1 minute.
6. Allow to air dry [or rinse [with potable water if necessary]].
7. Ventilate buildings and other closed spaces.
8. Do not house animals or employ equipment until treatment has been absorbed, set or dried.

NOTE: Not for use on handling and restraining equipment such as leashes, muzzles, halters or ropes, as well as forks, shovels, and scrapers used for removing litter and manure.

{Optional directions for refill pouch ONLY}

Refill Pouch

1. Open empty canister lid.
2. Open refill pouch [by cutting the top of the pouch] and remove saturated towelettes.
3. Place saturated towelettes in empty canister [invert the contents of the pouch into empty canister].
4. Replace canister lid.
5. [Remove the new lot number /[expiry date] from pouch and place it on top of the old lot number /[expiry date] sticker on the canister]

{Optional directions for refill pouch ONLY}

Re-sealable soft pack

1. Remove the seal (from the place shown on the label)
2. Remove one sheet at a time.
3. Close the seal (to prevent drying between uses)

ACCEPTED
with COMMENTS
EPA Letter Dated:

JAN 21 2010

Under the Federal Insecticide,
Fungicide, and Rodenticide Act as
amended, for the pesticide,
registered under EPA Reg. No. 74559-3

102209
Page 7 of 9
EPA Reg. No. 74559-3

10/11

{Administrative note: any language which appears between { } and is italicized will **NOT** appear on final printed labels. Any language that appears between [] is considered interchangeable language which may appear on final printed labeling}

STORAGE AND DISPOSAL:

Do not contaminate water, food or feed by storage and disposal.

Storage: Store in original container in areas inaccessible to small children. Keep securely closed.

Container Disposal: [Non-refillable container]. Do not reuse or refill this container. Rinse thoroughly. Offer for recycling if available or [wrap and] discard in trash.

{Optional for product supplied in refillable containers or refill package :}

Towelette Storage/Disposal: Dispose of [wipe] [towelette] [cloth] in trash after use. Do not flush. Reseal container label between uses to retain moisture. Refillable container. Refill this container with refill wipe only. Do not reuse this container for any other purpose. Discard in trash.

{Optional for product supplied in resealable packages :}

Towelette Storage/Disposal: Dispose of [wipe] [towelette] [cloth] in trash after use. Do not flush. Reseal container label between uses to retain moisture. Do not reuse empty [Pouch] [Container] [Package]. Discard empty [Pouch] [Container] [Package] in trash.

ENVIRONMENTAL HAZARDS {for labels with veterinary use patterns}

This product is toxic to birds, fish and aquatic invertebrates. Caution should be used when applying indoors because pets may be at risk.

Read the M.S.D.S. Before Using the Product.

M.S.D.S. is available on line at www.SciCan.com [www.virox.com]

[Si no puede leer en ingles, pregunte a se supervisor sobre las instrucciones de uso apropiadas anted de trabajar con este producto.]

EPA Reg. No. 74559-3

EPA Est. No. 74559-CAN-1

© 2008 Virox Technologies, Inc., 2770 Coventry Road, Oakville, ON L6H 6R1 Canada All Rights Reserved

AHP Logo**

** AHP technology and designs are trademarks of Virox Technologies, Inc. [used under license by SciCan Ltd.]

ACCEPTED
with COMMENTS
EPA Letter Dated:

JAN 21 2010

Under the Federal Insecticide,
Fungicide, and Rodenticide Act as
amended, for the pesticide,
registered under EPA Reg. No. 74559-3

{Administrative note: any language which appears between { } and is italicized, will NOT appear on final printed labels. Any language that appears between [] is considered interchangeable language which may appear on final printed labeling}

Virox Technologies, Inc.

Accel Tb Wipes Reference Sheet

HOSPITAL DISINFECTANT: BACTERICIDAL in the presence of 5% serum load and 1 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

Bactericidal Activity:

Pseudomonas aeruginosa [ATCC 15442] *Staphylococcus aureus* [ATCC 6538] *Escherichia coli* O157:H7 [ATCC 35150]
Salmonella enterica [ATCC 10708] *Acinetobacter baumannii* [ATCC 19606]
[formally known as *Salmonella choleraesuis*]

Antibiotic-Resistant [Strains of] Bacteria [Antibiotic-Resistant Bactericidal Activity]-

Staphylococcus aureus MRSA [ATCC 33592] *Enterococcus faecalis* VRE [ATCC 51575]

***VIRUCIDAL** in the presence of 5% serum load and 1 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

Viruses [*Virucidal Activity]:

Herpes Simplex Virus, Type 1 [HSV-1] Herpes Simplex Virus, Type 2 [HSV-2] Human Coronavirus [ATCC VR-740]
[ATCC VR-733] [ATCC VR-734] Rhinovirus Type 14 [ATCC VR-284]
Influenza A/Hong Kong [ATCC VR-544] Poliovirus Type 1, Strain Brunhilde Rotavirus WA
[ATCC VR-1000]

Kills HIV-1 [AIDS virus] [HTLV-III_B] when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.
Kills HBV & HCV when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

Kills Norovirus [Feline Calicivirus as the surrogate] [VR-782] when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

[Veterinary Viruses:]

Avian Influenza A, [H3N2] [ATCC VR-2072] Feline Calicivirus, Strain F9 [ATCC VR-782]

TUBERCULOCIDAL in the presence of 5% serum load and 5 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

Mycobacterium bovis [Tb] [OT 451C150]

FUNGICIDAL This product is fungicidal against the pathogenic fungi, *Trichophyton mentagrophytes* [ATCC 9533] [the fungus which causes Athlete's Foot] in the presence of 5% serum load and 10 minute contact time at 20°C (68°F) when used as directed on hard surfaces found in bathroom, shower stalls, locker rooms, or other clean, non-porous, hard surfaces commonly contacted by bare feet.

Trichophyton mentagrophytes [ATCC 9533]

Non-Food Contact Surface Sanitizer for hard non-porous surfaces - in the presence of 5% serum load and 30 second contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

Klebsiella pneumoniae [ATCC 4352] *Salmonella enterica* [ATCC 10708] *Pseudomonas aeruginosa* [ATCC 15442]
Staphylococcus aureus [ATCC 6538] [formally known as *Salmonella choleraesuis*] *Enterococcus faecalis* VRE [ATCC 51575]
Escherichia coli O157:H7 [ATCC 35150] *Staphylococcus aureus* MRSA [ATCC 33592]

Malodor[s] [Activity] [Counteractancy] - eliminates [destroys] odors and odor-causing bacteria in restroom areas, behind and under sinks and counters, and storage areas [and other places where bacterial growth can cause malodors].

EPA Reg. No. 74559-3

ACCEPTED
with COMMENTS
EPA Letter Dated:

JAN 21 2010

Under the Federal Insecticide,
Fungicide, and Rodenticide Act as
amended, for the pesticide,
registered under EPA Reg. No. 74559-3