

74559-3

6/7/2010

10/13



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

WASHINGTON, D.C. 20460

JUN - 7 2010

OFFICE OF  
PREVENTION, PESTICIDES  
AND TOXIC SUBSTANCES

Mr. David Swain  
Regulatory Agent for,  
Virox Technologies, Inc.  
2770 Coventry Road  
Oakville, Ontario L6H 6R1 Canada

Mail to: Attn: David Swain  
Scientific & Regulatory Consultants, Inc.  
P. O. Box 1014  
Columbia City, IN 46725

Subject: Accel Tb Wipes  
EPA Registration Number 74559-3  
Your Notification Dated May 10<sup>th</sup>, 2010  
EPA Received Date May 11<sup>th</sup>, 2010

The notification referred to above, submitted in connection with registration under the Federal Insecticide, Fungicide, and Rodenticide Act, FIFRA, as amended, to revise the storage and disposal section as per PR Notice 2007-4, is acceptable.

This notification has been made part of the permanent file.

If you have questions concerning this letter, please contact Karen M. Leavy at (703)-308-6237.

Sincerely,

A handwritten signature in dark ink, which appears to read "M. Swindell", is positioned above the printed name.

Marshall Swindell  
Product Manager 33  
Regulatory Management Branch I  
Antimicrobials Division(7510P)



Please read instructions on reverse before completing form.

Form Approved. OMB No. 2070-0060

	United States	<input type="checkbox"/> Registration <input type="checkbox"/> Amendment <input checked="" type="checkbox"/> Other	OPP Identifier Number
	Environmental Protection Agency		
	Washington, DC 20460		

Application for Pesticide - Section I

1. Company/Product Number 74559-3	2. EPA Product Manager M. Swindell	3. Proposed Classification <input checked="" type="checkbox"/> None <input type="checkbox"/> Restricted
4. Company/Product (Name) Accel Tb Wipes	PM# 33	
5. Name and Address of Applicant (Include ZIP Code) Virox Technologies, Inc. 2770 Coventry Road Oakville, ON L6H 6R1 Canada  <input type="checkbox"/> Check if this is a new address		6. Expedited Review. In accordance with FIFRA Section 3(c)(3) (b)(i), my product is similar or identical in composition and labeling to: EPA Reg. No. _____  Product Name _____

Section - II

<input type="checkbox"/> Amendment - Explain below.	<input type="checkbox"/> Final printed labels in response to Agency letter dated _____
<input type="checkbox"/> Resubmission in response to Agency letter dated _____	<input type="checkbox"/> "Me Too" Application.
<input checked="" type="checkbox"/> Notification - Explain below.	<input type="checkbox"/> Other - Explain below.

**Explanation:** Use additional page(s) if necessary. (For section I and Section II.)

Notification of new container disposal per PR Notice 2007-4. This notification is consistent with the provisions of PR Notice 98-10 and EPA regulations at 40 CFR 152.46, and no other changes have been made to the labeling or confidential statement of formula of this product. I understand that it is a violation of 18 U.S.C. Sec. 1001 to willfully make any false statement to EPA. I further understand that if this notification is not consistent with the terms of PR Notice 98-10 and 40 CFR 152.46, this product may be in violation of FIFRA and I may be subject to enforcement action and penalties under sections 12 and 14 of FIFRA. Please contact David Swain at 260-244-6273 (fax) with any questions.

Section - III

1. Material This Product Will Be Packaged In:				2. Type of Container	
Child-Resistant Packaging <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	Unit Packaging <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	Water Soluble Packaging <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No		<input checked="" type="checkbox"/> Metal	
				<input type="checkbox"/> Plastic	
				<input type="checkbox"/> Glass	
				<input type="checkbox"/> Paper	
				<input type="checkbox"/> Other (Specify) _____	
* Certification must be submitted					
3. Location of Net Contents Information <input checked="" type="checkbox"/> Label <input type="checkbox"/> Container		4. Size(s) Retail Container Various		5. Location of Label Directions <input type="checkbox"/> on label	
6. Manner in Which Label is Affixed to Product <input checked="" type="checkbox"/> Lithograph <input type="checkbox"/> Paper glued <input type="checkbox"/> Stenciled				<input type="checkbox"/> Other _____	

Section - IV

1. Contact Point (Complete items directly below for identification of individual to be contacted, if necessary, to process this application.)		
Name David Swain	Title Agent for Virox Technologies, Inc.	Telephone No. (Include Area Code) 250-244-6270
<b>Certification</b> I certify that the statements I have made on this form and all attachments thereto are true, accurate and complete. I acknowledge that any knowingly false or misleading statement may be punishable by fine or imprisonment or both under applicable law.		8. Date Application Received (Stamped)
2. Signature 	3. Title Agent for Virox Technologies, Inc.	
4. Typed Name David Swain	5. Date 5/10/2010	

35

30613

**SRC**  
**Scientific & Regulatory**  
Consultants, Inc.

May 10, 2010

Mr. Marshall Swindell, PM-33  
Document Processing Desk (NOTIF)  
Office of Pesticide Programs (7504P)  
U.S. Environmental Protection Agency  
Room S-4900, 4th Floor  
One Potomac Yard  
2777 S Crystal Drive  
Arlington, VA 22202

**SUBJECT: Notification (Container and Containment Regulations) per PR Notice 2007-4**

Accel Tb Wipes  
EPA Reg. No. 74559-3

Dear Mr. Swindell:

On behalf of our client, Virox Technologies, Inc., enclosed is a notification revising the storage and disposal section to comply with the Final Pesticide Container and Containment Regulations. This is a Notification per PR Notice 2007-4; therefore, no PRIA fee is required.

The new statements, as specified in the guidance, have been included as required text in the attached labeling. The new statements are highlighted in red text on the enclosed label. One other small addition has been made to comply with a requirement from California DPR. They required the addition of the phrase "hard, nonporous surfaces" to describe the surfaces behind sinks and counters in the malodors sections.

Accel Tb Wipes is a 100% repackaging of a registered end-use product. The product is identical in formulation to Oxivir Wipes (EPA Reg. No. 70627-60). The registrant of Oxivir Wipes, JohnsonDiversey, Inc., has already submitted these changes to federal EPA for review.

This notification, as stipulated by the guidance, is consistent with the provisions of PR Notice 98-10 and EPA regulations at 40 CFR 152.46, and no other changes have been made to the labeling or confidential statement of formula of this product. I understand that it is a violation of 18 U.S.C. Sec. 1001 to willfully make any false statement to EPA. I further understand that if this notification is not consistent with the terms of PR Notice 98-10 and 40 CFR 152.46, this product may be in violation of FIFRA and I may be subject to enforcement action, and penalties under sections 12 and 14 of FIFRA.


We look forward to a timely review of this submission. If you need any clarification or issues arise during the review of this information, please contact me immediately at 260-244-6270 for resolution.

Sincerely,



David Swain  
Agent for Virox Technologies, Inc.

Cc: Babak Givchchi



{Administrative note: any language which appears between { } and is italicized will **NOT** appear on final printed labels.  
Any language that appears between [ ] is considered interchangeable language which may appear on final printed labeling}

**Accel Tb Wipes**

[Ready to Use] [Easy to Use] [Towelette] [Wipe]  
One-Step Disinfectant and Deodorizing Wipes  
Sanitizer • Sanitizing Wipes • Cleaner • Disinfectant • Deodorizer  
\*Virucidal • Bactericidal • [Fungicidal] • Tuberculocidal  
Meets OSHA Bloodborne Pathogen Standard for HBV, HCV & HIV-1  
Kills Norovirus • Kills Rhinovirus [Type 14]  
Use on non-food contact surfaces, including: [Grocery Cart Handles] [Cash Registers]  
[Telephones] [Doorknobs]  
Just pull, wipe & toss! • Eliminates Buckets • No need to Spray!  
Pre-moistened • Pre-measured • Always Ready-to-Use  
No measuring or mixing required

[Hydrogen Peroxide]

For [Commercial] [Industrial] [&] [Institutional] Use [Only]

**ACTIVE INGREDIENT:**

Hydrogen Peroxide.....0.5%

**INERT [OTHER] INGREDIENTS:**.....99.5%


**TOTAL:**.....100.0%

**KEEP OUT OF REACH OF CHILDREN**

[See reference sheet [enclosed in each case] for [a complete list of organisms] [additional [features] [claims] [directions for use]] [for this product]].

**Net Contents:**

- 160 sheets (6" X 6.75") (15.24 cm X 17.14 cm)
- 160 sheets (8" X 9") (20.32 cm X 22.86 cm)
- 10 sheets (7.48" x 10") (19.0 cm x 25.4 cm) {for re-sealable package}



60 of 13

{Administrative note: any language which appears between { } and is italicized will **NOT** appear on final printed labels. Any language that appears between [ ] is considered interchangeable language which may appear on final printed labeling}

**{MARKETING CLAIMS}:**

A disinfectant cloth for use on telephone surfaces and other telecommunications equipment surfaces  
A Hydrogen Peroxide solution impregnated in a wiping cloth  
Aids in reducing [or in the reduction of] cross contamination between treated surfaces  
Always ready to use  
An effective counter cleaner  
An effective way to clean, deodorize against odor-causing organisms and disinfect hard, non-porous surfaces  
An effective way to disinfect hard, non-porous environmental surfaces  
Antibacterial  
Antibacterial [surface] Cleaner [cloth]  
Antibacterial Cell phone Wipe  
Antibacterial Surface Cleaner [Cloth]  
Anti-microbial cleaner [cloth]  
Anti-viral\* [Anti-Virus\*] [Surface] Cleaner  
Assures proper dilution strength for product effectiveness  
Bactericidal Cleaner  
Bactericidal  
Built-in Scrubbing Fibers  
Clean & disinfect any room, any time  
Clean & disinfect hard non-porous surfaces  
Clean & disinfect hard non-porous surfaces easily  
Cleans without scratching surface[s]  
Cleaner  
Cleaner and disinfectant towel [or cloth]  
Clean[s] and deodorize[s]  
Cleans and disinfects  
Cleans and disinfects [cell phone surfaces] [telephone equipment] [telephone equipment surfaces] [dictating equipment surfaces] [telecommunications equipment surfaces] [devices]  
Cleans and disinfects common hard, non-porous surfaces  
Cleans effectively  
Clean hard-to-remove stains  
Cleans, deodorizes and disinfects  
Cloth  
Convenient  
Convenient, quick & easy cleaning all around [the office] [business] [workplace]  
Cuts grease  
Cuts the toughest grime and grease without scratching surfaces  
Disinfectant [disinfecting] surface cleaner [cloth]  
Disinfects against *Staphylococcus aureus*, *Salmonella enterica*, and *Pseudomonas aeruginosa*  
Disinfects environmental surfaces  
Disinfects hard, non-porous environmental surfaces  
Disinfects headsets and telephones  
Disposable  
Each pre-moistened towelette [disposable wipe] cleans, deodorizes and kills listed bacteria all around you [office] [business] [workplace]  
Easy flip-top lid  
Effective against flu virus [Influenza A/Hong Kong] [Hong Kong Flu]  
Effective daily cleaner  
Effective in 1 minute against *Staphylococcus aureus*, *Salmonella enterica*, and *Pseudomonas aeruginosa*  
Effectively disinfects hard, non-porous surfaces  
Effectively disinfects in the presence of organic soil [5% blood serum]  
Everyday, every room disinfectant  
Exhibits effective virucidal activity: kills Influenza A/Hong Kong and Herpes Simplex Type 2 [cause of Genital Herpes]  
For hard, non-porous surfaces  
For use in critical care areas [where control of the hazards of cross-contamination between treated surfaces is of prime importance]  
For your toughest bathroom messes

{Administrative note: any language which appears between { } and is italicized will **NOT** appear on final printed labels. Any language that appears between [ ] is considered interchangeable language which may appear on final printed labeling}

Grease Cutting Formula

Great for grease cutting

Great for on-the-go cleaning

Great on-the-go disinfectant [cloth]

Hard non-porous surface [cleaning &] disinfecting made easy

Heavy-duty cleaner [or cleaning] [cloth]

Kills Avian Influenza A virus on pre-cleaned environmental surfaces

Kills *Staphylococcus aureus*, *Salmonella enterica*, and *Pseudomonas aeruginosa*

Kills the flu virus Influenza A

Kills Norovirus [Norwalk virus]

Kills Human Coronavirus

Kills antibiotic-resistant organisms

Kills HIV-1 [AIDS virus] [HTLV-III<sub>B</sub>]

Kills HCV

Kills HBV

Kills Rhinovirus [Type 14]

No mixing

No streaking

Non-woven disposable cloth

On-the-go disinfectant

Accel Tb Wipes control the spread of listed bacteria on treated hard, non-porous environmental surfaces

Accel Tb Wipes help control listed bacteria on treated surfaces while keeping your office [or workstation or workplace or business] clean and dust-free

Pre-measured

Pre-moistened [towelettes]

Pre-soaked [cloth or towelette]

Quick and easy cleaning everywhere

Quick drying

Removes tough stains

Restroom Disinfectant Wipes

Restroom Wipes [for hard, non-porous surfaces]

Scrubbing surface cleaner

Stay-moist package

Strong towelettes [or wipes]

Telephone [cleaner [and] disinfectant] wipes

These wipes are a convenient way to disinfect surfaces all around your [business] [office] [workplace]

Towel

Use on hard, non-porous surfaces only. [Test wipe on a small inconspicuous area first]

Use to disinfectant shopping [or grocery] cart handles

Versatile disinfectant [or disinfecting] cleaner

Virucidal\* [or] Kills viruses\* \*Influenza A/Hong Kong [Hong Kong flu], Herpes Simplex Type 2, HCV, HBV, HIV-1

Virus killing cleaner\*

Wipe down grocery [or shopping] cart handles

#### [FEATURES, CLAIMS & USES:]

[Notice to User: Accel TB may not be suitable for all items or surfaces. Prolonged exposure to copper, brass, anodized aluminum and some coatings should be avoided. If you are uncertain as to the material composition of your item, confirm with the manufacturer before proceeding [or visit the Compatibility section on [www.virox.com](http://www.virox.com) for a listing of compatible materials.]]

#### **For Use in Health Care, Veterinary, Animal Housing, Dental, [Life Science] Laboratory, Medical, Industrial & Institutional Facilities, Food Processing Establishments, and Hospitals**

This product [These] is [are] [is a] non-woven, disposable [towelette[s]] [wipe[s]], pre-saturated with Hydrogen Peroxide [based] disinfectant cleaner.

{Administrative note: any language which appears between { } and is italicized will **NOT** appear on final printed labels. Any language that appears between [ ] is considered interchangeable language which may appear on final printed labeling}

This product is a one-step [hospital-use] germicidal [disinfectant] cleaner and deodorant [odor-counteractant] [odor, neutralizer] designed for general cleaning, [and] disinfecting, [and] [deodorizing] [of] hard, non-porous environmental surfaces. It cleans quickly by removing dirt, grime, body oils and other common soils found in hospitals, nursing homes, [acute] [and] [long-term] care facilities, office buildings, schools and colleges [universities], [day care centers], [medical] [dental] offices, funeral homes, transportation facilities, train [depots] [stations] bus [depots] [stations], airports, casinos, cruise ships, correctional facilities, athletic facilities, athletic equipment, playgrounds, locker rooms, health clubs, exercise equipment, barber/beauty [salons] [shops], veterinary clinics, pet shops, [equine farms], animal life science laboratories, hotels, motels, public areas and restrooms, homes cleaned by [home cleaning services,] [professional cleaning services,] [medical assistance services,] and office buildings.

This product's [non-dulling] [advanced] [quick cleaning] [convenient] [easy-to-use] [highly effective] formula eliminates the time and labor normally required for rinsing. It is designed for use on [Use daily on] [the following] [hard, non-porous environmental [inanimate] surfaces]: vinyl, painted surfaces, plastic [surfaces], glass, mirrors, glazed ceramic, glazed porcelain, chrome, stainless steel, aluminum, laminated surfaces and baked enamel surfaces associated with floors, walls ceilings, tables, chairs, non-food contact countertops, telephones, fixtures, athletic equipment, exercise equipment, glazed tile, exterior surfaces of toilets, exterior surfaces of urinals, sinks, sink basins, faucets, tubs, chrome, stainless steel, shower rooms and locker room areas – any washable [non-food contact] surface where disinfection is required. This product's [non-dulling] [advanced] [quick cleaning] [convenient] [easy-to-use] [highly effective] formula eliminates the time and labor normally required for rinsing. Do not use on glasses, dishes and utensils. Not recommended for use on brass, copper or marble.

#### **[Hospitals/Healthcare Facilities/Dental Offices/Life Science Laboratories]**

This product is one-step [hospital-use] [laboratory-use] germicidal [disinfectant] cleaner and deodorant [odor-counteractant] [odor neutralizer] towelette designed for general cleaning, [and] disinfecting, [deodorizing] [of] hard, non-porous environmental surfaces. This product cleans quickly by removing dirt, grime, food residue, dead skin, blood, body soils and other organic matter commonly found in Dental Clinics, [Hospitals], [in Health Care Facilities], [laboratory settings]. It also eliminates odors, leaving [restroom] surfaces smelling clean and fresh. Use where odors are a problem.

This product cleans, disinfects and deodorizes hard, non-porous environmental [hospital] [medical] [dental] [laboratory] surfaces in one step [with no rinsing required]. Its non-abrasive formula is designed for use on [Use daily on] [the following] [hard, non-porous environmental surfaces]: [vinyl], [Naugahyde], painted surfaces, plastic [surfaces], glazed ceramic, glazed porcelain, chrome, stainless steel, laminated surfaces and baked enamel surfaces associated with floors, walls, ceilings, tables, chairs, non-food contact countertops, telephones, fixtures, glazed tile, toilet exteriors, urinal exteriors, sinks found in [health care facilities], [hospitals], [dental offices], [patient rooms], [operating rooms], [ICU areas], [shower rooms], and [locker rooms] [fume hoods and animal cages [found in laboratory settings]]. [It can be used to pre-clean and disinfect hospital items: wheelchairs, [hospital] [patient] bedrails and linings [patient bed linings] [vinyl mattress covers], wash basins, bed pans, non-critical medical equipment surfaces – almost any washable non-food contact surface [where disinfection is required].

#### **[[Bathroom] [Restroom] Disinfectant Cleaner:]**

This product is a one-step disinfectant cleaner and deodorant [odor-counteractant] [odor-neutralizer] designed for general cleaning, [and] disinfecting, [deodorizing] on hard, non-porous environmental surfaces.

This product cleans, disinfects and deodorizes surfaces by killing odor-causing microorganisms. Its [non-dulling] [advanced] [quick cleaning] [convenient] [easy-to-use] [highly effective] formula is designed for use on [Use it daily] on [hard, non-porous environmental] restroom surfaces]: glazed ceramic [restroom] tiles, glazed porcelain, chrome, stainless steel and plastic surfaces associated with floors, walls, fixtures, exterior [surfaces] of toilets, exterior [surfaces] of urinals, sinks, shower rooms and locker rooms. Not recommended for use on brass, copper or marble.

#### **[Animal Housing Facilities]**

This product cleans by removing dirt, grime, blood, urine, fecal matter and other common soils found in animal housing facilities, grooming facilities, farms, kennels, pet stores, veterinary clinics, laboratories or other small animal facilities. It [also] eliminates odors leaving surfaces smelling clean and fresh.



{Administrative note: any language which appears between { } and is italicized will **NOT** appear on final printed labels. Any language that appears between [ ] is considered interchangeable language which may appear on final printed labeling}

This product cleans, disinfects and deodorizes hard, non-porous environmental surfaces in one step. Its [non-dulling] [advanced] [quick cleaning] [convenient] [easy-to-use] [highly effective] formula is designed for use on [Use it daily] [Use daily to clean and disinfect] [hard, non-porous surfaces:] plated or stainless steel, chrome, glazed porcelain, glazed tile, laminated surfaces, [associated with floors, walls, non-food contact countertops, cages, kennels, animal equipment] found in [barns, pens and stalls] animal housing facilities.

#### [CLAIMS:]

[When used as directed, this product [is] highly effective against a wide variety [broad-spectrum] of pathogenic microorganisms [including bacteria, tuberculosis, antibiotic resistant bacteria, viruses and fungi. See reference sheet [enclosed in each case] for a complete list of organism].

**HOSPITAL DISINFECTANT: BACTERICIDAL** in the presence of 5% serum load and 1 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

#### Bactericidal Activity:

*Pseudomonas aeruginosa* [ATCC 15442]      *Staphylococcus aureus* [ATCC 6538]      *Escherichia coli* O157:H7 [ATCC 35150]  
*Salmonella enterica* [ATCC 10708]      *Acinetobacter baumannii* [ATCC 19606]  
 [formally know as *Salmonella choleraesuis*]

#### Antibiotic-Resistant [Strains of] Bacteria [Antibiotic-Resistant Bactericidal Activity]-

*Staphylococcus aureus* MRSA [ATCC 33592]      *Enterococcus faecalis* VRE [ATCC 51575]

**\*VIRUCIDAL** in the presence of 5% serum load and 1 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

#### Viruses [\*Virucidal Activity]:

Herpes Simplex Virus, Type 1 [HSV-1]      Herpes Simplex Virus, Type 2 [HSV-2]      Human Coronavirus [ATCC VR-740]  
 [ATCC VR-733]      [ATCC VR-734]      Rhinovirus, Type 14 [ATCC VR-284]  
*Influenza A/Hong Kong* [ATCC VR-544]      Poliovirus Type 1, Strain Brunhilde      Rotavirus WA  
 [ATCC VR-1000]

Kills HIV-1 [AIDS virus] [HTLV-III<sub>B</sub>] when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.  
 Kills HBV & HCV when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

Kills Norovirus [Feline Calicivirus as the surrogate] [VR-782] when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

#### [Veterinary Viruses:]

Avian Influenza A, [H3N2] [ATCC VR-2072]      Feline Calicivirus, Strain F9 [ATCC VR-782]

**TUBERCULOCIDAL** in the presence of 5% serum load and 5 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

*Mycobacterium bovis* [Tb] [OT 451C150]

**[FUNGICIDAL]** This product is fungicidal against the pathogenic fungi, *Trichophyton mentagrophytes* [ATCC 9533] [the fungus which causes Athlete's Foot] in the presence of 5% serum load and 10 minute contact time at 20°C (68°F) when used as directed on hard surfaces found in bathroom, shower stalls, locker rooms, or other clean, non-porous, hard surfaces commonly contacted by bare feet.]

[*Trichophyton mentagrophytes* [ATCC 9533]]

**Non-Food Contact Surface Sanitizer for hard non-porous surfaces** - in the presence of 5% serum load and 30 second contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

*Klebsiella pneumoniae* [ATCC 4352]      *Salmonella enterica* [ATCC 10708]      *Pseudomonas aeruginosa* [ATCC 15442]  
*Staphylococcus aureus* [ATCC 6538]      [formally know as *Salmonella choleraesuis*]      *Enterococcus faecalis* VRE [ATCC 51575]  
*Escherichia coli* O157:H7 [ATCC 35150]      *Staphylococcus aureus* MRSA [ATCC 33592]

100f13

{Administrative note: any language which appears between { } and is italicized will **NOT** appear on final printed labels. Any language that appears between [ ] is considered interchangeable language which may appear on final printed labeling}

**Malodor[s] [Activity] [Counteractancy]** – eliminates [destroys] odors and odor-causing bacteria on hard, non-porous surfaces in restroom areas, behind and under sinks and counters, and storage areas [and other places where bacterial growth can cause malodors].

# **DIRECTIONS FOR USE:**

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Not for personal cleansing. This is not a baby wipe. Do not flush in toilet.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

## **For Use as [Daily] [One-Step] Cleaner/Disinfectant:**

1. Pre-clean heavily soiled areas
2. Pull towelette from dispenser [canister] and wipe hard, non-porous environmental surfaces.
3. All surfaces must remain visibly wet for 1 minute. Use a 5 minute contact time for Tuberculosis. [Use a 10 minute contact time for fungi.]
4. Allow to air dry [Wipe surface dry] [or rinse] [with potable water if necessary].

Note: All food contact surfaces must be rinsed with potable water. Do not use on glassware, utensils, or dishes.

## **For Use as [Daily] [Bathroom] [Restroom] [Shower Rooms] [Locker Rooms] [Athletic Equipment] [One-Step] Cleaner/Disinfectant:**

1. Pre-clean heavily soiled areas
2. Pull towelette from dispenser [canister] and wipe hard, non-porous environmental surfaces.
3. To disinfect, all surfaces must remain wet for one [1] minutes. Use a 5 minute contact time for Tuberculosis. [Use a 10 minute contact time for fungi.]
4. Wipe surfaces dry [or rinse] [allow to air dry].

### **\*KILLS HIV AND HCV AND HBV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS**

**PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings [Hospitals, Nursing Homes] and other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/ body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS) Hepatitis B Virus (HBV), or Hepatitis C Virus (HCV).

### **SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV OR HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:**

**Personal Protection:** Disposable latex or vinyl gloves, gowns, face masks, or eye coverings as appropriate, must be worn during all cleaning of body fluids, blood, and decontamination procedures.

**Cleaning Procedures:** Blood/body fluids must be thoroughly cleaned from surfaces and objects before application of this product.

**Contact Time:** Allow surface to remain wet for 1 minute to kill HIV-1, HBV and HCV. Use a 5 minute contact time for Tuberculosis and 10 minute contact for fungi.

**Disposal of Infectious Material:** Blood/body fluids should be autoclaved and disposed of according to Federal, State, and local regulations for infectious waste disposal.

## **For Use as a Tuberculocide:**

1. Pre-clean heavily soiled areas.
2. Pull towelette from dispenser [canister] and wipe hard, non-porous environmental surfaces.
3. Allow the surface to remain wet for 5 minutes. Use additional towelettes if needed to ensure the correct wet contact time.
4. Wipe surfaces dry [or rinse] [or allow to air dry].
5. Change towelette after each use [application].

{Administrative note: any language which appears between { } and is italicized will **NOT** appear on final printed labels. Any language that appears between [ ] is considered interchangeable language which may appear on final printed labeling}

#### **To Kill Fungi:**

1. Pre-clean heavily soiled areas.
2. Pull towelette from dispenser [canister] and wipe hard, non-porous environmental surfaces.
3. All surfaces must remain visibly wet for 10 minutes. Use additional towelettes if needed to ensure the correct wet contact time.
4. Wipe surfaces dry [or rinse] [or allow to air dry].

#### **For use as a Non-Food Contact Sanitizer:**

1. Pre-clean heavily soiled areas.
2. Pull towelette from dispenser [canister] and wipe hard, non-porous environmental surfaces.
3. Allow the surface to remain wet for 30 seconds.
4. Wipe surfaces dry [or rinse] [or allow to air dry].
5. Change towelette after each use [application].

#### **For Use as a Cleaner and/or Deodorizer:**

Pull towelette from dispenser [canister] and wipe hard, non-porous environmental surfaces. Allow to air dry [or rinse [with potable water if necessary]].

#### **Animal Housing Facilities [including Veterinary Clinics, [Equine] Farms, Livestock, Swine and Poultry Houses]:**

**For cleaning and disinfecting the following hard, non-porous surfaces: feeding and watering equipment, utensils, instruments, cages, kennels, stables, catteries, etc.**

1. Remove all animals and feeds from premises, animal transportation vehicles, crates, etc.
2. Remove all litter, droppings and manure from floors, walls and surfaces of facilities [including pens, chutes, and barns] occupied or traversed by animals.
3. Empty all troughs, racks, and other feeding or watering appliances.
4. Thoroughly clean all surfaces with soap and detergent and rinse with potable water. Pull towelette from dispenser [canister] and wipe hard, non-porous environmental surfaces.
5. All surfaces must remain visibly wet for 1 minute.
6. Allow to air dry [or rinse [with potable water if necessary]].
7. Ventilate buildings and other closed spaces.
8. Do not house animals or employ equipment until treatment has been absorbed, set or dried.

NOTE: Not for use on handling and restraining equipment such as leashes, muzzles, halters or ropes, as well as forks, shovels, and scrapers used for removing litter and manure.

{Optional directions for refill pouch ONLY}

#### **Refill Pouch**

1. Open empty canister lid.
2. Open refill pouch [by cutting the top of the pouch] and remove saturated towelettes.
3. Place saturated towelettes in empty canister [invert the contents of the pouch into empty canister].
4. Replace canister lid.
5. [Remove the new lot number / [expiry date] from pouch and place it on top of the old lot number / [expiry date] sticker on the canister]

{Optional directions for refill pouch ONLY}

#### **Re-sealable soft pack**

1. Remove the seal (from the place shown on the label)
2. Remove one sheet at a time.
3. Close the seal (to prevent drying between uses)

{Administrative note: any language which appears between { } and is italicized will **NOT** appear on final printed labels. Any language that appears between [ ] is considered interchangeable language which may appear on final printed labeling}

## STORAGE AND DISPOSAL:

Do not contaminate water, food or feed by storage and disposal.

**Storage:** Store in original container in areas inaccessible to small children. Keep securely closed.

**Container Disposal:** Non-refillable container. Do not reuse or refill this container. Rinse thoroughly. Offer for recycling if available or [wrap and] discard in trash.

*{For product supplied in refillable containers or refill package :}*

**Towelette Storage/Disposal:** Dispose of [wipe] [towelette] [cloth] in trash after use. Do not flush. Reseal container label between uses to retain moisture. Refillable container. Refill this container with refill wipe only. Do not reuse this container for any other purpose. Discard in trash.

Cleaning the container before final disposal is the responsibility of the person disposing of the container. Cleaning before refilling is the responsibility of the refiller.

*{For product supplied in resealable non-refillable packages :}*

**Towelette Storage/Disposal:** Non-refillable container. Do not reuse or refill this container. Dispose of [wipe] [towelette] [cloth] in trash after use. Do not flush. Reseal container label between uses to retain moisture. Do not reuse empty [Pouch] [Container] [Package]. Discard empty [Pouch] [Container] [Package] in trash.

## ENVIRONMENTAL HAZARDS *{for labels with veterinary use patterns}*

This product is toxic to birds, fish and aquatic invertebrates. Caution should be used when applying indoors because pets may be at risk.

Read the M.S.D.S. Before Using the Product.

M.S.D.S. is available on line at [www.SciCan.com](http://www.SciCan.com) [www.virox.com]

[Si no puede leer en ingles, pregunte a su supervisor sobre las instrucciones de uso apropiadas antes de trabajar con este producto.]

EPA Reg. No. 74559-3

EPA Est. No. 74559-CAN-1

© 2008 Virox Technologies, Inc., 2770 Coventry Road, Oakville, ON L6H 6R1 Canada All Rights Reserved

AHP Logo\*\*

\*\* AHP technology and designs are trademarks of Virox Technologies, Inc. [used under license by SciCan Ltd.]

130413

{Administrative note: any language which appears between { } and is italicized will **NOT** appear on final printed labels. Any language that appears between [ ] is considered interchangeable language which may appear on final printed labeling}

**Virox Technologies, Inc.**

**Accel Tb Wipes  
Reference Sheet**

**HOSPITAL DISINFECTANT: BACTERICIDAL** in the presence of 5% serum load and 1 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

**Bactericidal Activity:**

<i>Pseudomonas aeruginosa</i> [ATCC 15442]	<i>Staphylococcus aureus</i> [ATCC 6538]	<i>Escherichia coli</i> O157:H7 [ATCC 35150]
<i>Salmonella enterica</i> [ATCC 10708]	<i>Acinetobacter baumannii</i> [ATCC 19606]	
[formally know as <i>Salmonella choleraesuis</i> ]		

**Antibiotic-Resistant [Strains of] Bacteria [Antibiotic-Resistant Bactericidal Activity]-**

*Staphylococcus aureus* MRSA [ATCC 33592]      *Enterococcus faecalis* VRE [ATCC 51575]

**\*VIRUCIDAL** in the presence of 5% serum load and 1 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

**Viruses [\*Virucidal Activity]:**

Herpes Simplex Virus, Type 1 [HSV-1] [ATCC VR-733]	Herpes Simplex Virus, Type 2 [HSV-2] [ATCC VR-734]	Human Coronavirus [ATCC VR-740]
Influenza A/Hong Kong [ATCC VR-544]	Poliovirus Type 1, Strain Brunhilde [ATCC VR-1000]	Rhinovirus Type 14 [ATCC VR-284] Rotavirus WA

Kills HIV-1 [AIDS virus] [HTLV-III<sub>B</sub>] when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.  
Kills HBV & HCV when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

Kills Norovirus [Feline Calicivirus as the surrogate] [VR-782] when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

**[Veterinary Viruses:]**

Avian Influenza A, [H3N2] [ATCC VR-2072]      Feline Calicivirus, Strain F9 [ATCC VR-782]

**TUBERCULOCIDAL** in the presence of 5% serum load and 5 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

*Mycobacterium bovis* [Tb] [OT 451C150]

**FUNGICIDAL** This product is fungicidal against the pathogenic fungi, *Trichophyton mentagrophytes* [ATCC 9533] [the fungus which causes Athlete's Foot] in the presence of 5% serum load and 10 minute contact time at 20°C (68°F) when used as directed on hard surfaces found in bathroom, shower stalls, locker rooms, or other clean, non-porous, hard surfaces commonly contacted by bare feet.

*Trichophyton mentagrophytes* [ATCC 9533]

**Non-Food Contact Surface Sanitizer for hard non-porous surfaces** - in the presence of 5% serum load and 30 second contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

<i>Klebsiella pneumoniae</i> [ATCC 4352]	<i>Salmonella enterica</i> [ATCC 10708]	<i>Pseudomonas aeruginosa</i> [ATCC 15442]
<i>Staphylococcus aureus</i> [ATCC 6538]	[formally know as <i>Salmonella choleraesuis</i> ]	<i>Enterococcus faecalis</i> VRE [ATCC 51575]
<i>Escherichia coli</i> O157:H7 [ATCC 35150]	<i>Staphylococcus aureus</i> MRSA [ATCC 33592]	

**Malodor[s] [Activity] [Counteractancy]** – eliminates [destroys] odors and odor-causing bacteria on hard, non-porous surfaces in restroom areas, behind and under sinks and counters, and storage areas [and other places where bacterial growth can cause malodors].

EPA Reg. No. 74559-3