



**UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
WASHINGTON, DC 20460**

**OFFICE OF CHEMICAL SAFETY
AND POLLUTION PREVENTION**

April 23, 2021

Lok Chum
Regulatory Affairs Manager
Virox Technologies Inc.
2770 Coventry Road
Oakville Ontario, L6H 6R1
Canada

Subject: PRIA Label Amendment – Acceptable Addition to Add SARS-COV-2 Claims
Product Name: OXY-1 Wipes
EPA Registration Number: 74559-10
Received Date: December 14, 2020
OPP Case Number: 00216929

Dear Mr. Chum:

The amended label referred to above, submitted in connection with registration under the Federal Insecticide, Fungicide and Rodenticide Act, as amended, is acceptable. This approval does not affect any conditions that were previously imposed on this registration. You continue to be subject to existing conditions on your registration and any deadlines connected with them.

A stamped copy of your labeling is enclosed for your records. This labeling supersedes all previously accepted labeling. You must submit one copy of the final printed labeling before you release the product for shipment with the new labeling. In accordance with 40 CFR 152.130(c), you may distribute or sell this product under the previously approved labeling for 18 months from the date of this letter. After 18 months, you may only distribute or sell this product if it bears this new revised labeling or subsequently approved labeling. "To distribute or sell" is defined under FIFRA section 2(gg) and its implementing regulation at 40 CFR 152.3.

Should you wish to add/retain a reference to the company's website on your label, then please be aware that the website becomes labeling under the Federal Insecticide Fungicide and Rodenticide Act and is subject to review by the Agency. If the website is false or misleading, the product would be misbranded and unlawful to sell or distribute under FIFRA section 12(a)(1)(E). 40 CFR 156.10(a)(5) list examples of statements EPA may consider false or misleading. In addition, regardless of whether a website is referenced on your product's label, claims made on the website may not substantially differ from those claims approved through the registration process. Therefore, should the Agency find or if it is brought to our attention that a website contains false or misleading statements or claims substantially differing from the EPA approved registration, the website will be referred to the EPA's Office of Enforcement and Compliance.

Your release for shipment of the product constitutes acceptance of these conditions. If these conditions are not complied with, the registration will be subject to cancellation in accordance

Page 2 of 2
EPA Reg. No. 74559-10
Decision No. 00216929

with FIFRA section 6. If you have any questions, please contact Alex Horansky by phone at 703-347-0128, or via email at Horansky.alex@epa.gov.

Sincerely,

A handwritten signature in dark ink, appearing to read 'DF', with a long horizontal stroke extending to the right.


Demson Fuller, Product Manager 32
Regulatory Management Branch I
Antimicrobials Division (7510P)
Office of Pesticide Programs

Enclosure

{*Administrative note:* Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike through~~ text means it has been removed.

OXY-1 Wipes

{Alternate Brand Names: PREempt Wipes, TRIumph Wipes, Avert Wipes, PEROXIgard Wipes, Rescue; Rescue Wipes One Step Disinfectant Cleaner & Deodorizer; Intervention Farm Animal Care Disinfectant Cleaner & Deodorizer; Intervention Animal Husbandry Disinfectant Cleaner & Deodorizer; Optim 1 Wipes; Rejuvenate Wipes; Restorox Wipes; Preempt Plus Wipes}



[USA Patent #: 7,632,523 & 7,354,604 {or insert new patent # when available}]
Accelerated Hydrogen Peroxide® and AHP® are owned by Virox Technologies Inc. and used under license

[Ready to Use] [Easy to Use] [Towelette] [Multi-Surface] [Wipe]
One-Step³ [Animal Premises] [Animal Husbandry] [Hospital] Disinfectant [Cleaner] [and Deodorizing] Wipes
Sanitizer • Sanitizing Wipes • Cleaner • Disinfectant • Deodorizer • Deodorant
[General] *Virucidal • Bactericidal • [Fungicidal] • [Tuberculocidal]
[General] *Virucide • Bactericide • [Fungicide] • [Tuberculocide]
Meets surface disinfection recommendations under the OSHA Bloodborne Pathogen Standard for *HBV, *HCV & *HIV-1
Kills Norovirus • Kills Rhinovirus [Type 37] • Kills Influenza • Kills Parvovirus[◇]
Use on non-food contact surfaces, including: [Grocery Cart Handles] [Cash Registers]
[Telephones] [Doorknobs]
Just pull, wipe & toss!¹ • Eliminates Buckets • No need to Spray!
Pre-moistened • Pre-measured • Always Ready-to-Use
No [measuring] [or] [mixing] required

[Hydrogen Peroxide]

[Intended For] [Hospital,] [Health Care,] [Commercial,] [Industrial] & [Institutional] [Residential] [Animal Health] [Animal Premises] Use[±] [Only]
[[Animal Health] Premises Disinfectant Cleaner³]
[For] Residential Use
[For] Home[±] Use
For Homecare[±] Use
For Home [and] [Healthcare] Use
For Office Use

For Veterinary[±] Use
For Aquaculture[±] Use
For Barn Use

{**Administrative note:** Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike through~~ text means it has been removed.

For Farm Use
For Dental Clinic Use
[For Spa and Salon Use[±]] [For Use in Spas and Salons[±]] [and Clinics] [ideal for Pedicure Foot (Spas) (Basins) (Baths)
(pedicure foot spa equipment) [±].]
Mold and Remediation
For Use in Personal Services
For Use in Alternative Care Clinics
For Use in Respiratory Therapy Clinics
For Use in Sleep Therapy Clinics
For Use in Allied Health Clinics
[For] [use in] Life Sciences [Research Facilities]
For Use in {insert institutional use site from table 8}
Hospital Disinfectant

No Added Fragrance(s) {Note to reviewer: only for products using non-fragranced CSFs}
(Fragrance Free Formula) (Unscented) {Note to reviewer: only for products using non-fragranced CSFs}

ACTIVE INGREDIENT:

Hydrogen Peroxide..... 0.5%

INERT [OTHER] INGREDIENTS:..... 99.5%

TOTAL:..... 100.0%

KEEP OUT OF REACH OF CHILDREN

[Efficacy claims backed by third party validation.]

[See reference sheet [enclosed in each case] for [a complete list of pathogenic organisms] [additional features, claims]
[Directions for Use] [claimed for this product] [or product name]].

{NOTE to Reviewer: Due to space constraints, the Container Label may not include all of the organism claims and marketing claims approved on the Master Label. Therefore, a separate reference sheet containing additional organism claims and marketing claims is added at the end under title "Reference Sheet". This reference sheet may be inserted into the shipping case and/or distributed separately or be available at the point of sale. When a separate reference sheet is used, it will either include the use direction, precautions, etc. or contain the instructions to the user to see the product/container label for use directions, precautions, etc.}

Net Contents: (Expressed as avoirdupois pounds and ounces)

Minimum Wt.: xx pounds and yy ounces
[optional: wipe dimensions:] X x X [inches] [in.]

EPA Reg. No. 74559-10

EPA Est. No. 74559-CAN-1

Manufactured By (For):
Virox Technologies Inc.
2770 Coventry Road
Oakville, ON, L6H 6R1
1 800 387-7578 [www.virox.com] [Insert URL]

*{Administrative note: Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. Strike-through text means it has been removed.*

{WARRANTY}

{Note to the reviewer: The warranty statement below will appear on select products}

"WARRANTY: Seller warrants that this product conforms to its chemical description and is reasonably fit for the purposes stated on the label when used in accordance with label directions under normal conditions of use, but to the extent consistent with applicable law, neither this warranty nor any other warranty of MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, expressed or implied, extends to the use of this product contrary to label instructions, or under abnormal conditions, or under conditions not reasonably foreseeable to Seller, and Buyer assumes the risk of any such use."

{FRAGRANCE CLAIMS}: {Note to reviewer: only for products using fragranced CSFs}

Fruity Floral Musk [Scent]

[Helps] eliminate odors with [Febreze] [fresh scent]

Leave nothing but a [fresh scent] [light scent] [clean scent] behind

Long lasting [Febreze] fresh scent

Scent [fragrance] free

The smell [that says] [of] clean

With Febreze [fresh scents] [freshness]

Natural Scent

Natural Almond Scent

Natural Cherry Scent

Sweet Almond Scent

Dark Cherry Scent

Cherry Scent

Cherry [almond] [crush] [sparkle] [spritz]

Citrus [mandarin] [shine] [fizz] [sparkle] [renewal] [spring], [crisp] [outdoor] [fresh orange] [citrus]

Refreshing [citrus] [mandarin] [spring] [lime] [breeze] [rain] [waterfall] [sparkle] [fizz]

Sparkle & Shine, [lime] [shine] [fizz] [sparkle], [green clean], [rushing waterfall], [spring clean], [clean spring] [waterfall]

Sparkling Sky, fresh [clean] [rain] [breeze] [sky] [renewal] [herbal rain]

Yuzuteria, [citrus spring], [lemon-lime], [grapefruit surge] [rush]

{SHORT GENERAL MARKETING CLAIM}: {Note to reviewer - This text will only appear on the appropriate container.}

Acute Oral Toxicity [LD50>5,000 mg/kg] [Category IV]

Acute Inhalation Toxicity [LC50>2.08 mg/L] [Category IV]

Acute Dermal Toxicity [LD50>5,050 mg/kg] [Category IV]

An effective One-step sanitizer-cleaner⁴ wipe (towelette) for use on hard non-food contact surfaces

[A disinfectant and cleaner [heavy duty wipe] that] cleans, disinfects and deodorizes in one labor saving step³

A disinfectant cloth for use on telephone surfaces and other telecommunications equipment surfaces

A Hydrogen Peroxide solution impregnated in a wiping cloth

Aids in reducing [or in the reduction of] cross contamination between treated surfaces[±]

Always ready to use

An effective counter cleaner

An effective way to clean, deodorize against odor-causing organisms and disinfect hard, non-food contact surfaces[±]

Antibacterial [surface] Cleaner [cloth] [wipe]

Antimicrobial cleaner [cloth] [wipe]

[Surface[±]] Cleaner

Assures proper dilution strength for product effectiveness

Bactericidal Cleaner

Bactericidal

Bactericidal, *Virucidal in 1 minute[¶]

Broad-spectrum [Multi-surface[±]] [Bactericidal] [*Virucidal] [Tuberculocidal] [Fungicidal][sanitizer] [Disinfectant] [activity]

Built-in Scrubbing Fibers

Cleans [and disinfects] hard, non-porous surfaces [in kitchen[±]] [and bathrooms]

Clean & disinfect any room, any time[±]

Cleans without scratching surface[s]

{*Administrative note:* Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

Cleans [and] [deodorizes] [disinfects] hard, non-porous surfaces [all over the house[±]] [kitchens[±] and bathrooms] [easily]
Cleans [fast] [all over the house[±]] [in kitchen[±]] [and bathrooms]
Cleans [removes] [eliminates] [destroys] baked-on food
Cleans [removes] [eliminates] [destroys] burnt-on food
Cleans [tough] [grease and grime] [fast]
Cuts [cleans] kitchen[±] [and] [bathroom] [and] [household grime and soils] [fast]
Cleans [and disinfects] [hard, non-porous surfaces] {or *insert any hard non-porous surface from table 7 and 11*} [found]
[in] [hospitals] [health care facilities] [offices] [colleges] {or *insert any institutional use site from table 8 or household use site from table 12*}
Cleaner and disinfectant towelette [or cloth]
Cleans effectively
Clean you expect. Ingredients you can trust
Cleans [removes] [eliminates] [destroys] stains
Cleans [tough] grease and grime [and] [institutional] [soils] [fast] [without scratching surface]
Cleans [tough] soap scum
Cleans without lightening or discoloration
Cleans and disinfects³ non-food contact hard non-porous surfaces (all over the house) (kitchens[±] and) (bathrooms)
(restrooms)
Compatible with [equipment[±]] and [hard] [non-porous] surfaces
Controls (kitchen[±] and) odor-causing bacteria between treated surfaces
Compatible with [Suitable for] [Will not harm] most hard non-porous [house hold] surfaces, [no bleach (or chlorine) stains]
Convenient
Convenient, quick & easy cleaning all around [the office[±]] [business[±]] [workplace[±]]
Destroys bacteria, *viruses, fungi, mycobacteria[†] and mildew on hard non-porous surfaces found in swine, poultry, dairy, equine farms and animal holding and research facilities
Destroys bacteria, fungi and *viruses on hard non-porous surfaces found in [shelters] [and] [or] [vet] [clinics] [practices]
Deodorizes odor-causing bacteria
Disinfects surfaces by killing [99.9%] SARS-Related Coronavirus 2[±] [SARS-CoV2] [Strain USA-WA1/2020] [in One Step³]
Disinfects (Deodorizes) (and Sanitizes) as it cleans³
Disinfects (Deodorizes) hard, non-porous surfaces³
Disinfects as it cleans³
Disinfects - kills germs for a deep down clean³
Disinfects headsets and telephones
Disinfects in 1 minute[¶]
Disposable
Dual action formulation providing broad spectrum activity against; fungi, mycobacteria[†], *viruses, and bacteria
Each pre-moistened towelette [disposable wipe] cleans, deodorizes and kills listed bacteria all around your [home[±]] [office[±]] [business[±]] [workplace[±]]
Easy flip-top lid
Effective against flu virus [Influenza A (Strain A/Hong Kong/8/68)]
Effective daily cleaner
Effective [shelter] [and] [or] [vet] [clinic] [practice] disinfectant
Effective in 1 minute against *Staphylococcus aureus*, *Salmonella enterica*, and *Pseudomonas aeruginosa*
Effectively disinfects hard, non-porous surfaces
Effectively disinfects in the presence of organic soil [5% blood serum]
Everyday, every room disinfectant
Exhibits effective virucidal activity: kills Influenza A(Strain A/Hong Kong/8/68) and Herpes Simplex Type 2
Easy to use – just wipe [and] [walk away] For non-antimicrobial uses only.
Effective against SARS-Related Coronavirus 2[±] [SARS-CoV-2] [Strain USA-WA1/2020] [in 1 minute]
Effective against SARS-CoV-2 virus[±] (agent of the respiratory illness COVID-19) responsible for the 2019 severe acute respiratory syndrome/disease/COVID-19 pandemic
Effective against 2019 SARS-Cov-2[±] (COVID 19 virus)
Effective – Cleans soap scum
Effective against [Kills] {insert any Bacteria claims from table 1} {OR} {insert any Virus claims from table 2} {OR} {insert any Fungi claims from table 4}

{*Administrative note*: Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike through~~ text means it has been removed.

Kills odor-causing bacteria
Eliminates odors [caused by bacteria], [mold], [mildew]
For hard, non-porous surfaces
For multiple cleaning tasks
For use as a mold preventative or remediation treatment in areas where mold [can be] [or has been] a problem.
For use [in] [Hospitals], [Health Care], [Daycare] settings [as a multi-purpose disinfectant cleaner].
Formulated with patented hydrogen peroxide technology
For everyday use [in institutional and (industrial) (institutional) facilities]
For a cleaner, fresher office (school) (college)
For any non-food contact hard, non-porous surface where germs are a concern.
For use [in] [on] [the] {*insert any institutional use site from table 8*} {*insert any residential use sites from table 12*} [as a multi-purpose disinfectant cleaner]
For non-scratch cleaning of showers and tubs, shower doors, fixtures and exterior surfaces of toilet bowls
For use [in] [Household], [Homecare], [Daycare] settings [as a multi-purpose disinfectant cleaner].
For use [in] [Animal] [Health] [premises] [settings] [facilities] [Animal Research Facilities]
For everyday use [in kitchens[±] and] [bathrooms[±]] [restrooms[±]] [all over the house[±]].
For a cleaner, fresher household
For everyday use
For use [in] [on] [the] {*insert any residential/household use site from table 12*} [as a multi-purpose disinfectant cleaner]
For your toughest bathroom messes
Good for most indoor cleaning tasks
Grease Cutting Formula
Great for grease cutting
Great for on-the-go cleaning
Great on-the-go disinfectant [cloth]
Hard non-porous surface [cleaning &] disinfecting made easy
Heavy-duty cleaner [or cleaning] [cloth]
Help make your {*insert any use site from table 8 or table 12*} [*insert product name*] clean
Intense cleaning action
Just pull, wipe & toss!¹
Kills [Eliminates 99.9%] SARS-Related Coronavirus 2[±] [SARS-CoV-2] [Strain USA-WA1/2020] [causative agent of COVID-19] [the virus that may cause COVID-19]
Kills [99.9% of] SARS-CoV-2[±], which causes COVID-19⁶
Kills [Effective against] the virus[±] that causes COVID-19⁶
Kills SARS-CoV-2 [COVID 19 virus] [virus] on hard non-porous surfaces
Kills [COVID-19 Pandemic] SARS-CoV-2 Virus[±] [formerly called SARS-nCoV]
Kills Bacteria in 1 minute
Kills *Viruses in 30 seconds, unless otherwise noted
Kills [99.9% of] germs [in 1 minute[¶]]
Kills [99.9% of] [the] [cold and flu viruses[~]] in 30 seconds
Kills antibiotic resistant bacteria²
Kills bacteria
**Kills bacteria commonly found in the [kitchen[±]] [and] [bathroom[±]]
Kills [Effective against] cold flu virus[~] [on hard, non-porous surfaces]
Kills germs
Kills germs while it cleans³
Kills household bacteria
Kills [household] bacteria [and sanitizes]
**Kills kitchen[±] and bathroom bacteria
Kills odor-causing bacteria
Kills odor-causing (household) bacteria
Kills (Destroys) 99.9% of germs
Kills (Destroys) 99.9% of germs while it cleans³
**Kills (Eliminates) 99.9% of odor-causing bacteria in the [kitchen[±]] [bathroom[±]] [restroom[±]] [all around the house[±]]
***Kills [kitchen[±]] [bathroom] [household] germs

{*Administrative note*: Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike through~~ text means it has been removed.

****Kills bacteria commonly found in kitchens[±] [in 1 minute]**

*****Kills bacteria and *viruses on hard, non-porous [kitchen[±]] [bathroom] [restroom] [household] surfaces in 1 minute[†]**

Kills Avian Influenza A virus on pre-cleaned environmental surfaces

Kills *Staphylococcus aureus*, *Salmonella enterica*, and *Pseudomonas aeruginosa*

Kills the flu virus Influenza A

Kills Influenza Virus

Kills Influenza A Virus

Kills Parvovirus[◊]

Kills Norovirus [Norwalk virus]

Kills Human Coronavirus

Kills antibiotic-resistant organisms²

Kills HIV-1 [AIDS virus] [HTLV-III_B]

Kills HCV[†]

Kills HBV^{††}

Kills Rhinovirus [Type 37]

Leaves [a fresh, clean scent] [no] [dull] [sticky] residue].

Leaves [office] [washroom] [restroom] smelling clean with a clean fresh scent

Leaves surfaces virtually streak-free

Leaves your [home] [house] [washroom] [restroom] smelling clean with a clean fresh scent

Leaves your kitchen[±] and bathroom with a fresh, clean scent

Long shelf life: 2 years

No streaking [filming] [sticky residue]

Non-tainting, leaves no unsightly residues

[No rinse required]. No streaking [or filming] formula

Non-corrosive [to eyes] [and] [skin]

No mixing, measuring or rinsing required

Not a skin [Dermal] sensitizer [non-sensitizer]

Non-woven disposable cloth

On-the-go disinfectant

Oxy-1 Wipes controls listed bacteria between treated surfaces[±]

Oxy-1 Wipes help control listed bacteria between treated surfaces[±] while keeping your office [or workstation or workplace or business] clean and dust-free

Powerful cleaning

Powers away tough stains [stuck on dirt]

Powerful bubbles blast away [dried] [stuck on] dirt [grime] [food] [spills] [grease] [greasy messes]

Powerful odor control [odor neutralizer] [deodorizer] [odor-counteractant] formula [on hard non-porous surfaces]

Primary Skin Irritation Category IV

Primary Eye Irritation Category IV

Pre-measured

Pre-moistened [towelettes]

Pre-soaked [cloth or towelette]

Quickly [eliminates] [removes] [food] odors

Quick and easy cleaning everywhere

Quick drying

Removes tough stains

Removes [the most common] [household] allergens [around your home] [around your kitchen[±]] [around your bathroom] [on your floors]⁺⁺

{*Note: the sentence above will be used with this claim qualifier*}: ⁺⁺This product [or product name] removes the following [non-living] allergens: [cat dander], [cockroach matter], [cockroach debris], [cockroach particles], [dog dander], [dust mite debris], [pet dander], [pollen particles]

Restroom[±] Disinfectant Wipes

Restroom Wipes [for hard, non-porous surfaces]

Sanitizes in 10 seconds

Specially designed to dissolve [penetrate] the toughest [stubborn] soils¹

Suitable for everyday use

{*Administrative note:* Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

Surface[±] cleaner (with antibacterial formula)
Suitable for use on [floors], equipment, cages and kennels, when used as directed
Step-saving, Labor-saving, Time-saving, Cost-saving formula
Smell of clean [is no smell]
Stay-moist package
Strong towelettes [or wipes]
Telephone [cleaner [and] disinfectant] wipes
These wipes are a convenient way to disinfect surfaces all around your business[±] [office[±]] [workplace[±]]
To be used in sanitizing and disinfecting of hard, non-porous non-food surfaces and non-critical equipment which may harbor germs
This product [*or product name*] is a broad spectrum, hard non-porous surface disinfectant that has been shown to be effective against SARS-CoV-2 [COVID-19 virus tested and listed on the label]
This product [*or product name*] has demonstrated effectiveness against SARS-CoV-2[±] [Human Coronavirus SARS-CoV-2 (the virus that may cause COVID-19)]
[Tough on germs.]
Towelette
Use only one product to get a streak-free shine plus antimicrobial action
Use this product [or product name] throughout your home to clean [beverage stains] [dirt] [food stains] [grime] [pet stains] [soap scum]
[Uses] [harnesses] the cleaning power of hydrogen peroxide
Use on hard, non-porous surfaces only. [Test wipe on a small inconspicuous area first]
Use to disinfectant shopping [or grocery] cart handles
Virucidal* [or] Kills viruses*
Will not leave grit or soap scum
Wipe down grocery [or shopping] cart handles

{GENERAL CLAIMS FOR A FLOOR CLEANER}:

[Acceptable for] [great for] [Perfect for] use on all finished floor types⁵
Great for use on [glazed tile] [glazed ceramic] [vinyl] [laminated] [linoleum] [synthetic (sealed) marble] [sealed stone] wood⁵
/ Try it on wood floors⁵
Long lasting freshness
New {*Note to reviewer: This statement will only be used for 6 months after the launch of a new product*}
New & Improved {*Note to reviewer: This statement will only be used for 6 months after the launch of a new product*}
Refills
Starter Kit
Swiffer Sweeper [antibacterial]
Toss away the dirt and grime
[Wet] Mopping clothes

{PRODUCT SPECIFIC CLAIMS}:

Cleans and removes bacteria³
Cleans and sanitizes⁴ [floors] [wall to wall]
Dual Action Cleaner³ – kills germs and cleans dirt
[Green] Trap & lock strip
Long lasting freshness
Formulated with [hydrogen] peroxide (Hispanic Translation)
[Quickly] cleans floors to a shine
Scrubbing strip powers through [tough] [dried-on] [sticky] [messes] [dirt]
[Textured wet cloths] [cleansers] Traps & Locks Dirt [& Grime] [deep in cloth]
Textured Wet Cloths Cut Through Dirt [& Grime]
With the oxygenated lift-action of Hydrogen Peroxide

{BONUS PACK CLAIMS}:

Bonus Pack

{*Administrative note:* Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

Bonus Size
Bundle Pack
Club Size
Combo Pack
Spring Cleaning [Bundle] Pack
Value Pack
Value Pack, 3 Pack, Three Pack, 2 Pack, Two Pack, Large Pack, Large Size
Wet & Dry [Bundle] Pack

{*Note: the following section on summary of claims is optional and might be inserted in a band at the bottom of the label*}

***VIRUCIDAL: 30 SECONDS**
ANIMAL PREMISES *VIRUCIDAL: 1 MINUTE
BACTERICIDAL: 1 MINUTE
NON-FOOD CONTACT SANITIZER: 10 SECONDS
FUNGICIDAL: 30 SECONDS
TUBERCULOCIDAL: 1 MINUTE

{*Note to the reviewer: the following short claims might be inserted on the container label only for the period of 6 months after the launch*}

[Better] [Improved] [New] [New and improved] [formula] [packaging]
[New and improved]
[New Claims]
[New Kill Claims]
[New Claims - Kills TB⁺ in 1 Minute]
[New - Saves Time]
[New] - [Kills TB⁺ in 1 Minute]

[FEATURES, CLAIMS & USES:]

[Notice to User: This product {*or product name*} may not be suitable for all items or surfaces. Prolonged exposure to copper, brass, anodized aluminum and some coatings should be avoided. If you are uncertain as to the material composition of your item, confirm with the manufacturer before proceeding [or visit the Compatibility section on www.virox.com for a listing of compatible materials.]

[For General Uses – Medical and Dental]

For Use in Health Care [facilities] [clinics], {*or insert any of the Medical/Dental Institutional Use sites from table 8*}.

This product [These] is [are] [is a] non-woven, disposable [towelette[s]] [wipe[s]], pre-saturated with Hydrogen Peroxide [based] disinfectant cleaner.

This product [*or product name*] is a one-step [hospital-use] [germicidal] [disinfectant] [cleaner] [and] [deodorant] [odor-counteractant] [odor neutralizer] [wipe] designed for [general] cleaning, [and] disinfecting, [deodorizing] [of] hard, non-porous [environmental] [surfaces].

It cleans quickly by removing {*insert any soil from table 5*} found in {*insert any of the Medical/Dental Institutional Use sites from table 8*}.

This product [*or product name*] is designed for use on hard, non-porous [environmental] surfaces: {*insert any surface material from table 6*} associated with [high] [hand] touch [environmental] surfaces [including] [such as] [floors], [walls], [ceilings] - {*or insert any Medical/Dental hard non-porous surface from table 7*} – or – any hard non-porous washable (non-food contact) surface where disinfection is required. [This product's non-dulling formula eliminates the time and labor normally required for rinsing]. [It [also] reduces odor-causing bacteria leaving surfaces smelling clean and fresh].

[For Veterinary Use]

{**Administrative note:** Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike through~~ text means it has been removed.

For Use in Veterinary clinics {or insert any of the Veterinary Institutional Use sites from table 8}

This product [or product name] is a one-step [germicidal] [disinfectant] [cleaner] [and] [deodorant] [odor-counteractant] [odor neutralizer] [wipe] designed for [general] cleaning, [and] disinfecting, (deodorizing) [of] hard, non-porous [environmental] surfaces.

It cleans quickly by removing {use any soil from table 5} found in [Animal][Premises][facilities][such][as] {insert any of the Veterinary Institutional Use sites from table 8}.

This product [or product name] is designed for use on hard, non-porous [environmental] surfaces: {insert any surface material from table 6} associated with [high] [hand] touch [environmental] surfaces, [including] [such as] [floors], [walls], [ceilings] – {or insert any Veterinary hard non-porous surface from table 7} – or – any hard non-porous washable (non-food contact) surface where disinfection is required. [This product's non-dulling formula eliminates the time and labor normally required for rinsing]. [It [also] reduces odor-causing bacteria leaving surfaces smelling clean and fresh].

[For Miscellaneous Use]

For Use [in] [on] [at] {insert any of the Miscellaneous Institutional Use site from table 8}

This product [or product name] is a one-step germicidal [disinfectant] [cleaner] [and] [deodorant] [odor-counteractant] [odor neutralizer] [wipe] designed for general cleaning, [and] disinfecting, (deodorizing) [of] hard, non-porous environmental surface.

It cleans quickly by removing {use any soil from table 5} found in {insert any of the Miscellaneous Institutional Use sites from table 8}.

This product [or product name] is designed for use on hard, non-porous [environmental] surfaces: {insert any surface material from table 6} associated with [high] [hand] touch [environmental] surfaces [such as] [including], [floors], [walls], [ceilings] – {or insert any Veterinary hard non-porous surface from table 7} – or – any hard non-porous washable (non-food contact) surface where disinfection is required. [This product's non-dulling formula eliminates the time and labor normally required for rinsing]. [It [also] reduces odor-causing bacteria leaving surfaces smelling clean and fresh].

(For Use in Spas, (Salons) & (Clinics))

(Surfaces) (and) (Footbaths) (pedicure foot spa equipment) in (health clubs), ((medical) spas), (tanning salons), (tanning), (hair) (nail) (pedicure) (salons), (barber) (beauty) (shops), (salons), (tattoo parlors) (massage therapy (clinics)):

This product [These] is [are] [is a] non-woven, disposable [towelette[s]] [wipe[s]], pre-saturated with Hydrogen Peroxide [based] disinfectant cleaner.

This product [or product name] is a one-step [hospital-use] [germicidal] [disinfectant] [cleaner], [*virucide] [and] [deodorant] [odor-counteractant] [odor neutralizer] [wipe] designed for cleaning, [and] disinfecting, [and] (deodorizing) [on] [of] hard, non-porous surfaces. It cleans quickly by removing body oils, dead tissue, soil and all other buildups or organic matter [found] on [hard non-porous (inanimate) [environmental]] [footbath] surfaces.

This product (Product Name) (is) (are) designed for cleaning and disinfection of hard non-porous surfaces in {insert professional beauty use sites from table 8}

It is designed for [daily] use on [hard, non-porous environmental (inanimate) surfaces] [such as] {insert professional beauty surfaces from table 7}

[For] [Use on] [Disinfection of] [Surfaces] [in] [Pharmaceutical] [and] [Compounding] [Processing] Facilities {or insert institutional use site from table 8}:

[This product [or product name] is for use on hard, non-porous [environmental] surfaces such as [floors], [walls] [and] [processing equipment] in [pharmaceutical] [and] [cosmetic] [compounding] [processing] facilities {or insert miscellaneous use site from table 8}.

{*Administrative note*: Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

[For Use in] Mold & Mildew Remediation [& Restoration]:

This product is ideal for use as a [cleaner], [disinfectant] and [deodorizer] in {*insert mold/mildew remediation use site from table 8*}. [This product was designed to leave surfaces clean and odor free.

This product [*or product name*] kills the growth of mold and mildew and the odors caused by [associated to] them on hard, non-porous surfaces such as [sealed] [concrete] [floors], [sealed walls]; [sealed drywall] [and] [gypsum board], [tabletops[±]], [bench-tops[±]], [drain troughs], [drip pans], [desks], [locker room benches], [and] [hard non-porous bath mats (e.g., rubber or plastic).] {*or insert institutional/miscellaneous surface from table 7*}

[For Use in] [Water/Flood Damage]:

This product kills the growth of mold and mildew, fungi, and bacteria, and [reduces] [removes] the odors caused by them on hard, non-porous environmental surfaces as a result of water and flood damage.

[For Use in] [Sewage Backflow]:

This product is a broad-spectrum [cleaner] [and] [deodorizer] that kills [reduces] [removes] (non-public health) (odor causing) organisms left on surfaces previously contaminated with raw sewage and black water contamination.

[For Use in] [Trauma, Crime Scene, and Fire Clean Up]:

OSHA Bloodborne Pathogen Standard Compliance. This product effectively cleans and disinfects³ surfaces and equipment that come into contact with blood and other potentially infectious materials (OPIM).

[For Use in] Equipment and PPE Disinfection (Non-Critical Devices):

To prevent cross contamination between treated surfaces[±], this product is effective for cleaning and disinfecting reusable person protective equipment (PPE), [gear], [tools], and [restoration] [equipment] {*or insert surface from table 7*} after removal from containment areas.

[For Use on] HVAC [hard non-porous surfaces of HVAC/R equipment]:

This product is suitable for hard, non-porous (exterior) surfaces of HVAC/R equipment, air conditioning equipment, refrigeration equipment and heat pumps. Allow surfaces to come to room temperature [prior to treatment]

[For] [Residential Use] [Home][Household]

This product [*or product name*] is an effective [cleaner] [disinfectant] [sanitizer] [wipe] on non-food contact surfaces in [around] {*insert residential use surfaces from table 11 or residential/household use site from table 12*}.

This product [*or product name*] is a[n] [Antibacterial] [Multi-Purpose] [Bathroom] [Floor] [Glass & Surface] [Kitchen[±]] [Restroom] [Cleaner] [Cleaning] [Disinfectant] [Disinfecting] [Sanitizer] [wipe]. This product is for use [in] [on] [the] {*insert any residential/household use sites from table 12*} [common hard non-porous household surfaces].

For everyday cleaning [and disinfecting] [use daily to clean [and disinfect]] [on] hard, non-porous surfaces [where germs^Δ are a concern] [made from] {*insert any surface material from table 10*} associated with kitchen[±] surfaces, [high] [hand] touch surfaces [all over the house], {*or – insert any of the hard non-porous surfaces listed in table 11*} and other hard [non-food contact] surfaces around the house [hold] [home]. [Not recommended for use on brass], [polished wood], [granite], [marble], [and] [leather].

[For] [Use in] [Greenhouse] [and] [Horticulture] [Facilities] [Use]

For use to disinfect (and clean) (environmental) hard non-porous (inanimate) surfaces [:] [such as] [including] (floors), (walls), (glasshouse structures), (ventilation¹) and (other) (equipment), (utensils), (trays), and (other) (containers), (water systems¹), (evaporative coolers), (storage rooms), and (vehicles) in (greenhouses) and (other horticultural settings) {*or insert horticultural use site/surface from table 8 or table 7*} (prior) (to) (introduction) or (reintroduction) of (plants), (seeds), or (soil).

This product [*or product name*] is a one-step [germicidal] disinfectant [cleaner] [and deodorant] designed for [general] cleaning, [and] disinfecting, (deodorizing) of hard, non-porous [environmental] surfaces³.

{**Administrative note:** Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

It cleans quickly by removing {*use any soil from table 5*} found in {*insert use sites from table 8*}. This product [*or product name*] is designed for use on hard, non-porous [environmental] surfaces [such as] [including] {*insert hard non-porous surface use site from table 7 or institutional surfaces from table 6*}

[For Use in (Laboratory Settings) (and) (Alternate Care Clinics)]

This product is ideal for use as a cleaner and disinfectant³ in Laboratory Settings such as (Life Sciences), (CAT Laboratories), (Pharmacies), (Animal Life Science Laboratories), (Vivariums), (Cleanrooms), (Critical Environments), (GMP facilities), (Compound pharmacy labs), (Biopharmaceutical), (Medical device facilities).

This product is ideal for use as a cleaner and disinfectant in Alternate Care Settings, Clinics & Offices such as (Sleep Therapy), (Respiratory Therapy), (Ophthalmic), (Massage), (Chiropractic), (Acupuncture), (Naturopathic), (Homeopathic), (Home Healthcare), (Otolaryngology (ear, nose and throat)), (podiatry), (foot care nursing), (midwifery). This product can also be used in (Emergency Care) (EMS) (Fire / Police / Ambulance), (Transportation Facilities) and (Airports).

This product is a one-step germicidal disinfectant cleaner and odor neutralizer designed for cleaning and disinfecting of hard, non-porous environmental surfaces. This product cleans by removing dirt, grime, food residue, dead skin, blood, biological waste (feces, urine, etc), and other organic matter commonly found on surfaces. This product was designed to leave surfaces clean, disinfected and odor free.

This product cleans, disinfects, and deodorizes hard, non-porous environmental surfaces in one step³. Its nonabrasive formula is designed for use on the following non-porous environmental surfaces (including) (such as) {*insert any surface material from table 6*} associated with floors, walls, ceilings, tables, chairs, countertops, telephones, fixtures, (fume hoods), (vivariums), (cages), (carts), (workstations), and (sinks), (including) (such as) those found in (laboratory settings).

[For use in {*insert use sites from table 8*}]

[This product [*or product name*] is a one-step³ [germicidal] disinfectant [cleaner] [and deodorant] designed for [general] cleaning, [and] disinfecting, (deodorizing) of hard, non-porous [environmental] surfaces.]

[It cleans quickly by removing {*use any soil from table 5*} found in {*insert use sites from table 8*}. This product [*or product name*] is designed for use on hard, non-porous [environmental] surfaces [such as] [including] {*insert hard non-porous surface use site from table 7*}

{Claims against Emerging Viral Pathogens}

{*This product qualifies for emerging viral pathogen claims per the EPA's 'Guidance to Registrants: Process for Making Claims Against Emerging Viral Pathogens not on EPA-Registered Disinfectant Labels' when used in accordance with the appropriate use directions indicated below.*

This product meets the criteria to make claims against certain emerging viral pathogens from the following viral category[ies]:

- Enveloped Viruses
- Large, Non-enveloped viruses
- Small, Non-enveloped viruses

For an emerging viral pathogen that is a/an...	...follow the directions for use for the following organisms on the label:
Enveloped virus	Poliovirus type 1
Large non-enveloped virus	Poliovirus type 1
Small non-enveloped virus	Canine Parvovirus

The statements shall be made only through the following communications outlets: technical literature distributed exclusively to health care facilities, physicians, nurses and public health officials, "1-800" consumer information services,

{*Administrative note:* Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. Strike-through text means it has been removed.

social media sites and company websites (non-label related). These statements shall not appear on marketed (final print) product labels.

[Product name] has demonstrated effectiveness against viruses similar to [name of emerging virus] on hard, non-porous surfaces. Therefore, [product name] can be used against [name of emerging virus] when used in accordance with the directions for use against {*Select one: For Small Non-enveloped Virus*} Canine Parvovirus {*For Enveloped or Large non-enveloped virus*} Poliovirus type 1} on hard non-porous surfaces. Refer to the [CDC or OIE] website at [pathogen-specific website address] for additional information.

{OR}

[Name of illness/outbreak] is caused by [name of emerging virus]. [Product name] kills similar viruses and therefore can be used against [name of emerging virus] when used in accordance with the directions for use against {*Select one: For Small Non-enveloped Virus*} Canine Parvovirus {*For Enveloped or Large non-enveloped virus*} Poliovirus type 1} on hard non-porous surfaces. Refer to the [CDC or OIE] website at [website address] for additional information.

{LIST OF APPROVED ORGANISM} Claims

[See Reference sheet for a complete list of [pathogenic] organisms] [this product [or product name] is effective against]

{Note to reviewer: ATTC numbers are optional and might not appear on the final label}

{Table 1 – Bacteria Claims}

BACTERICIDAL [Kills] [Effective Against] [Bacteria] [in 1 minute] [against]:

[Bacteria][Bactericidal] (Activity)	Contact Time
<i>Acinetobacter baumannii</i> (Multi-Drug Resistant (MDR))[(ATCC 19606)]	1 minute
<i>Bordetella bronchiseptica</i> [(ATCC 10580)]	1 minute
<i>Bordetella pertussis</i> [(ATCC 12743)]	1 minute
<i>Burkholderia cepacia</i> [(ATCC 25416)]	1 minute
<i>Chlamydia psittaci</i> (Strain 6 BC) [ATCC VR-125]	1 minute
<i>Enterobacter aerogenes</i> [(ATCC 13048)]	1 minute
<i>Enterobacter cloacae</i> (New Delhi metallo- β -lactamase 1[producing] [positive] (NDM -1)) [(CDC 100654)]	1 minute
<i>Enterococcus faecalis</i> (Vancomycin Resistant (VRE)) [(ATCC 51575)]	1 minute
<i>Enterococcus faecium</i> [(ATCC 51559)]	1 minute
<i>Escherichia coli</i> [(ATCC 11229)]	1 minute
<i>Escherichia coli</i> O157:H7 [(ATCC 35150)]	1 minute
<i>Escherichia coli</i> (Extended Spectrum Beta-lactamase (ESBL) [producing][positive]) [(ATCC BAA-196)]	1 minute
<i>Escherichia coli</i> (New Delhi metallo- β -lactamase 1(NDM-1) [producing][positive]) [(CDC 1001728)]	1 minute
<i>Klebsiella pneumoniae</i> [(ATCC 4352)]	1 minute
<i>Klebsiella pneumoniae</i> (Carbapenem Resistant (CRE)) [(ATCC BAA 1705)]	1 minute
<i>Klebsiella pneumoniae</i> (Extended Spectrum Beta-lactamase (ESBL) [producing][positive]) [(ATCC700603)]	1 minute
<i>Klebsiella pneumoniae</i> (New Delhi metallo- β -lactamase 1(NDM-1) [producing][positive]) [(CDC 1000527)]	1 minute
<i>Legionella pneumoniae</i> [(ATCC 33153)]	1 minute
<i>Listeria monocytogenes</i> [(ATCC 19117)]	1 minute
<i>Pasteurella multocida</i> [(ATCC 12947)]	1 minute
<i>Pseudomonas aeruginosa</i> [(ATCC 15442)]	1 minute
<i>Pseudomonas aeruginosa</i> – (Metallo-beta lactamase (MBL) positive) [(CDC 2012059)]	1 minute
<i>Salmonella enterica</i> (formerly known as <i>Salmonella choleraesuis</i>) [(ATCC 10708)]	1 minute
<i>Serratia marcescens</i> [(ATCC 14756)]	1 minute
<i>Shigella dysenteriae</i> [(ATCC 11835)]	1 minute

{*Administrative note*: Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike through~~ text means it has been removed.

<i>Staphylococcus aureus</i> [(ATCC 6538)]	1 minute
<i>Staphylococcus aureus</i> [Community Acquired Methicillin Resistant Genotype USA 400 (CA-MRSA 400)]	1 minute
<i>Staphylococcus aureus</i> (Methicillin Resistant (MRSA)) [(ATCC 33592)]	1 minute
<i>Staphylococcus aureus</i> (Vancomycin Intermediate Resistant (VISA)) [(HIP 5836)]	1 minute
<i>Staphylococcus aureus</i> (Vancomycin Resistant (VRSA)) [(VRS1)]	1 minute
<i>Staphylococcus epidermidis</i> (Methicillin Resistant (MRSE)) [(ATCC 51625)]	1 minute
<i>Staphylococcus pseudintermedius</i> [(ATCC 49444)]	1 minute
<i>Stenotrophomonas (Pseudomonas) maltophilia</i> [(ATCC 25596)]	1 minute
<i>Streptococcus pneumoniae</i> (Penicillin Resistant (PRSP)) [(ATCC 700677)]	1 minute
<i>Streptococcus pyogenes</i> [(ATCC19615)]	1 minute

[commonly found] [on hard non-porous environmental (inanimate) surfaces] [in {insert one or more use site from table 8}]
[using approved AOAC test methods] [under Good Laboratory Practices, [GLP's]], [in the presence of 5% serum load]
[and] [at] a 1 minute contact time [in 1 minute].

[Bactericidal Activity was determined using the AOAC Germicidal Spray Test Method]

TUBERCULOCIDAL [Kills] [Effective Against] [in 1 minute] [against]: {*Note to reviewer: Tb claim is not to be used in any label with residential language*}

[Tuberculocidal] (Activity)	Contact Time
[†] <i>Mycobacterium bovis</i> (BCG)	1 minute

[using approved AOAC test methods] [under Good Laboratory Practices, [GLP's]], [in the presence of 5% serum load]
[and] [at] a 1 minute contact time [in 1 minute].

{Table 2 – Virus Claims}

***VIRUCIDAL [Kills] [(is) Effective Against] [enveloped and non-enveloped] [*Viruses] [in 1 minute [unless otherwise noted]] [in 30 seconds]:** {*Note to reviewer: if viruses are NOT listed under the *Virucidal or *Virus heading, then each individual virus will be qualified by an **}

[*Viruses] [*Virucidal] (Activity):	Contact Time
Adenovirus Type 8 [(ATCC VR-1368)]	30 seconds
Avian Influenza A Virus (H7N9)	30 seconds
Avian Influenza A Virus (H3N2)	30 seconds
Cytomegalovirus [(ATCC VR-538)]	30 seconds
Feline Calicivirus [(ATCC VR-782)]	30 seconds
Hepatitis B Virus (HBV) (Duck Hepatitis B Virus as the surrogate)	30 seconds
Hepatitis C Virus (HCV) (Bovine Viral Diarrhea Virus as the surrogate)	30 seconds
Herpes Simplex Virus Type 1 [(ATCC VR-733)]	30 seconds
Herpes Simplex Virus Type 2 [(ATCC VR-734)]	30 seconds
HIV Type 1 (Strain HTLV III _B)	30 seconds
Human Coronavirus [(ATCC VR-740)]	30 seconds
Influenza A Virus (Strain A/Hong Kong/8/68) [(ATCC VR-544)]	30 seconds
Influenza B Virus [(ATCC VR-823)]	30 seconds
Norovirus (Feline Calicivirus as the surrogate)	30 seconds
Pandemic Influenza A Virus H1N1	30 seconds
Parainfluenza Virus Type 3 [(ATCC VR-93)]	30 seconds
Poliovirus Type 1 (Strain Chat) [ATCC VR-1562]	30 seconds
Respiratory Syncytial Virus (RSV) [(ATCC VR-26)]	30 seconds
Rhinovirus Type-37 [(ATCC VR-1147)]	30 seconds
Rotavirus (WA strain)	30 seconds

{**Administrative note:** Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. Strike-through text means it has been removed.

⁶ SARS-Related Coronavirus 2 [Human Coronavirus SARS-CoV-2 (the virus that may cause COVID-19*)] [SARS-CoV-2] [Strain USA-WA1/2020] [BEI NR-52281]	1 minute
Vaccinia Virus [(ATCC VR-119)]	30 seconds

[using approved ASTM test methods] [under Good Laboratory Practices, [GLP's]], [in the presence of 5% serum load]
[and] [in] 30 seconds [contact time].

[Kills] [(is) Effective Against] [Following] [Veterinary *Viruses] [Animal premises *Viruses] [Animal premises *Virucidal [Activity]] [in 1 Minute, unless otherwise noted]

[Veterinary] [Animal Premise] [*Viruses] [*Virucidal] (Activity):	Contact Time
Avian Reovirus (Strain 2177) [ATCC VR-2449]	1 minute
Bovine Viral Diarrhea Virus	30 seconds
Canine Distemper Virus [ATCC VR-128]	1 minute
Canine Influenza A Virus (H3N2) (USDA NVSL 004-IDV)	1 minute
Canine Parvovirus (CPV) (Strain Cornell) [ATCC VR-953]	5 minutes
Duck Hepatitis B virus	30 seconds
Equine Herpes virus type 1 (Strain EHv1)	1 minute
Equine Influenza A Virus (H3N8) (EIV)	30 seconds
Feline Calicivirus	30 seconds
Feline Herpes Virus (Feline Rhinotracheitis) [(ATCC VR-636)]	1 minute
Feline Infectious Peritonitis Virus (FIPV) [(ATCC VR-2202)]	1 minute
Feline Leukemia Virus (Strain CT600) [ATCC VR-1373] (FeLV)	1 minute
Feline Panleukopenia Virus [ATCC VR-648]	1 minute
Foot and Mouth Disease virus (FMDv) (Type O1) (BFS 1860)	3 minutes
Marek's Disease virus (Strain SB) (Charles River Laboratories)	1 minute
Minute Virus of Mice (Murine Parvovirus) [(ATCC VR-1346)]	1 minute
Mouse Hepatitis Virus [(ATCC VR-765)]	1 minute
Murine Norovirus (Strain MNV-1.CW1)	30 seconds
Newcastle disease virus (Strain B1) (Hitchner or Blacksburg) [ATCC VR-108]	1 minute
Porcine Epidemic Diarrhea virus (Strain Colorado (2013 Isolate)) (PEDv)	1 minute
Porcine Respiratory & Reproductive Syndrome virus (Strain NVSL) (PRRS)	1 minute
Porcine Rotavirus (Strain OSU) [ATCC VR-893]	1 minute
Pseudorabies virus (Strain Aujeszky) [ATCC VR-135]	1 minute
Vesicular Stomatitis virus (Strain Indiana Lab) [ATCC VR-158]	1 minute

{Note to reviewer: if viruses are NOT listed under the *virucidal or *virus heading then each individual virus will be qualified by an *}

{Table 3 – Sanitization Claims}

NON-FOOD CONTACT SANITIZER [in 10 seconds] against:

[Sanitizing] [Non-Food Contact Surface Sanitizer]:	Contact Time
<i>Acinetobacter baumannii</i> [(ATCC 19606)]	10 seconds
<i>Campylobacter jejuni</i> [(ATCC 29428)]	10 seconds
<i>Enterobacter aerogenes</i> [(ATCC 13048)]	10 seconds
(Vancomycin Resistant) <i>Enterococcus faecalis</i> (VRE) [(ATCC 51575)]	10 seconds
<i>Escherichia coli</i> O157:H7 [(ATCC 35150)]	10 seconds
<i>Listeria monocytogenes</i> [(ATCC 19117)]	10 seconds
<i>Pseudomonas aeruginosa</i> [(ATCC 15442)]	10 seconds
<i>Salmonella enterica</i> (formerly known as <i>Salmonella choleraesuis</i>) [(ATCC 10708)]	10 seconds
<i>Staphylococcus aureus</i> [(ATCC 6538)]	10 seconds
(Methicillin Resistant) <i>Staphylococcus aureus</i> (MRSA) [(ATCC 33592)]	10 seconds

{*Administrative note*: Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

{Table 4 – Fungi Claims}

FUNGICIDAL [Kills] [Effective Against] [FUNGI] [in 1 minute] [in 30 seconds, unless otherwise noted] [against]:

[Fungi] [Fungicidal] (Activity):	Contact Time
<i>Candida albicans</i> [(ATCC 10231)]	1 minute
<i>Microsporum canis</i> [(ATCC 10214)]	1 minute
<i>Trichophyton interdigitale</i> (formerly known as <i>Trichophyton mentagrophytes</i>) [(athlete's foot fungus)][(ATCC 9533)]	30 seconds

Mold/Mildew (Mildewcidal Activity) – [kills the growth of mold and mildew (at a 4 minute contact time): *Aspergillus brasiliensis* (formerly known as *Aspergillus niger*) [(ATCC 16404)] (and the odors caused by them) when applied to hard, non-porous environmental surfaces.]

[Mold/Mildew] (Mildewcidal Activity)	Contact Time
<i>Aspergillus brasiliensis</i> (formerly known as <i>Aspergillus niger</i>)	4 minutes

Mold/Mildew (Mildewstatic Activity) – [controls and prevents the growth of mold and mildew (at a 4 minute contact time): *Aspergillus brasiliensis* (formerly known as *Aspergillus niger*) [(ATCC 16404)] (and the odors caused by them) when applied to hard, non-porous environmental surfaces.]

[Mold/Mildew] (Mildewstatic Activity)	Contact Time
<i>Aspergillus brasiliensis</i> (formerly known as <i>Aspergillus niger</i>)	4 minutes

{OR}

[CLAIMS] [ORGANISMS]:

[Bacteria] [Bactericidal]:

Kills {*insert one or more bacteria from table 1*} on hard, non-porous [kitchen‡] [bathroom] [restroom] [household] [environmental] surfaces in 1 minute.

*[Viruses] [*Virucidal]:

*Kills {*insert one or more virus from table 2*} when used as directed on hard, non-porous (inanimate) [kitchen‡] [bathroom] [restroom] [household] [environmental] surfaces in 30 seconds, unless otherwise noted.

[Tuberculocidal] in the presence of 5% serum load and 1 minute contact time at 20°C (68 F) on hard, non-porous environmental surfaces

† *Mycobacterium bovis* (BCG) [(OT 451C150)]

Tuberculocidal activity was determined by the AOAC Tuberculocidal Activity of Disinfectant spray Products

Malodor[s] [Activity] [Counteractancy] – eliminates [destroys] odors and odor-causing bacteria on hard, non-porous surfaces in restroom areas, behind and under sinks and counters, and storage areas [and other places where bacterial growth can cause malodors].

[Fungicidal]:

Kills {*insert one or more fungi from table 4*} on hard, non-porous [kitchen‡] [bathroom] [restroom] [household] [environmental] surfaces in 30 seconds, unless otherwise noted.

{**Administrative note:** Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

Non-Food Contact Surface Sanitizer for hard non-porous surfaces – [in the presence of 5% serum load] [and] [in] 10 seconds [contact time] [on hard, non-porous environmental surfaces]: {*insert one or more organism from table 3*}

DIRECTIONS FOR USE:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.
A potable water rinse is required for food contact surfaces. [For all others, no rinsing is required].

Not for personal cleansing. This is not a baby wipe. Do not flush in toilet.

{*Note to reviewer: The following paragraph will be used on labels that list medical devices as defined by the FDA.*}

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

For Use as [Daily] [One-Step] Cleaner/Disinfectant [product]:

1. Pre-clean visibly soiled areas
2. [Pull] [Remove] towelette [cloth] from dispenser [canister] [tub] and Wipe hard, non-porous [environmental] surface(s) [or] [floors⁵].
3. Allow surface(s) [Allow this product [or product name] to penetrate and] [to] remain visibly wet for 1 minute^{fl}. [Use additional towelettes [cloths] if needed to ensure the correct wet contact time.]
4. [Wipe] [surfaces dry] [or rinse]. [or] [Allow to air dry]. [or] [Let air dry]. [If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]]

[For Use to Clean and Disinfect] [For Use on] Life Science [Laboratory] [Surfaces], [Instruments], [Equipment] [and] [Glassware]:

1. Pre-clean visibly soiled areas.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe hard non-porous surfaces.
3. Allow surface(s) [Allow this product [or product name] to penetrate and] to remain visibly wet for 1 minute^{fl}. [Use additional towelettes if needed to ensure the correct wet contact time.]
4. [For glassware/instrument: Wipe surface thoroughly and let air dry before reuse.] [For surfaces: Wipe [surfaces dry] [or rinse]]. [If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]].

For Use as a Cleaner/Disinfectant in non-food handling areas of Food Establishments [Food Processing Plants]:

1. Pre-clean visibly soiled areas.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe hard, non-porous surfaces.
3. Allow surface(s) [Allow this product [or product name] to penetrate and] to remain visibly wet for 1 minute^{fl}. [Use additional towelettes if needed to ensure the correct wet contact time.]
4. Wipe [surfaces dry] [or rinse]. [If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]].

For Use as (Daily) Bathroom (Restroom) (One-Step) Cleaner/Disinfectant:

1. Pre-clean visibly soiled areas.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe hard, non-porous surfaces.
3. Allow surface(s) [Allow this product [or product name] to penetrate and] to remain visibly wet for 1 minute^{fl}. [Use additional towelettes if needed to ensure the correct wet contact time.]
4. Wipe [surfaces dry] [or rinse]. [If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]].

For Use as a Bowl Cleaner/Disinfectant in Exterior Surfaces of Toilet Bowls [Dry/Vacuum Flush Toilets] [Dry/waterless Urinals] [Toilet Bowls Above the Waterline]:

1. Pre-clean visibly soiled areas.

{*Administrative note:* Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike through~~ text means it has been removed.

2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe exterior surfaces of toilet bowls, [Dry/Vacuum Flush Toilets], [Dry/waterless Urinals], [Toilet Bowls Above the Waterline].
3. Swab entire surface area especially under the rim. Allow entire surface to remain visibly wet for one (1) minute[¶].
4. Rinse swab applicator thoroughly.

For Use to Clean and Disinfect Bathtub and shower stalls

1. Pre-clean visibly soiled areas.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe hard non-porous surfaces.
3. Allow surface(s) [Allow this product [or product name] to penetrate and] to remain visibly wet for 1 minute[¶]. [Use additional towelettes if needed to ensure the correct wet contact time.]
4. Rinse surfaces thoroughly and let air dry.

For Use to Clean and Disinfect Mannequins Used in CPR Training:

Manikins should be cleaned as soon as possible at the end of each class to avoid drying of contaminants on surfaces. Disassemble the manikin as directed by manufacturer's instructions. Wipe all surfaces (of manikins) with this product [or product name] until visibly wet. Let stand for 1 minute[¶] and rinse with potable water.

To Disinfect Surfaces in Pharmaceutical [and Cosmetic] Processing Facilities:

1. Pre-clean visibly soiled areas.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe hard non-porous surfaces.
3. Allow surface(s) [Allow this product [or product name] to penetrate and] to remain visibly wet for 1 minute[¶]. [Use additional towelettes if needed to ensure the correct wet contact time.]
4. Rinse with sterile water.
5. [Wipe surfaces dry (or rinse). (If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]. (Allow to air dry [Let air dry]).]

To Disinfect Non-Critical[¶], Pre-cleaned Instruments:

1. Instruments must be thoroughly pre-cleaned to remove excess organic debris, rinsed, and rough dried.
2. (Clean and rinse lumens of hollow instruments before wiping with this product [or product name]).
3. Wipe all hard non-porous surfaces of instruments with this product {or product name} until thoroughly wet.
4. Let stand for 1 minute[¶] [at room temperature].
5. [Wipe surfaces dry (or rinse). (If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]. (Allow to air dry [Let air dry]).]

[¶]Non-Critical medical devices are items that come in contact only with intact skin.

For Use as a Tuberculocide:

1. Pre-clean visibly soiled areas.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe hard, non-porous environmental surfaces.
3. Allow the surface(s) [Allow this product [or product name] to penetrate and] to remain visibly wet for 1 minute. [Use additional towelettes if needed to ensure the correct wet contact time.]
4. [Wipe surfaces dry (or rinse). (If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]. (Allow to air dry [Let air dry]).]
5. Change towelette after each use [application].

{*Administrative note:* Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

***KILLS HIV AND HCV AND HBV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in [health care settings] [Hospitals, Nursing Homes, Dental Clinics] [and] [other] [settings] in which there is an expected likelihood of soiling of (inanimate) hard, non-porous surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/ body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS) Hepatitis B Virus (HBV), or Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV OR HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

Personal Protection: Disposable latex or vinyl gloves, gowns, face masks, and eye coverings must be worn during all cleaning of body fluids, blood, and decontamination procedures.

Cleaning Procedures: Blood/body fluids must be thoroughly cleaned from surfaces and objects before application of this product.

Contact Time: Allow surface to remain visibly wet for 30 seconds [to kill HIV-1, HBV and HCV] [use a 1 minute contact time to kill all [other] organisms listed on the label [Reference Sheet]].

Disposal of Infectious Material: Cleaning materials used that may contain blood/body fluids must be autoclaved and disposed of according to Federal, State, and local regulations for infectious waste disposal.

To Kill Fungi [Yeast]:

1. Pre-clean visibly soiled areas.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe hard, non-porous environmental surfaces.
3. Allow surface[s] [Allow this product [or product name] to penetrate and] to remain visibly wet for 1 minute. [Use additional towelettes if needed to ensure the correct wet contact time.]
4. [Wipe surfaces dry (or rinse). (If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]. (Allow to air dry [Let air dry]).]

To Kill *Trichophyton interdigitale* [(Athlete's foot fungus)]:

1. Pre-clean visibly soiled areas.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe hard, non-porous environmental surfaces.
3. Allow surface[s] [Allow this product [or product name] to penetrate and] to remain visibly wet for 30 seconds. [Use additional towelettes if needed to ensure the correct wet contact time.]
4. [Wipe surfaces dry (or rinse). (If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]. (Allow to air dry [Let air dry]).]

To Kill Canine Parvovirus [CPV]:

1. Pre-clean surfaces.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe hard, non-porous environmental surfaces.
3. Allow surface[s] [Allow this product [or product name] to penetrate and] to remain visibly wet for 5 minutes. (Use additional towelettes if needed to ensure the correct wet contact time.)
4. [Wipe surfaces dry (or rinse). (If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]. (Allow to air dry [Let air dry]).]

To Kill Foot and Mouth Disease virus [FMDv]:

1. Pre-clean visibly soiled areas.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe hard, non-porous environmental surfaces.
3. Allow surface[s] [Allow this product [or product name] to penetrate and] to remain visibly wet for 3 minutes. (Use additional towelettes if needed to ensure the correct wet contact time.)
4. [Wipe surfaces dry (or rinse). (If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]. (Allow to air dry [Let air dry]).]

To Kill Mold and Mildew:

1. Pre-clean surfaces.

{*Administrative note*: Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe product onto hard, non-porous surfaces until completely wet.
3. Allow the surface[s] to remain visibly wet for 4 (four) minutes. Allow [Let] to air dry [or wipe surface dry] [or rinse]. [If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]].

To Control Mold and Mildew:

1. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe product onto hard, non-porous surfaces until completely wet.
2. Allow [Let] to air dry [or wipe surface dry] [or rinse].
3. [If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]].
4. Repeat application weekly or when growth reappears.

To Kill [or] [control] Mold and Mildew [for Remediation and Restoration]:

1. Pre-clean surfaces.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe product onto hard, non-porous surfaces until completely wet.
3. Allow the surface to remain visibly wet for 4 (four) minutes.
4. Allow [Let] to air dry [or wipe surface dry] [or rinse]. [If streaking is observed wipe with a clean damp [cloth or] [microfiber cloth or] [paper towel]].
5. Repeat application weekly or when growth reappears.

[[For] Use in] Mold Remediation:

1. Repair the source of moisture (stop leaks, etc) to eliminate conditions favorable to mold growth such as moisture ingress and high relative humidity and to prevent further mold growth and/or re-growth.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe the wall with a clean cloth wetted with clean water or with this product [or product name]. Dispose of the cleaning cloth in a well-sealed plastic bag.
3. Wipe the contaminated hard, non-porous surfaces with this product [or product name].
4. Allow to remain visibly wet for 4 minutes.

For use as a [Non-Food Contact Sanitizer]:

1. Pre-clean visibly soiled areas.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe hard, non-porous environmental surfaces.
3. Allow the surface to remain visibly wet for 10 seconds.
4. [Wipe surfaces dry (or rinse). (If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]. (Allow to air dry [Let air dry]).]
5. [Change towelette after each use [application].]

For Use as a Cleaner and/or Deodorizer:

[Pull] [Remove] towelette from dispenser [canister] [tub] and wipe hard, non-porous environmental surfaces. Allow to air dry [or rinse].

For Use as a [General] Cleaner/Degreaser:

Wipe surfaces (and let air dry).

GENERAL CLEANING: [Instructions will be shown with icons depicting the use of Oxy-1 (or alternate brand name) system].

1. Wipe floor with cloth and allow floor to dry
2. Cleaning steps: [Attach], [Mop], [Toss]
3. Scrubbing Strip usage: [Flip to scrub]
4. Use scrubbing strip to remove tough spots from dried-on food^a to scuffs.

^aOn non-food contact surfaces only.

[For Use as a Cleaner/Disinfectant in] Animal Housing Facilities [in] [Including] [such as] {insert veterinary use site from Table 8}:

{Administrative note: Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike through~~ text means it has been removed.

For cleaning and disinfecting the following hard, non-porous surfaces: [feeding and watering equipment, utensils, instruments, cages, kennels, [stables], [catteries], etc.]

1. Remove all animals and feeds from premises, animal transportation vehicles, crates, etc.
2. Remove all litter, [droppings] and [manure] [fecal matter] from floors, walls and surfaces of facilities [including pens, stalls, chutes, runs and barns] occupied [or traversed] by animals.
3. Empty all [troughs], [racks], [and] [other] feeding or watering appliances.
4. Thoroughly clean all surfaces with soap and detergent and rinse with potable water. [Pull] [Remove] towelette from [dispenser] [canister] [tub] and wipe hard, non-porous [environmental] surfaces.
5. All surfaces must remain visibly wet for 1 minute[¶].
6. Allow to air dry [or rinse [with potable water if necessary]].
7. [Ventilate buildings and other closed spaces.]
8. [Do not house animals or employ equipment until treatment has been absorbed, set or dried.]
9. Thoroughly scrub all [treated] [feed racks,] [mangers,] [troughs,] [automatic feeders,] [fountains,] and [waterers] with soap or detergent, and rinse with potable water before reuse.

Poultry Houses and Poultry Hatcheries:

For cleaning and disinfecting the following hard non-porous surfaces: {insert veterinary use site from Table 7 or Table 8}

1. Remove all poultry and feeds from premises, trucks, coops, and crates.
2. Remove all litter and droppings from floors, walls and surfaces of facilities occupied (or traversed) by poultry.
3. Empty all troughs, racks, and other feeding and watering appliances.
4. Thoroughly clean all surfaces with soap and detergent and rinse with potable water. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe hard, non-porous environmental surfaces.
5. All surfaces must remain visibly wet for 1 minute[¶].
6. Allow to air dry [or rinse [with potable water if necessary]].
7. [Ventilate buildings, coops, and other closed spaces.]
8. [Do not house poultry or employ equipment until treatment has been absorbed, set, or dried.]
9. Thoroughly scrub treated feed racks, troughs, automatic feeders, fountains, and waterers with soap or detergent, and rinse with potable water before reuse.

NOTE: Not for use on handling and restraining equipment such as leashes, muzzles, halters or ropes, as well as forks, shovels, and scrapers used for removing litter and manure.

For Cleaning and Disinfection of [Veterinary Clinics] {or insert veterinary use site from table 8}:

[For cleaning and disinfecting of hard non-porous surfaces (and equipment)]

1. Pre-clean visibly soiled areas.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe hard, non-porous environmental surfaces.
3. Allow surface[s] [Allow this product [or product name] to penetrate and] to remain visibly wet for 1 minute[¶]. (Use additional towelettes if needed to ensure the correct wet contact time.)
4. (Wipe surfaces dry) (or rinse). (If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]. (Allow to air dry [Let air dry])).
5. Rinse waterers and feeders with potable water before reuse.

For Use in Facilities used for Temporary Confinement of Animals:

[Use this product {or product name} to clean and disinfect (inanimate) hard, non-porous surfaces associated with facilities used for the temporary confinement of animals. Sites may include, but are not limited to, barns, sheds, stables, pens, cages, and associated access alleys or walkways. The product {or product name} may also be used to clean and disinfect equipment related to the maintenance of animals found at fairs, exhibitions, animal auction yards, animal show/boarding facilities, or other similar agricultural facilities designed for temporary housing of animals.]

1. Pre-clean visibly soiled surfaces and rinse with water.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe all hard, non-porous surfaces with this product {or product name}.
3. Allow product {or product name} to remain in contact with surfaces for 1 minute[¶].

For Cleaning and Disinfection of Vehicle Cabins [Cabs]:

{**Administrative note:** Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike through~~ text means it has been removed.

To disinfect the non-porous, hard surfaces of vehicle cabins [cabs]: Pre-clean visibly soiled areas. Pull [Remove] towelette [cloth] from dispenser [canister] [tub] and wipe hard, non-porous environmental surfaces. Allow surface(s) [Allow this product [or product name] to penetrate and] [to] remain visibly wet for 1 minute. [Use additional towelettes [cloths] if needed to ensure the correct wet contact time.] [Wipe] [surfaces dry] [or rinse]. [Allow to air dry] [Let air dry]. [If streaking is observed, wipe with a clean, damp [cloth or] [microfiber cloth or] [paper towel]]

To Clean and Disinfect Barber and Beauty/Manicure Instruments and Tools:

1. Pre-clean visibly soiled items.
2. [Pull] [Remove] towelette from dispenser [canister] [tub] and wipe all hard, non-porous surfaces of pre-cleaned combs, brushes, scissors, clipper blades, razors, manicure implements and other non-porous instruments with this product [or product name] until thoroughly wet.
3. Let stand for 1 minute[¶].
4. Rinse the surfaces thoroughly with water.
5. Let air dry before reuse.

For Use to Disinfect (inanimate) (environmental) hard, non-porous surfaces associated with Aquaculture [including] [vehicles], [nets], [boots], [brushes] [and] [other similar equipment]:

1. Remove all feeds from premises or vehicles.
2. Remove excess debris
3. [Pull towelette from dispenser [canister] [tub] and] Wipe all surfaces of washable hard, non-porous environmental surfaces until (thoroughly) wet. To disinfect, all surfaces must remain visibly wet for one (1) minute[¶].
4. Rinse the surfaces thoroughly with water.
5. Let air dry [or wipe dry].

For use as a disinfectant, (*virucide) on hard, non-porous (inanimate) (environmental) surfaces of [non-circulating] [and] [circulating] footbaths and (pedicure foot spa equipment):

1. **Drain:** Drain the water (from) (the unit) (foot spa) (basin) (bowl) and remove any visible debris.
2. **Clean:** Thoroughly clean all hard, non-porous (inanimate) (environmental) surfaces with this product [*or product name*], then rinse with water and drain.
3. **Apply:** [Pull towelettes from dispenser [canister] [tub] [and] Wipe all hard non-porous (inanimate) (environmental) surfaces [with the wipe].
Allow surface(s) to remain visibly wet for one [1] minute[¶].
4. **Rinse:** [(cleaned) (treated)] Surfaces with fresh water.

Note (important): Change towelettes frequently to avoid re-deposition of soil.

{Optional directions for refill pouch ONLY}

Refill Pouch

1. Open empty canister lid.
2. Open refill pouch [by cutting the top of the pouch] and remove saturated towelettes.
3. Place saturated towelettes in empty canister [tub] [invert the contents of the pouch into empty canister [tub]].
4. Replace canister [tub] lid.
5. [Remove the new lot number/[expiry date] from pouch and place it on top of the old lot number /[expiry date] sticker on the canister [tub]]

{Optional directions for starter kit ONLY}

Starter Kit:

1. Connect Poles
2. Attach
3. Clean
4. Remove when dirty

{Packaging Related Statement – Appropriate language will be selected from this list, depending on packaging.}

{For canisters and buckets}

To open the package [To pull wipes through lid]: Remove lid. [Then remove and discard induction seal [from the canister]]. Pull out corner of first sheet [wipe] from centre of roll and push ½" to ¾" (1 to 2 cm) through "X" opening on lid. DO NOT

{*Administrative note*: Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

push finger through opening. Replace lid and pull wipe at an angle. The next sheet pops up automatically. When finished, close [lid] [flap] to retain moisture.

{*For re-fill pouch*}

Refill [Pouch]:

[[Look for] Product expiration date on back of refill pouch.] Open [Remove lid from] empty canister [bucket]. [Place] [Put] [Insert] refill [pouch] into bucket [canister]. [Use scissors] [to] [Cut] [Rip] [Tear] [Pull] [Open] refill pouch. DO NOT remove wipe roll from refill pouch. Pull out corner of first sheet [wipe] from centre of roll and push ½" to ¾" (1 to 2 cm) through "X" opening on lid. DO NOT push finger through opening. Replace lid and pull wipe [towel] at an angle. The next [sheet] [wipe] [towel] pops up automatically. When finished, close [lid] [flap] to retain moisture. If bucket becomes contaminated [or] [damaged] replace immediately. Replace bucket immediately if tabs below [lid] [flap] become damaged. {or}

Refill [Pouch]:

Open [Remove lid from] empty canister [bucket]. [Use scissors] [to] [Cut] [Rip] [Tear] [Pull] [Open] refill pouch. [Remove the new lot number / [expiry date] from pouch and place it on top of the old lot number / [expiry date] sticker on the canister]. Place saturated towelettes in empty canister [by] [invert[ing] contents of the pouch into empty canister]. Pull out corner of first sheet [wipe] from centre of roll and push ½" to ¾" (1 to 2 cm) through "X" opening on lid. DO NOT push finger through opening. Replace lid and pull wipe [towel] at an angle. The next [sheet] [wipe] [towel] pops up automatically. When finished, close [lid] [flap] to retain moisture. If bucket becomes contaminated [or] [damaged] replace immediately. Replace bucket immediately if tabs below [lid] [flap] become damaged.

{*For re-sealable flatpack, transport pack or softpack*}

To Open Package:

1. Open [click-lock] [closure] and peel back -or- remove seal.
2. Pull out [one] wipe [sheet] [at a time] [, use as directed].
3. [When finished,] [firmly] close [click-lock] [closure] to retain moisture [to prevent drying between uses].

{*For single count pack*}

To Open Package: Tear open packet. Unfold and use as directed. Discard after single use.

Not for personal cleansing. This is not a baby wipe. [Do not flush in toilet].

{Optional directions for refill pouch ONLY}

Re-sealable soft pack

1. Remove the seal (from the place shown on the label)
2. Remove one sheet at a time.
3. Close the seal (to prevent drying between uses)

STORAGE AND DISPOSAL:

PESTICIDE STORAGE: Store this product in a cool [dry] area, [away from direct sunlight and heat].

Do not contaminate water, food or feed by storage and disposal. {*Note to reviewer: this statement will be included for non-residential use products only.*}

For [canisters], [buckets], [tubs] [flat-packs], [transport-packs], [soft-packs], [lids], [flaps] [or] [pouches] that are **not intended to be refilled** by consumer.

CONTAINER HANDLING: Dispose of wipe in trash after use. Do not flush. Tightly close lid [flap] between uses to retain moisture. Non-refillable container [dispenser]. Do not reuse or refill this container [dispenser]. Offer empty container [dispenser] for recycling. If recycling is not available, discard container [dispenser] in trash.

For [canisters], [buckets], [flat-packs], [transport-packs], [soft-packs], [lids], [flaps] or [pouches] **intended to be refilled** by the consumer

{**Administrative note:** Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

CONTAINER HANDLING: Dispose of wipe [cloth] in trash after use. Do not flush [in toilet]. Tightly reseal [canister] [bucket] [pouch] [soft pack] [flat pack] [label] [lid] [flap] between uses to retain moisture. Do not reuse or refill except as described in the directions for use. Refill only with this product [-or- product name]. Offer empty container for recycling. If recycling is not available, [wrap and] discard [of empty container] in trash.

{Reference to SDS will not be on the residential use products}

[USER(S) SHOULD READ S.D.S. BEFORE USING PRODUCT]

[SDS is available on line at {the URL link for SDS}]

[SDS Ref. No. xxxxxxxxxx]

{*Worker Protection Standard: Spanish Statement*}

Si no puede leer en ingles, pregunte a su supervisor sobre las instrucciones de uso apropiadas antes de trabajar con este producto.(If you cannot read English, ask your supervisor to explain the appropriate directions for use before working with this product)

†AHP Technology and Designs are trademarks of Diversey, Inc. [used under license by {the name of company}]



{**QUALIFIER Table**}

{*Note to reviewer: Qualifiers connect to the corresponding symbol. All organisms do not have to be listed on the printed labels, but the listing must be consistent with EPA guidance requirements for germs claims.*}

†	†Hepatitis C Virus
††	††Hepatitis B Virus
∞	{ <i>Note to reviewer: This qualifier leads to a list of cold and flu viruses claimed on product label</i> } [∞] Parainfluenza Virus Type 3 [(ATCC VR-93)] [∞] Respiratory Syncytial Virus (RSV) [(ATCC VR-26)] [∞] Rhinovirus Type-37 [(ATCC VR-1147)] [∞] Human Coronavirus [(ATCC VR-740)] [∞] Influenza A Virus (Hong Kong Strain) [(ATCC VR-544)] [∞] Influenza B Virus [(ATCC VR-823)] [∞] Pandemic Influenza A Virus H1N1
‡	‡[on]Non-Food contact surfaces [of]
±	±[on]Hard, non-porous surfaces [of]
*	*{ <i>List of Viruses from Table 2: Virucidal Claims</i> }
**	**Bactericidal claims listed on (product) label
***	***Bactericidal, Fungicidal and Virucidal claims listed on (product) label
†	[† <i>Mycobacterium bovis</i>] [Tuberculocidal [BCG]] (TB [†] stands for BCG)
◇	◇{ <i>List of *Parvoviruses from Table 2 Virucidal Claims</i> }
¶	¶[Use a 3 minute contact time against Foot and Mouth Disease virus [FMDv].] [Use a 5 minute contact time against Canine Parvovirus [CPV].] { <i>Note to reviewer: This qualifier will only be used when FMDv or Canine Parvovirus is claimed on the label</i> }
++	++This product [or product name] removes the following [non-living] allergens: [cat dander], [cockroach matter], [cockroach debris], [cockroach particles], [dog dander], [dust mite debris], [pet dander], [pollen particles]
1	¹ When used according to cleaning use directions [instructions]
2	² { <i>List of antibiotic resistant bacteria from Table 1: Bactericidal Claims</i> }
3	³ When used according to disinfection use directions [instructions]
4	⁴ When used according to sanitization use directions [instructions]

{**Administrative note:** Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. Strike-through text means it has been removed.

5	⁵ Disinfects hard non-porous floor surfaces, including sealed wood floors, sealed tile and sealed grout (DO NOT use on unfinished, oiled or waxed wooden boards, non-sealed tiles or carpeted floors because they may be water sensitive) (for use on finished wood)
6	⁶ SARS-Related Coronavirus 2 [Human Coronavirus SARS-CoV-2 (the virus that causes COVID-19 [‡])] [SARS-CoV-2] [Strain USA-WA1/2020] [BEI NR-52281]

{List of soils, surface materials, hard non-porous surfaces and use site for use in Institutional products}

{Table 5 – Institutional Soils}

- Blood
- Body oils
- Dead skin
- Dirt
- Feces [Fecal matter]
- Finger prints
- Food residue[s] –or-soil[s]
- Grime
- Laboratory stains
- Mold stains
- Other common soils –and/or- stains
- Other organic matter
- Plant based resins
- Resin(s)
- Soils
- Urine

{Table 6 – Institutional Surface Materials}

- Aluminum
- Chrome
- Crystal (non-food contact items)
- [Baked] Enamel [Surfaces]
- Fiberglass fixtures
- Glass
- Glazed ceramic
- Glazed porcelain
- Glazed tile
- Laminated surfaces
- Linoleum
- Metals [except copper and brass]
- Mirrors
- Painted surfaces
- Plastic
- Plastic laminate
- Plastic surfaces
- Sealed fiberglass [polyester, vinylester and epoxy resins]
- Stainless steel
- Synthetic marble
- Vinyl

{Table 7 – Institutional Hard Non-porous surfaces}

Medical and Dental:

Hard non-porous surfaces of

- amalgamators -and/or- dental curing lights
- anesthesia machines
- autoclaves
- baskets
- bedpans
- bedrails
- [bedside] commodes
- bedside tables
- blood pressure cuffs
- blood pressure (BP) monitors
- BSL hoods
- cabinet [handles]
- call boxes
- CAT -or- Computerized Axial Tomography equipment
- carts
- centrifuges
- chairs
- clipper [handles]
- closet [handles]
- computer peripherals
- computer screens
- computer tables
- cords
- counters
- CPAP (Continuous Positive Airway Pressure) mask
- [crash] [emergency] carts
- dental countertops
- dental operator surfaces
- dentists' -or- dental chairs
- dental unit instrument trays

{*Administrative note:* Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

- desks
- desktops
- diagnostic equipment
- [exterior of] dialysis machines
- docking stations
- [exam -or- examination] tables
- external surfaces of [medical] equipment -or- [medical] equipment surfaces
- [external] [surfaces of] ultrasound transducers [- and/or- probes]
- exterior of pipes
- footboards
- fume Hoods
- glucometers
- gurneys
- handheld [electronic] devices
- headboards
- hose ends
- [hospital - or – patient] bed[s] [springs][railings][frames][lining]
- Infant incubators and care cribs
- IV poles [stands]
- IV pumps (exterior surface)
- keyboards
- laptops
- light lens covers
- loupes
- Manikins used in CPR training^Δ
- mammography equipment
- [Mayo] [instrument] stands
- medication carts
- microscope knobs
- [computer] mouse
- mobile devices
- mobile electronic equipment
- mobile workstations
- mouse pads
- [Surfaces of] MRI -or- Magnetic Resonance Imaging equipment
- nurse-call [devices] [buttons] [and cords]
- operating room tables and lights
- operatory light switches
- otoscope handles
- oxygen hoods
- overbed tables
- paddles
- pagers
- patient chairs
- patient monitoring equipment
- patient support and delivery equipment
- phlebotomy trays
- physical therapy (PT) equipment surfaces
- phone cradle
- plastic mattress covers
- power cords
- pulse oximeters
- PVC tubing
- Reception counters –or- desks –or- areas
- remote controls
- respirators
- respiratory therapy equipment
- scales
- sequential compression devices
- shower fixtures
- side rails
- slit lamps
- small surfaces
- spine backboards
- stethoscopes
- stools
- stretchers
- support bars
- tablet PCs
- tack
- toilet handholds
- traction devices
- ultrasound transducers
- walls [around toilet] [in patient rooms]
- wash basins
- wheelchairs
- work stations
- X-ray equipment

Veterinary:

Hard non-porous surfaces of

- animal equipment
- animal crates
- automatic feeders^Δ
- boots
- buckets
- brushes
- cages
- cage rack
- carts
- diagnostic equipment
- dog runs
- examination tables
- exterior surfaces of fountains
- external surfaces of [veterinary] equipment
- exterior surfaces of watering appliances
- farm vehicles (trucks)

{*Administrative note*: Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

- feed racks^Δ
- fencing
- foot baths
- gates
- gathering [guiding] equipment
- grooming equipment
- handling equipment
- harnesses & leashes
- imaging equipment

- pens
- reception counters -or- desks -or- areas
- rehabilitation equipment
- scales
- [wash] stalls
- transportation crates
- troughs^Δ
- veterinary care surfaces
- underwater treadmill

Professional Beauty:

Hard non-porous surfaces of

- [non-circulating] [circulating] footbaths [surfaces]
- countertops
- esthetic tools and equipment
- facial tables

- implements
- massage tables
- pedicure foot [spas] [basins] [baths]
- spa equipment
- workstations

Greenhouse (Horticulture):

Hard non-porous surfaces of

- belly pans
- ceilings
- drying tables
- [extraction] [trimming] equipment
- floors
- greenhouses
- handheld tools

- lighting
- processing tables
- racking
- scales
- trimmers
- walls

Miscellaneous

Hard non-porous surfaces of

- bath mats (e.g., rubber or plastic)
- [bathroom] fixtures
- [bathroom] faucet[s]
- biosafety cabinets
- blinds
- BP Monitors
- cabinets
- ceilings
- cell[ular] -or- mobile phones
- chairs
- changing tables
- chemotherapy hoods
- [children's] furniture
- closets
- computer keyboards & mice
- counters -or- countertops
- crawl spaces
- dashboard
- desk[s] [tops]
- [diaper -or- infant] changing [tables]
- drain troughs
- drip pans
- [restaurant^Δ] tables
- door[s] [and drawer] [handle[s]] [frames]

- empty garbage bins
- [filing] [medicine] cabinets
- [finished -or- painted] walls
- finished basements
- floors [around toilets]
- garbage storage
- [glass] [stainless steel] equipment
- furniture
- [hand]railings -or- rails
- areas -or- surfaces behind and under sinks
- HVAC Systems
- [instrument] stands
- lamps
- light fixtures -or- switches -or- panels
- lockers
- [medicine] cabinets
- non-circulating footbath surfaces
- phones
- rear view mirror
- room benches
- sealed [concrete] floors
- sealed drywall [and gypsum board]
- sealed walls
- shelves [and drawers]

{*Administrative note:* Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

- steering wheel
- stick shift
- tables
- tabletops
- trusses
- wall cavities
- wallboard
- [washable] walls
- windowsills
- Telephones

^ΔRinse with potable water after disinfecting.

{Table 8 - Institutional use sites}

Medical and Dental:

Hard non-porous surfaces of

- Acupuncture Clinics [Facilities]
- Ambulances -or- [Emergency Medical] Transport Vehicles
- Ambulatory Surgical Centers (ASC)
- Anesthesia Rooms or Areas
- [Assisted Living -or- Full Care] Nursing Homes
- Biopharmaceutical Areas
- Carts
- CAT Lab[oratories]
- Central Service Areas
- Central Supply Rooms -or- Areas
- Chiropractic Clinics
- Cleanrooms
- Colleges
- Containment Areas [(CDC) Biosafety labels 1 to 4]
- Critical Care Units -or- CCUs
- Critical Environments
- Dental Operatories
- Dental -or- Dentist's [Facilities] [Offices] [Clinics]
- Dialysis Clinics [Facilities]
- Doctors' Offices
- Donation Centers [blood] [plasma] [semen] [milk] [apheresis]
- Emergency Rooms -or- ERs
- Emergency Care Service (EMS) Areas (including Fire, Police and Ambulance)
- Examination Rooms
- Eye Surgical Centers
- FDA Regulated Sites
- Foot Care Nursing
- GMP Facilities
- Health Care Settings -or- Facilities
- Homecare
- Hospices
- Intensive Care Units -or- ICU[s] [areas]
- Isolation Areas
- Laundry Rooms
- Compound Pharmacy Laboratories
- [Life Sciences] [Laboratories] [Facilities]
- Long Term Care Facilities
- Massage Clinics [Facilities]
- [Medical] Clinics
- Medical Devices
- Medical Facilities [Offices] [Clinics]
- Medical -or- Physician's -or- Doctor's Offices
- Midwifery [Facilities] [Offices] [Clinics]
- Naturopathic -or- Homeopathic Clinics
- Newborn -or- Neonatal [Nurseries] [Intensive Care] Units [NICU]
- Nursing -or- Nurses' Stations
- Operating Rooms
- Ophthalmic Offices
- Ophthalmology Clinics
- Orthopedic Offices [Clinics]
- Otolaryngology Offices [Clinics]
- Out-Patient [Surgical Centers (OPSC)] Clinics [Facilities]
- Patient Areas
- Patient Restrooms
- Patient Rooms
- [Pediatric] Examination Rooms -or- Areas
- Pediatric Intensive Care Units [PICU]
- Pharmacies
- Physical Therapy Rooms -or- Areas
- Physicians' Offices
- Podiatry Offices [Clinics]
- Psychiatric Facilities
- General Admission Areas
- Radiology -or- X-Ray Rooms -or- Areas
- Recovery Rooms
- Rehabilitation Centers
- Respiratory Therapy Rooms [Areas] [Clinics] [Facilities]
- Schools
- Sleep Therapy Clinics [Facilities]
- Stomatology Hospitals

{*Administrative note:* Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike through~~ text means it has been removed.

- Surgery Rooms -or- Operating Rooms - or-ORs

- Waiting Rooms -or- Waiting Areas

Veterinary:

Hard non-porous surfaces of

- Agriculture farms and barns
- Agriculture & Livestock Buildings
- Animal [Containment] [Enclosures]
- Animal housing facilities and areas
- Animal shelters and kennels
- Animal [Pet] Housing [Kennels] [Facilities]
- Animal Examination Rooms
- Animal Holding Areas
- Animal Hospitals
- Animal Life Science Laboratories
- [Animal -or- Pet] Grooming Facilities
- Animal rescue Shelters [Facilities]
- Animal Transportation Vehicles
- Animal Quarters
- Adoption Facilities
- Barns [premises]
- Breeding Establishments
- Breeder [houses] [facilities] [premises]
- Broiler [houses] [facilities] [premises]
- Boarding Shelters
- Bovine Farms
- Catteries
- Cages
- Central Animal Facilities
- Dairy Farms
- Diagnostic Laboratories
- Equine Farms [Facilities]
- Exam Rooms
- [Exotic] [Veterinary] [Veterinarian] [Animal] [Clinics] [Import] [Export] [Facilities]
- Farms
- Farm premises
- Fish Farms [Facilities]
- Hatcheries
- Imaging Suites
- Kennels
- Laboratories
- Lab Animal Research Facilities
- Large Animal Facilities
- Life Science [Research] Facilities [Laboratories]
- Livestock -and/or- Swine -and/or- Poultry Facilities
- [Non-human] Primate Areas
- Pet [Areas] [Quarters]
- [Pet] Boarding [Facilities]
- Pet Shops -or- Stores
- Pet [Animal] Care Facilities
- Pet [Animal] Housing
- Poultry [Hatcheries] [Broiler Houses] [Breeder Houses] [Facilities] [Premises] [Houses]
- Poultry [Boot Wash] [Vehicle Wash]
- Practices
- Quarantine [Facilities]
- Radiology Clinics [Facilities]
- Reception Areas
- Rodent Facilities
- Small Animal Facilities
- Surgical Suites
- Swine [Farrowing] Nursery [Finisher]
- Swine [Boot Wash]
- Swine [Vehicle Wash]
- Tack Shops
- Therapeutic Laser Clinics [Facilities]
- Transportation vehicles
- Treatment facilities
- Truck Wash [Facility]
- Ultrasound Clinics [Facilities]
- Vehicle Wash
- Veterinary Clinics [Facilities] [practices] [hospitals] [Offices] [Waiting Rooms]
- Veterinary Examination Rooms
- Veterinary X-ray Rooms
- Veterinary Operating Rooms
- Veterinary Laboratories
- Vivariums
- Zoos

Professional Beauty:

Hard non-porous surfaces of

- Barbershops
- Beauty [Salons] [Shops]
- Clinics
- Facial Spas
- Footbaths
- Hair Salons
- Hairstyling Establishments
- Health Club[s] [Facilities]
- Lash bars
- Massage Therapy Clinics
- Medesthetic Clinics [Treatment Room]
- Nail Salons
- Pedicure Salons
- Spas [and] [Massage/Facial] [Hair/Nail] [Salons]
- Tanning [Facilities] [Spas]
- Tattoo Parlors

{*Administrative note*: Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

- Treatment rooms

- Wellness centers

Mold and Mildew Remediation and Restoration:

Hard non-porous surfaces of

- Crime scene
- Crisis (Crime) (and) (Trauma) (Fire) Scene
- Fire clean-up
- Grey or Black Water [Damage][Sewage] Areas
- HVAC Systems
- Mold and mildew remediation [sites]
- Professional restoration [sites]

- Restoration Industries [including water, fire, and mold restoration]
- Roofs
- Sewage backflow sites [areas]
- Trauma
- Water/flood damaged [restoration] areas
- Water and flood damage sites

Greenhouse (Horticulture):

Hard non-porous surfaces of

- Clean Rooms
- [Drying] [Curing] Rooms
- Extraction Rooms
- Greenhouses
- Growing Facilities
- Hydroponic Growing Areas

- Indoor Grow Rooms
- Licensed Producers
- Production Facilities
- Packaging [Processing] Areas
- Storage Areas

Miscellaneous

Hard non-porous surfaces of

- Airplanes [Airports] [Aviation]
- Amusement parks
- Attics
- Basements
- Bowling Alleys
- Buses
- Non-food handling areas of Cafeterias
- Campers
- Campgrounds
- Cellars
- Churches
- Colleges
- Convenience Stores
- Convention Centers
- Correctional Facilities
- Crawl Spaces
- Cruise Ships
- Crisis [Crime and Trauma] [Fire] Scenes
- [Home] Day Care Centers
- Dorms
- Dormitories
- Educational Institutions
- Elevators
- Factories
- Non-food handling areas of Fast Food Restaurants
- Non-food handling areas of Food Service Establishments [Restaurants]
- Funeral Homes

- Gas Stations
- Greenhouses
- Grey or Black Water [Damage] [Sewage] Areas
- Health Club Facilities
- Hotels
- HVAC Systems
- Industrial Facilities
- Institutional [Bathrooms]
- [Institutional] Laundromats
- Institutions
- Libraries
- Living Rooms
- Locker Rooms
- Lodging Establishments
- Lounges
- Malls
- Manufacturing Facilities
- General Merchandise Stores, Discount Retailers
- Military Installations
- Mold and Remediation Sites
- Mortuaries
- Motels
- Office Buildings
- [Police -and/or- Fire] Vehicles
- Public Areas
- Public Facilities
- Public Restrooms
- Public Telephone[s] [Booths]
- Public Transportation Systems

{Administrative note: Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike-through~~ text means it has been removed.

- Recreational [Athletic] Centers -or- Facilities
- Rest Stops
- Restoration Industries [Including Water, Fire and Mold Restoration]
- Restrooms -or- Restroom Areas
- Retail and Wholesale Establishments
- Retail Businesses
- [Retail] [Wholesale] Florists
- RVs
- School Buses
- Schools
- Sewage Backflow Areas
- School Buses
- Schools
- Shelters
- Shippers
- Shopping Centers
- Shops
- Shower Rooms
- Sports Arenas
- Sports Stadiums
- Subways
- Non-food handling areas of Supermarkets
- Transportation Terminals
- Transport vehicles
- Trailers
- Trains
- Universities
- Warehouse Clubs
- Water [/Flood] Damaged [Restoration] Areas

{List of surface materials, hard non-porous surfaces and use site for use in the Residential Use products}

{Table 9 – Residential Soils}

- Baby food
- Bacon fat
- Bathtub rings
- Beverage stains
- Baked/burnt on food
- Coffee
- Cosmetic stains
- Crayon / pencil marks
- Food stains
- Greasy [oily] food spills [splatters]
- Fingerprints
- Grime
- Hair dye
- Handprints
- Heel/Scuff marks
- Juice stains
- Pet stains
- Smudges
- Soap scum
- Sticky Food Messes – or – spills
- Tough soils
- Wine stains

{Table 10 – Residential Surface Material}

- Chrome
- Crystal (non-food contact items)
- Enamel
- Fiberglass fixtures
- Formica[^]
- Glass
- Glazed ceramic
- Glazed porcelain
- Glazed tile
- Linoleum
- Metals [except copper and brass]
- Mirrors
- Plastic laminate
- Plastic surfaces
- Sealed fiberglass [polyester, vinylester and epoxy resins]
- Sealed granite
- Sealed grout
- Stainless steel
- Synthetic marble
- Vinyl

[^]Formica is a registered trademark of the Formica Corporation

{Table 11 – Residential Hard Non-porous surfaces}

Hard, non-porous surfaces of

- Outer surfaces of air vents
- Appliances exterior
- Athletic equipment
- Bathrooms
- Bathroom fixtures
- [Bathroom] [kitchen†] faucet
- Bathroom doorknobs
- Bathtubs
- Bed frames^{±±}
- [Behind and under] counters

*{Administrative note: Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. Strike through text means it has been removed.*

- [Behind and under] sinks
- Blood Pressure monitors
- Burner trays
- Cabinets
- Cages
- Cell[ular] – or digital phones
- [Hard surfaces of] Chairs (not including high chairs)
- Changing tables⁺⁺
- Closets [handles] ⁺⁺
- Computers
- Computer keyboards
- Computer mouse
- Computer peripherals
- Computer screen
- cords
- Countertop⁺⁺
- Cupboard
- Desk[s][tops]
- Diaper pails
- Door[±] [knobs] [handles] [frames] ⁺⁺
- Exercise equipment
- Exterior of toilet surfaces
- Exterior of pipes
- Footboards
- Finished hardwood
- Fixture [(chrome) (stainless steel)]
- Floor (glazed tile)
- Floor (vinyl)
- Furniture [hard surfaces]⁺⁺
- Garbage [trash] can [bin]
- Glazed ceramic tile
- Glazed porcelain tile
- [Hand] railings
- [Hard[±] surfaces of] car [boat] [trailer] [interior]
- Infant cribs⁺⁺
- [Inside of] drawers
- Grill exterior
- Laminated surfaces
- Lamps
- Light fixture
- Medicine cabinet
- Metal blinds
- Metal surfaces
- Microwave exterior
- [Hard[±], surfaces of] [outdoor]Patio furniture
- Oven exterior
- Packing benches⁺⁺
- Phone cradles
- Piano keys
- Picnic table⁺⁺
- Plastic laundry basket
- Recycling bins
- Refrigerator / Freezer exterior
- Remote controls
- Sealed granite
- Shelves⁺⁺
- Shower[door][curtain]
- Showers [stalls]
- Sink [basin]
- Stainless steel surfaces
- Stereo consoles
- Storage bins
- Stool⁺⁺
- Strollers
- Tables⁺⁺
- Tabletop⁺⁺
- Telephone
- Tiled walls
- Toilet handholds
- Toilet rims [seats]
- Toys[±]
- Training toilets
- Tools
- [Bathroom] vanity
- Vinyl
- [Washable] walls
- Wash basins
- Walkers
- Washing machine exterior
- Windows
- Work bench

{Note to Review: This qualifier is not necessary if the phrase preceding the surface specifies hard, non-porous surfaces} [±]Hard, non-porous surfaces of

⁺⁺ Not for use on unfinished wood surfaces

{Table 12 - Residential/Household use sites}

- Basement
- Bathroom
- Bedroom
- Daycare settings
- Dorm rooms
- Dining room
- Exercise room
- Family room



{Administrative note: Note: Parenthesis () or [bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike through~~ text means it has been removed.

- Garage
- Garbage room
- Homecare settings
- Household settings
- Kitchen[‡]
- Laundry room
- Living room
- Nursery [room]
- Nursing homes
- Outdoor [patio] furniture
- Pet areas
- Playroom
- Residential settings
- Restroom
- Storage room
- Tool shed
- Washrooms


[‡]Non-food contact surfaces

{*Administrative note:* Note: [Bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Parenthesis () is for alternate text or text that will appear in parenthesis. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike through~~ text means it has been removed.


{*Graphic options*}



{*Administrative note:* Note: [Bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Parenthesis () is for alternate text or text that will appear in parenthesis. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike through~~ text means it has been removed.



{*Administrative note:* Note: [Bracketed text] is optional text, it might not appear on the printed label. Text appearing in {braces} is for administrative purposes only and will not appear on the printed label. Parenthesis () is for alternate text or text that will appear in parenthesis. Variables such as CAPITALIZATION, plurals or singular, **bold** and *italics* are optional unless specified in 40 CFR. Listed uses and use sites can be used partially or in their entirety. Underlined text is new. ~~Strike through~~ text means it has been removed.



1 MINUTE

{Note to reviewer: this graphic will only be used on product labels where all organisms are claimed at a contact time of 1 minute or less}

