

UNITED STATES ENVIRONMENTAL PROTECTION AGENCY WASHINGTON, DC 20460

December 21, 2017

OFFICE OF CHEMICAL SAFETY AND POLLUTION PREVENTION

Georgia Anastasiou Hu-Friedy Manufacturing Co., LLC c/o Lewis & Harrison 122 C St. NW, Suite 505 Washington, DC 20001

Subject: Notification per PRN 98-10 – Minor Label Revisions Product Name: HF ANTIMICROBIAL DUWL IODINE RESIN UNIT EPA Registration Number: 75372-1 Application Date: 11/20/2017 Decision Number: 536685

Dear Ms. Anastasiou:

The Agency is in receipt of your Application for Pesticide Notification under Pesticide Registration Notice (PRN) 98-10 for the above referenced product. The Antimicrobials Division (AD) has conducted a review of this request for its applicability under PRN 98-10 and finds that the action requested falls within the scope of PRN 98-10.

The label submitted with the application has been stamped "Notification" and will be placed in our records.

Should you wish to add/retain a reference to the company's website on your label, then please be aware that the website becomes labeling under the Federal Insecticide Fungicide and Rodenticide Act and is subject to review by the Agency. If the website is false or misleading, the product would be misbranded and unlawful to sell or distribute under FIFRA section 12(a)(1)(E). 40 CFR 156.10(a)(5) list examples of statements EPA may consider false or misleading. In addition, regardless of whether a website is referenced on your product's label, claims made on the website may not substantially differ from those claims approved through the registration process. Therefore, should the Agency find or if it is brought to our attention that a website contains false or misleading statements or claims substantially differing from the EPA approved registration, the website will be referred to the EPA's Office of Enforcement and Compliance.

If you have any questions, you may contact Donna Kamarei at (703)347-0443 or via email at Kamarei.donna@epa.gov.

Sincerely,

Warda

Demson Fuller, Product Manager 32 Regulatory Management Branch II Antimicrobials Division (7510P) Office of Pesticide Programs

HF Antimicrobial DUWL Iodine Resin Unit

For Use As An Antimicrobial Agent in Dental Unit Water Lines

Active Ingredient:	
lodine [*]	46. <mark>0</mark> %
Other Ingredients	54. <mark>0</mark> %
Total	
*Bound to polystyrene divinyl benzene quaternary ammonium	
anion exchange resin.	
This product contains 0.5 lb of iodine (I_2) p	er pound

KEEP OUT OF REACH OF CHILDREN CAUTION

FIRST AID

If in eyes:

- Hold eye open and rinse slowly and gently with water for 15-20 minutes.
- Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.
- Call a poison control center or doctor for treatment advice.

If on skin or clothing:

- Take off contaminated clothing.
- Rinse skin immediately with plenty of water for 15-20 minutes.
- Call a poison control center or doctor for treatment advice.

If swallowed:

- Call a poison control center or doctor immediately for treatment advice.
- Have person sip a glass of water if able to swallow.
- Do not induce vomiting unless told to do so by a poison control center or doctor.
- Do not give anything by mouth to an unconscious person.

Have the product container or label with you when calling poison control center or doctor, or going for treatment.

[See [side][back][panel][insert][accompanying instructions] for additional precautionary statements, first aid and directions for use.]

NOTIFICATION

75372-1

The applicant has certified that no changes, other than those reported to the Agency have been made to the labeling. The Agency acknowledges this notification by letter dated:

12/21/2017

EPA Label 10/25/2017 Page **1** of **4**

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

CAUTION. Causes moderate eye irritation. Avoid contact with eyes, skin and clothing. Wear protective eyewear (goggles) and gloves when handling. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse.

ENVIRONMENTAL HAZARDS

The exposed iodine resin pellets are toxic to fish and aquatic organisms. Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge.

DIRECTIONS FOR USE

It is a violation of federal law to use this product in any manner inconsistent with its labeling.

HF UWL lodine Resin Unit is designed for use on dental unit water lines attached to dynamic dental instruments. The product reduces microbial contamination and odor-causing bacteria in dental unit water lines.

Prior to installing HF DUWL lodine Resin Unit, thoroughly clean dental unit water lines following the manufacturer's instructions or with an appropriate cleaning agent. Insert HF DUWL lodine Resin Unit into the inlet and outlet fittings of dental unit water lines. Once the product is in place, the supply water will pass through the iodinated resin bed, which will continuously release low levels of antimicrobial iodine. After installation, check water flow and verify that the unit is operational. Create a check-off sheet to monitor volume of water used through unit. Shut off water to the dental unit when not in use.

The lifetime of this product is 240 liters or 365 calendar days (if water usage records are not monitored) or if iodine levels falls to less than 0.5ppm. After this period or level is reached, immediately replace the spent HF Antimicrobial DUWL lodine Resin Unit product/unit with a new or replacement product/unit. To install the new unit, follow the instructions above for installing a new product.

[If HF Antimicrobial DUWL lodine Resin Unit is used in conjunction with a bottle water system] The bottle in water bottle systems must be emptied each night and the empty bottle replaced on the manifold. [Prior to] [or] [Before] daily use, wash the inside of the bottle and rinse thoroughly. Wipe down the outside of HF Antimicrobial DUWL lodine Resin Unit and make sure the intake end is not obstructed. Fill with fresh water each morning before its first use. Replace HF Antimicrobial DUWL lodine Resin Unit after 240 liters of water usage, if iodine level falls to less than 0.5ppm, or if records of usage are not kept, after 365 calendar days of installation. Note: Adherence to CDC/ADA guidelines for aseptic procedures including a 2 minute morning flush, a 20-30 second flush between patients, and hand piece sterilization must be continued after the installation of this product.

If package seal is broken or product appears dry, wrap in newspaper and dispose of it in the trash.

Troubleshooting

If the water flow after product installation is slow or inhibited, check to ensure that the product was installed properly or if the product has been plugged by debris. For bottle installation ensure that the air supply is turned on at the maximum level (45 psi). In addition, ensure that the product is at the bottom of the water bottle. If there is an air leak from fittings, contact your dental supplier for replacement parts or service. Discoloration of tubing from yellow to blue to green to black is normal and is not a concern. Initial discoloration of water is caused by release initial flush of excess iodine. Operate water until flow appears clear.

STORAGE & DISPOSAL

STORAGE: Store at room temperature in original packaging, away from children, pets and direct sunlight. Avoid excessive heat. Do not allow cartridge to freeze. Do not store in open or unlabeled containers.

DISPOSAL: Securely wrap cartridge in newspaper and discard in trash.

EPA Reg. No. 75372-EPA Est No. Net Contents:

[Manufactured] [Distributed][by][for] Hu-Friedy Manufacturing Co., LLC D/B/A Hu-Friedy 3232 N. Rockwell, St Chicago, IL 60618

Filter Label 4.75" x 0.375"

Product Name	KEEP OUT OF REACH OF CHILDREN	EPA Reg. No. 75372-1
Active Ingredient: Iodine 46%	CAUTION: See insert for directions & precaut	ions EPA Est. No.:

OPTIONAL MARKETING LANGUAGE

- No risk of allergic reaction (isotopic iodine contains no protein)
- Container has no silver or other hazardous waste disposal as easy as throwing it in the trash
- EPA registered as an antimicrobial dental unit water line product
- Eliminates the need for manual addition of biocides into water systems. This product does it automatically.
- Consistent and accurate dosing.
- [Reduces][Controls][Inhibits] bacterial contamination (odor and fouling bacteria) in dental unit water lines.
- Meets or exceeds ADA guidelines for dental unit water lines.
- Meets or exceeds [ADA recommendation of] < 200 CFU/ml for dental unit water lines.
- One year life with minimal human interaction
- Easily installed in just minutes
- Ability to be used with either municipal or independent reservoir
- No corrosion or etching of metal parts
- No effect on bonding agents
- Ability to be used with current dental unit water bottles
- Shocking not required on monthly (or weekly) basis
- No daily, weekly or monthly protocols
- Works with any size bottle
- Non-corrosive
- Shocking not required on first use (or optional)
- Same product can be used for inline water or independent reservoir
- Less than 5-minute change out time.
- Antimicrobial Cleaner
- Cleans and maintains dental unit waterlines
- Inhibits the growth of microbial contamination in dental waterlines
- Kills odor causing bacteria
- For the cleaning of microbial contamination in dental until waterlines
- You're on your way to cleaner water!

OPTIONAL REMINDER HANG TAG LANGUAGE

- Equipped With [product name]
- Reminder. Change your filter annually. The lifetime of this product is 240 liters or 365 calendar days (if water usage records are not monitored). Replace [product name] after 240 liters of water usage, or if records of usage not kept, after 365 calendar days of installation.
- Date Installed:
- Installed By:
- Install New [product name]. Date for New Filter:
- Re-order Information: Product Name: Part Code:
- To reorder, call 1-800-HU-Friedy or order through your dental distributor