UNITED STATES ENVIRONMENTAL PROTECTION AGENCY WASHINGTON, D.C. 20460

03/30/2011

GEPA United States Environmental Protection Office of Pesticide Programs

WAR 30 2011

A-dec, Inc 2601 Crestview Drive Newberg, OR 97132

79662 - 1

Attention: Georgia Anastasiou

Subject: ICX EPA Registration No: 79662-1 Amendment Dated: December 29, 2010

This will acknowledge receipt of your Amendment to update the container handling statement per PR Notice 2007-4 and make other minor changes to the product label, submitted under the provision of FIFRA Section 3(c)(9).

Proposed Amendment:

- To modify the "Note" section in the directions of use.

General Comments:

Based on a review of the submitted material, the following comments apply: Under the Direction in the "Note" section advise users to "Conduct a flush treatment before using the ICX". The Amendment dated December 29, 2010 is acceptable, and is in compliance with the PR Notice 98-10. A stamped copy of the label has been added to your files

If you have any questions concerning this letter, please contact Zebora Johnson at (703) 308-7080

Sincerely "

Marshall Swindell Product Manger (33) Regulatory Management Branch I Antimicrobial's Division (7501P)

affdec°

ICX®

A-dec[®] ICX[®] tablets are specially formulated to maintain dental unit waterlines and prevent accumulation of odor and foul taste bacteria.

ACTIVE INGREDIENTS:

Sodium Percarbonate	6.96%
n-Alkyl(60% C14 30% C16 5% C12 5% C18) dimethyl benzyl ammonium chloride	
n-Alkyl(68% C ₁₂ 32% C ₁₄) dimethyl ethylbenzyl ammonium chloride	
Silver Nitrate	
Other Ingredients	
Total	

[Note to reviewer: EPA approved the location of the active ingredients statement on the back panel of the label on 12/15/09 in order to accommodate the small container size]

Net weight of tablet: 212.4 mg / 0.008 oz (2L application), 80.5 mg / 0.003 oz (0.7L application) Net weight of product: 10.62 g / 0.375 oz (2L application), 4.03 g / 0.142 oz (0.7L application) – 50 tablets per box

KEEP OUT OF REACH OF CHILDREN

WARNING

PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals

Causes moderate eye irritation, wear protective eyewear. Avoid contact with eyes, skin or clothing. Wash hands thoroughly with soap and water after handling and before eating, drinking, chewing gum, or using tobacco.

DIRECTIONS FOR USE:

 ICX^{\oplus} is intended for use in general dental applications by providing water to the dental syringe, various handpiece ancillary devices, and cuspidor cupfill. When used according to these instructions for every bottle refill, optimum efficacy for control of odor and foul taste bacteria can be achieved.

1. Empty any remaining water left in the bottle.

2. Drop tablet into an empty dental unit water bottle (0.7L tablet in 0.7 liter bottle, 2 L tablet in 2 liter). Avoid touching the tablet with unprotected skin.

3. Fill bottle with water and connect to the dental unit.

4. Wait 120 seconds for the tablet to fully dissolve.

NOTE -This product is intended for use with potable water and in conjunction with regular dental unit water testing, following flush (shock) treatment with an antimicrobial product registered with the U.S.EPA and recommended by the manufacturer

FIRST AID:

If solid product gets in eyes, hold eye open and rinse slowly and gently with water 15 to 20 minutes. Remove contact lenses, if present, after the first five minutes, then continue rinsing.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage disposal. Store at room temperature out of direct sunlight; DO NOT store above 30°C (86°F). Nonrefillable container. Do not reuse or refill this container. Offer for recycling if available or dispose of empty outer foil pouch in the trash.

A-dec Inc., 2601 Crestview Drive, Newberg, Oregon 97132, USA Phone: 1-800-547-1883 Fax: 503-537-2843 EPA Establishment Number 84100 Made in USA EPA

ACCEPTED with COMMENTS EPA Registration States -1

MAR 30 2011

2010 - 12 86.0613.00

Under the Federal Insecticide, Funcicide, and Rodenticide Act as amended for the pesticide, registered under EPA Reg. No. 79662-1