



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
WASHINGTON, DC 20460

OFFICE OF CHEMICAL SAFETY
AND POLLUTION PREVENTION

May 31, 2022

Dr. Thad Overmyer
CEO of Anodia Systems, LLC
109 Larrimore Lane
Danville, Kentucky 40422
Electronic Transmittal: drovermyer@anodiasystems.com

Subject: Label Amendment – Add label use directions for Cleaning Dental Water Bottle and Waterlines
Product Name: Mint-A-Kleen
EPA Registration Number: 85298-2
Received Date: 1/7/2022
Action Case Number: 00337882

Dear Dr. Thad Overmyer:


The amended label referred to above, submitted in connection with registration under the Federal Insecticide, Fungicide and Rodenticide Act (FIFRA), as amended, is acceptable. This approval does not affect any conditions that were previously imposed on this registration. You continue to be subject to existing conditions on your registration and any deadlines connected with them.

A stamped copy of your labeling is enclosed for your records. This labeling supersedes all previously accepted labeling. Pursuant to 40 CFR 156.10(a)(6) you must submit one copy of the final printed labeling before you release the product for shipment with the new labeling. In accordance with 40 CFR 152.130(c), you may distribute or sell this product under the previously approved labeling for 18 months from the date of this letter. After 18 months, you may only distribute or sell this product if it bears this new revised labeling or subsequently approved labeling. "To distribute or sell" is defined under FIFRA section 2(gg) and its implementing regulation at 40 CFR 152.3.

Should you wish to add/retain a reference to the company's website on your label, then please be aware that the website becomes labeling under FIFRA and is subject to review by the Agency. See FIFRA section 2(p)(2). If the website is false or misleading, the product would be misbranded and unlawful to sell or distribute under FIFRA section 12(a)(1)(E). 40 CFR 156.10(a)(5) lists examples of statements EPA may consider false or misleading. In addition, regardless of whether a website is referenced on your product's label, claims made on the website may not substantially differ from those claims approved through the registration process, FIFRA section 12(a)(1)(B). Therefore, should the Agency find or if it is brought to our attention that a website contains false or misleading statements or claims substantially differing from the EPA approved registration, the website will be referred to the EPA's Office of Enforcement and Assurance.

Your release for shipment of the product constitutes acceptance of these conditions. If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA section 6. If you have any questions, please contact Jamie Gobreski by phone at 202-566-0748 or via email at gobreski.jamie@epa.gov.

Sincerely,

A handwritten signature in black ink that reads "Steven Snyderman". The signature is written in a cursive style with a horizontal line underneath the name.

Steven Snyderman, Product Manager 33
Regulatory Management Branch II
Antimicrobials Division (7510P)
Office of Pesticide Programs

Enclosure: Stamped Final Label

**MINT-A-KLEEN®
EPA FINAL MASTER LABEL**

Front Panel

MINT-A-KLEEN®

[PROVEN PATENTED FORMULA & PROTOCOL FOR DENTAL UNIT WATER BOTTLE & WATERLINE TREATMENT BY ANODIA SYSTEMS, LLC.]

[For Use in any Dental Unit & Water Bottle.]

DENTAL [UNIT] [WATERLINE] & [WATER] [BOTTLE] CLEANER

ACTIVE INGREDIENTS

Chlorhexidine Gluconate 0.12%

OTHER INGREDIENTS 99.88%

TOTAL 100.00%

See carton insert for additional [Precautionary Statements], [First Aid], [Directions for Use], [and] [Container Handling].

READY TO USE | DO NOT DILUTE [MINT-A-KLEEN®]

[DO NOT ATTACH MINT-A-KLEEN® BOTTLE DIRECTLY TO DENTAL UNIT.]

Net Contents: _____ OZ (ML)

[4 FL OZ (118ML)]; [8 FL OZ (236 ML)]; [12 FL OZ (355 ML)]; [16 FL OZ (473 ML)];
[20 FL OZ (591 ML)]; [32 FL OZ (946 ML)]; [64 FL OZ (1.9 L)]; 1 GALLON (3.8 L)]

**[CAUTION]
[KEEP OUT OF REACH OF CHILDREN]**

[8 54542 00800 6] [8 54542 00831 0] [8 54542 00813 6] [8 54542 00812 9]

[Mint-A-Kleen 16 oz Bottle Barcode] [Mint-A-Kleen 1 Gallon Jug Barcode] [ABN: BottleX 32 oz Bottle Barcode] [ABN: Shocking 32 oz Bottle Barcode]


Watch our Online Instructional Video at www.AnodiaSystems.com
PATENTED Product By Anodia Systems, LLC, Protected by U.S. Patent: 10,785,998. [Mint-A-Kleen] is Developed, Patented, Bottled & Owned by Anodia Systems, LLC.

MADE IN USA

Questions? Please contact us directly:

Anodia Systems, LLC
109 Larrimore Lane
Danville, KY 40422
P: 866.246.2548
F: 866.926.8246
E: info@AnodiaSystems.com
www.AnodiaSystems.com

[EPA Reg. No. 85298-2 EPA Establishment No. 85298-KY-001]



MINT-A-KLEEN®**BACK PANEL [Master Label Directions for Use Set 1 of 3]****EPA FINAL MASTER LABEL****DENTAL [UNIT] [WATERLINE] & [WATER] [BOTTLE] CLEANER**

For the cleaning and control of microbial contamination in dental unit waterlines and bottles. For use in any dental unit.

DIRECTIONS FOR USE:

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING. **[FOLLOW ONLY MINT-A-KLEEN'S DENTAL UNIT WATER BOTTLE & WATERLINE CLEANING PROTOCOL INSTRUCTIONS, STEP BY STEP.]**

INITIAL CLEANING

*Note: Each week allow Mint-A-Kleen® to remain in the dental unit waterlines and bottle overnight and over the weekend.

Week 1: Use Mint-A-Kleen® at the end of every business day.


Weeks 2-5: Use Mint-A-Kleen® twice per week at the end of the business day.

After week 5: Use Mint-A-Kleen® once per week.

INITIAL CLEANING PURPOSE: Mint-A-Kleen® preps the dental unit waterlines by gradually penetrating the microbial contamination during this 5 week period. This timely initial cleaning prevents copious amounts of contamination from clogging the waterlines.

ROUTINE WEEKLY CLEANING


- Use Mint-A-Kleen® once per week.
- Use Mint-A-Kleen® at the end of the last business day of the week.
- Leave Mint-A-Kleen® in dental unit waterlines and bottle overnight and over the weekend.
- On Monday morning before dental procedures begin, Air Purge Waterlines again then Flush Waterlines with Water.

DIRECTIONS FOR WATERLINE AND BOTTLE CLEANING:


STEP 1:
AIR PURGE WATERLINES



STEP 2:
ADD MINT-A-KLEEN®



STEP 3:
AIR PURGE WATERLINES



STEP 4:
FLUSH WATERLINES
WITH WATER

STEP 1: AIR PURGE: REMOVE DENTAL UNIT WATER FROM DENTAL UNIT WATERLINES AND BOTTLE

- After the last patient leaves at the end of the business day, turn dental unit OFF.
- Empty water bottle. Now attach this empty water bottle to dental unit.
- Turn dental unit ON.
- Air purge waterlines by activating handpiece hoses and air/water syringes until water is eliminated.

STEP 2: ADD MINT-A-KLEEN®

- Turn dental unit OFF. Remove water bottle, pour 4 ounces of Mint-A-Kleen® into dental unit's empty water bottle.
- Place sterile 2x2 gauze on dental water bottle's opening. Turn dental water bottle upside down 3 times with gauze in place. This allows Mint-A-Kleen® to coat this dental water bottle's interior walls.

**MINT-A-KLEEN®
EPA FINAL MASTER LABEL**

BACK PANEL [Master Label Directions for Use Set 1 of 3 (continued)]

- Using the same 2x2 gauze, wipe down the pick-up tube.
- Attach dental water bottle to dental unit.
- Turn dental unit ON. Activate handpiece hoses and air/water syringes until Mint-A-Kleen® exits.
- Turn dental unit OFF. Allow Mint-A-Kleen® to treat waterlines and bottle for at least 12 hours before proceeding to the next step.

STEP 3: (REPEAT STEP 1: AIR PURGE DENTAL UNIT WATERLINES)

- Turn dental unit OFF. Remove water bottle, dispose of any excess Mint-A-Kleen®. Empty water bottle. Now attach this empty water bottle to dental unit.
- Turn dental unit ON.
- Air purge waterlines by activating handpiece hoses and air/water syringes until NO liquid comes out.

STEP 4: FLUSH DENTAL UNIT WATERLINES AND WATER BOTTLE

- Turn dental unit OFF. Remove water bottle from dental unit, dispose of any excess Mint-A-Kleen®.
- Rinse dental water bottle with water. Fill water bottle with water source.
- Turn dental unit ON. Flush waterlines with water for 30 seconds to remove any excess Mint-A-Kleen®.
- Begin dental procedures.

*Note: For optimal cleaning, air purge and clean all instruments and tubing that contain water using the aforementioned protocol, including but not limited to: handpiece, air/water syringe, scaler, quick disconnect line. Allow Mint-A-Kleen® to remain in the dental waterlines while the dental office is closed for the weekend, holiday break or vacation. Dental Schools that use Mint-A-Kleen® leave Mint-A-Kleen® in their dental waterlines for up to 3 months during their summer breaks.

[The following is the Bottle Cleaning Split Label using the same intent and language from Mint-A-Kleen's Set 1 of 3 EPA Approved Directions for Use for Cleaning Dental Unit Waterlines & Dental Water Bottle (Master Label Directions for Use Set 2 of 3)]

BACK PANEL [Master Label Directions for Use Set 2 of 3]

[Mint-A-Kleen® DENTAL WATER BOTTLE CLEANER

For the Cleaning and Control of Microbial Contamination in Dental Water Bottles. For Use in Dental Water Bottle.

DIRECTIONS FOR USE:

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING.

ROUTINE WEEKLY CLEANING:

- Use Mint-A-Kleen® once per week.
- Use Mint-A-Kleen® at the end of the last business day.
- Leave Mint-A-Kleen® in dental water bottle overnight and over the weekend.
- Monday morning before dental procedures begin, rinse dental water bottle with water source.

DIRECTIONS FOR BOTTLE CLEANING:

STEP 1: REMOVE WATER FROM DENTAL WATER BOTTLE

- After the Last Patient Leaves at the End of the Business Day, Turn Dental Unit **OFF**.
- Empty Water Bottle.

MINT-A-KLEEN®
EPA FINAL MASTER LABEL

BACK PANEL [Master Label Directions for Use Set 2 of 3 (continued)]

STEP 2: ADD Mint-A-Kleen®

- Pour 4 ounces of Mint-A-Kleen® into Empty Dental Water Bottle.
- Place Sterile 2x2 Gauze on Dental Water Bottle's Opening. Turn Dental Water Bottle Upside Down 3 Times with Gauze in Place. This Allows Mint-A-Kleen® to Coat this Dental Water Bottle's Interior Walls.
- Using the Same 2x2 Gauze Soaked in Mint-A-Kleen®, Wipe Down the Pick-Up Tube.
- Attach Dental Water Bottle to Dental Unit.
- Allow Mint-A-Kleen® to Treat Bottle for at Least 12 Hours Before Proceeding to the Next Step.

STEP 3: RINSE DENTAL WATER BOTTLE

- Remove Water Bottle. Empty Water Bottle and Dispose of Any Excess Mint-A-Kleen®.
- Rinse Dental Water Bottle with Water. Fill Water Bottle with Water Source.
- Attach Dental Water Bottle to Dental Unit. Turn Dental Unit **ON**.
- Begin Dental Procedures.

BACK PANEL [Master Label Directions for Use Set 3 of 3]

[The following is the Dental Unit Water line Cleaning Split Label using the same intent and language from Mint-A-Kleen's Set 1 of 3 EPA Approved Directions for Use for Cleaning Dental Unit Water lines & Dental Water Bottle (Set 3 of 3 Master Label Directions for Use)]

[MINT-A-KLEEN® DENTAL UNIT WATER LINE CLEANER

For the Cleaning and Control of Microbial Contamination in Dental Unit Water Lines. For Use in any Dental Unit.

DIRECTIONS FOR USE:

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING.

ROUTINE WEEKLY CLEANING:

- Use Mint-A-Kleen® once per week.
- Use Mint-A-Kleen® at the end of the last business day.
- Leave Mint-A-Kleen® in dental unit water lines and bottle overnight and over the weekend.
- On Monday morning before dental procedures begin, Air Purge Waterlines again then Flush Waterlines with Water.

DIRECTIONS FOR WATER LINE CLEANING:

STEP 1: AIR PURGE: REMOVE WATER FROM DENTAL UNIT WATERLINES AND BOTTLE

- After the Last Patient Leaves at the End of the Business Day, Turn Dental Unit **OFF**.
- Empty Water Bottle. Now Attach this Empty Water Bottle to Dental Unit. Turn Dental Unit **ON**.
- Air Purge Water Lines By Activating Handpiece Hoses and Air/Water Syringes until Water is Eliminated.

STEP 2: ADD MINT-A-KLEEN®

- Turn Dental Unit **OFF**. Remove Water Bottle, Pour 4 ounces of Mint-A-Kleen® into Empty Dental Water Bottle.
- Place Sterile 2x2 Gauze on Dental Water Bottle's Opening. Turn Dental Water Bottle Upside Down 3 Times with Gauze in Place. This Allows Mint-A-Kleen® to Coat this Dental Water Bottle's Interior Walls.
- Using the Same 2x2 Gauze Soaked with Mint-A-Kleen®, Wipe Down the Pick-Up Tube.

**MINT-A-KLEEN®
EPA FINAL MASTER LABEL**

BACK PANEL [Master Label Directions for Use Set 3 of 3 (continued)]

- Attach Dental Water Bottle to Dental Unit. Turn Dental Unit **ON**.
- Activate Handpiece Hoses and Air/Water Syringes until Mint-A-Kleen® Exits. Turn Dental Unit **OFF**.
- Allow Mint-A-Kleen® to Treat Water Lines and Bottle for at Least 12 Hours Before Proceeding to the Next Step.

STEP 3: (REPEAT STEP 1: AIR PURGE DENTAL UNIT WATER LINES)

- Turn Dental Unit **OFF**. Remove Water Bottle, Dispose of Any Excess Mint-A-Kleen®. Empty Water Bottle. Now Attach this Empty Water Bottle to Dental Unit. Turn Dental Unit **ON**.
- Air Purge Water Lines by Activating Handpiece Hoses and Air/Water Syringes until NO Liquid comes out.

STEP 4: FLUSH DENTAL UNIT WATERLINES

- Turn Dental Unit **OFF**. Remove Water Bottle from Dental Unit. Rinse Dental Water Bottle with Water.
- Fill Water Bottle with Water. Attach Dental Water Bottle to Dental Unit. Turn Dental Unit **ON**.
- Flush Water Lines with Water for 30 Seconds to Remove any Excess Mint-A-Kleen®. Begin Dental Procedures.

***NOTE:** Allow Mint-A-Kleen® to Remain in the Dental Water Lines while the Dental Office is Closed for the Weekend, Holiday Break or Vacation. Dental Schools that use Mint-A-Kleen® Leave Mint-A-Kleen® in their Dental Water Lines for up to 3 Months During their Summer Breaks.

CONTAINER HANDLING: Do not contaminate water, food or feed by storage and disposal.

Storage: store in a cool dark room away from heat or open flame. **Disposal:** wastes resulting from the use of this product must be disposed of on site or at an approved waste disposal facility.

Container: non-refillable container. Do not reuse or refill container. Rinse with water and place in trash or offer for recycling if available.

[PRECAUTIONARY STATEMENTS: Hazards to Humans and Domestic Animals. **CAUTION:** Causes moderate eye irritation. Harmful if inhaled or swallowed. Avoid contact with eyes or clothing or breathing vapor. Remove and wash contaminated clothing before reuse. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. **FIRST AID:** Have the product container or label with you when calling a poison control center or doctor, or going for treatment. For non-emergency information on this product, call Anodia Systems, LLC (1-866-246-2548), Monday through Friday, 9 a.m. to 5 p.m. EST. For medical emergencies, call the poison control center at 1-800-222-1222. **IF IN EYES:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice. **IF INHALED:** Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably mouth-to-mouth if possible. Call a poison control center or doctor for further treatment advice. **IF SWALLOWED:** Call a poison control center or doctor immediately for treatment advice. Have a person sip a glass of water if able to swallow. Do not induce vomiting unless told to by a poison control center or doctor. Do not give anything to an unconscious person. **PHYSICAL OR CHEMICAL HAZARDS.** Combustible. Do not use or store near heat or open flame.]

**[Anodia Systems, LLC
Phone:866.246.2548
www.AnodiaSystems.com]**

[Developed, Patented & Bottled
By Anodia Systems LLC]

EPA Reg. No. 85298-2 EPA Establishment No. 85298-KY-001

FINAL MASTER CARTON INSERT UPDATES**MINT-A-KLEEN®**

[PROVEN PATENTED FORMULA & PROTOCOL FOR DENTAL UNIT WATER BOTTLE & [WATERLINE] TREATMENT BY ANODIA SYSTEMS, LLC.]

[For Use in any Dental Unit & Water Bottle.]

DENTAL [UNIT] [WATERLINE] & [WATER] [BOTTLE] CLEANER

For the cleaning and control of microbial contamination in dental unit waterlines and bottles. For use in any dental unit & [water bottle].

ACTIVE INGREDIENTS

Chlorhexidine Gluconate 0.12%

OTHER INGREDIENTS 99.88%

TOTAL 100.00%

READY TO USE | DO NOT DILUTE [MINT-A-KLEEN®]

[DO NOT ATTACH MINT-A-KLEEN® BOTTLE DIRECTLY TO DENTAL UNIT.]

Net Contents: _____ OZ (ML)

[4 FL OZ (118ML)]; [8 FL OZ (236 ML)]; [12 FL OZ (355 ML)]; [16 FL OZ (473 ML)];
[20 FL OZ (591 ML)]; [32 FL OZ (946 ML)]; [64 FL OZ (1.9 L)]; 1 GALLON (3.8 L)]

[CAUTION]

[KEEP OUT OF REACH OF CHILDREN]

Watch our Online Instructional Video at www.AnodiaSystems.com

PATENTED Product By Anodia Systems, LLC, Protected by U.S. Patent: 10,785,998. [Mint-A-Kleen] is Developed, Patented, Bottled & Owned by Anodia Systems, LLC.

DIRECTIONS FOR USE:

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING. **[FOLLOW ONLY MINT-A-KLEEN'S DENTAL UNIT WATER BOTTLE & WATERLINE CLEANING PROTOCOL INSTRUCTIONS, STEP BY STEP.]**

[MINT-A-KLEEN® DENTAL [UNIT] [WATERLINE] & [WATER] [BOTTLE] CLEANER

For the cleaning and control of microbial contamination in dental unit waterlines and bottles. For use in any dental [unit] [&] [water] [bottle].

INITIAL CLEANING

*Note: Each week allow Mint-A-Kleen® to remain in the dental unit waterlines and bottle overnight and over the weekend.

Week 1: Use Mint-A-Kleen® at the end of every business day.

Weeks 2-5: Use Mint-A-Kleen® twice per week at the end of the business day.

After week 5: Use Mint-A-Kleen® once per week.

INITIAL CLEANING PURPOSE: Mint-A-Kleen® preps the dental unit waterlines by gradually penetrating the microbial contamination during this 5 week period. This timely initial cleaning prevents copious amounts of contamination from clogging the waterlines.


ROUTINE WEEKLY CLEANING

- Use Mint-A-Kleen® once per week.

FINAL MASTER CARTON INSERT UPDATES
[Directions for Use Set 1 of 3 Continued]

- Use Mint-A-Kleen® at the end of the last business day.
- Leave Mint-A-Kleen® in dental unit waterlines and bottle overnight and over the weekend.
- On Monday morning before dental procedures begin, Air Purge Waterlines again then Flush Waterlines with Water.


DIRECTIONS FOR WATERLINE AND BOTTLE CLEANING:



STEP 1:
AIR PURGE WATERLINES



STEP 2:
ADD MINT-A-KLEEN®



STEP 3:
AIR PURGE WATERLINES



STEP 4:
FLUSH WATERLINES
WITH WATER

STEP 1: AIR PURGE: REMOVE DENTAL UNIT WATER FROM DENTAL UNIT WATERLINES AND BOTTLE

- After the last patient leaves at the end of the business day, turn dental unit OFF.
- Empty water bottle. Now attach this empty water bottle to dental unit.
- Turn dental unit ON.
- Air purge waterlines by activating handpiece hoses and air/water syringes until water is eliminated.

STEP 2: ADD MINT-A-KLEEN®

- Turn dental unit OFF. Remove water bottle, pour 4 ounces of Mint-A-Kleen® into dental unit's empty water bottle.
- Place sterile 2x2 gauze on dental water bottle's opening. Turn dental water bottle upside down 3 times with gauze in place. This allows Mint-A-Kleen® to coat this dental water bottle's interior walls.
- Using the same 2x2 gauze, wipe down the pick-up tube.
- Attach dental water bottle to dental unit.
- Turn dental unit ON. Activate handpiece hoses and air/water syringes until Mint-A-Kleen®

[Master Carton Insert Directions for Use Set 1 of 3]

- exits.
- Turn dental unit OFF. Allow Mint-A-Kleen® to treat waterlines and bottle for at least 12 hours before proceeding to the next step.

STEP 3: (REPEAT STEP 1: AIR PURGE DENTAL UNIT WATERLINES)

- Turn dental unit OFF. Remove water bottle, dispose of any excess Mint-A-Kleen®. Empty water bottle. Now attach this empty water bottle to dental unit.
- Turn dental unit ON.
- Air purge waterlines by activating handpiece hoses and air/water syringes until NO liquid comes out.

STEP 4: FLUSH DENTAL UNIT WATERLINES AND WATER BOTTLE

- Turn dental unit OFF. Remove water bottle from dental unit, dispose of any excess Mint-A-Kleen®.
- Rinse dental water bottle with water. Fill water bottle with water source.
- Turn dental unit ON. Flush waterlines with water for 30 seconds to remove any excess Mint-A-Kleen®.
- Begin dental procedures.]

FINAL MASTER CARTON INSERT UPDATES [Directions for Use Set 2 of 3]

[The following is the Bottle Cleaning Directions For Use using the same intent and language from Mint-A-Kleen's Set 1 of 3 EPA Approved Directions for Use for Cleaning Dental Unit Waterlines & Dental Water Bottle (Master Carton Insert Directions for Use Set 2 of 3)]

[Mint-A-Kleen® DENTAL WATER BOTTLE CLEANER

For the Cleaning and Control of Microbial Contamination in Dental Water Bottles. For Use in Dental Water Bottle.

DIRECTIONS FOR USE:

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING.

ROUTINE WEEKLY CLEANING:

- Use Mint-A-Kleen® once per week.
- Use Mint-A-Kleen® at the end of the last business day.
- Leave Mint-A-Kleen® in dental water bottle overnight and over the weekend.
- Monday morning before dental procedures begin, rinse dental water bottle with water source.

DIRECTIONS FOR BOTTLE CLEANING:

STEP 1: REMOVE WATER FROM DENTAL WATER BOTTLE

- After the Last Patient Leaves at the End of the Business Day, Turn Dental Unit **OFF**.
- Empty Water Bottle.

STEP 2: ADD Mint-A-Kleen®

- Pour 4 ounces of Mint-A-Kleen® into Empty Dental Water Bottle.
- Place Sterile 2x2 Gauze on Dental Water Bottle's Opening. Turn Dental Water Bottle Upside Down 3 Times with Gauze in Place. This Allows Mint-A-Kleen® to Coat this Dental Water Bottle's Interior Walls.
- Using the Same 2x2 Gauze Soaked in Mint-A-Kleen®, Wipe Down the Pick-Up Tube.
- Attach Dental Water Bottle to Dental Unit.
- Allow Mint-A-Kleen® to Treat Bottle for at Least 12 Hours Before Proceeding to the Next Step.

STEP 3: RINSE DENTAL WATER BOTTLE

- Remove Water Bottle. Empty Water Bottle and Dispose of Any Excess Mint-A-Kleen®.
- Rinse Dental Water Bottle with Water. Fill Water Bottle with Water Source.
- Attach Dental Water Bottle to Dental Unit. Turn Dental Unit **ON**.
- Begin Dental Procedures.]

[The following is the Dental Unit Water line Cleaning Directions For Use using the same intent and language from Mint-A-Kleen's Set 1 of 3 EPA Approved Directions for Use for Cleaning Dental Unit Water lines & Dental Water Bottle (Set 3 of 3 Master Label Directions for Use)]

FINAL MASTER CARTON INSERT UPDATES [Directions for Use Set 3 of 3]

[**MINT-A-KLEEN® DENTAL UNIT WATER LINE CLEANER**

For the Cleaning and Control of Microbial Contamination in Dental Unit Water Lines. For Use in any Dental Unit.

DIRECTIONS FOR USE:

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING.

ROUTINE WEEKLY CLEANING:

- Use Mint-A-Kleen® once per week.
- Use Mint-A-Kleen® at the end of the last business day.
- Leave Mint-A-Kleen® in dental unit water lines and bottle overnight and over the weekend.
- On Monday morning before dental procedures begin, Air Purge Waterlines again then Flush Waterlines with Water.

DIRECTIONS FOR WATER LINE CLEANING:

STEP 1: AIR PURGE: REMOVE WATER FROM DENTAL UNIT WATERLINES AND BOTTLE

- After the Last Patient Leaves at the End of the Business Day, Turn Dental Unit **OFF**.
- Empty Water Bottle. Now Attach this Empty Water Bottle to Dental Unit. Turn Dental Unit **ON**.
- Air Purge Water Lines By Activating Handpiece Hoses and Air/Water Syringes until Water is Eliminated.

STEP 2: ADD MINT-A-KLEEN®

- Turn Dental Unit **OFF**. Remove Water Bottle, Pour 4 ounces of Mint-A-Kleen® into Empty Dental Water Bottle.
- Place Sterile 2x2 Gauze on Dental Water Bottle's Opening. Turn Dental Water Bottle Upside Down 3 Times with Gauze in Place. This Allows Mint-A-Kleen® to Coat this Dental Water Bottle's Interior Walls.
- Using the Same 2x2 Gauze Soaked with Mint-A-Kleen®, Wipe Down the Pick-Up Tube.
- Attach Dental Water Bottle to Dental Unit. Turn Dental Unit **ON**.
- Activate Handpiece Hoses and Air/Water Syringes until Mint-A-Kleen® Exits. Turn Dental Unit **OFF**.
- Allow Mint-A-Kleen® to Treat Water Lines and Bottle for at Least 12 Hours Before Proceeding to the Next Step.

STEP 3: (REPEAT STEP 1: AIR PURGE DENTAL UNIT WATER LINES)

- Turn Dental Unit **OFF**. Remove Water Bottle, Dispose of Any Excess Mint-A-Kleen®. Empty Water Bottle. Now Attach this Empty Water Bottle to Dental Unit. Turn Dental Unit **ON**.
- Air Purge Water Lines by Activating Handpiece Hoses and Air/Water Syringes until NO Liquid comes out.

STEP 4: FLUSH DENTAL UNIT WATERLINES

- Turn Dental Unit **OFF**. Remove Water Bottle from Dental Unit. Rinse Dental Water Bottle with Water.
- Fill Water Bottle with Water. Attach Dental Water Bottle to Dental Unit. Turn Dental Unit **ON**.
- Flush Water Lines with Water for 30 Seconds to Remove any Excess Mint-A-Kleen®. Begin Dental Procedures.

**MINT-A-KLEEN®
EPA FINAL MASTER LABEL UPDATES**

[*Note: For optimal cleaning, air purge and clean all instruments and tubing that contain water using the aforementioned protocol, including but not limited to: handpiece, air/water syringe, scaler, quick disconnect line.] Allow Mint-A-Kleen® to remain in the dental waterlines while the dental office is closed for the weekend, holiday break or vacation. Dental Schools that use Mint-A-Kleen® leave Mint-A-Kleen® in their dental waterlines for up to 3 months during their summer breaks.

CONTAINER HANDLING: Do not contaminate water, food or feed by storage and disposal.

Storage: store in a cool dark room away from heat or open flame.

Disposal: wastes resulting from the use of this product must be disposed of on site or at an approved waste disposal facility.

Container: non-refillable container. Do not reuse or refill container. Rinse with water and place in trash or offer for recycling if available.

PRECAUTIONARY STATEMENTS:

Hazards to Humans and Domestic Animals.

CAUTION: Causes moderate eye irritation. Harmful if inhaled or swallowed. Avoid contact with eyes or clothing or breathing vapor. Remove and wash contaminated clothing before reuse. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

FIRST AID:

Have the product container or label with you when calling a poison control center or doctor, or going for treatment. For non-emergency information on this product, call Anodia Systems, LLC (1-866-246-2548), Monday through Friday, 9 a.m. to 5 p.m. EST. For medical emergencies, call the poison control center at 1-800-222-1222.

IF IN EYES: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.

IF INHALED: Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably mouth-to-mouth if possible. Call a poison control center or doctor for further treatment advice.

IF SWALLOWED: Call a poison control center or doctor immediately for treatment advice. Have a person sip a glass of water if able to swallow. Do not induce vomiting unless told to by a poison control center or doctor. Do not give anything to an unconscious person. Have the product container or label with you when calling a poison control center or doctor or going for treatment.

PHYSICAL OR CHEMICAL HAZARDS. Combustible. Do not use or store near heat or open flame.



Watch our Online Instructional Video at www.AnodiaSystems.com

If you have any questions, please contact us directly.

Anodia Systems, LLC

109 Larrimore Lane | Danville, KY 40422 | P: 866.246.2548 | F : 866.926.8246

E : info@AnodiaSystems.com | www.AnodiaSystems.com