



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
WASHINGTON, DC 20460

OFFICE OF CHEMICAL SAFETY
AND POLLUTION PREVENTION

February 2, 2015

Ms. Susan Stracquatano
Agent For: Process Cleaning Solutions, Ltd.
Delta Analytical Corporation
12510 Prosperity Drive; Suite 160
Silver Spring, Maryland 20904

Subject: Label Notification per PRN 98-10 – Alternate Brand Name
Product Name: **PCS 7000**
EPA Registration Number: **90271-1**
Application Date: November 20, 2014
Decision Number: 499494

Dear Ms. Stracquatano:

The Agency is in receipt of your Application for Pesticide Notification under Pesticide Registration Notice (PRN) 98-10 for the above referenced product. The Antimicrobials Division has conducted a review of this request for its applicability under PRN 98-10 and finds that the action requested falls within the scope of PRN 98-10.

The label submitted with the application has been stamped "Notification" and will be placed in our records.

The following alternate brand names have been added to the product record:

- PCS 7000 Disinfectant Cleaner

If you have any questions, you may contact Killian Swift at 703-308-6346 or via email at Swift.Killian@epa.gov.

Sincerely yours,

A handwritten signature in black ink, appearing to read "DF", written over a light blue horizontal line.

Demson Fuller
EPA Product Manager 32
Regulatory Management Branch II
Antimicrobials Division 7510P

NOTIFICATION

90271-1

The applicant has certified that no changes, other than those reported to the Agency have been made to the labeling. The Agency acknowledges this notification by letter dated:

02/02/2015

PCS 7000

[Alternate Brand Name: ECOGent 7000]

[Alternate Brand Name: PCS 7000 Disinfectant Cleaner]

[Alternate Brand Name: ECOGent 7000 Disinfectant Cleaner]

For use on pre-cleaned hard non-porous environmental surfaces in health care facilities, institutions, schools and hospitality industries, where organic soils may be present.

Kills Bacteria (Bactericide)

Kills *Pseudomonas aeruginosa*, *Salmonella enterica* and *Staphylococcus aureus*

Kills *Clostridium difficile** (spores)

ACTIVE INGREDIENT:

Sodium Hypochlorite0.735%

OTHER INGREDIENTS:..... 99.265%

TOTAL 100.000%

Available Chlorine 0.7%

KEEP OUT OF REACH OF CHILDREN

CAUTION

See back panel for complete precautions

FIRST AID	
If in Eyes:	<ul style="list-style-type: none"> • Hold eye open and rinse slowly and gently with water for 15-20 minutes. • Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. • Call a poison control center or doctor for treatment advice.
If Swallowed:	<ul style="list-style-type: none"> • Call a poison control center or doctor immediately for treatment advice. • Have person sip a glass of water if able to swallow. • Do not induce vomiting unless told to do so by a poison control center or doctor. • Do not give anything by mouth to an unconscious person.
If on Skin or Clothing:	<ul style="list-style-type: none"> • Take off contaminated clothing. • Rinse skin immediately with plenty of water for 15-20 minutes. • Call a poison control center or doctor for treatment advice.
Have the product container or label with you when calling a poison control center or doctor, or going for treatment. In case of emergency, call poison control center 1-800-222-1222 for treatment advice.	

Manufactured for:
Process Cleaning Solutions Ltd
Peterborough, Ontario K9J8N4

EPA Reg. No. 90271-1
EPA Est. No. 67900-CAN-001

Net Contents:

PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals

CAUTION. Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

Physical or Chemical Hazards

Do not use this product with ammonia or acids such as vinegar, rust removers or toilet bowl cleaners. This product is not intended for use on fabrics.

[In accordance with PR notice 95-1, use the following complete Environmental Hazards statement for containers 5 gallons and larger]

Environmental Hazards

This product is toxic to fish and aquatic organisms. Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans or public waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office of the EPA.

[For containers smaller than 5 gallons use the following:]

Environmental Hazards

This product is toxic to fish and aquatic organisms.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

DISINFECTION

To clean and disinfect walls, tables, counter tops or other hard nonporous surfaces:

- Remove gross soil prior to disinfecting surfaces.
- Apply undiluted with a clean mop, cloth or sponge to pre cleaned surfaces. (When applying with a saturated cloth or sponge, gloves should be worn.)
- Completely wet surface and keep wet for five *or* [5] minutes.
- Wipe with a clean damp cloth or paper towel or allow to air dry.
- Rinse with potable water surfaces that come in contact with food and surfaces or objects that come in contact with children at the mouthing stage of development. For all other surfaces, no rinsing is required.

Avoid contamination of food.

Special Instructions for Cleaning Prior to Disinfection against *Clostridium difficile* spores

Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks or eye covering.

Cleaning Procedure: Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop and/or sponge saturated with this product. Cleaning is to include vigorous wiping and/or scrubbing until all visible soil is removed. Special attention is needed for high touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner such as from right to left or left to right on horizontal surfaces and top to bottom on vertical surfaces to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.

Infectious Materials Disposal: Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

Disinfection Procedure: Apply with a clean cloth, mop and/or sponge saturated with this product until surface is completely wet. Let stand for five *or* [5] minutes. Wipe with a clean damp cloth or paper towel or allow to air dry.

To Disinfect Non Critical Pre Cleaned Instruments - Instruments must be thoroughly pre cleaned to remove excess organic debris, rinsed, and rough dried. (Clean and rinse lumens of hollow instruments before filling with this product or before immersion.)

Immersion method: using a soaking tray, immerse instruments into this product and let stand for five *or* [5] minutes.

To kill *Clostridium difficile** *or C difficile** *or C diff** (spores) let stand for five *or* [5] minutes. Wipe with a clean damp cloth or paper towel and allow to air dry.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre clean or decontaminate critical or semi critical medical devices prior to sterilization or high level disinfection.

SANITIZATION OF NONPOROUS FOOD CONTACT SURFACES

Rinse Method: Prepare a sanitizing solution by thoroughly mixing 1 part of this product with 32 parts of water to provide approximately 200 ppm available chlorine by weight. Clean equipment surfaces in normal manner. Prior to use, rinse all surfaces thoroughly with the sanitizing solution, maintaining contact with the sanitizer for at least 2 minutes. Do not rinse equipment with water after treatment and do not soak equipment overnight.

Immersion Method: Prepare a sanitizing solution by thoroughly mixing 1 part of this product with 32 parts of water to provide approximately 200 ppm available chlorine by weight. Clean equipment in the normal manner. Prior to use, immerse equipment in the sanitizing solution for at least 2 minutes and allow the sanitizer to drain. Do not rinse equipment with water after treatment.

*Follow [the] Special Instructions for Cleaning Prior to Disinfection [listed on label].

STORAGE AND DISPOSAL

Do not contaminate food or feed by storage, disposal or cleaning of equipment.

Storage: Store in a cool dry area, away from direct sunlight and heat to avoid deterioration. In case of spill, flood areas with large quantities of water. **Product Disposal:** Product or rinsates that cannot be used must be diluted with water before disposal in a sanitary sewer. **Container Handling:** Nonrefillable container. Do not reuse or refill this container. Clean container promptly after emptying. Triple rinse as follows: Fill the container $\frac{1}{4}$ full with water and recap. Shake for 10 seconds. Drain for 10 seconds after the flow begins to drip and dispose of rinsate in a sanitary sewer. Repeat procedure two more times. Offer container for recycling if available or reconditioning if appropriate or place in trash.

[batch code]